



Symbiq Quick Tips-Alarms/Messages

To clear persistent **PROXIMAL OCCLUSION** alarm (after troubleshooting the line for occlusion):

- Press **Start** to clear alarm
- Press **Stop**
- Press **Load/Eject** to open and close the cassette door (to reset pressures)
- Resume infusion

If the alarm persists:

- Stop the infusion, **turn the pump off and back on**
- When the new patient screen is displayed, Change *New Patient?* to **no** (and resume the infusion)

If the alarm persists:

- Change the tubing

CHECK FLOWSTOP alarm when **Opening** Cassette Door:

- Press Silence. Flowstop may have been bumped when removing it from door - check cassette and **ensure flowstop is CLOSED**

CHECK FLOWSTOP alarm when **Closing** Cassette Door:

- Press Silence. Press **Load/Eject** button to open & close door
- Check cassette and **ensure flowstop is CLOSED**
- Check for closed clamps and open clamps, reload tubing (wait for door to close completely)
- Start infusion

CHECK CASSETTE alarm:

- Press **Alarm Tab** at the bottom of the screen to clear the alarm message
- Press **Load/Eject** to open cassette door
- Ensure cassette is “seated” properly - tug on tubing below cassette.
- Allow door to close completely before starting infusion

CHECK CASSETTE/Cassette Required Message:

- Typically happens when Start is pressed before door has completed closing cycle. **Wait until LED light (above channel) has gone off before pushing Start.**
- Press **OK** to message
- Press **Start**

AIR-IN-LINE alarm:

- Press Silence. Disconnect from patient
- Remove cassette, open Flowstop, invert and tap to dislodge air bubbles
- Reprime set to remove air
- Close Flowstop & insert cassette back in pump, reconnect set to patient
- Restart infusion



Symbiq Helpful Hints:

Forced Shutdown Message:

- May be caused by ...Trying to turn pump off before door is closed
- Please ...Select "OK" and allow door to shut completely then turn pump off.

***DO NOT FORCE SHUTDOWN UNLESS ABSOLUTELY NECESSARY

Secondary Line Maximum Rate of 500 mL/hr:

- If an infusion was programmed for >500mL/hr on a secondary line, it would begin pulling from primary container (due to backcheck valve). The pump is preventing this situation from happening.

PRIMING/CONNECTIONS:

- **Ensure that cassette is inverted, flowstop (white piece) is open, tap, then close flowstop when finished.**
- **At same time you invert the cassette, INVERT BACKCHECK valve and tap out air**; do the same with Clave ports.
- Priming set with vent for Glass bottles:
 - Ensure flowstop is closed.
 - Insert spike into bottle
 - Squeeze drip chamber
 - Open vent on piercing pin
 - **Open flowstop to prime**
- Changing glass bottles:
 - **CLOSE VENT WHEN RE-SPIKING BOTTLE**
 - Spike new bottle
 - Squeeze drip chamber
 - Open vent on piercing pin
- **When connecting male tip to Clave port, pull back collar, insert tip, push and turn.** After the tip is securely seated, then luer on collar. If you are experiencing distal occlusions, check your connection. If a secondary is not running, check your connection.

QUICKER SET UP:

- Press and hold On/Off button while bringing pump to patient room
- Can begin programming while pump is finishing self test (and opening and closing of doors)
- Load primed cassette. Once door has finished opening cycle (listen), can press Load/Eject button to close the door more quickly (quicker than waiting the 6 seconds for it to automatically close)
- FYI...Pump calibrates with cassette each time the door has been opened and closed (to ensure accurate delivery)

STANDBY:

Please use with caution: *Pump will not alarm to remind you it's in standby for 24 hours.*

Alternative to Standby: You can pause the infusion for a blood draw by pressing the Stop key and then silence the callback alarm if you get it.



Symbiq Helpful Hints:

HOW TO CLEAR A PROGRAM OR TURN 'A' CHANNEL OFF:

- Stop infusion
- Press Basic (if it's a program under the basic tab you want to clear)
- Press Clear button, press Continue to finish clearing program

HOW TO CLEAR A PIGGYBACK PROGRAM AND RESUME PRIMARY INFUSION:

- Stop Piggyback infusion
- Must clear piggyback program to resume basic infusion
- Press Piggyback tab
- Press Clear button, press Continue to finish clearing program
- Press Basic tab to get back to Basic program
- Press Start

***CHANGING A PROGRAM PARAMETER

IT IS IMPORTANT TO REMEMBER:

- If VTBI is changed, then delivery time will be changed
- If delivery time is changed, then VTBI will be changed
- The rate or dose will not change unless you access those fields and change them.
- Please look at Confirmation Screen parameters before starting an infusion.