



REIXS Beamline User Manual

(Work in Progress)

Feizhou He


Updated: November 7, 2011

TABLE OF CONTENTS

1.0	Introduction	1
2.0	Beamline User Controls	1
2.1	User Interface	1
2.1.1	Shutters, Valves and Vacuum Pressure.....	2
2.1.2	Beamline Parameters Monitoring.....	3
2.1.3	Additional Controls	5
2.2	Setup Endstation Branch.....	5
2.3	Selecting Energy Range and Polarizations.....	6
2.3.1	Selecting M2 Coatings and Gratings.....	7
2.3.2	Selecting EPU Harmonics and Photon Polarization	7
2.4	Setup Chopper	8
2.5	Opening the Beamline to the Ring.....	8
2.6	Getting the Beam into the Endstation	9
2.7	Beamline Monitoring and Diagnostics	10
2.8	Other.....	11
3.0	RSXS Endstation	12
3.1	Instruments	12
3.2	Motion Control.....	14
3.3	Data Collection Software.....	15
4.0	XES Endstation	15
4.1	Instruments	15
4.2	Motion Control.....	15
4.3	Data Collection.....	16

1.0 Introduction

The Resonant Elastic and Inelastic Soft X-ray Scattering (REIXS) Beamline is a soft X-ray beamline dedicated to soft X-ray scattering and soft X-ray spectroscopy experiments. The beamline is located at 10ID-2 port of Canadian Light Source.



Location of the REIXS Beamline on the experimental floor of CLS

This document describes the components and operations of REIXS Beamline.

2.0 Beamline User Controls


Two REIXS Operator Interface (OPI) computers (OPI1610-201 and OPI1610-202) are located inside the REIXS Cleanroom. The computer operating system is Scientific Linux. Users need to log in to access the beamline control.

2.1 User Interface

The REIXS Beamline User Interface can be accessed by typing

```
runREIXS_beamline &
```

from a command line window (xterm), or by double click the desktop icon "REIXS Beamline". The interface is written in the EPICS display-manager package EDM. It provides user with essential beamline controls. Important error conditions are also displayed in this window.




REIXS Beamline Main User Interface

2.1.1 Shutters, Valves and Vacuum Pressure

The upper half of the user interface includes the shutter/valve control, and shows the vacuum pressure along the beamline. The horizontal straight line represents the beam, from the frontend on the left side to the endstations on the right side.

The black letters below the beamline are the acronyms of each vacuum component:

PSH2	Photon Shutter 2
SSH	Safety Shutter
FM3	Fixed Mask 3
VA	Variable Aperture
M1A/B	M1A and M1B Mirrors
PSH3	Photon Shutter 3
4J1	4-Jaw #1
MONO	Monochromator
4J2	4-Jaw #2
CHP/ES	Chopper and Exit Slit
M3/M4	M3 and M4 mirrors
PSH4	Photon Shutter 4
M5	M5 Mirror
GC	Gas Cell
XES	X-Ray Emission Endstation (UofS)
RSXS	Resonant Soft X-Ray Scattering Endstation (UBC)



Shutters are shown as rectangles in or out of the beam, and the valves are double triangles. Click the [Open] or [Close] buttons above each shutter/valve to open or close corresponding shutter or valve. When closed, the shutter or valve is in red color. Open state is in green.

The numbers below the name of each vacuum section are the vacuum pressure, in Torr. The number is green when the pressure is less than 10×10^{-10} Torr, and yellow when the pressure is higher than 1×10^{-9} but less than 5×10^{-8} Torr. The number is red for pressure higher than 5×10^{-8} Torr. The vacuum pressure in any section shall be kept below 1×10^{-8} Torr. If vacuum pressure is too high, try using a smaller aperture size on Variable Aperture, or closing the Photon Shutter 3 until the pressure recovers.

2.1.2 Beamline Parameters Monitoring

The lower half of the user interface shows the status of the beamline components. The parameters are grouped into several panels.

EPU in 10ID




The 10ID panel shows the gap and phase information of both REIXS EPU and SM EPU. The polarization of REIXS EPU can be changed using the drop down menu button. The options include: Circular Left, Circular Right, Linear Horizontal, Linear Vertical -, Linear Vertical +, and Linear Inclined.

The small circles indicate the status of the respective motors, green means "Move Done", yellow means "Moving".

The EPU Harmonics selection button can be used to select which harmonics shall be used for the energy range.

Apertures




The Aperture panel shows the gap and position of each aperture. H – horizontal, V – vertical.

The setpoints for the Variable Apertures can be changed from this panel. The setpoints are in purple color, and the feedback values are black numbers above each input boxes.

The small circles indicate the status of the respective motors, green means "Move Done", yellow means "Moving".

Chopper



The speed and position of the chopper can be set in the Chopper panel. Use ONLY positive numbers for the "Speed" and "Move Relative" input. One revolution per second equals to 12 Hz switching. The blue numbers shows the current speed and position of the chopper.

Click [Stop] button will stop the chopper rotation. Once chopper is stopped, use "Move Relative" to adjust chopper position.

Monochromator

Monochromator	
Energy Setpoint	950.000 eV
Stop	Feedback 949.999 eV
M2 Mirror	1.698931 deg ●
Difference	0.000008
Translation	96.0310 mm ●
Grating	3.587376 deg ●
Difference	0.000004
Translation	139.5190 mm ●

The status of the monochromator is displayed here. The energy setpoint can be changed. The angles and translation positions of the mirror / grating are displayed, as well as the error (difference) between the real angles and the setpoints. The error is usually limited to the last two digits of the black numbers. The [Stop] button will stop the monochromator and EPU motions immediately.

Clicking the [Monochromator] button will bring up the "REIXS Monochromator" control interface.

Meters

Meters	
4-Jaw #1	Avg
Upper	2.9854 uA
Lower	3.8113 uA
ESlit	0.004 uA 5
M4 Au	0.435 nA 45
M5 Au	0.004 nA 8
M5 PD	0.000 uA 5
GC Au	0.000 uA 8
GC PD	0.038 uA 8
GCell	0.000 nA 5
Voltage	10.000 V OFF
EM	1.771e-12 A 45

The Meters panel shows the readings from meters along the beamline.

The realtime readings from each meter are in black color. The numbers on the right side shows the number of samples for averaging. When the number is green, the averaging is on. The number is red if averaging is off.

PD means Photodiode, Au means Gold Mesh. EM is the Keithley Electrometer connected to the Photodiode in the RSXS Endstation.

Clicking the [Meters] button will bring up the "REIXS Meters Setup" interface.

Scaler and RSXS Cryostat

Scaler		
Scan	Continuous	Counting Time
Stopped	Normal	5000.0 mSec
REIXS		SM
Timing	50.000	0.000 MHz
ChT	0.000	0.000
TEY	248388.000	0.000
MCP	194359.000	
RSXS Sample Temperature		
Setpoint	100 K	Heater Range
		High (25 W)
Sample	100.000 K	Heater Output
Cryostat	13.067 K	56 %

The Scaler panel shows the status of the SIS3820 scaler. The counting time can be set here, in milli-second. [Scan] button is used to trigger a single counting, while [Continuous] button will trigger continuous counting with the interval set in the Counting Time box.

The lower panel shows the results of last counting. The Timing value represents the actual counting time, 10MHz per second. The REIXS column and SM column show the counting results for the beam from REIXS EPU and SM EPU respectively, in the REIXS 2-in-1 operation mode.

ChT means Channeltron, TEY means Total Electron Yield, and MCP mean microchannel plate.

The RSXS Sample temperature panel is used to display and change the sample temperature in the RSXS chamber. The status of the cryostat and the heater is also shown.

2.1.3 Additional Controls

The buttons at the bottom of the user interface will call up additional control panels. Some of them will be explained in the following sections.

2.2 Setup Endstation Branch

The REIXS beamline has two endstations: XES and RSXS. At one time the beamline is setup for one of the endstation. Please ask beamline staff if you need to switch endstations.

Switching endstation involves two steps:

Move M3 mirror in or out of the beam

Translate XES endstation in or out of the beam

To use RSXS Endstation:


M3 in the beam: lift M3 tank until the encoder reach 0.8660mm.

XES out of the beam: Loosen the locking nuts on the switch yard on M5 tank. Move XES Endstation backward to limit switch. Hand tighten the locking nuts.

To use XES Endstation:

M3 out of the beam: lower the M3 tank until the encoder reach -0.XXXXmm.

XES in the beam: Loosen the locking nuts on the switch yard on M5 tank. Move XES Endstation forward until encoder reads 0.XXXXmm. Hand tighten the locking nuts.




Four locking nuts

M5 switch yard

2.3 Selecting Energy Range and Polarizations

The "REIXS Monochromator" window is used to setup monochromator and REIXS EPU.



In this window, you will be able to set the photon energy and select the mirror coatings and gratings.

The [STOP] button on top will stop REIXS EPU and Monochromator immediately.


The "Mono Calib" button is used to fine tune the energy calibration. Please refer to Section **Error! Reference source not found.** for details.

polarization of the photons from the EPU.

The REIXS EPU panel allows you to select the harmonic of the EPU output, and the

2.3.1 Selecting M2 Coatings and Gratings

The REIXS monochromator has three gratings, and there are four coatings on the M2 mirror.




In the "M2 Mirror" and "Grating" panel, user can select which grating and coating are used. The setpoint and feedback values are also displayed.

The table below lists the maximum usable range for each combination of the M2 Coating / Grating. For best flux and good higher order harmonics suppression, the optimized range shall be used.

Grating	M2 coating	Overall Range (eV)	Optimized Range (eV)
Ni-LEG	Si (SiO ₂)	167-525	230-525
	Ni	167-800	380-800
Au-LEG	Graphite	80-260	80-260
	Au	80-1000	
Au-HEG	Au	400-2000	800-2000

2.3.2 Selecting EPU Harmonics and Photon Polarization

Depending on the desired photon energy range, appropriate EPU harmonics shall be selected for optimum flux.



For circular polarization, only the first harmonics can be used. The energy range is from 100 eV to 1000 eV.

For linear polarization, use the first harmonics for energy between 80 eV and 1000 eV. Use the third harmonics for energy between 800 eV to 2500 eV. Use the fifth harmonics for energy above

1500 eV.

When the EPU control is established, the "Access" field shows "ENABLED"

The "Polarization" selection button can be used to select photon polarization. Available options are:

- Circular Left
- Circular Right
- Linear Horizontal
- Linear Vartical -
- Linear Vartical +
- Linear Inclined

For "Linear Inclined" polarization, the "Inclined Angle" must be set. The angle range is from -90° to $+90^\circ$.

2.4 Setup Chopper

Chopper must rotate in the positive direction for correct operation. If the chopper has moved in negative direction, it must rotate in positive direction for one revolution, so that the encoder can pass the home position reference mark and establish the correct angle readout.



LED indicators on the Chopper Controller in the NIM Crate show the status of the chopper system.

Once chopper is stopped, use "Move Relative" to adjust chpper position.


2.5 Opening the Beamline to the Ring

Once the setup of EPU and monochromator is done, follow steps below to open the beamline to the ring:

Open all vacuum valves

Go to 10ID POE ACIS panel. After injection, all five lights to the left shall be green. If any of the middle three lights is red, follow the "10ID POE Lockup Procedure" in Appendix.

Open the "Safety Shutter" by pushing the big green [Safety Shutters] button. The green light will be on.



From the REIXS user interface, open photon Shutter 2 (PSH2)

Open Photon Shutter 3 (PSH3)

Open Photon Shutter 4 (PSH4)

If any photon shutter cannot be opened, there could be some errors that trigger the interlock system.

Possible causes are: vacuum, cooling water, etc.

2.6 Getting the Beam into the Endstation

Verify apertures and slits along the beamline. Using the numbers below as a reference.

Variable Aperture:

Normal Mode: Horizontal Center 14.2 mm, Vertical Center -1.80 mm

Horizontal Gap 0 ~ 5 mm, Vertical Gap 0 ~ 5 mm

2-in-1 Mode: Horizontal Center 11.2 mm, Vertical Center -1.80 mm

Horizontal Gap 0 ~ 11.6 mm, Vertical Gap 0 ~ 5 mm

SM side Variable Aperture needs to be moved out, with
Horizontal Center -12.0 mm, Horizontal Gap 0.0 mm

4-Jaw #1:

Horizontal Center 0 mm, Vertical Center -1.15 mm


Horizontal Gap 20 mm, Vertical Gap 10 mm

4-Jaw #2:

Horizontal Center 0 mm, Vertical Center 0 mm

Horizontal Gap 20 mm, Vertical Gap 20 mm

Exit Slit:



Vertical Gap: usually 10 μm or 25 μm , range is 0 ~ 500 μm

Horizontal Center 0 mm, Horizontal Gap 20 mm

If the vertical gap of the exit slit needs to be changed, always approach the target value when closing the gap (gap decreasing).

2.7 Beamline Monitoring and Diagnostics

Several video cameras are installed along the beamline. The video can be accessed from any computer within the CLS intranet.

To access the camera inside 10ID POE Hutch, use the address below in a web browser:

<http://ccd1610-101.cs.csls.ca/>


To access the cameras along the beamline, use:

<http://v2e1610-101.cs.csls.ca/>

To access the video of endstations, use

<http://v2e1610-401.cs.csls.ca/>

There are multiple Beam Diagnostics Assemblies (BD1 to BD5) along the beamline to assist the beam alignment and to diagnose of the problem.



Each assembly may have several detectors attached, including photodiode for measuring the flux, gold mesh for measuring the drain current, and YAG Crystal: for visualizing the beam. The Beam Diagnostics can be inserted into or withdrawn from the beam by motorized linear feed-through.

To setup meters in the beamline, click the [Meters] button in the main user interface. The REIXS Meters Setup panel can be used to change most often used parameters for each meter.




2.8 Other

3.0 RSXS Endstation

3.1 Instruments

Please refer to the RSXS Endstation User Manual for details.

The picture below shows the instruments inside the rack for RSXS endstation.



The Ion Gauge controllers on top display the vacuum pressure in the scattering chamber and the load lock.

The NIM crate houses the electronics for pulse counting.

The SIS3820 Scaler is in the VME crate.

Verify the NIM crate and VME crate are working correctly.

Below the VME crate is the IOC1610-403 computer.

There are several Detectors in the RSXS chamber. To use each detector, following devices are necessary:


Chennaltron

High Voltage unit *NHQ203M*
Preamp Power Output unit *ORTEC 4003*
Discriminator unit *ORTEC 584*
Log/Lin Ratemeter *ORTEC 449-2*
Scaler *SIS3820*

Micro Channelplate

Micro Channelplate is 3300 Series from Quantar Technology

High Voltage unit *NHQ203M*
Position Analyzer *Model 2502A*
including Preamp / Processor Module and Rack Module
NIM-TTL-NIM Adapter *CAEN N89* or 4-8 Logic Fan-in Fan-out *CAEN N454*
Scaler *SIS3820*



TEY



Total Electron Yield is by measuring sample current.

- Current Preamplifier *SRS SR570*
- Voltage to Frequency Converter *NOVA N101VTF*
- NIM-TTL-NIM Adapter *CAEN N89* (optional)
- Scaler *SIS3820*

Photodiode



The photodiode in the scattering chamber is AXUV100EUT Photodiode from IRD. It can be used for reflectivity measurements, scattering measurements and fluorescence measurements. It is also very sensitive to visible light so the chamber has to be covered carefully when using photodiode. The HCG and RGA must be off as well.

Electrometer *Keithley 6514*

Sample Temperature Control

When the cryostat is on, the range of sample temperature in RSXS chamber is from 18K to 400K.

ARS Cryostat
Temperature Controller *Lakeshore 325*

Pumps

Scattering Chamber: turbo pump, cryopump, ion pump and rough pump

Loadlock: Turbo pump and rough pump

3.2 Motion Control

The diffractometer in the Scattering Chamber is driven by nine UHV stepper motors from AML. The motors are controlled by *spec* software.

The "RSXS Motors" panel can be used to setup and adjust parameters for all the motors in the Scattering Chamber.

	Two Theta th	Theta th	Chi chi	Phi phi	X x	Y y	Z z	Detector Z detz	Slit Wheel slit
	<input type="button" value="Stop"/>	<input type="button" value="Stop"/>	<input type="button" value="Stop"/>	<input type="button" value="Stop"/>	<input type="button" value="Stop"/>	<input type="button" value="Stop"/>	<input type="button" value="Stop"/>	<input type="button" value="Stop"/>	<input type="button" value="Stop"/>
MoveAbs (deg)	117.044	29.556	90.000	0.999	-0.310	-0.550	-0.300	0.000	252.000
Feedback (deg)	117.044	29.556	90.000	0.999	-0.310	-0.550	-0.300	0.000	252.000
MoveRel (deg)	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.000
HiLim (deg)	205.342	186.537	95.897	4.500	8.000	8.000	6.500	50.000	400.000
LoLim (deg)	-24.658	-23.463	85.897	-4.000	-8.000	-8.000	-7.250	-45.000	-400.000
	<input type="button" value="OFF"/>	<input type="button" value="OFF"/>	<input type="button" value="OFF"/>	<input type="button" value="OFF"/>	<input type="button" value="OFF"/>	<input type="button" value="OFF"/>	<input type="button" value="OFF"/>	<input type="button" value="OFF"/>	<input type="button" value="OFF"/>
Motor	<input type="button" value="ON"/>	<input type="button" value="ON"/>	<input type="button" value="ON"/>	<input type="button" value="ON"/>	<input type="button" value="ON"/>	<input type="button" value="ON"/>	<input type="button" value="ON"/>	<input type="button" value="ON"/>	<input type="button" value="ON"/>
Power	<input type="button" value="HARDW"/>	<input type="button" value="HARDW"/>	<input type="button" value="HARDW"/>	<input type="button" value="HARDW"/>	<input type="button" value="HARDW"/>	<input type="button" value="HARDW"/>	<input type="button" value="HARDW"/>	<input type="button" value="HARDW"/>	<input type="button" value="HARDW"/>
	<input type="button" value="SOFTW"/>	<input type="button" value="SOFTW"/>	<input type="button" value="SOFTW"/>	<input type="button" value="SOFTW"/>	<input type="button" value="SOFTW"/>	<input type="button" value="SOFTW"/>	<input type="button" value="SOFTW"/>	<input type="button" value="SOFTW"/>	<input type="button" value="SOFTW"/>
MaxSpd (deg/s)	3.000	2.000	0.895	0.382	7.500	7.500	7.500	7.500	30.000
Speed (deg/s)	2.000	1.000	0.358	0.286	5.000	5.000	5.000	5.000	30.000
BaseSpd (deg/s)	0.000	0.000	0.269	0.286	1.500	1.500	1.500	1.500	9.000
Accel. (sec)	2.000	2.000	1.000	1.000	1.000	1.000	1.000	1.000	1.000
Steps/Rev	1600	1600	200	200	200	200	200	200	200
EGU/Rev	0.200	0.200	8.952e-1	9.549e-1	0.500	0.500	0.500	0.500	3.000
Encoder (deg)	1.250e-1	1.250e-1	4.476e-1	4.775e-1	2.500e-1	2.500e-1	2.500e-1	0.002	0.015
Readback (deg)	1.250e-1	1.250e-1	4.476e-1	4.775e-1	2.500e-1	2.500e-1	2.500e-1	0.002	0.015
Deadband (deg)	0.001	0.001	0.001	0.001	0.005	0.005	0.005	0.005	0.015
Backlash	<input type="button" value="FAVORED"/>	<input type="button" value="FAVORED"/>	<input type="button" value="FAVORED"/>	<input type="button" value="FAVORED"/>	<input type="button" value="FAVORED"/>	<input type="button" value="FAVORED"/>	<input type="button" value="NO"/>	<input type="button" value="NO"/>	<input type="button" value="FAVORED"/>
Direction	<input type="button" value="DECREASE"/>	<input type="button" value="DECREASE"/>	<input type="button" value="INCREASE"/>	<input type="button" value="INCREASE"/>	<input type="button" value="INCREASE"/>	<input type="button" value="INCREASE"/>	<input type="button" value="INCREASE"/>	<input type="button" value="INCREASE"/>	<input type="button" value="INCREASE"/>
Backlash (step)	400	400	100	100	100	100	0	0	200

3.3 Data Collection Software

The data acquisition is through the *spec* software from *Certified Scientific Software*.

Please refer to *spec* manual and "*SPEC Macros for REIXS RSXS Endstation*" for more details.

4.0 XES Endstation

4.1 Instruments

Channeltron

MCP 2400

4.2 Motion Control

Sample X, Y, Z

Detector Rotation

Detector Tilt

Detector Translation
Spectrometer Rotation
Hex-Pot

4.3 Data Collection

Custom software