

## **Neuron Sensors Network Technology**

# **User Guide**

4\_10

**CAB** revision

© All rights reserved to EDU-LAB LTD.

The material in this book may not be copied, duplicated, printed, translated, re-edited or broadcast without prior agreement in writing from EDU-LAB LTD.

Edu-Lab Ltd, Karoo Close, Bexwell Business Park, Bexwell, Norfolk, PE38 9GA, UK Tel: +44 (0)1366 385 777 Fax: +44 (0)1366 386 535 Email: sales@edulab.co.uk Web: www.edulab.co.uk

EDU-LAB LTD

## Contents

Chapter	r 1 – Introduction	1
1.1 1.2 1.3 1.4 1.5 1.6 1.7 1.8 1.9	About data-loggers and sensors About logger sensors Plug and play system Modes of operation Data display and analysis About this User Guide Computer requirements Safety instructions Environmental conditions	1 2 3 4 5 5 6 6
Chapter	r 2 – Basic Setup	7
2.1 2.2 2.3 2.4 2.5 2.6	Installation Logger sensors main screen Connecting a sensor Sensor module box Color Temperature sensor module setup	7 7 8 9 9 9
Chapter	r 3 – On-line Experiment Mode	. 12
3.1 3.2 3.3 3.4 3.5 3.6 3.7	Setting up an on-line experiment Running an on-line experiment Running an on-line experiment with trigger Plotting an XY graph Single step mode Using photogates Video – record, play back and save with data	. 16 . 21 . 32 . 40 . 50 . 75 . 90
Chapter	r 4 – Off-line Experiment Mode	. 95
4.1 4.2 4.3	The off-line experiment mode menu Off-line experiment with battery module Off-line experiment mode with a PC	. 95 . 96 115
Chapter	r 5 – Using the RF Communication Modules	130
5.1 5.2 5.3	RF With a PC RF With a PC and Groups of sensors RF with the Monitor Display Unit (MDU)	130 130 131
Chapter	r 6 – Saving, Loading, Printing & Exporting to a Spreadsheet	132
6.1 6.2 6.3 6.4	Saving Loading Printing Exporting to a spreadsheet	132 132 133 133
Chapter	r 7 – Edu-Logger® <sup>™</sup> Tools	136
7.1 7.2 7.3 7.4 7.5 7.6 7.7	Set Sensor ID Number Set RF Module ID Number Restoring sensors' factory defaults Producing and saving instructions with and without the software configuration . Set significant figures Restore USB Language Select	137 137 138 140 148 148 148

Chapter	8 – Logger Sensors Modules	14	<b>19</b>
8.1	USB bridge module USB-100	14	49
8.2	Battery module Bat-100	14	49
8.3	RF Communication module RF-100	15	50
8.4	Digital display module VIEW-100	15	50
8.5	Graphic display module VIEW-101	15	50
8.6	Monitor display unit SX0	15	52
8.7	Voltage logger sensor NUL-101	15	53
8.8	Current logger sensor NUL-102	15	54
8.9	Temperature logger sensor NUL-103	15	54
8.10	Light logger sensor NUL-104	1	55
8.11	Oxygen logger sensor NUL-105	1	56
8.12	pH logger sensor NUL-106	1	57
8.13	Relative humidity logger sensor NUL-107	15	58
8.14	Heart rate & pulse logger sensor NUL-108	1	59
8.15	Photo gate logger sensor NUL-109	16	60
8.16	Pressure logger sensor NUL-110	16	64
8.17	Force logger sensor NUL-111	16	65
8.18	Sound logger sensor NUL-112	16	66
8.19	Motion logger sensor NUL-113	16	67
8.20	Magnetic logger sensor NUL-114	16	68
8.21	Conductivity logger sensor NUL-115	16	59
8.22	Spirometer logger sensor NUL-116	17	70
8.23	GSR logger sensor NUL-117	17	71
8.24	Colorimeter logger sensor NUL-119	17	72
8.25	CO <sub>2</sub> logger sensor NUL-120	17	73
8.26	Barometer logger sensor NUL-121	17	74
8.27	Blood pressure logger sensor NUL-122	17	75
8.28	Flow logger sensor NUL-124	17	76
8.29	Force plate logger sensor NUL-125	17	77
8.30	Rotary motion logger sensor NUL-126	17	78
8.31	Acceleration logger sensor NUL-127	17	79
8.32	Soil moisture logger sensor NUL-129	18	30
8.33	UVB logger sensor NUL-130	18	31
8.34	I urbidity logger sensor NUL-131	18	32
8.35	UVA logger sensor NUL-132	18	33
8.36	Surface temperature logger sensor NUL-133	18	34
8.37	Wide Range logger sensor NUL-134	18	35
Append	ix A – Edu-Logger® Modules	18	36
A.1	Accessory modules	18	36
A.2	Logger sensors modules	18	37

## Chapter 1 – Introduction

## 1.1 About data-loggers and sensors

In science experiments various values are measured, such as: temperature, light, sound, voltage, current, distance, velocity, acceleration, oxygen percentage in air or in solution, pH etc.

Most of the measured values are analogue, which change successively. The device that measures them is electronic, so the value is changed into an electrical signal. This device is called a sensor or transducer. There is a specific sensor for each value.

The electrical signal is then converted into digital form in order to be displayed or computerized by a component called ADC (Analogue to Digital Converter).

Digital numbers are changed in steps. The ADC resolution determines the size of the steps. The higher the resolution, the smaller the steps while the accuracy is higher.

A data-logger is a device that includes an ADC for reading values from several sensors and records this data in its internal memory. Usually, it has a display for viewing the measured values graphically. The data-logger can be connected to a PC for exporting the recorded data to Excel files or for saving it in files and viewing it on the screen.

The data-logger's capability is characterized by the number of sensors that can be connected to it simultaneously, its memory size, its ADC resolution level and its sampling rate. Its inputs are suited to all available sensors and also to sensors that will be developed in the future.

The data-logger software analyses-lab Ltd all available sensors.

#### 1.2 About logger sensors

The logger sensors system is rather different from almost all other educational data-loggers in that its sensors incorporate their own individually programmable microprocessors and have memory. Hence each sensor can be viewed as a data-logger in its own right, recording and storing data independently and can be referred to as logger sensors but within this guide they will just be called sensors.





Figure 1-1 Voltage-Current Measurement Experiment

To program the sensors they need to be connected to either a PC or to a logger sensors' Monitor Display Unit (MDU). They can be connected singly or in a chain. The order of connection does not matter and they can be added or removed from the chain without affecting the others.

Depending on their use, sensors can be powered directly from a PC or from the MDU. However, they can also be powered individually or in chains by a Battery Unit.

One chain of sensors can be divided into more chains by adding RF communication modules to all the sensors chains allowing remote connection of up to 30m (in open space).

10 bit resolution is provided for most sensors and sampling rates vary from 10000 per second to 1 per hour, depending on the sensor concerned. The rates available for each sensor match well to their likely use. Experiment durations are from 25ms to 31 days, depending on the sensor and the sampling rate. A Trigger setting with Pre-trigger is available through which sampling starts and data is displayed from just before the Trigger operated.

## 1.3 Plug and play system

The philosophy behind logger sensors is a "plug and play" system. Connect the experiment's required sensors to the PC (through a USB adaptor) or to a MDU (Monitor Display Unit) and perform the experiment with an intuitive software.

There is no need to study data-logger operation. All the measurements, the recording and analysis are done by the sensor itself according to its functions. The data sent to the PC or to the MDU is processed by the sensor.

The total memory of the system is increased, because each sensor has its own memory.

Each sensor has its own micro-controller (tiny computer), so it can control and adapt even the hardware to different functions. This is why many of the sensors have several ranges or different types of measurements that are usually done by more than one sensor.

The built-in software in the logger sensor can be upgraded at any time using software (without opening the sensor module).

#### **1.4** Modes of operation

The Logger Sensors system has two modes of operation: **On-line experiment** and **Off-line experiment**.

**On-line experiment mode** is where the sensors are connected to a PC or the MDU, programmed, and remain connected as the data is fed back continuously in real-time. All sensors gather data at the same rate, are triggered together (a choice has to be made of which sensor to trigger from) and run for the same time.

On-line experiment mode enables also collection of data at particular stages of the experiment instead of continuously. The **Single-step** mode (sometimes known as 'snapshot') can be used when measurements are taken while some of the experiment variables are changed manually, like changing the source voltage in an electrical circuit, volume in a pressure-volume experiment or adding drops into solutions, etc.

**Off-line experiment mode** is where the sensors have been programmed by connection to a PC or the MDU with different Sampling rates and Experiment durations. The sensors are then disconnected from the PC or MDU and put into battery-powered chains (or singly) to collect their data on the pressing of the sensors' Start/stop buttons. However, in this mode, each sensor can be programmed independently to be triggered to record data as well as having their Start/stop buttons pressed at different times. Five experiments can be stored in each sensor.

Reconnection of the sensors to a PC or the MDU is required to upload the data and analyse it. However, all graphs would be overlaid on the longest Time axis with t = 0s being where each sensor appears to have been triggered. Hence it does not take into account the time differences between the pressing of any sensors' Start/stop buttons and their independent triggering.

The sensors can remain connected to the PC at **Off-line experiment mode**. As before, the sensors can be programmed with different Sampling rates and Experiment durations. Again, each sensor can be programmed independently to be triggered to record data. Starting and stopping can be done either by clicking on relevant icons on the PC or by pressing the Start/Stop buttons on the sensors. Data is then uploaded in order to display and analyse it.

Collected data (in both modes) is stored in the sensors' internal memories, to be displayed as required on a PC or, in digital number format, on the MDU. Each sensor has an ID number that can be changed, when necessary, through software. This would be needed if you were using an array of the same type of sensor, i.e. several temperature sensors, and needed to see what was happening on each one. All sensors are automatically recognized by the system.

#### **1.5** Data display and analysis

The computer's display of data can be in the form of a graph, a table or both, plus a digital display of each sensor's current value e.g. 20 °C, 8.95V, 20.9%. The default graph display is of what the sensor measures (Y-axis) plotted against Time (X-axis) but, as mentioned in an earlier paragraph, it is possible to plot XY graphs in which one sensor's data is plotted against another's.

Graphs can have their axis scales pre-selected, can be zoomed to maximize their display in the Y-axis direction, or have small areas selected and zoomed to allow examination in more detail. The graphs of each sensor can be overlaid and their Y-axes moved to convenient positions on the screen.

A best-fit line/curve facility is provided to overlay the graphs, and extrapolation to zero of best-fit straight-line graphs is available. Areas under graphs can be easily calculated for use in determining such a quantity as Impulse from a Force-time graph. Graph plotting can be in 'points only' or in 'a joined up line' and a set of grid-lines can optionally be added.

A number of mathematical functions [log(A), ln(A), sqrt(A),  $A^2$ ,  $1/A^2$ , (A+B),  $(A \cdot B)$ , (A/B),  $(1000 \cdot A \cdot B)$ ,  $(A \cdot K)$ ,  $(e^A)$ ,  $(10^A)$  and (A/K)] are available with which to convert data, where A and B are the variables and K a constant. The function A/B would be useful in generating data (and a graph) of the resistance of a filament light bulb, where A is the voltage across the bulb and B the current flowing through it.

These functions allow one to deal with most data processing. The data can also be exported into a spreadsheet for further manipulation and processing.

Triggering is available to start data-logging when a particular sensor's measured value falls or rises below or above a set level. When triggering has been selected, a Pre-trigger display is made available on the graphs and in the tables so that sensor values just before the triggering took place can also be seen. This is particularly useful when looking at how the voltage across a coil changes when a magnet falls through it.

A video with play-back facility is provided via a webcam so that data can be collected synchronized with a video of the process. These can be saved for future use.

Worksheets, incorporating setup details, photographs and diagrams, can be developed and displayed, and saved with or without a Setup Configuration.

#### Note:

Additional features will be incorporated in later versions of Logger Sensors and more sensors are being developed. Revised versions of the software will be provided free to download as they become available.

#### **1.6** About this User Guide

This User Guide was designed to go through one complete chapter at a time.

Instructions to do things, click on a button or an icon, plug in a sensor etc., are indicated by a bullet  $\bullet$ .

Logger Sensors software is a Windows<sup>™</sup>-oriented program with dynamic windows that may, on occasions, hide each other. Moving the overlaid window (displayed as a table, graph, digital meter etc.), and then clicking on the previously hidden window, will bring the latter to the front. Alternatively, clicking on a window's icon or button can bring it to the front.

The various windows can be moved around the screen by clicking on their headlines (upper blue bars), holding the cursor down, and then dragging the window to the required position.

These windows can also be resized by clicking on their edges or corners until ---, -, -, or

are displayed, then hold down the cursor and drag to enlarge or shrink as required. Columns in table windows can also be changed in width and this is done by clicking on the column header until is displayed, hold down the cursor and drag to enlarge or shrink as required.

This manual describes the Edu-Logger® software through few sensors.

Because the measured data is processed in the logger sensor, the software treats each sensor in the same manner except the photo-gate sensor.

The Edu-Logger® software is very rich and at the same time intuitive. It is very simple to use.

There is another manual that describes available sensors and their typical experiments.

## 1.7 Computer requirements

Logger Sensors can be used with a PC running Microsoft<sup>®</sup> Windows 2000, NT, XP or Vista operating systems. It will **not** run with Windows 95, 98, 98SE or ME operating systems. It requires 20MB of hard disc space. There are no limitations of processor speed or RAM needed.

At least one USB port is required, together with another if a webcam or USB video camera is to be connected. Similarly a further port will be needed if a printer is to be connected.

Whilst no doubt a number of webcams and video cameras can be used, the Viewflex or Extravalue 1.3 Megapixel Snake Webcam – USB – No driver install required, Plug and Play, works fine with both Windows<sup>®</sup> XP and Vista. It is inexpensive and can also provide additional lighting if required.

If use is to be made of the export to a spreadsheet facility, then the PC must have a software such as Microsoft Excel<sup>®</sup>, Lotus 123<sup>®</sup>, Softmaker Planmaker<sup>®</sup>, OpenOffice Calc<sup>®</sup>, or similar CSV (Comma Separated Value) application.

If you wish to incorporate and save instructions for the use of activities, then access will be needed to a word processor or similar application to prepare such a file.

## 1.8 Safety instructions

- Read and understand these entire instructions before proceeding.
- Keep these instructions.
- This equipment is <u>not</u> designed to perform in an environment where failure may result in accident or injury.
- This equipment is designed to operate and perform in an educational environment. It is not designed to be used in a medical, marine or industrial environment.
- Disconnect module before servicing for any reason. Servicing should be performed by QUALIFIED PERSONNEL ONLY!
- Do not operate a module if it is broken or if components are revealed for any reason.
- Do not operate module with damaged cord, wires or electrical parts. Use only cords supplied with the module.
- Use only power source approved by EDU-LAB LTD or enclosed with the training system.
- The USB module is also an approved power source for the Edu-Logger® modules.
- Use modules only to measure the parameters for which they are designed as specified on the module label.
- Do not attempt to measure values exceeding those specified on the module label, particularly voltage and current.
- Clean only with a dry cloth.
- Do not install near any heat sources such as radiators, heat registers or other apparatus that produce heat.
- Unplug the module during lighting storms or when unused for long periods of time.

#### 1.9 Environmental conditions

- Do not expose modules to any kind of liquid.
- Operating temperature: 0°C to +40°C.
- Humidity: up to 95% at 35°C.



#### **Power requirements:**

Battery module (4 AA batteries) or computer USB outlet.

## Chapter 2 – Basic Setup

#### 2.1 Installation

The software and drivers must be installed before connecting any modules to the PC.

- Open the setup file on the CD you received with the system.
- Follow the instructions on the screen. The installation process is straightforward and the required drivers are installed automatically.

The installation is composed of two parts: Edu-Logger® software installation and USB driver installation. After the installation process is completed, the Logger Sensors software is ready to use.

#### Notes:

Upgrading the software can be done at any time. Installing the upgraded software just replaces the relevant files, so uninstalling the software before upgrading is not needed.

During upgrading the software the USB driver installation can be skipped by clicking the 'Cancel' button.

The Edu-Logger® shortcut icon is should appear on the PC desktop.

## 2.2 Logger sensors main screen

Double-click the Edu-Logger <sup>™</sup> shortcut icon <sup>™</sup> to display the following.



Note that the functions of each icon are indicated above but can also be seen by moving the cursor over each on-screen icon.

There are two icon bars on the main screen. The upper one, the **Main-icon bar**, includes the main Edu-Logger® functions:

Search for sensors On-Line experiment Off-Line experiment Tools About

The lower bar is the **Sub-icon bar**. Clicking on some of the Main-icon bar icons (**not** Search for sensors or Help) displays sub-icons relating to their associated functions.

#### 2.3 Connecting a sensor

- Connect the USB Bridge module *to a USB port on the PC.*
- Connect the Temperature sensor to the USB Bridge module using a Edu-Logger® short cable (the same procedure should be followed for any other sensor).

Each module has two identical cable sockets. You can use either. The sockets enable you to connect the sensor modules in a chain.

Click on the Search for sensors icon in the Main-icon bar. The program will scan and display the connected sensor's Module box, in this case that of the Temperature sensor, automatically in the Module window on the left-hand side of the screen as shown below.



The program starts automatically in the **Search for sensors** mode so, if you run the program whilst modules are connected to the PC, their presence will be scanned for and their Sensor Module boxes will be displayed in the Module window.

If you add or remove sensors from a chain of sensors you must rescan by clicking on the

## Search for sensors icon

Any newly detected Sensor Module boxes are displayed vertically in the Module window.

Each Sensor Module box is characterized by several parameters as shown in the following Section 2.4.

This is the "plug and play" method. Connect the sensors to the PC through the USB module, click 'Search' and the system is ready to use.

## 2.4 Sensor module box



- 1. **ID** Displays the sensor's ID. Up to 9 sensors of the same type can be connected in a chain.
- 2. **Value** Displays the sensor's numerical value.
- 3. **Module Setup** Opens the sensor's setup window.
- 4. **Type** Displays the sensor's type (Light, Temperature, etc.).
- 5. **Units** Displays appropriate units of measurement (Ix for a Light sensor, °F or °C for a temperature sensor, etc.).
- 6. **Colour** Shows the graph line colour of the sensor; this colour can be easily changed (see below). The numbers and units of the Y-axis will also be displayed in this colour.

### 2.5 Color

 Click on the Temperature Sensor Module box's colour button
 The following Colour box then appears:

Basic colors:
Custom colors:
Define Custom Colors >>
OK Cancel

• Select one of the colours by clicking the cursor on it and then click the **OK** button.

The Colour window disappears and you will see that the Colour box of the Temperature Sensor Module box has changed to your selected colour. This will also have changed the color of the graph line and Y-axis numbers' and unit's color, although you will not be aware of it at this stage.

 Click the cursor back on the Temperature Sensor Module box's color button, select the original red color, and then again click the OK button to return the system to its original state.

Note: The selected colour is for that particular sensor.

## 2.6 Temperature sensor module setup

 Click on the Module Setup button and of the Temperature Sensor Module box to display the Temperature 1 – Options tab window.

Temperature 1 - Options	×
Options Graph	
<u>Display</u> ☐ Table ☐ Digital ☑ Graph ☐ Grid	<ul> <li>Celsius</li> <li>Fahrenheit</li> </ul>

Two selection tabs – **Options** and **Graph** – are shown here, but an additional one – **Trigger** – will be available on some occasions, together with **Experiment duration** and **Sampling rate** when Off-line experiment mode is being used instead of On-line experiment mode.

#### 2.6.1 Options tab setup

The Options are:

#### <u>Display</u>

- Table When this box is ticked a table of data from the experiment will be displayed.
- **Digital** When this box is ticked the sensor's value will be displayed in a large digital form.
- **Graph** When this box is ticked data from the experiment will be displayed in graph form (default setting).
- **Grid** When this box is ticked a series of grid-lines will appear on the graph.

These four options appear for the vast majority of sensors.

On the right-hand side are two radio buttons which, for the Temperature sensor, allow the units of °C or °F to be selected by clicking on them. Other sensors will have radio buttons from which different ranges and functions can be selected.

Many of the sensors have several ranges that can be selected similarly. Some of the range selecting changes the internal electronics like amplification, filtering etc. This can done by the computer abilities of each sensor.

Click the cursor in the box alongside Digital. A ✓ sign appears in the box and a Digital window opens as shown below.



As before, this Digital window can be dragged to any desired position on the screen by clicking the left cursor on its headline, holding it down and dragging as required. It can also be changed in size and indeed in some circumstances it may be useful to have it almost fill the whole screen. You can also close it by clicking on its **Close button** but again **do not** do so at this stage.

- Now change the temperature units from °C to °F by clicking on the Fahrenheit radio button and see how it affects the value displayed on the Digital window and on the Sensor Module box.
- Change the temperature units back to °C by clicking on the **Celsius radio button**.

Such a change of units will be stored in the sensor's built-in memory and so whatever units are chosen, they will remain until a subsequent change is made.

- Warm the metal rod of the temperature sensor and see how that affects the value displayed on the **Digital** window and on the Sensor Module box.
- Now click on the Close buttons of each of the Table, the Temperature 1 Options tab and Digital windows, to close them.

#### Note: At this stage no data from the sensor will have been recorded.

- Now click the Close buttons at the right-hand ends of both the Sub-icon bar and the Main-icon bar.
- Finally disconnect the Temperature sensor sensor with its short cable from the USB Bridge module and the USB Bridge module from the PC.

## **Chapter 3 – On-line Experiment Mode**

In this chapter you will run an experiment with the temperature sensor in order to see and use the features of the **On-Line Experiment** mode. In this mode, the experiment is controlled by the PC. Results are recorded and displayed in real-time.

- Connect the USB Bridge module 🌋 to a USB port on the PC.
- Connect the Temperature sensor short cable.
- Double-click the Edu-Logger® shortcut icon <sup>22</sup> to display the following.



 Click on the On-Line Experiment icon
 to display the additional Sub-icon bar labelled below.

Note that the functions of each icon on the Sub-icon bar are indicated here but can also be seen by moving the cursor over each on-screen icon.



EDU-LAB LTD

 Click on the Module Setup button to display the Temperature 1 – Options tab window and ensure that the units of temperature radio button is set to Celsius and that the box alongside Graph is ticked. It should appear as below.

Temperature 1 - Options	×
<u>Display</u> ☐ Table ☐ Digital ☑ Graph ☐ Grid	<ul> <li>Celsius</li> <li>Fahrenheit</li> </ul>

The digital function was explained in the previous chapter.

■ Click the cursor in the box alongside **Table**. A ✓ sign appears in the box and a Table window opens as shown below.

🎲 Ter	nperature 1	
Time	Temperature	

Throughout the experiment this sensor's values **alone** will be shown in this table. If other sensors had been connected, their **individual** tables could be displayed or, as you will see later, a table can be displayed with all the sensors' values.

Columns in Table windows can also be changed in width and this is done by clicking on the column header until  $\neg$  is displayed, holding down the cursor and dragging to enlarge or shrink as required.

• Close this table by clicking again in the box alongside **Table**.

 Click in the box alongside Grid and note how a grid of feint lines is superimposed on the graph window as shown below.

Module window Temperature 1 19.3 °C	NeuLog <sup>14</sup> , Ver: 2.43	euron Sensors Network Technolog	, ý († († († († († († († († († († († († (†	2. % * • * *	= ** **	
	Next.60		Temperatur Opione G Tispla Tia Di Di Di Gri	e 1 - Options aph Y ble gifal aph d	Coliss Fatenhol	
	10 0 -10 -20 0	1 2 3	4 5 Time [s]	6 7	8 9	10

• Click on the **Graph tab** to display the Temperature 1 – Graph tab window shown below.

Т	emperature 1 - Graph		X
	Options Graph		
	Y max.	110 🚖	
	Ymin.	-25 🚖	
	Y-axis position	<b>1</b>	
	Ø** 😽		
	200m nt Color		

The **Y** max and **Y** min fields set the maximum and minimum Y-axis values of a graph display. Initially they are set for the total range of each sensor.

- Click on the Y max down-arrow ▼ and change the value to 50.
- Click on the Y min up-arrow ▲ and change the value to 0.

Note how the range of the Y-axis (Temperature) now goes from 0 ℃ to 50 ℃.

By clicking on the **Y-axis position arrows** I the position of the Y-axis can be moved.

Click on the Y-axis position arrow and change the value to 30.

Note how this has moved the Y-axis quite a way to the right. This is useful when there are many sensors in use as each will have its own Y-axis and it would not be helpful to have them overlap each other.

Clicking on the **Colour icon** copens the Colour window where the sensor's graph line colour can be changed. It does the same job as the **Colour button** on the Sensor Module box that you will have dealt with in Chapter 2, Section 2.5.

Click on the Colour icon 💁 . Select a green colour and then click the OK button.

You should now have in view something much like the screen displayed below with a number of icons on the Sub-icon bar and some more on the graph itself.



You will have noted the **Zoom-fit icon**  $p^{\text{res}}$  in the Temperature 1 – Graph tab window. This operates **only** when a graph is displayed, and so cannot be used here. Clicking on it adjusts just the Y axis s so that the graph's minimum and maximum values just fit within it. This **Zoom-fit** relates only to the sensor that you are dealing with – in this case that would be Temperature 1.

You will have noted other icons  $\boxed{2}$   $\boxed{2}$   $\boxed{2}$   $\boxed{2}$   $\boxed{2}$  just above the graph. The first three of these also deal with Zoom functions and they, along with the other three, will be dealt with in Chapter 3 - Section 3.2.

- Again click on the Colour icon Markov, select a red colour like what you saw originally, and then click the OK button.
- Return the Y-axis to –7 position.
- Now close the Temperature 1 Options tab window by clicking on the Close button in its top right-hand corner.

## 3.1 Setting up an on-line experiment

#### 3.1.1 Options for on-line experiment setup

In this mode the sensors always start with a series of default settings and a cleared dataset regardless of what was previously in their memories. This is very different to the situation with Off-line mode (see Chapter 4) where the sensors retain their last settings.

Before running an experiment you will need to set the **Experiment duration**, **Sampling rate** and **Participating sensor(s)**.

 Click on the Experiment setup icon in the Sub-icon bar to open the Experiment setup – Options tab window within the Graph window as shown below.

Experiment setup - Options		×		
Options Graph Trigger XY	Graph			
<u>Display</u>				
□ Table ☑ Graph				
Experiment duration	🗸 10 seconds	•		
Sampling rate	🗸 10 per second	•		

You will be making use of the Graph option and, by default, it displays as enabled by a tick  $(\checkmark)$  in the box alongside **Graph**.

 Click the cursor in the box alongside **Table** – a tick will appear in it and the Table window shown below also opens within the Graph window.

Time	Temperature 1	Manual values	
rime	Tremperature 1	Mailuai values	

The data from the Temperature 1 sensor will be displayed in such a table. If other sensors were connected, columns would be generated for their data too. This is very different from the table created for each separate sensor which would be produced via the **Module Setup** 

#### button 🧐

Note also a 'Manual values' column. When relevant, values not obtainable from connected sensors can be typed in here. If you were conducting an experiment investigating how the Volume of a fixed mass of air changed with Pressure, then the values of Volume would be entered into the table and the title 'Volume' and appropriate units would be typed into the column's header.

The columns in Table windows can also be changed in width and this is done by clicking on the column header until  $\checkmark$  is displayed, then hold down the cursor and drag to enlarge or shrink as required. The Table window, and indeed any other window, can be moved about on the screen by clicking the cursor on its upper blue headline, holding it down, and then dragging to the required position. So if there is an overlap with the Experiment setup – Options tab window, move one or other so that both are visible.

Click on the **Close button** in the top right-hand corner of the Table window to close it.

The **Experiment duration** defines the time from when sampling (data logging) starts until sampling ends. This duration can only be selected from a pre-selected list.

 Click on the down-arrow alongside Experiment duration to display the available durations and select '30 seconds' by clicking on it.

Experiment setup - Options		×
Options Graph Trigger XY	Graph	
<u>Display</u>		
□ Table ☑ Graph		
Experiment duration	🗸 10 seconds	•
Sampling rate	<ul> <li>10 seconds</li> <li>30 seconds</li> <li>1 minute</li> </ul>	<ul> <li>Image: A second s</li></ul>
	- 2 minutes	
	<ul> <li>3 minutes</li> <li>4 minutes</li> </ul>	
	- 5 minutes	
	🗕 6 minutes	*

The **Sampling rate** defines how many measurements per second (per minute or per hour) are taken. The program may select a sampling rate automatically according to the experiment duration. Like the Experiment duration, the Sampling rate can only be chosen from a pre-selected list. Indeed, some of the sampling rates are blocked, they being dependent on the Experiment duration or the sensor in use. For example, the Temperature sensor cannot sample at a faster rate than 100 samples per second (there is no reason to go higher) so the higher rates are blocked even if another sensor is connected which can sample more quickly.

 Click on the down-arrow alongside Sampling rate to display the available rates and select, if not already the case, '10 per second' by clicking on it.

Some of the higher sampling rates may be blocked when experiment durations are long. This is because the maximum number of samples that can be taken and stored per module is 64,000 samples. High rates and long durations would obviously require high numbers of samples.

#### 3.1.2 Graph setup for on-line experiment

• Click on the **Graph** tab in the **Experiment setup** window. The following window appears.

	×
aph	
Temperature 1	•
<u>.</u>	
	aph Temperature 1

The **Participating sensor(s)** box displays a list of the sensors connected to the PC – in this case just one. If the data from a sensor is to be used in an experiment, then the box alongside its name will need to be ticked, as indeed it is here for Temperature 1. Un-ticking a participating sensor's box will result in its data not being plotted on the graph.

In some experiments more than two sensors can be in use but one may wish to only display the Y-axes of two of them. Then the box alongside **Display 2 axes only** would need to be ticked. The two Y-axes you wish to be displayed would then be selected by the down-arrows alongside. With just the one sensor in use you can obviously only have one Y-axis displayed so you might as well disable the selection of **Display 2 axes only**. Regardless of ticking the box alongside **Display 2 axes only**, the graph lines of **all** the participating sensors will be plotted.

 Click the cursor in the box alongside **Display 2 axes only** to disable its selection and change the Experiment setup – Graph tab window to that shown below.



#### 3.1.3 Trigger setup for on-line experiment

Click on the Trigger tab in the Experiment setup window to display the Experiment setup

 Trigger tab window shown below.

Experiment setup - Trigger	×
Options Graph Trigger XY Graph	
<b>Trigger</b>	

There are occasions when you want sampling to commence **only** after a certain event has occurred. With a Force sensor this might be when the force exceeds 2.00N or with a Light sensor when the Illuminance falls below 1000 lx. With a Temperature sensor it would be when the temperature has risen above or fallen below a particular value.

It is the **Trigger** that enables you to start the data logging **only** after the event occurs.

 Click in the box alongside the word **Trigger** so that the Experiment setup – Trigger tab window displays as below.

Options Graph Trigger XY Graph
<b>∨</b> Trigger
Sensor Temperature 1 💌
Level 0 *C
Mode Rise 🔽

If left as displayed in the above window, the sampling would start when the temperature rises above  $0^{\circ}$ . The **Level** in this case sets the temperature at which triggering occurs. **Mode** allows you to select whether it is a '**Fall** below' or a '**Rise** above'.

You can Trigger the start of data logging with any of the sensors that have been connected, selecting as appropriate via the down-arrow alongside the sensor's name. Here, of course you only have a temperature sensor connected.

- Click in the box alongside Level, delete the '0' and type in '30'.
- If the box alongside Mode does not already indicate 'Rise', click on the down-arrow and select 'Rise'.

Then, (**do not do it**) following the **Run experiment icon** \$ being clicked, data logging would be triggered when the Temperature 1 sensor registers a rise of temperature taking it above 30 °C. A triggered setup will be used in running an experiment with a Force sensor in Section 3.4.

Click in the box alongside **Trigger** to disable this mode.

#### 3.1.4 XY Graph, and selecting the units of time on the X-axis

In the simplest situations you will use just one sensor and the graph will then just be of what it senses and records against Time. In the experiment which follows this section you will have a graph of Temperature displayed on the Y-axis and Time on the X-axis. However, there are many occasions when you will be using two (or more) sensors such as Pressure and Temperature together, or Current and Voltage together. Then the interest would be in seeing how the data they record relate to each other, i.e. might they be directly proportional to each other, or may they not?

As only one sensor has been connected at the moment, the **XY Graph** mode cannot be put to use. However, you will see below how it is set up and you will make use of it in Sections 3.4 and 3.5.

Click on the XY Graph tab to display the Experiment setup – XY Graph tab window below.

xperiment setup - XY Graph	X
Options Graph Trigger XY Graph	
X-axis	
• Time	
Seconds 📃	
C Sensor	

In this situation you will see that the X-axis radio button is set to **Time. Temperature** then, by default, would be plotted on the Y-axis.

It is here that you can change how the X-axis records Time, regardless of whether or not you wish to plot an XY graph. It can be set to always display in Seconds or, by clicking the **down-arrow** in the box below the **Time** radio button, it can be displayed in Days, Hours, Minutes and Seconds. Be careful on interpreting **Time** such as '1:20 (Minutes)'. It means '1 minute and 20 seconds' and **not** '1 minute and 20/100ths of a minute'. You can swop between these styles of display after the graph(s) have been plotted. Leave it set on 'Seconds'.

If a Pressure sensor had also been connected, you could choose which axis on which to plot Pressure or Temperature. Although you cannot do it at this stage, if a Pressure sensor had been connected, and you had clicked the radio button alongside **Sensor**, the following window would be displayed.



This would indicate that the Temperature would be plotted on the X-axis with, by default, the Pressure along the Y-axis. However, by clicking on the **down-arrow** alongside Temperature 1, Pressure 1 would be displayed and that could be chosen instead for the X-axis and, by default, Temperature 1 would then be plotted on the Y-axis.

 Click on the Close button in the top right-hand corner of the Experiment setup – XY Graph tab window to close it.

## 3.2 Running an on-line experiment

NOTE: Take care. In this experiment you will be using hot water.

- Pour about 200ml of hot water at around 60 ℃ into a 250ml beaker.
- Pour about 200ml of cold water at a temperature below that of the laboratory/room into another 250ml beaker.

You have already setup the **Experiment duration** to 30 seconds and the **Sampling rate** to 10 per second so all is now ready to go and you should see the following screen displayed.



21

At this stage the X-axis only displays as far as 10 seconds but, on running the experiment, this axis will 'expand' to accommodate the Experiment duration that you had set.

- Click on the **Run experiment icon** in the Sub-icon bar to start the datalogging.
- Insert the Temperature sensor rod into the hot water, watch the graph, note the increasing value of Temperature on the Temperature Sensor Module box and that data is being recorded in the Table.
- After about 15 seconds, take the Temperature sensor rod out of the hot water and insert it into the cold water.

After 30 seconds a message **'Experiment Completed'** appears and you should get a graph similar to the one shown below.

**NOTE:** If you had wished to stop the data logging before the 30 seconds had elapsed, you could have clicked on the **Stop experiment icon** in the Sub-icon bar.



Click on the Experiment setup icon in the Sub-icon bar to display the Experiment setup – Options tab window below.



Click in the box alongside **Table** to tick (✓) it and you should then see a Table of data from the experiment similar to that shown below.

🦻 Tab	le	
Time [s]	Temperature 1 [°C]	Manual values
0.00	10.8	
0.10	10.8	
0.20	10.8	
0.30	10.8	
0.40	10.8	
0.50	10.8	
0.60	10.9	
0.70	11.2	
0.80	11.6	
0.90	11.9	
1.00	12.5	
1.10	13.1	
	10.0	

If you wished to see what the temperature was at a particular time you would either click on the **Up/Down scroll buttons**  $\checkmark$  **v** or click and drag the **Movable scroll button**  $\blacksquare$  of the Scroll bar upwards or downwards.

 Click on the Close button in the top right-hand corners of both the Table window and the Experiment setup – Options tab window to close them.

You will have noticed in Chapter 2 that the Graph window has a number of icons, pictured and labelled below, just above the graph. You will now explore their use. Their purpose can also be seen by passing the cursor over each icon on the screen.



Zoom window	<ul> <li>Enables one to zoom part of the graph. The zoom functions change the Y maximum and the Y minimum of the displayed graph(s) together with the X-axis.</li> </ul>
Unzoom-all	<ul> <li>Returns the graphs to maximum zoom-out (the original view).</li> </ul>
Zoom-fit	<ul> <li>Adapts the Y maximum and Y minimum of the Y-axes only to just above and below the displayed graph's maximum and minimum values respectively. It has no effect on the X-axis</li> </ul>
Show cursors	<ul> <li>Displays two cursors that enable analysis of the graph at certain points or in a range.</li> </ul>
Show functions	<ul> <li>Enables one to view and operate several functions on a graph.</li> </ul>
Change to point or line graph	<ul> <li>Toggles between a points plot and a line graph</li> </ul>

The following steps will show you what these icons allow you to do.

#### 3.2.1 Zoom functions

Click on the Zoom window icon . Locate the mouse cursor at a point above the graph, press the left mouse button and keep it pressed, move the mouse cursor to the right and down below the graph. A rectangle should appear as shown below.



 Now release the mouse button. The part of the graph trapped in the rectangle should now be displayed in the Graph window as shown below.



Note that both the X-axis and Y-axis have been changed to the extent of the rectangle.

- Click on the Change to point/line icon and note its effect. Now click again on the icon to return the graph to its original state.
- Click on the Unzoom-all icon . This should have returned the graph to the original view.
- Click on the Zoom-fit icon to obtain a graph similar to that on the next page.



Note that this has expanded the graph vertically so that the Y-axis scale maximum is just above the highest value of temperature obtained, and the Y-axis scale minimum is just below the lowest value of temperature obtained.

• Again click on the **Unzoom-all icon** (a) to return the graph to the original view.

When more than one graph is displayed (i.e. when using several sensors) **Zoom-fit** for a particular sensor can also be performed.

 Click on the Module Setup button of the Temperature Sensor Module box to display the Temperature 1 – Options tab window shown below.

Temperature 1 - Options	×
Options Graph	
<u>Display</u> ☐ Table ☐ Digital ☑ Graph ☐ Grid	<ul> <li>Celsius</li> <li>Fahrenheit</li> </ul>

Click on the Graph tab to display the Temperature 1 – Graph tab window shown below.



- Now click on its Zoom-fit icon prime and you should see that this has performed the same task as the other Zoom-fit icon and .
- Close the Temperature 1 Graph tab window by clicking on the Close button in its top right-hand corner.
- Again click on the **Unzoom-all icon** (a) to return the graph to the original view.

#### 3.2.2 Using cursors

The **Show cursors icon** enables you to put two movable (horizontally) vertical cursors (lines) onto the Graph window.

Click on the Show cursors icon to display the graph with the two cursors shown below.



The values at the bottom left-hand corner of the graph indicate (for the graph above):

X1=7.28s and Y1=53.56°C – values where the first purple cursor intersects the graph line. X2=17.27s and Y2=59.3°C – values where the second green cursor intersects the graph line.  $\Delta$ X=9.99s – the difference between X1 and X2.  $\Delta$ Y=5.73°C – the difference between Y1 and Y2.

If you had more than one graph line (resulting from using more than one sensor) then their cursor values and differences would be obtained by clicking on the down-arrow alongside **Temperature 1** and selecting accordingly.

Each cursor can be moved to the left or right by clicking on them with the left-mouse button, holding the button down, and dragging to the position required.

Whilst this is useful to determine values and differences between values on graphs, you can also use it to find the **Gradient** at different positions on graphs.

 Move the cursors to a position similar to those on the graph displayed below, such that ΔX=1s.



Noting that, on the example graph above,  $\Delta Y=2.25$  °C it is easy to see that the Gradient ( $\Delta Y/\Delta X$ ) at that position is 2.25 °C/s. More precise values around a point would be obtained by making  $\Delta X$  smaller. With  $\Delta X=1$ s it is an easy calculation for you to make. However, if  $\Delta X$  had been 0.27s, the calculation would have been more difficult and may have had you reaching for a calculator. Edu-Logger®<sup>TM</sup>, however, provides a series of useful functions, one of which is to automatically calculate a Gradient.

## 3.2.3 Using functions: linear-fit, calculate function, statistics and math

 Move the first purple cursor to a position of interest on the graph and then click on the Show functions icon 4 to display the Functions tab window below.



The functions available are:

Linear-fit	-	This overlays a best-fit straight line on the selected graph and provides its equation in the form of $Y = mX + c$ . Extrapolation to zero is also available
Linear-fit	_	This overlays a best-fit straight line on the section of the selected graph
(between cursors)		between the cursors and provides its equation in the form of $Y = mX + c$ .
		Extrapolation to zero is also available
Area	_	This calculates and displays the value of the area under the selected graph
		between the cursors and the line Y = 0 (Not necessarily the X-axis).
Gradient [cursor	_	This calculates and displays the value of the gradient of the selected graph
1]		at the point of intersection with the purple cursor 1.
Polynomial	_	This overlays a polynomial-fit of selectable order between the cursors of the
[between cursors]		selected graph, and displays its equation

 Click on the down-arrow alongside 'Linear-fit' and then select 'Gradient [cursor 1]' from the drop-down menu shown below.

Functions		×
Functions	Statistics Math	
Tempera	ature 1	•
Li	near-fit	-
	near-fit	
	near fit (between cursors)	
	rea radient (cursor 1)	
P A	olynomial (between cursor:	s)

You will then see the Functions tab window displayed with a **Calculate function icon** () as shown below.



28

Click on the Calculate function icon Point Chosen (the intersection of the purple cursor with the graph line) is indicated as follows:



Linear-fit, Linear-fit [between cursors] and Area, are not Functions appropriate to use here but will be dealt with in Section 3.3 (Area) and 3.5 (Linear-fit).

Click on the Statistics tab and then the Calculate statistics icon and you will have displayed, as below, the Maximum, Minimum and Average temperatures recorded, together with their Standard deviation.



Click on the Math tab and then on the down-arrow alongside 'log[A]' to view the range of operations available as shown below.

Math				×
Functions	Statistics	Math		
А				
Tempera	ature 1		-	
В				
Tempera	ature 1		-	
log(A)			-	
log(A)			^	
In(A) sart(A)				
A^2			-	
1/A 1/A^2				
A+B			(278)	
A-B			~	

• Use the scroll bar to see the rest of the operations available.

With the **Temperature-Time** graph that you have on screen there is no useful mathematical operation that can be performed on its data, so do **not** implement a choice at this stage. However, as you will note from the range of operations available, it could be used to multiply the data from a Voltage sensor (A) by that from a Current sensor (B) so that a graph can be plotted of Electrical Power against Time or, with those same two sensors, to divide the data from a Voltage sensor (A) by that from a Current sensor (B) to provide a graph of Electrical Resistance against Time. This will be looked at later in Section 3.4.

- Click twice on the Close button in the top right-hand corner of the Math tab window to close it.
- Click on the **Show cursors icon** to remove the cursors from the graph.

There are occasions when it is helpful to retain the graph from an experiment on the screen, re-run the experiment, and have its graph displayed also. To do this you make use of the

Freeze current graph/s facility accessed via its icon

- First click on the Freeze current graph/s icon is on the Sub-icon bar.
- Place the Temperature sensor's rod back into the beaker of cold water.
- Click the on the Run experiment icon sin the Sub-icon bar.



After 30 seconds you should see two graphs similar to those displayed below.

Such overlays of graphs can be repeated as often as required.

If you wished to erase an existing graph, or graphs, before re-running an experiment, then you would click on the **Clear experiment results icon** in the Sub-icon bar.

- Click on the Clear experiment results icon and note that both of the graphs have been erased.
- Click on the Close button in the top right-hand corner of the Graph window to close it.
- Now click the Close buttons at the right-hand ends of both the Sub-icon bar and the Main-icon bar.
- Unplug the Temperature sensor I from the USB Bridge module but leave the latter plugged into the PC.

## 3.3 Running an on-line experiment with trigger

In the experiment so far you have used the **Run experiment icon** for start the data logging. However, as you will have noted in Section 3.1.3, there are occasions when you want sampling to commence only after a certain event has occurred.

Here you will be using a Force sensor. If you do not have a Force sensor available, but wish to see how to operate the Trigger, plug the Temperature sensor back into the USB Bridge module and follow the instructions in Section 3.1.3.

- Arrange a runway (shelf) at an angle of around 20° and place a brick at its lower end.
- Place the Force sensor flat up against the brick so that its hook is facing up the runway. You will need to move the brick into the middle of the runway to do this.
- Connect the Force sensor module Section to one of the sockets of the USB Bridge module
   Ising a Edu-Logger®<sup>™</sup> short cable.
- Double-click the Edu-Logger® shortcut icon to display the window shown below.



Click on the Module Setup button force Sensor Module box to display the Force 1 – Options tab window shown below.



• Check that **just** the box alongside **Graph** is ticked.

- If not already selected, click on the **10N** radio button.
- Click on the 'Upwards' force icon so that an impact force will register as positive. You will see an OK briefly displayed.

These settings will remain in the sensor's built-in memory until they are changed again.

Click on the Calibrate and zero icon
 Again an OK will be briefly displayed. Note that the value of the force displayed on the Force Sensor's Module box is now 0.00N or near.

Again this setting will remain in the sensor's built-in memory until it is changed again.

- Click on the Close button in the top right-hand corner of the Force 1 Options tab window to close it.
- Now click on the On-line experiment icon in the Main-icon bar to display the Graph window shown below.

Module window Force 1 0.00N	NeuLo	OG <sup>TU</sup> Neur	on Sensors Netwo	ork Technology	ي مرا ا	(2) ( () () () () () () () () () () () () ()		🧇 🚰 🔐	- 57		
	Next op			, 3		Ś Time [s]	Б	17		, 9	10

Click on the Experiment setup icon icon icon to open the Experiment setup - Options tab window displayed below.

xperiment setup - Options					
Options Graph Trigger XY	Graph				
Display					
□ Table ☑ Graph					
Experiment duration	🖌 10 seconds 📃				
Sampling rate	🖌 10 per second 📃				

- Click on the down-arrow alongside Experiment duration and select '1 second'.
- Click on the down-arrow alongside Sampling rate and select '100 per second' the fastest allowable rate with that duration.
- 34
- Now click on the **Trigger tab** to open the Experiment setup Trigger tab window below.



 Click in the box alongside the word **Trigger** so that the Experiment setup – Trigger tab window displays as shown below.

Experiment setup - Trigger	×
Options Graph Trigger XY Graph	
✓ Trigger	
Sensor Force 1	
Level 0 N	
Mode Rise 💌	
2	

 Click the cursor in the box alongside Level, delete '0' and type in '1'. Leave the Mode as 'Rise' so that it displays as shown below.

Experiment setup - Trigger 🛛 🛛 🔀
Options Graph Trigger XY Graph
✓ Trigger
Sensor Force 1
Level 1 N
Mode Rise 💌

 Click on the Close button in the top right-hand corner of the Experiment setup – Trigger tab window to close it.

It is now set up so that, after initially clicking on the **Run experiment icon** for a data logging will only begin once the force rises above 1N.

Click on the Run experiment icon in the Sub-icon bar.

The message 'Waiting for Trigger' will appear on the screen.

• Position a model car about 30cm up the runway. Then release it letting it collide with the Force sensor at the foot of the runway, but catching it as it bounces off.

You should now have displayed a Force against Time graph similar to that shown below.



You will see that there is a new button, the **Pre-Trigger button**  $\leq$  just below the '0' of the X-axis. Clicking on this button will allow you to see what was 'happening' at, and just before, the Trigger was activated.

Click on the **Pre-Trigger button** *H* to display a slightly extended graph as shown below.



**Do not do it**, but clicking on the **Unzoom-all icon** would return you to the original graph.

Click on the Zoom-fit icon to enlarge the graph as shown on the following page.

ieuLog 10 20 1 14 - 3 -												
2 - 1 -												
a -	-0.1	ů	0.1	0.2	0.3	0.4 Time	0.5	0.6	0.7	0.8	0.9	1.

### 3.3.1 Area under graph

With such a Force against Time graph the **Impulse** on impact can be calculated from the area under the section of graph (that shows the impact) down to the line where the Force = 0.00N.

Click on the Show cursors icon and move the two cursors so that the purple cursor is positioned where the force just starts to rise, and the green cursor where the force has returned to zero, as shown in the graph below.



• Now click on the **Show functions icon**  $\overset{\textcircled{}}{43}$  to display the Functions tab window below.



Click on the down-arrow alongside 'Linear-fit', select 'Area' from the drop-down menu, and then click the Calculate function icon to display the value of the area as shown below.

Functions				X
Functions	Statistics	Math	Ĩ	
Force 1				-
				-
			- 1	
Area= 0.0	09 [Ns]			
	9			

- Click on the Close button in the top right-hand corner of the Functions tab window to close it.
- Click on the **Show cursors icon** again to remove the cursors from the graph.
- Click on the **Unzoom-all icon** or return the graph to the original view.
- Click again on the **Pre-Trigger button** <sup>1</sup> to get the whole of the 'blip' on the graph.

Rather than just enlarge all of the graph in the Y-axis direction, which is what the **Zoom-fit** icon does, you can also enlarge just the section of interest using the **Zoom-window** icon does.



As before click on the Show cursors icon and move the two cursors so that the purple cursor is positioned where the force just starts to rise, and the green cursor where the force has returned to zero, as shown in the graph below.



As before, you could now click on the **Show functions icon** 3, select 'Area', click on the **Calculate function icon** and again would be displayed the value of the area under the graph between the cursors and the line Force = 0.00N.

This value represents the sum of work done by the force pulse  $-W = F \cdot dt$ 

Click on the Show cursors icon to remove the cursors.

Both Pre-trigger and Post-trigger data can be examined in their Table.

 Now click on the Module Setup button of the Force Sensor Module box to display the Force 1 – Options tab window shown below.



## 3.3.2 Table of data

 Click in the box alongside **Table** to place a tick in it and a Table of data similar to that displayed below will appear.

🍄 Force	1	
Time [s]	Force 1 [N]	^
-0.05	0.07	
-0.04	0.07	
-0.03	0.03	
-0.02	0.03	
-0.01	0.03	
0	3.43	
0.01	3.43	
0.02	3.43	
0.03	3.43	
0.04	0.10	
0.05	0.10	
0.06	0.10	
0.07	0.10	
0.08	0.07	
0.09	0.07	
0.1	0.07	*

Use of the scroll bar will allow you to examine both Pre-trigger and Post-trigger data.

- Click on the Close buttons in the top right-hand corners of both the Force 1 –Options window and the Force 1 Table window to close them.
- Click on the **Unzoom-all icon** (1) to return the graph to the original view.
- Now click on the Close button in the top right-hand corner of the Graph window, and then on the Close buttons in the top right-hand corners of both the Sub-icon and Main-icon bars to close them.
- Disconnect the Force sensor from the USB Bridge Module , but leave the latter connected to the PC.

# 3.4 Plotting an XY graph

There are a number of experiments in which two sensors are used and their data plotted, not against Time, but with what the two sensors have measured plotted against each other. Examples of this are graphs of (i) the Voltage across a component (e.g. resistor, light emitting diode or bulb) against the current flowing through it and (ii) the Pressure of a fixed mass of gas (at constant volume) against the Temperature of that gas. In the following experiment you will be looking at the first of these two examples as a small filament light bulb is switched on.

You will be connecting up the circuit shown below.



Circuit diagram

This circuit is ideally mounted on a box so that there are sockets available for the sensors to be plugged into. The ideal switch is a microswitch with lever and the bulb an m.e.s. tubular 6.5V 300mA. The battery should be a 6V one and here it has been placed inside the box. A photograph of such a setup is shown below.



- Connect a Voltage sensor <sup>±20V</sup>/<sub>±20V</sub> to one of the sockets of the USB Bridge module using a Edu-Logger®<sup>™</sup> short cable.
- Connect a Current sensor ±2.5A to either the other socket of the USB Bridge module
   If or to that of the Voltage sensor ±20V using a Edu-Logger®<sup>™</sup> short cable.
- Plug the Voltage sensor's red and black 4mm plugs into the sockets across the bulb, red plug to red socket and black plug to black socket.
- Plug the Current sensor's red and black 4mm plugs into the other two sockets, red plug to red socket and black plug to black socket.

Double-click the Edu-Logger® shortcut icon sensor modules display in the Module window as shown below.

Module window	NeuLog <sup>™</sup> Neuron Sensors Network Technology	😥 🛃 🗒 🛠 🍉	
voltage 1 0.00 v	Ver: 2.44		
Current			
0mA			

Click on the **On-line Experiment icon** in the Main-icon bar to display the Graph window below.

Module window	NeuLog	Neuron Sensor	s Network Technolo	•y 🧯		<u> </u>	>	8	
Voltage 1 0.00v	Ver: 2.44	4	se se (	unta ( <mark>1</mark> 744)	Ø 🕑	i 📰 🎸	a 🖓	8	
	NeuLog	<b>A</b> (43) <b>—</b>							
	V         mA           20         2,400           18         2,200           16         2,000           14         1,800           12         1,400           14         1,800           12         1,400           10         1,000           6         800           2         200           98         0           10         200           4         -600           -4         -600           -5         -600           -6         -1,000           -1         -1,000								
	-12 - 1,600 -14 - 1,800 -162,000 -182,200 - 2,000								
	0	0.1	0.2 0.3	0.4	0.5 Time [s]	0.6 0.7	0.8	0.9	1.0

Click on the Experiment setup icon in the Sub-icon bar to display the Experiment setup – Options tab window below.

Experiment setup - Options	×
Options Graph Trigger XY (	Graph
<u>Display</u>	
☐ Table ▼ Graph	
Experiment duration	🖌 10 seconds 📃
Sampling rate	🖌 10 per second 📃

The switching on of the bulb is a very fast event, so click on the down-arrow alongside Experiment duration and select '150 milliseconds'. Now click on the down-arrow alongside Sampling rate and select '1000 per second' if it is not already automatically selected for you. The Experiment setup – Options tab window should now appear as shown below.

Experiment setup - Options	×
Options Graph Trigger XY	Graph
<u>Display</u>	
└ Table ☑ Graph	
Experiment duration	🖌 150 milliseconds 📃
Sampling rate	🖌 1000 per second 📃

Click on the Graph tab to display the Experiment setup – Graph tab window displayed below.

Experiment setup - Graph		×
Options Graph Trigger XYG	raph	
Display 2 axes only	Voltage 1 Current 1	•
Participating sensor(s)	<u>):</u>	

It should appear exactly as shown above with the box alongside **Display 2 axes only** ticked and those below **Participating sensor(s)** also ticked for 'Voltage 1' and 'Current 1'. Reset if necessary.

With such a fast event you will need to use the **Trigger** to begin the data logging.

 Click on the Trigger tab to display the Experiment setup – Trigger tab window and then click in the box alongside Trigger to display the following window.

Experiment setup - Trigger	×
Options Graph Trigger XY Graph	
✓ Trigger	
Sensor Voltage 1 💌	
Level 0 V	
Mode Rise 💌	

To trigger the start of sampling you will make use of the rise in voltage across the bulb as the switch is closed.

• Check first that the sensor indicated is 'Voltage 1' and that the **Mode** is 'Rise'. Now click in the box alongside **Level**, delete '0' and type in '0.1' as shown below.

Options Graph Trigger XY Graph
✓ Trigger
Sensor Voltage 1 💌
Level 0.1
Mode Rise 💌

Alternatively you could have set the Trigger so that it is based on the Current sensor, Current 1, rising to a level of say 5mA. Clicking on the down-arrow in the box alongside **Sensor** would have displayed 'Current 1' and the **Level** would then indicate in units of mA (milliamperes).

At this stage it is useful to first see each of the graphs plotted against Time.

- Click on the Close button in the top right-hand corner of the Experiment setup Trigger tab window.
- Click on the Run experiment icon series in the Sub-icon bar. The message 'Waiting for Trigger ' will be displayed on the screen.

 Now close the switch and hold it closed for a second or so and graphs similar to those below should be plotted.



Click on the Pre-trigger button does not be build by the second the graphs so that they display what was happening to the Current through the bulb, and the Voltage across it, just before and at the moment, the switch was closed.



The interesting section of the above graphs is obviously just before the bulb was switched **ON** to when the bulb was fully lit where the two graph-lines level off. So it is useful to be able to examine this section in more detail and you will see how to do this on the following page.

Click on the Zoom window icon . Locate the mouse cursor at a point above the 'spike' on the graph and to its left, press the mouse left button and keep it pressed, move the mouse cursor to the right and down a little below the graphs to where they have both levelled off. A rectangle should appear as shown below.



 Release the mouse button and the section of graph within the zoomed box will be displayed greatly enlarged as shown below.



Up until this point, running the experiment has resulted in graphs of 'Voltage across the bulb against Time' and 'Current through the bulb against Time' being plotted separately. This done, you can now plot what is known as an **XY graph** with the 'Voltage across the bulb' on the X-axis and the 'Current through the bulb' on the Y-axis.

Click again on the Experiment setup icon in the Sub-icon bar to open the Experiment setup window again. Click its XY Graph tab to display the Experiment setup – XY Graph tab window overleaf.

Experiment setup - XY Graph	
Options Graph Trigger XY Graph	
X-axis	
• Time	
Seconds	
C Sensor	

• Click on the radio button alongside **Sensor** to display the following window.

Experiment setup - XY Graph	×
Options Graph Trigger XY Graph	
X-axis	
C Time	
<ul> <li>Sensor</li> </ul>	
Voltage 1	
1	

- Click on the down-arrow in the box beneath the Sensor radio button and select 'Voltage 1' from the drop-down menu if it is not already selected.
- Now click on the Close button in the top right-hand corner of the Experiment setup XY Graph tab window to close it.

You should now see that you have a graph window displayed with **Voltage** along the X-axis and **Current** along the Y-axis.



Alternatively you could have plotted the graph with Voltage along the Y-axis and Current along the X-axis.

This graph shows us that until 4V the bulb acts as a constant resistance.

The bulb resistance increases and the current goes down when the bulb turns on.

### 3.4.1 Using the Functions – Math tab

With both voltage and current data available you can now obtain a graph of how the resistance of the bulb changed with time by generating a graph of (voltage across the bulb  $\div$  current flowing through the bulb) plotted against Time. The math tab window of the functions window provides a number of mathematical functions with which to generate new data from the existing data. In this case, to calculate the resistance of the bulb.

Click on the Experiment setup icon in the Sub-icon bar to open the Experiment setup – Options tab window again. Click its XY Graph tab to display the Experiment setup – XY Graph tab window below.

	3
Options Graph Trigger XY Graph	
X-axis	I
💭 Time	
<ul> <li>Sensor</li> </ul>	
Voltage 1	

 Click on the radio button alongside **Time** to again display the Experiment setup – XY Graph tab window as shown below.

Experiment setup - XY Graph		×
0	Options Graph Trigger XY Graph	
	X-axis	
	• Time	
	Seconds 🗾	
	🔿 Sensor	
L		

 Click on the Close button in the top right-hand corner of the Experiment setup – XY Graph tab window to close it again. Click on the Show functions icon 3 to display the Functions tab window. Click on its Math tab to display the Math tab window and, using the down-arrows and scroll bar as necessary, select as shown on the next page.

Math	×
Functions Statistics	Math
А	
Voltage 1	-
В	
Current 1	-
А/В	<b></b>
K	
0	

'K' would provide the value of a constant, negative or positive, if one had been required (it isn't).

Now click the Calculate math icon is to display a graph of how the Resistance of the bulb changed with time and the slightly changed Functions – Math tab window shown below.

Math	X
Functions Statistics	Math
А	
Voltage 1	•
В	
Current 1	-
A/B	
К	
0	<b>F</b>

 Click on the Close button in the top right-hand corner of the Functions – Math tab window to close it and you should now see displayed the original graphs of Voltage against Time and Current against Time, and an additional graph (in blue) of Voltage/Current against Time as shown below.



This last graph is of course the resistance of the bulb plotted against Time. It has slightly odd units of [V/mA] or 'Volts  $\div$  milliamperes' which equates to 'kilohms'.

We can see how the resistance of the bulb goes up when the bulb turns on.

Click again on the Show functions icon 3 to display the Functions Math tab window.
 Click on the Erase icon 3 to erase the newly calculated graph.

**NOTE:** Be careful to not confuse it with the near identical **Clear experiment results icon** in the Sub-icon bar.

- Now click on the Close buttons in the top right-hand corners of the Functions Math tab and Graph windows to close them.
- Click on the Close buttons in the top right-hand corner of the Sub-icon and Mainicon bars to close them.
- Unplug the Voltage ±20V and Current ±2.5A sensors from both the Light bulb circuit and from the USB Bridge module 1, but leave the USB Bridge module connected to the PC.

# 3.5 Single step mode

An experiment can be run taking and recording measurements only when you wish to do so. This is known as **Single Step mode** and is accessed from the **Single step icon**  $\boxed{free}$  on the Sub-icon bar. In this mode data will only be collected from the sensor(s) on each click of the **Single step icon**  $\boxed{free}$ .

# 3.5.1 Introduction

This mode will be used when what is, or are being measured, are discontinuous or do not change as a function of Time. It may be that you wish to record the temperatures of a range of different soil or sand samples that have been exposed to the Sun for an hour. You would

simply insert a Temperature sensor  $\blacksquare$  into each sample in turn and record their temperatures into a table by repeatedly clicking on the **Single step icon**  $\blacksquare$ . A graph of the data would probably not be useful, though if one was plotted its X-axis would display as a Counter (1, 2, 3, 4 etc.). You can try this by following the instructions below.

- Fill three 250ml beakers with different soils or sands of differing colours and leave them in the Sun for about an hour. If it is not a sunny day, pour warm water at three different temperatures into the beakers instead.
- Connect the Temperature sensor to the USB Bridge module using a Edu-Logger®<sup>™</sup> short cable.
- Double-click the Edu-Logger® shortcut icon <sup>22</sup> to display the window shown below.

If the temperature is **not** displayed in °C, click on its **Module setup button** and reset to the Celsius radio button.



Click on the On-line experiment icon in the Main-icon bar to display the Sub-icon bar and Graph window shown below.

Module window	NeuLog <sup>11</sup> Neuron Sensors Network Technology
Temperature 1 18.1 °C	Ver 2.44 🛃 🚛 🐜 🐜 🎮 🏧 🖓 💷
	110
	80 70
	60
	28
	10- 0- 
	-20 0 1 2 3 4 5 6 7 8 9 10 Time

Now click on the Experiment setup icon is to display the Experiment setup – Options tab window shown below.

Experiment setup - Options		X
Options Graph Trigger XY (	Graph	
<u>Display</u>		
☐ Table ☑ Graph		
Experiment duration	🖌 10 seconds	•
Sampling rate	🖌 10 per second	-

- Click the cursor in the box alongside **Graph** to untick it.
- Click in the box alongside **Table** to place a tick in it.
- Now click the Close button in the top right-hand corner of the Experiment setup -Options tab window to close it.

You will now see that the Graph window has been removed but a Table window has appeared as shown below.

- Place the Temperature sensor's rod into the first beaker, wait for a few seconds and then click on the Single step icon in the Sub-icon bar.
- Now move the Temperature sensor's rod into the second beaker, again wait for a few seconds and then click on the Single step icon in the Sub-icon bar.
- Now move the Temperature sensor's rod into the third beaker, again wait for a few seconds and then click on the Single step icon in the Sub-icon bar.

You should now see a table with the three temperatures displayed as below.

🦈 Table		
Samples	Temperature 1 [°C]	Manual values
1	22.8	
2	29.7	
3	34.3	
	1	

- Click the **Close button** in the top right-hand corner of the Table window to close it.
- Click on the Close buttons in the top right-hand corners of both the Sub-icon and Main-icon bars to close them.
- Disconnect the Temperature sensor from the USB Bridge module but leave the USB Bridge module connected to the PC.

### 3.5.2 Single step with a linear-fit line

**Single step mode** can be used when you need to collect sensor readings at specific values. You could investigate how the current through a  $100\Omega$  resistor, a diode, a 6V m.e.s. filament bulb and a light emitting diode (led), vary with applied voltages of approximately 0V,  $\pm 1.5V$ ,  $\pm 3V$ ,  $\pm 4.5V$  and  $\pm 6V$ , provided simply by two sets of cell holders. Such a setup does not require a continuously variable voltage supply and so is easier to provide multiple sets for class use. By finally setting up **XY graphs** with the data and overlaying these with a (best-fit) **Linear-fit line**, it can be seen which of the components obeyed **Ohm's Law**. The following instructions are for investigating a  $100\Omega$  resistor.

- Connect a Voltage sensor ±200 into one of the sockets of the USB Bridge module using a Edu-Logger®<sup>™</sup> short cable.
- Connect a Current sensor ±2.5A into either the other socket of the USB Bridge module or into the unused socket of the Voltage sensor ±200 using a Edu-Logger®™ short cable.

 Construct the circuit shown in the diagram and photograph below, but leaving the red lead (dotted line) from the Current sensor to the battery disconnected.



Double-click the Edu-Logger® shortcut icon below.



Click on the Voltage sensor's Module Setup button to display the Voltage 1 – Options tab window shown below.



• Click the cursor in the box alongside **Grid** so that both **Graph** and **Grid** are ticked.

 Now click on the Graph tab so that the Voltage 1 – Graph tab window is shown as below.



- Use the down-arrow alongside Y max and reset its value to '10'.
- Use the up-arrow all alongside Y min and reset its value to '-10'.
- Click the Close button in the top right-hand corner of the Voltage 1 Graph tab window to close it.
- Click on the Current sensor's Module Setup button to display the Current 1 Options tab window below.

Current 1	
Options Graph	
Display	
☐ Table ☐ Digital ☑ Graph ☐ Grid	

- Click the cursor in the box alongside Grid so that both Graph and Grid are ticked.
- Now click on the Graph tab so that the Current 1 Graph tab window is shown as below.



- Use the **down-arrow** I alongside **Y** max and reset its value to '100'.
- Use the **up-arrow a**longside **Y** min and reset its value to '-100'.
- Click the Close button in the top right-hand corner of the Current 1 Graph tab window to close it.
- Click on the On-line Experiment icon in the Main-icon bar to display the Graph window below.



 Connect the Voltage sensor's red plug into the 0V (really -6V) socket of the lower cell holder.

The Voltage Sensor's Module box should now display a value around '-6V' and the Current Sensor's Module box a value around '-50mA'. If the Current Sensor's Module box is displaying a value of around '+50mA' instead, swop the Current sensor's connections over.

 Click on the Single step icon in the Sub-icon bar and the first plot will displayed as shown in the Graph window below.



- Now move the Voltage sensor's red plug into the 1.5V (really -4.5V) socket of the lower cell holder.
- Click on the Single step icon <sup>main</sup> in the Sub-icon bar.
- Repeat the process of moving the Voltage sensor's red plug into the -3V, -1.5V, 0V (join of the two cell holders), +1.5V, +3V, +4.5V and +6V sockets of the cell holders, clicking on the Single step icon in the Sub-icon bar after each move.
- You should have a graph similar to that shown below.



- Remove the Voltage sensor's red plug from the 6V socket of the upper cell holder.
- Now click on the Experiment setup icon icon icon icon bar to display the Experiment setup – Options tab window below.



Click on the XY Graph tab to display the Experiment setup – XY Graph tab window below.



 Click on the radio button alongside Sensor to display the following Experiment setup – XY Graph tab window shown below.

Experiment setup - XY Graph	×
Options Graph Trigger XY Graph	
X-axis	
🔘 Time	
• Sensor	
Voltage 1	

- Check that the box below the Sensor radio button displays 'Voltage 1' as this is, as the Independent variable, what is to be plotted along the X-axis. If 'Current 1' had been the Independent variable and so should be plotted along the X-axis, then you would click on the down-arrow and select 'Current 1'. Do not do it.
- Click the Close button in the top right-hand corner of the Experiment setup XY Graph tab window to close it. This should leave on the screen a graph of 'Current through the resistor' plotted against the 'Voltage across the resistor', similar to that shown below.



 Now click on the Show functions icon 4 to display the Functions tab window with 'Linear-fit' as shown below.

Functions	×
Functions Statistics Math	
Voltage 1	-
Linear-fit	-
12	

- Click on the down-arrow alongside 'Voltage 1' and select 'Current 1' instead.
- Now click the Calculate function icon 
   Icon to display the Functions tab window with the equation of the Linear-fit line similar to that below.

Functions		X
Functions	Statistics Math	
Current	1 🔹	
Li	near-fit 📃 💌	
Linear-fit:	: Y=6.69X -3.39	
	🛃 📵 😰	

 Click the Close button in the top right-hand corner of the Functions tab window to close it. This should leave on the screen a graph of 'Current through the resistor' plotted against the 'Voltage across the resistor', together with the overlaid Linear-fit line, as shown below.



If you wished to remove the Linear-fit line then click the **Show functions icon**  $\textcircled{3}{3}$  again and then click the **Erase icon** in the Functions tab window. Whilst the **Erase icon** looks much like the **Clear experiment results icon** in the Sub-icon bar, clicking on it **only** removes the **Linear-fit line**.

Additionally you should have noticed the **Extrapolation to zero icon** is which, if you had a graph line that needed to be extended back to zero, would be one to click on.

You may also have noticed that the Functions tab window also has a 'Linear-fit [between cursors]' facility as shown below. This allows you to select just a section of a graph **between cursors** in which to calculate and display a Linear-fit line.

Math	×
Functions Statistics Math	
Current 1	
Linear-fit	-
Linear-fit	
Area	
Gradient (cursor 1)	
Polynomial (between cursors)	

- Now click on the **Close button** in the top right-hand corner of the Graph window.
- Click on the Close buttons in the top right-hand corners of both the Sub-icon and Main-icon bars to close them.
- Unplug the Voltage ±20V and Current ±2.5A sensors from both the circuit and the USB
   Bridge module sourced to the PC.

# 3.5.3 Extrapolation to zero

This same **Single step mode** can be used in an investigation to predict the value of Absolute Zero by recording the variation of Pressure of a fixed mass of gas/air (at constant volume) with Temperature. The data collected can then be put into an XY graph with the Temperature plotted along the X-axis from +100 °C to -400 °C and the Pressure along the Y-axis. With a best-fit (linear-fit) line applied to the plots and extrapolated back to where the Pressure is zero, the predicted Absolute Zero of temperature can be determined. The value obtained in this simplified experiment rarely matches that in the textbooks as the tubing connected to the Pressure sensor will contain gas/air at a different (lower) temperature that that recorded by the Temperature sensor. You can conduct the experiment by following the instructions below.

#### NOTE: Take care. In this experiment you will be using very hot water.

- Connect a Pressure sensor into one of the sockets of the USB Bridge module using a Edu-Logger®<sup>™</sup> short cable.
- Connect a Temperature sensor into either the other socket of the USB Bridge module into either the other socket of the Bridge into either the other socket of the Bridge module into either the other sensor into a Edu-Logger®™ short cable.
- Fit the Pressure sensor's adapter/disconnector to that on the end of the tube connected to the boiling tube.
- Attach the Temperature sensor sensor to the side of the boiling tube with a couple of rubber bands.

The photograph below shows the setup so far.



- Fill an electric jug kettle with cold water and ice.
- Place the boiling tube with Temperature sensor attached into the jug of cold water so that the water and ice just cover it. Anchor it in place using a retort stand, bosshead and clamp.
- Double-click the Edu-Logger® shortcut icon below.



- Leave it to cool down to near 0 °C keeping the water well stirred.
- Click the On-line experiment icon in the Main-menu bar to display the Graph window below.



 Click on the Single step icon in the Sub-icon bar and you should then see the first point plots displayed as shown on the following page.

Veul.	.00	6											8
~			-	A+3 /	2								
1	ĥ	-	kPa 700		122.00								
1	00		658										
	90		600										
	80	-	550										
	70	111	500										
	60		450										
	58		400										
	40	ssure	350										
	30	Pre	300										
	20	1.11	250										
	10		200										
			150	•									
	10	1	100	•									
	20		50										
		<u>80</u>	•	1	2	3	4	5 Samples	6	7	8	9	10

- Switch the kettle ON and allow it to heat up for a minute or so to get to a temperature of around 15°C. Switch OFF and stir the water until the temperature and pressure readings have stabilised.
- Click on the Single step icon in the Sub-icon bar again to obtain the next plots on the graph.
- Repeat the switching ON and OFF of the kettle to raise the temperature by around 20 °C each time until the water is near or at its boiling point, stirring until the readings have stabilised, then clicking each time on the Single step icon in the Sub-icon bar to record and display the plots.

You should now see a set of plots displayed on the graph window similar to those shown below.

Neu	Log									×
E		0	-	🖉 🥵 🔨						
	<u>kPa</u> 700		110 1	Seconderio de company						
-1	658		100							
1	600		90							
1	550	1	80							
-2	500		70							
2	450	100	60							
9	400	ture	50							
ressu	358	npera	40							
٩.	258	Ter	30							
	200		20							
	150	1.1.1	10							
- 5	100									1
	50	1.1	-10							
		÷	-20			,	,			
			0	1	2	3 Sem	4 ples	5	6	7

Click on the Experiment setup icon in the Sub-icon bar to display the Experiment setup – Options tab window shown below.

Experiment setup - Options	xperiment setup - Options 🛛 🔀							
Options Graph Trigger XY	Graph							
<u>Display</u>								
└ Table ✔ Graph								
Experiment duration	🖌 10 seconds 📃							
Sampling rate	🖌 10 per second 📃							

Click on the XY Graph tab to display the Experiment setup – XY Graph tab window displayed below.

xperiment setup - XY Graph	X
Options Graph Trigger XY Graph	
X-axis	
<ul> <li>Time</li> </ul>	
Seconds 📃	
Sensor	

 Click on the radio button alongside Sensor to display the following Experiment setup – XY Graph tab window shown below.

Experiment setup - XY Graph 🛛 💈	3
Options Graph Trigger XY Graph	
X-axis	
💭 Time	
Sensor	
Temperature 1	
1	

G 💽 🐜 🔿 🎮 ns Graph Trigger XY Graph perature 1 -20 -10 ó 40 50 Temperature 1 [\*C] 60 70 90 100 10 30 80

You should now see the points plotted on an XY graph similar to that displayed below.

- Click on the Close button in the top right-hand corner of the Experiment setup XY Graph tab window to close it.
- Click on the Show function icon (4) to display the Functions tab window shown below.

Functions Functions Statistics Math	×
Temperature 1	-
Linear-fit	•

■ Click the **down-arrow** I alongside 'Temperature 1' and select 'Pressure 1' instead.

Click on the Calculate function icon (Interpretent to display the graph below with its best-fit (linear-fit) line and note that there is quite a good linear-fit.

Neu		0.		R+3	1											8
Pressure	kPa kPa 650 600 550 500 450 450 350 350 350 250 200 150			<u> </u>			Fun	ctions Inctions S Pressure 1 Pressure 1 Linear-fit Y	italistics 1 ar-fit =0.28x +3	Math 3.49						
	100 50 0		-20	-10	Ó	10	20	30	40	50	60	70	80	90	100	110

This graph can now be extrapolated back to where the Pressure would be zero.

Click on the Extrapolation to zero icon to display a graph similar to that below.



You can now determine what the Temperature would be if the Pressure had been reduced to zero. You can either see this direct from the graph itself or from the information given on the Functions tab window 'X(Y=0) = -334.63' indicating a temperature of -334.63°C. Not quite the -273.15°C for Absolute Zero that the textbooks all quote, but not far off.

- Now click on the Close buttons in the top right-hand corners of both the Functions window and the Graph window to close them.
- Then click on the Close buttons in the top right-hand corners of both the Sub-icon and Main-icon bars to close them.

- Carefully (it will be hot) unclamp and then remove the boiling tube and Temperature sensor from the hot water.
- Disconnect the Pressure sensor from the boiling tube using the adapter/disconnector.
- Remove the Temperature sensor strapping it to the boiling tube.
- Pour away the hot water.
- Disconnect both the Pressure sensor and Temperature sensor from the USB
   Bridge module but leave the USB Bridge module connected to the PC.

### 3.5.4 Inputting data into the manual values column

**Single step mode** can also be used when one of the quantities you wish to measure and record has no sensor available at the moment. An investigation of how the 'Pressure of a fixed mass of gas/air (at constant temperature)' varies with 'Volume occupied by the gas/air' can be dealt with in this way to investigate Boyle's Law. The Pressure would be recorded when the Volume of the air was near 60ml, 55ml, 50ml, 45ml, 40ml, 35ml and 30ml, the Volumes being arranged via the graduations on a gas syringe connected to the Pressure sensor but adding an extra few millilitres to allow for the volume of air in the connecting tube. Whilst the Volumes would show initially as just 1, 2, 3, 4, 5, 6 and 7, you can then enter their actual values into the 'Manual values' column that appears in tables as an extra. Then, with the aid of the Functions – Math tab, values of '1/Volume occupied by the gas/air' can be calculated and used in plotting a graph of 'Pressure of a fixed mass of gas/air at constant temperature' against 1/Volume occupied by the gas/air' and seeing if Boyle's Law is obeyed. The following instructions will show how to use the **Single step mode** in this way.

- Plug the Pressure sensor into either of the two sockets on the USB Bridge Module
   using a Edu-Logger®<sup>TM</sup> short cable.
- Take a 60ml syringe and remove its plunger. Put a little silicon grease around the plunger's seal and then insert it back into the syringe to the 60ml graduation.
- Connect the syringe to the Pressure sensor via the adapter/disconnector.

The basic setup is as in the photograph below.



 Double-click the Edu-Logger® shortcut icon so that the Pressure Sensor's Module box displays in the Module window as shown overleaf.

Module window	NeuLog <sup>™</sup> Neuron Sensors Network Technology	😥 🛃 🗒 🛠 🍉	
Pressure 1 98.0kPa			
<b>6</b>			

Click the On-line experiment icon in the Main-menu bar to display the Graph window below.

	Module window Pressure 1 98.0kPa	NeuLog <sup>™</sup> Neuron Sensors Network Technology Q Q Q Q Q Q Q Q Q Q Q Q Q Q Q Q Q Q Q	
Tine		NeuLog File 50 50 50 50 50 50 50 50 50 50	

Click on the Single step icon in the Sub-icon bar and you should then see the first point plot displayed as shown below.

Ne	uLog												×
	0	0		A-3	7								
ľ	kPa 700												_
	650	2											
	000	2											
	550	-											
1	500	-											
1	450	1											
	400	-											
sure	350	2											
Pres	300	-											
	250	-											
	200	2											
L	150	1											
	100												
	100												
	50												
	0	-	0	i	2	3	4	5 Samples	6	7	8	9	10

- Move the plunger to the 55ml graduation and click the Single step icon in the Subicon bar.
- Move the plunger to the 50ml graduation and click the Single step icon K
- Move the plunger to the 45ml graduation and click the Single step icon K
- Move the plunger to the 40ml graduation and click the Single step icon Tell again.
- Move the plunger to the 35ml graduation and click the Single step icon the again.

You should now see that you have a graph similar to that shown below.

Ne	uLog							E
	0	A (+3)						
	<u>kPa</u> 700 -	[						
L	650							
L	600 -							
L	550							
	500							
	450 -							
æ	400							
nssa,	350 -							
ď.	300 -							
	250							
	200							
	100				200			
	50 -		1.1710					
	0 -							
		0	1	2	3 Samples	4	5	6

Click on the Experiment setup icon in the Sub-icon bar to display the Experiment setup – Options tab window shown below.

Experiment setup - Options	2
Options Graph Trigger XY (	Graph
<u>Display</u>	
☐ Table ✔ Graph	
Experiment duration	🖌 10 seconds 📃 💌
Sampling rate	🖌 10 per second 📃

• Click the cursor in the box alongside **Table** to reveal the Table window with an empty 'Manual values' column as shown below.

Samples	Pressure 1 [kPa]	Manual values
1	97.9	
2	106.2	
3	116.6	
4	128.6	
5	143.8	
6	160.7	

• First click the cursor in the header marked 'Manual values' and type in 'Volume [ml]'.

**NOTE:** It is essential that this column header is dealt with at this stage and that the units are within square brackets [].
Now enter the volumes of air in the syringe that corresponded to the recorded pressures. Add on the extra volume of air that is in the connecting tubing, probably about 3ml. A Table window similar to that shown below should now appear.

Samples	Pressure 1 [kPa]	Volume [ml]	
1	97.9	63	
2	106.2	58	
3	116.6	53	
4	128.6	48	
5	143.8	43	
6	160.7	38	

You can move (click on and drag) the various windows about if sections are being obscured.

Click on the Show function icon icon to display the Functions tab window shown below.

Functions	×
Functions Statistics Math	
Pressure 1	•
_⊿ <sup></sup> Linear-fit	-

At this stage you will need to ensure that the calculations that are about to be made are dealt with to a sensible precision. The measurements of Volume were to two significant figures, so those of any calculations using them should also be to two significant figures.

 Click on the **Tools-icon** in the Main-icon bar to display the Tools-icon bar shown below. Ignore the other boxes and icons here. They will be dealt with in Chapter 8.



If necessary, click on the up-arrow or down-arrow to set the significant figures to
 '2' in the Significant figures field indicated above.

### 3.5.5 Using the functions math tab to operate on data

 Now click the Math tab of the Functions window to display the Functions – Math tab window shown below.

Math				×
Functions	Statistics	Math	]	
А				
Pressure	1		•	
В				
Pressure	1		•	
log(A)			•	
K				
0			P	

- In box A click on the down-arrow alongside ' Pressure 1' and select 'Volume'.
- In the box alongside 'log(A) click on the **down-arrow** and select '1/A'.

The Functions – Math tab window should now be displayed as shown below.

Math				3
Functions	Statistics	Math		
A				
Volume			-	
В				
Pressure	1		•	
1/A			•	
К				
0		(		

Now click on the Calculate math icon to display the combination of Experiment setup – Options, Table, and Graph and Functions – Math tab, windows shown below. Click and drag the windows as required.

Experiment setup - Options	X	to a Table	1			X
Options Graph Trigger XY Graph	Bors network rectinitionary	Samples	Pressure 1 [kPa]	Volume [ml]	1/(Volume ) [1/ml]	
Display	the second se	1 64 1	97.9	63	0.016	
✓ Table		2 VII 2	106.2	58	0.017	
<ul> <li>Grahu</li> </ul>		3	116.6	53	0.019	
Experiment duration 10 seconds		4	128.6	48	0.021	
		5	143.8	43	0.023	
Sampling rate	<u> </u>	6	160.7	38	0.026	
			L.			×
Math	A+3 🔼					
Functions Statistics Math					[1/m]]	
A						
v diane						
В						
B Pressure 1						8
B Pressure 1					0.03 -	: :
B Pressure 1					0.03 -	
B Pressure 1					0.03 -	
B Pressure 1 1/A K 0					0.03 -	
B Pressure 1 17/4 K 0					0.03 -	
B Pressue 1 T/A K D					0.03 -	lume )
B Pressue 1 7/A K 0 0 8 8 9 8 9 8 9 8 9 8 9 8 9 8 9 8 9 8					0.03 -	(Volume )
B Pressure 1 ▼ 17/A K 0 0 0 0 0 0 0 0 0 0 0 0 0					0.03 -	1/(Volume )
B Pressure 1 ▼ 17/A K 0 200 200					0.03 -	1/(Volume )
B Pressure 1 17/A ▼ K 0 250 260 260 260					0.03 -	1/(Volume )
B Pressure 1 7/A K 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0					0.03	1/(Volume)
B Pressure 1 17/A K 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0					0.03 -	1/(Volume )
B Pressure 1 17/A K 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0					0.02	1/(Volume )
B Pressure 1 17/A K 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	•	· · ·			0.03 - 0.02 - 0.01 -	o 1/(Volume)

 Click on the XY Graph tab of the Experiment setup – Options window to display its XY Graph tab window shown below.

Experiment setup - XY Graph	×
Options Graph Trigger XY Graph	
X-axis	
<ul> <li>Time</li> </ul>	
Seconds 🗨	
C Sensor	

• Now click on the **radio button** alongside 'Sensor' 'to display the following window.

 Click the down-arrow alongside 'Pressure 1' and select '1/Volume [1/ml]' to reveal the following Graph window.

NeuLog						8
- Q - Q	A+3					
kPa						
700 -						
650 -						
600 -						
550 -						
500 -						
450 -						
400						
2						
360 - S						
2 300 -						
250 -						
200 -						
150 -			100			
100 -			1.0			
50						
	0.01	0.015	0.02 1/(Volume ) [1/ml]	0.025	0.03	0.035

 Now click back on the Functions tab of the Functions – Math tab window to display the Functions – Functions window shown below.

Functions	X
Functions Statistics Math	
Pressure 1	-
_⊿ <sup></sup> Linear-fit	•

Click on the Calculate function icon and the best-fit (linear-fit) graph will be revealed as shown below.

NeuLog			×
Fo	0		
kP 700			
650		Functions Calification Math	
600 550		Pressure 1	
500			
450		Linear-fit: Y=6254.206X -1.536	
350 300			
250			
200			
100			
60			
		0.01 0.015 0.02 0.025 0.03 1/(Volume) [1/ml]	0.035

The final task is to now extrapolate that graph back to 0,0 to see if it passes through it, or nearly so.

Click on the Extrapolation to zero icon to display a graph similar to that below.



With a straight-line graph passing through 0,0 it shows that the 'Pressure of a fixed mass of air (at constant temperature)' is directly proportional to the 'Volume occupied by that air' – an alternative way of expressing Boyle's Law :

# Pressure of a fixed mass of gas volume occupied by that gas = a constant value (at constant temperature)

You could have seen whether or not the law was verified without getting the above graph plotted. Instead of calculating '1/Volume occupied by the air', you could have simply used the Functions – Math tab facility to calculate 'Pr essure 1. Volume' ( $A \cdot B$ ) and noted in the Table that the values of this product were nearly the same each time.

- Now click on the Close buttons in the top right-hand corners of Experiment setup XY graph tab, Functions – Functions tab, Table and Graph windows to close them.
- Then click on the Close buttons in the top right-hand corners of both the Tools-icon and Main-icon bars to close them.
- Disconnect the **Pressure sensor** from the syringe using the adapter/disconnector.
- Disconnect the Pressure sensor if from the USB Bridge module is but leave the USB Bridge module connected to the PC.

# 3.6 Using photogates

Photogates are used to measure the Time(s) for which their infrared beams are interrupted. By inputting the length(s) of Timing-cards passing through photogates, both Velocity and Acceleration can be calculated too.

Cut out, as accurately as possible, two Timing-cards from black plastic card, one with a single 50mm 'flag' and another with two 50mm 'flags' separated by a 40mm gap as shown below, and another with at least three 'flags' of the same or different size and spacings.



Single 'flag' timing-card



Double 'flag' timing-card

### 3.6.1 Time and velocity mode



- Attach a Single 'flag' timing-card to a model car.
- Place the model car on an approximately 20° sloping runway and temporarily secure it in position.
- Fix a photogate sensor on a retort stand using a bosshead and position it half-way down the runway and such that the timing card will interrupt its beam.
- Connect the Photogate sensor Iso one of the sockets of the USB Bridge Module using a Edu-Logger®<sup>™</sup> short cable.

The arrangement of the equipment should be much like that in the photograph above.

Double-click the Edu-Logger® shortcut icon so that the Photogate sensor module displays in the Module window as shown below.

Module window	NeuLog <sup>111</sup> Neuran Sensors Network Technology	😥 🛃 🗒 🛠 🍉	5 2 X
Photogate 2 1	Ver: 2.40		
8 <b>–</b>			

Click on the On-line Experiment icon in the Main-icon bar to display the Table shown below.

🦈 Photoga	ite		
Samples	Time [s]	Velocity [m/s]	
			1

 Click on the Experiment setup icon in the Sub-icon bar to display additionally the Photogate – Time and Velocity mode window shown on the following page. Note the purpose of each section.



Time and Velocity mode:Setup with a single 'flag' timing-card and one photogate to<br/>measure either Time or VelocityAcceleration mode 1:Setup with a double 'flag' timing card and one photogate to<br/>measure Acceleration.Acceleration mode 2:Setup with a single 'flag' timing-card and two photogates to<br/>measure Acceleration.Status graph mode:Setup with single, double, triple 'flags' or more on timing-<br/>card and one photogate to display graphically the digital<br/>status, 1 or 0, of the photogate's output as the 'flags'<br/>passed through it.

- Click in the box alongside **Digital** to place a tick in it to open the Digital meter window.
- Click on, and drag, the Table, Digital meter, and the Photogate Time and Velocity mode window, to arrange as shown below. It may be helpful to resize some of these windows too.



78

- Use the **down-arrow** i just above the image of the Timing-card to change the value to '50' as 50mm (5cm) is the width of its 'flag'. Leave the other setting on **Time**.
- Click on the Run experiment icon state in the Sub-icon bar. Mark the release position, let the model car run through the photogate to display the first measurement in the Table and a Digital meter display of how long the Timing-card took to pass through it. Resize the Digital meter window as required to accommodate the value and its unit.
- Repeat twice more, releasing the model car from the same position, to obtain a set of measurements similar to those shown below. Note that an average value is provided too.

Module window NeuLog <sup>™</sup> Neuron Sensors Network Technology	Q 🛃 📅 🛠 🤛	
Photogate	🖗 🕶 🖙 🌄 🕷	2 2
	Photogate	
	Samples Time [s] Velocity [m/s]	
	1 0.1281 0.78	
	2 0.1280 0.78	
	3 0.1229 0.81	
Tme S		
Photogate		
<b>0.1229</b>	s]	

Note that the Table also has velocities recorded as well as times. If you had clicked on the **down-arrow** alongside **Time** and selected **Velocity**, the Table window would have appeared the same but the Digital meter would have shown velocities instead.

- Now click on the **Stop experiment icon** in the Sub-icon bar.
- Click on the Clear experiment results icon sin the Sub-icon bar.

### 3.6.2 Acceleration mode 1

The apparatus used here is almost the same as in Time and Velocity mode, with the exception of the use of a Double 'flag' timing-card instead of a Single 'flag' version. A photograph of the setup is shown below.



- Attach a Double 'flag' timing-card to the model car.
- Now click on the **Acceleration mode 1** image to display the image displayed below together with changes to the Table and Digital meter windows.



 Use the up-arrows i just above the image of the Double timing-card to change the values to '50' as 50mm is the width of both of the 'flags'.

- Click on the Run experiment icon in the Sub-icon bar. Let the model car run through the photogate to display the first measurement in the Table and Digital meter display of how great its acceleration was. Widen the Digital meter window as required to accommodate the value.
- Repeat twice more, releasing the model car from the same position, to obtain a set of measurements similar to those shown below.



- Now click on the Stop experiment icon in the Sub-icon bar.
- Click on the Clear experiment results icon 4 in the Sub-icon bar.

### 3.6.3 Acceleration mode 2



- Replace the Double 'flag' timing-card by the Single timing-card again.
- Connect a further Photogate to the USB Bridge module using an additional Edu-Logger®<sup>™</sup> short cable and mount it slightly lower down the runway. Again check that its height is correct for the Timing-card being able to break its beam.

The apparatus should appear much as in the photograph in the previous page.

- Click on the **Search for sensors icon** <sup>1</sup> in the Main-icon bar.
- Click on the On-line Experiment icon in the Main-icon bar to display the Table shown below.

amples	Time [s]	Valacity Im (a)	
		velocity [m/s]	

 Click on the Experiment setup icon in the Sub-icon bar to display additionally the Photogate – Time and Velocity mode window below.



 Click on the Acceleration mode 2 image of a Single-timing card with two photogates to display the following.



**NOTE:** When using two identical sensors, such as here with two photogates, it is essential that they have differing Sensor IDs (above they are '1' and '2'. If you find that they have the same ID Number you will need to change one of them as shown in Chapter 8, Section 8.1. Note also that you do need to identify and select in 'Choose photogates' sequence' which photogate is passed through first and which second. This is most easily done by noting the result of an obvious acceleration which should have a positive value – swop in Sensor IDs if it comes up negative.

- Click in the box alongside **Digital** to place a tick in it to open the Digital meter window.
- Click on, and drag, the Table, Digital meter, and the Photogate Acceleration mode 2 window so that all are viewable.
- If necessary, use the **up-arrow a** just above the image of the Single timing-card to change the value to '50' as 50mm (5cm) is the width of the 'flag'.
- Click on the Run experiment icon in the Sub-icon bar. Let the model car run through the photogates to display the first measurement in the Table and on the Digital meter of how great its acceleration was. Widen the Digital meter window as required to accommodate the value.

 Repeat twice more, releasing the model car from the same position, to obtain a set of measurements similar to those shown below.



- Now click on the **Stop experiment icon** in the Sub-icon bar.
- Click on the Clear experiment results icon in the Sub-icon bar.
- Disconnect only the lower photogate <a>[</a>
   If the lower photogate

This mode can be used to measure velocities and momentum in collision.

 Use the Scroll down arrow and change 'acceleration' to 'velocity'. The table will show four columns:



Fields for the carts' mass are opened too.

- Write the carts' mass in the fields.
- Push two carts against each other so they will collide after they pass through the gates and then return again through the gates.



The table will show the velocity and the momentum of each cart before and after the collision.

### 3.6.4 Status graph mode

This mode produces a graph showing how the digital status (0 or 1) of the photogate changes with time as a Timing card passes through.



The apparatus setup required now is much like that used for the Time and Velocity mode and is shown in the photograph in the previous page.

- Attach a Triple 'flag' timing-card to the model car.
- Click on the Search for sensors icon <sup>1</sup>/<sub>2</sub> in the Main-icon bar to display the following.



Click on the On-line Experiment icon in the Main-icon bar to display the Table window shown below.

ne		
Time [s]	Status	
	Time [s]	Time [s] Status

 Click on the Experiment setup icon in the Sub-icon bar to display additionally the Photogate – Time and Velocity mode window displayed below.



• Click on the **Status graph mode** image of the single photogate with the Triple-timing card (far right) to reveal the Status graph mode window below.

Photogate		×
Display	Experiment duration	
⊽ Table □ Graph	25 milliseconds	

- Click on the down-arrow beneath Experiment duration and select '1 second'.
- Click in the box alongside **Graph** to put a tick in it.

• Size the various windows and then click on, and drag, the Table, Graph and the Photogate – Status graph mode windows, so that all are viewable as shown below.

Madule window		ý 🖉 🖳 🛠 🍉	
Photogate 1 0	- 🥹 / 🚁 🐅 🔛	🎫 🚧 🚰 🛋 😽 🌆	St 🖉 🔹
Photogate		Photogate	
		Samples Time (s) Status	
<u>Display</u> ♥ Table ♥ Graph	Experiment duration 1 second		
	NeuLog		
	g		
	•		
	E	\$	

 Click on the Run experiment icon in the Sub-icon bar and pass run the model car with the multi-flagged timing-card through the photogate to obtain results similar to those shown below.

Module window		nology 🧔 🦉	2%	< 🔶	3	_	
Photogate 1 0		i 🔤 🔤 🍻	<b>P</b>	<b>~</b>		5	
			🦈 Photogate				
Photogate			Samples	Time [s]	Status		^
			28	0.0360	1		
			29	0.0380	1		_
			30	0.0400	1		
			31	0.0420	1		
Dicolay			32	0.0440	0		
Table	Experiment duration		33	0.0460	0		
Graph	1 second 📃		34	0.0480	0		
			35	0.0500	0		
			36	0.0520	0		
			37	0.0540	0		
			38	0.0560	0		~
		47 0.56 0.65 s	0.74	0.83	0.92	2 1	

Click on the Stop experiment icon in the Sub-icon bar.

You may have noticed that the Graph window has a **Pre-trigger button**  $\leq$  in its lower left-hand corner.

Click on this Pre-trigger button does not below, as below, that a change of Status from Digital 0 to Digital 1 started the data logging.



Without the timing-card interrupting the photogate beam the Status is 0 and when it interrupts the beam it is 1. The 0 is indicated by a voltage output from the photogate a little above 0V and the 1 by a voltage output of near 5V.

- Now click on the Close buttons in the top right-hand corners of the Status graph, Table and Photogate mode windows to close them.
- Then click on the Close buttons in the top right-hand corners of both the Tools-icon and Main-icon bars to close them.
- Disconnect the photogate from the USB Bridge module but leave the USB
   Bridge module connected to the PC.

## 3.7 Video – record, play back and save with data

The Logger Sensors software allows you to record a video of the experiment using a video camera or webcam connected to the PC. The video can then be played back to review the experiment.

Before proceeding with the video configuration, you must make sure that the camera is connected and installed on your PC. It may well need a Driver file to be installed which would have been provided with the camera. Updates of such Driver files are usually obtainable via the Internet.

In this example you will just use a single Temperature sensor, but it could just as well be a chain of sensors.

### 3.7.1 Configuring the video camera or webcam

NOTE: You may well find that your video camera or webcam is already configured for use with Edu-Logger®<sup>™</sup> and that no adjustments need to be made. However, if that is not the case, the details below show how to access the setup for configuring (setting up the various properties e.g. brightness, contrast etc.) of a video camera or webcam.

- Connect a compatible webcam or video camera.
- Check that the USB Bridge module signal plugged to a USB port on the PC.
- Connect the Temperature sensor is to the USB Bridge module is using a Edu-Logger® short cable.
- Double-click the Edu-Logger® shortcut icon <sup>22</sup> to display the following.



Click on the On-line experiment icon in the Main-icon bar so that the Graph window is displayed as shown below.

Module window	NeuLog <sup>11</sup> <sub>Neuron Sensors Network Technology</sub> 😥 🚑 🚉 🛠 🍉 🗉 🗉
1 1 18.3 °C	Ver 2.41
<b>0</b>	
	90
	70
	e
	90 T
	10 10
	-10
	0 1 2 3 4 5 6 7 8 9 10 Time

Click on the **Display Video icon** on the On-Line Experiment Sub-icon bar to display additionally the Video Camera Frame in the closed state shown below.



Note that the functions of each icon are indicated above but can also be seen by moving the cursor over each on-screen icon.

Click on the Open camera icon and adjust the camera's focusing device to give a sharp image within the Video Camera Frame of you holding the Temperature sensor.

**NOTE:** If you need to reconfigure the camera or webcam, click on the **Set Camera options** icon icon A window will open similar to that shown below which is from an Viewflex or Extravalue 1.3 Megapixel Snake Webcam – USB – No driver install required, Plug and Play that works fine on both Windows<sup>®</sup> XP and Vista. This will be specific to your video camera or webcam, so you would need to refer to its reference manual.

Video Source	? 🗙
Capture Source Device Settings Camera Controls	
Select a Video Device:	
USB Video Device	•
OK Cancel 4	
	-PPIY

 Configure if required and then click on its Apply and/or OK button(s) to close the window.

# 3.7.2 Recording, playing back, and saving experiments with video

The next stage would be to make all the key settings required for the experiment e.g. Experiment duration, Sampling rate, Trigger setting, whether you want a Table and/or Graph etc. Here you will use the simplest default settings with the Temperature sensor.

 Click on the Experiment setup icon in the Sub-icon bar to open the Experiment setup – Options tab window within the Graph window as shown below.

Experiment setup - Options		×
Options Graph Trigger XY	Graph	
<u>Display</u>		
∟ Table ✓ Graph		
Experiment duration	🖌 10 seconds 📃	
Sampling rate	🖌 10 per second 📃	

 This default setting is fine so click on the Experiment setup – Options tab window's Close button in its top right-hand corner to close it.

**NOTE:** To see both the graph being plotted and the video, you will need to reduce the size of the Video frame by clicking the cursor on one of its **corners** and dragging towards the **centre** of the frame. You will also need to drag the reduced sized frame to a convenient position on the PC screen so that the graph is not totally obscured.

- Click on the Run experiment icon in the Sub-icon bar. The video recording will start as the first item of data is recorded and will stop automatically after 10 seconds, or sooner if you click on the Stop experiment icon in the Sub-icon bar.
- View the recorded experiment (video and graph) by clicking on the Play video icon .
   Note that this icon now incorporates a Pause button . You use this button if you wished to pause the video and graph plotting at some point. Clicking again on the Play video icon .
   video icon .

If you wished to end the video and graph plotting before it had finished, you would click on the **Stop video icon** 

Save the video and the associated experiment's results by clicking on the Save Experiment icon . Initially the Edu-Logger® Docs folder will open. Open the Experiments folder within it (create this folder if it does not already exist). Type a suitable title (e.g. Temperature) into the box alongside File name, check that the 'Save as type' is Experiment File, and then click the OK button.

- Now erase both the experiment's results and video (but not the saved file) by clicking on the Clear experiment results icon in the Sub-icon bar.
- Click on the Close button in the top right-hand corner of the Video Camera Frame to close it.
- Now rerun the saved video and experiment by clicking on the Open experiment icon
   , navigate to the Experiments folder, select the 'Temperature' file (click on it) and then click the Open button.

The experiment data and video recording will load and the camera frame will open automatically.

 Click the Play video icon vito to watch the graph of the experiment being plotted in synchronization with the video of it.

Clicking on the **Single Frame icon** puts a copy of the image in the video frame into the **Clipboard**. From there it can be pasted temporarily into an open word processor or similar file for incorporation into a report at a later time.

- Now click on the Close buttons in the top right-hand corners of the Video Camera Frame and the Graph window to close them.
- Then click on the **Close buttons** in the top right-hand corners of both the On-line Experiment icon-bar and Main-icon bar to close them.
- Disconnect the Temperature sensor sensor from the USB Bridge module and also the USB Bridge module from the PC.
- Unplug the webcam from the PC.

# **Chapter 4 – Off-line Experiment Mode**

In this mode the experimental results are not displayed in real-time. The sensors are preprogrammed to perform measurements. The experiment setup, together with the last data set, is saved in each sensor's internal non-volatile memory.

**Off-line experiment mode** is where the sensors have been programmed by connection to a PC or the MDU with different Sampling rates and Experiment durations. The sensors are then **disconnected** from the PC or MDU and put into battery-powered chains (or singly) to collect their data on the pressing of the sensors' **Start/stop** buttons. However, in this mode, each sensor can be programmed independently to be triggered to record data as well as having their **Start/stop** buttons pressed at different times.

Reconnection of the sensors to a PC or the MDU is required to upload the data and analyse it. However, all graphs would be overlaid on the longest Time axis with t=0s being where each sensor appears to have been triggered. Hence it does not take into account the time differences between the pressing of any sensors' **Start/stop** buttons and their independent triggering.

At **Off-line experiment mode** the sensors could **remain** connected to the PC. As before, the sensors can be programmed with different Sampling rates and Experiment durations. Again, each sensor can be programmed independently to be triggered to record data.

Starting is enabled by clicking the **Run sensors icon** described by pressing the **Start/Stop** buttons on the sensors, so all the sensors start at the same time. Stop the

experiment by clicking on the **Stop measurement icon I** on the PC, by pressing the **Start/Stop** buttons on the sensors, or wait for the experiment duration to be completed, so the experiment ends automatically.

Data is then uploaded in order to display and analyse it. However, all graphs would be overlaid on the longest Time axis with t = 0s being where each sensor appears to have been triggered.

The activities in this chapter have been designed merely to show how to operate in this mode and its special features.

## 4.1 The off-line experiment mode menu

Double-click the Edu-Logger® shortcut icon is to display the following.



EDU-LAB LTD

Click on the Off-line experiment icon on the Main-icon bar. A new Sub-icon bar will then be displayed as shown below.



As with the On-line Experiment mode, note that the functions of each icon are indicated as above but can also be seen by moving the cursor over each on-screen icon.

- Click on the Close button in the top right-hand corner of the Graph window to close it.
- Also click on the Close button in the top right-hand corner of the Sub-icon bar to close it.

## 4.2 Off-line experiment with battery module

In this Off-line experiment mode you will be using sensor(s) whilst it/they are connected and

powered by a Battery module . Each sensor can store up to 5 different experiments data. Before connection, check that this module's batteries are OK by pressing its button and seeing that the red LED alongside it lights. If the LED does not light then you will need to replace the batteries.

### 4.2.1 Setting up a single sensor without a trigger

In this mode each sensor retains in its memory the settings with which it was previously programmed so, when reprogramming, it is best to first reset sensors to their default settings unless you know that the settings are those you know you wish to use. This is done using the

**Restore sensor factory defaults icon** whe Tools-icon bar shown below.



Connect the Temperature sensor 3 to the USB Bridge module 3 using a short cable.

Each module has two identical cable sockets. As before, you can use either. The sockets enable you to connect the sensor modules in a chain.

Click on the Search for sensors icon in the Main-icon bar. The program will scan and display the connected sensor's Module box, in this case that of the Temperature sensor, automatically in the Module window on the left-hand side of the screen as shown below.



The **Temperature Sensor Module** box labelled below again appears in the Module window on the far left of the screen.



 Click on the Tools icon in the Main-icon bar to display the Tools-icon bar shown below.



- Now click on the Restore sensor factory defaults icon to reset the Temperature sensor to its default settings.
- Now click the Off-line experiment icon window shown below.



 Click the Module Setup button of the Temperature Sensor Module box to display the Temperature 1 – Options tab window shown below.

Temperature 1 - Options	×
Options Graph Trigger	
<u>Display</u>	
☐ Table ☐ Digital ✔ Graph ☐ Grid	C Fahrenheit
Experiment duration	🖌 10 seconds 📃
Sampling rate	🖌 10 per second 📃

Note that this window is different from the one displayed in Section 2.6 for the On-line Experiment mode. The **Experiment duration** and **Sampling rate** have to be set up here, as does the use of a **Trigge**r, although the Trigger is not being dealt with until Section 4.2.5.

NOTE: It is essential that at least one change is made to the settings as only then are they downloaded into the Temperature sensor module's memory. If you only want the default settings you would need to change one of them and then go back to it. As you make changes an OK message will be displayed indicating that the change has been downloaded into the sensor's module memory.

- Check that the radio button is set to 'Celsius'. Reset if necessary.
- Click on the down-arrow alongside Experiment duration and set it to '30 seconds'.
- Leave the **Sampling rate** set to '10 per second'.
- Ignore the Graph tab as its settings are not stored in the sensor's memory.
- Click on the **Trigger** tab to display the Temperature 1 Trigger tab window below.

Temperature 1 - Tr	igger 🛛 🔀
Options Graph T	rigger
Trigger	
Level	0 • <b>c</b>
Mode	Rise

 In this first example the Trigger is **not** being used so check that the box alongside **Trigger** is **not** ticked. If necessary, untick it by clicking in the box.

- Close the Temperature 1 Trigger window by clicking on the Close button in its top right-hand corner.
- Now disconnect the short cable, with the Temperature sensor still attached, from the USB Bridge module and plug it into the Battery module
- If you have a digital display module (VIEW 100), you can connect it to the chain.

### 4.2.2 Running a single sensor without a trigger

#### NOTE: Take care. In this experiment you will be using hot water.

- Pour about 200ml of hot water at about 60 °C into a 250ml beaker.
- Pour about 200ml of cold water at a temperature below that of the laboratory/room into another 250ml beaker.
- Press the Start/Stop button on the Temperature sensor. Its red light emitting diode (LED) will turn ON. The sensor will have started sampling the temperature 10 times per second for 30 seconds and is recording the data in its internal memory.
- After about 5 seconds insert the Temperature sensor rod into hot water.
- After about a further 10 seconds take the Temperature sensor rod out of the hot water and insert it into the cold water.

After 30 seconds, the red LED will turn **OFF**. The data logging is complete. Pressing the **Start/Stop** button sooner would have also stopped the data logging.

- Remove the Temperature sensor from the cold water.
- Unplug the Temperature sensor switch its Edu-Logger® short cable from the Battery
   Unit and reconnect it to the USB Bridge module
- Check that the PC is still in Off-line Experiment mode and click on the Load data from sensors icon in the Sub-icon bar and note the following message:

#### Temperature 1 Experiments

Clicking on Experiments will show the list of stored experiments in the sensor's memory (up

to 5). With the **Load data icon** select **Experiments** to upload the required experiment data.

• Click on **Temperature 1**, on the message, to upload the last stored experiment data.

 The sampled data will be uploaded to the PC and a graph similar to that below will be displayed.



As in On-line mode, all the various Zoom, cursors, functions, point/line graph, clear experiment results, freeze graph and export to spreadsheet, facilities, are available via their icons.

You can upload an experiment, freeze it and upload another one.

- Click the Clear experiment results icon in the Sub-icon bar.
- Leave the Temperature sensor and USB Bridge module connected to the PC.

### 4.2.3 Setting up two sensors without triggers

Additionally plug the Light sensor into either the socket on the Temperature sensor

or the unused one on the USB Bridge module 🌋 using a further short cable.

- Click on the Search for sensors icon in the Main-icon bar so that both sensors' Module boxes are displayed in the Module window.
- Click on the Tools icon in the Main-icon bar to display the Tools-icon bar shown below.



 Click on the Restore sensor factory defaults icon in the Tool-icon bar to reset both sensors to their default settings. Click the Off-line experiment icon in the Main-icon bar to display the Graph window shown below.

Module window	NeuLog <sup>11</sup> Neuron Sensors Network Technology	
Temperature 1 15.5 °C	Ver 2.44 🔄 🖗 🎄 🌲 🥮 🏁 🖓 🚱 💷 🖓 🔐 🎘	
Light 1		
	950 900 90 850 80 750	
	70 - 700 60 - 650 600 말 50 - 550	
	500 450 30 400 20 350	
	- 300 10 - 250 0 - 200 - 155 - 10 - 100	
	20 50 0 1 2 3 4 5 6 7 8 9 Time [s]	10

 Click on the Temperature sensor's Module setup button to display the Temperature 1 – Options window below.

Temperature 1 - Options	×
Options Graph Trigger	
<u>Display</u>	<ul> <li>Celsius</li> </ul>
☐ Table ☐ Digital ☑ Graph ☐ Grid	C Fahrenheit
Experiment duration	🖌 10 seconds 📃
Sampling rate	🖌 10 per second 📃

- Check that the radio button alongside 'Celsius' is checked.
- Use the down-arrow alongside Experiment duration and set this at '2 minutes'.
- Use the down-arrow alongside sampling rate and set this at '5 per second'.
- Again ignore the Graph tab as in this mode the sensor does not store such a setting in its memory.
- No new setting is needed for the Trigger tab, so close the Temperature 1 Options tab window by clicking on the **Close button** ⊠ in its top right-hand corner.

Click on the Light sensor's Module Setup button to display the Light 1 – Options tab window shown on the following page.

ight 1 - Options	X
Options Graph Trigger	
<u>Display</u> ☐ Table ☐ Digital ☑ Graph ☐ Grid	<ul> <li>● 1,000 lx</li> <li>● 6,000 lx</li> <li>● 150,000 lx</li> </ul>
Experiment duration	🖌 10 seconds 📃
Sampling rate	🖌 10 per second 📃

- Check that the radio button is set to '1000 lx'.
- Use the down-arrow alongside Experiment duration to set it to '1 minute'.
- Use the down-arrow alongside Sampling rate to set it to '50 per second'.
- Ignore the **Graph** tab as the sensor does not store this setting in its memory.
- Click on the **Trigger** tab to display the Light 1 Trigger tab window displayed below.

Light 1 - Trigger	8
Options Graph	rigger
Trigger	
Level	0 lx
Mode	Rise 🔽

The Trigger is not being used here so the box alongside Trigger should **not** have been ticked. Untick if necessary.

■ Now close the Light 1 – Trigger tab window by clicking on the Close button in its top right-hand corner.

Both the Temperature and Light sensors have now been setup.

Disconnect the short cable with the two sensors attached from the USB Bridge module
 and connect to the Battery module

### 4.2.4 Running two sensors without triggers

#### NOTE: Take care. In this experiment you will be using hot water.

- If not already available, pour about 200ml of hot water at about 60°C into a 250ml beaker.
- Place the Temperature sensor's rod into the hot water and press the Start/stop button on the Temperature sensor.
- After about 1 minute press the Start/stop button on the Light sensor and move it so that different levels of Illumination are sensed.
- When both LEDs have gone out, disconnect the Edu-Logger®<sup>™</sup> short cable with the sensors attached from the Battery module and reconnect to the USB Bridge module
- Check that the PC is still in Off-line Experiment mode, click on the Load data from sensors icon in the Sub-icon bar and note the following message.

#### All (last experiments) Temperature 1 Light 1 Experiments

When **Experiments** is selected, the list of each sensor appears so one experiment of each sensor can be uploaded.

• Click on **All** and you should have displayed graphs similar to those shown below.



Note that **both graphs** begin at **Time=0.00 minutes** even though the Light sensor did not have its **Start/stop** button pressed until 1 minute after the Temperature sensor had its **Start/stop** button pressed. In this mode, the graphs do not take into account such time differences.

If you had clicked on **Temperature 1** just the Temperature graph would have been displayed. Similarly if you had clicked on **Light 1** then just the Illuminance graph would have been displayed. Although not a sensible one to plot here, you could also have plotted an XY

graph of Temperature against Illumination by clicking on the **Experiment setup icon** in the Sub-icon bar and using one of its Graph tabs to arrange this.

- Click on the Clear experiment results icon keiling in the Sub-icon bar.
- Disconnect the Temperature sensor sensor from the chain but leave the Light sensor and USB Bridge module still connected to the PC.

**NOTE:** The two sensors could also have been programmed and operated separately from each other. Uploading their data subsequently could then have been done (i) separately, in which case only that sensor's data would be displayed on a graph or (ii) both together, in which case both sets of data could be overlaid on the same graph if wished (or indeed uploaded separately).

### 4.2.5 Setting up a single sensor with a trigger

- Click on the Search for sensors icon <sup>1</sup>/<sub>2</sub> in the Main-icon bar.
- Click on the Tools icon in the Main-icon bar to display the Tools-icon bar shown below.



- Click on the Restore sensor factory defaults icon in the Tool-icon bar to reset the sensor to its default settings.
- Click on the Off-line experiment icon in the Main-icon bar to display the window shown below.


Click on the Light sensor's Module Setup button to display the Light 1 – Options tab window shown below.

Light 1 - Options	×
Options Graph Trigger	
<u>Display</u> ☐ Table ☐ Digital ☑ Graph ☐ Grid	<ul> <li>● 1,000 lx</li> <li>● 6,000 lx</li> <li>● 150,000 lx</li> </ul>
Experiment duration	🖌 10 seconds 📃 💌
Sampling rate	🖌 10 per second 📃

- Check that the radio button is set to '1000 lx'.
- Use the down-arrow alongside Experiment duration to set it to '1 second'.
- Use the down-arrow alongside Sampling rate to set it to '100 per second'.
- Ignore the **Graph** tab as the sensor's memory does not store this setting.
- Click on the **Trigger** tab to display the Light 1 Trigger tab window shown below.

Light 1 - Trigger	
Options Graph T	igger
Trigger	
Level	0 lx
Mode	Rise

- Click in the box alongside **Trigger** to put a tick into it.
- Click in the box alongside Level, delete the '0' and type in '5'.

Check that the box alongside Mode indicates 'Rise'. Use the down-arrow to change this if necessary. The Light 1 –Trigger tab window should now appear as shown below.

Light 1 - Trigger	×
Options Graph Trigger	
<b>✓</b> Trigger	
Level 5 Ix	
Mode Rise 👤	

- Now close the Light 1 Trigger tab window by clicking on the Close button in its top right-hand corner.
- Disconnect the Edu-Logger®<sup>™</sup> short cable with the Light sensor we attached from the USB Bridge module and then reconnect to the Battery module .

## 4.2.6 Running a single sensor with a trigger

- Place a finger over the hole on the side of the Light sensor to stop any light from getting in.
- Push the Start/stop button on the Light sensor. Point the sensor's access hole towards the light and then remove your finger from the hole.
- Disconnect the short cable with the Light sensor if from the Battery module and reconnect it to the USB Bridge module
- Check that the PC is still in Off-line Experiment mode and then click on the Load data from sensors icon in the Sub-icon bar to obtain a graph similar to that shown below.



Click on the **Zoom-fit icon** to expand the graph in the Y-axis direction and obtain a graph like the one below.



- Click on the Pre-trigger button 
   <u>H</u> in the bottom left-hand corner of the graph to see how
   quickly the sensor responded after being triggered. A second click on the Pre-trigger
   button will slightly extend the time.
- Click on the Clear experiment results icon in the Sub-icon bar.
- Leave the Light sensor and USB Bridge module science connected to the PC.

### 4.2.7 Setting up two sensors with triggers

- Connect a Temperature sensor sensor to a socket on either the Light sensor or the USB
   Bridge module vising an additional short cable.
- Click on the Search for sensors icon <sup>1</sup>/<sub>2</sub> in the Main-icon bar.
- Click on the Tools icon in the Main-icon bar to display the Tools-icon bar shown below.



- Click on the Restore sensor factory defaults icon in the Tool-icon bar to reset both sensors to their default settings.
- Click on the Off-line experiment icon in the Main-icon bar to display the window shown below.



 Click the Module Setup button of the Temperature Sensor Module box to display the Temperature 1 – Options tab window shown below.

Temperature 1 - Options	×
Options Graph Trigger	
<u>Display</u> ☐ Table ☐ Digital ☑ Graph ☐ Grid	<ul> <li>Celsius</li> <li>Fahrenheit</li> </ul>
Experiment duration Sampling rate	<ul> <li>✓ 10 seconds</li> <li>✓ 10 per second</li> </ul>

- Check that its radio button is set to 'Celsius'.
- Use the down-arrow alongside Experiment duration and set this at '2 minutes'.
- Use the down-arrow alongside Sampling rate and set this at '5 per second'.
- Ignore the **Graph** tab as its settings are not stored in the sensor's memory.
- Click on the Trigger tab to display the Temperature 1 Trigger tab window shown below.

Temperature 1 - Trigger	×
Options Graph Trigger	
<b>Trigger</b>	
Level 0 •C	
Mode Rise 💌	

- Click in the box alongside **Trigger** to put a tick into it.
- Click in the box alongside Level, delete the '0' and type in '30'

Check that the box alongside Mode indicates 'Rise'. Use the down-arrow it to change this if necessary. The Temperature 1 – Trigger tab window should now appear as shown below.

Temperature 1 - Tri	igger	×
Options Graph Tr	rigger	
<mark>▼</mark> <u>Trigger</u>		
Level	30 <b>  •C</b>	
Mode	Rise 💌	

- Now close the Temperature 1 Trigger tab window by clicking on the Close button in its top right-hand corner.
- Now click on the Light sensor's Module Setup button to display the Light 1 Options tab window shown below.

ight 1 - Options	×
Options Graph Trigger	
<u>Display</u> ☐ Table ☐ Digital ☑ Graph ☐ Grid	<ul> <li>1,000 lx</li> <li>6,000 lx</li> <li>150,000 lx</li> </ul>
Experiment duration Sampling rate	<ul> <li>✓ 10 seconds</li> <li>✓ 10 per second</li> </ul>

- Check that the radio button is set to '1000 lx'.
- Use the down-arrow alongside Experiment duration to set it to '1 second'.
- Use the down-arrow alongside Sampling rate to set it to '100 per second'.
- Ignore the Graph tab as its settings are not stored in the sensor's memory.

• Click on the **Trigger** tab to display the Light 1 – Trigger tab window shown below.



- Click in the box alongside **Trigger** to put a tick into it.
- Click in the box alongside Level, delete the '0' and type in '5'.
- Check that the box alongside Mode indicates 'Rise'. Use the down-arrow it to change this if necessary. The Light 1 –Trigger tab window should now appear as shown below.

Light 1 - Trigger			×
Options Graph T	rigger		
✓ Trigger			
Level	5	lx	
Mode	Rise	-	

- Now close the Light 1 Trigger tab window by clicking on the **Close button** in its top right-hand corner.
- Disconnect the short cable with the Temperature sensor and Light sensor attached from the USB Bridge module and connect it to a Battery module

## 4.2.8 Running two sensors with triggers

#### NOTE: Take care. In this experiment you will be using hot water.

- If not already available, pour about 200ml of hot water at about 60°C into a 250ml beaker.
- Also pour about 200ml of cold water into another 250ml beaker.
- First place the Temperature sensor's rod into the cold water and press **the Start/stop** button on the Temperature sensor.
- After about 1 minute press the Start/stop button on the Light sensor and move it so that different levels of Illumination are sensed.
- When the LED on the Light sensor has gone out move the Temperature sensor's rod into the hot water.
- When both LEDs have gone out, disconnect the Edu-Logger®<sup>™</sup> short cable with the sensors attached from the Battery module and reconnect it to the USB Bridge module
- Check that the PC is still in Off-line Experiment mode and then click on the Load data from sensors icon in the Sub-icon bar. Note the following message.

All (last experiments) Temperature 1 Light 1 Experiments

 Click on All to display graphs similar to those shown below. Clicking instead on Temperature 1 or Light 1 would have resulted in just their graph being displayed.



Click on the **Zoom-fit icon** to expand the graph in the Y-axis direction and obtain a graph like the one below.



 Click on the **Pre-trigger button** in the bottom left-hand corner of the graph above to see how quickly the Light sensor responded after being triggered. Its graph will probably be somewhat similar to that shown below.



- Click on the Clear experiment results icon sin the Sub-icon bar.
- Disconnect the Temperature sensor but leave the Light sensor and USB Bridge module connected to the PC.

## 4.3 Off-line experiment mode with a PC

## 4.3.1 Setting up a single sensor without a trigger

- Click on the Search for sensors icon <sup>1</sup>/<sub>2</sub> in the Main-icon bar so that the Light sensor's Module box is displayed in the Module window.
- Click on the Tools icon in the Main-icon bar to display the Tools-icon bar shown below.



- Click on the Restore sensor factory defaults icon in the Tool-icon bar to reset the Light sensor to its default settings.
- Click the Off-line experiment icon in the Main-icon bar to display the following Graph window.



Click on the Light sensor's Module Setup button to display the Light 1 – Options tab window shown below.

ight 1 - Options	X
Options Graph Trigger	
<u>Display</u>	<ul> <li>1,000 k</li> </ul>
Table	C 6,000 lx C 150,000 lx
Graph Grid	
Experiment duration	🖌 10 seconds 📃
Sampling rate	🖌 10 per second 📃

- Check that the radio button is set to '1000 lx'.
- Use the down-arrow alongside Experiment duration to set it to '30 seconds'.
- Use the down-arrow alongside Sampling rate to set it to '5' per second'.

Ignore the **Graph** tab as its settings are not stored in the sensor's memory.

Also ignore the **Trigger** tab as no Trigger is to be set here.

 Close the Light 1 – Options tab window by clicking on the Close button in its top right-hand corner.

The Light sensor has now been setup.

NOTE: Leave the Light sensor Science of the PC via the USB Bridge 📌 .

#### 4.3.2 Running a single sensor without a trigger

 Click the Run sensors icon in the Sub-icon bar or press the Start/Stop button on the Light sensor

Note that you have **two** ways of starting the data logging in this mode.

Move the Light sensor about so that different levels of Illumination are sensed.

 When the LED has gone out click on the Load data from sensors icon in the Subicon bar to display a Graph window similar to the one below.

1,000 -	
950	
900	
850	
800	
750	
/00	
500 T	
600	
± 500 -	
450	
400	
350	
300 -	
250	
200 -	-
150 AMMAA MA	
	$\sim$
	/
	30
Time [s]	

Click on the Clear experiment results icon in the Sub-icon bar.

NOTE: Leave the Light sensor key connected to the PC via the USB Bridge module 1

#### 4.3.3 Setting up two sensors without triggers

- Plug the Temperature sensor sensor into either the socket on the Light sensor is or the unused one on the USB Bridge module sensor unused an additional Edu-Logger®<sup>™</sup> short cable.
- Click on the Search for sensors icon <sup>1</sup>/<sub>2</sub> in the Main-icon bar so that both sensors' Module boxes are displayed in the Module window.
- Click on the Tools icon in the Main-icon bar to display the Tools-icon bar shown below.



 Click on the Restore sensor factory defaults icon in the Tool-icon bar to reset both the Light sensor and the Temperature sensor to their default settings. Click the Off-line experiment icon in the Main-icon bar to display the Graph window displayed below.

Module window	NeuLo	<b>g</b> <sup>TU</sup> Neuron Sense	ors Network Technolo	wy 🌾	) 🔁 j	<u> </u>			
Temperature 1 15.5 °C	Ver: 2.44		æ 🤹 .	<u>i</u> e	Ø 🗗		2~ <u>8</u> ~	8	
Light	NeuLog								
25k	100 - 950 100 - 900 90 - 850 800								
	80 - 750 70 - 700 60 - 650								
	enteradua 20 - 550								
	20 350 - 300 10 250								
	0 - 200 - 150 - -10 - 100 - -20 - 50 -								
		i i	2 3	4	5 Time (s)	6	7	8	9 10

 Click on the Temperature sensor's Module setup button to display the Temperature 1 – Options window below.

Temperature 1 - Options	X
Options Graph Trigger	
<u>Display</u>	<ul> <li>Celsius</li> </ul>
Table Digital	C Fahrenheit
Graph	
Experiment duration	🖌 10 seconds 📃
Sampling rate	🖌 10 per second 📃

- Check that the radio button is set to 'Celsius".
- Use the down-arrow alongside Experiment duration and set this at '1 minute'.
- Use the down-arrow alongside Sampling rate and set this at '30 per minute'.

Ignore the **Graph** tab as its settings are not stored in the sensor's memory.

Also ignore the **Trigger** tab as no Trigger is to be set here.

Close the Temperature 1 – Options tab window by clicking on the Close button in its top right-hand corner.

Click on the Light sensor's Module Setup button to display the Light 1 – Options tab window shown below.

Light 1 - Options	X
Options Graph Trigger	
<u>Display</u> ☐ Table ☐ Digital ☑ Graph ☐ Grid	<ul> <li>1,000 lx</li> <li>6,000 lx</li> <li>150,000 lx</li> </ul>
Experiment duration Sampling rate	✓ 10 seconds ✓ 10 per second

- Check that the radio button is set to '1000 lx'.
- Use the down-arrow alongside Experiment duration to set it to '5 seconds'.
- Use the down-arrow alongside Sampling rate to set it to '50 per second'.

Again ignore the **Graph** tab as its settings are not stored in the sensor's memory, and also ignore the **Trigger** tab as no Trigger is to be set here either.

 Close the Light 1 – Options tab window by clicking on the Close button in its top right-hand corner.

Both the Temperature and Light sensors have now been setup.

NOTE: Leave both of the sensors connected to the PC via the USB Bridge **1**.

#### 4.3.4 Running two sensors without triggers

#### NOTE: Take care. In this experiment you will be using hot water.

- If not already available. pour about 200ml of hot water at about 60 ℃ into a 250ml beaker.
- Similarly, if not already available, pour about 200ml of cold water at a temperature below that of the laboratory/room into another 250ml beaker.
- Place the Temperature sensor's rod into the hot water.

Click the Run sensors icon 4 in the Sub-icon bar and note the following message.

All (last experiments) Temperature 1 Light 1 Experiments

- Click on Light 1 to start the Light sensor's data logging and move the Light sensor about so that different levels of Illumination are sensed.
- When the Light sensor's LED has gone out, again click the Run sensors icon in the Sub-icon bar and, in response click on Temperature 1 to start its data logging.
- After about 30 seconds move the Temperature sensor's rod into the cold water.
- Now click on the Load data from sensors icon in the Sub-icon bar and note the following message.

All (last experiments) Temperature 1 Light 1 Experiments

Click on All to display graphs similar to those shown below.



Note that **both graphs** begin at **Time = 0s** even though the Temperature sensor did not start its data logging until after the Light sensor had finished its data logging. The graphs do no take into account such time differences in this mode.

If you had clicked on **Temperature 1** just the Temperature graph would have been displayed. Similarly if you had clicked on **Light 1** then just the Illuminance graph would have been displayed. Although not a sensible one to plot here, you could also have plotted an XY graph of Temperature against Illumination as outlined in Section 3.5.

121

As in On-line experiment mode you also have access to Zoom-fit, Functions etc.

- Click on the Clear experiment results icon in the Sub-icon bar.
- Unplug the Light sensor but leave the Temperature sensor connected to the USB
   Bridge module and the latter still plugged into the PC.

### 4.3.5 Setting up a single sensor with trigger

- Click on the Search for sensors icon in the Main-icon bar so that the Temperature sensor's Module box is displayed in the Module window.
- Click on the Tools icon in the Main-icon bar to display the Tools-icon bar shown below.



- Click on the Restore sensor factory defaults icon in the Tool-icon bar to reset the Temperature sensor to its default settings.
- Click the Off-line experiment icon with in the Main-icon bar to display the following Graph window.



 Click on the Temperature sensor's Module Setup button to display the Light 1 – Options tab window shown on the following page.

Cemperature 1 - Options	X
Options Graph Trigger	
<u>Display</u> ☐ Table ☐ Digital ☑ Graph ☐ Grid	<ul> <li>Celsius</li> <li>Fahrenheit</li> </ul>
Experiment duration Sampling rate	<ul> <li>✓ 10 seconds</li> <li>✓ 10 per second</li> </ul>

- Check that the radio button is set to 'Celsius'.
- Use the down-arrow alongside Experiment duration to set it to '1 minute'.
- Use the down-arrow alongside Sampling rate to set it to '60 per minute'

Ignore the **Graph** tab as its settings are not stored in the sensor's memory.

 Now click on the Trigger tab to open the Temperature 1 - Trigger tab window shown below.

Temperature 1 - Trigger	×
Options Graph Trigger	
Level C	
Mode Rise 💌	

- Click in the box alongside **Trigger** to put a tick in it.
- Now click in the box alongside Level, delete '0' and replace it by '30'.
- Leave the Mode as 'Rise'.

This will have set up the system so that the Temperature will start to be logged once it has risen past 30 °C and will be logged every second for one minute.

Close the Temperature 1 - Trigger tab window by clicking on the Close button in its top right-hand corner.

NOTE: Leave the sensor connected to the PC via the USB Bridge **1**.

## 4.3.6 Running a single sensor with a trigger

#### NOTE: Take care. In this experiment you will be using hot water.

- If not already available, pour about 200ml of hot water at about 60°C into a 250ml beaker.
- Similarly, if not already available, pour about 200ml of cold water at a temperature below that of the laboratory/room into another 250ml beaker.
- Place the Temperature sensor's rod into the cold water.
- Click the Run sensors icon in the Sub-icon bar or press the Start/Stop button on the Temperature sensor

Note that you have **two** ways of starting the data logging in this mode.

- After about 10 seconds move the Temperature sensor's rod into the hot water.
- When the LED has stopped flashing on the Temperature sensor click on the Load data from sensors icon in the Sub-icon bar to display a Graph window similar to that shown below.



 Click on the Pre-trigger button ≤ near the bottom left-hand corner of the graph to show what the temperature sensor was recording just before and at the time data logging was triggered to start. The Graph window should look similar to that displayed on the following page.



Click on the Clear experiment results icon ker in the Sub-icon bar.

NOTE: Leave the Temperature sensor connected to the PC via the USB Bridge module.

#### 4.3.7 Setting up two sensors with triggers

- Plug the Light sensor into either the socket on the Temperature sensor sont or the unused one on the USB Bridge module solution using an additional Edu-Logger®<sup>™</sup> short cable.
- Click on the Search for sensors icon in the Main-icon bar so that both sensors' Module boxes are displayed in the Module window.
- Click on the Tools icon in the Main-icon bar to display the Tools-icon bar shown below.



 Click on the Restore sensor factory defaults icon in the Tool-icon bar to reset both the Light sensor and the Temperature sensor to their default settings. Click the Off-line experiment icon in the Main-icon bar to display the Graph window displayed below.

Module window	NeuLo	9 Neuron Sens	iors Network Techno	ology 💃		1.%	٠		3	
Temperature 1 15.5 ℃	Ver: 2.44		æ ź.	1. C	Ø 🖗		<u>8</u>	- 🎸	3	
ight Light 25kk	Secul. ag           C         C         E           110         -1,000         900           90         850         900           90         850         800           90         60         650           90         500         500           90         850         900           90         850         900           90         850         900           90         850         900           90         850         900           90         900         900           90         800         900           90         800         900           90         800         900           90         800         900           90         800         900           90         800         900           90         800         900           90         900         900           90         900         900           90         900         900           90         900         900           90         900         900           900         900									
	20 350 20 300 10 250 0 200 10									
	-10 100 -20 50 0	0 1	2 3	4	5 Time (s)	6	7	8	9	10

Click on the Temperature sensor's Module setup button to display the Temperature 1 – Options window below.

Temperature 1 - Options	
Options Graph Trigger	
<u>Display</u> ☐ Table ☐ Digital ✔ Graph ☐ Grid	<ul> <li>Celsius</li> <li>Fahrenheit</li> </ul>
Experiment duration Sampling rate	<ul> <li>✓ 10 seconds</li> <li>✓ 10 per second</li> </ul>

- Check that the radio button is set to 'Celsius".
- Use the down-arrow alongside Experiment duration and set this at '30 seconds'.
- Use the down-arrow alongside Sampling rate and set this at '5 per second'.

Ignore the **Graph** tab as its settings are not stored in the sensor's memory.

• Now click on the **Trigger** tab to open the Temperature 1 - Trigger tab window shown overleaf.



- Click in the box alongside **Trigger** to put a tick in it.
- Now click in the box alongside Level, delete '0' and replace it by '30'.
- Leave the Mode as 'Rise'.
- Close the Temperature 1 Trigger tab window by clicking on the Close button in its top right-hand corner.
- Click on the Light sensor's Module Setup button to display the Light 1 Options tab window shown overleaf.

Light 1 - Options	×
Options Graph Trigger	
<u>Display</u> ☐ Table ☐ Digital ☑ Graph ☐ Grid	<ul> <li>● 1,000 lx</li> <li>● 6,000 lx</li> <li>● 150,000 lx</li> </ul>
Experiment duration Sampling rate	<ul> <li>✓ 10 seconds</li> <li>✓ 10 per second</li> </ul>

- Check that the radio button is set to '1000 lx'.
- Use the down-arrow alongside Experiment duration to set it to '5 seconds'.
- Use the down-arrow alongside Sampling rate to set it to '20 per second'.

Ignore the **Graph** tab as its settings are not stored in the sensor's memory.

 Click on the Trigger tab to display the Light 1 – Trigger tab window shown on the following page.



- Click in the box alongside **Trigger** to put a tick in it.
- Click in the box alongside Level, delete '0' and replace it by '20'.
- Leave the **Mode** as 'Rise'.
- Close the Light 1 Options tab window by clicking on the Close button X in its top right-hand corner.

This will have set up the system so that the Temperature will start to be logged once it has risen past  $30^{\circ}$ C and will be logged five times each second for thirty seconds. Similarly, Illuminance will be logged as soon as its level rises past 20 lx and will then be logged at a rate of twenty samples per second for five seconds.

#### NOTE: Leave both of the sensors connected to the PC via the USB Bridge

#### 4.3.8 Running two sensors with triggers

#### NOTE: Take care. In this experiment you will be using hot water.

- If not already available, pour about 200ml of hot water at about 60°C into a 250ml beaker.
- Similarly, if not already available, pour about 200ml of cold water at a temperature below that of the laboratory/room into another 250ml beaker.
- Place the Temperature sensor's rod into the cold water.
- Click the Run sensors icon 3/2 in the Sub-icon bar and note the following message.

All (last experiments) Temperature 1 Light 1 Experiments

- Arrange the Light sensor is so that the value displayed in the Temperature Sensor Module box is lower than 20 lx.
- Click on Light 1, wait for a couple of seconds, and then move the Light sensor into much brighter light.
- When the Light sensor's LED has gone out, again click the Run sensors icon in the Sub-icon bar and, in response, click on Temperature 1.
- After about 10 seconds move the Temperature sensor's rod into the hot water to start its data logging.

**NOTE:** The sensors' data logging could also have been started by pressing their **Start/Stop** buttons.

 When the LED on the Temperature sensor has stopped flashing, click on the Load data from sensors icon in the Sub-icon bar and note the following message.

> All (last experiments) Temperature 1 Light 1 Experiments

• Click on **All** to display graphs similar to those shown below.



 Click on the Pre-trigger button ≤ near the bottom left-hand corner of the graph to show what the temperature sensor was recording just before, and at the time, data logging was triggered to start. The Graph window should look similar to that displayed on the following page.



Note that **both graphs** show the sensors' Trigger points ( $30^{\circ}$ C and 20 lx) at **Time = 0s** and their values for a short time before the sensors were triggered.

If you had clicked on **Temperature 1** just the Temperature graph would have been displayed. Similarly if you had clicked on **Light 1** then just the Illuminance graph would have been displayed. Although not a sensible one to plot here, you could also have plotted an XY graph of Temperature against Illumination as outlined in Section 3.5.

As in On-line experiment mode you also have access to Zoom-fit, Functions etc.

- Now click on the **Close button** ⊠ in the top right-hand corner of the Graph window to close it.
- Then click on the **Close buttons** in the top right-hand corners of both the Off-line experiment-icon and Main-icon bars to close them.
- Disconnect both the Light sensor and the Temperature sensor from the USB
   Bridge module and unplug the latter from the PC.

## Chapter 5 – Using the RF Communication Modules

## 5.1 **RF** With a PC



For communication between a single sensor, or chain of sensors, and a PC, simply plug the sensor(s) and RF Communication module into a Battery module via its USB (B) socket. Then, at the PC end, plug another RF Communication module into the USB Bridge module's USB (B) socket with the latter plugged into a PC. Communication is then as if the sensor(s) were connected directly.

The RF module has also an ID number. Two RF modules that communicate with each other should have the same ID number.

How to setup the RF module ID number is explained in section 7.3.

Notes:

#### Only two RF modules with the same ID number should be used in a classroom.

The PC does not require any Wi-Fi<sup>™</sup> or Bluetooth<sup>™</sup> installation.

## 5.2 RF With a PC and Groups of sensors

We can use one PC with up to nine groups of sensors operated as described in section 5.1.

Each group of sensors and battery module have an RF module with a different ID number connected to them.

The PC can address each group at a time by setting the ID number of the PC RF module to the required group RF module ID number.

## 5.3 **RF** with the Monitor Display Unit (MDU)



For communication between a single sensor, or chain of sensors, and the Monitor Display Unit (MDU), simply plug the sensor(s) and RF Communication module into a Battery module via its USB (B) socket. Then, at the MDU end, plug another RF Communication module into its USB (B) socket. Communication is then as if the sensor(s) were connected directly.

Whilst connection can be made in this way over up to 30m, within a building, or around buildings, it is likely to be much less.

## Chapter 6 – Saving, Loading, Printing & Exporting to a Spreadsheet

## 6.1 Saving

If an 'Experiments' folder has not already been created to both Save and Load files from, it is essential that you make one. It would best be located in the 'Edu-Logger®, Docs' folder.

Saving the files of data from experiments in both On-line and Off-line experiment modes

simply involves clicking on the **Save experiment icon** [Image: in their respective Sub-icon bars to open a folder similar to those shown below (Microsoft<sup>®</sup> XP left, Vista<sup>®</sup> right).

Save As			? 🛛	Save As						X
Save in: My Recent Documents Desktop My Documents My Computer	Experiments 100 Ohm resistor Boyles law bub switch on the tand cold Gilyah bub offlineB light and temp pressure temperature		•	Save in:	Beperimer Name Light bu	nts Date taken Ib.exp ture.exp	Tags	Size	€ ㎡ 囲•	
My Network Places	File name:     I       Save as type:     Experiment File	•	Save Cancel		File name: Save as type:	Experiment F	File		•	Save Cancel

Into the box alongside **File name:** type an appropriate name which will enable you to recognise what it contains. The **Save as type:** should automatically be 'Experiment File'.

## 6.2 Loading

Loading data files from experiments in both On-line and Off-line experiment modes involves first clicking on the **Open experiment icon** in their respective Sub-icon bars to open a folder similar to those shown below (Microsoft<sup>®</sup> XP left, Vista<sup>®</sup> right).

Open	🕐 🔀			
Look in: Experiments Constant and cold Documents Desktop My Documents My Computer	E - Look in: Recent Paces Desitop Oris Computer	Experiments     Name     Date taken     Light bulb.exp     Temperature.exp	Tags Size	⊨ to of Er- Rating
My Network Places File name: Files of type: Experiment File	Open Cancel	File name: Files of type: Experiment	File	Open     Cancel

The appropriate file is then clicked on to highlight it and then the **Open** button clicked. The **Files of type:** should be displayed as 'Experiment file'. The highlighted file of data will then be loaded into Edu-Logger®.

The same routine is used for saving and loading video files that are associated with experiments, though information on that is also provided in Section 3.6.2 Recording, playing back and saving experiments with video.

## 6.3 Printing

Obviously a printer must be connected.

To print a graph from an experiment in either On-line or Off-line experiment mode the **Print** icon on the Sub-icon bar needs to be clicked. This will **only** print the **graph**.

If other windows within Edu-Logger®<sup>™</sup> require printing, press the **PrtSc** (Print Screen) key on the PC's keyboard to put an image of the whole of the screen into the Clipboard. Then Paste into a graphics package such as Microsoft Windows Paint<sup>®</sup>, OpenOffice Draw<sup>®</sup>, Adobe Photoshop<sup>®</sup>, Corel Paint Shop Pro<sup>®</sup>, GNU Image Manipulator Program (GIMP)<sup>®</sup> or Serif PhotoPlus<sup>®</sup>, cut out whatever you wish (e.g.Table of data), save and then open in, for example, Microsoft WordPad<sup>®</sup>, Microsoft Word<sup>®</sup>, Lotus WordPro<sup>®</sup>, OpenOffice Writer<sup>®</sup> or Softmaker TextMaker<sup>®</sup>.

If it is a large Table of data that requires printing then it is best to Export (See Section 7.4 below on Exporting to a spreadsheet) the data into a spreadsheet such as Microsoft Excel<sup>®</sup>, Lotus 123<sup>®</sup>, Softmaker PlanMaker<sup>®</sup> or OpenOffice Calc<sup>®</sup> and print it directly from it.

## 6.4 Exporting to a spreadsheet

Edu-Logger's® Functions – Math tab enables the processing of data in many ways (e.g. multiplying Voltage by Current (A \* B) to give Power, dividing Voltage by Current (A/B) to give Resistance, calculating the reciprocal of Volume (1/A) or 1/Distance<sup>2</sup> (1/A<sup>2</sup>). However, there are occasions when other mathematical functions are required that have not been incorporated into Edu-Logger®<sup>TM</sup> (e.g. sin A, cos A, tan A, sin<sup>2</sup> A, cos<sup>2</sup> A). An example of this would be with an investigation of Malus' Law in which the Intensity of polarised light passing through two polarisers is proportional to cos<sup>2</sup>  $\theta$  where  $\theta$  is the Angle between the polarisers. So in this circumstance it is useful to be able to export data to a spreadsheet and do further processing there.

Although the icon in the Sub-icon bar indicates the use of Microsoft Excel<sup>®</sup>, you can export to any spreadsheets that make use of CSV (Comma Separated Values) data format such as Lotus 123<sup>®</sup>, Softmaker Planmaker<sup>®</sup> or OpenOffice<sup>®</sup>.

To do this you would click on the **Export to Excel file icon** in the Sub-icon bar to open folders such as those in the next page (Microsoft<sup>®</sup> XP left, Vista<sup>®</sup> right) in which to save the data file. It may have other files already saved in it.

Save As				? 🔀	Save As					×
Save in	Experiments	<u>-</u>	- 🗈 💣 📰	•	Save in:	Bxperimer	nts	•	+ 🗈 💣	•
My Recent Documents Desktop My Documents My Computer	Boyle's Law				Recert Places Desktop Onts Computer Network	Name Light bu	Date modified Ib resistance	Туре	Size	Tags
My Network Places	File name:			Save		File name:	I		•	Save
	Save as type: Ex	cel File	<b>_</b>	Cancel	1	Save as type:	Excel File		•	Cancel

These folders additionally allow you to save in the CSV file format too as you would see if you had clicked on the **down-arrow** salar alongside 'Save as type' to display the following, again Microsoft<sup>®</sup> XP left and Vista<sup>®</sup> right.

Save As				? 🛛	Save As					×
Save in: 📔	Experiments		· + • •		Save in:	Biperiment	ts	-	+ 🗈 💣	<b></b>
My Recent Documents Deskop My Documents My Computer	Boyle's Law				Recent Places	Name	Date modified b resistance	Туре	Size	Tags
My Network File	ile name:		•	Save		File name:			•	Save
Sa	ave as type:	Excel File	-	Cancel		Save as type:	Excel File		•	Cancel
	Ċ	SV File			C.		CSV File			14

If you have Microsoft Excel<sup>®</sup> on your PC you would then select 'Excel File', but otherwise select 'CSV File. The following screen-dump reflects the use of Microsoft Excel<sup>®</sup> but it is a similar process with other spreadsheets.

You would then type a meaningful name (e.g. in this example it was **Malus Law**) into the Filename box and then click the **Save button**. Then, on opening the spreadsheet, the table of data would be displayed as below and further processing done before plotting a graph of Corrected Light Intensity against cos<sup>2</sup> Angle between polarisers.

	Aicrosoft E	xcel - Malı	us Law												- a 🗙
:0	Ele Edit	View Ins	sert Format <u>T</u> ools	Data	Window Help							Ту	pe a questio	n for help	# ×
ID		3017	1 1 Ca 10 1 -0	- (Y -	· . Σ · 2	1.00	Arial		• 10	- 👚 B	IU		at 😐	10 . 31	· A · [
	035	+	15											a Fille and Fille	_
	A	B	C	D	E	F	G	н	I IS .	J	K	L	M	N	0
1	Samples	Light [lx]	Angle [degrees]												
2	1	26	90		_										
3	2	37	80											_	
4	3	64	70											2	
5	4	100	60		_										
6	5	136	50		-									-	-
1	6	180	40												
8	/	203	30		-						-			-	
9	8	256	20												
11	10	337	10											-	
12	10	341					-								-
13											-				_
14							-							-	
15															
16														_	-
17															
18															
19															
20															
21															
22															
23															_
24					_										_
25															_
20					_									_	
21					-										_
20					-									-	
20	-														_
31														-	-
32														-	
33															
	+ H\St	eet1 / She	et2 / Sheet3 /		-				<						>1
Dea	- 1 A	utoShapes *	1100	40	🔹 🛃 👌 -	2 · A ·									

## Chapter 7 – Edu-Logger®<sup>™</sup> Tools

- Connect the USB Bridge module **\*** to a USB port on the PC.
- Connect the Temperature sensor 📃 to the USB Bridge module **\*** using a short cable.
- Double-click on the **shortcut icon** to display the following.

Module window	Noul og	👘 🖓 🖳 🖉 🛷 👞	
	Neuron Sensors Network Technology	The last was in the	
Temperature			
1 10.000	Ver: 2.44		
10.0 %			
S 📕			

The Logger Sensors software has several tools through which a sensor's ID and software language can be changed.

Another important tool saves an experiment's configuration. For each experiment setup an instruction document can be prepared. When the document is opened (doc, docx, wps, odt, tmd, rtf, pdf etc.) the experiment setup is automatically loaded.

Click on the **Tools icon** to reveal the Tools-icon bar shown below. As with all the icon bars, moving the cursor over the icons reveals their purpose as shown below.



## 7.1 Set Sensor ID Number

Setting a number in the box alongside the ID icon and then clicking on the **ID** icon will change all the connected sensors' IDs to this number. Here you will just deal with one sensor.

This function allows the use of up to 9 sensors of the same type.

- With just the Temperature sensor sensor connected, click on the up-arrow all alongside the ID icon and change the number in the box to '2'.
- Click on the ID icon . The sensor will get the instruction to change its ID to '2'.

The **Search** function will run automatically and you will see that the sensor's Module box in the Module window now shows '2'.

Now return the sensor's ID number back to '1' by clicking on the down-arrow alongside the ID icon , change the ID to '1' and then click on the ID icon .

## 7.2 Set RF Module ID Number

Setting a number in the box alongside the RF icon and then clicking on the **RF** icon will change all the connected RF module IDs to this number. Here you will just deal with one sensor.

This function allows the use of up to 9 RF modules.

- Connect one RF module to the USB module and check that no other RF module is connected to a battery module.
- Click on the up-arrow alongside the RF icon and change the number in the box to '2'.
- Click on the RF icon 
   The RF module will get the instruction to change its ID to '2'.
- Repeat the above steps and change another RF module ID to '2'.

## 7.3 Restoring sensors' factory defaults

When sensors have been in use in **Off-line** mode they retain whatever settings and data they last had. Clicking on the Restore sensors' factory defaults icon resets all sensors connected (by cable or via RF Communication modules) to the same defaults that would automatically be set in On-line mode. Here you will see the effect with just the Temperature sensor.

Click on the Off-line experiment icon and graph window will then be displayed as shown below.

Module window	NeuLog <sup>11</sup> Neuron Sensors Network Technology	
Temperature 1 19.7 ℃		
	Yeal og       Image: Second secon	
	-20 0 1 2 3 4 5 6 7 8 Time[s]	9 10

Click on the Module Setup button (1) to display the Temperature 1 – Options window below. It may not look quite like this as all will depend on its last usage and setup.

Femperature 1 - Options					
Options Graph Trigger					
<u>Display</u> ☐ Table ☐ Digital ☑ Graph ☐ Grid	<ul> <li>Celsius</li> <li>Fahrenheit</li> </ul>				
Experiment duration Sampling rate	<ul> <li>✓ 1 minute</li> <li>✓ 30 per minute</li> </ul>				

- Click on the down-arrow alongside Experiment duration and select a value of '2 minutes'.
- Click on the down-arrow alongside Sampling rate and select a value of '60 per minute'.

 Now click on the Trigger tab to display the Temperature 1 – Trigger tab window shown below. Again it may well not appear quite as shown as all will depend on its last usage and setup.

Temperature 1 - Trigger 🔀				
Options Graph Trigger				
✓ Trigger				
Level	20 <b>*C</b>			
Mode	Rise 💌			

- If the box alongside **Trigger** is not ticked, click the cursor in it to tick it.
- Click the cursor in the box alongside **Level** and enter a value of '40'.
- Click the cursor on the down-arrow alongside Mode and select 'Fall'.
- Close the Temperature 1 Trigger tab window by clicking on the Close button in its top right-hand corner.

All would now have been set up to collect data for two minutes at a rate of 60 per minute once the temperature had fallen through the 40 °C point. Assume that this was how it was set up.

It is useful to always have a default setup -a series of settings that the sensor can always be set to and the next few instructions will show you how to return a sensor (or sensors if more than one had been connected) to the default setting.

- Click on the Tools icon in the Main-icon bar and then on the Restore sensors' factory defaults icon .
- Click again on the Off-line experiment icon <sup>1</sup>/<sub>24</sub> in the Main icon bar.

Click again on the Module Setup button is to display the Temperature 1 – Options window. It should now display its default settings as shown below.

Temperature 1 - Options				
Options Graph Trigger				
<u>Display</u>	<ul> <li>Celsius</li> </ul>			
☐ Table ☐ Digital ☑ Graph ☐ Grid	🔿 Fahrenheit			
Experiment duration	🖌 10 seconds 📃			
Sampling rate	🖌 10 per second 📃			

Again click on the Trigger tab to display the Temperature 1 – Trigger tab window. It should now display its default settings as shown below.

Temperature 1 - Trigger			
Options Graph Trigg	er		
Trigger			
Level	0 <b>*C</b>		
Mode	Rise		

■ Click on the **Close button** log of the Temperature 1 – Trigger tab window to close it.

# 7.4 Producing and saving instructions with and without the software configuration

Whilst very often it is easier to produce instructions on paper-based worksheets, there are occasions when it is helpful to be able to display on a PC in conjunction with Edu-Logger®<sup>™</sup>. This can be done in two ways: (i) brief or detailed instructions but without automatically saving the configuration of the setup, and (ii) to save the current configuration of the software including the experiment setup, the used windows and their location, together with information for the user. The latter is useful when you would want to repeat an experiment and do not want to have to set everything up from the beginning again.

## 7.4.1 Producing and saving with the software configuration

Before using this function a document file (doc, docx, wps, odt, tmd, rtf, pdf etc.) with information for the user must be prepared using a word processor or similar application. An example is provided on the following two pages.

# Investigating the Respiration of Maggots Introduction Aerobic respiration is the process in which glucose and oxygen, through the action of enzymes in mitochondria, produce carbon dioxide, water and lots of energy. It can be summarised by the equation: $C_6H_{12}O_6 + 6O_2 = 6CO_2 + 6H_2O + Energy$ (Glucose + Oxygen = Carbon dioxide + water + energy) In this experiment the use of oxygen in the aerobic respiration of maggots in a sealed container can be observed over time with the Edu-Logger®™ Oxygen sensor **Procedure** USB Bri

Place about thirty maggots into a 100 ml beaker, place a piece of gauze fabric over its open end and anchor in place with a rubber band. Pour water at around room temperature into a trough to a depth of about 3cm. Thread the Oxygen sensor through the hole in the centre of a bee-hive shelf so that the sensor is pointing upwards. Place the beaker of maggots on the bee-hive shelf and anchor the Oxygen sensor to it, again with a rubber band. Now place the beaker with Oxygen sensor attached and the bee-hive shelf into the water in the trough. Position a cut-off fizzy drink bottle over the whole lot and into the water so that the air is sealed inside. The photograph above shows how it might appear.


Plug the USB Bridge module into a USB port on the PC and connect the Oxygen sensor is to one of the USB Bridge module's sockets using a short cable.

Double-click the **Edu-Logger® shortcut icon** to display the window shown below.



Click on the On-line experiment icon in the Main-icon bar to display the following Graph window.



Click on Oxygen sensor's Module Setup button to display the Oxygen 1 – Options window shown below.



- Check first that the radio button is indicating '% in air' as shown above.
- Wait for two minutes or so to allow the Oxygen sensor's readings to stabilise. If the Oxygen Sensor's Module box is not displaying a value near 20.9% click on the Calibration icon which should reset it to that value.
- Check that the box alongside Graph is ticked. The Graph tab can be left on its default setting so close the Oxygen 1 – Options window by clicking on its Close button 3.
- Click on the Experiment setup icon is to display the Experiment setup Options window shown below.



- Use the down-arrow alongside Experiment duration, then the scroll-bar on the drop-down menu, and select '1 day'.
- Use the down-arrow alongside Sampling rate, then the scroll-bar on the drop-down menu, and select '60 per hour'. No changes are needed to the other tabs so close the Experiment setup Options window by clicking on its Close button 2.

Click on the Tools icon in the Main-icon bar to reveal the Tools-icon bar shown below.



- Click on the Save configuration icon in the Tools-icon bar. Navigate to the 'Experiments' folder and open it.
- Into the box alongside File name type 'Respiration of Maggots'. The Save as type should be 'Experiment File'. It should appear similar to that shown below.

Save As					? 🔀	Save As					×
Save in	Experiments		• ÷ €	d* 💷 •		Save in:	Experime	nts	- + 6		•
My Recent Documents Desktop My Documents My Computer	100 0hm resis Boyles law Boulb switch or hot and cold Bight bulb Bight variation Bight variation Bight variation Bight and offineBight a pressure temp	nd temp erature				Recent Places Desktop Orris Computer Network	Name	Date modif Type Ib.exp ature.exp	Size	Tags	
My Network Places	File name: Save as tupe:	Respiration of Maggots		3	Save		File name:	Respiration of Maggots		•	Save

- Now click the **Save** button.
- Close the Graph window by clicking on the Close button K in its top right-hand corner.
- Now click the Close buttons at the right-hand ends of both the Tools-icon bar and the Main-icon bar.
- Finally disconnect the Oxygen sensor from the USB Bridge module and the USB Bridge module from the PC.

This would now be set up ready for the user to load the activity by first clicking on the **Tools** 

icon in the Main-icon bar and then on the Load activity icon in the Tools-icon bar. He/she then locates the 'Respiration of Maggots' file in the Experiments folder, opens it and follows the on-screen instructions.

Whilst this illustrates how to do this for an experiment using an Oxygen sensor, it is applicable to any experiment setup with any sensors.

#### 7.4.2 Producing and saving without the software configuration

As in Section 7.3.1 a document file with information for the user must be prepared using a word processor or similar application. The example shown on the following page is one where the software configuration has not been done in advance and only rudimentary details are provided as to how to conduct the experiments. Obviously these could have been very detailed 'hold by the hand' instructions, the choice will be dictated by circumstances – the ability and knowledge of the users.



Setup of pH sensor experiment

The stomach's mucosa has cells within it which secrete hydrochloric acid of between pH 2 and pH 3. This acid activates pepsinogen, secreted by other cells in the mucosa, converting it into the protein-digesting enzyme pepsin. Additionally this acid kills off many of the bacteria that accompany food. However, this acid can also bring about indigestion.

Indigestion is the pain or discomfort in the chest or abdomen that one might get after a rather heavy meal. Its technical name is dyspepsia. Some people refer to it as heartburn owing to the burning feeling they have in their chest. There are many possible causes and sometimes no reason for it at all. The most common cause is when the stomach acid flows back up into the oesophagus in a process known as acid reflux. Alternatively it may be the result of the protective lining of the stomach allowing the acid to come into contact with the wall of the stomach.

Very occasional indigestion is not likely to be a great problem, but repeated indigestion needs consultation with a doctor in case it is the result of a medical condition or problem. To deal with occasional bouts of indigestion some people take indigestion tablets such as *BisoDol<sup>®</sup>*, *Setlers<sup>®</sup>* and *Rennie<sup>®</sup>*, but how effective are they?

#### Your task

In this experiment you will be adding indigestion tablets to some cola which has a pH around that of stomach acid. Using a pH sensor you will need to monitor, and report on, how the pH changes with time in order to (i) compare the effectiveness of various crushed indigestion tablets and (ii) to see how crushing/not crushing them affects their operation.

If you have not already done so, connect the pH sensor to the USB Bridge module using a short cable. Now plug the USB Bridge module into a USB port on the PC.

Click on previous on the Task bar at the foot of the screen to display Edu-Logger®.

Switch between Dindugestion Tablets -... and Preulog on the Task bar as and when necessary.

Click on the Search for sensors icon 💯 in the Main-icon bar.

Place the pH sensor into a beaker of buffer solution of pH 7.0 and then click on the pH sensor's **Module Setup button** to display the pH 1 – Options tab window shown below.

pH 1	×
PH 1 Options Graph Display □ Table □ Digital ✓ Graph □ Grid	×

Click on the **Calibration icon** to set the pH sensor to pH 7.0. The rest of its pH settings will have been calibrated too. You should see an **OK** message displayed briefly as you click on the icon. Close the pH 1 – Options tab window by clicking on the **Close button** in its top right-hand corner.

Now click on the **On-line experiment icon** in the Main-icon bar and you are **almost** ready to start.

At this stage, if you have not already done so, you need to sort out how you are going to conduct your experiments. Discussion with your colleagues would be useful. You will need to consider (or find out in a pilot experiment) what Experiment duration and Sampling time might be most suitable, and then set up the system with them. To compare the graphs that you obtain you may be best to use the Freeze current graph(s) facility to overlay them all. Alternatively you could save each one separately with a unique name (e.g. that of the tablet, crushed or uncrushed) and then Open their files one by one when you have completed all the experiments.

After use, wash the pH sensor in either distilled or tap water and replace in its 'pot'.

With this example you would then save the document with the file name 'Indigestion Tablets' with whatever file extension (doc, docx, wps, odt, tmd, rtf, pdf etc.) was appropriate in an 'Experiments' folder (as a sub-folder of Edu-Logger®, Docs), creating this if necessary.

All would then be set up ready for the user to load the activity by first clicking on the **Tools** 

icon  $\bigotimes$  in the Main-icon bar and then on the Load activity icon  $\bigotimes$  in the Tools-icon bar. He/she then locates the 'Indigestion Tablets' file in the Experiments folder, opens it, and follows the on-screen instructions.

Again, whilst this illustrates how to do this for an experiment using a pH sensor, it is applicable to any experiment setup with any sensors.

### 7.5 Set significant figures

When a sensor is not available to make a measurement (e.g. Volume of gas in a syringe) you can enter the required data through the 'Manual values' column of a Table window. If this data is to then be used in calculations it is essential that those calculations are made to the correct number of significant figures. For example, if the Volumes of gas were 63, 58, 53, 48, 43 and 38ml, they have been measured and recorded to two significant figures. Therefore any calculations using them cannot be made to more than two significant figures. So, before any Functions – Math tab calculations are made, it is essential that you first select the number of significant figures that the calculations are to be made to.

The number of significant figures is set with the **up-arrow** and **down-arrow** of the **Significant figures field** on the Tools-icon bar.

### 7.6 Restore USB

The Logger Sensors software searches automatically for the first free USB port. The software is limited to a USB port no larger than 8.

If a USB port was used in the past, the PC sees it as still in use, even if it is not used.

The **Restore USB field** enables selection of a USB port that was previously used.

Usually you will not need to use this function.

### 7.7 Language Select

The language used within the Edu-Logger®™ software can be changed by clicking on the

**Help icon** in the Main-icon bar to reveal the Tools-icon bar shown below.



By clicking on the **down-arrow** alongside the national flags of Spain and the United Kingdom in the Language select field Exercise T, English or Spanish can be selected. All the words and phrases are then put into the selected language.

## Chapter 8 – Logger Sensors Modules

This chapter deals with details of the Logger Sensors and associated modules and their special features.

A sensor's Operating Range is denoted by the Y max and Y min values specified.

The **Sampling Rates** of the sensors cover two ranges: (i) from 10,000 samples per second to 1 per hour and (ii) from 100 samples per second to 1 per hour. Not all sensors can be sampled at fast rates as some do not respond quickly. In **On-line Experiment** mode all the modules connected together will automatically run at the same rate, but in **Off-line Experiment** mode they can operate at differing rates.

The **Experiment Duration** should be set by the user; whilst these are most often quite short, they could also be long (many days).

**Sampling Rates** and **Experiment Durations** are interdependent and so, for very fast rates, only short durations are available. The combination of fast rates and long durations is limited by the memory storage capacity of the modules.

### 8.1 USB bridge module 🌋 USB-100

This Edu-Logger®<sup>™</sup> module enables a fast connection of the sensors to the PC. It operates with Windows 98(SE), XP and Vista.

The USB module is the first in a chain of sensors connected to a USB port of the PC. It provides both power from the PC to the sensors and communication between the PC and the sensors.

This module has a USB cable on one side and two Logger Sensor sockets on the other side.

# 8.2 Battery module 🗐 Bat-100

The Battery module supplies power to a sensor, or chain of sensors, operating in the Off-Line mode, and to sensors connected to an RF Communication module.

This module can be checked for the 'goodness' of its internal batteries by pressing its pushbutton. This will turn **ON** its light emitting diode (LED) when the batteries are O.K.

#### **Specifications:**

4 AA batteries Battery test button

### 8.3 RF Communication module RF-100

The RF Communication module allows remote operation of a single sensor or a chain of sensors. The remote single sensor or chain is connected to an RF Communication module and to a Battery Module which powers them all. Another RF Communication module should be connected directly to the Monitor Display Unit, or to the PC via the USB Bridge module. More than two units can be used to cater for more chains of sensors or more independent sensors.

The PC does not need to have either Bluetooth<sup>™</sup> or Wi-fi<sup>™</sup> incorporated. Everything required is built into the RF Communication Module.

#### **Specifications:**

Frequency: 2.4GHz DSSS (Direct-sequence spread spectrum). Bit Rate: 1Mbps. Maximum distance of use in open space: 30m.

### 8.4 Digital display module



The VIEW-100 is a small display module that can be connected to any chain of logger sensors working off-line with a Battery Module.

The VIEW-100 automatically searches for the connected sensors and displays one of them digitally. Scrolling for displaying the reading of another sensor is done by pressing the pushbutton switch on the module.

### 8.5 Graphic display module



The Graphic Display Unit (GDU) is used to run experiments without a PC. The GDU displays the sensor's measurements in digital and graphical forms. It can also be used to program the sensor's experiment setup as well as viewing the input from up to five sensors at a time.

This unit has a user-friendly design with a colour graphic display and touch screen.

The Graphic Display Unit is used when a PC is not available for each group. It can work with up to 5 sensors in parallel.

Some of the unit's features are:

- Automatic recognition of sensors.
- Uses preset experiment parameters for easy initiation.
- Communicates with all the sensors or one at a time.
- Controls each sensor's range and measurement units.
- Internal charging circuit is incorporated.
- Mode to view sensor values in real time up to five at a time.
- Automatic Power-Off for longer battery life.

Sensors are connected to the GDU via its USB (A) socket. Remote connection is also possible by plugging an RF Communication module into it and another into the sensor or sensor chain. This enables both setting up and analysing the collected data.

The VIEW-101 can be connected to a chain of sensors ending with the battery module connected to the last sensor of the chain.

When the GDU VIEW-101 receives power, it starts scanning and identifies the connected sensors. The located sensors are displayed on the left side of the screen.

The viewer uses the same icons as the Edu-Logger <sup>™</sup> software:



- Search for connected sensors.



- Run experiment while displaying the results. The data is also saved in the in the modules' internal memory and can be uploaded to the viewer at any time.



- Stop experiment run.



 Upload experiment data from the sensors. This function acts also as zoom Out.



- Zoom fit.



- Erase screen.



- Experiment setup.
- Tools: Set sensor ID, Set RF ID, Set screen shut down, Turn the viewer off. The last two functions are aimed for saving battery.

### Monitor display unit SX0 8.6



The Monitor Display Unit (MDU) is used to run experiments without a PC. The MDU displays the sensor's measurements in digital form only. It can also be used to program the sensor's experiment setup as well as viewing the input from up to two sensors at a time.

This unit has a user-friendly design with a liquid crystal display (LCD) screen and easy to use three-button keyboard. It does not display graphs.

When used with a Photogate, **only** Time can be measured.

The Monitor Display Unit is used when a PC is not available for each group.

Some of the unit's features are:

- Automatic recognition of sensors.
- Power supply for sensors.
- Uses preset experiment parameters for easy initiation.
- Communicates with all the sensors or one at a time.
- Controls each sensor's range and measurement units.
- Works with internal rechargeable batteries (2 NiMH AA cells 2300mAH). •
- Internal charging circuit is incorporated.
- LED charging indicator.
- Battery status check.
- Mode to view sensor values in real time up to two at a time.
- Can freeze readings to view values at a specific time.
- Automatic Power-Off for longer battery life.

Its batteries are charged by connection to a charging unit.

Sensors are connected to the MDU via its USB (B) socket. Remote connection is also possible by plugging an RF Communication module into it and another into the sensor or sensor chain. This enables both setting up and analysing the collected data.

### 8.7 Voltage logger sensor 201 NUL-101

Voltages can be measured across various resistive, capacitive and inductive components, as well as those of photovoltaic cells, batteries and power supplies. This sensor can also be used to measure electrode potentials in Redox reactions and to investigate the charging and discharging of capacitors.

When used in conjunction with the Current sensor the dependence of the current flowing on the applied voltage can be studied in various electric circuits.

This sensor can be used to measure low voltage AC and DC circuits. With its 4mm plugs it can easily be connected into electric circuits.

It can also measure, **using a step-down transformer**, the AC voltage of the Main supply and check its frequency 50/60 Hz (**the input is limited to 60Hz and ±20V maximum**).



Experiment Duration: 50 milliseconds to 31 days.

Range and operation modes	ADC resolution	Accuracy	Resolution	Max Sample Rate (S/sec)
+20)/	10 hit	10/	0.011/	100 Slow
±20V	TU DIL	170	0.01V	3000 Fast

### 8.8 Current logger sensor **±2.5A** NUL-102

This sensor can be used to measure the current in parallel or series low voltage a.c. and d.c. circuits and also investigate the dependence of the current flow through components on the voltage across them.

With its 4mm plugs it can easily be connected into electric circuits.

	Current 1	C	urrent 1 - Graph	
	Options Graph		Options Graph	
Current 1 40mA	<u>Display</u> □ Table □ Digital ☑ Graph □ Grid		Y max. Y min. Y-axis position Zoom fit Color 2500 € -2500 €	

Experiment Duration: 50 milliseconds to 31 days.

#### **Specifications:**

Range and operation modes	ADC resolution	Accuracy	Resolution	Max Sample Rate (S/sec)
+0500mA	10 hit	10/	10~	100 Slow
±200011A	TU DIL	170	TUTTA	3000 Fast

# 8.9 Temperature logger sensor 📓 NUL-103

This is one of the most versatile sensors. It can be used in Biology to monitor ecological systems, to study photosynthesis or to study the effect of temperature on enzymes; in Chemistry, to study exothermic or endothermic reactions, and in Physics to study heat/energy transfer.



#### Experiment Duration: 1 second to 31 days

The sensitive element is within a 180 mm long, 3.2 mm diameter stainless steel tube. This sensor can be used for temperature measurements in solids, liquids or gases.

Range and operation modes	ADC resolution	Accuracy	Resolution	Max Sample Rate (S/sec)
–35°C to 150°C	10 hit	±1°C	0.1°C	100
–31°F to 23 <sup>8</sup> °F		±2°F	0.1°F	100

## 8.10 Light logger sensor **S** NUL-104

This sensor is very versatile with applications in many areas of the natural sciences. In Biology it can be used to study photosynthesis. In Chemistry to study light-emitting chemical reactions. In Physics it can be used to study the effect of changing voltage on a light-bulb's output.

With three ranges, it can be used in low light environments such as a classroom, or high light environments as in daylight outdoors. It measures Illumination.

With both fast and slow modes, it can be used to measure fast light changes such as those produced by light bulbs connected to an AC supply, as well as the near steady levels outside on a sunny day.



Experiment Duration: 50 milliseconds to 31 days.

The light sensor in located in a plastic box just behind an access hole.

Range and	ADC	Resolution	Max Sample
operation modes	resolution		Rate (S/sec)
Illumination: 0 to 1,000 lx	12 bit	1 lx	
Illumination: 0 to 6,000 lx		6 lx	100
Illumination: 0 to 150,000 lx		150 lx	
Signal: 0 to 1,000 lx		1 lx	
Signal: 0 to 6,000 lx		6 lx	3000
Signal: 0 to 150,000 lx		150 lx	

### 8.11 Oxygen logger sensor 🎴 NUL-105

This sensor can be used to make measurements of the level of free-oxygen in air or dissolved oxygen in water.

The Free-oxygen in air mode is used to measure changes in oxygen levels during combustion or in reactions that produce oxygen (hydrogen peroxide decomposition). The Dissolved-oxygen mode is useful in the study of photosynthesis.

To change mode, click the Module Setup icon <sup>1</sup> on the Sensor Module box in Edu-Logger®<sup>™</sup> software or use the change range option on the Monitor Display Unit.



#### Experiment Duration: 1 second to 31 days.

The oxygen sensor is designed for use both in the school laboratory and in the field. It employs easy-to-use polarographic (Clark) technology and replaceable membranes are available for it. The electrode itself is constructed of Delrin<sup>®</sup> for durability.

With its integral thermistor, it provides dependable temperature-compensated measurements. The thermistor is housed in stainless steel and sealed on the electrode's outer wall providing fast, accurate readings.

The installation and replacement of the membrane is quick and easy. Simply fill the membrane cap assembly with DO electrolyte and screw it into place. Two membrane cap assemblies are included with each sensor. Store in de-ionised water between measurements and overnight. Long-term, dissemble, rinse in de-ionised water and store dry.

#### Sensor Calibration:

Calibration of the probe is simply achieved in open air, taking this as a standard level of 20.9%. First connect the sensor to a voltage source (the USB Bridge plugged into a PC, Edu-Logger's® Monitor Display Unit or Battery Unit) and wait for stabilization of the readings (about 2 minutes). Press the push-button on the sensor's box for about 3 seconds when the readings are stable. The sensor will then be calibrated at 20.9%.\* Alternatively the sensor can be connected to a PC running the Logger Sensors software via the USB module. First click on the Module setup icon on the Oxygen sensor's Module box to open its Module setup window. Then click on the **Calibration icon**.

\*This is an assumed stable level in the Earth's atmosphere at sea-level.

#### **Specifications:**

Range and operation modes	ADC resolution	Resolution	Max Sample Rate (S/sec)
0 to 25% (in air)	10 bit	0.1%	100
0 to 100% (dissolved)		0.1%	
0 to 12.5 mg/L (dissolved)		0.1 mg/L	

#### Specifications for the electrode are as follows:

Body diameter:	12mm
Overall length:	150mm
Cap:	16mm OD ±30mm Long
Construction:	Polarographic (Clark Type) design with Silver Anode/Gold Cathode/anode
	system, Delrin body and PTFE membrane
Range/Output:	<ul> <li>0-20 ppm Dissolved Oxygen (0-200% saturation)</li> </ul>
	<ul> <li>Output: 0-40 mA / 0-400 mA</li> </ul>
Response time:	98% of full response in 60 seconds at 25 $^{\circ}$ C

### 8.12 pH logger sensor SNUL-106

This sensor can be used to measure the static pH values of common liquids (water, milk, soft drinks, vinegar, etc.) as well as the changing values in titrations or experiments such as those looking at the effect of antacids.

	pH 1	×	pH 1 - Grap	h	×
рН 1 4.0рН	Öptions Graph Display Table Digital ✓ Graph Grid		Options G Y max Y min. Y-axis Zoom fit	araph c. 14 € 0 € s position 7 ↓ >	

#### Experiment Duration: 1 second to 31 days.

The pH sensor is designed for long life in a variety of general purpose situations. Its sealed reference system and gel fill make it easy to use and maintain. With an epoxy body it is a durable electrode for use both in the laboratory and in the field.

#### Sensor Calibration:

This sensor gives a fast response across the full pH range and can be calibrated with any standard buffer solution.

Connect the sensor to a voltage source (the USB Bridge plugged into a PC , Edu-Logger's® Monitor Display Unit or Battery Unit). Insert the sensor into a pH = 7 buffer and press the sensor's push button for about 3 seconds. The reading is calibrated to 7. Alternatively the sensor can be connected to a PC running the Logger Sensors software via the USB Bridge.

First click on the Module setup button 🧱 on the pH sensor's Module box to open its Module

setup window. Then click on the Calibration icon 2.

#### **Specifications:**

Range and operation modes	ADC resolution	Resolution	Max Sample Rate (S/sec)
0 to 14	16 bit	0.01pH	100

#### Specifications for the electrode are as follows:

Body diameter:	12mm
Overall length:	150mm
Cap:	16mm OD ±30mm Long
Construction:	Epoxy Body, Round Bulb ASG VIII pH Glass, Sealed, Gel-filled Single-
	Junction Reference with fibre frit, Ag/AgCI wire, ATC
Range/Output:	■ 0-14 pH
_	<ul> <li>mV output with isopotential point at 0±20mv at pH 7</li> </ul>
Response time:	98% of full response in 30 seconds at 25 $^{\circ}$ C

### Relative humidity logger sensor **NUL-107** 8.13



Measuring Relative Humidity, this sensor can find use in recording variations with weather conditions and the biological effect on such organisms as seedlings and insects.

	Rel. Humidity 1	×	Rel. Humidity 1 - Graph	
Rel. Humidity 1 53.7%	Options       Graph         Display       Table         Digital       Graph         Grid       Grid		Options       Graph         Y max.       95 €         Y min.       0 €         Y-axis position       -7 ↓         Joint Color       Zoom fit	

#### Experiment Duration: 1 second to 31 days.

It is located in a plastic box with exposure of the sensor being through a hole in the side.

Range and operation modes	ADC resolution	Accuracy	Resolution	Max Sample Rate (S/sec)
0 to 95% RH	16 bit digital	±5% RH	0.1%	100

#### Heart rate & pulse logger sensor 🜌 NUL-108 8.14



This sensor can be used to monitor and compare pulse rates under various exercise and rest conditions and to compare the "normal' and "after exercise" pulse rates. Additionally, it can show how blood volume/flow rates in the finger or ear lobe vary with time.

1 0100	iraph 🛛
1     280 Arb     Options Graph     Options     Ymax       Display     BPM     Wave     Ymax       Pulse     Graph     Graph       1     64 BPM	Graph Graph x. 1000 ÷ n. 0 ÷ s position 100 ↓ S position 100 ↓

Experiment Duration: 1 second to 31 days.

To operate, connect the clip to a finger or ear lobe and start measuring either connected via the USB Bridge to a PC, or to the Monitor Display Unit.

On the PC you can choose to see the pulse wave showing changes of blood volume/flow in the finger or earlobe with time (and calculate the pulse) or get the value of the pulse rate directly via the software. The operating mode is changed by clicking on the sensor's Module Setup icon 🚳 on the Sensor Module box to display the Heart and Pulse sensor module setup window, and selecting the mode as required.

For best results, the sensor should be kept away from direct sunlight and high intensity lights.





Specifications (BPM stands for beats per minute):

Range and operation modes	ADC resolution	Accuracy	Resolution	Max Sample Rate (S/sec)
30 to 240 BPM	10 bit		2	100
0-1024 Analogue Values	TO DIL	1	1	10,000

#### Specifications for the electrode are as follows:

They are both plethysmograph-based and so record changes in blood volume/flow. The sensors consist of an infrared LED transmitter and a matched infrared phototransistor receiver.

# 8.15 Photo gate logger sensor 🌆 NUL-109

This sensor can be used to study various kinds of motion. With four modes of operation, time, velocity or acceleration can be measured with one or two photo gates and associated timing-cards, as well as showing pictorially the status (digital 1 or 0) of the voltage output of the photo gate as timing cards pass through it.



The modes of operation are selected by clicking on the relevant picture.

#### 8.15.1 Digital meter mode

As well as recording Time, Velocity and Acceleration in tables, this mode allows individual measurements to be displayed together with their units and average values.

0.59	[m/s^2]
Average = 0.63	[m/s^2]

#### 8.15.2 Time and Velocity mode



The displayed value can be either the Time it takes for the timing-card to pass through the photo gate, or its Velocity. The timing-card's (flag) length must be entered in the box above its image. Repeated passes can be made and their values recorded, together with an average value.

#### 8.15.3 Acceleration mode 1



The displayed result will be the acceleration of the double timing-card passing through the photo gate. The timing-cards' (flag) lengths must be entered in the box above their images. Repeated passes can be made and values recorded, together with an average value.

#### 8.15.4 Acceleration mode 2



The displayed result is the acceleration of the timing-card passing through the two photo gates. The timing-card's flag length must be entered in the box above its image.

The ID numbers of the two photo gates used must be selected and entered into the two boxes indicating which is passed through first and which second.

Repeated passes can be made and values recorded, together with an average value.

This mode can be used to measure velocities and momentum in collision as described in the following screen.

Module window	NeuLog <sup>TU</sup> Neuron Sensors Network Technology
Photogate 1 0	
Photogate 2 0	Photogate
	Lisplay Table Metocaty ■ m/s Choose photogates' sequence First Second 100   gram 100   gram
	Samples     Velocity ID1 [m/s] Velocity ID2 [m/s]     Momentum ID1 [kg/m/s]

Fields for the carts' mass are opened too.

The two carts are pushed against each other so they collide after passing through the gates and then return again through the gates.





The table will show the velocity and the momentum of each cart before and after the collision.

#### 8.15.5 Status graph mode



As a timing-card passes through the photo gate, a graph of its Digital Status (1 or 0) against Time is produced. When the beam is interrupted, it displays as **1** indicating a voltage output from the photo gate of near 5V. When the beam is uninterrupted it displays as **0** indicating a voltage output from the photo gate of near 0V.

The photo gate is contained in a strong plastic frame with an infrared light emitting diode (LED) on one side and an infrared-sensitive phototransistor on the other side.

Range and operation modes	ADC resolution	Accuracy	Resolution	Max Sample Rate (S/sec)
Single timing-card with one photogate	16 bit	100 μS	100 μS	10,000
Double timing-card with one photogate	digital			
Single timing-card with two photogates				
Digital status one photogate				

### 8.16 Pressure logger sensor 🥝 NUL-110

This sensor can be used to monitor chemical reactions that involve gases and to investigate both Boyle's Law and the Pressure Law for ideal gases. It can also prove useful in studies of weather phenomena.

	Pressure 1	×	Pressure 1 - Graph	×
	Options Graph Display © kPa Table © psi Display Provide		Options Graph	
Pressure 1 128.3kPa	_ Digital ✓ Graph ✓ Grid	Atm	Y min. Y-axis position 7 X Zoom fit Color	

Experiment Duration: 1 second to 31 days.

The pressure sensor is located in a plastic box. The sensing part is connected to a small stainless steel tube for connection to pressure sources such as a syringe via an adapter.

#### **Specifications:**

Range and operation modes	ADC resolution	Accuracy	Resolution	Max Sample Rate (S/sec)
0 to 7 atm	10 bit	±1%	0.01 atm	100
0 to 100 p.s.i.		20°C-30°C	0.1 p.s.i.	
0 to 700 kPa			0.1 kPa	
0 to 7 bar			0.01 bar	

#### 164

#### 165

### 8.17 Force logger sensor 😵 NUL-111

This sensor can measure the mass to weight relationship and study how different pulley systems affect the effort needed to lift weights. It can also be used to measure push/pull forces and impacts.



Experiment Duration: 50ms to 31 days.

The force sensor is contained in a metal box. There is a hook at the bottom of the box that can be connected to various pulling loads. A simple bumper (for push/impact measurements) could be made and attached using a bolt which is placed through a length of plastic tubing.

The sensor can be hung from a universal laboratory stand via a rod through the hole in its box.

This sensor can be operated either facing upwards, downwards and any intermediate (including horizontal) position.

#### Sensor zeroing

Connect the sensor to a voltage source (Edu-Logger's<sup>®</sup> Bridge, Monitor Display Unit or Battery Unit). To zero the reading, just press the sensor's push-button for about 3 seconds. Alternatively the sensor can be connected to a PC running the Logger Sensors software via

the USB module. First click on the **Module setup** icon **3** on the Force sensor's Module box to open its Module setup window. Then click on the **Calibration icon**.

Range and operation modes	ADC resolution	Resolution	Max Sample Rate (S/sec)
±10 N ±50 N	10 bit	0.01 N	3000

## 8.18 Sound logger sensor NUL-112

This sensor has two modes of operation. In slow mode it can be used to measure Soundpressure level in decibels. In fast mode it can be used to compare different sources of sound, their waveforms can also be displayed. The frequencies of tuning forks and wind-chimes could be determined and simple electronic signal generators calibrated using it. With two sound sensors the velocity of propagation of sound in various media could be determined by timing a pulse travelling between them.

Experiment Duration: 25 milliseconds to 31 days.

The sound sensor is located in a plastic box accessible to the atmosphere via a hole in its side.

Range and operation modes	ADC resolution	Accuracy	Resolution	Max Sample Rate (S/sec)
Level: 40 to 110 dB	10 hit	±2 dB	0.1 dB	100
Signal: 0-1024		1	1	10,000

## 8.19 Motion logger sensor MUL-113

This sensor uses an ultrasonic transducer to transmit an ultrasonic wave and measure the time of the echo return. In this way, the sensor measures the distance to an article located against it.

Using the module software, it is able to calculate also the item velocity and acceleration.

The sensor has three modes of operation.

Motion 1 0,16m/s^	·2			
Motion 1 0.382 m/s Motion 1 0.229 m	Motion 1 - Options Options Graph Display Table Digital Graph Grid	™ m/s m/s^2	Motion 1 - Graph           Options         Graph           Y max.         Y min.           Y-axis position         Zoom fit	

Experiment Duration: 1 second to 31 days.

#### **Specifications:**

Range operation	e and n modes	ADC resolution	Resolution	Max Sample Rate (S/sec)
Distance:	0.25-6m	10 bit	2mm	100
Velocity:	± 10 m/s		m/s	
Acceleration:	± 100 m/s2		0.01 m/s2	

#### Note:

The measurement is based on ultrasonic waves, a sound wave humans cannot hear. The wave is not narrow. Echo can be received from bodies near the line between the motion sensor and the measured article.

The size of the measured article must be at least 10 x 10 cm.



# 8.20 Magnetic logger sensor **NUL-114**

This is a very sensitive magnetic fields sensor. It can measure a very low level of magnetic fields such as the magnetic field of Earth.

The logger sensor has only one range and measures the magnetic fields in milli Tesla (mT)

	Magnetic 1 - Options	×	Magnetic 1 - Graph
Magnetic 1 -0.029mT	Options Graph Display □ Table □ Digital ☑ Graph □ Grid		Options       Graph         Y max.       10 €         Y min.       -10 €         Y-axis position       -7 €         Zoom fit       Color

Experiment Duration: 25 milliseconds to 31 days.

Range and operation modes	ADC resolution	Resolution	Max Sample Rate (S/sec)
±10 mT	10 bit	0.001mT	3000

## 8.21 Conductivity logger sensor **1** NUL-115

This logger sensor is based on a probe with two flat electrodes with known surface area and distance between them. A signal is supplied to the electrodes and by testing the signal behaviour, the conductivity of the solution is calculated.

The logger sensor has three ranges for displaying the solution conductivity:

 $\mu$ s/cm – micro Siemens per centimetre mg/L – milli gram per Liter ppm – part per million

Conductivity 1 0.0 ppm				
Conductivity 1 0.0 mg/L © Conductivity 1 0.0 us/cm	Conductivity 1 - Options Options Graph Display Table Digital Graph Grid	© us/cm © mg/L ⊙ ppm	Conductivity 1 - Graph Options Graph Y max. Y min. Y-axis position Zoom fit Color	20000 € -300 € -7 ()

Experiment Duration: 1 second to 31 days.

Range and operation modes	ADC resolution	Resolution	Max Sample Rate (S/sec)
0-20000 μs/cm	10 bit	0-2000 μs/cm - 0.1 μs/cm Over 2000 - 1 μs/cm	100
0-18000 mg/L		0-1000 mg/L - 0.1 mg/L Over 1000 - 1 mg/L	
0-18000 ppm		0-1000 ppm -0.1 ppm Over 1000 -1 ppm	

# 8.22 Spirometer logger sensor MUL-116

The Spirometer enables measuring the volume of our lungs. The sensor includes a tube and it measures the air flow that passes through it. The volume (in litres) is calculated by the software area calculation function.

The tube has a narrow part in its center and it measures the flow rate by measuring the pressure difference between the two ports of the tube.



Experiment Duration: 1 second to 31 days.

#### **Specifications:**

Range and operation modes	ADC resolution	Resolution	Max Sample Rate (S/sec)
±10 L/s	10 bit	0.05 L/s	100

#### Note:

Disposable laminated papers are included with the sensor.

A paper must be rolled and put in the tube before using it and exhaling through it.

# 8.23 GSR logger sensor **NUL-117**

The Galvanic Skin Response (GSR) logger sensor measures the conductivity of our skin, especially between our hand fingers.

The conductivity of our skin changes according to unconscious emotion effects such as sudden noise, smell, touch, pain or view.

This sensor has two ranges: conductivity in micro Siemens and arbitrary numbers.

GSR				
1 0.0000 us (%)	GSR 1 - Options Options Graph Display Table Digital Options Graph Options Graph Options Graph	gnal Y	- Graph ns Graph max. 65500 min. 0	× •
1 Garb	✓ Graph ☐ Grid	Y-	exis position -7	••

Experiment Duration: 1 second to 31 days.

#### **Specifications:**

Range and operation modes	ADC resolution	Resolution	Max Sample Rate (S/sec)
0-65279 Arbitrary units 0-10 μS	10 bit	10nS	100

#### Note:

The skin response time from the sudden effect is between 0.1 to 0.5 seconds. The level of the response changes dramatically from one person to another.

The user must put his hand on a table, chair or on his lap and be still.

# 8.24 Colorimeter logger sensor **Z** NUL-119

The colorimeter measures the levels of colour components (red, green and blue) of a solution.

Another option is to measure the solution absorbance of each colour component.

The colorimeter logger sensor has an opening for a special square solution tube. The colorimeter turns on 3 different lights (RGB) in known values and measures the received light that passes through the solution.

This sensor has two modes of operation.

Colorimeter			
1	Colorimeter 1 - Options	×	🛛 Colorimeter 1 - Graph 🛛 🕅
0.57 Absr	Options Graph		Options Graph
Colorimeter 1 27.04%T	<u>Display</u> ☐ Table ☐ Digital ✔ Graph ☐ Grid	<ul> <li>● Red ● %T</li> <li>● Green ● Absorbance</li> <li>● Blue</li> </ul>	Y max. Y min. Y-axis position Zoom fit Color

Experiment Duration: 1 second to 31 days.

#### **Specifications:**

Range and operation modes	ADC resolution	Resolution	Max Sample Rate (S/sec)
R.G.B. 0-4 absorbance	10 bit	0.01 % 0.01 abs	100

#### Note:

The tube must be put in a certain way inside the colorimeter. The tube has one rough side. This side should not be in the line of the light.

# 8.25 CO₂ logger sensor NUL-120

This logger sensor is based on an electromechanical reaction between  $\text{CO}_2$  gas and the sensor.

The result of the electromechanical reaction is voltage, measured by the logger sensor.

The measure units are ppm (parts of CO<sub>2</sub> per million of air).



Experiment Duration: 1 second to 31 days.

#### **Specifications:**

Range and operation modes	ADC resolution	Resolution	Max Sample Rate (S/sec)
350-10,000 ppm	10 bit	1 ppm	100

#### Note:

The sensor must be calibrated before every use according to the following simple procedure:

Connect the sensor to the PC through a USB module or to a battery module.

Wait for two minutes to warm the sensor and press the pushbutton switch on the CO2 logger sensor for 3 seconds.

# 8.26 Barometer logger sensor MUL-121

This sensor measures atmospheric barometric pressure using a barometer pressure sensor in it.

This sensor has four common ranges for displaying atmospheric pressure:

KPa – Kilo Pascal Atm – Atmospheres in/hg – Mercury in inches

The fourth range is the altitude (in meters).

The highest barometric pressure is on the sea level – when we go up, the pressure goes down.

Barometer 1 61.3m				
Barometer 1 29.67 in/hg				
Barometer 1 0,99 atm	Barometer 1 - Options Options Graph Display Table Digital Graph	● kPa ● atm ● in/hg ● Altitude	Barometer 1 - Graph       Options     Graph       Y max.     106 ÷       Y min.     80 ÷	
1 100.6 kPa	Grid		Zoom fit Color	

Experiment Duration: 1 second to 31 days.

Range and operation modes	ADC resolution	Resolution	Max Sample Rate (S/sec)	
80-106 KPa	10 bit	0.1 KPa	100	
0.80-1.00 Atm		0.01 Atm		
23.60-31.30 in/hg		0.01 in/hg		
-380-1950 m		0.1 m		

# 8.27 Blood pressure logger sensor W NUL-122

The sensor measures the pressure in the air pillow wrapped on the tested person's arm.

Heartbeats affect the blood pressure. That is the reason for the difference between the systolic and diastolic pressure of the tested person.

The sensor has three ranges:

- 1. The average pressure in the air pillow in mm/Hg.
- 2. The pressure beats with no units.
- 3. The sum of the two above signals.

The two parameters of the blood pressure are the average pressure where the pressure beats are 10% of their maximum value.

Blood pressure 1 0.31 Atmosta 0.31						
Blood pressure	Blood pressure 1 - Options	X	Blood pre	ssure 1 - Granh I		×
100.04	Options Graph			Graph		
Blood pressure	Display ☐ Table ☐ Digital ✔ Graph ☐ Grid	<ul> <li>mmHg</li> <li>Arb</li> <li>Arb+mmHg</li> </ul>	Y m Y m Y-az Zoom	ax. in. is position * fit Color	250 € 0 € 8 • ►	

Experiment Duration: 1 second to 31 days.

Range and operation modes	ADC resolution	Resolution	Max Sample Rate (S/sec)
0-250 mm/Hg	10 bit	0.30 mm/Hg	
0-250 Arb		0.30 mm/Hg	
0-250 mm/Hg + Arb		0.30 mm/Hg	

## 8.28 Flow logger sensor MUL-124

This sensor measures water flow. It includes a rotation wheel that rotates when water flows through it. The sensor has inlet and outlet pipes.

The wheel is floating on a bearing and is not connected mechanically to anything else. Its speed is measured by magnetic field changes.

	Flow rate 1 - Options	Flow rate 1 - Graph	×
	Options Graph	Options Graph	
Flow rate 1 0.0000 m/s	<u>Display</u> ☐ Table ☐ Digital ✔ Graph ☐ Grid	Y max. Y min. Y-axis position Zoom fit Color	

Experiment Duration: 1 second to 31 days.

Range and operation modes	ADC resolution	Resolution	Max Sample Rate (S/sec)	
0-4.7 m/s	10 bit	0.0001m/s	100	

# 8.29 Force plate logger sensor NUL-125

This sensor measures heavy weight or forces. We can stand or jump on it.

Handles can be connected to its plates for measuring pulling forces.

	Force plate 1 - Options	×	Force plate 1 - Graph	×
Force plate 1 510N	Dptions] Graph Display ☐ Table ☐ Digital ☑ Graph ☐ Grid		Options       Graph         Y max.       3500 €         Y min.       -800 €         Y-axis position       -7 ↓         ✓       ✓         Zoom fit       Color	

Experiment Duration: 1 second to 31 days.

Range and operation modes	ADC resolution	Resolution	Max Sample Rate (S/sec)	
–800-(+3500) N	10 bit	0.3N	100	


This sensor measures angles, rotation speed or rotation acceleration.

The sensor has a pulley connected to its shaft and the pulley rotation is measured.

This sensor has three modes of operation.

Rotary 1 0,00rev/s					
🎨 💻					
Rotary 1 O <sup>rad/s^2</sup>	t				
۲					
Rotary 1	Rotary 1 - Options		Rotary 1 - Grap	h	×
0.0rad/s	Options Graph		Options Graph		1
Rotary	<u>Display</u> ☐ Table ☐ Digital ☑ Graph ☐ Grid	Angle     Rad/s     Rad/s <sup>2</sup> Rev/s	Y max. Y min. Y-axis po	370 € -10 € sition -7 < ►	
61.71 °			Zoom fit	Color	

Experiment Duration: 1 second to 31 days.

#### **Specifications:**

Range and operation modes	ADC resolution	Resolution	Max Sample Rate (S/sec)
0°-360°	10 bit	0.08°	100
±345 Rad/s		0.6 Rad/s	
±32,222 Rad/s <sup>2</sup>		11 Rad/s <sup>2</sup>	

#### Note:

To zero the measured angle, click on the **Calibration icon** 2.

# 8.31 Acceleration logger sensor SNUL-127

This logger sensor includes a 3D (three dimensions) acceleration sensor, but only one dimension acceleration can be displayed at a time.

This sensor has three modes of operation.

	Acceleration 1 - Options	l	×	Acceleration 1 - Graph	×
	Options Graph			Options Graph	
Acceleration 1 -0.45 <sup>m/s^2</sup>	<u>Display</u> ☐ Table ☐ Digital ✔ Graph ☐ Grid	● X-axis ● Y-axis ● Z-axis		Y max. Y min. Y-axis position Zoom fit Color	

Experiment Duration: 25 milliseconds to 31 days.

#### **Specifications:**

Range and operation modes	ADC resolution	Resolution	Max Sample Rate (S/sec)
Acc-x -80 a 80 m/s <sup>2</sup> Acc-y -80 a 80 m/s <sup>2</sup> Acc-z -80 a 80 m/s <sup>2</sup>	10 bit	0.15 m/s²	3000

#### Note:

The logger sensor measures all the three accelerations at the same time and stores it.

We can run an experiment on-line, upload one dimension acceleration, freeze it and upload another, using the off-line method.



This logger sensor is based on measuring the vacuum pressure in a tensiometer.

A tensiometer is a closed tube with a special ceramic part in its end.

The tensiometer is filled with water and put in the soil. If the soil is dry, water goes out by diffusion through the ceramic holes and vacuum pressure is created in the tensiometer.

When we wet the soil, the vacuum in the tensiometer pulls water into the tensiometer and the vacuum decreases.

This is why soil moisture is measured in pressure level.

	Soli moisture					
	2	Soil moisture 2 - Options	×	🛛 Soil moisture 2 - Graph		X
	-0.02kPa	Options Graph		Options Graph		
	ی ک	Display	⊙ cbar ⊙ kPa	Y max.	50 主	
Ī	Soil moisture	☐ Digital ✓ Graph		Y min.		
	2	🗖 Grid		Y-axis position		
	-0.04cbar			🚺 🔊 🐆		
	۵ 🛑			Zoom fit Color		

Experiment Duration: 1 second to 31 days.

#### **Specifications:**

0.1.1.1.1

Range and operation modes	ADC resolution	Resolution	Max Sample Rate (S/sec)
0-50 cBar 0-50 KPascal	10 bit	0.02 cBar 0.02 KPascal	100

# 8.33 UVB logger sensor NUL-130

There are different wave lengths of ultra violet light.

The UVB wave length range is 280-320 nm, which is 2% of the total UV radiation.

The UVB radiation affects the vitamin generation in the human body, the immune suppression, skin cancer and cataract.

The intensity of this light is measured in mW/m<sup>2</sup> (milli Watt per square meter).



Experiment Duration: 1 second to 31 days.

Range and operation modes	ADC resolution	Resolution	Max Sample Rate (S/sec)
0-1500 mW/m <sup>2</sup>	10 bit	0.8 mW/m <sup>2</sup>	100

### 8.34 Turbidity logger sensor SNUL-131

This sensor measures the reflected light that enters into a tube containing a solution. As solution turbidity is higher, more light is reflected and measured by the light sensor, located perpendicular to the tube.

The turbidity measurement units are Nephelometric Turbidity Unit (NTU).

	Turbidity 9 - Options	×	Turbidity 9 - Graph	X
	Options Graph		Options Graph	
Turbidity 9 29,47 NTU	<u>Display</u> ☐ Table ☐ Digital ✔ Graph ☐ Grid		Y max. Y min. Y-axis position Zoom fit Color	

Experiment Duration: 1 second to 31 days.

Range and operation modes	ADC resolution	Resolution	Max Sample Rate (S/sec)
0-200 NTU	10 bit	0.08NTU	100



There are different wave lengths of ultra violet light.

The UVA wave length range is 320-370 nm, which is 98% of the total UV radiation.

The UVA radiation affects the photo aging and photo chemical smog.

The intensity of this light is measured in mW/m<sup>2</sup> (milli Watt per square meter).



Experiment Duration: 1 second to 31 days.

Range and ADC operation modes		Resolution	Max Sample Rate (S/sec)
0-50,000 mW/m <sup>2</sup>	10 bit	74 mW/m²	100

### 8.36 Surface temperature logger sensor MUL-133

This logger sensor is very similar to the temperature sensor NUL-103, without the stainless steel tube.

The sensor can be put on any surface and even in water.

This sensor has two modes of operation.

	Temperature 1	×	Temperature	1 - Graph		×
	Options Graph	19	Options Gra	aph		
Temperature 1 17.8 °C	Display Table Digital ✓ Graph Grid	Celsius	Y max. Y min. Y-axis p Zoom fit	position Second Color	110 <b>÷</b> 25 <b>÷</b> 7 <b>↓</b>	

Experiment Duration: 1 second to 31 days.

Range and operation modes	ADC resolution	Resolution	Max Sample Rate (S/sec)
–35°C a 150°C –31°F a 238°F	10 bit	0.1°C 0.1°F	100



This sensor is based on a thermocouple sensor that enables measuring very high level and even flame temperature.

This sensor has two modes of operation.

W. R. T.				M
1 0.0652 ∘F ₩. R. T. 1 7318 °C	W.R.T. 1 - Options Options Graph Display Table Digital Graph Grid	<ul> <li>Celsius</li> <li>Fahrenheit</li> </ul>	W. R. T. 1 - Graph Options Graph Y max. Y min. Y-axis position Zoom fit Color	¥ 1200 € -200 € -7 ↓ ►

Experiment Duration: 1 second to 31 days.

Range and operation modes	ADC resolution	Resolution	Max Sample Rate (S/sec)
–200°C a 1200°C –328°F a 2200°F	10 bit	0.1°C 0.1°F	100

## Appendix A – Edu-Logger® Modules

### A.1 Accessory modules

Catalog No.	Module	Module Symbol	Photographic Image	Purpose
900-115	Battery Module			Powers logger sensors when not connected to PC or Monitor Display Unit
900-116	RF Communication Module	(((		Enables wireless connection of logger sensors to PC and Monitor Display Unit
900-114	Monitor Display Unit			Enables programming of logger sensors without the use of a PC, together with digital display of sensors' data
900-150	Graphic Display Unit			Allows Visual experiments
900-113	USB Bridge Module	¥		Connects logger sensors to PC

### A.2 Logger sensors modules

Catalog No.	Module	Module Symbol	Photographic Image	Purpose
900-101	Voltage logger sensor	±20V	AD A	Measures voltage in DC and AC circuits
900-102	Current logger sensor	±2.5A		Measures current in DC and AC circuits
900-103	Temperature logger sensor			Measures temperature
900-104	Light logger sensor	C	Heating .	Measures level of Illumination
900-105	Oxygen logger sensor	02		Measures % oxygen in air and dissolved in water
900-106	pH logger sensor			Measures relative pH
900-107	Relative Humidity logger sensor	o		Measures humidity
900-108	Heart Rate & Pulse logger sensor		6	Measures pulse rate and blood flow
900-109	Photo/Lightgate logger sensor	Ũ		Measures time and, indirectly, speed/velocity, acceleration and momentum.

Catalog No.	Module	Module Symbol	Photographic Image	Purpose
900-110	Pressure logger sensor		Headage	Measures gas or air pressure
900-111	Force logger sensor			Measures forces, both push and pull
900-112	Sound logger sensor			Measures sound level and displays waveforms
900-123	Motion logger sensor		Nedday 	Measures distance, velocity and acceleration
900-124	Magnetic logger sensor	Q	e montante	Measures magnetic field intensity
900-125	Conductivity logger sensor	+ -	-0	Measures solution conductivity
900-126	Spirometer logger sensor			Measures lung air flow and volume
900-127	GSR logger sensor	GSR		Measures Galvanic Skin Response
900-128	Logger EKG sensor			
900-129	Colorimeter logger sensor			Measures solution RGB color transfer and absorbance
900-130	CO₂ logger sensor			Measures CO₂ in air
900-131	Barometer logger sensor			Measures air pressure and altitude

Catalog	Module	Module Symbol	Photographic Image	Purpose
900-132	Blood pressure logger sensor			Measures blood pressure
900-134	Flow Rate logger sensor		6	Measures water flow
900-135	Force plate logger sensor			Measures high weight
900-136	Rotary motion logger sensor			Measures rotary speed, acceleration and rounds. Can also be used as a pulley sensor.
900-137	Acceleration logger sensor		Next.org	Measures 3D acceleration
900-138	Salinity Logger sensor			Measures solution Turbidity
900-139	Soil moisture logger sensor		R	Measures soil moisture
900-140	UVB logger sensor	UVB	Test.by	Measures UVB radiation
900-141	Turbidity logger sensor		Houldong Parkay Carlay	Measures solution turbidity
900-142	UVA logger sensor	UVA	Heatog	Measures UVA radiation
900-143	Surface temperature logger sensor			Measures surface temperature

Comprehensive details of each of these modules is provided in Chapter 8.