

Troubleshooting guide (Biomedicine T5, MolMed HT11)

Problem	Possible causes	Probable solutions
Very few cells in the culture, thus no RNA & poor FACS data	Wrong cell count due to the difficulty in distinguishing RBC from leukocytes	Add RBC lysis buffer step
	Cells were washed away and lost	Culture more cells Use 14mL tubes, which gives firmer and “protected” pellets
No or very low CD93 mRNA expression	Failed reverse transcription (No cDNA)	Check reference gene
	Failed amplification on CD93 (bad assay)	Check known positive controls
	High apoptosis and low neutrophils population (time frame)	Count live and dead cells when we harvest them and compare
	Mixed cell population and different cell types have different CD93 mRNA expression levels	Cell-type-specific culture Grow more cells
One of the SNP assay did not work	Bad assay	Run the products on an agarose gel Positive control from known samples
	No DNA	Check the other assay, it worked.
Bacterial contamination	Not working in sterile lab	Count the cells to see
		Work in sterile area (hood or similar)
		Put the lid on as soon as possible
		Increase antibiotics concentration in the media

	Some indicators (eg pH) to monitor	
	Working with a bacteria incubator	Don't do that.
	Colleagues' extensive use of bacteria culture than past years	Clean it well... or well...
No or degraded RNA	Presence of RNase	Use RNA stabilizing agent or RNase inhibitor at all times
	Not low enough temperature	
	Long processing time	
Low intensity for CD93 in ELISA	Removal of "jelly" may have removed some cells and soluble CD93	Work sterile
	Soluble CD93 stuck on bacteria	Culture more cells
	24h lifespan of neutrophils	Decrease grow time, faster handling time