

User Manual - Manuel d'utilisation
Gebrauchsanweisung - Gebruikershandleiding
Användarhandbok - Brugermanual
Manuale d'uso - Manual de Usuario
Manual do Utilizador

P280 Overlay
P005856 / P005858


162867
Rev.001

User Manual

P280 Overlay
P005856 / P005858

162867
Rev.001

EN

Hill-Rom Industries S.A.
188, Rue du Caducée
34195 MONTPELLIER Cedex 5 - FRANCE
Tel: +33 (0)4 64 04 64 04
Fax: +33 (0) 4 67 04 64 00
www.hill-rom.com

First edition, first printing: June 2010

162867, Rev.001/June 2010

The information contained in this manual is confidential and may not be reproduced or divulged in any form or by any means without the prior written permission of Hill-Rom.

Hill-Rom® is a registered trademark of Hill-Rom Inc.

Velcro® is a registered trademark of Velcro Industries BV.

Terralin® is a registered trademark of Schülke & Mayr GmbH.

Dismozon® is a registered trademark of Bode Chemie GmbH.

Sanivap® is a registered trademark of Sanivap S.A.

Hill-Rom reserves the right to make changes to the design, characteristics and models without prior notice. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

To order copies of this manual, refer to the last page, identify your national Hill-Rom representative and order the article with the part number 162867.

© 2010 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

Summary

Introduction, specifications

The structure of the manual	1
Symbol definitions	2
Applications	3
Description	3
Indication	3
Contraindication	3
Safety and Usage Tips	4
Putting into service	4
Risk prevention	4
Electrical safety	5
Technical specifications	6
Control unit	6
Overlay	7
Conditions for transport, storage and use	7
Product overview	8
General symbols on the top cover	9
General symbols on the control unit	10
Keys and indicator lights	10
Identification label	11
Product identification labels	12
Control unit	12
Upper and lower covers	13
Packaging	13

Installing the patient

Before placing the patient on the P280 Overlay	15
List of compatible bed chassis and mattresses	15
Installation	16

Mobilizing the patient

Functions	19
Operating modes	19


Securing the patient

CPR	21
Safety lock-out	22
Alarms	22
Alarm disable	22
Mains power fault	22
Low pressure	22

Malfunction	22
Moving the bed	
Transport mode	23
Transferring the overlay	25
Packing instructions	26
Disinfection, maintenance	
Disinfection	27
Safety Recommendations	27
Infection control	27
Recommendations for cleaning and disinfection	27
Cleaning and disinfection after the departure of a patient or whenever the linen is changed	29
Cleaning and disinfection at regular intervals or when a serious risk of contamination occurs	30
Maintenance	32
Safety recommendations	32
Preventive Maintenance	32
Troubleshooting	33
De-commissioning	34
Appendix	
Warranty and after sales service conditions	35
Compliance	35
Electromagnetic conformance	36
Complies with electromagnetic emission standards	36
Compliance with electromagnetic immunity	36
Recommended separation distances	38

The structure of the manual


For every type of use, Hill-Rom mattresses provide patients with optimal comfort and greater independence for a feeling of well-being that is conducive to a swift recovery. They are also easy to use for caregivers.

Symbol definitions

This manual contains different typefaces and icons designed to improve readability and increase understanding of its content. Note the following examples:

- Standard text: used for regular information.
- **Bold text:** emphasizes a word or phrase.
- ⓘ highlights special information or explains very important instructions.

The symbols below represent different risks or hazards:

Symbol	Description
	WARNING • This symbol indicates that the failure to follow the associated recommendation can put the patient or the user in danger, or damage the equipment.
	CAUTION • This symbol indicates that the failure to follow the associated recommendation can result in damage to the equipment.
	Tip
	Risk of falling
	Caught hazard warning
	Risk of crushing an upper limb
	Chemical Hazard Warning
	Electric Shock Hazard

Applications

Description

The P280 Overlay is a therapeutic overlay system designed by Hill-Rom to meet the following needs:

- the efficient automatic redistribution of pressure, especially in the head, sacrum and heel zones;
- greater patient safety.

It must be installed on a standard foam type mattress.

Indication

It helps to prevent and treat pressure ulcers in high-risk adult patients, with the recommended weight limits of 30 and 140 kg. It is not suited to patients with atypical anatomies.

It can be used as an overlay in the following environments:

- in hospitals or other medical infrastructures administering short-term care;
- in long-term care establishments;
- at home.

 *In accordance with the NPUAP/EPUAP directives¹, we recommend that the condition of each patient be regularly checked. For patients with special needs, we recommend the use of the more suitable Immersion™ Therapy system. Caregivers are responsible for taking this decision, in accordance with modern care practices.*

Contraindication

The P280 Overlay must not be used with patients suffering from:

- an unstable fracture of the spine.
For any other unstable fractures, a medical examination is necessary to decide whether the use of the P280 Overlay is appropriate;
- cervical or skeletal traction.

1. NPUAP / EPUAP - Pressure Ulcer Prevention - Quick Reference Guide, January 2010
NPUAP / EPUAP - Pressure Ulcer Treatment - Quick Reference Guide, December 2009

Safety and Usage Tips

Putting into service

Before using the mattress, it is essential to have a thorough understanding of this manual. This manual contains instructions for general use and maintenance and guarantees improved safety. Caregivers must have access to this manual.

Caregivers must be informed of the risks that may be encountered in the use of electric devices.

When using the mattress with medical devices (accessories), the user must ensure that security and conformity requirements are met.

Before using the mattress for the first time or when removing from storage:

- check the condition and conformity of the electrical system with the applicable safety standards;
- connect the P280 Overlay system to the mains power supply (See “Electrical safety” page 5);
- make sure that all the functions of the system are in good working order (See “Installation” page 16);
- make sure that the device and the care environment are in a good state of hygiene (See “Disinfection” page 27).

Risk prevention

Improper use of the P280 Overlay system can result in risks for the patient or the user. The following recommendations must be read and followed.

In view of the multitude of models of frames and siderails, and for safety reasons, Hill-Rom advises that all necessary precautions must be taken, especially with regard to the height of the siderails and the dimensions of the sleep surface. If this overlay is used on a bed with siderails that are less than 22 cm above the mattress, patients must not be left unattended.

For safety reasons, it is advisable to use the lock functions:

- whenever intervening on the patient or the device (e.g.: examinations, transfers, maintenance);
- when the patient is in an unusual condition or behaves abnormally (e.g.: excited, confused, disoriented, obsessive, old or of a weak constitution).

It is up to the suitably qualified medical staff to determine the conditions of use of the device and the necessary level of observation or immobilization.

It is imperative to respect all precautions related to the safety of the caregivers. Special attention must be taken when redistributing the load application points, since there is a danger that the bed will tip when the chassis is moved.

Never block the air vents on the base and the back of the control unit. Doing so would prevent efficient ventilation.

Keep the control unit away from liquids.


Check that the control unit is securely attached to the bed panel. Avoid all mechanical impacts.

Electrical safety

The mains power supply for the bed must comply with following standards:

- NF C 15-100 and NF C 15-211 (France);
- International Electrotechnical Commission (IEC) 60364 for other countries.

Check that the supply voltage of the P280 Overlay system shown on the identification label corresponds to the power supply voltage in the establishment (See "Product identification labels" page 12).

The P280 Overlay system should be connected to a power system equipped with a maximum 30 mA earth leakage circuit breaker, in compliance with IEC 60364-5-53.

Connect the device to the nearest wall power outlet so that the shortest possible length of power cord is on the floor.

In compliance with standards relating to electromagnetic interference for medical equipment, this product does not interfere with other medical devices or is not susceptible to interference when combined with other medical devices that also comply with the electromagnetic standards in force.

Some devices, particularly older ones that do not comply with the electromagnetic compatibility standards, may however undergo interference or may themselves interfere with the working of this product.

The users of such devices are responsible for ensuring that any malfunctions will not endanger the patient or any other person.

When direct intravascular or intracardiac connections are in use, the electric potentials of all the unprotected metal parts need to be equalized.

This label indicates that oxygen tents must never be used and that only the use of nasal tubes and oxygen masks is authorized. For reasons of safety, masks and tubes should always be kept at a higher level than the sleep surface.


Technical specifications

i *Hill-Rom has an continuous improvement policy for its products.
Therefore specifications are liable to be altered without notice.*

Control unit

Characteristics	Values	
Model	P005856	P005858
Power supply	220-240 V~, 50 Hz, 0.07 A	120 V AC, 50-60 Hz, 0.17 A
Nominal power	9.64 W	10.28 W at 50 Hz 8.46 W at 60 Hz
Dimensions	29.1 x 20 x 11.7 cm	29.1 x 20 x 11.7 cm
Weight	2.2 kg	2.2 kg
Sleeve material	ABS PC	ABS PC
Fuse	T1A, 250 V	T1A, 250 V
IEC 60601-1 classification	Class II with functional earth terminal	Class II with functional earth terminal
Degree of protection against electric shock	Type BF	Type BF
Degree of protection provided by the sleeve (IEC 60529)	IP20	IP20
Protection against inflammable anesthetic mixtures	Not for use with flammable anesthetics.	Not for use with flammable anesthetics.
Operation of the appliance	Continuous/alternating	Continuous/alternating
Noise level (ISO 3744)	<45 dB(A)	<45 dB(A)
Compression	0-105 mbar	0-105 mbar
Compressor flow rate	10 l/min	10 l/min

Overlay

Characteristics	Values
Dimensions	 203 x 90 x 10 cm
Weight	About 7.8 kg
Cover material	Polyurethane
Air bladder material	Polyurethane
 Therapeutic weight limit	The maximum and minimum permissible patient weight for guaranteed clinical efficiency in all the usual positions of the adjustable head section, excluding the essential bedding accessories.
Degree of protection provided by the sleeve (IEC 60529)	IP24
Self-sufficient	10 hours (transport mode)

Do not use the device outside the therapeutic weight range.

Conditions for transport, storage and use

Symbol	Features	Use	Transport, storage ^a
	Temperature	+10°C - +40°C	-18°C - +60°C
	Humidity	10% - 90%	10% - 95%
	Atmospheric pressure	900 mbar - 1,060 mbar	700 mbar - 1,060 mbar

a. Applicable only if the P280 Overlay is stored in its original packaging.


The device must be stored in its original packaging:

- *protected against light and damp;*
- *at least 10 cm from the floor to avoid any penetration of water;*
- *away from dust;*
- *clear of passages.*


Never stack more than 5 overlays.

Product overview


Item	Name
A	Control unit
B	Detachable upper cover
C	Pneumatic mattress
D	Foam base of the foot section
E	Lower cover
F	Power cord
G	Hose

General symbols on the top cover

	Do not walk on or run over the power cord.
	Foot end
	Annotation zone
	Top of the mattress
	Install on a standard foam type mattress.
	Adjust the straps at the head and foot ends.

General symbols on the control unit

Keys and indicator lights

Keys		
	Symbol	Description
A		Alternating pressure cycles (every 5 minutes)
B		Comfort setting: reduce the pressure
C		Comfort setting: increase the pressure
D		Selection of the operating mode
E		Security lock
F		On/Off
G		Alarm disable

Indicator lights		
	Symbol	Description
1		Maximum inflation (P-max)
2		Alternating low pressure mode
3		Continuous low pressure mode
4		Over-inflation
5		Mains power fault
6		Low pressure alarm
7		Service


Identification label

	Functional earth terminal		Warning: high voltage!
	Manufacturer		Refer to the user manual.
	Compliant with Directive 93/42/EEC		Attention, read the safety instructions carefully
	Alternating Current		DO NOT DISCARD. Abide by the local recovery rules.
	Fuse		Temperature limits
	Class II device		Atmospheric pressure limits
	Type BF device		Hygrometry limits
IP20	Protection code of the control unit sleeve		Catalog part number

Product identification labels

Control unit


Refer to the identification label on the back of the control unit for the device model and its serial number.

P005856	The identification label for P005856 contains the following information: - Hill-Rom Industries S.A. 188, Rue du Caducée - Parc Euromédecine 34195 Montpellier Cedex 5 - France 220-240V ~, 50 Hz, 0.07A - T1A 250V P280 Overlay - REF: P005856 - Modification: - CE 0459 - IP20 1060 mbar 900 mbar 10% 90% +10°C +40°C 15
P005858	The identification label for P005858 contains the following information: - Hill-Rom Industries S.A. 188, Rue du Caducée - Parc Euromédecine 34195 Montpellier Cedex 5 - France 120V ~, 50/60 Hz, 0.17A - T1A 250V P280 Overlay - REF: P005858 - Modification: - CE 0459 - IP20 1060 mbar 900 mbar 10% 90% +10°C +40°C 15

Upper and lower covers

Open the zip fastener of the overlay to access the identification labels of the upper and lower covers.

Packaging


To identify the model on the packaging label:

P005856	P005858
REF P005856 Quantity: 1 (Q) P280 Overlay Hill-Rom CE 0459 SN L123AV12345 COUNTRY OF ORIGIN / MADE IN FRANCE Hill-Rom Industries S.A. 18, RUE DU CADUZE 34116 MONTPELLIER CEDEX 5 FRANCE	REF P005858 Quantity: 1 (Q) P280 Overlay Hill-Rom CE 0459 SN L123AV12345 COUNTRY OF ORIGIN / MADE IN FRANCE Hill-Rom Industries S.A. 18, RUE DU CADUZE 34116 MONTPELLIER CEDEX 5 FRANCE

Installing the patient

Before placing the patient on the P280 Overlay

Assess the various risks, and the following aspects in particular (incomplete list covering risks due to reasonably foreseeable misuse):

- *risks of entrapment;*
- *potential falls from the bed;*
- *patient in state of confusion;*
- *patient's learning ability;*
- *children (aged less than 12, or less than 1.46 m tall);*
- *persons lacking the mental capacity to recognize risky actions;*
- *unauthorized persons.*

List of compatible bed chassis and mattresses

Check that the P280 Overlay is compatible with the bed/mattress combination, and their respective dimensions in particular.

i *The dimensions of the P280 Overlay are contained in the Technical Specifications of the Overlay page 7.*

In order to meet the safety demands defined in the “Hospital Bed Safety Workgroup” guide and IEC standard 60601-2-52, the P280 Overlay can be used in combination with the following bed chassis and foam mattresses:


- Hill-Rom AvantGuard 800;
- Hill-Rom AvantGuard 1200;
- Hill-Rom AvantGuard 1600;
- Hill-Rom NP50, single-density foam mattress;
- Hill-Rom NP100, dual-density foam mattress;
- Hill-Rom NP150, visco-elastic foam mattress.

Installation

1. Unpack the control unit and the overlay.

 **Take care not to damage them when unpacking.**

2. Check that all the components are present and correct and that the power cord is not damaged.
3. Place the rolled up mattress on the top of the bed and unroll it.
4. Check that the symbol  on the cover is at the foot of the bed.

Make sure that the overlay is not used directly on the bed's sleep surface. This device is designed to be used on top of a mattress.

5. Attach the overlay using the straps at the head and foot ends.
6. Adjust the length of the straps to safely secure the overlay.

Make sure that the device is correctly installed and attached, well centered on the sleep surface and securely fixed to the board at the foot end to avoid the danger of entrapment.

Make sure that the attachments do not become entangled in the mobile parts of the bed chassis, such as actuators, CPR handles, etc. With articulated chassis, make sure that the straps of the overlay are only attached to the mobile head and foot sections and NOT to the main fixed chassis.

7. Attach the control unit to the foot board of the bed using the universal attachment support.

8. Connect the hose to the control unit.
The connector clicks when correctly in position.

9. Connect the power cord to the control unit.

10. Use the strap loops (on the right of the overlay) and the Velcro® straps provided to keep the entire length of the overlay in place.

11. Connect the power cord to the wall outlet.

With articulated chassis, make sure that the power cord is not attached to any mobile parts or pinched, in order to avoid the risk of cutting the cord.


12. Switch on the main switch on the right of the control unit.

13. Press  to start the system.

The light comes on.

When connected to the mains power for the first time, the overlay will automatically inflate to maximum inflation (P-max) to fully expand the air bladders. It then automatically returns to the default mode: comfort setting "5", "Alternating low pressure" mode, 10 minute cycle.

Do not place the patient on the mattress while the red indicator light is on during the initial inflation phase.

Note the position of the CPR valve. The CPR valve must remain accessible for use in the event of emergency deflation.

-  ***Placing a cotton sheet on the overlay will improve patient comfort and facilitate care.***

Mobilizing the patient

i This overlay is designed for optimal therapeutic benefits when the AHS is inclined between 0° and 45°. Performance will diminish at angles greater than 45°.

Functions

Operating modes

Select the operating mode by pressing  until the light corresponding to the required mode comes on.

Alternating low pressure

In alternating low pressure mode, the internal pressure of the different zones of the overlay varies, creating zones of high and low pressure under the patient's body.

When this operating mode is selected, the light near the  symbol comes on.

Press  to select an alternating low pressure cycle. Select a 10, 15, 20 or 25 minute cycle for optimal patient comfort and according to the required result.

Continuous low pressure

In continuous low pressure mode, the pressure at the interface with the patient's body is uniform.

When this operating mode is selected, the light near the  symbol comes on.

Comfort settings

The pressure in the air bladders of the P280 Overlay can be adjusted in all the operating modes using the comfort settings keys:  decreases the pressure and  increases the pressure.

Adjust the comfort settings regularly according to the condition of the patient.


To find the right comfort setting:

1. Fully inflate the overlay in the required operating mode.
2. Gradually decrease the pressure until the patient feels comfortable, but without reaching the lowest pressure level.

The pressure of the overlay can also be checked manually:

1. Slide your hand under the air cushions beneath the patient's sacrum, with the palm facing upwards.
2. Raise your hand: there must be a 2 or 3 cm gap under the sacrum.
3. If this is not the case, increase or decrease the pressure accordingly.

i For reference only, in theory, the level 1 comfort setting corresponds to 30 kg and level 10 corresponds to 140 kg.

Over-inflation

In order to provide added support for the patient when seated, the pressure can be increased slightly above the selected level in the alternating low pressure and continuous low pressure modes.

The light near the  symbol comes on, in addition to the light of the selected operating mode.

Maximum inflation (P-max)

When connected to the mains power supply for the first time, the P280 Overlay automatically inflates to P-Max and the light near the  symbol comes on.

It then automatically returns to the default operating mode (comfort setting: 5, alternating low pressure, 10 minute cycle).

i The overlay can be fully inflated manually when the patient is being moved or repositioned or wants to leave the bed. The overlay then returns to the previously selected operating mode after 20 minutes.

Securing the patient

CPR

Never allow a non-qualified person to operate this function and check that no obstacles (e.g., limbs, accessories, objects, power cables) or persons are under the head section.

- i** The CPR valve can be easily identified by the red and yellow hose connector on the left of the control unit.

1. Disconnect the hose connector from the control unit.

The overlay takes about 15 seconds to deflate.

2. Lower the head section of the bed, or place the bed chassis in CPR position (refer to the bed chassis manufacturer's instructions).
3. Place a CPR panel under the patient, or follow the CPR function protocols.

Cancel the CPR

1. Connect the hose connector to the control unit.

The overlay returns to the preceding comfort settings.

2. Return the bed chassis to the appropriate position (refer to the bed chassis manufacturer's instructions).

Safety lock-out

This function locks the user interface to prevent the operating mode from being changed unintentionally.

1. Press  until the light comes on to lock the user interface.
2. Press  until the light goes out to unlock the user interface.

 *If no functions are used for 5 minutes, the user interface is automatically locked.*

Alarms

Alarm disable

To stop the “Mains fault” and “Low pressure” alarms from sounding, press .

 *If the problem continues, the alarm sounds again after 3 minutes.*

Mains power fault

If the P280 Overlay is disconnected or there is a mains power failure, an audible alarm sounds and the light near the  comes on.

 *The audible alarm and the light are active when the P280 Overlay is being transported (See “Transport mode” page 23).*

Low pressure

When the pressure in the air bladders is too low, an audible alarm sounds and the light near the symbol  comes on (refer to the Service Manual P/N: 162879).

Malfunction

In the event of a malfunction of the P280 Overlay system, an audible alarm sounds and the light near the symbol  comes on (refer to Service Manual P/N: 162879).

Moving the bed

Transport mode

The following operations must be completed **before** disconnecting the P280 Overlay from the mains and transporting the patient on the bed:

1. Press  to select Continuous low pressure mode.

The light of the symbol  comes on.

2. Wait until the pressure is perfectly stable and the compressor stops. This must not take any longer than 5 minutes.

These two steps guarantee that the system deflates in a uniform manner.

3. Switch off the control unit by pressing .

The light of the symbol  goes out.

If an electronic module is disconnected when the system is switched on, the module may be destroyed.

The Mains fault alarm sounds and the  light comes on. The other lights are off.

4. Deactivate the Mains fault audible alarm by pressing  (See “Alarm disable” page 22).
5. Disconnect the P280 Overlay from the mains by removing the plug from the wall socket.


Never pull on the power cord or you may damage it. A damaged power cord is an electric shock hazard.

The P280 Overlay automatically switches to transport mode and remains inflated.

6. Safely stow the power cord.

Make sure that the cord does not drag on the ground, so that the bed chassis does not run over it and so that no-one trips up on it. Where necessary, use the attachments provided with the overlay.



When there is a risk of the control unit being damaged, or when there is a mains power failure:

1. disconnect the hose from the control unit;
2. place the transport cap on the hose as quickly as possible to prevent the air escaping from the overlay.

The cap clicks into place.

The P280 Overlay remains inflated for 10 hours when the transport cap is in place.


Make sure that the overlay does not remain disconnected from the mains for more than 10 hours when a patient is installed.

Reconnect the P280 Overlay as soon as the destination is reached. It returns to the preceding mode.

- ① *The CPR function remains available in transport mode. If the hose is connected, simply disconnect it from the control unit. Otherwise, remove the transport cap to allow the air to escape.*

Transferring the overlay

To transfer the P280 Overlay from one bed to another:

Never leave the patient on the overlay.

Make sure that the bed chassis brakes are applied to prevent any accidental movements.

1. To prepare the P280 Overlay for transfer, refer to steps 3 to 5 of Transport mode (page 23).
 2. If necessary, detach the Velcro® straps securing the power cord to the bed chassis.
 3. Undo the two straps securing the overlay to the bed and the mattress (at the head and foot ends).
 4. Place the control unit on the overlay to prevent it from being damaged.
- ⓘ *If there is no mattress on the other bed, transfer the base mattress with the P280 Overlay.*
5. To install the P280 Overlay on the other bed, refer to steps 2 to 13 of Installation (page 16).

Packing instructions

Make sure that the bed chassis brakes are applied to prevent any accidental movements.

1. Switch off the P280 Overlay by pressing .

The light goes out.

2. Disconnect the hose from the control unit.
3. Disconnect the P280 Overlay from the mains by removing the plug from the wall socket.


Never pull on the power cord or you may damage it. A damaged power cable is an electric shock hazard.

The Mains fault alarm sounds and the  light comes on. The other lights are off.

4. Deactivate the Mains fault audible alarm by pressing  (See “Alarm disable” page 22).
5. Detach the Velcro® straps securing the power cord to bed.
6. Remove the cord from the retaining loops of the straps.
7. Undo the two straps securing the overlay to the bed and the mattress (at the head and foot ends).
8. Starting at the head end, slowly roll up the overlay to allow the remaining air to escape.

Never fold the mattress lengthways in order to conserve the life time of the footrest foam.

9. Store the P280 Overlay (See “Conditions for transport, storage and use” page 7).


Disinfection, maintenance

Disinfection

Safety Recommendations

- *Check that the brakes are applied on the bed on which the overlay is installed.*
- *Lock out all electrical functions.*
- *Disconnect the overlay and stow the power cord.*
- *Check that all the connectors are tight to prevent any water ingress.*
- *Do not use too much water or high-pressure jets to wash the overlay.*
- *Never use water at a temperature of more than 60°C.*
- *Avoid excess water on the connector jacks.*
- *Refer to the recommendations of the cleaning and disinfecting product manufacturer.*
- *Thoroughly dry before reusing.*

Failure to implement one or more of these recommendations may lead to damage or deterioration, preventing use of the overlay and rendering the warranty void.

Infection control

Insufficient cleaning = Risk of infection (biological danger)!

All parts must be kept clean at all times in order to avoid the risk of infection. All necessary precautions must be taken to remove all visible stains and soiling.

Recommendations for cleaning and disinfection

The following recommendations are not designed to replace existing cleaning protocols drawn up by the hygiene officer or by other bodies for your hospital.

The cleaning and disinfecting method described below applies specifically to the overlay and its accessories and is designed to save time and to help combat nosocomial infection more effectively.

Hill-Rom recommends that the overlay should be disinfected before first use and before scrapping, in accordance with applicable local regulations.

Regularly check the cover for cuts, tears, cracks or snags. Never use an overlay with a damaged cover.

If the internal foam in the heel zone is contaminated and presents a risk of contamination, it must be replaced.

Suitable individual protective gear must be worn during the phases of the cleaning operations (blouse, gloves, eye protection, etc.). Refer to the instructions of the cleaning and disinfecting products.

i *Cleaning and disinfection are two separate procedures.*

The following products should not be used.


Never use cleaning agents, detergents, degreasing agents or industrial solvents containing any of the following products in order to avoid damaging the overlay:

	Phenol		Hydrochloric, nitric or sulfuric acids
	Cresol		Soda

Never use highly alkaline products or highly concentrated detergents.

Never use abrasive cleaning materials or products, such as scouring pads.

Recommended cleaning and disinfection

Parts to be kept clean at all times.

- After the departure of a patient
- Whenever the linen is changed
- If inspection reveals ingress of fluids

Cleaning and disinfection after the departure of a patient or whenever the linen is changed

Recommended products

Hill-Rom recommends the following product:


Terralin® Protect (Schülke)

Concentration: 5 ml/g (=0.5%), contact time: 1h

Cleaning

Clean the overlay with a cloth dampened with hot water and a cleaning/disinfecting product approved by the establishment.
Do not use excessively liquid, acid or abrasive products.

Check that the hose connector is securely connected to avoid any damage to the PCB by ingress of liquid.

Cleaning tough stains


i Quickly wipe away any traces of pharmaceutical solutions used for the patients to avoid damage to the surface.

- Remove tough stains with a standard textile stain remover or Javel water in a concentration of less than or equal to 1,000 ppm and a soft brush.
- To loosen heavy, dried-on soil or excreta, soften by saturating.

Disinfection

When there is visible soiling, Hill-Rom recommends that the overlay be disinfected with an intermediate disinfectant (tuberculocidal) in accordance with the applicable regulations (e.g. the requirements of directive 93/42/EEC).

Follow the recommendations below for all other disinfectants:

	It is possible to use chlorine-based solutions: concentration must be equal to or less than 1,000 ppm.
	Ethanol (alcohol) can be used: concentration must not exceed ¼ of ethanol for ¾ of water.

Cleaning and disinfection at regular intervals or when a serious risk of contamination occurs

Follow the same instructions as above (See “Cleaning and disinfection after the departure of a patient or whenever the linen is changed” page 29), but with the products shown below.

Recommended products

After a patient with an infectious disease:

- **Dismozon® Pur (Bode)**

Concentration: 40g/l (= 4%), contact time: 1h

After a patient with Clostridium difficile:

- **Dismozon® Pur**

Concentration: 15g/l (= 1.5%), contact time: 2h

- **Anioxy-Spray WS (Anios)**

Ready-to-use, contact time: 30 min


- **Javel water**

Concentration $\leq$ 1,000 ppm, contact time: 10 min

Dry steam cleaning

This overlay can be dry steam cleaned. However, in order to avoid any damage or deterioration caused by high pressure or abnormal surface temperature, the following precautions should be taken:

- Avoid excess water and use low-pressure steam on the electric parts.
- Do not use accessories such as high-pressure hoses (A). It is preferable to use soft non-metallic brushes (B) and wipers (C) in such a way as to reduce the pressure to an acceptable level.


- Avoid water and steam getting into connectors that are not in use.
- Do not brush.
- Use low pressure on the labels and markings.
- Carefully dry and test the overlay before reuse.

Hill-Rom recommends using the Sanivap® method (values recommended by the supplier):

- application time: about 1 second;
- distance between the source and the surface to be disinfected: about 3 mm;
- application: exterior parts of the mattress (covers);
- frequency: between each patient.

Cleaning the upper cover

The upper cover of the overlay can be machine-washed. Follow the instructions below to prepare for machine washing.

Do not always machine-wash, as doing so will reduce the life time of the components. The overlay should be machine-washed if particular risks of infection occur.


	Maximum temperature: 70°C. Delicates.
	Only clean with oxygenated agents.
	Do not iron.
	Do not dry clean.
	Can be tumble-dried. Moderate temperature: 50°C.

Final steps

- *Always remove all traces of the products used when washing or disinfecting from the overlay.*
- *Make sure that all the parts of the overlay are perfectly dry before installing.*

Maintenance

Safety recommendations

Only facility-authorized personnel should perform maintenance.

Before maintenance or servicing works:

- Check that the brakes are applied on the bed on which the overlay is installed.
- Lock out all electrical functions.
- Unplug the overlay.
- Secure the sleep surface and take whatever steps are necessary to prevent any movement.

For the complete maintenance procedure, refer to the P280 Overlay Service Manual. Contact your local Hill-Rom representative (see back of manual).

Preventive Maintenance

NOTE:

The frequency of inspections must be adapted to the general condition of the product and its use, for example, if the overlay is used by heavy patients. It is the responsibility of the facility to implement a preventive maintenance program for the overlay's functions under its conditions of use.

The overlay and accessories should be inspected at least once a year to keep it in good condition and working properly.

Every three years, it is preferable to ask Hill-Rom After-Sales Service or a Hill-Rom approved supplier to inspect the overlay in order to keep it in safe and good working order over time. Depending on the maintenance operations and observations, the date of the next inspection must be recommended every time the bed is serviced.

In order to benefit from optimal and rapid service when calling Hill-Rom about your P280 Overlay, provide the serial number of the P280 Overlay about which you are calling (See "Product identification labels" page 12).

Troubleshooting

The P280 Overlay is designed to function automatically. Consequently, troubleshooting of the control unit is limited to a few checks.

Always disconnect the overlay before troubleshooting.

If an electric power supply problem occurs, check the following:

- the connection to the mains is OK;
- the power cord is not damaged;
- the fuse is in good condition.

If a pressure problem occurs, check the following:

- the hose is intact and securely connected to the control unit and the air bladders;
- the overlay is not damaged (tears or cuts);
- the air filter on the rear of the control unit is clean.

i Refer to the *P280 Overlay Service Manual* for detailed troubleshooting instructions.


Replacing the fuse

1. Remove the cover of the fuse box on the rear of the control unit using a small screwdriver.

Check that the new fuse is type T1A, 250 V and meets the VDE standard.

2. Replace the old fuse.
3. Close the fuse box properly.

Cleaning and replacing the air filter

1. Remove the air filter cover on the rear of the control unit.

The air filter can be cleaned with water and a cleaning agent and re-used.

2. Thoroughly clean and dry the air filter.
3. Install the filter and its cover.

i The air filter can be ordered separately. For more information, refer to the list of spare parts of the P280 Overlay (P/N: 162883).

If the problem is not solved, contact your local Hill-Rom representative (see back of manual).

De-commissioning

The P280 Overlay and its accessories should be cleaned and disinfected before de-commissioning.

Decommissioned equipment materials (plastics, electrical components, etc.) must be recycled in accordance with local recycling regulations. Always meet the applicable demands and local rules relating to environmental protection, in particular for waste from medical devices.

Do not dispose of electric and electronic equipment in the waste bin (as per directive 2002/96/EC).

Never discard the batteries or accumulators of your device. They may contain substances and metals that are hazardous for the environment and health (as per Directive 2006/96/EC).

Warranty and after sales service conditions

The warranty for our overlays will be rendered null and void, in part or in total, in the event of:

- repairs, installations, assembly, modifications or checks and tests are not conducted by the manufacturer's maintenance personnel or by personnel authorized by the manufacturer;
- the electric system does not provide the conditions allowing for the use of medical appliances that comply with the standard EN 60601-1. In hospital wards in particular, the electric system does not meet the requirements applying to medical establishments;
- the device is not used in accordance with this manual;
- accessories are used that do not meet the requirements of this manual.

Contact details of the after-sales service departments in every country can be found on the back of this manual.

Compliance

The P280 Overlay is a class IIa medical device that complies with the requirements of Directive 93/42/EC and its amendments.

The P280 Overlay is designed and manufactured in accordance with the following standards and classifications:

	Description
Technical and quality assurance standards	EN IEC 60601-1
	EN IEC 60601-1-2
	EN ISO 13485

Electromagnetic conformance

Complies with electromagnetic emission standards

Manufacturer's recommendations and declaration – electromagnetic emissions		
The P280 Overlay is intended for use in the electromagnetic environment specified below. The customer or the user of the P280 Overlay should make sure it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment - Recommendations
RF emissions CISPR 11	Group 1	The P280 Overlay only uses RF energy internally. Consequently, it only produces very weak RF emissions that are unlikely to cause interference with nearby electronic equipment.
RF emissions CISPR 11	Class B	The P280 Overlay is suitable for use in all establishments, including domestic premises and premises directly connected to the public low-voltage power supply in residential buildings.
Harmonic emissions IEC 61000-3-2	Class A	
Emissions of voltage fluctuations/flicker IEC 61000-3-3	Compliant	

Compliance with electromagnetic immunity

Manufacturer's recommendations and declaration – electromagnetic immunity			
The P280 Overlay is intended for use in the electromagnetic environment specified below. The customer or the user of the P280 Overlay should make sure it is used in such an environment.			
Immunity test	Test level IEC 60601	Compliance	Electromagnetic environment - Recommendations
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV on contact ±8 kV in the air	±6 kV on contact ±8 kV in the air	Floors should be made of wood, concrete or ceramic tiles. If the floor is covered with a synthetic material, relative humidity must be at least 30%.
Fast electric transients in bursts IEC 61000-4-4	±2 kV on power supply lines ±1 kV for the input/output lines	±2 kV on power supply lines ±1 kV for the input/output lines	The quality of the electric power supply should be the equivalent of that found in typical commercial or hospital environments.
Shock waves IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±1 kV line(s) to line(s) ±2 kV line(s) to earth	The quality of the electric power supply should be the equivalent of that found in typical commercial or hospital environments.

Manufacturer's recommendations and declaration – electromagnetic immunity			
The P280 Overlay is intended for use in the electromagnetic environment specified below. The customer or the user of the P280 Overlay should make sure it is used in such an environment.			
Immunity test	Test level IEC 60601	Compliance	Electromagnetic environment - Recommendations
Voltage dips, brief outages and variations in voltage on the electric power supply lines IEC 61000-4-11	• < 5% U_T^* (dips >95% of U_T) for 0.5 cycle • 40 % U_T (dips = 60% of U_T) for 5 cycles • 70% U_T (dips = 30% of U_T) for 25 cycles • <5% U_T (dips >95% of U_T) for 5 seconds	• < 5% U_T (dips >95% of U_T) for 0.5 cycle • 40 % U_T (dips = 60% of U_T) for 5 cycles • 70% U_T (dips = 30% of U_T) for 25 cycles • <5% U_T (dips >95% of U_T) for 5 seconds	The quality of the electric power supply should be the equivalent of that found in typical commercial or hospital environments. If the user of the P280 Overlay requires the appliance to be in continuous use during mains power cuts, it is advisable to power the P280 Overlay from an uninterruptible power source or a battery.
Magnetic field at the frequency of the mains power supply (50/60 Hz) IEC 61000-4-8	3 A/m	3 A/m	The magnetic fields at the frequency of the mains power supply should have levels characteristic of a representative place located in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Veff From 150 kHz to 80 MHz 3 V/m From 80 MHz to 2.5 GHz	$[V_1] = 3 \text{ Veff}$ $[E_1] = 3 \text{ V/m}$	Portable and mobile RF communications appliances should not be used too close to any part of the P280 Overlay, including the cables. The recommended distance should be respected, which is calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \left[\frac{3,5}{\sqrt{f}} \right] \sqrt{P}$ $d = \left[\frac{3,5}{E_1} \right] \sqrt{P} \text{ From 80 MHz to 800 MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \text{ from 800 MHz to 2.5 GHz}$ where P is the maximum output power of the transmitter in Watts (W), according to the manufacturer of the transmitter and d is the recommended separation distance in meters (m). The field intensities of fixed RF transmitters, determined by an on-site electromagnetic investigation^o, should be below the level of compliance in each frequency range^o. Interference may occur close to devices identified with the following symbol: 

- a. U_T is the nominal value of the supply voltage applied during the test.
- b. The field levels of fixed transmitters, such as radio telephone bases (cell/wireless) and terrestrial mobile radios, amateur radios and AM, FM and TV communication radios cannot be theoretically evaluated precisely. Site measurements are required in order to obtain the electromagnetic environment due to fixed RF transmitters. If the field level measured in the working environment of the P280 Overlay is greater than the above applicable levels of compliance, the proper operation of the P280 Overlay must be checked. If any anomalies are detected, additional measures must be taken, such as redirecting or relocating the reference equipment.
- c. Above the frequency band 150 kHz to 80 MHz, the field level must be less than 3 V/m.

NOTE:

These recommendations may not apply to certain situations. Electromagnetic propagation is affected by absorption and by the reflections of structures, objects and persons.


Recommended separation distances

Recommended separation distances between portable and mobile RF communications devices and the P280 Overlay			
The P280 Overlay is designed for use in an electromagnetic environment in which interference due to radiated RF is monitored. The customer or user of the P280 Overlay can help to prevent electromagnetic interference by keeping a given distance between the portable and mobile RF communications appliances (transmitters) and the P280 Overlay, as recommended below, and according to the maximum power of the communications appliances.			
Maximum power assigned to the transmitter W	Separation distance versus the frequency of the transmitter m		
	From 150 kHz to 80 MHz	From 80 MHz to 800 MHz	From 800 MHz to 2.5 GHz
	$d = 1, 16\sqrt{P}$	$d = 1, 16\sqrt{P}$	$d = 2, 33\sqrt{P}$
0.01	0.12 m	0.12 m	0.24 m
0.1	0.37 m	0.37 m	0.74 m
1	1.16 m	1.16 m	2.33 m
10	3.67 m	3.67 m	7.37 m
100	11.60 m	11.60 m	23.30 m
For transmitters whose maximum transmission power is not shown above, the recommended separation distance d in meters (m) can be estimated using the equation of application to the frequency of the transmitter, where P is the maximum output power of the transmitter in Watts (W), according to the manufacturer.			
NOTE: At 80 MHz and at 800 MHz, the separation distance for the highest frequency range applies. NOTE: <i>These recommendations may not apply to certain situations. Electromagnetic propagation is affected by absorption and by the reflections of structures, objects and persons.</i>			

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P005856	P005858
REF P005856 Quantity: 1 Transport and Storage Requirements Conditions de Transport et stockage T°: -18°C to +60°C R.H: 10 to 95% P°: 700 to 1060 mBar P280 Overlay 0459 220-240V ~ / 50 Hz 0.07 A SN L123AV12345 COUNTRY OF ORIGIN / MADE IN FRANCE Hill-Rom Industries S.A. 180, RUE DU CADUCEE 34100 MONTPELLIER CEDEX 5 FRANCE	REF P005858 Quantity: 1 Transport and Storage Requirements Conditions de Transport et stockage T°: -18°C to +60°C R.H: 10 to 95% P°: 700 to 1060 mBar P280 Overlay 0459 120V ~ / 50/60 Hz 0.17 A SN L123AV12345 COUNTRY OF ORIGIN / MADE IN FRANCE Hill-Rom Industries S.A. 180, RUE DU CADUCEE 34100 MONTPELLIER CEDEX 5 FRANCE

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

P005856	This is the P280 Overlay label for model P005856. It includes the following information: - Manufacturer: Hill-Rom Industries S.A., 188, Rue du Caducée - Parc Euromédecine, 34195 Montpellier Cedex 5 - France - Electrical specifications: 220-240V ~, 50 Hz, 0.07A; T1A 250V - Model: P280 Overlay - Reference: REF: P005856 - Modification: Modification - CE mark with 0459 - IP20 protection rating - Temperature range: +10°C to +40°C - Pressure range: 900 mbar to 1060 mbar - Humidity: 10% to 90%
P005858	This is the P280 Overlay label for model P005858. It includes the following information: - Manufacturer: Hill-Rom Industries S.A., 188, Rue du Caducée - Parc Euromédecine, 34195 Montpellier Cedex 5 - France - Electrical specifications: 120V ~, 50/60 Hz, 0.17A; T1A 250V - Model: P280 Overlay - Reference: REF: P005858 - Modification: Modification - CE mark with 0459 - IP20 protection rating - Temperature range: +10°C to +40°C - Pressure range: 900 mbar to 1060 mbar - Humidity: 10% to 90%

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
P005856


Technical drawing of a P280 Overlay for P005856. The drawing includes the Hill-Rom Industries S.A. logo and address: 188, Rue du Caducée - Parc Euromédecine, 34195 Montpellier Cedex 5 - France. The voltage is 220-240V ~, 50 Hz, 0,07A. The drawing shows a T1A 250V terminal block. The overlay is labeled P280 Overlay. The reference is REF: P005856. The modification is Modification: 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15. The drawing includes a CE mark with the number 0459. The drawing shows a temperature range of +10°C to +40°C. The drawing shows a pressure range of 900 mbar to 1060 mbar. The drawing shows a humidity range of 10% to 90%.

P005858

Technical drawing of a P280 Overlay for P005858. The drawing includes the Hill-Rom Industries S.A. logo and address: 188, Rue du Caducée - Parc Euromédecine, 34195 Montpellier Cedex 5 - France. The voltage is 120V ~, 50/60 Hz, 0,17A. The drawing shows a T1A 250V terminal block. The overlay is labeled P280 Overlay. The reference is REF: P005858. The modification is Modification: 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15. The drawing includes a CE mark with the number 0459. The drawing shows a temperature range of +10°C to +40°C. The drawing shows a pressure range of 900 mbar to 1060 mbar. The drawing shows a humidity range of 10% to 90%.

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P005856	P005858
REF P005856 Quantity: 1 (Q) Transport and Storage Requirements Conditions de Transport et stockage T°: -18°C to +60°C R.H: 10 to 95% P°: 700 to 1060 mBar P280 Overlay Hill-Rom CE 0459 220-240V ~ / 50 Hz 0.07 A SN L123AV12345 COUNTRY OF ORIGIN / MADE IN FRANCE Hill-Rom Industries S.A. 18, RUE DU CADISSE 94160 MONTREUIL-CEDEX 5 FRANCE	REF P005858 Quantity: 1 (Q) Transport and Storage Requirements Conditions de Transport et stockage T°: -18°C to +60°C R.H: 10 to 95% P°: 700 to 1060 mBar P280 Overlay Hill-Rom CE 0459 120V ~ / 50/60 Hz 0.17 A SN L123AV12345 COUNTRY OF ORIGIN / MADE IN FRANCE Hill-Rom Industries S.A. 18, RUE DU CADISSE 94160 MONTREUIL-CEDEX 5 FRANCE

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