- DK Brugervejledning
- EN User manual
- DE Gebrauchsanleitung
- FR Manuel
- ES Manual
- RU Руководство

# Congratulations on your new coloQuick

We are so pleased that you bought our product

If the product is to do you good we advise you to read this manual, and, most of all, that you let all the persons who are to use the product, read it.

We hope you will keep the manual and use it as an inspiration to get hold of your colostrum protocols.

## **How to get a quick start**

Find a suitable freezer – preferably an upright model, in order to have a proper overview of all your colostrum portions.

Start by filling the freezer with colostrum, without regard for the quality.

You are not supposed to use the coloQuick water bath, until the colostrum bank has been totally loaded.

Only when you have 6-8 portions in your freezer you are ready to start with the coloQuick system. From now on you will estimate how much colostrum you need to have in stock in order not to run out.

Before long your freezer will be filled up - exclusively with colostrum of the best quality- and in this way you are able to secure the most flying and optimal start for all your calves.

No matter how your colostrum procedures were until now we are sure that the coloQuick system will give your calves a big boost.

Once more – good luck with your calves in future.

Best regards

Coloquick International

Visit us on [www.coloquick.com](http://www.coloquick.com)

Congratulations on your new coloQuick	20
Service	22
When you receive a coloQuick	23
Firstly the calf – secondly the cow	24
General information	26
Type specifications	28
Troubleshooting – coloQuick thawing unit	29
Troubleshooting – Pasteur	30
Installation of coloQuick thawing unit	32
Filling station	33
Fact sheet	34
Standard procedures:	
Colostrum bank	110
Milking the cow	112
Checking the quality of the colostrum	114
Placing the bag in the cartridge	116
Making notes on the cartridge	118
When a calf is born	120
While the colostrum is heated	122
When the colostrum is ready	124
Handling of the cartridge	126
Feeding by teat	128
Feeding by esophageal tube feeder	130
Cleaning	132
Hygiene in general	134
Pasteurisation	136

*Please note that pictorial materials are only instructive. It will not necessarily reflect the precise specifications and details of the model in question.*

## Service

When contacting Calvex please have the following data available.  
You will find the information on the data plate of the machine.

Type: \_\_\_\_\_

Year of production: \_\_\_\_\_

Machine serial number: \_\_\_\_\_

*Calvex A/S, Holstebrovej 104, 7800 Skive, Danmark · phone +45 9753 7333*

In this manual the following symbols were used:

*Caution – Important  
security instructions.*

*General information  
and advice.*

*Environmental information.*

*We reserve the right to amendments without further notice.*

## When you receive a coloQuick

### On supply

Supply, everything is in one box.  
Unwrap the box.

Avoid damaging the materials when unwrapping.  
Avoid to cut open the cardboard wrapping with a knife or the like.

On supply there is equipment inside the water bath.  
Exercise necessary caution so the equipment and the water bath will not be damaged.

Check that you received all the parts for the coloQuick-concept (see the below parts list).

### Parts list Thawing unit

- 1 Thawing unit
- 1 Filling station
- 10 Cartridges
- 35 Colostrum bags
- 1 Colostrometer
- 1 Set of teats
- 1 Esophageal tube feeder
- 2 Belts
- Hoses, fittings, etc.


### Parts list Pasteur

- 1 Pasteur
- 1 Filling station
- 10 Cartridges
- 35 Colostrum bags
- 1 Colostrometer
- 1 Set of teats
- 1 Esophageal tube feeder
- 2 Belts
- Hoses, fittings, etc.


# Firstly the calf – secondly the cow

It is much more important to feed the calf with colostrum quickly after calving than it is to milk the cow. Therefore, feeding of the calf has to be separated from milking of the cow. Only in that way you can secure a perfect start for all calves.

## coloQuick thawing unit

# Pasteur

ENGLISH

# General information

## Precautions

Please read this manual carefully before installation and putting into operation:

- for the safety of you and your surroundings
- for the correct operation and use of the machine.

Always keep the manual together with or nearby the machine and let it come with, if you move or give it to another person.

The producer is not responsible, if the machine is installed incorrectly, and the use leads to damage.

## Precautions for kids and exposed persons

- Please be aware not to let any persons, including children, with impaired physical senses, impaired mental abilities or lack of experience and knowledge operate the machine. They must be supervised by a person, who is responsible for their security, or they must be instructed how to use the machine.
- Don't let children play with the machine.
- Please store all packaging materials unreachable to children. There is danger of suffocation or personal injury.
- Please don't let children and pets come near the machine, while the lid is open.

## Precautions in general

- The specifications of the machine must not be changed. This could cause risk of personal injury and damage of the machine.
- Please don't use inflammable products or components wetted with inflammable products in the machine or close to or on it. There is danger of explosion or fire.

## Maintenance and cleaning

- Switch off the machine, pull out the plug.
- Don't use the machine without the appurtenant filters. Make sure they are mounted correctly. If not it may cause water leaks.
- Don't splash water directly on the machine.

## CAUTION

Please do not use this machine, if the lid to the electrical installations is open.

## Installation

- The machine is heavy and difficult to handle. Please be careful when moving it.
- The best way of transportation is on a pallet. The size of it fits on one half euro pallet.
- Don't install and connect a damaged machine.
- Make sure to remove all the packaging before start-up.
- Make sure the plug is not connected during installation.
- Electrical and plumber works and installation must be done only by skilled persons. Then personal or materials damage will be reduced.
- Don't place or use the machine in a place where the temperature falls below 0°C or rises above 30°C.
- Install the machine on a firm surface – place it horizontally.

## Connection, water

- Please don't connect the machine with old hoses, which have already been used. Use only new hoses.
- Be careful not to damage the hoses.
- Don't connect the machine to new pipes or pipes that have not been used for a long time. Let the water run for some minutes and then connect the inlet hose.
- Make sure the hoses and the couplings are tight before the machine is put into operation.

## Connection, power

- Please check that the machine has an earth connection.
- Make sure the electrical data on the data plate correspond to the power supply.
- Always use a correct and legally installed electricity network.
- Don't use connectors and extension cords. There is a risk of fire.
- Don't replace the electrical cable. Contact the service centre.
- Be careful not to damage the mains plugs and the cable.
- Don't put the mains plug in the socket until the installation has been completed. Make sure there is an access to the plug after the installation.
- Don't switch off the machine by pulling the cable. Always take hold on the mains plug itself.

### **Right to exchange – damaged bags**

Please note: to be sure of a quick solution you must save damaged bags for inspection.

Have you any questions as to the return, please contact Calvex on phone No. +45 97537333 or mail to: info@coloquick.com.

### **Disposal of the machine**

1. Pull out the plug
2. Snip off the power cable and discard it
3. Remove the lid lock and discard it. This will prevent children and animals from locking themselves up in the machine.

### **Danish environmental considerations / packaging materials**

Place the packaging material in suitable containers for recycling.

***All parts may be recycled!***

## Type specifications

<b>Model:</b>	coloQuick thawing unit	Pasteur
<b>Version:</b>	IV	III
<b>Type:</b>	Colostrum heater	Pasteuriser and colostrum heater
<b>Connection:</b>	230V/400V 3-phase (EU) 208V/240V 3-phase (US)	230V/400V 3-phase (EU) 208V/240V 3-phase (US)
<b>Power:</b>	6100W	6200W
<b>Water contents (litres):</b>	85	85
<b>Consumption (kWh):</b>	- 1.8 (thawing)	4.8 (pasteurisation) 1.8 (thawing)
<b>Dimensions, LxWxH:</b>	686 x 515 x 732 mm	708 x 525 x 735 mm
<b>Green cartridge Dimensions, LxWxH:</b>	417 x 378 x 62 mm	417 x 378 x 62 mm
<b>Weight water bath (kos):</b>	60	68
<b>Weight filling station (kos):</b>	8,9	8,9

## coloQuick thawing unit

Problem	What is happening	Solution	Comments
Motor runs, but the cartridge in the water tank is not rotating.	Split-pin has broken.	Mount a new bolt. Check that the motor and the tank shafts are end to end.	 *
Motor runs, but temperature is not rising.	<p>No voltage to heating element.</p> <p>Level sensor does not function correctly or is defective. <i>(NB: Level sensor has not been mounted in all models).</i></p> <p>Voltage to heating element.</p>	<p>Affect the button of the thermal fuse.</p> <p>Affect the movable part of the sensor. Replace the sensor, if necessary.</p> <p>Heating element defective – must be replaced.</p>	 
Motor is not running.	<p>Lid switch not activated.</p> <p>Fuse defective.</p>	<p>Check plate of the lid switch. Maybe it is skew. Lid switch defective – must be replaced.</p> <p>Replace fuse (max. 3A inert).</p>	 
Motor makes noise or vibrates slightly.	If the motor is worn it will make more noise and maybe vibrate slightly.	Replace motor.	
Thermostat reads incorrectly.	Very large temperature fluctuations.	Thermostat must be replaced. <i>(NB: Thermostat has not been mounted in all models).</i>	
Temperature does not rise, or heating takes too long.	Split-pin has broken.	Replace split-pin.	* Same picture as on top.

ENGLISH

## Pasteur

Problem	What is happening	Solution	Comments
<p>Display shows error code E01:</p> <p><i>Motor error</i></p>	<p>No feedback from motor for 20 seconds.</p>	<p>Check connections to motor. Cables have fallen off. Motor lacks voltage from power supply.</p>	
<p>Display shows error code E02:</p> <p><i>Water inlet – error</i></p>	<p>Appears if the water inlet valve has been open for more than 20 minutes.</p>	<p>Open the water inlet. Inlet choked.</p>	
<p>Display shows error code E05:</p> <p><i>Heat element error</i></p>	<p>Appears if the temperature for a period of more than 5 minutes will not rise by more than 0.2 degrees.</p>	<p>Heat element was not connected. Thermal fuse switched off.</p>	
<p>Display shows error code E08:</p> <p><i>Water level sensor error</i></p>	<p>Appears if there is an “illegal” feedback from the water level sensors. For instance if L1 (lower sensor) does not warn that it is covered - contemporary with the second sensor from the bottom breaking off. From the error parameter it appears which sensors were active when the error occurred, and the sensor sensitivity.</p>	<p>Dirt must be removed from the sensors. Check connection to sensors.</p>	
<p>Display shows error code E09:</p>	<p>Appears if the water level stands below L2 (second from the bottom) when running a program.</p>	<p>Dirt must be removed from discharge and discharge valve.</p>	

## Pasteur

Comments	What is happening	Solution	Comments
Motor runs, but cartridges in the water tank do not rotate. They are not rotating, if you cannot hear the water sloshing around.	Bolt has broken.	Mount a new bolt. Check that the motor and the tank shafts are end to end.	
Motor runs, but temperature is not rising.	No voltage to heating element.	Affect the button of the thermal fuse.	
Motor is not running.	Lid switch not activated. Lock light flashes.  Lid switch has been activated. Lock light flashes.  Fuse defective.	Check plate of the lid switch. Maybe it is skew.  Lid switch defective. Replace it.  Replace fuse (max. 3A inert).	 
Motor makes noise or vibrates slightly.	If the motor is worn it will make more noise and maybe vibrate slightly.	Replace motor.	

## Installation of thawing unit

Place the thawing unit on a firm, plane surface and in a place that fits into the work routines.

Attach handles to the lid of the thawing unit.

Fill water into the water bath up to the mark.

Be sure there are no foreign bodies in the water tank.

Check that the cartridges are properly fixed into the cartridge holder.

Get an authorised electrician to install the power supply to the water bath correctly as prescribed.

Never turn on the power with no water in the water bath!

Please note that the guarantee will not cover incorrect installation of power!

In case of 3-phase alternating current the phase sequence must be correct (L1-L2-L3).


The water tank must **always** be switched on – and have the same temperature all the time.

Change of temperature, if necessary:

Shortly press ▽. Then set the wanted temperature by means of the ▽ and △ buttons. After 20 seconds the changes are automatically stored. Do not use the U- and P-keys.

*When using the system for the first time it is important to manually check by at thermometer how much the colostrum temperature falls, from it is taken up from the water bath – and until it is fed. This must be repeated if the feeding procedure is changed. Please note that the large surface of the bag also means that the milk will be cooled quickly!*

ENGLISH

## Remember the 3 Qs

It is all important for the newborn calf's success as a calf, heifer / fat calf and milk cow that it gets sufficient colostrum of good quality quickly after it was born.

We are talking about the 3 Qs – Quick, Quantity and Quality:

**Quick** – is for colostrum quickly after birth. The coloQuick concept prescribes that the calf must have colostrum within 1 hour after birth. With the coloQuick this has now been feasible.

Time is the most important factor, because that is the one you can directly influence the most when using coloQuick.

**Quantity** – the quantity of colostrum. A normal calf of 40 kos must have 4 litres of colostrum (10 % of the calf's weight).

If you feed colostrum by teat the calf can normally not drink 4 litres, and then it is even more important that the calf has colostrum of first quality very quickly after birth.

**Quality** – Good quality means at least 50 grams of antibodies per litre of colostrum. But there are also many other valuable substances in the colostrum.

Quality also means top hygienic colostrum. Colostrum from mastitis cows and colostrum with visible changes must be discarded.

On the whole think in that way that the calf should not drink any milk you would not drink yourself.

## Good "immunization" of the calf

It is most important that the newborn calf has as many antibodies transferred to its blood as altogether possible.


If there is at least 15 grams of antibodies per litre of blood on the calf's 2nd day of life then we can talk about a "good immunization". The Americans talk about "Failure of Passive Transfer" as an expression of a bad immunization.

This figure can be measured or calculated. Please see the below table.

With the coloQuick concept you make sure that the calf reached a "good immunization" with more than 22 grams of antibodies per litre of blood.

By comparison a calf will obtain 3 grams of antibodies per litre of blood, if it has only 2 litres of colostrum with 40 grams of antibodies per litre 6 hours after birth. This calf will have a catastrophic start of life.

Quantity & Quality IgG

- Good immunization.
- Acceptable.
- Bad immunization.

## Quality – colostrometer

The contents of antibodies in the colostrum may vary a lot – from less than 10 to more than 200 grams of antibodies per litre. The average is typically approx. 40-50 grams of antibodies per litre. We define a good quality as more than 50 grams per litre.

Unfortunately, it is a widespread misunderstanding that you can estimate the contents of antibodies in the colostrum without measuring it. That is not possible!

A colostrometer will measure the contents of antibodies in the colostrum. The colostrometer is a hydrometer which is lowered into the colostrum.

If it floats at a high level so the surface of the colostrum intersects the colostrometer in the dark green area, the quality is good. If the colostrometer sinks deep down so the surface of the colostrum intersects in the red area, the quality is poor. The light green area is in the middle between them.

The colostrometer shows the quality on a colour scale and states a density – as an example a density of 1.045 corresponds to 50 grams of antibodies per litre of colostrum.

The colostrometer is estimated for a colostrum temperature of 20° C. If the temperature is higher the colostrometer will underestimate the quality – and the opposite at lower temperatures.

If the colostrum foams there is also an indication error.

## Quality – diseases


Make sure to minimize the risk of using colostrum with paratuberculosis and salmonella bacteria.

Never use colostrum from risk cows.


## Pasteurisation of colostrum

In the old days people believed you could not pasteurise the colostrum without destroying the antibodies. But modern research confirms that you really can.


The formula is heating of the colostrum to 60°C for 60 minutes. The secret is a very precise control of the temperature.

If the temperature is too high the wanted killing of bacteria will still go on, but the antibodies will also be destroyed.

If the temperature is not high enough, the antibodies are preserved, but the wanted killing of bacteria will not take place.

## Colostrum bank

**Freezer** – the coloQuick concept prefers a freezer rather than a fridge for the colostrum bank.

- It clearly gives the simplest “flow” through the colostrum bank
- The colostrum keeps in the freezer for one year
- It only takes a few minutes to heat the colostrum from frozen condition to drinking temperature
- Keeping the colostrum in the freezer ensures the quickest possible cooling.
- It is easier to have an overview of the colostrum cartridges in an upright freezer than it is in a chest freezer

**Fridge** – if you still prefer to keep your colostrum bank in a fridge, then please be aware:

- that the temperature must be max. 4-5°C
- that the storage time is max. 48 hours
- that “the flow” will not end in a mess
- that fresh colostrum which is placed in the fridge will cool down slowly, and bacteria will grow quickly, as long as the temperature is more than 10-15°C. Even down to 4-5°C there will still be a certain growth of bacteria.

## The colostrum temperature

The right temperature for milk for feeding is the body temperature of the calf, which means that the temperature of the colostrum must be approx. 39°C.

The colostrum temperature falls relatively quickly, when it is taken up from the water bath. This is due to the large surface of the colostrum bag.

If you feed the colostrum by teat then please be aware of this temperature drop.

When feeding by esophageal tube feeder it is important that the colostrum is not too hot. It is unpleasant to the calf, and it might be highly responsive to it.

## Feeding by esophageal tube feeder

*Also see page 130.*

Treat the calf gently and calmly and learn to feed by esophageal tube feeder with good skills. Remember, the calf is a newborn baby! If you feel uneasy then ask your general vet to instruct you in the feeding by esophageal tube feeder.

- Hygiene is most important in connection with all handling of the colostrum. Always wear disposable gloves when feeding the calf, or wash your hands carefully.
- Make sure the esophageal tube feeder has a proper and pleasant temperature. This is especially important, if the esophageal tube feeder is kept in the freezer while not used.
- Make sure the colostrum is not too hot. It is unpleasant to the calf, and it might be highly responsive to it.
- Place the colostrum cartridge at a convenient height – on your back, over your shoulders or in a cartridge holder on intermediate bars or the like. The higher it is placed the faster the colostrum bag will be emptied.
- “Grease” the tip of the esophageal tube feeder before use – for instance by dipping it in colostrum, water or something like that.
- Place the calf in a corner of the box – or the like – so you don’t have to put efforts into holding the calf. Avoid “fighting”. If the calf will not stand up then, if possible, make sure it lies “upright” on its chest.
- Stand astride with the head of the calf between your legs.
- Press your thumb into the corner of the calf’s mouth behind its teeth. Now you have a good grip on the head, and it is easy to insert the esophageal tube feeder.
- Move the tip of the esophageal tube feeder to the back of the tongue of the calf, until the sucking reflex is stimulated. Then gently and with feeling push the esophageal tube feeder down into the gullet of the calf. It is important to follow the calf’s signals and not to rush.
- Feel on the throat of the calf when the probe is led down through the gullet – especially the drop-shaped tip is clearly felt. If the probe is in the windpipe you cannot feel it on the calf’s throat.
- Once the esophageal tube feeder is down and in place, you open for the colostrum – not until then!
- Adapt the pace to the signals of the calf. Turn off the milk supply, if the calf is choking on the milk – that is, if some colostrum comes up again through the gutter and could end in the windpipe.
- When the calf got its colostrum then turn off the probe, and easily remove it.

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## Colostrum bank

Build up a colostrum bank of good and quality controlled colostrum.

We recommend that the colostrum is kept in an upright freezer, where the green cartridges can be placed, easy to survey and easy to take out.

Start by filling the freezer with colostrum, without regard for the quality.

### Before putting into operation:

You are not supposed to use the coloQuick water bath until you have 6-8 portions in your colostrum bank.

- Milk the cow
- Give the calf its colostrum
- Check the quality of the excess colostrum with a colostrometer
- All colostrum is filled into bags – without regard for the quality
- Make a note of the quality on the writing boards on the cartridges
- Place the cartridge in the freezer.

### Putting into operation

Only when you have 6-8 portions in your freezer you are ready to start with the coloQuick system.

In future you will estimate how much colostrum you need to have in stock in order not to run out.

Place the colostrum bank in a suitable place that fits into the everyday work routines.

The colostrum may be kept in the freezer for up to 1 year.

DK

EN

## Malke koen

Koen malkes, når det passer ind i bedriftens rutiner.

## Milking the cow

The cow is milked when it fits into the other routines of the farm.

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The hygiene must be first-class. The cow has not been milked for several weeks, and therefore udder and teats must be most carefully cleaned.

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The premilking must be extra careful.

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Milking equipment, etc. must be handled with at least the same cleaning and hygiene procedures as apply to the milking of dairy milk.

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The freshly milked colostrum must not be exposed to manure pollution.

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The colostrum must be cooled as soon as possible.

Brug altid engangshandsker.

Always remember to use disposable gloves.

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## Checking the quality of the colostrum

Cool the colostrum to 20°C.

Pour the colostrum in the  
hopper of the filling station.

Please note the litre scale.

Use only top hygienic milk  
without blood and other visible  
contamination.

Discard mastitis milk.

Check the quality of the  
colostrum with a colostrumeter.

Gently lower the colostrumeter  
into the colostrum in the  
hopper of the filling station.

The scale to be read at the  
surface of the milk.  
The deeper the colostrumeter  
falls – the poorer quality.

Good quality is poured into the  
colostrum bag.

Squeeze the cartridge to get  
the air out of the bag, before  
the lid is screwed on.

Poor quality is discarded, or it is  
used for second-time feeding.

*Important – all colostrum must  
be checked by a colostrumeter.  
You cannot see or feel, whether  
the quality is okay.*

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Pres de to perforeringer i  
posens bund ned over de to  
"pigge".

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## Placing the bag in the cartridge

The cartridge must be half-  
open (90°).

Place it with the two "spikes"  
to the left.

Press the two perforations at  
the bottom of the bag down  
onto the two "spikes".

With your right hand fix  
the connecting pipe of the  
bag onto the groove of the  
cartridge.

Hold the connecting pipe in  
place with your right hand and  
close the cartridge.

Lock the cartridge with the red  
"locking buttons".

Check that the bag is tight.

*Be sure always to have a number  
of cartridges with colostrum bags  
ready for use.*

DK

EN

## Notater på kassetter

Skriv notater på skrivepladen

- fx dato, kvalitet, konummer m.v.
- med whiteboard-pen
- eller med blyant, som kan slettes med viskelæder.

## Making notes on the cartridges

Write your notes on the writing board

- e.g. date, quality, cow number, etc.
- with a white board pen
- or with a pencil that can be erased.

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Fix the writing board onto the cartridge.

There is room for more than one writing board.

Kassetten placeres i råmælksbanken.

The cartridge is placed in the colostrum bank.

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## When the calf is born

Immediately take a portion of colostrum from the colostrum bank.

Put the cartridge with the disposable bag into the water bath. Make sure the cartridge is properly fixed.

The water must cover the colostrum bag, once the cartridge is in upright position.

Turn on the motor in the water bath and wait for 20 minutes.


Tip: Set the alarm, if any, on your mobile.

Please note that colostrum is a natural product – and therefore thawing time may vary. In most cases the colostrum will be ready in 10-15 minutes, but in some cases it may take up to 20 minutes.

If you add two cartridges at a time, you must increase time to 25 minutes.

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## While the colostrum is heated

As soon as possible separate the calf from its mother.


Remove mucous matter from the mouth and throat.

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
## When the colostrum is ready

Take the cartridge out of the water bath. Place it in upright position with the lid of the bag up.

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Unscrew the lid.

Spot check the temperature now and then just before feeding.

Remember that due to the large surface of the bag the colostrum will cool quickly.

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Connect the tubing – for teat or esophageal tube feeder.

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## Handling of the cartridge

Wear the cartridge over your shoulder by a belt.

Wear the cartridge as a rucksack by two belts.

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Hang up the cartridge on intermediate bars on a holder.

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The above are just examples – may be handled in other ways, too.

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## Feeding by teat

The calf's need to suck and its ability to utilise the colostrum is most distinct immediately after birth. So therefore: feed as soon as possible after birth.

Make sure the colostrum has the right temperature – and most important, not too hot.

Do not stress the calf.

Avoid any disturbance or interruption, once the calf is sucking.

Hold the teat so the neck of the calf is stretched out.

See to it that the colostrum is not getting too cold.

If necessary, repeat warming up.

Let the calf drink as much as possible.

A normal calf of 40 kos will normally drink approx. 2 ½ litres of colostrum.

If the calf has less than approx. 4 litres of colostrum, it must be fed with colostrum with even higher contents of antibodies.

**i**

*Teat, etc. must be cleaned up thoroughly – and scalded after use. Always wear disposable gloves.*

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## Feeding by esophageal tube feeder

- The calf must always be standing up.
- Place the esophageal tube feeder above the tongue, until the sucking reflex is stimulated.
- Slowly push the esophageal tube feeder down into the gullet.

You must be able to feel the tip of the probe at the left side of the throat.

**i**

Only when the esophageal tube feeder is placed correctly you may open to the milk.

It is important that no milk is spilt, which may end in the windpipe.


Be sure that the colostrum holds the right temperature – and important that it is not too hot.

Turn off the tap of the probe, before it is removed.

**i**

Remember always to wear disposable gloves.

Please also see the section "Fact sheet" previously in this manual.

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## Cleaning

### Filling station

- Clean the filling station after use every time.
- Start by rinsing with cold water.
- Then wash with hot water (min. 50°C) and detergent.
- Rinse with clean water.

### Thawing unit


- Change the water in the water bath – and rinse with water.
- Then wash with hot water (min 50°C) with detergent. Sweep away cleaning water.

### Probe / teat and tubing

- Be sure to thoroughly rinse probe/teat and tubing immediately after use.
- Probe and tubing are rinsed thoroughly with cold water.
- Clean the outside of probe/teat and tubing with a sponge or the like.
- Then scald probe/teat and tubing thoroughly with hot water – preferably 80°C.
- Empty probe and tubing.
- Keep probe/teat and tubing in the freezer to prevent bacteria growth.

Be sure to take the chill off probe/teat and tubing before use next time.

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## Hygiene in general

- Always wear disposable gloves
- when handling the colostrum
  - when milking etc. – and last, but not least
  - when feeding the calf.

The use of disposable gloves is very hygienic and helps you "think hygienically".

See to it that the colostrum is quickly cooled down and frozen.

The colostrum bags **cannot** be reused.

It is not possible to clean the bags for reuse. If they are scalded the bacterium barrier of the bags will be destroyed.

### Colostrum journal:

Instead use the bags as a colostrum journal or use the forms in the box with the colostrum bags.

**i**

Be very careful to avoid that the calf ingests any bacteria, which being in the guts will cause sick calves and a bad utilisation of the antibodies of the colostrum.

If for example the calf gets manure in its mouth the size of a pea, before it gets any colostrum – the utilisation of the colostrum is practically nil!!

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Once the colostrum has been checked for good quality – and has been poured into bags – place it in the Pasteur.

Immediately after the pasteurisation – quickly cool down the colostrum – e.g. in a bucket of cold water.

Dry the cartridge thoroughly – and add it to the deep-frozen colostrum bank.

*Pasteurisation is a heat treatment in order to eliminate as many harmful bacteria as possible – without damaging the important antibodies.*

*The colostrum is pasteurised at 60°C for 60 minutes.*