

User Manual

P03123 Rev 2.11.0_Au/Eu



www.signostics.com.au

CONTENTS

1. Introduction	1
About Signos RT	1
Patient Environment	2
Classifications	2
About the User Manual	2
Intended uses of Signos RT	4
Transducer applications	5
Training	7
Contact information	8
Trademarks and document copyright	8
Labeling	9
2. Getting Started	11
Signos RT components	11
Display unit components	11
Display unit controls and functions	12
Turning Signos RT on	12
The microSD card	13
Recovering security key information	13
Log in screen	13
Urgent Mode	14
How to insert the transducer	14
Connecting and disconnecting the Power Supply	14
3. Scanning	15
Setting Up	15
Default Settings	15
Entering exam and patient information	15
Imaging modes	15
Changing imaging modes	15
Scan button	15
B-mode Imaging	16
M-mode Imaging	16
PW Doppler Imaging	16
Scroll wheel adjustments	17
Recommended ultrasound transmission gel	17
Recommended ultrasound probe sterile sheaths	17
4. Signos RT screens and menus	18
Home screen	18
Ultrasound screen	18
Navigating the touch screen	19
Menus	19
Message boxes	19
Screens and windows	19
Tool Bar icons	20
The Title bar	21
Application icons	21
Screen title	21
Icon tray	21
System Menu	22
Home screen	22
Ultrasound	22
Printer	22
Preferences	23
Flashlight	27
Status icons	28

Log in/Log out/Remove Card	30
Sleep	30
Shutdown.....	30
About	31
Turning Signos RT off.....	31
5. Controls for the Ultrasound application	32
Display unit controls	32
Touch screen.....	32
Back button.....	32
Power button	33
Scan button	33
Scroll wheel	33
The Tool Bar	33
Main Menu	34
Exam management	34
Exam management screen.....	35
Starting a new exam	35
Editing patient data.....	36
Reviewing exams.....	36
Closing exams.....	36
Image/Exam Scroller	37
Movie Clip Controls	37
Saving a Movie	37
Scan Review Mode.....	37
Measurement.....	38
Using measurements.....	38
M-mode and PW Doppler image measurement.....	40
Fetal Biometry Measurements.....	40
Annotations.....	40
Text annotations	40
Voice recordings.....	41
Delete Image	41
Hide Tool Bar.....	41
6. The Main Menu	42
Scan settings	42
Changing and saving scan settings	43
Customizing scan presets.....	44
Loading a preset.....	44
System-defined presets.....	45
Sequences.....	46
Bladder Auto sequence	46
Bladder Calipers sequence.....	47
Grid and Rulers.....	48
Image.....	48
Print Images... ..	49
Print Screen	49
Preferences	49
Other Main Menu options	49
About.....	49
Home	49
7. Scanning shortcuts.....	50
Pop-up menus	50
Scan button.....	50
Scroll wheel	50
Touch screen.....	51
Panning and zooming	51
Keyboard gestures.....	51

8. Maintenance and Troubleshooting	52
Maintenance	52
Battery Maintenance.....	52
Transducer cleaning and disinfection.....	52
Cleaning	53
Disinfection	53
Troubleshooting.....	55
Signostics Customer Support	55
9. Connectivity and accessories	56
Accessories	56
System connectivity.....	56
Power sources.....	56
Internal battery.....	57
Mains charger	57
SigViewer software.....	57
10. Specifications.....	58
System dimensions	58
Temperature and humidity limits	58
Mains charger	58
Internal battery.....	59
MicroSD card	59
Measurement accuracy	59
Degree of Protection Against Ingress of Water	60
11. Safety	61
Ergonomics.....	61
External materials	61
Disposal	61
Electrical safety.....	61
Electromagnetic compatibility (EMC).....	62
Mains charger	64
Biological safety.....	64
Acoustic Output	64
References	66
Glossary.....	67
Key to icons	68
Appendix 1: Definitions and parameters of scan settings	69
Appendix 2: System-defined preset parameters.....	71
Index	73

1. INTRODUCTION

This **Signos RT** User Manual is applicable to customers in the Australia, Europe, and countries that recognize CE marking only. Customers in countries in the United States and Canada must refer to **Signos RT** User Manual P03622 or contact [Signostics](#) for further assistance (page 8).

About Signos RT

Signos RT Personal Ultrasound has been designed to deliver affordable and conventional handheld ultrasound diagnosis to the point-of-care. The **Signos RT** generates and transmits ultrasound energy in the form of pulses in the 3-5MHz range into a patient, detects the reflected pulses, and processes the information in order to generate ultrasound images.

The system is available with 3MHz sector and 3-5MHz sector transducers. Both sector transducers allow penetration in abdominal and cardiac applications, with the 3-5MHz sector transducer having better superficial resolution. **Signos RT** allows for image storage, review and measurement.

Signos RT operates in three imaging modes, B-mode (Fig 1.1), M-mode (Fig 1.2), and PW (or Pulsed Wave) Doppler mode. These imaging modes are further discussed in [Scan](#) (page 15).




Fig 1.1 B-mode image of liver and right kidney



Fig 1.2 M-mode image of lung slide




Fig 1.3 Pulsed Wave (PW) Doppler Image of arterial blood flow

Patient Environment

The **Signos RT** user is expected to bring the device into the patient environment for use. The Signostics power supply (P03198), any interconnected computers, and any other peripherals must be kept outside the patient environment. See the diagram below for a drawing of the patient environment. When the power supply (P03198) is connected to the device, ensure the connecting micro-USB lead can be easily disconnected.




Fig 1.4 Patient Environment

Warning

- ✓ Equipment not suitable for use in the presence of a **FLAMMABLE ANAESTHETIC MIXTURE WITH AIR OR WITH OXYGEN OR NITROUS OXIDE**.
- ✓ Do not use the system near high-frequency surgical equipment, as it could create a burn hazard.

Classifications

The **Signos RT** is internally battery powered during scanning.

The Signostics Power Supply P03198 classification for **Protection against electric shock: Class II equipment**.

The 3MHz Sector (P03479) and 3-5MHz Sector (P03611) transducers are: **Type BF Applied Parts**.

The **Signos RT** probe is classified as **IPX1**.

The **Signos RT** is not for use within an oxygen rich environment.

The **Signos RT** is suitable for continuous operation (see section 10).

About the User Manual

This User Manual is a reference tool for users of the **Signos RT** Personal Ultrasound; it does not constitute medical advice nor provide clinical training, instruction in exam protocols, or information on how to interpret scans.

This manual should be read before **Signos RT** is used in the clinical setting.

Warning and Caution conventions

A Warning describes precautions a user should take to prevent injury or loss of life.

A Caution describes precautions a user should take to prevent damage to the device.

Manual conventions

The following style conventions are used in this manual:

- Buttons found on your **Signos RT** are indicated in bold italics, e.g., ***Back button***, ***Scroll wheel***. This style is also used to describe areas of the touch screen, e.g., ***Image name***.
- ‘Tap’ refers to touching the screen quickly with the stylus.
- ‘Click’ refers to pressing and releasing a button on the display unit, such as the ***Back button*** or the ***Scroll wheel***.
- ‘Scroll’ refers to rotating the ***Scroll wheel***.
- Text on the touch screen, such as screen names or menu options, are indicated in bold, e.g., **Home** screen. Icon names also appear in bold, e.g., **Main Menu**.
- New terms that describe functions for the ultrasound application are introduced in italics, e.g., *exam*.
- Links to other sections within the manual appear colored and underlined, e.g., see [Contact information](#) (page 8).

Version information

This User Manual applies to **Signos RT** with software version V2.11.x

Select **Main Menu**>[About](#) (page 49) for details on what software version is installed on your **Signos RT**.

Contents of the Signos RT pack

Your **Signos RT** pack contains a **Signos RT** unit (consisting of an ultrasound probe, one or more ultrasound transducers, and a display unit with connecting cable) and the following accessories:

- A Quick Reference Guide.
- A microSD card.
- A mains charger.
- A USB 2.0 compatible connection cable.
- A DVD containing:
 - The User Manual.
 - **SigViewer** software and User Manual.
 - Training materials including:
 - A generic ALARA education program (see enclosed ISBN 1-93004 7-71-1, *Medical Ultrasound Safety*).
 - Terms and conditions of Warranty.
- A sample of Ultrasound Transmission Gel.



Fig 1.5 Signos RT Pack Contents

Intended uses of Signos RT

The **Signos RT** ultrasound system is for non-invasive imaging of the human body and is to be used as an adjunct to physical examination for abdominal, fetal, cardiac, peripheral vessel and musculoskeletal applications. Users must have ultrasound training before using the device.

Warning: Where **Signos RT** imaging findings are inconclusive, the patient should be managed by further appropriate clinical investigation.

Warning: The **Signos RT** may not be effective in imaging deep structures and therefore may not be effective in imaging obese patients.

This system can be used to obtain 2D, Pulsed Wave (PW) Doppler, or M-mode images for:

Abdominal applications

- ✓ Ruling-in the presence and extent of aortic aneurysms;
- ✓ Ruling-in the presence and extent of free peritoneal fluid;
- ✓ Ruling-in the presence and extent of gross gall bladder pathology; and
- ✓ Assessing the size of the urinary bladder.

Warning: Studies have shown the effectiveness of ultrasound in assessing free fluid is reduced for penetrating trauma.

Fetal applications

- ✓ Determining gross fetal presentation;
- ✓ Ruling in the presence of fetal heart beat; and
- ✓ Estimating gestational age.

Warning: As a rule-in imaging device, it is NOT intended that **Signos RT** will be used for assessment of fetal morphology or fetal abnormalities.

Cardiac applications

- ✓ Ruling-in the presence and location of fluid in the pericardial space;
- ✓ Ruling-in the presence and location of fluid in the pleural space; and
- ✓ Ruling-in the presence and extent of pneumothoraces.

Warning: As a rule-in imaging device, it is NOT intended that **Signos RT** will be used to provide detailed assessment of cardiac anatomy.

Musculoskeletal applications

- ✓ Observing muscle movement for use as a biofeedback tool.

Warning: As a rule-in imaging device, it is NOT intended that **Signos RT** will be used to provide detailed assessment of musculoskeletal structures.

Peripheral Vessel

- ✓ Locating and identifying large blood vessels.

Contraindications

- ✓ The **Signos RT** is designed for percutaneous scanning only. Do not attempt intracavity imaging; in particular, trans-esophageal, trans-vaginal and trans-rectal scans are contraindicated.

Transducer applications

The intended use applications depend on the attached transducer (see Tables 1.1 and 1.2). These applications are further defined in [Intended uses](#) of **Signos RT** (page 4).

Table 1.1 Signos RT applications for the 3MHz transducer (P03479)

Clinical Application		Mode of Operation		
General	Specific	B-mode	B-mode/ M-mode	B-mode/ PW Doppler
Ophthalmic	Ophthalmic			
Fetal imaging & other	Fetal	✓	✓	✓
	Abdominal	✓	✓	✓
	Intra-operative (specify)			
	Intra-operative (neuro)			
	Laparoscopic			
	Pediatric	✓	✓	✓
	Small organ (specify)			
	Neonatal cephalic			
	Adult cephalic			
	Trans-rectal			
	Trans-vaginal			
	Trans-urethral			
	Trans-esophageal (non-cardiac)			
	Musculoskeletal (conventional)	✓	✓	✓
	Musculoskeletal (superficial)			
	Intravascular			
	Other (specify)			
Cardiac	Cardiac adult	✓	✓	✓
	Cardiac pediatric	✓	✓	✓
	Intravascular (cardiac)			
	Trans-esophageal (cardiac)			
	Intra-cardiac			
	Other (specify)			
Peripheral vessel	Peripheral vessel	✓	✓	✓
	Other (specify)			

Table 1.2 Signos RT applications for the 3-5MHz transducer (P03611)

Clinical Application		Mode of Operation		
General	Specific	B-mode	B-mode/ M-mode	B-mode/ PW Doppler
Ophthalmic	Ophthalmic			
Fetal imaging & other	Fetal	✓	✓	✓
	Abdominal	✓	✓	✓
	Intra-operative (specify)			
	Intra-operative (neuro)			
	Laparoscopic			
	Pediatric	✓	✓	✓
	Small organ (specify)			
	Neonatal cephalic			
	Adult cephalic			
	Trans-rectal			
	Trans-vaginal			
	Trans-urethral			
	Trans-esophageal (non-cardiac)			
	Musculoskeletal (conventional)	✓	✓	✓
	Musculoskeletal (superficial)			
	Intravascular			
	Other (specify)			
Cardiac	Cardiac adult	✓	✓	✓
	Cardiac pediatric	✓	✓	✓
	Intravascular (cardiac)			
	Trans-esophageal (cardiac)			
	Intra-cardiac			
	Other (specify)			
Peripheral vessel	Peripheral vessel	✓	✓	✓
	Other (specify)			

Training

Warning

- ✓ You must have appropriate ultrasound training before using the system.

The Signos RT device is intended to be used by clinicians with appropriate professional qualifications and clinical training for their area of intended use. The list of clinicians and the intended areas of clinical use are:

- Qualified physicians – Clinical applications depending on their medical training
- Qualified nurses – Bladder assessment to determine the volume of the bladder
- Qualified midwives – Assessment of fetal presentation
- Physical therapists – Musculoskeletal (Conventional) for biofeedback of muscle contractions

All users should read the generic ALARA education program supplied with your **Signos RT** (see enclosed ISBN 1-93004 7-71-1, *Medical Ultrasound Safety*). This program outlines the guiding principle for diagnostic ultrasound, where the qualified user keeps ultrasound exposure to “as low as reasonably achievable” while performing a diagnostic examination.

In addition to the above, users must have appropriate training in ultrasound for their specific area of clinical application. It is recommended that users also acquire certification in clinical ultrasound within their area of expertise. Appropriate information on training may be obtained by contacting [Signostics](#) or your local professional body. Examples are:

Australia

- Australasian College of Emergency Medicine (ACEM), www.acem.org.au
- Royal Australian College of Surgeons (RACS), www.surgeons.org
- Royal Australian College of General Practitioners (RACGP), www.racgp.org.au
- Australasian Society for Ultrasound in Medicine (ASUM), www.asum.com.au

Europe

- European Federation of Societies for Ultrasound in Medicine and Biology (EFSUMB), www.efsumb.org

Contact information

Australia



Signostics Limited
 PO Box 1048 Pasadena
 1284 South Road
 Clovelly Park, SA 5042
 Australia
 Telephone +61 8 7424 0600
 Fax +61 8 7424 0601
 Freecall within Australia 1800 SIGNOS (1800 744 667)
 Email support@signostics.com.au
 Website www.signostics.com.au

European Economic Area



Authorized representative:
 Priory Analysts Ltd
 The Pinnacle
 160 Midsummer Boulevard
 Milton Keynes
 MK9 1FF
 United Kingdom
 CE 0805

Rest of the World

Signostics Limited
 PO Box 1048 Pasadena
 1284 South Road
 Clovelly Park, SA 5042
 Australia
 Telephone +61 8 7424 0600
 Fax +61 8 7424 0601
 Email support@signostics.com.au
Website www.signostics.com.au

Trademarks and document copyright

Signostics, **Signos RT** and **SigViewer** are registered trademarks or trademarks of Signostics Limited.

All other trademarks are the property of their respective owners.

© Signostics Limited, 2015.

Labeling

Symbol	Explanation
	Manufacturer
	Manufacturer's declaration of product compliance with applicable EEC directives and the Notified Body reference number
	Intertek nationally recognized testing laboratory certification mark for Canada and United States
	TUV Japanese power supply certification mark
	TUV Rheinland nationally recognized testing laboratory certification mark
	UL Recognized Component mark certifying compliance with Canadian and United States requirements
	Tested to comply with FCC standards
	Class II Equipment
	Type BF applied part (B – Body, F = floating applied part)
	Caution
	Consult instructions for use
	Do not dispose of this product in normal trash or landfill (see European Commission Directive 93/86/EEC). Refer to local regulations for disposal
IPX0 - D	Degree of protection from moisture for display unit. The display unit has no protection against moisture
IPX1 - P	Degree of protection from moisture for the ultrasound probe. The probe is protected against effects of vertically falling water only when the transducer is securely fitted
REF	Part Number
SN	Serial number
	Date of manufacture
	Temperature limitation
	Humidity limitation
	Atmospheric pressure limitation
	Stack box this way up
	Direct Current (DC)
	Alternating Current (AC)

Symbol	Explanation
	Dangerous voltage
	Freeze/Unfreeze Button
	Back Button Control
	Power Button
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician or any other practitioner licensed by law of the State in which he practices to use or order the device

2. GETTING STARTED

Signos RT components



Fig 2.1 Signos components

Display unit components




Fig 2.2 Display unit components

Caution

- ✓ Only use the stylus on the touch screen. Pens, paperclips and other metal objects will damage the screen.
- ✓ Avoid unnecessary bending or winding of the connecting cable.

Display unit controls and functions

Back button	Moves the Signos RT screen one step back towards the Home screen. On the Home and Ultrasound screens, it launches the System Menu. During duplex scanning, it changes modes.
Headphone connector	For connecting headphones.
Loud speaker	Plays sounds and replays voice recordings that may be attached to scanned images.
MicroSD card slot	Contains the microSD card.
Microphone	For voice recording.
Power button	Turns Signos RT on or off when held for 1 second. Will put Signos RT to sleep if pressed briefly. The Power Button also includes an indicator LED with the following possible states: ➤ Solid Green – Device powered on and ready for use or in use ➤ Blinking Green – Device is in low power standby mode and ready for use after touching any button ➤ Solid or blinking amber – Battery low. Device needs to be recharged
Probe LED	The Probe LED is an indicator LED with the following possible states: ➤ Solid Green – Probe is in use ➤ Solid amber – Probe is in use and the battery is low. ➤ Blinking Green – Probe is in low power mode and ready for use ➤ Blinking amber – Probe is in low power mode, and the battery is low.
Reset button	Causes the device to reboot, and will also reset some settings, e.g., time and date.
Scan button	Starts and stops (freezes) scanning.
Scroll wheel	Has separate scrolling and clicking actions and is used for user interface navigation and control.
Stylus holder	Holds the stylus when not in use.
Top panel cover	Protects the connections and controls accessible through the top panel. This cover is designed to fit snugly. To open the top panel cover, face the device toward you and push the top panel cover toward you from the back while lifting at the lip on the front.
Flashlight	Provides general purpose illumination.
Touch screen	Displays ultrasound images, controls and text, which may be controlled by tapping the screen with the stylus.
USB connector	For connecting Signos RT for charging, transferring data or printing.

Turning Signos RT on

Turn **Signos RT** on by pressing and holding the **Power button** for 1 second. A loading screen will show during start-up.

You will be prompted to set [Date and Time](#) and your Locale (or language) the first time you turn on your **Signos RT** (page 24).

The microSD card

Signos RT stores patient and image data on a microSD card. Insert the card into the microSD card slot (Fig 2.3).

You will be asked to format your microSD card the first time you insert it. The **Format Card** screen (Fig 2.4) will appear where you may set:

- The **Security Policy**. The options are **Password Protected** and **No password**.
- A **Card Name**.
- A **User name**.
- A **Password** (if required).

Enter text using the stylus on a keyboard that appears when you select a text box (see [Text box](#), page 20).

Password protection can be used to protect patient data from unauthorized access. Consult local laws or your institutions privacy policies to determine whether to use a password. If you require a password, select and confirm a password. Check the **Show Password Text** check box to view alphanumeric text in the password boxes rather than asterisks. You may change the password at any time (see [Username & Password](#), page 27).

If you choose not to have password-protection, your **Signos RT** will boot directly to the [Ultrasound](#) screen (page 18) if a functioning probe is detected. The information on your microSD card will still be encrypted and not readily available without a **Signos RT** or [SigViewer](#) software (page 57).

If you choose not to use the microSD card, you may still use the **Signos RT** in [Urgent Mode](#) (page 14).

If the inserted microSD card has previously been used in an earlier version of **Signos RT**, an **Updating Repository...** message may appear the first time the card is inserted into the newer **Signos RT**. This data update will begin and finish automatically.

To safely eject the microSD card, use your fingertip to gently depress then release the card, which will pop out of its slot for removal. Do not use any other object, such as the stylus, to depress the card as this may cause the card to eject very quickly out of its slot.

For information on:

- Loading and unloading the microSD card, see [Log in/Log out/Remove Card](#) (page 30).
- Card features, see [Storage Card Info](#) (page 29).
- Changing the **Security Policy**, see [Username & Password](#) (page 27).
- Card names and labels, see [Labels](#) (page 25).

Recovering security key information

The microSD card contains a critical file used to encrypt the data. Data on the card becomes inaccessible if this file is lost or corrupted. [SigViewer](#) software can backup these files to enable data recovery (page 57).

Log in screen

The **Log in** screen (Fig 2.5) appears after start-up if your microSD card is password-protected. The microSD card name and user name appear automatically.




Fig 2.3 microSD card insertion

Fig 2.4 Format Card

Fig 2.5 Log in

Select the **Password** text box to make a [keyboard](#) appear at the bottom of the screen (page 20). Enter the password by tapping the keyboard with your stylus. There is no limit to the number of times an incorrect password may be entered.

Urgent Mode

Press the **Urgent Mode** button to begin scanning immediately in an emergency. You can access **Urgent Mode** if:

- No microSD card is inserted.
- A password-protected microSD card is inserted. **Urgent Mode** bypasses the need for password entry.

Images scanned in **Urgent Mode** are temporarily saved to a small internal memory, which can store approximately 25 images. New images overwrite the oldest images in memory if you exceed this limit. To save these images, insert a microSD card and log in following the prompts. These images will be saved to an **Unidentified** patient and identified by exam date and time. Patient or other data on the microSD card cannot be accessed in **Urgent Mode**.

Save the data before shutting down by inserting a microSD card if none is already inserted, and log in following the prompts. Unsaved data will be lost if you shut down **Signos RT**.

How to insert the transducer

The detachable ultrasound transducer is the section at the end of the probe (Fig 2.6).

The transducer is press-fit, and can be inserted by sliding it into the cavity at the end of the ultrasound probe (Fig 2.7), ensuring the silver dot on the transducer is aligned with the silver dot on the ultrasound probe. The transducer snap-locks into position.

Transducers may be interchanged at any time. Remove by pulling firmly on the transducer. **Signos RT** will automatically detect a transducer change.

The transducer model number and part number are printed on the transducer chassis (S3 – P03479 or S3-5 – P03611).

Warning

- ✓ Check the transducer before use to ensure it is properly attached.

Connecting and disconnecting the Power Supply

Warning

- ✓ Recharge **Signos RT** only with the mains charger provided (**Signostics** Power Supply P03198).
- ✓ Only connect the **Signostics** Power Supply to a mains supply rated at 100-240V and 50-60Hz.
- ✓ Do not use the device or **Signostics** Power Supply if there are signs of damage.

The **Signostics** Power Supply (P03198) is used to provide power to the device for recharging the battery. Plug the power supply into an AC outlet, and the micro-USB connector on the power supply lead into the display unit USB connector under the top panel cover. The power supply will recharge both the display unit and ultrasound probe batteries in approximately 3 hours (see Fig 2.8). You should fully charge the device before using it for the first time. **To disconnect all mains voltages from the equipment remove the USB cable from the display unit USB socket.**



Fig 2.6 S3 or S3-5 transducer



Fig 2.7 Transducer insertion



Fig 2.8 Display Power Supply Connection

3. SCANNING

Setting Up

Before scanning, check you have access to everything you need and the device is ready to scan as follows:

- The **Signos RT** is powered on and has battery charge available to scan (the display unit and probe LEDs are green or flashing green).
- The **Signos RT** probe has been cleaned or disinfected as appropriate.
- You have suitable ultrasound transmission gel available.
- The patient is positioned appropriately and you have access to the anatomy to be scanned.

Default Settings

Each transducer defaults to a set of preset scan settings at initialization which ensure safe operation of the device. The settings are set to values which should be adequate for average patients, and the depth is set to the maximum depth. The user is able to modify the default preset and create new customized presets (see [Scan settings](#) page 42). Scan settings cannot increase acoustic output level, ensuring the device always operates with safe acoustic output levels.

Entering exam and patient information

Enter patient data via the  [Exam management](#) icon on the Tool Bar (page 34) or via **SigViewer** software. Patient data may be entered before beginning a scan, or added at any later time.

Signos RT will prompt you to enter details if you close an [active exam](#) without entering patient data. Scans will automatically be saved to an **Unidentified** patient when no patient data is entered.

Imaging modes

Warning

- ✓ Check the connecting cable, connectors, and system housings before use for cracks or fraying. Do not use if damaged.

Signos RT employs three separate imaging modes, known as B-mode, M-mode and PW Doppler mode.

Changing imaging modes

The current imaging mode is indicated by the yellow letter in the upper right corner of the touch screen. When not scanning, you can change the imaging mode by:

- Tapping the yellow icon on the touch screen and selecting the required imaging mode.
- Holding down the **Scan button** for more than 1 second and selecting the required imaging mode.

Scan button

Click the **Scan button** to begin scanning. The scanned image will constantly be updated on the screen.

When actively scanning, click this button to stop and store the scan. The last scanned image will be retained on the screen and in the exam (see [Exam management](#), page 34).

B-mode Imaging

Signos RT generates conventional sector-shaped B-mode images of underlying anatomy. B-mode scans are started when the *Mode Identifier* is **B**, and the **Scan button** is pressed. The image automatically generated is a sector or pie shape of approximately 90° (see Fig 3.2). When a scan is running, the **Scan button** is pressed to stop the scan and the final image is automatically saved.

For probe orientation, the **Scan Button** on the probe corresponds to the **Flip Icon** on the display. The convention is for the **Scan Button** to be on the patient's right side for transverse scans, and superior for sagittal or coronal scans (Fig 3.1), with a white dot (flip icon) on the image to denote scan button position.




Fig 3.1 B-mode scans with correct probe orientation




Fig 3.2 B-mode scan showing field of view and orientation

M-mode Imaging

M-mode images can record tissue motion over time

With M-mode selected, starting a scan will initiate a B-mode image with a line down the centre of the image showing the position of the M-mode line of sight (Fig 3.3). With the B-mode scan running, either pressing the **Back button** or tapping the **Duplex Thumbnail** will freeze the B-mode image and commence generating the M-mode trace (Fig 3.4). Changing back from the M-mode trace to the B-mode positioning is via pressing the **Back button** or tapping the **Duplex Thumbnail**.

M-mode images show a stationary interface as a straight line on the trace. Heart rates can be measured using the  [Measurement](#) icon (page 38).




Fig 3.3 Positioning scan for M-mode

PW Doppler Imaging

PW Doppler images can record fluid motion over time.

With PW Doppler mode selected, starting a scan will initiate a B-mode image with a line and a gate down the center of the image. The line and gate are a locating aid, assisting the user to align the tissue of interest (Fig 3.5). The gate position can be moved up and down and widened by using the scroll wheel. With the B-mode scan running, either pressing the **Back button** or tapping the **Duplex Thumbnail** will freeze the B-mode image and commence generating a PW Doppler trace (Fig 3.6). Changing from the PW Doppler trace back to the B-mode positioning scan is via pressing the **Back button** or tapping the **Duplex Thumbnail**. A trace depicting fluid velocity is displayed on the image.




Fig 3.4 M-mode trace

Note: PW Doppler does not provide absolute velocity measurements.




Fig 3.5 B-mode scan with PW Doppler line and gate




Fig 3.6 PW Doppler trace

Scroll wheel adjustments

When **Signos RT** is actively scanning, clicking the **Scroll wheel** provides a shortcut to adjust commonly altered [scan settings](#) depending on the mode.

In B-mode and M-mode the first click causes the **Gain** control to appear on the screen (Fig 3.7). The second click replaces this with the **Scan Depth** control. Subsequent clicks cycle through the gain and scan depth controls alternately.

In PW Doppler mode, during the B-mode positioning scan the **Doppler Gate Depth**, **Doppler Gate Width**, **Gain**, and **Scan Depth** can be adjusted. During the Doppler scan, **Doppler Gain**, **Baseline**, Pulse Repetition Frequency (**PRF**), **Gate Depth**, and **Gate Width** can be adjusted.

When an adjustable setting is displayed, scrolling the **Scroll wheel** will adjust the parameter. These controls will disappear after 4 seconds, or when the scan is stopped.




Fig 3.7 Gain control

Recommended ultrasound transmission gel

Signostics recommends the use of *Aquasonic 100* Ultrasound Transmission Gel manufactured by Parker Laboratories, INC. Fairfield, New Jersey 07004. A sample is provided with the system.

Recommended ultrasound probe sterile sheaths

Where fluid contamination is possible, such as during ultrasound guided fluid drainage of vascular access, the ultrasound probe shall be covered with an appropriate sterile sheath, which will promote asepsis and minimize cleaning. Signostics recommends the use of Eclipse Probe Cover – Latex Free 38-03 by Parker Laboratories, INC. Fairfield, New Jersey 07004.

Warning

- ✓ Be aware of latex allergy. Some commercially available probe covers contain latex. Refer to the US FDA alert title: 'Medical Alert: Allergic Reactions to Latex-Containing Medical Devices', issued March 29, 1991.

4. SIGNOS RT SCREENS AND MENUS

Home screen

Four applications may be accessed from the **Home** screen (Fig 4.1). The name of the highlighted icon is displayed in the blue bar at the bottom of the screen, e.g., reads **Ultrasound** in Fig 4.1.

The icons that appear are:



Ultrasound: launches the ultrasound application. This application will launch automatically on start-up if the probe is connected. See the next chapters for a complete description of its functionality.



Printer: opens the **Print Manager**. See [Printer](#) for more information (page 22).



Preferences: opens the **Preferences** menu. See [Preferences](#) for more information (page 23).



Flashlight: turns on the flashlight. See [Flashlight](#) for more information (page 27).

To return to this screen, launch the [System Menu](#) (page 22) by tapping the application icon or screen title in the top left corner of the screen, then selecting  **Home**.



Fig 4.1 Home

Ultrasound screen

The screen in Fig 4.2 appears after you launch the ultrasound application. This application will launch automatically when an ultrasound probe is detected.

- The *Title bar* appears at the top of the touch screen.
- The *Tool Bar* appears at the bottom of the touch screen.
- Ultrasound images are displayed centrally.



Fig 4.2 Ultrasound

Navigating the touch screen

Signos RT communicates with the user in a number of ways.

- Menus
- Message boxes
- Screens and windows
- Tool Bar icons

Most of these options can be navigated either by using the stylus to tap the touch screen, or by using both the scrolling and clicking functions of the **Scroll wheel**. The **Back button** also helps navigate various options.

Menus

Menus can be accessed and navigated in 2 ways.

Tapping the touch screen with the stylus will select menus and menu options. Tap the  **Back** icon that appears in the top right corner of a full screen menu to back out of the current menu.

The **Scroll wheel** can also be used to navigate through menu options. The scrolling action scrolls from one item to the next in a menu; the clicking action selects the highlighted item. Clicking the **Back button** will return you to the previous menu or screen.

Message boxes

System events are communicated via on-screen message boxes. Some message boxes require acknowledgement by pressing the **OK** button.

Screens and windows

Signos RT may open full-size screens or smaller windows that either collect information (e.g., password information on the **Log in** screen) or impart information, (e.g., the **Battery Info** window).

A number of different controls may be found.

- Check box
- Date and time box
- Drop-down box
- List box
- Slider
- Text box

These options can be manipulated using either the stylus or **Scroll wheel**. Tap the  **Close** icon in the top right corner to close a screen or window.

Check box

Select and de-select a check box by tapping either the box or associated text with the stylus; or by clicking the **Scroll wheel** on the highlighted option.

Date and time box

Tap the date or time component you wish to change (e.g., month, hour), then use the up and down arrows at the right end of the box to adjust it. Alternatively, click the **Scroll wheel** on the setting you wish to change, scroll to change the setting, and then click to confirm changes.

Drop-down box

Tap on either the down arrow or the text in the box to access options; select one by tapping its name. Alternatively, click the **Scroll wheel** to open the box, scroll to the desired option, and then click to select it.

List box

Tap the list to access options; select one by tapping its name. Alternatively, click the **Scroll wheel** to select the box, scroll to the desired option, and then click to select it. Perform the required action by pressing the appropriate button at the bottom of the screen or window, e.g., **Edit**, **Delete**, **OK**.

Slider

Move the pointer on the slider to vary a parameter through a range. To use the stylus, press and drag the pointer. To use the **Scroll wheel**, select the slider by clicking, which will highlight the pointer in blue. Scroll the wheel to move the pointer, and then click to accept the changes.

Text box

Text boxes allow you to enter alphanumeric text such as passwords, patient information and image annotations. Tap the text box to enter text. A keyboard will appear at the bottom of the screen. Three different keyboards (Fig 4.3) can appear for entering lowercase text, numbers, or uppercase text and punctuation (Caps).

- To open the numeric keyboard, tap the  icon in the top left corner of the alphabet keyboard.
- Open the Caps keyboard by tapping the **Caps** button on the alphabet keyboard. Tapping the **Shift** button will convert the lowercase keyboard to the Caps keyboard for one character.
- To return to the alphabet keyboard, tap the  icon in the top left corner of the numeric keyboard.

Improve efficiency with [keyboard](#) gestures (page 24).




Fig 4.3 Three keyboards
(left) Lowercase keyboard; (middle) numeric keyboard; (right) Caps keyboard


Tool Bar icons

The Tool Bar displays icons at the bottom of the touch screen; any icons that are not selectable are grayed out. Icons can be selected using the stylus or **Scroll wheel**. The Tool Bar may be hidden; to view the Tool Bar again, click the **Back button** or the **Scroll wheel**; or hold the stylus at the very bottom of the screen. See [The Tool Bar](#) for more information (page 33).

The Title bar

The Title bar (Fig 4.4) is at the top of the touch screen. This bar indicates which application is currently displayed and has status icons that can be activated with the stylus. When in ultrasound mode, the Title bar also indicates which frequency transducer is attached.

Fig 4.4 Title bar



Application icons

The icon displayed in the top left corner indicates which application **Signos RT** is currently displaying. The options are:

-  [Home](#)
-  [Ultrasound](#)
-  [Printer](#)
-  [Preferences](#)
-  [Flashlight](#)

Tapping on this icon always launches the [System Menu](#) (page 22).

Screen title

The screen title shows the name of the currently active application, menu, or screen or window. Tapping on this title always launches the [System Menu](#) (page 22).

Icon tray

A number of icons may appear in the icon tray.

-  **Back** appears on **Signos RT** menus. Tap this icon to revert to the previous menu or screen.
-  **Close** appears on **Signos RT** screens and windows. Tap this icon to close the screen or window.

The following icons may also appear in the Title bar. Tapping these icons will launch a window describing the current status of **Signos RT** settings.

-  [Probe Status](#)
-  [Battery Info](#)
-  [Storage Card Info](#)
-  [Volume Info](#)
-  [Print Manager](#)

When these icons are hidden due to space constraints, they can be accessed through System Menu>[Status icons](#) (page 28). However status icons that are displaying critical information, such as a low battery condition, will never be hidden.

System Menu

The System Menu (Fig 4.5) can be accessed by:

- Tapping the application icon or screen title in the Title bar.
- Clicking the **Back button** while on the **Home** or **Ultrasound** screen.

Home screen

Selecting  **Home** will take you to the [Home screen](#) (page 18).

Ultrasound

Selecting  **Ultrasound** will launch the main ultrasound application.




Fig 4.5 System Menu

Printer

Selecting  **Printer** will take you to the [Print Manager](#) screen (Fig 4.6 and page 30). Printing jobs can be set through **Main Menu>Print Images...** (page 49).

When connected with a PictBridge printer via a USB port, the printer details will appear in the top section of the screen, and print queue details in the bottom of the screen. If no printer is connected, the print queue will display jobs that have been created but not yet printed.

If there is a problem with the printer, such as no paper, printing will pause. Fix the printer problem, acknowledge the error message and press **Continue**.

Select a print job, then press **Job Details** (Fig 4.7) to show information about that job. Press **Delete Job** to delete the print job. You will not be asked to confirm your request.

If the printer does not support the set print parameters, the **Print Settings** window (Fig 4.8) will open to adjust printing preferences. The **Print Settings** window may be launched through **Main Menu>Print Images...**




Fig 4.6 Print Manager




Fig 4.7 Job Details




Fig 4.8 Print Settings

Preferences

Select  **Preferences** to go to the main **Preferences** menu (Fig 4.9). The following options may be adjusted:

- [Audio](#) (page 23)
- [Colors](#) (page 23)
- [Date and Time](#) (page 24)
- [Display](#) (page 24)
- [Keyboard](#) (page 24)
- [Labels](#) (page 25)
- [Locale](#) (page 25)
- [Power](#) (page 26)
- [Security](#) (page 26)
- [Touch Screen Calibration](#) (page 26)
- [USB](#) (page 27)
- [Username & Password](#) (page 27)




Fig 4.9 Preferences

These preferences are saved to the microSD card. In **Urgent Mode**, the preferences for the **Signos RT** are the same as for the last card inserted. Any changes made to the **Preferences** in **Urgent Mode** will not be saved.

Audio

Control the volume of different audio channels on the **Audio** screen (Fig 4.10).

- **Headphones Volume** sets the volume at which sounds are played through headphones.
- **Loud Speaker Volume** sets the volume at which sounds are played through the loudspeaker. Loud speaker output is always muted when headphones are connected.

Check the **Mute Audio Output** box to silence all sounds normally emitting from **Signos RT**.

The **Keyboard/Touch Clicks** box enables click noises when the touch screen and control buttons are used.

- **Record Volume** sets the volume at which sounds are recorded through the microphone.



Fig 4.10 Audio

Colors

Change between **Normal** (color) and **Monochrome** (grayscale) color palettes for the **Signos RT** interface on the **Colors** screen (Fig 4.11).




Fig 4.11 Colors

Date and Time

Configure the time and date settings linked to scanned images on the **Date and Time** screen (Fig 4.12).

- **Time Zone** sets your current time zone relative to UTC (Coordinated Universal Time).
- **Date** sets the date. Select the month/date/year, then use the up and down arrows to adjust settings. The day of the week will automatically appear. Date format can be set via **Preferences>Locale** (page 25).
- **Time** sets the time. Use the stylus to select the hour:minute:second (and AM/PM if appropriate), then use the up and down arrows to adjust settings. The time format can be set via **Preferences>Locale**.

You may also opt to update time and date settings when synchronizing your **Signos RT** with **SigViewer** software.

If date and time are not set when the **Signos RT** is turned on, e.g., if you have reset your device, a similar screen will show on start-up.




Fig 4.12 Date and Time

Display

A dimmed screen uses less power and can extend battery life. Control screen brightness under three discrete power scenarios on the **Display** screen (Fig 4.13).

Use the slider to adjust settings, with values ranging from 0–100.

- **Full Brightness (Battery)**: screen brightness when **Signos RT** is actively being used on its internal battery.
- **Dimmed Brightness (Battery)**: screen brightness after the **Signos RT** screen has dimmed to conserve power on internal battery (see [Power](#), page 26).
- **Brightness (AC)**: screen brightness when **Signos RT** is powered by a mains or USB connection.




Fig 4.13 Display

Keyboard

Keyboard gestures improve time efficiency when using the keyboard, and involve making a stylus stroke on the keyboard part of the screen. Enable keyboard gestures using the **Keyboard** screen (Fig 4.14).

- An upward stroke on a letter in the keyboard will generate the uppercase letter.
- A forward stroke will generate a space.
- A backward stroke will generate a backspace.
- A downward stroke equals Return or Enter.




Fig 4.14 Keyboard

Labels

The **Labels** screen (Fig 4.15) allows you to edit the card name and to add labels to the card and device, e.g., contact information in case you misplace your card or **Signos RT**.

- **Device Label** adds a label to the **Signos RT**. This label is visible only on this screen.
- **Card Name** allows you to edit or change the card name that appears on the **Log in** screen.
- **Card Label** adds a label to the microSD card. This text is saved in an unencrypted file called **CardLabel.txt** on the microSD card; if the card is found and inserted into a PC, the information on the label can be retrieved. Note that patient data on the microSD card is encrypted and cannot be retrieved without a **Signos RT** device or [SigViewer](#) software_ (page 57).



Fig 4.15 Labels

Locale

Set language, time and date preferences on the **Locale** screen (Fig 4.16).

You will need to reboot the system for locale changes to take effect.



Fig 4.16 Locale

Power

Configure power settings to optimize performance and battery life on the **Power** screen (Fig 4.17). Turn off power-saving options by choosing **Disabled** from the relevant drop-down list.

- Backlight dimming saves power by reducing the brightness of the screen after a period of inactivity when **Signos RT** is running on battery. The screen returns to normal brightness immediately on use. Set the time after which you would like **Signos RT** to dim using the first drop-down box.
- Sleep mode saves power by turning off the screen when **Signos RT** is not being used. Use the **Sleep Timeout** drop-down boxes to set the time after which you would like **Signos RT** to sleep. Different timeouts can be set for battery or AC power. **Signos RT** will also go to sleep if you depress the **Power button** briefly.

Signos RT wakes from sleep instantly when any of the following are clicked: **Back button**; **Scan button**; **Scroll wheel**; **Power button**. Tapping the screen or moving the **Scroll wheel** will not wake **Signos RT** from sleep to prevent accidental awakenings in transit. **Signos RT** will charge more quickly in sleep than if the screen is lit.

Note that the backlight and sleep timeouts start to countdown together; if the sleep timeout is shorter than the backlight timeout, **Signos RT** will go to sleep without dimming the screen.




Fig 4.17 Power

Security

Vary the security settings protecting the microSD card data on the **Security** screen (Fig 4.18). Turn off security options by selecting **Disabled** from the relevant drop-down list.

- Use the **Idle Timeout** drop-down box to set **Signos RT** to ask for your password after a certain time of inactivity has elapsed.
- Use the **Password Timeout** drop-down box to set **Signos RT** to ask for your password after a certain amount of uptime has elapsed.

You will need to enter the password before changes will take effect.




Fig 4.18 Security

Touch Screen Calibration

The **Touch Screen Calibration** screen (Fig 4.19) allows re-calibration of touch input with respect to the touch screen.

The cursor will appear as a + cross in the middle of the screen. Use the stylus to touch and hold the center of this cross until it moves. The cross will then move to each corner to enable calibration.

Exit the **Touch Screen Calibration** screen at any time by clicking the **Back button**.




Fig 4.19 Touch Screen Calibration

USB

Set the default function that occurs when **Signos RT** is connected to a PC or PictBridge printer on the **USB** screen (Fig 4.20).

- **Ask Each Time:** after connection to USB, the following three options will appear in a new screen (Fig 4.21).
 - **File Transfer:** data from the microSD card is made available to the PC. **Signos RT** will display a screen to show that it is being used by the connected computer. Disconnect **Signos RT** from the USB cable to return to normal function.
 - **Print:** connects **Signos RT** to a PictBridge printer, and any queued print jobs will begin to print. You may not perform scans while connected to a printer, but other **Signos RT** functions are retained.
 - **Remote Display:** will launch **SigViewer for Signos RT** on the connected PC. A new window will open displaying the contents of the **Signos RT** screen. Your **Signos RT** will not respond to any keyboard commands or mouse clicks made on the **PC**. See the **SigViewer** User Manual for more information. **Note:** this mode is primarily intended for reviewing images. Although a scan can be performed while connected to a PC in remote display mode, please be aware that Signos RT only conforms to CISPR11 RF Emissions Class A while scanning in this configuration. Refer to "Electromagnetic compatibility (EMC)" in Section 11.

You will need to exit this screen for changes to take effect.




Fig 4.20 USB




Fig 4.21 Select USB Function

Username & Password

Change the security policy, user name and password for the microSD card on the **Username & Password** screen (Fig 4.22).

Use the **Security Policy** drop-down box to change the requirement for a password on the microSD card.

Press the **Change** button to save changes, or **Cancel** to exit the screen. If the card is password protected, you will need to re-enter the current password for changes to take effect.

This **Preferences** option is not available in **Urgent Mode**.




Fig 4.22 Username & Password

Flashlight

Turn on the  **Flashlight** from either the **Home** screen or the System Menu. Turn off the flashlight by clicking the **Scroll wheel**; subsequent clicks will toggle the light on and off. An on-screen image will also indicate whether the light is on or off (Fig 4.23).




Fig 4.23 Flashlight

Status icons

Five windows can be accessed from this menu option (Fig 4.24).

-  [Probe Status](#) (page 28)
-  [Battery Info](#) (page 29)
-  [Storage Card Info](#) (page 29)
-  [Volume Info](#) (page 30)
-  [Print Manager](#) (page 30)




Fig 4.24 Status icons

Probe Status

Information on the  **Probe status** window (Fig 4.25) may need to be accessed if there is a fault with the **Signos RT**.

- **Serial #:** the unique serial number of the probe.
- **Probe Status:** reports whether the probe is functioning correctly or has experienced an error.
- **HW Revision #:** the revision of the probe electronics hardware.
- **SW Version #:** the version of the probe software.
- **BL Version #:** the version of the probe bootloader.
- **CPLD Version #:** the version of the probe programmable logic device.




Fig 4.25 Probe Status

A **Probe Status** icon will appear in the Title bar if **Signos RT** detects an error with the probe. The icons are:

-  Probe functioning normally.
-  (flashing) Error detected with probe (see [Troubleshooting](#), page 52).
-  (flashing) Probe missing or not detected (see [Troubleshooting](#), page 52).

If there is an error with the probe, contact [Signostics](#) (page 55), or your **Signos RT** distributor to arrange for repair or replacement.

Battery Info

The  **Battery Info** window (Fig 4.26) shows the current capacity of the device batteries; whether the device is running off battery, USB or mains; battery voltages and current; and whether the batteries are currently charging.

Press the **Preferences** button to open the [Power](#) screen (page 26).

There are five different battery icons that appear in the Title bar.

-  (bars moving up) **Signos RT**' internal batteries are being charged.
-  **Signos RT** is connected to mains power and the batteries are not charging.
-  **Signos RT** is connected via USB to a power supply, e.g., through a PC, and the batteries are not charging.
-  Error detected with the battery (see [Troubleshooting](#), page 52).
-  **Signos RT** is running on battery. The number of bars indicates charge level: a fully charged battery is represented by ; low battery level is represented by  (flashing).




Fig 4.26 Battery info

A message will appear to indicate battery level is critically low when the battery capacity falls to 8%.

If **Signos RT** is off when connected to a charger, a screen will briefly appear to indicate charging. The **Power button** will flash while charging continues.

Storage Card Info

The  **Storage Card Info** window (Fig 4.27) shows:

- Card status
 - The  **Secure Card** icon indicates files are being saved to a microSD card (with or without password-protection).
 - In **Urgent Mode**, a  gray padlock will appear, with the text **Urgent Mode**.
- **Total space** on the card, free and used.




Fig 4.27 Storage Card Info

The icons that may appear in the Title bar are:

-  Card loaded and 25% full. The amount of blue on the disk indicates used card space: an empty card is represented by ; a full card is represented by .
-  **Urgent Mode**.
-  (flashing) Storage card error.

Volume Info

The  **Volume Info** window (Fig 4.28) allows you to:

- Adjust the **Volume** of the loud speaker or headphones (whichever is being used) using the slider. You will hear tone volume as the slider moves.
- Mute **Signos RT**.

Press the **Preferences** button to launch the **Preferences**>[Audio](#) screen (page 23).

The various icons that may appear in the Title bar are:

-  **Signos RT** set to play tones.
-  **Signos RT** muted.




Fig 4.28 Volume Info

Print Manager

The  **Print Manager** screen (Fig 4.29) shows the status of current print jobs. For further description of this window, see [Printer](#) (page 22).

The various icons that may appear in the Title bar are:

-  (*flashing*) Printing.
-  Print jobs queued but no printer connected.
-  Error with printer, e.g., no paper.




Fig 4.29 Print Manager

Log in/Log out/Remove Card

Log in appears on the System Menu when the **Signos RT** is in **Urgent Mode**; use this option to load a password-protected card.

Log out appears when a password-protected microSD card is loaded; use this option to safely unload the card.

Remove Card appears on the System Menu when a card without a password is loaded; use this option to safely unload the card.

A message will appear when the card may safely be removed. If you eject a card without selecting **Log out** or **Remove Card**, **Signos RT** will prompt you to re-insert the microSD card; failure to do so will cause you to lose unsaved data.

Sleep

Select **Sleep** to immediately put **Signos RT** to sleep. For more information on Sleep mode, see [Power](#) (page 26).

Shutdown

Select **Shutdown** to turn **Signos RT** off. You will be asked to confirm the **Shutdown** request.

About

Select **About** to view the exact model specifications of your **Signos RT** and the unique hardware serial number of your device (Fig 4.30). The **Details** button will launch a **Help** page describing **Attributions**.




Fig 4.30 About hardware

Turning Signos RT off

To turn **Signos RT** off, either:

- Hold down the **Power button** for 1 second; or
- Select **Shutdown** from the [System Menu](#) (page 22).

5. CONTROLS FOR THE ULTRASOUND APPLICATION

Display unit controls

Touch screen

The controls and text displayed on the touch screen are selected either by tapping with the stylus, or by using the scrolling and clicking actions of the **Scroll wheel**.

In an [active exam](#) in the ultrasound application, the following screen will appear (Fig 5.1).




Fig 5.1 Key elements of the ultrasound screen

- The **Patient name** displays the patient name or if that is not set the patient EMR Id in brackets.
- The **Exam identifier** displays the time and: **[Active]** for an active exam that is currently recording scan images; or **[Review]** for a closed exam that has been re-opened for review. For more information on exams, see [Exam management](#) (page 34). Tap either the **Patient name** or **Exam identifier** to open the [Exam pop-up menu](#) (page 50).
- The **Flip** icon shows the left–right orientation of the image. Tap the **Flip** icon to change image orientation. The icon on the left side of the screen indicates the original orientation; the icon on the right side of the screen indicates flipped orientation.
- The **Image name** displays the text [annotation](#) associated with the image (page 40). When an image is first created, this annotation is set to a numbered default (Image 1, 2, 3 etc). Tap the **Image name** to open the [Image pop-up menu](#) (page 50).
- The **Preset name** displays the name of the current scan [preset](#) (page 42). If the current scan settings have not been saved into a preset, the name will display an asterisk. Tap the **Preset name** to open the [Settings pop-up menu](#) (page 50).
- The **Mode identifier** indicates whether the **Signos RT** is currently scanning in B-mode, M-mode, or PW Doppler. Tap the **Mode identifier** to open the [Mode pop-up menu](#) (page 50).

Back button

Click the **Back button** to move **Signos RT** one step back towards the **Home** screen. On the **Home** and **Ultrasound** screens, it launches the System Menu. While M-mode or PW Doppler scans are in progress, use the **Back Button** to toggle between B-mode scans and M-mode/PW traces.

Power button

Press the **Power button** for 1 second to switch **Signos RT** on and off. Press the **Power button** briefly to put **Signos RT** to sleep.

Scan button

Click the **Scan button** to begin scanning. The scanned image will constantly be updated on the screen.

When actively scanning, click this button to stop and store the scan. The last scanned image will be retained on the screen and in the exam (see [Exam management](#), page 34).

Scroll wheel

The **Scroll wheel** has separate scrolling and clicking actions. The function of the **Scroll wheel** varies depending on whether or not **Signos RT** is actively scanning.

When **Signos RT** is not actively scanning, the **Scroll wheel** enables general navigation and selection of on-screen options. The scrolling action scrolls you from one item to the next and the clicking action selects the highlighted item.

When **Signos RT** is actively scanning, clicking the **Scroll wheel** provides a shortcut to adjust commonly altered [scan settings](#) (see page 17).

The Tool Bar

The Tool Bar displays icons at the bottom of the touch screen; any icons that cannot be selected are grayed out. Icons can be selected using the stylus or **Scroll wheel**. The icons that appear when you first open the ultrasound application are shown in Fig 5.2.




Fig 5.2 Tool Bar icons

Some icons will launch a new screen or menu; others will replace the Tool Bar with new controls pertinent to their function.

The Tool Bar may be hidden; to view the Tool Bar again, click the **Back button** or the **Scroll wheel**; or hold the stylus at the very bottom of the screen.

Main Menu

Selecting the  **Main Menu** icon will launch the screen in Fig 5.3. The menu options are:

- **Scan Presets**
- **Start Sequence**
- **Grid and Rulers**
- **Image**
- **Print Images...**
- **Print Screen**
- **Review Last Scan**
- **Preferences**
- **About...**
- **Home**




Fig 5.3 Main Menu

Each of these options will be discussed in full in the next chapter, [The Main Menu](#) (page 42).

Exam management

Select  **Exam management** to manage information relating to patients and exams.

A **Signos RT exam** describes a single interaction (e.g., appointment) between a patient and **Signos RT** user. You may save up to 50 different ultrasound images in a single exam.

A *patient* is an individual whose details may be saved on the microSD card for linking with one or more exams. Fig 5.4 shows the relationship between patients, exams and images (see also [Starting a new exam](#), page 35).




Fig 5.4 Schematic relationship between patients, exams and images

Exams are identified by patient name; different exams linked with the same patient are identified by exam date and time. If you do not enter patient details for an exam, the exam will be linked with an **Unidentified** patient. **Unidentified** patient exams are identified by exam date and time.

The exam to which images are being saved is called the *active exam*. There can only be one active exam open at a time, and there is always an active exam open.

After an active exam is closed, no more images may be added but it may be opened for review and editing. Up to five exams may be open at one time. When multiple exams are open, the exam image displayed on the screen is called the *current exam*.

Exams may be set to be either *editable* or read only (see

Preferences, page 49). An active exam can always be edited. For exams:

- Images can be associated with [measurement](#) calipers or shapes (page 38), or with text/voice [recordings](#) (page 40).
- Images can be [deleted](#) (page 41).
- Patient details linked with that exam can be [changed](#) (page 36).

Clicking the **Scan button** will usually jump **Signos RT** to the [active exam](#). The exception arises when jumping to the active exam may cause you to lose unsaved data, e.g., if you are managing scan presets. You will be prompted to save data before beginning to scan.

Exam management screen

Select the  **Exam Management** icon to open the **Exam Management** screen (Fig 5.5). The six buttons are:

- [New Exam](#) (page 35)
- [Patient Details](#) (page 36)
- [Review Exam](#) (page 36)
- [Close Exam](#) (page 36)
- [Close All Exams](#) (page 36)
- Close




Fig 5.5 Exam Management

The first five options are described in the next sections. The last option, **Close**, will exit the **Exam Management** screen.

If you start scanning without going through the **Exam Management** screen, the saved images will be linked with an **Unidentified** patient.

Urgent Mode

In **Urgent Mode**, only three options are available.

- **New Exam**, which starts a new [active exam](#) for an **Unidentified** patient.
- **Switch to Secure Mode**, which will take you to the **Log in** screen in order to load a microSD card.
- **Close**.




Fig 5.6 Patient Search

Starting a new exam

Press the **New Exam** button to start a new [active exam](#). The **Patient Search** screen (Fig 5.6) will appear which allows you to:

- Choose a patient whose details are already stored.
- Select **New** to create a new patient.

To search for a patient start typing any part of the name in the **Patient Name Search** text box; patients whose name contains the typed text will appear (Fig 5.7). Tap the  button to the right of the patient's name to

view complete patient details, e.g., for two patients with the same name. Select the patient's name from the list box and press **OK**.

To create a new patient, tap the **New** button. The **Edit New Patient** screen (Fig 5.8) will appear which allows you to enter the details for a new patient.

To enter details for a new patient, use the text boxes, date box and drop-down boxes provided. Tap the **OK** button to accept the changes and begin scanning.

Press **Cancel** to exit without setting patient details. Images in the active exam will be saved to an **Unidentified** patient.

Editing patient data

Press the **Patient Details** button to edit patient information for the [current exam](#).

- If no patient details are set for this exam, the **Edit New Patient** screen will open (see Fig 5.8). Use this option to link **Unidentified** patients with new or existing patient details.
- If patient details are set for this exam, the **Edit Patient** screen will open (similar to Fig 5.8).

Reviewing exams

Press the **Review Exam** button to open an existing exam. This option will launch the **Exam Search** screen (Fig 5.9) from which you may:

- See all exams linked to a patient.
- See a list of all exams saved on the card.
- Delete exams.

To view all exams linked with one patient, select the desired patient name from the list box (see [Starting a new exam](#), page 35). The most recent linked exam will appear in the **Exams** box; tap this box to view a list of all linked exams.

Press the **All Exams...** button to view a list of all exams saved to the loaded card (Fig 5.10). Exams will be identified by date, time and patient name. To return to the patient list, press the **Lookup Patient...** button.

To delete exams, press the **Delete** button to open the **Delete Options** pop-up menu. Only closed exams may be deleted. You will be asked to confirm delete requests.

- Use **Delete Exam** to delete a single exam.
- Use **Delete Oldest** to create space on the microSD card. At least 10% of card space is freed by deleting the oldest exams, i.e., 200MB on a 2GB card.
- Use **Delete All** to delete all but the open exam(s).

Closing exams

Select the **Close Exam** button to close the [current exam](#).

Select the **Close All Exams** button to close all open exams. A new [active exam](#) will automatically start.




Fig 5.7 Patient Searching



Fig 5.8 Edit Patient




Fig 5.9 Exam Search




Fig 5.10 Exam Search: All Exams...

Image/Exam Scroller

This function allows you to:

- Scroll through images in the [current exam](#) (*image scroller*).
- Toggle between open exams (*exam scroller*). These exams may include the [active exam](#).

Select the  **Image/Exam Scroller** icon to launch the image scroller (Fig 5.11), which shows thumbnails of images in the current exam. Three thumbnails are shown at a time, with the current image in the middle. Scroll through the thumbnails of all images in the exam by selecting the arrows at either end of the scroller. The screen will automatically display the image in the center of the scroller. Alternatively, use the stylus to tap the image you wish to display.

To activate the exam scroller, either click the **Scroll wheel** or tap the small square to the left of the scroller. The exam scroller will then appear (Fig 5.12), with exams identified by patient surname and date of exam.

You can also toggle between images using the [Image pop-up menu](#), or between open exams using the [Exam pop-up menu](#) (page 50).

When the current image is a saved movie clip, the movie will begin playing as soon as the image is viewed.

Movie Clip Controls

Select the  **Movie Clip Controls** icon to save a movie or enter Scan Review Mode.

Saving a Movie

After selecting **Movie Clip Controls**, select “**Save as a movie clip**” (Fig 5.13).

Scan Review Mode

After selecting **Movie Clip Controls**, select “**Enter scan review**” to enable the **Scan Review Mode** Toolbar (Fig 5.14). This mode allows you to:

- Scroll through images in the movie buffer rotating the **Scroll wheel**.



- Click the **Back Button** to disable the image review mode. Scroll down to access other scan review functions.
- Save the currently displayed image by selecting **the Save Single Frame** icon .
- Save the movie clip by selecting the **Save Movie** icon .
- Exit Scan Review mode by selecting the **Exit scan Review** icon .




Fig 5.11 Image scroller




Fig 5.12 Exam scroller




Fig 5.13 Movie Control Selection




Fig 5.14 Movie Scan Review Toolbar

Measurement

Select the  **Measurement** icon to insert calipers or shapes to measure distance, area and volume on B-mode images (Fig 5.15).

- [Calipers](#) measure a straight line distance between two points (page 38).
- [Arrow](#) inserts an arrow marker that can be used to highlight part of the image, but does not actually perform a measurement (page 39).
- [Polygon](#) produces a shape with any number of corners that can measure both area and perimeter (page 39).
- [Ellipse](#) produces an elliptical shape that can measure both area and circumference (page 39).
- [Bladder Volume](#) places an outline around a bladder (or other anechoic region) that estimates volume from a single image (page 39).
- [Beats Per Minute \(BPM\)](#) places a straight line between two points on an M-mode, or PW Doppler trace and calculates the beats per minute based on the time represented by the straight line (page 39).
- [Fetal Biometry](#) places a special caliper, or ellipse, measurement on a B-mode image that estimates gestational age.




Fig 5.15 Measurement pop-up menu

Note that volume measurements are also possible using a bladder volume [sequence](#) (page 46). See also M-mode and PW Doppler image measurement (page 40).

You can insert four measurement shapes per image. Tap the image outside the shape to end editing that shape. Re-select a shape to open it for editing.

The units of measurement are cm for distance, cm² for area, and ml for volume; the values appear in the bottom left-hand corner of the screen. Measurement shapes are numbered in the order in which they were added to the image. The shape, its number and its value/comment are drawn in the same color.

Measurements are more [accurate](#) over a longer distance (page 59). Use the [pan/zoom](#) controls to ensure the object being measured fills as much of the screen as possible (page 37).

The other menu options:

- **Delete All Measurements** deletes all measurement shapes associated with the image.
- **Delete Measurement** deletes the selected shape from the image.

Using measurements

Calipers

A caliper measures linear distance on a B-mode image. Select **Calipers** from the [Measurement pop-up menu](#) to add a caliper (Fig 5.16).

The distance between the points is displayed in the bottom left corner. No measurement value will show until the shape has been moved or re-shaped. Move the caliper using the:

-  Cross-hair handle to position each end of the caliper in the correct place.
-  Square handle in the middle to move the entire caliper around the screen.
-  Double arrow handle to rotate the caliper around its middle.




Fig 5.16 Calipers

Arrow

An arrow is simply used to highlight part of the image, but does not actually perform a measurement. Move the arrow using the:

-  Square handle at the tip of the arrow to move the entire arrow around the screen.
-  Double arrow handle to rotate the arrow around its tip.

Polygon

A polygon shape measures area and perimeter on a B-mode image. Select **Polygon** from the [Measurement pop-up menu](#) to add a polygon shape (Fig 5.18).

The enclosed area (**A:**) and perimeter (**C:**) measurement is displayed in the bottom left corner. No measurement value will show until the shape has been moved or re-shaped from the original insert. Move the polygon using the:

-  Square handle in the middle to move the entire polygon around the screen.
-  Circular handle to position each corner of the polygon in the correct place.

Tap the side of the polygon to add a corner. Remove existing corners by tapping them. A polygon has a minimum of three corners.

Ellipse

An ellipse shape measures area and circumference on a B-mode image. Select **Ellipse** from the [Measurement pop-up menu](#) to add an ellipse shape (Fig 5.19).

The enclosed area (**A:**) and circumference (**C:**) measurement is displayed in the bottom left corner. No measurement value will show until the shape has been moved or re-shaped from the original insert. Move the ellipse using the:

-  Double arrow handles on the circumference to make the ellipse the correct shape.
-  Square handle in the middle to move the entire ellipse around the screen.
-  Double arrow handle inside the shape to rotate the ellipse around its middle.

Bladder Volume

A Bladder Volume shape estimates bladder volume on a B-mode image. Select **Bladder Volume** from the [Measurement pop-up menu](#) to insert a shape. You will need to tap the screen to indicate bladder position; try to tap as closely as possible to the centre of the dark, anechoic region. An outline bordering the bladder and an estimate of bladder volume will appear. Fig 5.20 shows a Bladder Volume shape outlining the enclosed area; the estimated volume measurement appears in the bottom left corner.

Change the Bladder Volume shape by dragging the outline with your stylus. Re-shaping works best by moving the outline in and out, rather than along the border.




Fig 5.17 Arrow




Fig 5.18 Polygon




Fig 5.19 Ellipse




Fig 5.20 Bladder Volume

The bladder volume calculated from a single scan is an estimation assuming the bladder shape is symmetrical. This measurement should only be used as an indicative value, as asymmetrical bladders could have actual volumes very different to the single plane estimation.

M-mode and PW Doppler image measurement

Circumference and area measurements are not possible on M-mode or PW Doppler images since the horizontal axis represents time. You can still insert polygon, ellipse and bladder volume measurements to highlight a region of interest, but the reading in the bottom left corner will read **No Measurements Possible**. Shapes may not be rotated on these images.

The calipers can be used to measure either distance (in cm) and time (in milliseconds) in M-mode, or just time (in milliseconds) for PW Doppler; or heart rate (in beats per minute (bpm)) (Fig 5.21)). Beats per minute measurements are added using the Beats Per Minute (BPM) menu. The following options are available:

- **BPM - single cycle** measures heart rate over one cardiac cycle.
- **BPM – two cycles** measures heart rate over two cardiac cycles.
- **BPM – four cycles** measures heart rate over four cardiac cycles.




Fig 5.21 M-mode caliper measurement

Fetal Biometry Measurements

Fetal biometry measurements allow estimation of the gestational age by performing measurements of fetal anatomy.

Measurements can be converted to gestational age using one of three predefined charts. The chart to be used must be selected in the Ultrasound preferences before fetal biometry measurements can be made (see [Preferences](#)).

The following fetal biometry measurements are provided:

- **Crown Rump Length (CRL)**
- **Biparietal Diameter (BPD)**
- **Head Circumference (HC)**
- **Femur Length (FL)**
- **Abdominal Circumference (AC)**




Fig 5.22 Fetal Biometry Measurement

Annotations

Annotations allow you to link text or voice recordings with an image.

Text annotations

Text annotations allow you to label images. Text appears in a single line at the bottom of the screen (Fig 5.23). This label becomes the image name when printing, toggling between images, etc. You may choose an annotation from the preset list, or you may enter your own label.

Select **A Image Annotations** to open the Annotations pop-up menu containing preset annotations. Modify the included text labels via **Main Menu>Scan Presets>Manage Presets...** (page 44).




Fig 5.23 Text annotation

Alternatively, select **Text Annotation...** at the bottom of the menu to open a keyboard to enter your own label. Multiple lines of text can be entered, but only one line of text will be displayed. An ellipsis will appear at the end of the line if on-screen text has been truncated. Select **Text Annotation...** to view and edit the complete text annotation again.

Every image must be named. If a text annotation is entirely deleted the image name return to its default, eg “Image 1”.

Tap the text annotation on the screen to launch the [Image pop-up menu](#) (see page 50).




Fig 5.24 Voice Recording window

Voice recordings

Voice recordings allow you to link a recorded message with an image. Images with linked voice recordings will display a  **Voice Recording** icon in the bottom left corner (Fig 5.25). Tap this icon to launch the **Voice Recording** window and play the existing voice recording.

To record a new voice recording, select **Voice Recording...** at the bottom of the Annotations pop-up menu and press the  **Record** button. Speak clearly into the **Signos RT** microphone. Press the  **Stop** button to end recording. The slider on the right of the pop-up menu will indicate the length of the recording. To add further voice recordings, simply press **Record** again; the new recording will be added to the existing recording.

To play a voice recording, press the  **Play** button or tap the  **Voice Recording** icon.

To delete all of the recording, press the  **Delete** button.




Fig 5.25 Voice recording

Delete Image

Select the  **Delete Image** icon to erase the current image. You will not be asked to confirm the delete request. Undo the deletion via **Main Menu>Image>Undo>Delete Image**. If all images in an exam are deleted, that exam will also be deleted; however, patient details will be retained.

Hide Tool Bar

Select the  **Hide Tool Bar** icon to minimize the Tool Bar. Display the Tool Bar again by:

- Clicking the **Back button**; or,
- Clicking the **Scroll wheel**; or,
- Tapping the stylus on the very narrow blue bar at the bottom of the screen.

Note that use of the **Back button** or **Scroll wheel** does not register either of these controls' normal functions, but will serve merely to unhide the Tool Bar.

6. THE MAIN MENU

Select  **Main Menu** to launch the **Main Menu** (Fig 6.1). The menu options are:

- [Scan Presets](#) (page 42)
- [Start Sequence](#) (page 46)
- [Grid and Rulers](#) (page 48)
- [Image](#) (page 48)
- [Print Images...](#) (page 49)
- [Print Screen](#) (page 49)
- [Preferences](#) (page 49)
- [About...](#) (page 49)
- [Home](#) (page 49)




Fig 6.1 Main Menu

Scan settings

You can adjust ultrasound scan settings for different scan types to achieve optimal image quality. See page 69 for definitions and parameters of adjustable [scan settings](#).

Several commonly used scan presets are hardwired into **Signos RT**; these are called *system-defined presets* and have a  padlock icon next to their name. [System-defined presets](#) may not be changed but can be copied and edited under a new name (page 45). You may also set and save your own scan settings.

The scan settings used by **Signos RT** at any time are described as *active settings*.

Main Menu>Scan Presets will open the **Scan Presets** menu (Fig 6.2).

- Choose [Edit Scan Settings...](#) to modify the active settings (page 43).
- Choose [Manage Presets](#) to create, edit or delete presets (page 44).
- Choose [Select Preset](#) to choose a preset to load (page 44).
- Choose **Default Presets...** to set the default preset for each transducer (Fig 6.3). Select the desired preset from the drop-down list provided. Press **OK** to apply changes.

The default preset automatically loads when a new exam starts. The default preset for the new transducer will also automatically load if the transducer is changed. Presets may not be deleted or have their target transducer edited if they are currently set as default presets.

- Choose **Diagnostic Area...** to define which [system-defined presets](#) appear on your **Signos RT** (Fig 6.4). You may select between **General Medical**, **Veterinary**, **Physiotherapy**, or **Obstetric**. Press **OK** to apply changes.

Changing the **Diagnostic Area** will restore all system-defined presets and default presets to the default configuration. User-defined presets will not be affected.

- Choose **Revert to Factory Presets** to restore any deleted system-




Fig 6.2 Scan Presets




Fig 6.3 Default Presets




Fig 6.4 Diagnostic Area...

defined presets and restore factory default presets for each transducer. User-defined presets will not be affected. You will be asked to confirm this choice via the **Factory Presets** window.

You can also modify active settings and load presets via the [Settings pop-up menu](#) that appears on the touch screen in an [active exam](#) (page 67).

Changing and saving scan settings

Main Menu>Scan Presets>Edit Scan Settings... will open the **Edit Settings** screen, where you can edit and save all adjustable scan settings (Fig 6.5). Note that when editing the active scan settings only settings that apply to the current scan mode will be available.

Use the scrollbar at the right side of the screen to scroll up and down the screen. The current setting is shown in large font at the right of the screen. Use the sliders or drop-down boxes to change the current settings.

Use the **Save As** button at the bottom of the screen to save your preset. You will be prompted to name your preset when you save it. Scan settings are saved to the microSD card.

If you use the **Back button** or **Close** icon to exit this screen, you will not be prompted to save your changes to a new preset. Instead, the **Preset name** in the top right corner of the screen will show an asterisk next to its name, which indicates unsaved changes have been made to the named preset. To save the active settings, go to **Edit Scan Settings>Save As**.

This screen can also be accessed:

- Through **Main Menu>Scan Presets>Manage Presets...>Edit**.
- Through the first option on the [Settings pop-up menu](#) (page 50).




Fig 6.5 Edit Settings

Customizing scan presets

Main Menu>Scan Presets>Manage Preset will open the **Manage Presets** screen, where you can customize presets (Fig 6.6).

To edit an existing scan preset, highlight its name and then press **Edit...**

[System-defined presets](#) may not be edited.

To create a new scan preset, highlight the name of an existing preset on which you wish to base the new preset, then press **Copy...**

Press **Delete** to delete a preset.

Re-order the presets in the list using the  and  arrow buttons.

Selecting either **Edit...** or **Copy...** will open the [Edit Settings](#) screen (Fig 6.5). Change the desired scan settings, then press:

- **Manage Image Annotations** to change the text annotations that appear when  **Image Annotations** is selected from the Tool Bar. Press:
 - **New...** to add a new annotation for this preset.
 - **Remove** to delete the selected entry.
 - **OK** to save changes and exit.
 - **Cancel** to exit without saving changes.

For more information on:

- [Annotations](#), see page 40.
- Annotations linked to [system-defined presets](#), see page 45.
- **Save As** to save the changes in a new preset. You will be prompted to name your new preset.
- **Save** to save your changes, over-writing the original preset. If you have changed the name of the preset in the text box at the top of the screen, the preset will be saved to the new name. This option is only available if you have accessed the **Edit Settings** screen using the **Edit...** button.

Scan button and the **Scroll wheel shortcuts** can be used to adjust imaging modes and some scan settings (page 50).

Loading a preset

To load a scan preset, you can:

- Navigate **Main Menu>Scan Presets>Select Preset** (Fig 6.7). Highlight the name of the desired preset.
- Open the [Settings pop-up menu](#), and then tap the name of the desired preset.

Alternatively, if you have a saved image with scan settings that you wish to re-load, bring up the required image on the touch screen using **Exam Management>Review Exam**. With this image on the screen, select **Main Menu>Image>Image Information...** (page 48) or open the [Image pop-up menu](#) (page 50). Scroll to the bottom of the page, and press **Use Image Scan Settings**. The word **Imported*** will appear in the **Preset name**.




Fig 6.6 Manage Presets




Fig 6.7 Select Preset

System-defined presets

System-defined presets, marked with a  padlock icon, are hardwired into **Signos RT** and may not be edited. Different presets appear depending on which [Diagnostic Area](#) you have selected (page 42).

- Twelve presets are available for **General Medical**.
- Ten presets are available for **Veterinary**.
- Seven presets are available for **Physiotherapy**. These are a subset of the **General Medical** presets.

For the scan settings and image annotations of each preset, see [Appendix 2](#) (page 71).

Delete system-defined presets using **Main Menu>Scan Presets>[Manage Presets](#)**. Restore any deleted system-defined presets using **Main Menu>Scan Presets>[Revert to Factory Presets](#)** (page 42).

Sequences

Two frequently used scan *sequences* have been pre-set into **Signos RT**. Each sequence consists of a number of *steps*.

- **Bladder Auto**: allows you to estimate bladder volume using two **Bladder Volume** shapes on two images (page 39).
- **Bladder Calipers**: allows you to estimate bladder volume using three **caliper** measurements from two images (page 38).

Choose **Main Menu>Start sequence** and select a sequence option to begin running the sequence (Fig 6.8).

When a sequence is running, new icons will appear in the Tool Bar that allow you to manipulate sequence steps.




Fig 6.8 Start Sequence

-  **Previous Sequence Step** moves back to the previous step in the sequence.
-  **Delete and Redo** deletes the image saved in the current sequence step and allows you to re-scan this image. Click the **Back button** to cancel.
-  **Next Sequence Step** skips to the next step in the sequence.
-  **End Sequence** ends the sequence at the current step.
-  **Repeat Sequence** repeats the sequence from the first step. All images from the original sequence are saved.

Only one image may be saved to each step, but each step may be repeated until a satisfactory image is obtained. Each image is automatically labeled with the sequence step. The required measurement shapes that need to be correctly positioned will appear as appropriate.

If completing the sequence would mean that the number of images in the **active exam** would exceed the maximum allowed (50), an error message will appear allowing you either to open a new exam or to cancel the sequence.

Bladder Auto sequence

Choose **Start Sequence>Bladder Auto** to estimate bladder volume using **Bladder Volume** shape measurements (page 39). The scan settings used for this sequence are the default **Bladder** settings (page 71).

Step 1: Sagittal scan

- Perform a sagittal scan with the transducer positioned just superior to the symphysis pubis in the abdominal midline. When a good image of the maximum longitudinal-section of the bladder is obtained, stop scanning to freeze the image.
- If the image is not satisfactory, select  **Delete and Redo** to repeat the scan.
- If the image is satisfactory, tap the screen in the center of the bladder to place the Bladder Volume shape.
- Select  **Next Sequence Step**.

Step 2: Transverse scan

- Rotate the probe clockwise through 90° to perform a transverse scan.

- When a good image of the maximum cross-section of the bladder is obtained, stop scanning to freeze the image.
- Tap the screen in the center of the bladder to place the Bladder Volume shape. An estimate of bladder volume will appear in the bottom left corner of the screen.
- Select  **Next Sequence Step**.

Step 3: Sequence complete

- The calculated bladder volume in ml will appear on the final screen.
- Exit the sequence by selecting  **End Sequence**.

Bladder Calipers sequence

Choose **Start Sequence>Bladder Calipers** to estimate bladder volume using [caliper](#) measurements (page 38). The scan settings used for this sequence are the default [Bladder](#) settings (page 71).

Step 1: Sagittal scan

- Perform a sagittal scan with the transducer positioned just superior to the symphysis pubis in the abdominal midline. When a good image of the maximum longitudinal-section of the bladder is obtained, stop scanning to freeze the image.
- If the image is not satisfactory, select  **Delete and Redo** to repeat the scan.
- If the image is satisfactory, use the caliper to measure maximum bladder length.
- Select  **Next Sequence Step**.

Step 2: Transverse scan

- Rotate the probe clockwise through 90° to perform a transverse scan.
- When a good image of the maximum cross-section of the bladder is obtained, stop scanning to freeze the image.
- Two calipers appear. Use them to measure the maximum height and depth of the bladder.
- Select  **Next Sequence Step**.

Step 3: Sequence complete

- The calculated bladder volume in ml will appear on the final screen.
- Exit the sequence by selecting  **End Sequence**.

Grid and Rulers

Adjust viewing options for the touch screen through **Main Menu>Grid and Rulers** (Fig 6.9).

You can turn on a 1cm grid that covers the entire screen (Fig 6.10), or rulers marked in 1cm increments that appear on any side of the screen (Fig 6.11). Hide the grid and ruler markings by de-selecting the options.

You may also choose to hide all on-screen text by checking the **Hide Screen Text** option. The on-screen text will re-appear if you de-select the check box, or when you open a new exam.




Fig 6.9 Grid and Rulers




Fig 6.10 1cm on-screen grid




Fig 6.11 Left and bottom 1cm rulers

Image

Main Menu>Image allows you to view image information and to manipulate images (Fig 6.13).

- Choose **Image Information** (Fig 6.12) to find:
 - A description of **Image Properties**, such as [imaging mode](#), size and resolution.
 - A description of **Image Scan Settings**, such as [preset](#) name, [gain controls](#), and [depth](#).
 - The **Use Image Scan Settings** button, which loads the [scan settings](#) used to generate the displayed image. Save these settings using the [Edit Settings...](#) page (page 43).

Images are named by their text annotation.

- **Flip Orientation** flips the image orientation left to right (page 32).
- **Undo** and **Redo** apply only to image deletion. You can undo an image deletion by selecting **Undo**; you can re-delete the same image by selecting **Redo** or the  **Delete Image** icon on the Tool Bar.




Fig 6.13 Image




Fig 6.12 Image Information

Items on this menu can also be accessed via the [Image pop-up menu](#) (page 50).

Print Images...

Choose **Main Menu>Print Images...** when an exam is open to select images for printing (Fig 6.14). You can print [annotations](#) and [measurements](#) using the **Print Measurements** check box. You can control the printing of duplex images using the **Duplex Image Options** drop down. Duplex images refer to M-mode and PW Doppler images that have an associated B-mode image.

Select images using the check box in the bottom right corner of each thumbnail, and then press **Print**. If a printer is not connected when print jobs are queued, commands will be remembered until either **Signos RT** is connected to a printer; or the print jobs are canceled.

Press **Print Settings** set printing options and defaults (Fig 6.15).




Fig 6.14 Print Images

Print Screen

Main Menu>Print Screen will print the current image displayed on the screen, at the current magnification, with all annotations and measurements.

Preferences

Main Menu>Preferences gives options for the following settings (Fig 6.16).

- **Hide Tool Bar on Scan Start** will minimize the Tool Bar when a scan starts.
- **Maximum Scan Duration** will cause a scan to stop when maximum scan time has elapsed if not stopped manually. Adjustable from 0:30–5:00 minutes.
- **Fetal Biometry Charts** allows the chart used for fetal biometry measurements to be set. The following charts are available:
 - **ASUM** (Australasian Society for Ultrasound in Medicine “Policies and Statements D7 – Statement on Normal Ultrasonic Fetal Measurements. Revised May 2001”).
 - **BMUS** (British Medical Ultrasound Society “Ultrasound 2009;17(3):161-167 Fetal size and dating: charts recommended for clinical obstetric practice.”).
 - **Hadlock** (Hadlock,F., et al. “Estimating Fetal Age: Computer-Assisted Analysis of Multiple Fetal Growth Parameters.” Radiology, 152: 1984, pp 497-501; Hadlock,F., et al. “Fetal Crown-Rump Length: Re-evaluation of Relation to Menstrual Age (5-18 weeks) with High-Resolution, Real-Time Ultrasound.” Radiology, 182: February 1992, pp 501-505.)




Fig 6.15 Print Settings




Fig 6.16 Main Menu Preferences

Other Main Menu options

About...

Main Menu>About launches a window giving version information about the ultrasound application.

Home

Main Menu>Home takes you back to the [Home screen](#) (page 18).

7. SCANNING SHORTCUTS

Pop-up menus

A screen similar to Fig 7.1 will appear when the image is frozen during an [active exam](#). Tap different areas with the stylus to reveal pop-up menus.




Fig 7.1 Pop-up menus

- The *Exam pop-up menu* is opened by tapping the **Patient name** or **Exam identifier**; use this menu to toggle between open exams.
- The *Image pop-up menu* is opened by tapping the **Image name**; use this menu to toggle between images in the exam. This menu contains options also available in **Main Menu**>[Image](#) (page 48).
- The *Mode pop-up menu* is opened by tapping the **Mode identifier** use this menu to switch between different scan modes (pages 15 and 16).
- The *Settings pop-up menu* is accessed by tapping the **Preset name**; use this menu to [edit scan settings](#) or to [load a new preset](#) (pages 43 and 44).
- Tap the [Flip icon](#) to flip the image left to right.

Mode and Settings pop-up menus are not available when an exam is opened for review.

Scan button

Hold down the **Scan button** for 1 second when not scanning to switch [imaging modes](#) (page 15).

Scroll wheel

Click the **Scroll wheel** while actively scanning to adjust common scan settings (see page 17).

Touch screen

Panning and zooming

Use the stylus and touchscreen to pan and zoom around an image:

- Pan by pointing the stylus at your region of interest and dragging it to the center of the screen.
- In B-mode, zoom in (X- and Y-scales equally) by drawing small clockwise circles with the stylus. Scrollbars will appear at the right and bottom edges of the screen. Draw anticlockwise circles to zoom out.
- In M-mode or PW Doppler, the vertical and horizontal scales can zoom independently. Draw your clockwise circles near the scale you wish to zoom:
 - Near the left or right edges of the screen to zoom the vertical scale.
 - Near the top and bottom edges of the screen to zoom the horizontal scale.
 - In the middle of the screen to zoom vertical and horizontal scales together.

Keyboard gestures

Use **Preferences**>[Keyboard](#) to enable keyboard shortcuts (page 24)

8. MAINTENANCE AND TROUBLESHOOTING

Maintenance

Warning

- ✓ No modification of this equipment is allowed.
- ✓ This device contains no user-serviceable parts. Please return your device to your dealer or to Signostics for maintenance or repair.
- ✓ The internal battery is not to be replaced by the user – doing so may create fire or explosion hazard.

Signos RT requires no routine or periodic maintenance on either hardware or software components. When software updates are available for download, the user should install them.

Inspect the probe and display unit housings, connectors and cables prior to each use for cracks; if cracks are found, use of the device should be discontinued immediately. Contact Signostics Customer Support or your **Signos RT** dealer to arrange repair or replacement.

If **Signos RT** behaves in an abnormal manner or when damage is suspected use of the device should be discontinued immediately. Contact Signostics Customer Support or your **Signos RT** dealer to arrange repair or replacement.

Signos RT must be returned to Signostics Customer Support or your **Signos RT** dealer to arrange battery replacement. Recycle batteries where possible and abide by local disposal guidelines and directives.

Caution

- ✓ Ultrasound transducer crystals are fragile and are easily damaged if knocked, dropped or excessively vibrated.

If **Signos RT** is dropped or suffers a heavy knock

- Inspect the device. If you discover major damage do not use the device; report the fault to your **Signos RT** dealer for repair.
- Switch on and check normal function of all controls and carefully analyze **Signos RT** images for degraded quality.

If there are any concerns about image quality of the device, report the faults to your **Signos RT** dealer for inspection and/or repair.

Battery Maintenance

The **Signos RT** batteries cannot be removed by the user. Precautions need to be taken if the device is not used for long periods of time. Before storage, fully charge the device. Turn the device off when placing in storage. Take the device out of storage at least every six months and recharge the device.

Caution

- ✓ The **Signos RT** batteries should be charged every six months at a minimum, even if you are not using your device.

Transducer cleaning and disinfection

Warning

- ✓ Remove all particles and other matter from crevices and surfaces when cleaning the system and components.
- ✓ The device is supplied unsterile.

- ✓ Clean and disinfect the transducer between patients
- ✓ Before cleaning or disinfection turn **Signos RT** off and disconnect from the power supply
- ✓ Do not submerge the transducer, the display unit or the mains charger as electric shock could result. The ultrasound probe is IPX1 which allows water falling vertically onto the probe. The display is IPX0 and has no protection against ingress of water.

Caution

- ✓ Use only recommended disinfection methods
- ✓ Before cleaning or disinfecting the ultrasound probe and transducer, ensure the transducer is inserted into the probe securely. Failure to insert the transducer into the probe before cleaning may result in damage to either or both of the transducer and probe.
- ✓ Do not use abrasive cleaners, isopropyl alcohol or solvents such as benzene on any part of **Signos RT**.
- ✓ Do not heat sterilize any part of the **Signos RT** device.

Cleaning

Signostics recommends the use of a soft cloth, lightly dampened with warm water or a mild soap solution for cleaning the display unit, probe, and transducers. Wipe surfaces with the moist cloth, and towel dry with a soft clean cloth.

Where a personal, laptop, or tablet computer is used, the user should refer to the user manual for that component for cleaning requirements.

Disinfection

The following disinfecting means can be used with the **Signos RT**:

- Pi-Spray II by *Pharmaceutical Innovations* of New Jersey
- Sani-Cloth® Active by *PDI* of Flint, UK
- Tuffie 5 Wipes by *VernaCare* of Lancashire, UK
- Protex Disinfectant Wipes by *Parker Laboratories* of New Jersey
- Any Isopropyl alcohol impregnated wipes with $\leq 70\%$ w/w alcohol

Contact [Signostics](#) for additional disinfection options.

Caution

- ✓ Minimize application of alcohol based disinfectant to colored overmold materials. Long term use may result in material degradation.

The ultrasound probe and transducer may be disinfected using a wipe or spray method of disinfection. First clean the ultrasound probe as above. Then wipe or spray the probe with disinfectant, and leave for the disinfectant manufacturers recommended contact time. Remove any residue with a moistened soft cloth. Do not allow any solutions to dry on the probe or transducer.

The display unit may only be disinfected by the wipe method of disinfection. First clean the display unit as above. Then wipe the display unit with disinfectant, and leave for the disinfectant manufacturers recommended contact time. Remove any residue with a moistened soft cloth. Do not allow any solutions to dry on the display unit.

Warning

- ✓ After cleaning or disinfection examine the ultrasound probe, transducer, and display unit as appropriate for cracks or leaks, and if damage exists discontinue use of the system and contact [Signostics](#)

Troubleshooting

Table 8.1 Troubleshooting guide

Symptom	Possible cause	Remedy
The Signos RT does not turn on	Flat battery	Connect to the mains charger (P03198) and turn on.
The Signos RT displays a “Transducer Error” message when starting or during a scan	The transducer is not connected correctly to the probe	Unplug the transducer and reinsert.
The Signos RT displays a “Probe Error” message when starting a scan or during scanning	Probe fault. Low Probe battery.	Unplug the USB connector from the probe, wait 10s, and reconnect. Check the battery status and connect the display unit to the mains charger if the display or probe battery is low
The Signos RT status bar shows the Probe is missing 	Low Probe battery	Disconnect the Probe from the USB cable. Plug the mains charger directly into the Probe USB connector and let charge for 15 minutes. Plug back into the display USB cable.
The Signos RT battery fails to recharge after a period of heavy use. The battery status indicates one or both batteries are not recharging despite the system being connected to the mains charger (P03198)	Battery temperature exceeds 45°C	Allow device to cool for 10–15 minutes, then try charging the battery again.
Storage card error 	microSD card corrupt or not formatted	Re-format microSD card or insert new card.
Signos RT fails to connect to a host PC USB port	Host PC USB port is not USB 2.0 compliant.	Some desktop and tower PCs have USB ports on the front of the case that are not compliant with the USB 2.0 standard. Connect the Signos RT device to the host PC via the USB 2.0 port on the rear of the desktop or tower PC.
Signos RT becomes non-responsive or will not power up after recharge.	Unknown system error	Shutdown Signos RT and reboot. If unable to shutdown and reboot press the reset button adjacent to the microSD card slot using a thin object like a paperclip.
Any error message not resolved by the actions above or the action suggested on screen.	Various errors	Ensure the latest revision of software is installed on the device. Contact Signostics customer support for help with any unresolved errors.

Signostics Customer Support

If problems persist or recur, contact Signostics Customer Support.

Free call within Australia 1800 SIGNOS (1800 744 667)

For further details, see [Contact information](#) (page 8).

9. CONNECTIVITY AND ACCESSORIES

Accessories

The following accessories are supported by **Signos RT**:

- S3 Transducer (P03479)
- S3-5 Transducer (P03611)
- Sigviewer software (P03630)
- **Signostics** Power Supply (P03198)
- MicroSD Card (P02857)

System connectivity

Warning

- ✓ The user must not touch the any device connectors while in physical contact with the patient.
- ✓ Where **Signos RT** is connected to a personal computer via USB, the personal computer and connected peripherals shall be outside the patient environment and must be certified to the respective IEC standards and particular national deviations (IEC 60950 or IEC 60601-1:2nd and/or 3rd edition). Furthermore, any user connecting **Signos RT** to other equipment is configuring a medical system, and therefore is responsible for ensuring the system complies with IEC 60601-1: 3rd Edition Clause 16, and/or the relevant IEC 60601-1-1 collateral standard for use with IEC 60601-1:2nd Edition. Extension cords or Multiple Socket Outlets shall not be connected to the system when the PC is only IEC 60950 compliant a medical grade isolation transformer complying with IEC 60601-1 2nd Edition and/or IEC 60601-1 3rd Edition shall be used.

Contact [Signostics](#) for guidance on connecting **Signos RT** to external computers and peripherals (page 8).

Failure to comply with these guidelines may result in electric shock.

- ✓ Only connect accessories or items that are specified as being compatible with **Signos RT**. Contact [Signostics](#) for information on compatible accessories and systems.

There are two main points for system connectivity.

- MicroSD card slot for image storage and transfer.
- USB connection port for direct connection to a PC loaded with **SigViewer** software or a PictBridge printer.

In **Urgent Mode**, internal memory can temporarily store images. Images cannot be permanently saved to this memory.

The Ultrasound probe can be connected directly to a PC loaded with SigViewer or a tablet computer, and provides full functionality. Refer to the Sigviewer User Manual (P03629).

Power sources

During operation, **Signos RT** is electrically powered by an internal battery. The battery can be recharged by a mains-powered supply via the provided mains charger (P03198).

When **Signos RT** is connected to a Personal Computer (see System Connectivity, page 56) via a USB cable, the device will slowly recharge.

Internal battery

When powered by a fully charged battery, **Signos RT** will run for more than 12 hours with an average operating regime of 4×1 minute scans per hour.

From a fully discharged state, while **Signos RT** is in sleep mode, the battery recharges to 80% capacity when connected to a:

- Mains charger for 2-3 hours.
- A 500mA USB host port for more than 12 hours.

Mains charger

Warning

- ✓ Recharge **Signos RT** only with the mains charger provided (**Signostics** Power Supply P03198).
- ✓ Do not open or modify the **Signostics** Power Supply P03198 – Risk of electric shock
- ✓ Only connect to a mains supply rated at 100-240V and 50-60Hz.
- ✓ Do not use the device or mains charger if there are signs of damage.

The power supply can deliver a charging power to **Signos RT** of between 0–10W.

From a fully discharged state, full charge will take approximately 3-4 hours.

Note: Scanning whilst charging from mains is not recommended. **Signos RT** conforms to CISPR 11 RF Emissions Class A in this configuration; refer to "Electromagnetic compatibility (EMC)" in Section 11.

SigViewer software

Information can be transferred between **Signos RT** and a PC loaded with **SigViewer** software (P03630). This software is contained on an enclosed DVD, and allows the user to:

- Control the **Signos RT** ultrasound probe and view images.
- Transfer patient information between your PC and **Signos RT**.
- Transfer exam information, images, annotations and measurements to your PC. You may also edit exam information on your PC and transfer the updates back to **Signos RT**.
- View and store patient, exam and image information on your PC.
- Add and edit measurements and annotations on your PC.
- Magnify and display the **Signos RT** screen on your PC.
- Update **Signos RT** software.

For further information on **SigViewer** software, please refer to the **SigViewer** manual (P03629).

10. SPECIFICATIONS

System dimensions

Table 10.1 Dimensions of probe, display unit and screen

Component	Height	Width	Depth	Weight
Ultrasound probe	110mm (4.3 inches) 150mm with transducer (5.9 inches)	50–60mm (2-2.4 inches)	20–32mm (0.8-1.3 inches)	~116g (4.1 oz) ~200g with transducer (7.1 oz)
Display unit	115mm (4.5 inches)	68mm (2.7 inches)	20mm (0.8 inches)	~188g (6.6 oz)
Touch screen	60mm (2.4 inches)	45mm (1.8 inches)	n/a	n/a

Table 10.2 Dimensions of stylus and cable

Component	Length	Diameter
Stylus	95mm (3.7 inches)	3mm (0.12 inches)
Connecting cable	800mm (31.5 inches)	5mm (0.20 inches)

Temperature and humidity limits

Table 10.3 Operating, charging, transport and storage condition ranges for Signos RT

	Operating and Charging	Transport and Storage
Temperature (°C)	0 (32°F) – 40°C (104°F) 0°C (32°F)  40 °C (104°F)	-20 (-4°F) – 45°C (113°F) -20°C (-4°F)  45 °C (113°F)
Relative humidity (non-condensing)	20–80%	10–90%
Pressure	70–106 kPa (10.3-15.4 psi)	50–106 kPa (7.3-15.4 psi)

Mode of Operation

CONTINUOUS – At and below 30°C

NON-CONTINUOUS – Above 30°C (50% duty cycle up to 10 minutes scanning/10 minutes not scanning)

Warnings

- ✓ After storage at extreme temperatures, check the transducer surface temperature before applying to a patient. A cold or hot surface may burn a patient.

Caution

- ✓ Only operate, charge and store **Signos RT** within the approved environmental parameters.

Mains charger

Signostics Power Supply P03198

Input: 100-240V AC supply 50-60Hz

Output: 5.0V ±0.25V at 2.0A maximum current

Internal battery

Display unit

Signostics Li-ion Polymer Cell

Description: Rechargeable lithium polymer battery

Rating: 3.7V 2000mAh

Ultrasound probe

Signostics Li-ion Polymer Cell

Description: Rechargeable lithium polymer battery

Rating: 3.7V 1350mAh

MicroSD card

SanDisk microSD™ card.

Measurement accuracy

Good scan technique and accurate placement of measurement shapes is essential. Measurements reflect a physical property such as distance for interpretation by a physician. They do not account for variations in acoustic speed of the body, which can vary significantly.

Table 10.4 Accuracy of measurements

Measurement	Units	Range	Accuracy	Transducers	Modes
Distance					
Axial	cm	Full scale	±5% or ± 0.1cm	3 MHz 3-5 MHz	B-mode M-mode
Lateral	cm	Full scale	±5% or ± 0.3cm	3 MHz 3-5 MHz	B-mode
Circumference					
Ellipse	cm	Full scale	±5% or ± 0.4cm	3 MHz 3-5 MHz	B-mode
Polygon	cm	Full scale	±5% or ± 0.4cm	3 MHz 3-5 MHz	B-mode
Area					
Ellipse	cm ²	Full scale	±7% or ± 0.1cm ²	3 MHz 3-5 MHz	B-mode
Polygon	cm ²	Full scale	±7% or ± 0.1cm ²	3 MHz 3-5 MHz	B-mode
Volume¹					
Single Plane	ml	Full scale	±15% or ± 15ml	3 MHz 3-5 MHz	B-mode
Dual Plane	ml	Full scale	±15% or ± 15ml	3 MHz 3-5 MHz	B-mode
Time	ms	Full scale	±5% or ± 10ms	3 MHz 3-5 MHz	M-mode

¹ When used as per instructions scanning a bladder in a tissue equivalent phantom.

Degree of Protection Against Ingress of Water

The display unit (P03477) is IPX0.

The ultrasound probe (P03508), when the transducer (P03479 or P03611) is securely inserted, is IPX1 below the level of the USB connection.

The transducer lens can be submersed in liquid for testing purposes.




Fig 10.1 Water Protection

11. SAFETY

Ergonomics

Warning

- ✓ Avoid musculoskeletal strain with prolonged use of **Signos RT**. Follow guidelines below.

The **Signos RT** is intended for quick-look applications by qualified health professionals. It is not intended for continual use in radiology or other departments. If you need to use the device for a continual period, take the following precautions:

- Position yourself comfortably, either with a chair with appropriate lower back support or by sitting or standing upright.
- Minimize twisting, relax your shoulders and support your arm with a cushion
- Hold the transducer lightly, keep your wrist straight and minimize the pressure applied to the patient
- Take regular breaks

External materials

The parts of the device that come into patient contact are biocompatible as assessed in ISO10993-1.

Disposal

Warning

- ✓ Do not incinerate or discard the device in general waste at end of life. The lithium battery is a potential environmental and fire safety hazard.

Signos RT contains lithium-polymer batteries and the system should be disposed of in an environmentally responsible manner in compliance with federal and local regulations. It is recommended the system is taken to a recycling centre which specializes in the recycling and disposal of electronic equipment.

Where the device has been exposed to biologically hazardous material it is recommended the device be disposed of using biohazard containers and in compliance with federal and local regulations. It is recommended the system is taken to a waste centre which specializes in the disposal of biohazard waste.

Electrical safety

Warning

- ✓ **Signos RT** complies with the requirements of IEC 60601-1:2nd Edition and IEC 60601-1:3rd Edition. To avoid the risk of injury or electrical shock, comply with all safety instruction and warnings.
- ✓ The user must not touch the any device connectors while in physical contact with the patient.
- ✓ Where **Signos RT** is connected to a personal computer via USB, the personal computer and connected peripherals shall be outside the patient environment and must be certified to the respective IEC standards and particular national deviations (IEC 60950 or IEC 60601-1:2nd and/or 3rd edition). Furthermore, any user connecting **Signos RT** to other equipment is configuring a medical system, and therefore is responsible for ensuring the system complies with IEC 60601-1: 3rd Edition Clause 16, and/or the relevant IEC 60601-1-1 collateral standard for use with IEC 60601-1:2nd Edition. Extension cords or Multiple Socket Outlets shall not be connected to the system when the PC is only IEC 60950 compliant a medical grade isolation transformer complying with IEC 60601-1 2nd Edition and/or IEC 60601-1 3rd Edition shall be used.

Contact [Signostics](#) for guidance on connecting **Signos RT** to external computers and peripherals (page 8). Failure to comply with these guidelines may result in electric shock.

Electromagnetic compatibility (EMC)

Caution

- ✓ **Signos RT** contains sensitive components and circuits. Failure to observe proper static control procedures may result in damage to the system. Any faults should be reported to **Signostics** or your **Signos RT** dealer for repair.

Warning

- ✓ **Signos RT** complies with the Electromagnetic Compatibility requirements of AS/NZ CISPR 11:2004 and IEC 60601-1-2:2004. However, electronic and mobile communications equipment may transmit electromagnetic energy through air and there is no guarantee that interference will not occur in a particular installation or environment. Interference may result in artifacts, distortion, or degradation of the ultrasound image. If **Signos RT** is found to cause or respond to interference, try re-orienting **Signos RT** or the affected device, or increasing the separation distance between the devices. Contact **Signostics** for further information.

The **Signos RT** is intended for use in the electromagnetic environment specified below. The user of the **Signos RT** should assure that it is used in such an environment.

Table 11.1 Guidance and manufacturer's declaration: electromagnetic emissions

Emissions test	Compliance	Electromagnetic environment: guidance
RF emissions CISPR 11	Group 1	Signos RT uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B ¹	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

¹Class B applies to **Signos RT** in all but the following configurations: (a) Scanning while charging, and (b) scanning while display unit is connected to a PC (e.g. for remote display). **Signos RT** has Class A compliance in these configurations, meaning it is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. If **Signos RT** is found to cause or respond to interference follow the guidelines in the warning section above

Table 11.2 Guidance and manufacturer's declaration: electromagnetic immunity

The user should assure the **Signos RT** is used in the electromagnetic environment specified below. The use of accessories, transducers and cables other than those sold by Signostics can lead to more emissions or decrease the electromagnetic immunity of this system

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment: guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8kV air	±6 kV contact ±8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. When ±4kV contact discharge is applied to the device, the scan may be interrupted. Normal scan operation will be restored by clicking the Scan button .
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1kV line(s) to line(s) ± 2kV line(s) to earth	± 1kV differential mode ± 2kV common mode	Mains power quality should be that of a typical commercial or hospital environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment: guidance
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T ¹ (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	<5% U_T ¹ (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
^{2,3} Conducted RF IEC 61000-4-6	3 Vrms 150kHz – 80MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of Signos RT , including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter Recommended separation distance $d = 1.2\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80MHz–2.5 GHz	3 V/m	$d = 1.2\sqrt{P}$ 80MHz to 800MHz $d = 2.3\sqrt{P}$ 800MHz to 2.5GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separations distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ⁴ , should be less than the compliance level in each frequency range ⁵ . Interference may occur in the vicinity of equipment marked with the following symbol. 

Notes for

Table 11.2

¹ U_T is the AC mains voltage prior to application of the test level

² At 80MHz and 800 MHz, the higher frequency range applies

³ These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

⁴ Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which **Signos RT** is used exceeds the applicable RF compliance level above, **Signos RT** should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorientating or relocating **Signos RT**.

⁵ Over the frequency range 150kHz to 80MHz, field strengths should be less than 3V/m

Table 11.3 Separation Distances

Recommended separation distances between portable and mobile RF communications equipment and the Signos RT			
The Signos RT is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Signos RT can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Signos RT as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Mains charger

The external power supply complies with the requirements of IEC 60601-1.

Warning

- ✓ Do not open or modify the **Signostics** Power Supply P03198 – Risk of Electric Shock

Biological safety

The device can be disinfected for re-use on immunocompromized and non-immunocompromized patients after:

- Frequent contact with healthy skin (frequent low risk).
- Occasional contact with mucous membranes or contamination with particularly virulent or readily transmissible organisms (occasional intermediate risk).

Acoustic Output

Warning

- ✓ The ALARA principle (As Low As Reasonably Achievable) should be employed for all medical ultrasound exposure

Ultrasound imaging has been in regular use for over 20 years and has an excellent safety record. Even though there are no known risks of ultrasound imaging, it heats the tissues slightly and may produce small pockets of gas in body fluids or tissues (cavitation). The long-term effects of tissue heating and cavitation are not known.

Given potential does exist for bioeffects with ultrasound exposure, any exposure should be kept As Low As Reasonably Achievable. Scans should only be performed when there is a medical benefit and in the shortest

time commensurate with obtaining an adequate study. As the acoustic output of the device does not exceed an MI or TI of 1.0, visual display of MI and TI values is not required. A generic ALARA education program is supplied with your **Signos RT** (see enclosed ISBN 1-93004 7-71-1, *Medical Ultrasound Safety*).

Table 11.4 Acoustic output and uncertainties for Signos RT

	$I_{SPTA,3}$ [mw/cm ²]	TI type	TI Value	MI	$I_{PA,3}@MI_{max}$ [W/cm ²]
P03479 3MHz Output	97.9	TIB	0.653	0.599	60.9
P03479 3MHz Uncertainties	±13.5%		±13.5%	±6.75%	±13.5%
P03611 3-5MHz Output	162	TIB	0.684	0.670	97.2
P03611 3-5MHz Uncertainties	±13.5%		±13.5%	±6.75%	±13.5%

$I_{SPTA,3}$ Derated spatial peak, temporal average intensity (milliwatts/cm²)

TI Thermal index

MI Mechanical index

$I_{PA,3}@MI_{max}$ Derated spatial-peak pulse-average intensity (watts/cm²)

Uncertainties are the cumulative effect of measurement uncertainties.

REFERENCES

- IEC 60601-1:1988 (2nd edition), Medical Electrical Equipment-Part 1. *General Requirements for Safety*.
- IEC 60601-1:2005 (3rd edition), Medical Electrical Equipment-Part 1. *General Requirements for Basic Safety and Essential Performance*.
- IEC 60601-1-1:2000 (2nd edition), Medical Electrical Equipment-Part 1. *General Requirements for Safety – Section 1.1. Collateral Standard. Safety Requirements for Medical Electrical Systems*.
- IEC 60601-1-2:2001 (2nd edition), Medical Electrical Equipment. *General Requirements for Safety. Collateral Standard. Electromagnetic compatibility requirements and tests*.
- IEC 60601-1-2:2007 (3rd edition), Medical Electrical Equipment. *General Requirements for Basic Safety and Essential Performance. Collateral Standard. Electromagnetic compatibility requirements and tests*.
- IEC 60601-2-37:2007 (2nd edition), *Particular requirements for the basic safety and essential performance of ultrasonic medical diagnostic and monitoring equipment*.
- ISO 10993-1: 2003, Biological Evaluation of medical devices
- ISBN 1-93004 7-71-1, *Medical Ultrasound Safety*.

Glossary

Active exam	The open exam that is currently recording images. Only active exams may have images added. One active exam is open at all times.
Active settings	The ultrasound scan settings that Signos RT is currently using.
Annotation	Any text or voice notes associated with a particular image.
Current exam	The exam whose images are displayed on the touch screen. The open exam may be the active exam or a previously recorded exam that has been opened for review.
Exam	A unique interaction between ultrasound user and patient. Up to 50 images may be saved to each exam. Exams are identified by patient name and/or by date and time of scan.
Patient	A single person, with unique details such as name, gender, date of birth and EMR (electronic medical record).
Scan settings	A combination of settings for ultrasound scanning, which includes parameters such as near field gain, far field gain, scan depth, gain rate and gain delay. A particular combination of scan settings may be saved for future use (see Preset).
Sequence	A pre-defined series of scans that can be used to efficiently progress between scan types in frequently-used multi-scan diagnoses or procedures.
Measurement	Any distance or area measurement shape added to an individual image. Shapes can be calipers (linear distance); polygon (area and perimeter); ellipse (area and circumference); or bladder volume (volume). Distance (cm), area (cm ²) and volume (ml) can be measured on B-mode images; heart rate (beats per minute) or time (ms) vs. distance (cm) can be measured on M-mode images.
Preset	A pre-saved and named combination of scan settings. Included system-defined presets may not be edited.
Storage card	The microSD card that can be inserted into Signos RT for saving ultrasound images and patient information.
Title bar	The top section of the Signos RT screen, which contains the screen title, icon tray and application icon.
Tool Bar	The bottom section of the Signos RT screen when Signos RT is in its ultrasound application. It contains icons such as Main Menu , Exam management and Measurement .

Key to icons

Application icons

-  Home
-  Ultrasound
-  Printer
-  Preferences
-  Flashlight

Title bar

-  Probe status ( )
-  Battery Info (    )
-  Storage Card Info (   )
-  Volume Info ()
-  Print Manager ( )
-  Back
-  Close

Sequence controls

-  Previous Sequence Step
-  Delete and Redo
-  Next Sequence Step
-  End Sequence
-  Repeat Sequence

Tool Bar

-  Main Menu
-  Exam Management
-  Movie Clip Icon
-  Measurement
-  Image Annotations
-  Delete Image
-  Hide Tool Bar

Movie controls

-  Play movie
-  Save single frame
-  Save movie
-  Exit scan review

Other icons

-  Scan Review mode
-  Voice recording

APPENDIX 1: DEFINITIONS AND PARAMETERS OF SCAN SETTINGS

Operating Frequency	Refers to the transmit frequency of the probe.
Scan Depth	Refers to the maximum depth of tissue displayed on the ultrasound image. Adjustable from 1–18cm.
Gain	Controls the brightness of the ultrasound image. Increasing gain amplifies returning echoes to make the image brighter, while decreasing gain makes the image darker. Adjustable from 0–80dB.
Far Field Gain	Controls the brightness at the bottom (deeper part) of the ultrasound image. Increasing gain amplifies returning echoes to make the image brighter, while decreasing gain makes the image darker. Adjustable from 0–80dB.
Gain Delay	Refers to the initial depth/delay before amplification of the returning signal is applied. It is sometimes useful not to amplify superficial echoes when you are interested in deeper anatomy, e.g., an obese patient. Adjustable from 0–18cm.
Gain Rate	Controls the rate of amplification increase to counter the effects of the attenuation (weakening) of the ultrasound signal with depth. Adjustable from 0–1.5dB/cm/MHz.
Transfer Curve	Adjusts the grayscale characteristics of the display. Adjustable from -3 (maximum grayscale expression) to 3 (maximum contrast).
Bidirectional Scanning	Enables scanning on transducer sweeps in forward and reverse directions, increasing the frame rate from 8 frames per second to 16 frames per second.
Image Post Processing	Enables an image post processing filter that reduces speckle noise artifacts.
Initial Scan Mode	Sets the scan mode when a preset is first selected. Note there is an option of None meaning that the scan mode is not changed.
Initial Image Flipped State	Sets the flipped state for scanning when a preset is first selected. Note there is an option of Leave As Is meaning that the flipped state is not changed.
M-mode scan settings	
Sweep Speed	Adjusts the sweep speed for M-mode scanning. Adjustable from 20-250 pixels per second, e.g., the screen is 240 pixels wide, so setting the sweep speed to the default of 80 pixels per second will result in the scan taking 3 seconds to fill the screen.
PW Doppler scan settings	
PW Doppler Gain	Adjusts the amplification of the PW Doppler signal.
PW Doppler Depth	Adjusts the depth at which PW Doppler measures velocity.
PW Doppler Gate Length	Adjusts the gate length at which PW Doppler measure velocity. For example, if the depth is set to 5mm and gate length to 3mm, then all the velocities between 5mm and 8mm will be detected and displayed.
Doppler Baseline Offset	

Adjusts the offset applied to the PW Doppler signal before display. For example, if a Doppler signal is aliasing, the Doppler Baseline Offset can be adjusted to increase the dynamic range available in one direction (while reducing the dynamic range in the other direction).

PW Doppler Audio

Enabled audio output of the PW Doppler signal.

APPENDIX 2: SYSTEM-DEFINED PRESET PARAMETERS

General Medical and Obstetric

Table A.0.1 Scan settings for medical and obstetric system-defined presets

Preset	Obs	Depth (cm)	Near field gain (dB)	Far field gain (dB)	Gain delay (cm)	Gain rate (dB/cm/MHz)	Transfer curve	Doppler Gain
S3 Default	✓	18	15	60	0	0.45	0	50
Abdomen		14	10	60	0	0.45	0	
Aorta		10	10	60	0	0.50	0	50
Bladder	✓	14	5	50	0	0.50	0	
Bladder Auto	✓	14	5	50	0	0.50	0	
Lung		8	0	45	0	0.45	0	
Obstetric early	✓	12	15	50	0	0.45	0	
Obstetric	✓	17	15	60	0	0.50	0	
Fetal Heart *	✓	17	15	60	0	0.50	0	
Cardiac*		16	10	60	0	0.45	-2	
S3-5 Default	✓	18	15	60	0	0.45	0	50
Large Vessel (3-5)		5	5	60	0	0.50	0	50

* Uses bidirectional scanning.

Table A.0.2 Scan settings for physiotherapy system-defined presets

Preset	Depth (cm)	Near field gain (dB)	Far field gain (dB)	Gain delay (cm)	Gain rate (dB/cm/MHz)	Transfer curve
S3 Default	18	15	60	0	0.45	0
Bladder Auto	14	5	50	0	0.50	0
TA	8	5	50	0	0.45	0
Pelvic Floor	14	10	45	0	0.45	0
S3-5 Default	18	15	60	0	0.45	0

Veterinary

Table A.0.3 Scan settings for veterinary system-defined presets

Preset	Depth (cm)	Near field gain (dB)	Far field gain (dB)	Gain delay (cm)	Gain rate (dB/cm/MHz)	Transfer curve	Doppler Gain
S3 Default	18	15	60	0	0.45	0	50
Bladder Auto	14	5	50	0	0.50	0	
Dog-large	13	15	50	0	0.45	0	
Dog-small	8	5	50	0	0.45	0	
Lung	8	0	45	0	0.45	0	
S3-5 Default	18	15	60	0	0.45	0	50
Large Vessel (3-5)	5	5	60	0	0.50	0	50

INDEX

A

about
 hardware 31
Signos RT 1
 software 50
 User Manual See User Manual
 accessories 3, 56
 acoustic output 64
 active
 exam See exams
 settings See scan settings
 ALARA principle 64
 alcohol 54
 annotations See **Image Annotations**
 applications See intended uses
 area See **Measurement**
Audio screen 23, 30

B

Back button 12, 19, 26, 32
Back icon 19, 21
 backlight dimming See Power
 battery 24, 26, 59
 disposal 61
Info window 29
 recharge 14, 57
 replacement 53
Bidirectional Scanning 68
 biological safety See safety
Bladder Auto sequence 46
Bladder Calipers sequence 47
 bladder volume 39, See also **Measurement**
 B-Mode Imaging 16
 brightness See display

C

calipers 38, See also **Measurement**
 card See microSD card
 cavitation 64
 charger See mains charger, See mains charger
 check box 19
 circumference See **Measurement**
Close icon 19, 21
Colors screen 23
 connecting cable 11, 58
 connectivity 56
 contact information 8, 55
 contraindications See safety
 controls See display unit controls
 copyright 8
 current exam See exams
 Customer Support 55

D

Date and Time screen 24

adjustment 19
Default Presets 42
Delete Image 41
 depth 17
Diagnostic Area 42
 dimensions 58
 dimming See Power
 display
 (brightness) screen 24
 unit 11, 58
components 12
 controls 12, 32
 disposal 61
 distance See **Measurement**
 Doppler Audio 69
 Doppler baseline 17
 Doppler gain 17
 Doppler gate depth 17
 Doppler gate width 17
 Doppler pulse repetition frequency 17
 drop-down box 20
 dynamic adjustment 17, 33

E

editable exam See exams
 electrical safety See safety
 electromagnetic compatibility 62
 ellipse 39, See also **Measurement**
 Exam identifier 32
 exams 34, 66
 active 35, 66
 current 35, 66
 delete 36
 editable 35, 49
 empty 41
Management screen 35
 open mode 49
 pop-up menu See pop-up menus
 read only 35, 49
 scroller See **Image/Exam Scroller**
 external materials 61

F

file transfer See USB
 Flashlight 18, 27
 Flip icon 32
 flip image See orientation

G

gain 17
 change 17
 delay 68
 rate 68

H

headphones 12, 23, 30

Hide Tool Bar	41
Home screen	18, 22, 50
I	
icons	18, 21, 28, 67
idle timeout	See Security screen
image	
information	48
pop-up menu	See pop-up menus
properties	48
quality	See scan settings
Image Annotations	40, 66
manage	44
text	40
voice	41
Image name	32
Image/Exam Scroller	37
imaging mode	15
intended uses	4, 5
internal	
battery	See battery
memory	14, 25, 56
IPX Rating	60
K	
keyboard	20, See also text box
gestures	24
screen	24
L	
Labels screen	25
language	25
latex	17
list box	20
Locale screen	24, 25
Log in	30
screen	13
Log out	30
loud speaker	23, 30
M	
Main Menu	34, 42
mains charger	56, 64
input	58
maintenance	53
Measurement	66
accuracy	59
B-mode	38
calipers	38
circumference	39
delete	38
distance	38
ellipse	39
heart rate	40
M-mode	40
numbering	38
polygon	39
PW Doppler mode	40
time	40
menus	19

pop-up	See pop-up menus
message box	19
microphone	41
microSD card	13, 29, 56
card label	25
card name	13, 25
eject	13, 30
insert	13
password	13, 27
recovering lost information	13
security	13, 14
Security Policy	13
user name	13, 27
mirror image	See orientation
M-mode Imaging	16
Mode identifier	32
Mode pop-up menu	See pop-up menus
Movie controls	37
musculoskeletal strain	61
mute	23, 30

N

navigate	See touch screen
----------------	------------------

O

operating frequency	68
orientation	32

P

Pan/Zoom Controls	37
password	
timeout	See Security screen
patients	15, 34, 35, 66
edit	36
polygon	39, See also Measurement
pop-up menus	
Exam	51
Image	49, 51
Mode	51
Settings	51
Power	
button	12, 26, 31, 33
screen	24, 26, 30
sources	56, 57

Preferences

Main Menu	49
System Menu	23, 30
preset	See scan preset
Preset name	32
asterisk	43
Imported*	44
print	27
Print Manager	See printing
Print Screen	See printing
printing	18, 22, 30
Print Screen	49
Probe status window	28
PW Doppler Imaging	16

R

read only exam	See exams
recommended gel	17
recording	41
references	65
Remove Card	30
reset button	12
Revert to Factory Presets.....	42
review exam	36
rulers	See Grid and Rulers

S

safety	61
biological	64
contraindications.....	4
electrical	61
saving a movie	37
scan	
button.....	12, 26, 35
shortcuts	51
depth	17
preset	66
create new	44
default	42
edit	44
factory	42
load.....	44
save	43
settings.....	42, 48
active.....	66
definitions	66
edit	43
pop-up menu	See pop-up menus
save	43
start	15, 33
screen title	21
screens.....	19
scroll wheel	33
shortcuts.....	17, 33, 51
Scroll wheel	12, 19, 26
SD card.....	See microSD card
Security Policy	13
Security screen	26, See also microSD card
sequences	46, 66
serial number	
hardware	31
probe	28
software	See version
service	See maintenance
setting up	15
settings	See scan settings
Show Password Text	13
shutdown	30, 31
SigViewer software	56, 57
sleep mode.....	See Power screen
slider.....	20
software	See also SigViewer software
update	57
version	See version
solvents.....	54

sterile sheaths.....	17
storage card	66, See also microSD card
Storage Card Info window.....	29
stylus	19, 58
shortcuts	24
support	55
sweep speed	68
System Menu	18, 22, 32

T

text box	20
time and date	See Date and Time screen
Title bar.....	21, 66
Tool Bar	33, 66
hide	41
unhide	41
un-hide.....	33
touch screen	19, 58
Calibration screen.....	26
navigating.....	19, 32
trademarks.....	8
transducer	1, 11
changing	14
Transfer Curve	68
transfer files	See USB
turn off	31
turn on	12

U

ultrasound	
application	18, 22
gel.....	17
probe.....	11, 58
risks	64
screen.....	18
Unidentified patient	34, 35
change patient details	36
Updating Repository	13
Urgent Mode	14
save.....	14
USB	
connection	56, 57
while scanning	56, 61
screen.....	27
User Manual.....	2
conventions.....	2
version	3
user name	See microSD card
Username & Password screen	27

V

version.....	50
User Manual	See User Manual
Volume Info window	30

W

windows	19
---------------	----

X

X-axis range	68
--------------------	----

Z

zoomSee **Pan/Zoom Controls**