

Product Note P9



ModuLight User Guide



The ModuLight-module is a programmable light source that has been designed to investigate photo-electrical devices, such as solar cells. The add-on module will operate in combination with Ivium potentiostats, through the peripheral I/O-port.

The ModuLight, by default, contains 7 LEDs with wavelengths ranging from 460-740 nm. On request LEDs can be exchanged for others from the same product range (see table below).

During operation an LED can be programmatically selected. The sinewave generator of the potentiostat can then be used to modulate the intensity of the LED with a frequency of 10 μ Hz-2MHz. The extensive Solar cell applications that are included in the Ivium software allow a full characterization of the solar cell. The functionality includes E/I curves as function of the light intensity, IMVS/IMPS, and solar cell modelling resulting in all characteristic values of the studied object.

1. Modulight specifications

Mode	Color	Wavelength [nm]	Power/Intensity [mW] / [lm]	Digital Code
1	Cool White	6500 K	250 lm	IIO
2	Blue	460	40 lm	IOI
3	Green	523	200 lm	IOO
4	Amber	590	105 lm	OII
5	Red	623	160 lm	OIO
6	Deep Red	660	900 mW	OOI
7	Far Red	740	705 mW	OOO
	<i>Cool White</i>	<i>5500 K</i>	<i>250 lm</i>	<i>Optional</i>
	<i>Violet</i>	<i>405</i>	<i>1400 mW</i>	<i>Optional</i>
	<i>UV</i>	<i>365</i>	<i>1470 mW</i>	<i>Optional</i>
	<i>IR</i>	<i>940</i>	<i>1150 mW</i>	<i>Optional</i>
	<i>IR</i>	<i>850</i>	<i>800 mW</i>	<i>Optional</i>
	<i>Dental Blue*</i>	<i>460</i>	<i>1100 mW</i>	<i>Optional</i>

** Dental Blue LED emitter provides superior radiometric power in the wavelength range specifically required for dental curing light applications.*

Table 1 Wavelength specification, specification subject to change as LED technology progresses.

Light intensity	can be modulated with the FRAoutput of the build-in sinewave generator of the IviumStat/CompactStat from 10 μ Hz to 2 MHz.
Bias resolution	16 bit, 0.0015%
Wavelength	set programmatically 460-740 nm in 7 steps
Power	0 - 1430 mW (depending on LED)
Bandwidth	0 => 2 MHz
Light aperture	circular, d = 34mm, 9.08cm ²
Power requirements	external adapter: 100-240 V, 45-65 Hz at DC-connector: 5 V, 1A
Size	w x d x h = 12 x 13 x 2.5 cm
Weight	0.5 kg
Interfacing/connectivity	DB37, connects to the potentiostat peripheral port D-Sub 15pins
Use	only i.c.w. Ivium potentiostats

2. Installation

Connect the ModuLight to the potentiostat using the cable that is included with the instrument. Insert the connection cable into the SD37connector at the back of the Modulight and the other end into the peripheral port connector of your potentiostat:

- The SD37-SD37 cable connects the ModuLight to the IviumStat and CompactStat;
- The SD37-SD15 cable connects the ModuLight to the Vertex.s and sModule.

The HD15 connector at the front of the Modulight is intended for the connection of the light-intensity meter. This Ivium product will be available soon.

Insert the external 5V power adapter. The ModuLight is now ready for operation, using the IviumSoft-software. The test object (solar cell) can be connected to the cell cable and illuminated with the ModuLight.


The operation of the ModuLight is independent of its position, i.e. the ModuLight can be placed under, on top, to the side, etc. of the test object. At the bottom of the ModuLight a screw socket is available to put the ModuLight on a common camera stand.

3. Direct operation

The ModuLight can emit a modulated light flux of 6 different wavelengths and 1 wavelength spread (white), using 7 different LEDs as light source. The present specifications of each wavelength are given in table 1 (above).


Wavelength Selection

The wavelength/LED can be selected by using the 3 digital outputs of the peripheral port. The digital codes are given in table 1. (Note that these codes are according to the electronic standards, i.e. counting from right to left). The digital outputs can be accessed in the IviumSoft from the "Direct" tab in the "Extern" tab below. As an example, to choose the color "white", the digital code is "IIO". This means that "Dig out 3" and "Dig out 2" need to be checked and "Dig out 1" is unchecked.




DC Intensity

The DC intensity of the light is controlled via the Analog output¹. This output has a range of 0 – 4V, which corresponds to a light output of 0–100% of the full output power. For example enter "1V" into the ch1 field and click on "ch1". This will result in a light output of the chosen colour corresponding with 25% of full-drive current².



AC Modulation

The AC modulation of the light is controlled via the AC-output of the peripheral port. The AC-output is accessed in IviumSoft from the "AC" tab which is found under the Direct control menu tab.



The rate of modulation is controlled by the Frequency setting. An AC amplitude setting of 0 – 1V corresponds to a modulation of 0 –100% of the full output power³, which is added to the selected dc power.

Examples

1-

Set DC Voltage (Analog output ch1) = 2V; click "Ch1"

Set AC Voltage (AC signal Amplitude) = 0.500V

Set Frequency (AC signal Frequency) = 2Hz; click "Apply"

Now the light will start modulating between 100% intensity and 0% intensity at the chosen frequency of 2 Hz. The 2V DC sets the light amplitude at 50% of maximum drive current³, then 0.5V AC adds a modulation intensity of 50% of full output current³, and so the combined signal oscillates between full output power and zero.

¹ The DC intensity/Analog output1 can also be controlled directly from within some methods (Mixed Mode) and as a Direct command from Batch Mode.

² See paragraph 7. Ivium Modulight Photometry

³ See paragraph 7. Ivium Modulight Photometry

2-

Set DC Voltage = 3.6V (90% of max power)
Set AC Voltage = 0.1V (10% of max power)
Set Frequency (AC signal Frequency) = 100Hz

This produces a modulation between 80% of full-on and 100% of full on with a dc offset power at 90% of full-on, at a frequency of 100Hz⁴.

3- Overdriving to non-linearity (not recommended):

DC Voltage = 4V (100% of max power)
AC Voltage = 0.05V (5% of max power)

This produces a modulation between 95% intensity and 100% intensity. Because more than 100% intensity is not possible, the first half of the sine wave will modulate between 95% and 100%, the second half of the sine wave will yield 100%. This could produce a non-linear modulation and is not recommended.

Note1: When controlling the ModuLight from the direct mode, it should be taken into account that communication between PC and Ivium device only takes place once per second. This means that there may be a delay of up to one second before a change takes effect.

Note2: The peripheral port signal output depends on the Ivium instrument type. I.e. when the "Plus-module" option is checked, the AC-output signal is a factor of 5.33 lower than the set value; when the "E extended range" option is checked, the AC-output signal is a factor of 2 lower than the set value. This means that in these cases the amplitude of the light intensity variation is equally less. If this is an issue, please contact Ivium Technologies for a solution.

Note3: Due to the non-linear nature of the LED light output as a function of the driving current⁴, the DC voltage control 0 - 4V also shows a non-linear control function.

⁴See paragraph 7. Ivium ModuLight Photometry

4. Running a Method

4.1 Single method

A single method can be run by first setting the desired parameters of the ModuLight, such as wavelength and intensity, in the "Direct" mode tab (see paragraph 3). Then in the "Method" mode tab, a method can be run. For example, to measure an E/I curve of a solar cell, set the desired wavelength (Digital outputs) and intensity (DC Voltage, Ch1), then the E/I curve can be measured using the linear sweep method⁵.

To run an impedance measurement, the ModuLight parameters such as wavelength and DC intensity need to be set first in the "Direct" mode tab. In the "Method" mode tab select the desired impedance technique (for example "Constant E"). Then activate the "Advanced" method, now the "MeasConfig" is available in the parameter list. The "MeasConfig" parameter determines how and which signals are recorded and shown in the graph (see also paragraph 5).

4.2 Automated methods

Batch mode

It is possible to run several methods successively in the "BatchMode". In between methods it is possible to change the wavelength (light colour) and DC intensity of the ModuLight. To do this add a "DirectCommand" line in the appropriate place in the batch program. When this line is selected in the Line-properties, the wavelength can be chosen by ticking the "SetDigOut" box. The boxes for the individual digital outputs become available and by ticking the relevant boxes the wavelength can be chosen. In the same way ticking the "SetDAC" box will make the analog outputs available. Inserting the desired value in the "DAC 1" field (=Analog output 1, ch1) will set the intensity of the light.

Mixed mode

In the transient technique "MixedMode" a sequence of operations can be programmed in the method parameter "Stages" (pop-up). In the "Properties for Level", at the bottom of the list, the following parameters allow operation of the ModuLight:

- AnOut1: when ticked, a value can be entered from 1 - 4 V which determines the DC intensity of the light from 0 - 100% (in the same way as in the direct mode operation, see paragraph 3).
- Digouts: when ticked, an integer value can be entered to select the wavelength. This integer corresponds to an 8 bit conversion for the digital outputs, i.e. 0 = 000 = Far Red, 6 = 110 = Cool White.
- The AC modulation of the light can be controlled by ticking the "Record ac" box. In the main method parameters the modulation Frequency (Hz) and AC amplitude can be set (amplitude: 0 - 1V corresponding with 0 - 100% intensity; in the same way as in the direct mode operation, see paragraph 3).

⁵ Modulation cannot be applied while performing a DC sweep – only during EIS. If an AC modulation was selected under the Direct Mode tab, it will simply switch off during the DC sweep. (However, note that sweep measurements incorporating pulsed modulation of the light amplitude can be achieved using Mixed Mode by segmenting the sweep into levels and having different values of the analogue output in each sweep level.)

5. Measurement/signal configuration (AC transfer function)

Standard electrochemical impedance is determined from measurements of the AC current and AC potential at the instrument's cell connectors. A standard impedance (Z) sweep is derived from E/I, where the Y-input/signal is E (voltage) and the X-input/signal is I_{we} (current at the working electrode). In the table 2 (below) this situation is given in line 0: MeasConfig = standard.

However, with Ivium potentiostats/IviumSoft it is possible to select other signals for the X- and Y-inputs. In the Impedance techniques the advanced method parameter "MeasConfig" is available. This parameter allows various signals to be used for the X- and Y-inputs. These signals can then be plotted in the impedance Result graph as if these were I_{we} and E respectively.

Table 2 below shows alternative choices which are accessible in IviumSoft from the "MeasConfig" parameter in the Impedance techniques (Advanced mode only). For the Modulight the relevant choices are lines 0, 1 and 2 (BOLD), corresponding with standard, INT_{ac} I and INT_{ac} E.

Table 2

	MeasConfig	X	Y	remarks
0	standard	I we	E	
1	INT_{ac} I	ac intern	E	Internal DSG not applied
2	INT_{ac} E	I we	ac intern	Internal DSG not applied
3	EXT _{ac} EI	X periph.	Y periph.	
4	EXT _{ac} I	X periph.	E	
5	EXT _{ac} E	I we	Y periph.	
6	EXT _I INT _E	X periph.	ac intern	
7	EXT _E INT _I	ac intern	Y periph.	
8	DirectE	I we	E	CE as reference, instead of (RE-S)
9	DirectE_INT _I	ac intern	E	CE as reference, instead of (RE-S)
10	DirectE_EXT _I	X periph.	E	
11	BiStat	I we	I we2	
12	EXT _I BiStat	X periph.	I we2	

In the table:

- I_{we} and E: Are generated or measured as usual via the potentiostat/cell connector;
- AC intern: Is the internally generated (FRA) ac signal applied via the peripheral port connector (=AC-out in the pin-out in the instrument manual). It drives the ModuLight AC intensity.
- Internal DSG not applied: Signifies that the Internal Direct Signal Generator is not applied, so no AC perturbation is carried through to the E and I_{we} signals; these are DC signals only.

When measuring impedance the IviumSoft always defines this as Y/X so, for standard impedance (Z), this would be E/I_{we} in the table above.

When other signals are used for Y and X the data is still displayed as impedance (Z=E/I) by the software however it is actually a transfer function which is defined by the MeasConfig selection and the user must interpret the data.

To invert the transfer function to X/Y, select Admittance (Y) instead of Impedance (Z) in the software.

Example 1: Solarcell IMPS - Intensity Modulated Photocurrent Spectroscopy

The IMPS technique measures the transfer function between modulated light intensity and the resulting AC current generated by the cell⁶. To run this technique:

- Select from the method tree: Impedance> Constant E.
- Activate advanced method parameters
- Set parameter "MeasConfig" to "INT_ac E"

With this parameter-setting, the internal applied AC signal is used as Y (light-intensity) and the photo current is used as X:

$$\text{IMPS} = Y/X = \text{light intensity}/\text{photo-Current} = \text{AC intern}/I_{we}.$$

This potentiostat's dc potential may be set to short circuit conditions ($E=0$), but it may also be made at, for example, the maximum power point or E_{oc} ($I=0$).

Note that the potentiostat's DC potential may be set to short circuit conditions ($E=0$), but it may also be operated at other settings, for example maximum power point or E_{oc} ($I=0$).

Example 2: Solarcell IMVS - Intensity Modulated Photovoltage Spectroscopy

The IMVS technique measures the transfer function between modulated light intensity and the resulting AC voltage generated by the cell³. To run this technique:

- Select from the method tree: Impedance> Constant I.
- Activate advanced method parameters
- Set parameter "MeasConfig" to "INT_ac I"

With this parameter-setting, the internal applied AC signal is used as Y (light-intensity) and the photo potential is used as X:

$$\text{IMPS} = Y/X = \text{photo potential}/\text{light intensity} = E/\text{AC intern}$$

Note that the galvanostat current may be set to OCP conditions ($I=0$), but it may also be made at, for example, the maximum power point.


⁶ There are no universal standards for display of IMPS and IMVS. If the inverse ratio is required then Admittance (Y) should be selected for display within IviumSoft.

6. Experimental example


Experiments were conducted with a 1cm² Dye Sensitized SolarCell (DSSC). The solarcell was connected to an Ivium CompactStat, with WE/S to the positive pole, and CE/RE to the negative pole (a negative current means that the solarcell is producing power).

As light source, the Ivium ModuLight module was used with the 635nm setting, connected to the CompactStat peripheral port's AnalogOut1 for DC intensity and ACout for AC intensity modulation.


- The DC LSV scans were done by first setting the illumination level with AnalogOut1 (in the Direct Mode) and thereafter performing a scan.
- The DC experiments with pulsed light were done with the Mixed Mode technique (the Anout1 parameter of the 2nd level was used to set the desired light intensity).
- The EIS scans were done by first setting the bias-illumination level with AnalogOut1 (in the Direct Mode tab), and thereafter performing an EIS scan with "MeasConfig" parameter to "Int_AC_E" for potentiostatic scans and "Int_AC_I" for galvanostatic scans. The applied frequency range was 10kHz to 0.1Hz.




E/I curves for DSSC at various light intensities: 15lm, 30lm, 60lm




Screenshot of the solar cell modelling circuit, using the 15lm intensity curve




Screenshot of the solarcell modelling circuit, using the 60lm intensity curve



Open Cell Potential of DSSC for variable intensity light pulses



Current response of DSSC at Estat 0V, for variable intensity light pulses



IMVS photo-electric-impedance of a DSSC at OCP, at various light intensities, left to right: 15lm, 18lm, 23lm, 30lm




Photo-electric-impedance of a DSSC at various constant currents, 0.1Hz to 1kHz, at fixed light intensities of 60lm, left to right: 0mA (OCP), 1mA, 3 mA




Photo-electric-admittance of DSSC at various DC bias, at fixed light intensity of 60lm, left to right: 0.7V, 0.6V, 0.0V (IMPS)




Photo-electric-admittance of DSSC at 600mV DC, at various light intensities, left to right: 15lm, 30lm, 60lm

7. Ivium ModuLight Photometry

Luminous Flux

The applied LED series is adapted to the so-called luminous efficiency function. This function describes the average sensitivity of the human eye to light at different wavelengths. This means that to obtain a relatively equal light output at different wavelengths the radiant power will vary.

The relationship between luminous and radiant intensity can be described as:


$$I_V = 683 \times I \times \bar{Y}(\lambda)$$

I_V Luminous Intensity in Lm/sr

I Radiant Intensity in W/sr

$\bar{Y}(\lambda)$ Luminous Function

sr Steradium



Luminous efficiency function

Luminous Intensity

In the graphs in paragraph 8 and 9, the relative light output (%) as a function of the driving current is given for each of the LEDs used in the ModuLight. 0 mA corresponds to 0% light, 1000mA corresponds with 100% relative light output.

In the ModuLight the LED intensity is controlled via the setting of Analog output 1, linear from 0 - 4V. The setting of 0V corresponds to a driving current of 0 mA, 4 V corresponds to a driving current of 1000mA.



Example

Using a the RED LED (623nm, 260lm/W, output 160lm at 100%):



- setting 0V will give 0 light
- setting 4V will give 100% light (1000mA), corresponding with 160lm or 615mW (= 160/260)
- For 50% light output, a driving current is needed of ca. 430mA (see also graphs in paragraph 8 and 9). To achieve 50% a setting of 1.72V ($4V * 0.430/1$) is required, corresponding with 80lm or 307.5mW or (= 80/260)

8. LED specifications (Default configuration)



Cool White 6500k 200lm @ 1A





Blue 460nm 40lm @ 1A





Green 523nm 200lm @ 1A





Amber 590nm 105lm @ 1A





Red 623nm 160lm @1A



Deep Red 660nm 900mW @1A





Far Red 740nm 705mW @1A





9. LED specifications (Optional)



Cool White 5500k 227lm @1A





Warm White 3000k 180lm @1A





Violet 400nm 1400mW @1A





UV 365nm 1120mW @1A





IR 940nm 1150mW @1A



IR 850nm 1150mW @1A




Dental Blue* 460nm 1100mW @1A



* Dental Blue LED emitter provides superior radiometric power in the wavelength range specifically required for dental curing light applications.

10 Luminous Efficacy Tables

Wavelength λ (nm)	Photopic Luminous Efficacy V_λ	Photopic Conversion lm/W	Scotopic Luminous Efficacy V'_λ	Scotopic Conversion lm/W
380	0.000039	0.027	0.000589	1.001
390	0.000120	0.082	0.002209	3.755
390	0.000120	0.082	0.002209	3.755
400	0.000396	0.270	0.009290	15.793
410	0.001210	0.826	0.034840	59.228
420	0.004000	2.732	0.096600	164.220
430	0.011600	7.923	0.199800	339.660
440	0.023000	15.709	0.328100	557.770
450	0.038000	25.954	0.455000	773.500
460	0.060000	40.980	0.567000	963.900
470	0.090980	62.139	0.676000	1,149.200
480	0.139020	94.951	0.793000	1,348.100
490	0.208020	142.078	0.904000	1,536.800
500	0.323000	220.609	0.982000	1,669.400
507	0.444310	303.464	1.000.000	1,700.000
510	0.503000	343.549	0.997000	1,694.900
520	0.710000	484.930	0.935000	1,589.500
530	0.862000	588.746	0.811000	1,378.700
540	0.954000	651.582	0.655000	1,105.000
550	0.994950	679.551	0.481000	817.700
555	1.000.000	683.000	0.402000	683.000
560	0.995000	679.585	0.328800	558.960
570	0.952000	650.216	0.207600	352.920
580	0.870000	594.210	0.121200	206.040
590	0.757000	517.031	0.065500	111.350
600	0.631000	430.973	0.033150	56.355
610	0.503000	343.549	0.015930	27.081
620	0.381000	260.223	0.007370	12.529
630	0.265000	180.995	0.003335	5.670
640	0.175000	119.525	0.001497	2.545
650	0.107000	73.081	0.000677	1.151
660	0.061000	41.663	0.000313	0.532
670	0.032000	21.856	0.000148	0.252
680	0.017000	11.611	0.000072	0.122
690	0.008210	5.607	0.000035	0.060
700	0.004102	2.802	0.000018	0.030
710	0.002091	1.428	0.000009	0.016
720	0.001047	0.715	0.000005	0.008
730	0.000520	0.355	0.000003	0.004
740	0.000249	0.170	0.000001	0.002
750	0.000120	0.082	0.000001	0.001
760	0.000060	0.041	0.000000	0.000
770	0.000030	0.020	0.000000	0.000



Source: Table 6-1 of Williamson & Cummins, *Light and Color in Nature and Art*, Wiley, 1983. The Photopic conversion (lm/W) is obtained by multiplying V_λ by 683 and the Scotopic conversion is obtained by multiplying V'_λ by 1700 as suggested by those authors.



www.ivium.com