

EMPLOYER MANUAL

—

How to get the best use out of me
Let's get to work!

—



Measurement system for the
Anterior surface of the eye



Introduction

Thank you for purchasing this Eye Surface Profiler (ESP) system and for your confidence in Eaglet Eye. The ESP gives you highly precise measurements of the front surface of the eye, not only of the cornea but also a large part of the sclera.

We are committed to assuring your complete satisfaction with the ESP.

At Eaglet Eye, we continue our research and development to further extend the capabilities of all our products and solutions. This may result in improvements of certain features of the ESP system that have not yet been updated in the user manual. We strive to minimize these differences and come with regular updates of the manual.

For your safety it is essential that you read this manual carefully and that you familiarize yourself with its content before you start using this device. In particular, pay close attention to the safety instructions in this manual.

If you have any questions or desire further information on this product, please contact us via the website or email us. Our team will be glad to help you.

The Eaglet Eye team

Description

The Eye Surface Profiler (ESP) is a system using Fourier Transform Profilometry that measures the anterior surface of the cornea, the limbus and large parts of the sclera. Up to 20mm diameter of the visible anterior surface can be measured. The ESP requires fluorescein instilled in the eye so that the projected light is scattered back from the eye front surface.

The measurement process consists of 3 single shots in rapid succession that last a few milliseconds. The ESP calculates a 3-dimensional model of the anterior surface of the eye with on average more than 250.000 measurement points. From these measurement points, maps and parameters are calculated automatically that assist in the evaluation of the shape of the cornea, limbus and sclera. Separately, the pupil size and iris size are determined automatically from a top view image.

The ESP measures the complete anterior front surface of the eye up to 20mm without extrapolation. The corneal limbus is the border of the cornea and the sclera. The average corneal diameter is 11.77 ± 0.37 mm for males and 11.64 ± 0.47 mm for females. The new software will detect the limbus automatically and will give you detailed information about the whole cornea.

The ESP is a Profilometry based system and has therefore no interference with the shadow of the nose.

Contents


01. Quick Reference Guide	03
02. Components	04
Important elements	04
View window	04
Connections	04
03. Safety	05
Setting up and installation of equipment	07
Technical details of the ESP	08
04. Operational instructions	09
Home screen overview	09
Database Location	09
Add a patient	11
New measurement	11
Finding an existing patient	15
Search results	15
History	15
Viewing results	16
Parameters	16
Sagittal height	20
Figures	22
Source	22
Grid	23
Bisphere elevation	23
Bisphere 3D	24
Tangent angles map	24
Tangent angles 3D map	25
Height	25

Height 3D	26
Fluorescein simulation map	26
Tangential curvature	27
Axial curvature	28
Refractive power	29
Corneal elevation	29
Quad map	30
Notes	30
Evaluating functionalities	31
05. Exporting, Printing and Comparing Results	32
06. Lens Design Module	36
07. Troubleshooting	38
08. Maintenance and service	39
Care and maintenance	39
Casing	39
Optical components	39
Software installation	39
Minimum PC requirements	40
Conditions of warranty	40
Liability for malfunction and damages	41
09. Compliance	42
EMC	42
CE	42
Index	43
Address of manufacturer and service department	45


01 Quick Reference Guide

Prepare


Check: patient should not wear contact lenses.



Adjust chin rest so that eye is at level of marker.




Apply enough fluorescen






Focus & Centration




Let patient look at the red cross.



Focus by adjusting the alignment for the Z axis


- ✗ Out of Focus: too Close 
- ✓ In focus: points straight above each other 
- ✗ Out of focus: too far away 

Centration by adjusting the alignment for the x and y axes

- ✗ Too far right 
- ✓ Cross hairs in the center of the pupil 
- ✗ Too far to the left 


Review

"Blink, eyes wide open."




Hit fire button.

Assess image in quality control screen and save.



Deselect the images you don't want to process.

Push the 'save and add more' bar if you want to make more measurements right away



03

02 Components

Important elements

1. Forehead support
2. Optimal eye height indicator
3. Chin rest
4. Chin rest height adjustment swivel
5. Joystick
6. Fire button
7. View window

View window

8. Digital camera
9. Focus help LEDs
10. Fixation cross
11. Measure pattern projectors

Connections

12. USB Port
13. DC power inlet
14. Power adapter



12. 13. 14.



03 Safety

To assure safe operation, it is imperative that the instrument is used according to the instructions in this manual. Therefore, you should become thoroughly familiar with the content of the user manual before putting the system into operation. The Eye Surface Profiler (ESP) system must only be operated by trained personnel capable of using it properly on the basis of their training, expertise and practical experience.

Please save this user manual and ensure that it is accessible to operating personnel at all times.

This instrument must only be used for purposes as laid down in this user manual. Before using this instrument for the first time, you must be instructed on how to use it by our staff or an authorized dealer.

Only operate the instrument with the original components supplied by Eaglet Eye and if the instrument is in good working condition. Should the instrument be defective, do not operate it, and contact the supplier.

Before maintenance and cleaning the instrument, always pull out the power plug of the ESP and all equipment attached to it, for example a printer. Do not connect cables if this proves difficult. If it is impossible to make a connection, check to see whether the plug fits into the socket. In case you ascertain a defect, contact your authorized dealer to have it repaired by our service team.

When disconnecting electrical connections, always pull on the plugs, not on the cables. Supplementary equipment which is connected to the analog or digital interfaces of the instrument

must comply with relevant specifications. For detailed instruction on maintenance and cleaning turn to the chapter on Maintenance and Service in this manual.

If the ESP is coupled with external electrical equipment (for ex. a printer) this must not result in a decrease of patient safety. If the coupling of equipment causes the tolerance levels for leakage current to be exceeded, safety features, which include a disconnection device, must be provided for.

Do not operate the delivered equipment

- in areas with risk of explosion,
- in the presence of flammable anesthetics or volatile solvents such as alcohol, benzene or similar substances.

Do not use or store the instrument in damp rooms. Avoid placing the instrument in the vicinity of dripping, running or spraying water and ensure that no moisture can penetrate the instrument. For this reason, do not place any containers filled with fluids near the instrument. When cleaning the instrument with a damp cloth, be sure that no moisture enters it.

This instrument is a high-quality technical product containing sophisticated optics. In order to assure faultless and safe operation, service inspections must be performed according to the service plan of Eaglet Eye. See also the chapter on Maintenance and Service in this manual.

Eaglet Eye emphasizes that the user bears the full responsibility for the correctness of data measured, calculated or displayed using the ESP. The manufacturer will not accept claims based on erroneous data.

Setting up and installation of equipment

Before initial operation, the ESP needs to be set up and connected by an authorized dealer or service department.

Before installing the ESP, please consider the transport, storage and current room temperature where the device is to be installed.

Make sure the device is at room temperature. The ESP must be placed at a place prevented from direct (day)light from influencing the components.

A reflection-free examination must be assured. Therefore, the ESP should be used in a darkened room.

When placing the measurement head on the base plate, make sure the wheels are exactly aligned in the rail. This can be easily checked when holding a card at the base plate and at the end of the rail, the wheels should touch the card at both sides. If this is not the case, adjust the wheels in the rail.

This is an optical device and should be handled with care. Do not subject it to vibrations, jolts or high temperatures.



Important elements

Model: ESP 20200

Seizes (Depth*Width*Height)= 35 cm*35 cm * 48 cm

Weight: 12 kg

Energy and power: PLUG POWER SUPPLY, 10W, EU, 7.5V

- Power Supply Output Type: Fixed
- Input Voltage AC Min: 90V
- Input Voltage AC Max: 264V
- Output Voltage Nom.: 7.5V
- Output Current Max: 1.14A
- Output Power Max: 8W
- Output Connector: Barrel Plug
- SVHC: No SVHC (16-Jun-2014)
- Applications: General Purpose
- No. of Outputs: 1
- Output Current: 1.14A
- Output Voltage: 7.5V
- Power Rating: 8W

Power consumption: 8W

In this screen (see picture below) it is also possible to select a database if the eye care professional saved several databases.



The Eye Surface Profiler software home screen consists of three columns:

Column	Description
Customer column	From this column recent customers can be selected as well as new patients that can be added
New measurement column	Enables access to new measurements if a patient is selected
Results column	Contains a list of recent measurements

Add a patient

Before starting a measurement the patient must be entered into Eaglet-Eye Surface Profiler.

To add a patient, take the steps listed below:

1. Go to the home screen;
2. Click on the plus sign in the 'Customer column';
'Customer column' changes to four input fields for the patient data
3. Enter the required data;
 - Given name
 - Family name
 - Birthdate
 - Patient ID
4. Click 'Apply'.


New measurement

A new measurement can be made by pressing the button 'New measurement'. This button is disabled if no customer is selected.

Making a new measurement is done by following the steps listed below:

1. Make sure the patient does not wear contact lenses;
2. Select the patient from the customers (if the patient is not available enter the patient following the 'add a patient' procedure);
The name of the patient appears in the customer column title.
The button 'New measurement' is enabled

3. Apply fluorescein;
The patients eye reflects light used by the ESP
Make sure that a large part of the sclera is yellow in order to get the best results
4. Tell patient to place head on chin rest;
5. Adjust chin rest with swivel so that the patient's eye is roughly at level of marker;
6. Select new measurement in the ESP-software;
A live image of the patients eye is shown in the ESP-Software
7. Focus device and ask patient to look at the red cross;



The focus lights help you to focus by aligning and positioning the lights on the eye as is shown on the right.

8. Hit fire button;
The patient sees two flashes and the ESP-software shows the images. If the images are unsuitable because of the patient blinking, click 'discard'
9. Select measurements by clicking a check mark in the box next to the images and click 'save'.

The ESP software starts processing the images


Instead of clicking the 'save' button after a measurement, the operator can choose the 'save and add more' button to expand a measurement session. The picture below displays a measurement session of both eyes. If the operator would click the 'save' button (the one in the right), all of these nine measurements will be saved.



In order to provide the highest quality of the measured results there are some tools that will help the operator to optimize the focus, centration and sharpness of the camera before making a measurement. In the left bar in the picture above the operator can choose the 'eye tracking' option.


The colour bar on the right gives live indication of the quality of focus, brightness and sharpness while the operator adjusts the position of the measurement head. The two dots on the right of the eye still have to be placed above each other. The green circle (with green plus) has to be placed between the white focus points.

Please instruct the patient to look straight ahead at the red cross in order to get an optimal centration.



The measurements are saved in the ESP software.

All three measurements are selected and they can be unselected by clicking the check mark. Below these three measurements you can see an indication of the quality of the measurements, 100% meaning enough light, fluorescein and excellent focusing by the operator.



Quality	
Overall	90%
Focus	90%
Centration	100%
Brightness	95%

Finding an existing patient

A patient can be found in the ESP-software by entering the given name or family name in the 'Name' input field.


Search results

The patients matching the criteria will appear in the search results while typing. You can select the patient by clicking the name from the results.

You can also search by added notes, for example 'keratoconus'.

History

When no data is entered in the 'Name' input field, the most recent selected patients are listed in the left column.



Viewing results

Results can be selected in the Results column. By default the ESP-Software lists the measurement sessions by date.

When a patient is selected from the left column, measurement sessions of this patient are listed only. A measurement session can be expanded.



Parameters

The following figures will display the functionalities and parameters of the ESP-software.

Display
Source
Grid
<input checked="" type="radio"/> Bisphere elevation
Bisphere 3D
Tangent angles
Tangent angles 3D
Height
Height 3D
Fluo simulation
Tangential curvature
Axial curvature
Refractive power
Corneal elevation
Quad map
Notes

Refractive		
Sphere	-42.80	
Astigmatism	-0.80	92°
Sims		
SimsKs	41.90	8.13mm
SimsKf	39.70	8.49mm
Astigmatism	-1.80	87°
Central Asphericity		
Ks	46.20	8.14mm
Kf	44.90	8.38mm
Qs, Qf	-0.07	0.26
Qavg	0.00	
Best-fit bisphere		
Inner radius	8.67mm	
Limbus	6.46mm	1.10mm
Outer radius	13.38mm	
Keratocornis		
KP:	26%	
Pupil		
Limbus	12.92mm	1.10mm
HVID	12.27mm	
Pupil p	5.13mm	
Pupil offset	-0.39mm	0.05mm
PAP+, PAP-	7.67mm	7.71mm
FAC+, FAC-	8.53mm	8.78mm

Refractive

Sphere The Maloney spherocylinder is used. Spherical sphere equivalent corresponds to the plus half of the cylinder (Astigmatism).

Astigmatism The refractive error of the astigmatic eye stems from a difference in degree of curvature refraction of the two different meridians.

SimK

SimK parameters (Anterior Curvature) are calculated with the so called keratometric index, which is known from Placido topographers, and is equal to 1.3375. This does not correspond to the actual index of refraction of the cornea (1.376) since it takes into account the posterior surface parameters calculated with keratometric index are called simulated (SimKf, SimKs, ...). Steep (SimKs) and SimK flat (SimKf) are calculated from the pair of meridians 90° apart with the greatest difference in average power, from 0.5 to 2.0 mm distance from the centre.

This maximizes the Astigmatism parameter, which is the difference between SimKf and SimKs.

Astigmatism As mentioned above, the astigmatism value is the difference between SimKs and SimKf.

Corneal

Asphericity

Ks K values are calculated with actual refractive index of the cornea (1.376). Ks stands for steep curve.

Kf AK values are calculated with actual refractive index of the cornea (1.376). Kf stands for flat curve.

Qs (steep)

Qf (flat)

Qavg (average)

e is the symbol for eccentricity

Q is the symbol for asphericity

p is the symbol for shape factor

The three are related through the following formula:


$$Q = p - 1 = -e^2$$

Table conic sections

	e	p	Q
Oblate ellipse	$0 < e < 1$	$P > 1$	$Q > 0$
Sphere	$e = 0$	$P = 1$	$Q = 0$
Prolate ellipse	$0 < e < 1$	$1 > p > 0$	$0 > Q > -1$
Parabola	$E = 1$	$P = 0$	$Q = -1$
Hyperbola	$P > 1$	$P < 0$	$Q < -1$

Note that eccentricity (e) cannot distinguish between oblate and prolate ellipsoid. However e^2 is negative for an oblate ellipsoid and e^2 is positive for a prolate ellipsoid.

Best-fit Bisphere



Inner radius After subtracting a best fit sphere from the corneal part, the inner radius gives information about the radius of the inner sphere (the corneal sphere).

Limbus The limbus is indicated from the differences between the corneal and the scleral best fit spheres. There are two spots in the cross section where the two best fit spheres interfere (see picture above, the squares surround these two spots). The first number displays the mean distance from apex to the center of the interfering part of the two best fit spheres

	(corneal and scleral). The second number represents the mean length of this interfering area. This is still an indication.
Outer radius	After subtracting a best fit sphere from the scleral part, the outer radius gives information about the radius of the outer sphere (the scleral sphere).
Keratoconus	
KPI	The keratoconus prediction index is a compilation index of DSI, OSI, CSI, SAI, SimK1, SimK2, IAI, and AA (footnote 1). The mathematical formula for the simulated KPI, using simulated component indices, is described below: $\text{KPI} = 4.62 * [0.30 + 0.01 (-41.23 - 0.15 * \text{DSI} + 1.18 * \text{OSI} + 1.49 * \text{CSI} + 4.13 * \text{SAI} - 0.56 * \text{SimKs} + 1.08 * \text{SimKf} - 3.74 * \text{IAI} + 0.10 * \text{AA})] - 60.25$
Warning:	Keratoconus indices are designed to be an adjunct to a clinical examination, and should never be used as a sole assessment tool in the diagnosis of keratoconus (footnote 2).
Pupil	
Limbus	Indicates the start of the limbus area measured from the apex, the second value corresponds to the width of the limbus.
HVID	Horizontal Visual Iris Diameter
Pupil Ø	Pupil diameter
Pupil offset	Distance in x and y from the apex to the pupil centre.
PAP+, PAP-	The palpebral aperture relative to the pupil, these values are the length of the eye above and below the pupil.
PAC+, PAC-	The palpebral aperture relative to the corneal centre (apex), these values are the length of the eye above and below the apex.

Footnotes.

1. See Index, Page 43.


2. (Maeda N, Klyce SD, Smolek MK, Thompson HW. Automated keratoconus screening with corneal topography analysis. *Invest Ophthalmol Vis Sci* 1994; 35: 2749–57).

Sagittal height

Tsag-Nsag is the new parameter added to the ESP software that will quantify scleral asymmetry and give eye care professionals better control over lens fits.

Literature describes the flattening of the sclera on the nasal side of the eye, compared to the temporal side which shows no such behaviour, which introduces an asymmetry to the eye surface.

A pronounced asymmetry between the nasal and temporal scleral curves will impact lens fit.



The image above indicates a horizontal cross section of a typical eye surface. The sagittal height is calculated over a total chord of 16mm in this example. Tsag and Nsag each get attributed exactly half of the total chord length and the respective sagittal heights are determined from the apex outwards.

Other sagittal height parameters displayed in the ESP software are:

360° This is sagittal height calculated over as many meridians as possible, building towards the average sagittal height over the full 360 degrees of the eye.

Meridional This is the sagittal height for one single meridian.

MINsag/MAXsag Similar to the 360 degrees sagittal height is that it includes all available meridians at the given chord length. Out of the available sagittal heights it will then select the maximum and its angle and also the minimum and its angle.


Sagittal height		
Chord length	19.00	<input checked="" type="checkbox"/>
Meridional	5.87mm	32.8°
		48.7°
360°	5.88mm	44.0°
Tsag	6.50mm	30.0°
Nsag	5.46mm	40.6°
Tsag - Nsag	1.05mm	
MINsag	5.80mm	@169.0°
MAXsag	5.95mm	@8.0°

Figures

The following figures will show the different functionalities and results after the measurement has been made.


The colour scale of the tangent angles maps can be adjusted and therefore customized to the individual preference of the eye care practitioner.

Source




This picture of the source shows an excellent centration and focus. The white dots on the horizontal red line are an indicator for the centration; the two white dots on the right part of the iris (placed above each other) indicate the focus of the camera and therefore the sharpness of the measurement.

Grid




The grid map visualizes the vertical lines, as projected on the eye and reflected back into the camera. The clear contrast between the lines indicates that there was enough use of fluorescein during this measurement, which optimizes the quality. To read the quality of the measurement, there is a 'Quality bar' on the right, which indicates the centration, focus and brightness of the measurement.

Bisphere elevation




The bisphere elevation map allows the operator to simulate two best fit spheres on the eye: a corneal sphere and a scleral sphere. To get the most detailed version of the measurement you need to simulate two spheres due to the difference in steepness of the corneal curvature compared to the scleral curvature.

Bisphere 3D




The bisphere 3D map enables the operator to look at the surface of the eye in 3D in micron precision.

Tangent angles map



This map shows the surface of the eye in tangent angles and is able to point out the limbal area. In the parameter list the inner edge of the limbus and limbus width is indicated.


Tangent angles 3D map



The tangent angles map is also offered in 3D version. The colour scale of the tangent angles maps can be adjusted and therefore customized to the preferences of the eye care practitioner.


Tangent angles are important to establish the exact landing zone of the lens. The tangent angles are measured at every meridian and at every chord length.

Height




In the height map the details of the sagittal heights are displayed in a 2D map. These are the raw data.

Height 3D




In the height 3D map the details of the sagittal heights are displayed in a 3D version.

Fluorescein simulation map



The Fluorescein simulation map indicates a simulation of fluorescein patterns of a contact lens on the eye. It takes as input the front surface of the cornea and the back surface of a contact lens, and the specified relative positioning of the lens on the cornea. The simulation allows the eye care practitioner to assess the fit of the contact lens on the eye. The simulated lens is spherical in shape.

Axial curvature




The axial curvature map, also called 'sagittal' or 'power' map, is the most straight forward of all the topographical maps. It shows variations in corneal curvature as projections and uses colours to represent dioptric values.

Warm colours such as orange and red show steeper areas; cool colours such as blue and green denote the flatter areas. To get a global view of the corneal curvature as a whole, the axial curvature map is the best choice.


The downside to the axial curvature map is its tendency to ignore minor variations in curvature.

Refractive power




The refractive power map utilizes the measured dioptric power and applies Snell's law to describe the cornea's actual refractive power. The aspheric contour of the cornea and spherical aberrations can be compensated by a refractive map. The most important part of the refractive map is the centre. Abberations on the pupil almost invariably impact the visual performance.

Corneal elevation




The ESP is a proflometry based system and provides a qualitative and quantitative evaluation of corneal curvature of the whole cornea. An elevation map shows the measured corneal height from which a best fitted reference sphere has been subtracted. This map shows the true shape of the eye.

Quad map



The quad map displays the grid map, tangential curvature map, corneal elevation map and bisphere elevation map at the same time to compare those four.

Notes




The notes map can be used to describe the measurement or to report on the patient.

Evaluating functionalities

By right-mouse clicking on the picture different evaluating functionalities appear.


- Show eye allows to view the measured surface on the eye with the translucency of your choice.
- HVID shows the Horizontal Visual Iris Diameter.
- None deselected all previous clicked options.
- Limbus draws a circle around the calculated limbus.
- Pupil indicates the location of the pupil on the map.
- Conic indicates the flat and steep curves of the eye.
- Grid shows the compass rose angles.
- Values shows a grid of spot values on the map.
- Scleral angles shows the scleral angles of the eye measurement.
- The asphericity is calculated at 15 degrees intervals.

** Any map or cross section can be minimized in order to get a bigger view of either the cross section or the map.*




05 Exporting, Printing and Comparing Results

When a patient is selected in the home screen, it is possible to print, export, remove or compare any results simply by clicking the check mark on the right.




The same options are also to be found in the screen after opening a specific measurement.


Just select one or more measurements in the right bar (see right picture below) and then choose between the three options below: print, export or discard. The measurements you have selected will appear in the small box in the left of the screen.



The printing option enables the eye care practitioner to print both OD and OS (see printing example below).




The export option enables the eye care practitioner to export measurements in different formats and extensions.




To compare files just select two measurements from the right bar and click the comparing files option on the left box below the different maps bar. It is possible to compare two measurements from any day of one patient. It is not possible to compare measurements from several patients to each other.

There are three options to compare results.


1. Difference



2. Side by side (see picture below)



3. Overlap (see picture below.)




In the right of the screen you can find the parameters of the measurement in order to help the eye care practitioner to choose the right values for the lens design.

There are three options to display a lens design for the measurement: fit, lens and fluorescein simulation.

The screenshot shows a software interface with the following sections:

- Display:** Fit, Lens
- Fluor simulation:** (Selected)
- General:**
 - Thickness (mm): 0.100
 - Producer: Sumipro
 - Material: LDF90 (n=1.463)
- Shape:**

Sphere (mm)	Sagittal (mm)	Diameter (mm)
8.24	2.57	11.96
9.90	3.44	13.97
9.88	3.99	14.98
9.95	4.61	15.98
- Correction:**
 - Sphere (D): 0
 - Cylinder (D): 0
 - Axis (°): 0
- Calculate



07 Troubleshooting

- When starting 'New measurement' no image appears

To connect the ESP hardware a USB 2.0 connection is required.

Reassure if the USB port

to which the hardware is connected complies with the USB 2.0 standard.

- The camera image freezes Check if the camera image is frozen by moving your hand in front of the ESP camera. Take the steps listed below:

1. Close ESP software
2. Pull out USB from base of the measurement head
3. Insert USB back into base of the measurement head
4. Start ESP software

- Software is unable to process the image

1. Close the error dialog
2. Close the ESP software
3. Start the ESP software

- To recover the measurements that were lost during the process error, follow the next steps

1. Open file browser
2. Go to C:\Users\COMPUTERNAME\AppData\Roaming\Eaglet Eye\Eye Surface Profiler
3. Click measurement.zip



08 Maintenance and Service

Care and maintenance

Casing

When cleaning always pull out the plug! Do not use any aggressive cleansers, which contain chlorine, solvents, abrasive or caustic substances!

It is best to clean the surface of the casing with a soft cloth and an anti-static cleansing agent. Otherwise, wipe off outer surfaces with a damp cloth. In case of residue, remove with a mixture of equal parts of spirit and distilled water with a small dose of conventional detergent.

Optical components

There are no parts inside the instrument which should be cleaned by the user. Opening the casing without explicit authorization by Eaglet Eye will void the warranty.

Software installation

Eaglet Eye service department will provide you with the necessary software installation file(s) and thereafter with the required software update file(s). These will be in the form of .msi files and will come with complete step-by-step instructions.

Before installing a new version, please ensure that you de-install the older version first through the Windows Control Panel. To remove the older version please follow these steps: Open Programs and Features by clicking the Start button, then clicking Control Panel, next clicking Programs, and then clicking Programs and Features. Now select a program, and then click Uninstall.

Minimum PC requirements

ESP software will only run on Windows 7 or higher. The computer system should have at least 4 GB of RAM. Also, the computer must have at least one available USB 2.0 port.

Conditions of warranty

This Eaglet Eye product is a high-quality instrument. It was manufactured with great care using high-quality materials and modern production technology. It is important that you read the user manual before operation and observe the safety precautions.

You have a warranty on this instrument in accordance with legal warranty regulations starting on the day of purchase. This warranty covers all kinds of malfunction caused by faulty material and manufacture.

Malfunction caused by improper operation and external influences is not covered. Should you nevertheless have reason for legitimate complaint during the period of warranty, the malfunction will be eliminated at no cost to you.

These warranty claims can be made if a receipt documenting date of purchase is presented. If unauthorized persons manipulate the instrument, all warranty claims become void, for improper alterations and handling can cause considerable damage to the user and patients.

We request that complaints concerning damages which occur during transport or delivery be made to the shipping firm immediately and that the damage be documented on the consignment note so that the claim can be settled properly. Our terms of trade and delivery are valid as laid down at the time of purchase.

Liability for malfunction and damages

Eaglet Eye only considers itself responsible for the safety, reliability and proper function of the instrument if it is used in compliance with this user manual.

There are no parts attached to or inside the instrument that should be serviced or repaired by the user. If assembly work, expansions, adjustments, maintenance work, alterations or repairs are carried out by non-authorized personnel or if the unit is handled improperly, Eaglet Eye is not liable for any damages that might occur.

If authorized persons perform the kind of work described above, documentation of the nature and extent of the repair work must be demanded of them, if pertinent with a statement regarding changes made in specifications or application. The document must include date and particulars of services performed as well as name and address of firm and signature.

If desired, Eaglet Eye will put wiring diagrams, lists of replacement parts, additional descriptions and instructions for installation at the disposal of such authorized persons.

In case of maintenance and repairs, only original Eaglet Eye parts are to be used.

09 Compliance

EMC



CE
0344

Index

Area Analyzed (AA)

The simulated index is the ratio of the actual data area to the area circumscribed by a circle 4.5 mm in radius.

Average Central Power (ACP)

ACP is simulated as the average dioptric power of all points within the central 3 mm.

Central/Surround Index (CSI)

The simulated center/surround index reports the difference between the average area corrected power between the central area (3-mm diameter) and an annulus surrounding the central area (3–6 mm).

Differential Sector Index (DSI)

The simulated differential sector index reports the greatest difference in average area-corrected power between any two 45° sectors.

Irregular Astigmatism Index (IAI)

Irregular astigmatism index is simulated as the average summation of area-corrected dioptric variations along every semi-meridian for the entire analyzed surface and normalized by the average corneal power and number of all measured points.

IAI

i: semi-meridional position; j: point location relative to center; $P_{i,j}$: corneal power on the point (i,j); A: area which corresponds to power $P_{i,j}$

Inferior-Superior (IS)

Inferior-Superior index is simulated as the difference between the inferior and superior average dioptric values approximately 3 mm peripheral to the corneal vertex as defined by the center of the map. Opposite Sector Index (OSI) The simulated opposite sector index reports the greatest difference in average area-corrected power between opposite 45° sectors Surface Asymmetry Index (SAI) The surface asymmetry index is simulated as the centrally weighted average of the summation of differences in corneal power between corresponding points 180° apart on all meridians.

Standard Deviation Powers (SDP)

SDP is simulated as the standard deviation of all measured powers present on the map. Surface Regularity Index (SRI) The surface regularity index is simulated as a difference in power gradient between successive points on a hemi-meridian that is assigned a positive value and added to the running sum. This sum is divided by the number of points that went into the sum and then scaled. The mathematical formula of the simulated SRI index is:

Address of manufacturer and service department

You can receive additional information from our service department and authorized dealers.

Eaglet Eye B.V.
Jules Breukersstraat 21
6041 BP Roermond
The Netherlands
E-Mail: support@eaglet-eye.com
Web contact: www.eaglet-eye.com

