

ElePhor 24S

AUTOMATIC ELECTROPHORESIS ANALYSER

USER MANUAL

INDEX

	Page
1.0 DESCRIPTION OF THE INSTRUMENT	4
1.1 Standard components	4
1.2 Technical characteristics	5
1.3 Host computer requirements	5
1.4 Installation	6
2.0 TURNING ON	7
2.1 Instrument preparation	8
2.2 ELEPHOR 24Smanagement software start-up	9
3.0 'ELEPHOR 24S ' MANAGEMENT SOFTWARE	10
3.1 Main Panel	10
3.2 Starting Work	12
3.3 Patient Lists	14
3.4 Archive Management	16
- Storage	16
- Search by date	17
- Search by name	18
- Data backup	19
- Restore Data	20
- Deletion	21
3.5 Display graphs	22
- Delete areas	24
- Insert minimums	24
- Remove minimums	24
- Edit baseline	24
- Calibrate Fraction	25
- Highlight fraction	25
- Set fraction name	25
- Cancel Changes	25
- Print current graph	25
3.6 Configuration	26
- EXAM parameters	26
- VARIOUS	29
- REAGENTS	30
3.7 Printout	31
- Print graphs in archive	31
- Print graphs of the day	32
4.0 MAINTENANCE	33
4.1 Daily maintenance	33
4.2 Routine maintenance	33

5.0	QUALITY CONTROL QC	34
5.1	New Control	35
5.2	Modify Control	36
5.3	Delete Control	37
5.4	Exec a QC statistics	38
A	REMOTE HOST	40
A.1	Setting Remote Host Info	40
A.2	Start – Selecting data from Remote Host	41
A.3	Sending data to the Remote Host	42
A.4	Send Database Items to Remote Host	43
B	VETERINARY SCIENCE MANAGEMENT	44
B.1	Configuring species	44
B.2	Managing species	46
C	TROUBLESHOOTING	47
	REMOTE HOST INTERFACE SPECS	50

1.0 DESCRIPTION OF THE INSTRUMENT

ELEPHOR 24S is able to perform the electrophoretic procedures of the following clinical tests automatically:

- Serum proteins
- Haemoglobins
- Lipo proteins
- Urinary proteins

furthermore additional tests may be created and customized, such as Serum protein multifractionation.

ELEPHOR 24S uses supported cellulose acetate strips which are managed completely automatically. Simply insert them in the appropriate slots and tell the machine which strips to use.

The patient archive is stored on an external computer (HOST) connected to the machine through a standard USB cable. The number of patients which can be stored in the archive (never less than 100,000 in any case) therefore depends on the capacity of computer's magnetic storage medium (HARD DISK).

1.1 STANDARD COMPONENTS

ELEPHOR 24S comes complete with the following standard components:

- Eight-thin plate serum applicators for MICRO methods
- 24-place sampler (3 supports x 8 samples)
- 3 supports for cellulose acetate strips
- 1 reagent tank with two compartments for automatic reagent handling
- Usb cable
- Power cord
- Destainer, washing and waste tank connection tube (external diameter = 5 mm; internal diameter = 3 mm)
- One reagent tank (buffer)
- Manuals
- **ELEPHOR 24S** software CD-ROM

1.2 TECHNICAL CHARACTERISTICS

- Power supply voltage: 220 V 50Hz
- Power consumption: 250 W max
- Automatic reagent loading and handling
- Sample refrigeration by a Peltier cell
- Reagent (buffer and stain) recycling upon job completion
- Fluid circuit washing routine upon job completion
- Electrophoresis power supply unit with voltage (100V-300V) and current regulation (7mA-30mA)
- Pneumatic support handling
- Operating capacity of 16 tests per hour
- Measurement with non cleared supports
- Dimensions (WxHxP): 56x46x54 cm
- Weight : 30 Kg

1.3 HOST COMPUTER REQUIREMENTS

The minimum configuration for correctly interfacing with the **ELEPHOR 24S** system is:

- 700 MHz PENTIUM III PC (or equivalent)
- 128MB RAM
- 1.2 Gbyte hard disk
- 800x600 resolution 65000 colour video adapter
- 15” monitor
- 2X CD-ROM drive
- Windows XP SP1 operating system

Recommended configuration:

- 1200 MHz PENTIUM III PC
- 256MB RAM
- 1.2 Gbyte hard disk or better
- 800x600 resolution 16 million colours video adapter
- 17” monitor
- 2X CD-ROM drive or better
- Windows XP SP3 operating system or better

About 30MByte of free space is needed on the hard disk for installation. The remaining space may be used for the patient archive until full.

Note : to avoid visual artifacts set Windows ® Desktop theme in “Windows classic” mode and set numeric format in “0.00” notation using ‘.’ (dot) as decimal separator.

1.4 INSTALLATION

The instrument must be placed on a perfectly flat surface in a room with a temperature of between 15°C and 30°C and humidity of 40/80% max. Make sure the ventilation slits are not obstructed and leave at least 10 cm of space behind the back panel.

- Remove the **lock screw from the mechanical arm before use**
- Connect the USB cable provided to 'USB' connector on the left panel of the instrument (fig. 1.4.1) and connect the other end to the computer USB hub





fig. 1.4.1

- Connect the power cord to the Power Plug
- Connect the five connectors on the back of the instrument (fig. 1.4.2) using the tube provided



Where :

(1) =stainer tank, (2) = destainer tank, (3) = buffer tank, (4) = washing solution tank and (5) waste tank

2.0 *TURNING ON*

To turn the instrument on press the power switch on the left panel of the device; the system starts to move back to the home position.

Note: during the first activation, be sure to remove the lock screw from the mechanical arm !

2.1 INSTRUMENT PREPARATION

Before starting an electrophoresis process, prepare the electrophoretic chamber in the following way:

- 1 Pull the electrophoretic chamber out from its slot by the sides
- 2 Place the two floats in the electrophoretic chamber and fasten them using the support with holes provided
- 3 Fill the two compartments with 150 ml of buffer solution until the compartment level reaches the reference line; level off the quantity between the two compartments.
- 4 Put the electrophoretic chamber back into its slot and push home.
- 5 Pull out the sampler and pipette 25-30 μl of sample into each well (1)
fig. 2.1.1




Fig. 2.1.1

- 6 Fill the depositor washing compartment (2) with 1 ml of bidistilled water
- 7 Insert the special strip of absorbent paper in the compartment (3) in position E on the treatment rack
- 8 Put the sampler back in place and push home
- 9 Place the applicator on the support plate in position D (Fig. 2.1.2)
- 10 Place the supports to be prepared on the support plate (1-3 in Fig. 2.1.2)
- 11 Run the **ELEPHOR 24S** software to start the operating cycle.





Fig. 2.1.2

2.2 ELEPHOR 24S MANAGEMENT SOFTWARE START-UP

- Turn the computer on and let the operating system (Windows 9x) start up
- When the Windows Desktop appears (fig 2.2.1), click the start button and choose **ELEPHOR 24S** from the menu, or point the mouse at the **ELEPHOR 24S** icon on the Windows Desktop and double-click.

Note : to avoid visuals artifacts set Windows ® Desktop theme in “Windows classic” mode.




Follow the instructions in section 3.2 to start working.

3.0 'ELEPHOR 24S' MANAGEMENT SOFTWARE

The Management Software provided with the **ELEPHOR 24S** analyser allows electrophoresis traces to be selected, customized, corrected, stored and printed through an intuitive hierarchical organization. Every function or group of functions is brought together on a single screen, which can be selected from the main panel by simply clicking the corresponding button

3.1 MAIN PANEL



- (1) Archive Management Button
- (2) Patient List Button
- (3) Start Button
- (4) Graph Display Button
- (5) Configuration Button
- (6) Print Management Button
- (7) Quality Control Button

3.2 STARTING WORK


Before starting a “work routine” you must Log in into the software (default user name ADMIN with password admin) and choose the type of test and the number of supports the sample to be analysed is to be deposited on. The position of the support corresponds to the position of the samples contained in the sample holding base.

To select a test and the samples to prepare and/or measure, carry out the following procedure on the Management Software (hereinafter referred to as MS) main screen:

- point the mouse at the “Start” button and left-click:




the following screen appears:



the screen is split into two areas. The left part shows the sample grid (seen from above) and the right consists of buttons for the methods and functions which may be carried out (support preparation, support measurement or both).

The sample grid is arranged in rows from top to bottom. Therefore row 1 corresponds to support 1, row 2 to support 2, etc.; the position of the first sample in each row is indicated by a yellow dot.

Therefore by selecting a row from the sample base, you also select the position of the support to be used:



All you have to do after that is select the type of test required by clicking the corresponding button.

Note that by default the MS is set for support preparation and its subsequent measurement (the two checkboxes are selected). Therefore, if you only wish to do one of the two functions, deselect the option concerned by pointing the mouse at the checkbox and left-clicking (the box is no longer selected):

- Once the supports/samples to be used and the test to be performed have been selected a button marked “Start” appears:

By clicking this button, the chosen process starts and a status screen appears which shows the progress of the various stages of the method for each support, as well as data regarding certain parts of **ELEPHOR 24S** (migration voltage and current) and the levels of the reagents in the tanks.

Note: it is advisable to wait until ELEPHOR 24S has carried out all the initialization operations before leaving the MS workstation unmanned (if necessary) since error messages may appear which must be corrected by the operator so that the processes may be restarted correctly.

One typical error which may occur is insufficient reagent in the storage tank. In the case of a method with automatic reagent level handling, this results in failure to reach the operating level and leads to the corresponding error message.

Once the ELEPHOR 24S has started up correctly, the MS workstation may be left unmanned.

In order to correct an error, the operator must fill the corresponding reagent storage tank, click the “Continue” button on the error screen, and click the “Start” button on the main screen again.

Note: button named “Remote Host” will be explained into appendix A

3.3 PATIENT LISTS

This button opens a screen with information relating to the patients to be tested. Up to 24 patients may be managed at a time (8 patients x 3 supports) plus eight other patients which may be retrieved from the archive using search functions.

PATIENT LISTS - PATIENT ARCHIVE

Name	<input checked="" type="checkbox"/> Age	<input checked="" type="checkbox"/> Sex	<input checked="" type="checkbox"/> Date	<input checked="" type="checkbox"/> ID Code	<input checked="" type="checkbox"/> T.P.	<input checked="" type="checkbox"/> Division	Q.C.
1	<input type="text"/>	<input type="text"/>	29/04/10	<input type="text"/>	<input type="text"/>	<input type="text"/>	-----
2	<input type="text"/>	<input type="text"/>	29/04/10	<input type="text"/>	<input type="text"/>	<input type="text"/>	-----
3	<input type="text"/>	<input type="text"/>	29/04/10	<input type="text"/>	<input type="text"/>	<input type="text"/>	-----
4	<input type="text"/>	<input type="text"/>	29/04/10	<input type="text"/>	<input type="text"/>	<input type="text"/>	-----
5	<input type="text"/>	<input type="text"/>	29/04/10	<input type="text"/>	<input type="text"/>	<input type="text"/>	-----
6	<input type="text"/>	<input type="text"/>	29/04/10	<input type="text"/>	<input type="text"/>	<input type="text"/>	-----
7	<input type="text"/>	<input type="text"/>	29/04/10	<input type="text"/>	<input type="text"/>	<input type="text"/>	-----
8	<input type="text"/>	<input type="text"/>	29/04/10	<input type="text"/>	<input type="text"/>	<input type="text"/>	-----

Electrophoresis of SEROPROTEINE

The "**PATIENT ARCHIVE**" menu contains the data relating to the support selected from the archive.

The patient's data (name, surname, age, etc.) and total proteins value may be edited. The "**SHOW**" button displays the graphs which correspond to the selected support. Any corrections or changes to the traces displayed may be made on this screen.

PATIENT LISTS - CURRENT DATA

	Name	<input checked="" type="checkbox"/> Age	<input checked="" type="checkbox"/> Sex	<input checked="" type="checkbox"/> Date	<input checked="" type="checkbox"/> ID Code	<input checked="" type="checkbox"/> T.P.	<input checked="" type="checkbox"/> Division	Q.C.
1		<input type="text"/>	<input type="text"/>	29/04/10	<input type="text"/>	<input type="text"/>	<input type="text"/>	-----
2		<input type="text"/>	<input type="text"/>	29/04/10	<input type="text"/>	<input type="text"/>	<input type="text"/>	-----
3		<input type="text"/>	<input type="text"/>	29/04/10	<input type="text"/>	<input type="text"/>	<input type="text"/>	-----
4		<input type="text"/>	<input type="text"/>	29/04/10	<input type="text"/>	<input type="text"/>	<input type="text"/>	-----
5		<input type="text"/>	<input type="text"/>	29/04/10	<input type="text"/>	<input type="text"/>	<input type="text"/>	-----
6		<input type="text"/>	<input type="text"/>	29/04/10	<input type="text"/>	<input type="text"/>	<input type="text"/>	-----
7		<input type="text"/>	<input type="text"/>	29/04/10	<input type="text"/>	<input type="text"/>	<input type="text"/>	-----
8		<input type="text"/>	<input type="text"/>	29/04/10	<input type="text"/>	<input type="text"/>	<input type="text"/>	-----

Electrophoresis of SEROPROTEINE

Select "**CURRENT DATA**" from the "**PATIENT MANAGEMENT**" menu to input the patient data corresponding to the pherogram to be measured.

Once this function is chosen, the set of selected supports appears and it is possible to input patient data for all the supports displayed.

Simply type the surname, name, sex, sampling date (the default is always the current data), patient ID, total proteins and department into the **Patient Data** table


The Q.C. field is used to select a previously defined Control.

3.4 ARCHIVE MANAGEMENT

This screen provides complete control over the report archive. Patients may be searched for by name or by storage date, reports which have not yet been stored or edited may be saved, and data may be backed up on removable media or deleted.

ARCHIVE MANAGEMENT - STORE

Select "**STORE**" from the "*ARCHIVE MANAGEMENT*" menu to save the supports measured along with the corresponding patient data in the archive.




When this option is selected, all valid supports are selected for storage. Deselect any supports you do not wish to save by left-clicking beside them. Once selection/deselection is complete, click "**DONE**" to start saving the graphs in the archive.

ARCHIVE MANAGEMENT - SEARCH BY DATE

Select "**SEARCH BY DATE**" from the "*ARCHIVE MANAGEMENT*" menu to search for a support in the archive by its storage date.

Select the **year**, then the **month** then the **day** concerned: a list of supports stored on the date input appears. Scroll through the list to find the support required.



Once selected, click "**DONE**" to load the support.

The patient data for the chosen support may be displayed by choosing "**PATIENT ARCHIVE**" from the "**PATIENT LISTS**" menu.

ARCHIVE MANAGEMENT - SEARCH BY NAME


NAME	ID	EXAM	DATE
USER1	123456789	SIEROPROTEINE	12/1/04
USER2	AAA123456	SIEROPROTEINE	12/1/04
USER3	LOL123455	SIEROPROTEINE	12/1/04
USER4		SIEROPROTEINE	12/1/04
USER5		SIEROPROTEINE	12/1/04
USER6	ACABBBACC	SIEROPROTEINE	12/1/04
USER7		SIEROPROTEINE	12/1/04
USER8		SIEROPROTEINE	12/1/04

Select "**SEARCH BY NAME**" from the "*ARCHIVE MANAGEMENT*" menu to search for a support in the archive by the patient's name, the patient's exam or id code (any combination is valid).

Type the patient's name into the appropriate field to start searching: the result is highlighted in blue. Left-click the result to select the patient found along with the types of tests concerned and their storage date and press Confirm to display it.

The patient data relating to the graph and support selected may be displayed by choosing "**PATIENT ARCHIVE**" from the "**PATIENT LISTS**" menu.

ARCHIVE MANAGEMENT - BACKUP



It is advisable to save the archive data periodically on a floppy disk or backup unit (for example IOMEGA ZIP).

Select "**BACKUP**" from the "**ARCHIVE MANAGEMENT**" menu to save all or part of the data archive.

In order to make a backup:

1. Select the backup unit (floppy disk drive A: by default)
2. Select the type of backup
 - Whole Database (default)
 - One particular year only
 - One particular month only.


Then input the year and/or month to be saved.

Once done, click the **DONE** button to save the data on the backup media selected (for example floppy disk A:)

In the case of floppy disks, the backup procedure may require more than one disk; it is advisable to number each disk used during the procedure.

The last disk must be labelled "**Restore Disk**" and may be used to restore the data saved if necessary.

ARCHIVE MANAGEMENT - RESTORE



Select "**RESTORE**" from the "**ARCHIVE MANAGEMENT**" menu to restore data previously stored on BACKUP media to the archive.

In order to restore data:

1. Select the backup unit (floppy disk drive A: by default)
2. Select what to restore
 - Whole archive (default)
 - One particular year only
 - One particular month only.

Then select the year and/or month to be restored and insert the disk labelled "**Restore Disk**" in the chosen unit (for example floppy drive A:).

Once done, click the "**DONE**" button to start restoring the data. If necessary, insert the other numbered disks in the order in which they are requested.

Note: button "Remote Host" will be explained into appendix A.

ARCHIVE MANAGEMENT - DELETION



Select "**DELETE**" from the "**ARCHIVE MANAGEMENT**" menu to delete all or part of the archive.

First of all, select the type of deletion required:

- Delete Year
- Delete Month
- Delete Support.

Once the support has been selected, you may also choose to delete only certain graphs.


By default all graphs are selected for deletion. To prevent certain graphs from being deleted, deselect them by left-clicking.

3.5 DISPLAY GRAPHS

After preparing/measuring a support, its graphs may be displayed.

The "**DISPLAY**" menu shows a list of available supports (prepared and measured).

By selecting a support with the mouse, a preview appears where it is possible to display the contents of



the whole support and immediately see if any errors have occurred during the preparation and measurement stages.


In order to leave the preview of the entire support and display the individual graphs for editing or correcting, simply click the mouse on one of the traces displayed.

Note:

In order to display traces retrieved from the archive, select '**DISPLAY**' again from the '**PATIENT LISTS**' -> '**PATIENT ARCHIVE**' screen.

It is possible to display both the patient data associated with a graph and the percentage values of its individual fractions from the individual graph **correction** and display menu.

Every graph has a correction menu with icons placed along side it.



From the menu it is possible to:

- Delete a trace area
- Insert a minimum between two fractions
- Delete a minimum between two fractions
- Correct the baseline
- Calibrate the value of a fraction
- Highlight a fraction by dividing it
- Change the fraction label (name)
- Restore a graph
- Print a graph

Left-click one of the icons on the correction menu to select the type of correction required. Right-click to deselect it.

DELETE TRACE AREAS

When this button is clicked, a vertical light blue line appears which may be positioned at the point where the area to be deleted is found using the mouse.

A yellow square appears next to the displayed graph showing the optical density (O.D.) of the trace. This parameter changes with the movement of the cursor line.

The O.D. is useful when determining the position where deletion is necessary.

By left-clicking the first half of the graph, everything between the beginning of the graph and the cursor line is deleted; by left-clicking the latter half of the graph, everything between the cursor line and the end of the graph is deleted.

Right-clicking or selecting another button from the correction menu cancels the deletion function.

INSERTING MINIMUMS

When this button is clicked, a vertical light blue line appears which may be positioned at the point where the minimum is to be inserted using the mouse.

A yellow square appears next to the displayed graph showing the optical density (O.D.) of the trace. This parameter changes with the movement of the cursor line.

The O.D. is useful when determining the position where the minimum is to be inserted between two graph fractions.

Left-click to insert the minimum.

Right-clicking or selecting another icon from the correction menu cancels the minimum insertion function.

REMOVE MINIMUMS

When this button is clicked, a vertical light blue line appears which may be positioned on the minimum to be removed using the mouse. Left-click to remove the selected minimum.

Right-clicking or selecting another icon from the correction menu cancels the minimum removal function.

A yellow square appears next to the displayed graph showing the optical density (O.D.) of the trace. This parameter changes with the movement of the cursor line.

The O.D. is useful when checking the position of a minimum between two graph fractions.

EDIT BASELINE

This function is useful for correcting supports which have undergone uneven destaining, and therefore have ODs above the baseline value at the beginning or towards the end of the graph.

In this case the graph must be corrected so that its beginning and/or end have a uniform O.D. value equal to zero.

When this button is clicked, two adjustment bars appear along side the graph which may be used to move the baseline where necessary.

Select one of the two bars by left-clicking, and adjust the corresponding cursor either using the cursor keys on the keyboard or by dragging with the mouse (hold the left button pressed). Adjust the other cursor if necessary in the same way.

Once cursor adjustment is complete, click the button again and correct the graph baseline obtained by moving the end points of the trace onto the zero O.D. line.

CALIBRATION FRACTION

The % value of a single measured fraction may be calibrated applying a calibration factor. When this button is pressed a grid of values will be displayed. Set the single calibration factor by changing the value and pressing the 'Enter' key to validate it. After you changed the values click on "Confirm" to apply it.

SET FRACTION NAME

Click the grid containing the fraction names to display an input box where a new name for the fraction or sub-fraction may be input.

INSERT NEW FRACTION

Right clicking on the grid containing fraction names it will be displayed an input box where setting a new fraction name and a new (optional) range of values in pct.

The new added fraction will be placed where the right mouse button was pressed.

HIGHLIGHT FRACTION

This function highlights certain fraction characteristics by dividing it into several areas (maximum of 3). In this way precise information may be obtained regarding the concentrations of certain pathological components (e.g. monoclonal peaks).

CANCEL CHANGES

If you make irreversible mistakes while correcting a graph (for example when correcting its baseline or deleting an area), it is possible to restore the graph as measured. Click this button to restore the graph to its original state.


The graph is restored when it is displayed. Once a new support has been displayed, it is no longer possible to go back and restore the original state of a previously displayed graph.

PRINT GRAPH

Click this button to print the graph displayed.

Printing takes place according to the settings made by choosing "VARIOUS" from the "**CONFIGURATIONS**" menu".

3.6 CONFIGURATIONS - EXAM



Select "**EXAM**" from the "**CONFIGURATIONS**" menu to configure the parameters for individual tests.

The following parameters may be configured on this screen for each test:

- Test name
- Number of fractions and their names (max 10)
- % reference tables
- g/dl reference table for total proteins
- Deposition step (MICRO and SEMIMICRO)
- Number of points to eliminate in the deposition area
- Migration voltage
- Migration current
- Multifractionation technique (serum proteins only)
- Method operating values

Note: you can change the test name in order to rename a pre-defined exam.

A new name is automatically selected for a fraction (or edited) after the name of the fraction itself is written. If the FRACTIONS field is not completed, the fraction itself is not included in the graph minimums calculation.

Numerical percentage values are expressed as ##.# (for example 12.7), i.e. with only one decimal place and a decimal point separator. The program beeps if characters or numbers with a different format are typed in.

After inputting the min-max % reference table, a congruity test may be carried out by clicking the "T" key found beside the table.

Any errors found are automatically highlighted in sequence by a box which appears on the cell(s) containing the incorrect value(s).

The min-max total proteins reference values in g/dl together with the min-max % reference table determines the min-max reference table in g/dl for the single fractions;

For example:

$$\text{Min. ALBUMIN value in g/dl} = \text{Min. g/dl value for TOTAL PROTEINS} * \text{Min. \% value for ALBUMIN}$$

If a SEMI-MICRO deposition step is chosen, only four traces may be selected (choose "**CURRENT DATA**" from the "**PATIENT LISTS**" menu).

The odd traces are available for measurement, whereas the even ones are automatically defined as NOT VALID.

On this screen, you can choose whether to carry out constant voltage or constant current migration for the type of test concerned and set the appropriate value.


Clicking the "METHOD" button displays a new screen where the time of operations performed by ELEPHOR 24S may be edited. All times are in seconds, changing these times changes the method and affects the quality of the traces produced. The default method parameters are a good compromise between process speed and result quality. It is therefore inadvisable to change these settings unless you are fully aware of what you are doing.



Parameters modification related to tracks editing are shown on this screen:

- Application Line : by decreasing this value more points at the end of the track will be displayed. (needs membrane re-scanning)
- Filter Calibrator : allows to change the filter parameter used by the instrument to compensate variation density of the acetate membrane (default value is 23).
Change of this parameter is password protected (PW : device) and should be performed only by authorized technical service personnel .
It is highly recommended not to modify this parameter unless to be able to prove its correctness by means of a proper serum calibrator.
- Calibration Factors : allows to set a calibration factor for each fraction defined into the exam

CONFIGURATIONS - VARIOUS



Select the "**VARIOUS**" button from the "**CONFIGURATIONS**" menu to edit the printing and instrument connection parameters.

REPORT PRINTOUT CONFIGURATION PARAMETERS

"PAGE HEADER": this is the space reserved for the header of the printed report, it's possible to modify font, size and text formatting.

"Notes at Bottom of Page": this space is used to print fixed text at the bottom of values

"Additional Line at the Bottom Of Page": a single line can be typed in which is printed under the footnotes with the same layout as used when writing.

"Print Report in Colour": it is possible to choose either a colour or black and white trace printout. This selection only refers to the trace printout. Depositions are normally printed in red.

"Mark Values out of Range": this option prints any fraction values which are out of the normal range in bold.

"Kit Management": it's reserved for analyzers with smart card reader.



It's not a real configuration, but it offers a way to manage some operations on Reagents individually (Fill, Flush, Recycle). In this screen it is also possible to Recycle all Reagents (Recycle button) or to Recycle Reagents and Clean Sides (Recycle + Clean button).

VERY IMPORTANT NOTE:

YOU MUST BE SURE THAT ALL CONNECTORS OF SIDES ARE PLUGGED IN !
ANYWAY IT MAY RESULT IN DAMAGES TO THE CIRCUIT OF REAGENTS.

3.7 PRINTING

The whole printing process may be set up in one place on this screen, where the graphs to be included in the report may be selected.

PRINT - PATIENT ARCHIVE




Select "**PATIENT ARCHIVE**" from the "**PRINT**" menu to print the support selected from the patient archive.

When this option is selected, all the valid graphs for the support are selected for printing.

Graphs may be deselected by left-clicking beside the traces which you do not wish to print. Once selection is complete, click "**DONE**" to start printing.

Printing may be interrupted at any time by clicking the "**CANCEL PRINT**" button on the top right of the screen.

PRINT - CURRENT DATA



Select "**CURRENT DATA**" from the "**PRINT**" menu to print the supports produced during the operating cycle.

When this option is selected, all valid supports and graphs are selected for printing.

Supports and/or graphs may be deselected by left-clicking beside the supports or graphs you do not wish to print.

Once deselection is complete, click "**DONE**" to start printing the graphs.

Printing may be interrupted at any time by clicking the "**CANCEL PRINT**" button on the top right of the screen.

4.0 MAINTENANCE

4.1 DAILY MAINTENANCE

At the end of each work session, it is advisable to carry out the following cleaning operations:

- Wash the serum plate under running water and dry with absorbent paper
- Empty the electrophoretic chamber and recycle the buffer; remove the floats from their housings, rinse them with deionized water and leave them to dry.
- Remove the acetate supports from the strip racks.
- Clear the depositor by a soft bristle brush.

4.2 ROUTINE MAINTENANCE

It is advisable to wash the following parts more carefully at regular intervals:

- Migration chamber: fill the chamber, fill it with distilled water and leave it for an hour
- Clear the buffer side:
Remove the buffer tank and put a tank with a 10% hypochlorite solution. From the Reagents Menu, click on the buffer 'Fill' button and after a couple of seconds press the 'Recycle' button.
- Immerse the serum plate in 10% hypochlorite solution for an hour and rinse with distilled water.

5.0 QUALITY CONTROL

Here must be defined all Samples used for control purpose. It's possible define one or more samples for each exam indicating the average value attended for every fraction.

When a Control is defined it's possible to select it by the Q.C. dropdown menu into the "List of Patient" screen.

The screenshot shows a software window titled "Quality Control". The window has a standard Windows-style title bar with a close button (X) in the top right corner. The main area is divided into two sections. On the left, there are two labels: "Control :" followed by an empty text input field with a dropdown arrow, and "Exam :" followed by a dropdown menu currently displaying "SIEROPROTEINE". On the right side, there is a vertical panel containing two buttons: "New Control" at the top and "Search" below it. The central area below the input fields is a large, empty rectangular box, likely intended for a list of controls or data.

5.1 NEW CONTROL

The first thing to do to properly use the QC management is to select the type of exam from the drop-down "Exam:" In this menu we offer exams currently configured in the software, click the mouse on the exam for which you want to make a QC or where you want to add a new control.

By default the software comes with no parameters related 'standard' samples, as this is made by "known sample" for the laboratory or using special kits.

The screenshot shows a software window titled "Quality Control". At the top left, there is a "Control:" label followed by an empty text input field. Below it is an "Exam:" label followed by a dropdown menu showing "SIEROPROTEINE". To the right of these fields is a "New Control" button. Below the "Exam:" dropdown is a "Reference Values" section containing a table with 10 columns and 2 rows. The first cell of the second row is highlighted in black. To the right of the table is a "Search" button. At the bottom of the window are "Confirm" and "Close" buttons.

First create a control sample pressing the button "New Control":


after selecting the exam by the drop-down menu and fill in the text box "Control" with the name that identify that control.

Fill in the bottom grid from left to right with the theoretical values of the control and click on "Confirm".

5.2 MODIFYING CONTROL DATA


Once you create the control is always possible to change the values expected in the following way:

- Select the exam
- From the control box, select the name of the control to change
- Once you pop the grid values, select the "Edit" button and proceed to the modification of values in the grid



The dialog box titled "Quality Control" has a title bar with a close button. It contains two dropdown menus: "Control:" with "TEST CONTROL1" selected and "Exam:" with "SIEROPROTEINE" selected. To the right are three buttons: "New Control", "Delete Control", and "Search". Below the dropdowns is a "Reference Values" table with five columns and one row containing the values 60.2, 2.0, 7.8, 12.4, and 17.6. A "Modify" button is located below the table.

60.2	2.0	7.8	12.4	17.6



The dialog box titled "Quality Control" is identical to the previous one, but the "Reference Values" table now has ten columns. The first five columns contain the values 60.2, 2.0, 7.8, 12.4, and 17.6, while the remaining five columns are empty. The "Modify" button has been replaced by "Confirm" and "Close" buttons.

60.2	2.0	7.8	12.4	17.6					

Click on "Confirm" to store the modified data.

5.3 DELETING A CONTROL

To delete a control from the list of available controls, select the Exam and the control name as done in the previous paragraph, once the pop-up screen:

The image shows a 'Quality Control' dialog box. At the top, there are two dropdown menus: 'Control : TEST CONTROL1' and 'Exam : SIEROPROTEINE'. Below these is a table titled 'Reference Values' with 10 columns and 2 rows. The first row contains the values 60.2, 2.0, 7.8, 12.4, and 17.6, followed by five empty cells. The second row is empty. To the right of the dialog box are three buttons: 'New Control', 'Delete Control', and 'Search'. At the bottom of the dialog box are two buttons: 'Confirm' and 'Close'.

Reference Values									
60.2	2.0	7.8	12.4	17.6					

Click on “Delete Control” and press “Confirm” to submit the operation.

5.4 EXECUTION OF A QC STATISTICS

To perform a useful statistical QC for a monthly or annual review and select the control you want to make statistics and press "Search":




It will appear a box with a list of all the years of work in which the reference control was used, clicking the right mouse button on the year with a series of options appear:

- Run QC: QC statistics starts and displays the result
- Delete: delete the selected item (Year / Month) from the QC
- Abort: Closes the menu

If the entry for the year you left click mouse, a screen appears containing the list of months, the year selected, which was used as the reference standard.



Now right-click mouse menu appears as described above from which to start statistics on the option with the click left button:



The screen displays the average values for each reference sample, relative to the day of the month in which it was performed for each portion of the examination in which the QC will appear:

- Reference value
- Median
- Variance
- Diversion
- CV%
- The total number of points on which it is performed to Statistical


Values of samples are shown on the graph in the range of confidence-2s, +2 s that is in the range of 2 times the standard deviation.

Pressing the button press produces printed copy of the QC for all these villages, regardless of the currently displayed.

APPENDIX A *REMOTE HOST*

ELEPHOR 24Ssoftware is capable of exchange data with a remote host in order to manage exams Work list automatically.

A.1 *SETTING REMOTE HOST INFO*




Pressing the "**REMOTE HOST**" button you'll be able to indicate where remote host data are located. You must first indicate the name of input file without extension (this file contains the work lists) after this you can use the browse button "... " to locate your file folder. The same is for the output file (this file contains the results of exams).

A.2 START – SELECTING DATA FROM REMOTE HOST

In order to select data from a remote host you must:

- Press the button "**START**" from the main screen of program
- Press the button "**REMOTE HOST**"



Now you'll see all available exams queued and for each exam the work list.

Selections are possible into the right box clicking on the little box on the left of the patient name by the left button of the mouse. Press the right button to select all patient listed into the box.


When you have selected some patients press the "**Select Items**" button to add this selection to the sample grid.

Note: you can select only a single exam list per operative cycle.

Pressing the button "**Delete All Selections**" all elements selected for that exam will be deselected, pressing the button "**Delete Last Selection**" the last elements queued are deselected.

To obtain a printed work list press the button "**Print All Selections**".

A.3 SENDING DATA TO THE REMOTE HOST



To send the results of an exam to the remote host you must:

- Press the button Start from the main screen an then select "**REMOTE HOST**"
- Press thebutton "**Items To Send**"
- Select items into the box
- Press the button "**Send Items**"

It's also possible to send data from the local database to the remote host (see next session).

A.4 SEND DATABASE ITEMS TO REMOTE HOST

NAME	AGE	SEX	ID	EXAM	DATE
<input type="checkbox"/> USER1	30	M	123456789	SIEROPROTEINE	12/1/04
<input checked="" type="checkbox"/> USER2	22	F	AAA123456	SIEROPROTEINE	12/1/04
<input type="checkbox"/> USER3	29	F	LOL123455	SIEROPROTEINE	12/1/04
<input checked="" type="checkbox"/> USER4	66	F		SIEROPROTEINE	12/1/04
<input checked="" type="checkbox"/> USER5	50	F		SIEROPROTEINE	12/1/04
<input checked="" type="checkbox"/> USER6	0	M	AC&BBBACC	SIEROPROTEINE	12/1/04
<input type="checkbox"/> USER7	0	M		SIEROPROTEINE	12/1/04
<input type="checkbox"/> USER8	0			SIEROPROTEINE	12/1/04

To send Database items to a remote host you must:

- Select the button "**Database**" from the main screen
- Press the button "**Remote Host**"
- Search the items to send using the name, exam and data filter
- Select items you want to send
- Press the button "**Send**"

APPENDIX B VETERINARY SCIENCE MANAGEMENT

This appendix describe program items to manage animal species, by using a dedicated software program. (version 'V')

All things described for the standard version are valid too.

B.1 CONFIGURING SPECIES

Using the "Settings" window you can create for each exam one or more species with a particular range of normal values and a particular number of fractions(max 10 fractions).

Note: Veterinary program version doesn't allow the creation of a user defined exam.



Pressing the button "SPECIES" will appear a drop-down menu where you can:

1. Create a new species
2. Select a species
3. Delete a species

You can create a species by pressing the button "SPECIES" and choosing "Create Species", it will be created a new table of normal ranges and a new label where you can input the name of the created species.

To delete a species press the button "SPECIES" and choose "Delete Species", it will appear the list of species available where you can select the item to delete.

Note : the default species cannot be deleted, but can only be modified in all its parts.

B.2 MANAGING SPECIES

Species management can be done for each sample by selecting the appropriate species item into the "Work Lists" window.

Where not stated, default species parameters are used.

Samples recalled from the database not allow species to be changed.

	Nominativo	<input checked="" type="checkbox"/> Et�	<input checked="" type="checkbox"/> Sex	<input checked="" type="checkbox"/> Data	<input checked="" type="checkbox"/> Cod. ID	<input checked="" type="checkbox"/> P.T.	<input checked="" type="checkbox"/> Reparto	SPECIE
1	TEST3_1	31_	M	14/03/05	331	3.1	3333331	CANE
2	TEST3_2	32_	M	14/03/05	332	3.2	3333332	CANE
3	TEST3_3	33_	M	14/03/05	333	3.3	3333333	CANE
4	TEST3_4	34_	M	14/03/05	334	3.4	3333334	CANE
5	TEST3_5	35_	F	14/03/05	335	3.5	3333335	CANE
6	TEST3_6	36_	F	14/03/05	336	3.6	3333336	CANE
7	TEST3_7	37_	F	14/03/05	337	3.7	3333337	CANE
8	TEST3_8	38_	F	14/03/05	338	3.8	3333337	CANE

ELETTROFORESI DELLE SIEROPROTEINE

To change a species on a sample click on the drop-down list button into the field "SPECIES" and select the species you want. To change the species all eight samples in one pass click on the button "SPECIES" and make your selection.

APPENDIX C TROUBLESHOOTING

Problem		Probable cause	Solution
1	Analyzer not linked to the P.C. Host	<p>The analyzer is turned off</p> <p>The analyzer is on and the Busy indicator on the status window is active</p> <p>USB cable is not connected to the P.C. correctly</p>	<p>Turn on the analyzer</p> <p>Reset the analyzer by pressing the Reset button and re-run the software</p> <p>Control the connection of the cable between the analyzer and the P.C. host</p>
2	Samples are not applied correctly on the acetate plate	<p>Serum sample is not in appropriate quantity or it is not uniformly dispensed in the serum plate</p> <p>Tips of applicator are not independent the one from the others</p> <p>Damaged applicator tips</p> <p>Acetate membrane is too wet before the application of samples</p>	<p>Remove the serum samples from the plate and re-add, well distributed, a sufficient quantity (25 microliters each)</p> <p>Clean the applicator assuring the independent moving of each tip</p> <p>Replace damaged tips with new one</p> <p>Increase the drying 1 time on the method parameters configuration (up to 120 sec.); increase the application time (up to 25 sec.)</p>
3	Insufficient migration current	<p>One or both wicks are glided down respect their correct position</p> <p>Insufficient buffer solution level in the electrophoretic chamber</p> <p>Acetate membrane is bad mounted on the plastic frame</p>	<p>Replace the wicks assuring the correct position</p> <p>Refill the buffer solution in both sides of the chamber, up to reach the contact with the wicks (total buffer volume: 150 ml, equally distributed in the sides of the electrophoretic chamber)</p> <p>Assure the correct mounting of the membrane on the plastic frame</p>

Problem		Probable cause	Soluzione
4	Too high migration current (over 10 mA) (Migration at constant voltage)	Wrong dilution of the buffer (too much concentrated) Acetate membrane too much wet before migration	Replace the buffer solution correctly diluted (1:10) Increase the drying time 1 in the method parameter configuration (up to 120 sec.)
5	Too high migration voltage (over 250 volts) (Migration at constant current)	Too dry acetate membrane One or both wicks are glided down respect their correct position	Decrease the drying 1 time on the method parameters configuration (up to 30 sec.); Replace the wicks assuring the correct position
6	Distortion on migrations	One or both wicks are glided down respect their correct position Too high buffer solution level in the electrophoretic chamber Acetate membrane is too wet before the migration	Replace the wicks assuring the correct position Recheck the buffer solution level in both sides of the chamber, up to reach the contact with the wicks (total buffer volume: 150 ml, equally distributed in the sides of the electrophoretic chamber) Increase the drying 1 time on the method parameters configuration (up to 120 sec.)
7	Insufficient staining of the fractions (too low alpha1 fraction)	Exhsausted staining solution Destaining solution too diluted Too low quantity of samples in the wells	Replace the staining solution Replace the destaining solution with well diluted one (dilution is = 1:20) Remove the serum samples from the plate and re-add, a sufficient quantity (25 microliters each) , well distributed
8	Inhomogeneous destaining of the acetate membrane background	Destaining tank empty during the analytical cycle	Check, before the startup of the machine, the presence of enough destaining solution in the external tank (for a complete cycle of 3 membrane, the minimum quantity of 500 ml is required)

9	Truncation of the gamma fraction in the graph	Back migration of the gamma zone fraction during the migration	Decrease, in the params configuration, the deposition line value (up to 130) and re-scan the membrane (before to re-scan, please assure the complete drying of the membrane)
10	Gamma fraction is moving forward respect the application point	Too high migration Tension/Currente at the startup	See points 4 and 5

REMOTE HOST INTERFACE SPECS

Remote Host Specifications

The data relating to a test and to the patient related thereto are sent by a remote computer (host) by sending a text file to the NetDir folder of ELEPHOR 24S management program, located on the local PC connected to the instrument. The coding of the data contained in such file will be explained later in this section.

Requirements: in order to properly interface ELEPHOR 24S analyser with a remote host, first you have to network the PC (by a LAN network, serial connection, etc...) physically connected to the analyser and the remote host. Such operation does not need the information provided herein and is pertaining exclusively to the final user and/or the host program providers. Interfacing the two computers requires the compliance with the standard rules of LAN connections or other management specifications allowing text file exchange and folder sharing by two or more computers (please note that if you have computers running under Windows 2000/XP or Linux, you may need to configurate the user accounts with administrator rights).

In general, you just need to share the NetDir folder of the computer the analyser is connected to and a folder of the computer the host is part of, such folder will be indicated in the ELEPHOR 24Snet program configurations.

Coding of acceptance data:

Acceptance takes place through a text file in the NetDir folder sent by the remote computer, such file will have a fixed name followed by the extension ".txt", typical of text files. (Eg.: Acceptance.txt)

The acceptance/host software shall take exclusive care of the removal of such file from the NetDir folder, the updating of the data contained therein shall be made concurrently by the two software (ELEPHOR 24S and the host program), consequently we recommend the use of methods allowing to synchronise the file access by the acceptance/host software.

Acceptance file data format :

The acceptance file is an unformatted text file the lines of which are the acceptance data of an individual test in the following form:

```
<Name>,<age>,<gender>,<date>,<Idcode>,<TP>,<department>,<Testcode>,<Acceptedflag>,<Sentflag><CR><NL>
```

The following rules are applicable:

- Each line consists of a set of 10 (ten) string type fields separated by a comma
- Each line ends with a character <CR> (Carriage Return) and a character <NL> (New Line)
- Each field must be indicated, for the fields that are not desired just key in a character <SP> (SPace) and maintain the correct sequence as above specified.

Field meaning:

<Name>=Name of the Patient the test is relating to (max 30 characters)

<age>=age (max 3 characters)

<gender>=gender (1 character M/F)

<date>=date related to the test in dd/mm/yy format (this field must be compulsory)

<Idcode>=patient identity code (max 9 characters)

<TotalProteins>=float type value, format 0.00 accounting for total proteins

<department>=name of relevant department (max 25 characters)

<Testcode>=test code (1=Seroproteins,2=Haemoglobins,3=Lipoproteins). This field is COMPULSORY.

<Acceptedflag>=boolean flag indicating the occurred acceptance of the test relating to that patient. It is set at <F> (False) by the host the last time the file is sent and is modified into <T> (True) by ELEPHOR 24S software when such patient is accepted.

<Sentflag>=boolean flag indicating the occurred dispatch by ELEPHOR 24S software of the report data relating to that patient. It is set at <F> (False) by the host the first time the acceptance file is sent and is modified into <T> (True) by ELEPHOR 24S software when the report data are sent.

Typical line examples of an acceptance file, please note that the use of characters <> is only the result of an agreement in order to highlight the used fields and characters, they are not shown in the text file:

Acceptance of Seroproteins

MARIO ROSSI,33,M,15/09/02,0001,6.0,HAEMATOLOGY,1,F,F<CR><NL>

when the test is accepted the line becomes:

MARIO ROSSI,33,M,15/09/02,0001,6.0,HAEMATOLOGY,1,T,F<CR><NL>

when the test is sent to the remote host, the line becomes:

MARIO ROSSI,33,M,15/09/02,0001,6.0,HAEMATOLOGY,1,T,T<CR><NL>

Now, the report data of the patient MARIO ROSSI are available and the remote host may retrieve them from the report file.

For partial acceptance, for example, with no information about age, gender, Total Proteins and department you would obtain the following:

MARIO ROSSI,<SP>,<SP>,15/09/02,0001,<SP>,<SP>,1,F,F<CR><NL>

When the field <Sentflag> becomes <T> the data have been keyed in the report file and may be decoded and processed by the host.

Report file specifications

The report file is organised very similarly to the acceptance file: an unformatted text file containing the patient's data and the report data including the scan data.

Individual report line data format:

```
<Name>,<age>,<gender>,<data>,<Idcode>,<TP>,<department>,<Testcode>,<Acceptedflag>,<Sens  
tflag>,<Testname>,<RefMinPT>,<RefMaxPT>,  
<NumRef>,<RefMinPct1>,...,<RefMinPct NumRef>,<RefMaxPct1>,...,<RefMaxPct NumRef>,  
<NumVal>,<ValPct 1>,...,<ValPct NumVal>,  
<Label1>,...,<Label NumVal>,  
<NumMarks>,<Mark 1>,...,<Mark NumMarks>,  
<NumPoints>,<Point1>,...,<Point numPoints><CR><NL>
```

The same rules of the acceptance file are applicable and the initial fields are identical to those of the acceptance file. Now let us examine the report specific fields:

<Testname>=string indicating the test name as it is shown on the analyser

<RefMinPT>=float string type min reference value for Total Proteins in g/dL (format 0.00)

<RefMaxPT>=float string type max reference value for Total Proteins in g/dL (format 0.00)

<NumRef>=this value indicates the number of elements forming the test reference table, followed by NumRef values being the percentage minimum reference values, followed by NumRef values being the percentage maximum reference values. All values of reference tables are of float string type in format 0.0, so following the <NumRef> field there will be 2*NumRef values forming the percentage reference tables. To retrieve the reference tables expressed in g/dL use the parameters RefMinPT and RefMaxPT with the percentage Min/Max reference tables.

<NumVal>=Number of elements forming the table of the read concentration values, such value takes also into account a fraction partitioning values. Then there are NumVal values being the percentage concentration values.

All concentration values are of float string type, format 0.0

These are followed by NumVal strings <Label1>,...,<Label NumVal> representing the names of fractions and partitioning; please consider that the partitioning names are always preceded by a character <@> (at).

In order to obtain the concentration values expressed in g/dL use the TP value (if available) and the percentage concentration values.

<NumMarks>=Number of elements being fraction marks, followed by NumMarks values of full type representing the abscissa where the fraction marks are shown, such values are associated with the scan data values

<NumPoints>=Number of points being the scan, followed by NumPoints values of full type being the scan curve.

The report file name reflects the rules of the acceptance file, it will be of 'Report.txt' type. The host management program takes care of the removal of such file from the folder.