



## Intraoral Camera



0086

**Gendex Dental Systems**  
901 West Oakton Street  
Des Plaines, IL 60018 U.S.A.  
Tel 1-888-275-5286  
Fax 1-847-550-1322  
[www.gendex.com](http://www.gendex.com)

Authorized Representative:  
**KaVo Dental GmbH**  
Bismarckring 39  
88400 Biberach, Germany  
[www.gendex-dental.com](http://www.gendex-dental.com)  
[www.kavo.com](http://www.kavo.com)

**GENDEX®**  
Imaging Excellence Since 1893

**GENDEX®**  
Imaging Excellence Since 1893

©2009 Gendex Dental Systems, 032-0233P1Rev2

# GXC-300 Intraoral Camera Features

1

The Gendex GXC-300 Intraoral Camera provides live video display, as well as freeze and capture capabilities, for dental imaging.

- Proprietary, seven-element glass optical lens produces exceptional image quality
- Uses a 1/4" Sony CCD sensor with resolution of 410,000 pixels
- Capture control capability directly on the handpiece
- High-definition imaging
- Connects to the PC via a USB port
- Camera disconnects from the PC, allowing portability from PC to PC
- Uses a bright LED light source
- Allows intraoral and extraoral imaging without a lens change
- Lightweight, compact, and simple in design, allowing ease of operation and installation



## Gendex GXC-300 Intraoral Camera Specifications

Name	Gendex GXC-300
Power	DC 5V
Image Sensor	1/4" Sony Super HAD CCD 410,000 pixels
Effective Pixel	768(H) X 494(L)
Horizontal Resolution	480 TV-lines
AGC/ White balance	Auto
Lamp	High luminance white chip LED: 4EA
Cable	Length: 10 ft. (3 m) Connector: USB

## Package Contents



GXC-300 Intraoral Camera



Camera Holster



User Manual




USB Driver  
Installation CD-ROM




Foot Switch  
(Optional)

## Names and Functions of the Camera Parts

This device is an intraoral camera system that provides the dentist and the patient with intraoral and extraoral images. This assists the dentist in describing the dental procedure being performed, as well as allowing the dentist to show the results to the patient.



	Parts	Functions
1	<b>Viewing Head</b> 	Seven-element glass optical lenses
2	<b>Slide Switch</b> 	Slide switch offers a sliding continuous focus mechanism to view intraoral and extraoral images without a lens change. Slide to the tooth icon to capture up to three teeth. Slide to the face icon to capture all other images, including smile, arch, and portrait.
3	<b>Capture Button</b> 	Freeze live images with a click of the Capture button. - Single-touch: Freeze image - Press and hold for 2 seconds or more: Freeze image and transfer it into the imaging software
4	<b>Light Source</b> 	Four bright white LEDs

This page intentionally left blank

## GXC-300 Installation

### Minimum system requirements:

- Microsoft Windows XP SP2 or higher
- DirectX 9.0b or higher
- Pentium 4 processor at 2 Ghz or higher
- 1 GB RAM
- Monitor with minimum 1024 x 768 resolution capability
- Graphics/Display card capable of minimum 1024 x 768 resolution at 24-bit color and 32 MB video RAM
- CD-ROM drive
- USB 2.0 port
- GxPicture 2.62 or later

The installation procedure has two parts:

First - install the GxPicture software

Second - install the intraoral camera

**IMPORTANT!** Be sure to follow the instructions to install the software before connecting the camera. Otherwise the camera will not function properly.

### Install the GxPicture software

1. Note the Gendex products currently installed and using GxPicture:


- Gendex Digital Cameras
- GxTwain
- DenOptix
- DenOptix QST
- VisualiX, GX-S USB, HDI
- VisualiX eHD
- KaVo Dig eXam
- Orthoralix DPI series
- Orthoralix DDE series

2. Remove any earlier versions of GxPicture.
3. Insert the enclosed Installation CD-ROM into the CD-ROM drive.

- On the Welcome screen, click "Next". Carefully read the Software License Agreement and click "Yes" to accept the Software License and continue. You can click "No" to decline and exit the installation.



- On the features setup screen, select "Gendex Digital Cameras" and the other Gendex product noted in step 1.




- Click "Next".


- Select the GXC-300 from the product list and click "Next".



- Select "Do Not Use Filter Engine (for users of VixWin and GXTwain)" and click "Next".



- Select "Yes, I want to restart my computer now" and click "Finish" to complete the software installation.




## Install the Intraoral Camera

1. Connect the USB cable of the camera to an available USB port on your computer.
2. When the Found New Hardware Wizard screen appears, select "No, not this time" and click "Next".
3. Select "Install the software automatically (recommended)" and click "Next".


*Note: If a Windows Logo Testing window appears, click "Continue Anyway".*

4. After the camera is installed, click "Finish".



5. If you purchased the optional Foot Switch, you can now connect it to another USB port on your computer. Refer to the adapter's documentation for more information.
6. Install the camera holster in a place that is convenient and does not exceed the maximum length of the USB cable, which is 10 ft. (3 m). The camera holster should also be positioned so that the cable from the camera does not create a tripping hazard. To install the holster, clean the mounting surface, remove the backing from the double sided-tape on the holster, then firmly press the holster against the mounting surface. After approximately two hours, the holster can be used to store the camera. 

*Note: Remove the protective film covering the glass of the viewing head before using the camera.*



## Intraoral Camera Sheath

Refer to the intraoral camera sheath instructions for installation and replacement procedures.

To order more sheaths, please contact your local sales representative.

**⚠ Caution:** The sheath must be removed after each patient. The intraoral camera must be cleaned and disinfected after each patient by rubbing with isopropyl alcohol (70%) for at least one minute. After cleaning the camera and handpiece, install a new hygienic cover.

After each replacement, check the hygienic sheath for any tears or rips. Replace any damaged sheaths before use.

Avoid contact between the hygienic sheath and any sharp objects, such as the patient's teeth, orthodontia, etc. If there is any form of contact made, inspect the sheath for any tears or rips and replace it if necessary.

The hygienic cover is disposable and should be replaced before each use of the camera.


The hygienic cover is designed only for intraoral use. If you use it for taking extraoral images, the images will not be clear as compared to intraoral use because the transparent film reflects light and changes the focus.

This page intentionally left blank

# 3

## GXC-300 Operation

1. To use the camera, remove it from the holster and it will automatically turn on.



2. Remove the protective film covering the glass of the viewing head.


3. For first-time use, right-click the GxVideoApp icon in the system tray. If this is not the first use of the camera, proceed to step 6.



4. Click "Select Device".



5. Click "Gendex GXC-300 Camera" and click "OK".



6. Left-click the GxVideoApp icon one time to open the Video Application.




Note: By default the image is displayed as a mirrored image. To flip the display, refer to the procedure "Set Capture Properties."

7. Position the camera and then press and release the Capture button to freeze the image on your viewing screen. If you want to save this image, press and hold the Capture button for at least two seconds. This will send the image to the imaging software on the computer.




8. If you have purchased the optional Foot Switch, you can use it in place of the Capture button on the camera. When the Foot Switch is plugged in and one of the two foot pedals is pressed, the Capture button on the camera will be disabled until the Foot Switch is removed and is no longer detected. The blue C pedal is used exactly like the Capture button on the camera. The red M pedal toggles between live video and previously stored images.



**⚠ Caution:** Ensure that the camera is properly placed into the holster with the cord end facing down, as shown. The camera will fail to turn OFF if placed in the holster with the cord end facing up.

9. Place the camera in the camera holster to turn the camera off. Refer to the diagram on the right for proper placement of the camera in the holster.




## Set Capture Properties


1. Select "Preferences" from the Setup menu.




2. Click "Configuration" in the window that appears.



3. Double-click "Capture" to display the capture properties.



4. Click the Video Image tab and select "Flip Vertical" to change the display from a mirrored image to a non-mirrored image.



5. Click OK.
6. Click OK

## Important Safety Precautions

Warning: Failure to follow instructions may result in serious bodily injury or death.

- If the camera gets wet, disconnect it, wipe it down, and then allow it to air dry. Do not attempt to dry with a heater or microwave.
- During severe thunderstorms, stop using the camera. Failure to do so may cause damage to the product.
- Do not disassemble, drop, or otherwise subject the camera to hard shocks; the warranty is null and void in such cases. There is a chance of shock or fire hazard.

## Handling

 Caution: Failure to follow instructions may result in light bodily injury or damage to the product.

- Do not drop or strike the camera since the precision optics could be damaged.
- If the camera is left operational but not in the holster for an extended period of time, the handset will be very warm to the touch.
- When the camera is not in use, place the camera in the camera holster, which will turn it off.
- The intraoral camera system should only be used for dental purposes.
- The camera is not waterproof. Do not immerse in water and avoid situations where water or other foreign matter could get on or inside the product. Failure to do so may damage the product.
- If the camera is wet, or if your hands are wet, do not use the product.
- Avoid subjecting the camera to vibrations or shocks. When the product is not in use, make sure it is placed in the camera holster.
- Be careful not to use damaged products. Doing so may cause injury.
- Avoid using the camera in a dusty environment.
- Remove the protective film covering the glass of the viewing head.

Attention: Federal law restricts the sale of this device to licensed dental practitioners. Consult accompanying documents for more information.



## **Specifications**

### **USB Compatibility**

USB 2.0 and USB 1.1 autoswitching control

### **Electrical Ratings**

5VDC, 500mA

### **Requirements for safety approvals**

This device must be connected to a personal computer that meets the standard specified in IEC 60950, revision 2003.

### **Ambient conditions**

Transportation and storage temperature: -40°C to +70°C (-40°F to +158°F)

Operating temperature: +10°C to +40°C (+50°F to +104°F)

Relative humidity: 10% to 90%, non-condensing

### **Standard Conformity**

The GXC-300 camera complies with the following standards:

UL 60601-1 1st Edition (Medical Electrical Equipment, Part 1: General Requirements for Safety)

CAN/CSA-C22.2 No. 601.1-M90 (Medical Electrical Equipment - Part 1: General Requirements for Safety)

EMC: EN 60601-1-2: Medical electrical equipment. General requirements for safety. Collateral standard. Electromagnetic compatibility.

## Manufacturer's Declaration

This camera, classified as Medical Electrical Equipment, requires special precautions regarding EMC and must be installed and put into service according to the EMC information provided in the accompanying product documentation. Portable and mobile RF communications equipment can effect Medical Electrical Equipment. This camera complies with EMC requirements when used with the cables and accessories supplied with the product. The use of accessories and cables other than those sold by Gendex Dental Systems and specified as replacement parts for internal components, may result in increased emissions or decreased immunity of this camera. This camera should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, this camera should be observed to verify normal operation in the configuration in which it will be used.

Table 201 - Guidance and manufacturer's declaration - electromagnetic emissions for all EQUIPMENT and SYSTEMS

Guidance and manufacturer's declaration – electromagnetic emissions		
This camera is intended for use in the electromagnetic environment specified below. The customer or the user of this camera should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	This camera uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	This camera is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Table 202 – Guidance and manufacturer's declaration – electromagnetic immunity for all EQUIPMENT and SYSTEMS

Guidance and manufacturer's declaration – electromagnetic immunity			
This camera is intended for use in the electromagnetic environment specified below. The customer or the user of this camera should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/-6 kV contact +/-8 kV air	+/-2, 4 & 6 kV contact +/-2, 4 & 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	+/-2 kV for power supply lines +/-1 kV for input/output lines	+/-2 kV for power supply lines Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+/-1 kV differential mode +/-2 kV common mode	+/-0.5 & 1 kV differential mode Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % $U_T$ (>95 % dip in $U_T$ ) for 0,5 cycle  40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles  70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles  <5 % $U_T$ (>95 % dip in $U_T$ ) for 5 sec	Not applicable  Not applicable  Not applicable  Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of this camera requires continued operation during power mains interruptions, it is recommended that this camera be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3A/m	Not applicable	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE $U_T$ is the a.c. mains voltage prior to application of the test level.			

Table 204 – Guidance and manufacturer’s declaration – electromagnetic immunity for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

Guidance and manufacturer’s declaration – electromagnetic immunity			
This camera is intended for use in the electromagnetic environment specified below. The customer or the user of this camera should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2,5 GHz</p>	<p>3 V</p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of this camera, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1,2\sqrt{P}$ $d = 1,2\sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2,3\sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
<p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p><sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which this camera is used exceeds the applicable RF compliance level above, this camera should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating this camera.</p> <p><sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

Table 206 - Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and this camera			
This camera is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this camera can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this camera as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter M		
	150 kHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.  NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			