
CIBS OFFICE

USER MANUAL

TANGO VII

BLUES II

Norsk

Gratulerer.....	5
Om CIBS Office.....	5
Blues II (halsstryngen).....	6
Tango VII (trådløs sender for fasttelefon).....	6
CIBS Office funksjoner.....	8
Lading av Blues II.....	8
Montere Tango VII.....	8
Konfigurere Tango VII.....	9
Slå På/Av Blues II.....	10
Paring av enheter.....	10
Hvordan.....	11
justere volum (Blues II).....	11
justere volum (Tango VII).....	11
justere tonenivået (Tango VII).....	11
svare på en telefonsamtale (Tango VII).....	11
avslutte en telefonsamtale (Tango VII).....	12
ringe fra fasttelefon (Tango VII).....	12
programmere (Blues II).....	12
Betydning av lysdioder og farger.....	13
Fabrikkinnstilling.....	13
Feilsøking og ofte stilte spørsmål.....	14
Oppbevaring og vedlikehold.....	14

Gratulerer

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English

Congratulations.....	5
About CIBS Office.....	5
Blues II (Neck Loop).....	6
Tango VII ((wireless transmitter for desk telephones).....	6
The functions of CIBS Office.....	8
Charging of Blues II.....	8
Installing the Tango VII.....	8
Configure Tango VII.....	9
Switch Blues II on/off.....	10
Pairing the units.....	10
How to.....	11
adjust the volume (Blues II).....	11
adjust the volume (Tango VII).....	11
adjusting the tone level (Tango VII).....	11
answer a telephone call (Tango VII).....	11
end a telephone call (Tango VII).....	12
make a call from your desk phone (Tango VII).....	12
programme (Blues II).....	13
Description of light diodes and colours.....	13
Factory default settings.....	14
Trouble shooting and FAQ.....	14
Storage and maintenance.....	15

Congratulations

Congratulations with your new CIBS Office. We hope you get a lot of pleasure from this solution which gives you better sound from your landline phone.

This manual will help you to get started and ensure that you get the most out of your new solution. Your local dealer offers a range of accessories which will further enhance the functionality of this product.

Before use, the unit must be fully charged (see page 8 for instructions). Charging time is approximately 2 hours.

Be aware!

When connected to a headphones/earbuds the Neck Loop (Blues II) can emit audio of a very high volume. Exposure to this may cause permanent hearing damage. Be careful when increasing the volume to avoid such injuries.

Be careful when using the devices while engaged in activities which demand your full attention.

Never attempt to dismantle the product yourself. Only approved local dealers or service centres may open the products. Contact your local dealer if any parts of the unit need replacing.

Keep out of reach of children:

The plastic packaging of the product and accessories are not children's toys. The packaging or the small parts it contains may cause choking if swallowed.

Ignition of flammable atmospheres

Do not use CIBS Office in environments where there is a danger of igniting combustible gasses.

About CIBS Office

CIBS Office consists of two devices: Blues II (Neck Loop) and Tango VII (wireless transmitter for desk telephones).

Blues II (Neck Loop)

See the illustrations on page 2 for an explanation of the points below

- **Answer button ①:**
Used to answer incoming phone calls.
- **Socket for charging and accessories ②:**
The accompanying power supply cable is connected here, as well as other accessories with 2.5 mm jack.
- **USB socket ③:**
To charge the neck loop by connecting it to a computer
- **Headphone socket ④:**
Socket to connect earbuds with microphone, or standard headphones (accessories).
Used when you are not using a hearing aid.
- **Microphone ⑤**
- **Battery indicator light ⑥:**
Indicates the battery level of the Blues II.
See the table on page 13 for an explanation of the light signals.
- **On/Off button ⑦:**
To switch the unit on and off.
- **Volume control ⑧:**
Push upwards for higher volume and downwards for lower volume.
- **Pairing button ⑨:**
Press to set Blues II in pairing mode.
- **Speed-Dial button ⑩:**
This function is not available when Blues II is connected to Tango VII.

Tango VII (wireless transmitter for desk telephones)

See the illustrations on page 1 for an explanation of the points below

- **Selector for the telephone handset ⑪:**
Press to receive the phone call through the telephone handset.
- **Selector for dialogue mode ⑫:**
This function is not available on Tango VII.
- **Selector for Blues II ⑬:**
Press to receive the phone call through the Blues II.

- **Volume control for the telephone handset 14 a/b:**
Press to adjust the volume of the conversation through the telephone handset.
- **Volume indicator for the telephone handset 14 c:**
Indicates the volume of the conversation through the telephone handset.
- **Configuration button 14 d:**
Press to configure Tango VII, in order to get optimal communication with your telephone.
- **Pairing button 14 e:**
Press to set the Tango VII in pairing mode.
- **Tone adjustment for the telephone handset 15:**
Press to adjust the tone level of the conversation through the telephone handset.
- **Volume control for Blues II 16 a/b:**
Press to adjust the volume of the conversation through Blues II.
- **Volume indicator for Blues II 16 c:**
Indicates the volume of the conversation through Blues II.
- **Microphone sensitivity button 16 d:**
Press to adjust the microphone sensitivity in Tango VII.
- **Tone adjustment for Blues II 17:**
Press to adjust the tone level of the conversation through Blues II.
- **Communication light 18**
Indicates connection between Tango VII and Blues III
- **Input for the telephone handset 19:**
This is where you connect the telephone handset of your desk phone.
- **Input for the desk phone 20:**
This is where you connect your desk phone with the accompanied cable.
- **Input for handset lifter (accessory) 21:**
This is where you connect the handset lifter, which enables you to answer the phone call remotely via your Blues II.
- **Audio output 22:**
This is where you connect e.g. external telecoil amplifier etc.
- **Dialogue input 23:**
Not available on Tango VII

- **Power input 24:**
Connection for the accompanied power supply.

In order to receive sound from the Blues II, you need to put your hearing aid/CI in T or MT position and wear the neck loop around your neck. Alternatively, you can connect earbuds/headphones to the Blues II.

The functions of CIBS Office

CIBS Office offers a wireless solution for your desk phone. By connecting Tango VII to your desk phone and wearing Blues II neck loop, you can enjoy a wireless conversation on the phone, with a range up to 40 meters. Tango VII also amplifies the sound in the telephone handset.

Charging of Blues II

Make sure that Blues II is fully charged before using it.

You can charge Blues II the following ways:

- Using the charger for wall outlet.
- Directly from a PC by USB cable (accessory)
- Using the Twin Charger Cradle (accessory)

When charging, the Battery indicator lamp 6 will flash red.

When Blues II is charged, the Battery indicator lamp 6 will illuminate solid blue.

When the battery is low, the Battery indicator lamp 6 will flash red and at the same time you will hear a sound signal from the Blues II.

NOTE: The lifetime of the battery will be significantly reduced if the device is left uncharged for long periods. We advise you to charge the device at least once a month.

When the operation time of the battery begins to decrease and you wish to replace it, please send it to your local dealer.

The article number of the battery is: 956-001

Installing the Tango VII

Connect the accompanying power supply to the power socket 24. Tango VII will switch itself on automatically. This will be indicated by all the lights illuminating for some seconds. Tango VII does not have an on/off switch, but will automatically go into stand-by mode when not in use.

- EN**
1. Find the cable that connects the telephone receiver to the desk phone. Detach it from the desk phone and connect it to the input for telephone receiver on the Tango VII **19**.
 2. Then, connect the accompanied cable (60cm) to the desk phone where the above cable was detached (1), and to the input for desk phone on Tango VII **20**.
 3. Handset lifter (accessory) can be connected to the input for handset lifter **21**.

Configuration Tango VII

Configuration means adapting the Tango VII to your desk phone. See the Configuration list on the last page.

If your desk phone is listed on the attached Configuration list, apply the configuration digit that corresponds with your phone, and continue to alternative a.

If your desk phone is not listed, continue to alternative b.

Alternative a.

Use a pen or similar pointed object to press and hold configuration button **12d**. The configuration number is indicated by a light in the indicator bar **12c**, with configuration number 01 at the bottom and 08 at the top. Press **12a** or **12b** repeatedly to select the configuration digit of your phone. When the correct digit is selected, release the button **12d**. Make sure that you can hear the dialling tone in your phone and check the microphone by blowing on it. Tango VII is now successfully configured.

Alternative b.

Press button **11** to activate the telephone receiver. Use a pen or similar pointed object to press and hold configuration button **12d**. The configuration number is indicated by a light in the indicator bar **12c**, with configuration number 01 at the bottom and 08 at the top. Lift the telephone handset and press **12a** or **12b** repeatedly until you hear the dialling tone clearly through the receiver. Check that the microphone works by blowing on it. If not, step to the next configuration number and repeat the microphone test. When you both hear the dialling tone and the microphone works, your Tango VII is successfully configured. The configuration button **12d** can then be released. If too much time elapses, the dialling tone will cease and you must replace the handset and lift it again in order to continue configuration.

Configuration of microphone sensitivity (loudness of your voice)

Use a pen or similar pointed object to press and hold button **16d**, while adjusting the sensitivity up or down by pushing button **16a** or **16b**. The light bar **10c** indicated the level of microphone sensitivity.

Switch Blues II on/off

How to switch the Blues II on:

- Press and hold the On/off button **7** until you see the Battery indicator lamp **6** turn blue. Simultaneously you will hear a signal via the T-Loop or via connected headphones/earbuds.

How to switch the Blues II off:

- Press and hold the On/off button **7** until you see the Battery indicator lamp **6** turn red. Simultaneously you will hear a signal via the T-Loop or via connected headphones/earbuds.

Pairing the units

The units must be programmed to be able to communicate, with each other. This is done by connecting these together wirelessly (pairing).

NOTE! It is important that Blues II and Tango VII are within a distance of less than 1 meter during pairing and that Blues II are fully charged.

Pair Blues II with the Tango VII

- Press and hold the pairing button **9** until the Battery indicator light **6** flashes alternately red and blue and you hear a double tone.
- The Blues II is now ready to pair for approximately 3 minutes.
- Press and hold the Tango VII's pairing button **14e** until the communication light **13** flashes alternately red and blue.
- The Tango VII is now ready to pair.
- When Blues II and Tango VII are paired together the communication light **13** will light stable blue.
- If the pairing was unsuccessful, repeat the procedure.

How to...

... adjust the volume (Blues II)

- The volume can be adjusted gradually (press once for each volume level) or by pressing the volume button **8** in the desired direction until the desired volume has been reached.
- You will hear a signal when the volume has reached the highest or lowest level.
- The Blues II remembers its last volume position, so that when switched off and on it will start at the same volume position.

... adjust the volume (Tango VII)

- For the telephone handset
 - Make sure that Tango VII is lit on the volume indicator for the telephone handset **14c**. If not, press the selector for telephone handset **11**
 - Adjust the volume by pressing the volume control for the telephone handset **14a/b**
 - You will see the sound level in the volume indicator for the telephone handset **14c**.
- For the neck loop (Blues II)
 - Make sure that Tango VII is lit on the volume indicator for Blues II **13c**. If not, press the selector for Blues II **13**
 - Adjust the volume by pressing the volume control for Blues II **13a/b**
 - You will see the sound level in the volume indicator for Blues II **13c**.

... adjust the tone level (Tango VII)

- For the telephone handset
 - Make sure that Tango VII is lit for the volume indicator for the telephone handset **14c**. If not, press the selector for telephone handset **11**
 - Adjust the tone level by pressing the tone adjustment for the telephone handset **15**
- For the neck loop (Blues II)
 - Make sure that Tango VII is lit on the volume indicator for Blues II **13c**. If not, press the selector for Blues II **13**
 - Adjust the tone level by pressing the tone adjustment for Blues II **17**

... answer a telephone call (Tango VII)

- Through the telephone handset
 - Make sure that Tango VII is lit for the volume indicator for the telephone handset **14c**. If not, press the selector for telephone handset **11**

- Lift the handset from your telephone.
- You will then hear the call through the telephone handset.

• Through neck loop (Blues II)

- Make sure that Tango VII is lit on the volume indicator for Blues II **13c**. If not, press the the selector for Blues II **13**.
- Press the answer button **1**.
- Lift the handset from your telephone.
- You will then hear the call in Blues II.
- Note! By using a handset lifter (optional), you will hear the ring signal directly in the neck loop, and you can easily answer the call by pressing the answer button **1**. Then you do not need to lift the handset.

... end a telephone call (Tango VII)

- Through the telephone handset
 - Replace the handset of your telephone.
- Through neck loop (Blues II)
 - Replace the handset of your telephone.
 - Press the answer button **1**.
 - Note! By using a handset lifter (optional), you will be able to close the call by pressing the answer button **1**. Then you do not need to replace the handset.

... make a call from your desk phone (Tango VII)

- Through the telephone handset
 - Make sure that Tango VII is lit for the volume indicator for the telephone handset **14c**. If not, press the selector for telephone handset **11**.
 - Lift the handset of your telephone, and dial the desired telephone number.
 - You will then hear the call through the telephone handset.
- Through neck loop (Blues II)
 - Make sure that Tango VII is lit on the volume indicator for Blues II **13c**. If not, press the the selector for Blues II **13**.
 - Press the answer button **1**.
 - Lift the handset of your telephone, and dial the desired telephone number.
 - You will then hear the call in Blues II.
 - Note! By using a handset lifter (optional), the handset lifter will automatically lift the handset when you press the answer button **1**. Then you do not need to lift the handset.

... programme (Blues II)

- Go to www.commidt.com to download both the operating instructions and the programme (CIBS Interface).
- The CIBS Interface programme offers the following options:
 - Turn on or off monitoring voice (enables you to hear your own voice)
 - Locate the serial number

Useful tips:

Using the accompanying shirt/blouse clip can reduce noise caused by the device rubbing off clothing. The clip can be attached to the rear of the Blues II, in the slot on the desired side

By using the external microphone (accessory, article no: 901-001) you can use the Blues II under your clothing, offering even better sound to your audience. These and other microphones are connected at the socket

Description of light diodes and colours

Battery indicator light 6:	
What you see	What this implies for the Blues II
Solid blue light	The battery level is between 100% - 20%
Flashing blue light	Incoming call (only when using the handset lifter, accessory)
Flashing red light	The battery level is below 20%
Solid blue light (while charging)	Charging complete
Flashing red light (while charging)	The Blues II is charging
Alternate blue and red light	The Blues II is ready to pair

Factory default settings

- This function will delete the stored pairing to the Blues II. Any settings made after you received the device will also be deleted.
- Press and hold the pairing button 7. First you will hear a double tone which indicates the device is set to pairing mode. You must keep the button pressed until you hear the second double tone.
- The flashing red and blue light will still flash
- Switch off the Blues II. When you switch the device on again, it will automatically go into pairing mode (the battery indicator light 6 will flash alternately red and blue). This indicates that the device has been reset to factory settings.

Trouble shooting and FAQ

The Blues II will not switch on!

- Charge the battery

The Blues II will not charge!

- Make sure the charge plug is connected correctly in the Blues II and that the power supply is working.

I hear crackling!

- For best possible quality the Blues II must be within range of the Tango VII.
- If you are using headphones, make sure the plug to the headphones is correctly connected.

I cannot hear anything from the Blues II!

- Increase the volume of the Blues II.
- Check that the Blues II is paired to the Tango VII.
- If you are using headphones, make sure the plug to the headphones is correctly connected. If you are using the hearing device, make sure the coil is lying correctly around the neck and that the hearing device is set to T or MT.

Will the Blues II and Tango VII work with other Bluetooth devices?

- The Blues II and Tango VII is made to work with most Bluetooth products. This depends on the Bluetooth product being compatible with Bluetooth version 1.1, 1.2 or 2.0, supporting headsets, handsfree and/or A2DP profile.

Storage and maintenance

- Always store the Blues II in a safe place with the power switched off.
- Do not store the Tango VII or Blues II in extreme temperatures (over 45°C/113°F, in direct sunlight, or below -10°C/14°F).
- Do not expose the Tango VII or Blues II to rain or other liquids.
- Clean the devices with a dry cloth.

Certification and safety approvals


CE

This product is CE marked accordingly to the provisions of the R & TTE Directive (99/55/EC), COMMidt AS, declares that this product is in compliance with the essential requirements and other relevant provisions of Directive 1995/5/EC.

For further information, please consult: <http://www.commidt.com>

Bluetooth

The Bluetooth word® and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by COMMidt AS is under licence. Other trademarks and trade names are those of their respective owners.

Electrical and electronic equipment (EEE) contains materials, components and substances which can be dangerous and detrimental to human health and the environment if waste electrical and electronic equipment (WEEE) is not disposed of correctly. Products marked with the “crossed-out waste bin” symbol shown below are electrical and electronic equipment.

This symbol indicates that waste electrical and electronic equipment must not be discarded along with unsorted household waste, but must be collected separately. In the development and manufacture of your COMMidt product, high-quality materials and components have been used which can be recycled. When disposing of electrical and electronic equipment you should therefore make use of a collection arrangement intended for this purpose, which may involve either delivery at your local collection point or recycling station or collection at your home. Further information should be obtained from your local dealer.

Batteries contain substances that may be harmful to human health and the environment if not handled properly.

Batteries are labelled with a crossed dustbin as shown.

It symbolizes that waste batteries must not be disposed of as unsorted municipal waste, but must be collected separately.

Some batteries are also labelled with the chemical symbol Hg (mercury), Cd (cadmium) or Pb(lead). This is particularly harmful substances, and therefore it is important that these batteries will be collected.

It is important that you hand your waste batteries to the collection systems which are established. On this way you help to ensure that batteries are recycled in accordance with the law and will not unduly impact the environment.

When disposing of waste batteries you must use a collection system established for that purpose whether it be handed in at your local collection point or recycling or collected directly from the household. Detailed information must be obtained from the local authorities.

Removal of the battery:

1. Unscrew the screw on the rear of the unit
2. Remove the front cover
3. Unscrew the screw on the circuit board.
4. Remove the circuit board
5. Unplug the battery

Declaration of conformity

We COMMIdt, Kirkegata 57-59, 7600 Levanger, Norway declare under our own responsibility that the Bluetooth products to which this declaration refers conforms with the relevant standards or other standardised documents.

- 93/42/ EØF
- NS-EN 12182
- ETSI EN 301 489-1 V1.4.1 (2002-08)
- ETSI EN 301 489-17 V1.2.1 (2002.08)
- EN 300 328 V1.6.1 (2004-11)
- EN 61000-4-3, 5.1
- EN 61000-4-2, 5.2

According to the regulations in Radio Equipment and Telecommunications Terminal Equipment (R&TTE) Directive 1999/5/EC

This device complies with Part 15 of the FCC Rules.

Operation is subject to the following two conditions:

- (1) this device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

This device is a class 2 Bluetooth device (max 2 mW EIRP). Due to the low level of transmit output power the device is not subject to the routine RF exposure valuation as per Section 2.1093 of the FCC rules.

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