

# ACCU-CHEK<sup>®</sup> *Spirit*


**INSULIN PUMP**

**User Guide**



---

## The Accu-Chek Spirit – At a Glance



This user's guide is valid for version 2.XX of the Accu-Chek Spirit.

---

---

Dear Accu-Chek Spirit insulin pump user,

Congratulations on your new Accu-Chek Spirit insulin pump. You and your Accu-Chek Spirit insulin pump are unique. Therefore, it is important that you, your doctor or healthcare team fill in the following information:

Serial number of your Accu-Chek Spirit insulin pump \_\_\_\_\_

Insulin (name/type) \_\_\_\_\_

Accu-Chek Spirit Training Place \_\_\_\_\_

Date(s) \_\_\_\_\_

\_\_\_\_\_

Trainer \_\_\_\_\_

**Useful addresses**

Your contact for medical emergencies and/or information about insulin pump therapy (doctor/healthcare team)

\_\_\_\_\_

\_\_\_\_\_

Phone \_\_\_\_\_

---

Your contact for assistance with and/or additional information about your Accu-Chek Spirit insulin pump:

**United Kingdom:**

Roche Diagnostics Ltd.  
Charles Avenue  
Burgess Hill  
West Sussex, RH15 9 RY  
United Kingdom  
Phone +44 1444 256 000  
lewes.insulinpumps@roche.com  
www.accu-chek.co.uk

**Australia:**

Roche Diagnostics Australia Pty. Ltd.  
31 Victoria Avenue  
Castle Hill, NSW 2154  
Australia  
Phone +61 (0) 2 98602100 (Pump Hotline)  
1800 802 409 (Orders)  
Fax 1800 066 598 (Orders)  
sydney.accu-chekextra@roche.com  
www.accu-chek.com.au

**New Zealand:**

Roche Diagnostics N.Z. Ltd.  
P.O. Box 62-089  
15 Rakino Way  
Mt. Wellington, Auckland  
New Zealand  
Accu-Chek Enquiry Line: 0800 80 22 99  
Phone +64 9 276 4157  
Fax +64 9 276 8917  
info@accu-chek.co.nz  
www.accuchek.co.nz

**South Africa:**

Roche Products (Pty) Ltd.  
Diagnostics Division  
P.O. Box 1927, Randburg 2125  
9 Will Scarlet Road  
Randburg  
South Africa  
Phone +27 11 504 4600  
Fax +27 11 781 0269  
www.diabetes.co.za

---

<b>1</b>	<b><i>Introduction</i></b> . . . . .	<b>9</b>
1.1	<b>Who this user guide is written for</b> . . . . .	<b>9</b>
1.2	<b>How to use this user guide</b> . . . . .	<b>10</b>
1.3	<b>Working safely with your Accu-Chek Spirit insulin pump</b> . . . . .	<b>10</b>
1.4	<b>Intended use of your Accu-Chek Spirit insulin pump</b> . . . . .	<b>11</b>
<b>2</b>	<b><i>Before you begin</i></b> . . . . .	<b>13</b>
2.1	<b>Warnings and cautions for proper use and your safety</b> . . . . .	<b>13</b>
2.1.1	Set-up warnings and cautions . . . . .	13
2.1.2	Warnings and cautions regarding hygiene . . . . .	15
2.1.3	Warnings and cautions on proper use . . . . .	16
2.1.4	Warnings and cautions for secure use . . . . .	18
2.1.5	Warnings and cautions about environmental hazards . . . . .	20
2.1.6	Warnings and cautions regarding safe battery use . . . . .	21
2.1.7	Warnings and cautions about wearing the insulin pump . . . . .	22
2.2	<b>Guarantee</b> . . . . .	<b>23</b>
<b>3</b>	<b><i>The Accu-Chek Spirit insulin pump</i></b> . . . . .	<b>25</b>
3.1	<b>Display</b> . . . . .	<b>25</b>
3.2	<b>Backlight</b> . . . . .	<b>26</b>
3.3	<b>Keys and key combinations</b> . . . . .	<b>26</b>
3.4	<b>KeyLock</b> . . . . .	<b>29</b>
3.5	<b>Scrolling</b> . . . . .	<b>30</b>
3.6	<b>STATUS screen</b> . . . . .	<b>31</b>
3.7	<b>Alarm signals</b> . . . . .	<b>31</b>
3.8	<b>RUN and STOP modes</b> . . . . .	<b>32</b>
3.9	<b>STOP-Warning</b> . . . . .	<b>33</b>
<b>4</b>	<b><i>Accessories, disposables and software</i></b> . . . . .	<b>35</b>
4.1	<b>Sterile products</b> . . . . .	<b>36</b>
4.1.1	Cartridge . . . . .	37
4.1.2	Infusion sets . . . . .	38
4.2	<b>Adapter</b> . . . . .	<b>39</b>
4.3	<b>Battery</b> . . . . .	<b>40</b>
4.4	<b>Battery cover</b> . . . . .	<b>41</b>
4.5	<b>Battery key</b> . . . . .	<b>42</b>
4.6	<b>Emergency kit</b> . . . . .	<b>43</b>
4.7	<b>Software</b> . . . . .	<b>44</b>

---

<b>5</b>	<b><i>Prepare the pump for use</i></b> . . . . .	<b>45</b>
5.1	<b>Insert and change the battery</b> . . . . .	<b>45</b>
5.2	<b>Start-Up procedure</b> . . . . .	<b>48</b>
5.2.1	Prepare the cartridge . . . . .	52
5.2.2	Fill the cartridge . . . . .	54
5.3	<b>Insert the cartridge, adapter and infusion set</b> . . . . .	<b>58</b>
5.3.1	Connect the cartridge, adapter and infusion set . . . . .	58
5.3.2	Insert the cartridge . . . . .	60
5.4	<b>Prime the infusion set</b> . . . . .	<b>64</b>
5.5	<b>Prepare the infusion site</b> . . . . .	<b>66</b>
5.6	<b>Change the infusion set</b> . . . . .	<b>69</b>
5.7	<b>Change the cartridge and infusion set</b> . . . . .	<b>71</b>
5.8	<b>Change the adapter</b> . . . . .	<b>73</b>
5.9	<b>How to wear your pump</b> . . . . .	<b>73</b>
<b>6</b>	<b><i>Select a user menu (STANDARD, ADVANCED or CUSTOM)</i></b> . . . . .	<b>75</b>
6.1	<b>The three levels of the Accu-Chek Spirit</b> . . . . .	<b>75</b>
6.2	<b>Navigate through the menus</b> . . . . .	<b>76</b>
6.2.1	Scrolling . . . . .	77
6.2.2	Looping . . . . .	77
6.2.3	Backward loop . . . . .	77
6.2.4	Exit options . . . . .	77
6.2.5	Select a user menu . . . . .	78
<b>7</b>	<b><i>STANDARD user menu</i></b> . . . . .	<b>83</b>
7.1	<b>Set time and date</b> . . . . .	<b>83</b>
7.2	<b>Basal rate</b> . . . . .	<b>87</b>
7.2.1	Program a basal rate profile . . . . .	88
7.2.2	Copy an hourly basal rate . . . . .	90
7.2.3	Basal rate and bolus . . . . .	91
7.3	<b>Start insulin delivery</b> . . . . .	<b>92</b>
7.4	<b>Stop insulin delivery</b> . . . . .	<b>93</b>
7.5	<b>Program a bolus</b> . . . . .	<b>95</b>
7.5.1	Standard bolus . . . . .	96
7.5.2	Program a “Quick” Standard bolus . . . . .	97
7.5.3	Cancel a “Quick” Standard bolus . . . . .	99
7.5.4	Program a menu-guided “Scroll” Standard bolus . . . . .	100
7.5.5	Cancel a menu-guided “Scroll” Standard bolus . . . . .	101

---

<b>7.6</b>	<b>Temporary basal rate (TBR)</b> . . . . .	<b>102</b>
7.6.1	Program a temporary basal rate . . . . .	103
7.6.2	Cancel a temporary basal rate . . . . .	105
<b>7.7</b>	<b>Setup menu standard</b> . . . . .	<b>107</b>
7.7.1	KeyLock function . . . . .	107
7.7.2	Adjust the beep volume . . . . .	109
7.7.3	Alarm signals . . . . .	110
7.7.4	Automatic off . . . . .	111
7.7.5	Battery type . . . . .	112
7.7.6	Display orientation . . . . .	113
<b>7.8</b>	<b>Review data memory (INFORMATION)</b> . . . . .	<b>114</b>
7.8.1	Review the bolus history . . . . .	114
7.8.2	Review the alarm history . . . . .	116
7.8.3	Review the daily insulin totals history . . . . .	116
7.8.4	Review the temporary basal rate history . . . . .	117
7.8.5	Review the time remaining . . . . .	118
7.8.6	Quick info screen . . . . .	119
<b>8</b>	<b>ADVANCED user menu</b> . . . . .	<b>121</b>
<b>8.1</b>	<b>Extended bolus</b> . . . . .	<b>122</b>
8.1.1	Program an Extended bolus . . . . .	123
8.1.2	Cancel an Extended bolus . . . . .	124
<b>8.2</b>	<b>MultiWave bolus</b> . . . . .	<b>126</b>
8.2.1	Program a MultiWave bolus . . . . .	127
8.2.2	Cancel a MultiWave bolus . . . . .	129
<b>8.3</b>	<b>Basal rate profiles</b> . . . . .	<b>131</b>
8.3.1	Set an additional basal rate profile . . . . .	131
8.3.2	Select a basal rate profile . . . . .	132
<b>8.4</b>	<b>Alarm clock</b> . . . . .	<b>133</b>
<b>8.5</b>	<b>Setup menu advanced</b> . . . . .	<b>134</b>
8.5.1	Time format . . . . .	134
8.5.2	Date format . . . . .	135
8.5.3	Bolus increment . . . . .	136
8.5.4	Prime quantity . . . . .	137
8.5.5	Lock basal rate profiles . . . . .	138
8.5.6	Language . . . . .	139
8.5.7	Display contrast . . . . .	140

---

- 9 Data transfer . . . . . 141**
- 10 Alerts and errors . . . . . 145**
  - 10.1 Alerts . . . . . 147**
    - 10.1.1 Alert A1: CARTRIDGE LOW . . . . . 147
    - 10.1.2 Alert A2: BATTERY LOW . . . . . 148
    - 10.1.3 Alert A3: REVIEW TIME AND DATE . . . . . 148
    - 10.1.4 Alert A4: ALARM CLOCK . . . . . 149
    - 10.1.5 Alert A5: PUMP TIMER . . . . . 149
    - 10.1.6 Alert A6: TBR CANCELLED  
(Temporary Basal Rate cancelled) . . . . . 150
    - 10.1.7 Alert A7: TBR OVER (Temporary Basal Rate over) . . . . . 151
    - 10.1.8 Alert A8: BOLUS CANCELLED . . . . . 151
  - 10.2 Errors . . . . . 152**
    - 10.2.1 Error E1: CARTRIDGE EMPTY . . . . . 152
    - 10.2.2 Error E2: BATTERY DEPLETED . . . . . 153
    - 10.2.3 Error E3: AUTOMATIC OFF . . . . . 153
    - 10.2.4 Error E4: OCCLUSION . . . . . 154
    - 10.2.5 Error E5: END OF OPERATION . . . . . 157
    - 10.2.6 Error E6: MECHANICAL ERROR . . . . . 158
    - 10.2.7 Error E7: ELECTRONIC ERROR . . . . . 159
    - 10.2.8 Error E8: POWER INTERRUPT . . . . . 160
    - 10.2.9 Error E10: CARTRIDGE ERROR . . . . . 161
    - 10.2.10 Error E11: SET NOT PRIMED . . . . . 162
    - 10.2.11 Error E12: DATA INTERRUPTED . . . . . 163
    - 10.2.12 Error E13: LANGUAGE ERROR . . . . . 164
- 11 Troubleshooting . . . . . 165**
- 12 Life with your Accu-Chek Spirit insulin pump . . . . . 169**
  - 12.1 Recommendations for daily use . . . . . 169**
    - 12.1.1 Short interruption of insulin pump therapy . . . . . 169
    - 12.1.2 Continuation of insulin pump therapy . . . . . 170
    - 12.1.3 Interruption of insulin pump therapy for  
a longer period of time . . . . . 170
  - 12.2 Weather conditions . . . . . 171**



---

<b>12.3 Your insulin pump and water</b> . . . . .	<b>172</b>
12.3.1 Daily situations. . . . .	172
12.3.2 Accidental contact with water. . . . .	173
12.3.3 What to do after water contact . . . . .	173
12.3.4 Other liquids. . . . .	174
<b>13 On the go with your Accu-Chek Spirit insulin pump</b> . . . . .	<b>175</b>
<b>13.1 Electromagnetic fields and hazardous areas</b> . . . . .	<b>175</b>
<b>13.2 Sports</b> . . . . .	<b>176</b>
<b>13.3 Travelling</b> . . . . .	<b>176</b>
<b>14 Care of your Accu-Chek Spirit insulin pump</b> . . . . .	<b>177</b>
<b>14.1 System check</b> . . . . .	<b>177</b>
<b>14.2 Maintenance and cleaning</b> . . . . .	<b>179</b>
14.2.1 Cleaning your Accu-Chek Spirit insulin pump . . . . .	179
14.2.2 Battery information. . . . .	179
<b>14.3 Storing your Accu-Chek Spirit insulin pump</b> . . . . .	<b>180</b>
<b>14.4 When your insulin pump is dropped</b> . . . . .	<b>181</b>
<b>14.5 Repair</b> . . . . .	<b>182</b>
<b>14.6 Disposal</b> . . . . .	<b>183</b>
<b>15 Technical data</b> . . . . .	<b>185</b>
<b>15.1 General technical data</b> . . . . .	<b>185</b>
<b>15.2 Technical standards on electromagnetic emissions</b> . . . . .	<b>190</b>
<b>15.3 Technical standards on electromagnetic immunity</b> . . . . .	<b>191</b>
<b>15.4 Configuration parameters</b> . . . . .	<b>196</b>
<b>16 Annexes</b> . . . . .	<b>199</b>
<b>16.1 Abbreviations</b> . . . . .	<b>199</b>
<b>16.2 Beeps and melodies</b> . . . . .	<b>201</b>
<b>16.3 Vibrations</b> . . . . .	<b>202</b>
<b>16.4 Symbols</b> . . . . .	<b>203</b>
16.4.1 General symbols. . . . .	203
16.4.2 Display symbols . . . . .	205
<b>16.5 Sterile products and accessories</b> . . . . .	<b>208</b>
16.5.1 Sterile products . . . . .	208
16.5.2 Accessories . . . . .	209
<b>17 Glossary</b> . . . . .	<b>211</b>
<b>18 Index</b> . . . . .	<b>217</b>
<b>19 Alert overview</b> . . . . .	<b>223</b>

---

---

# 1 Introduction

## 1.1 Who this user guide is written for

You are a highly motivated diabetes patient who is willing to work with your doctor and/or healthcare team on insulin therapy. From the start it is very important that you

- ▶ take the time to train with a healthcare professional on the insulin pump
- ▶ take the time for frequent testing of your blood glucose levels, for example using an Accu-Chek blood glucose meter in order to help you find your personal settings quickly and avoid incorrect dosages at an early stage

This user guide is written for you as well as for healthcare professionals. It is designed to provide you with the necessary information for safe and efficient use of your Accu-Chek Spirit insulin pump. Regardless of your experience with pump therapy, **please read this user guide carefully before using your Accu-Chek Spirit insulin pump.** Always refer to this user guide for questions about operating and troubleshooting your Accu-Chek Spirit insulin pump.

Roche Diagnostics is proud to work with you in the management of your diabetes therapy. However, remember that Roche Diagnostics does not make any recommendations about your diabetes therapy, such as how to set up your personal programming or define what features are appropriate. Always follow the instructions given by your doctor or healthcare team. Before using your Accu-Chek Spirit insulin pump, you must be trained by your doctor and/or healthcare team on insulin pump therapy as well as on your Accu-Chek Spirit insulin pump.

We also recommend that you have someone around you (e.g. family, friends) who understands diabetes and pump therapy with your Accu-Chek Spirit insulin pump to help you in case of an emergency.

Should you have any questions, we urge you to contact your local Customer Care Area, your doctor or healthcare team. In this way you can be assured of a trouble-free experience with the Accu-Chek Spirit insulin pump.

## 1.2 How to use this user guide

**Note** The Accu-Chek Spirit user guide presents you with sample screens only. The screens in this user guide may vary slightly from your actual pump screens.

At the back of this user guide you will find various references in the form of tables and definitions.

## 1.3 Working safely with your Accu-Chek Spirit insulin pump

Your Accu-Chek Spirit insulin pump meets the safety requirements of your national medical device legislation and complies with or exceeds the international standards for Electromagnetic Compatibility regarding its use. Section 1.4 introduces you to the intended use of the Accu-Chek Spirit insulin pump, its safety features, and to warnings and precautions you should be aware of while using your insulin pump.

When your insulin pump is in use, its two microprocessor safety systems continuously monitor and control the insulin pumps functioning. Over 9 million safety checks are performed daily. If your insulin pump detects a deviation of its normal state, an alert (warning instruction) or an error (error message) will occur.

See the sections 10 “Alerts and errors” and 15 “Technical data” for further information on how to address these issues.

**Note** Your insulin pump is a valuable medical device. We strongly encourage you to add your insulin pump to your homeowner’s insurance policy to protect you in case of accidental damage or loss. Ask your insurance agent for details.

## **1.4 Intended use of your Accu-Chek Spirit insulin pump**

Your Accu-Chek Spirit insulin pump is a prescription device and has been developed exclusively for the continuous subcutaneous delivery of U100 short-acting insulin or fast-acting insulin analogue. The insulin type for the treatment of your insulin dependent diabetes mellitus will be prescribed by your doctor or healthcare team.

Do not use your insulin pump for the delivery of medications other than U100 short-acting insulin or fast-acting insulin analogue.

---

## 2 Before you begin

### 2.1 Warnings and cautions for proper use and your safety









To ensure that there is no risk to you or your health, you should be aware of the following warnings regarding set-up, hygiene and secure use prior to using your Accu-Chek Spirit insulin pump.

Please read through these warnings and cautions before starting. Warnings, cautions and other important safety information can be found in this section and are clearly visible throughout the user guide.

#### 2.1.1 Set-up warnings and cautions

- ➔ Before using a new Accu-Chek Spirit insulin pump, **always check your personal settings before use in order to avoid incorrect insulin delivery**. Record your personal settings from your current insulin pump and ensure that they are correctly programmed into a new Accu-Chek Spirit insulin pump. Make sure that the time and date are set correctly. If you are not sure about your personal settings or you are less experienced, contact your doctor or healthcare team to check your personal settings.
- ➔ **Do not interrupt the start-up procedure by pressing keys or by any handling on the insulin pump**. Interruption of the start-up procedure may lead to malfunctions of your insulin pump.
- ➔ **Incorrect programming** of your insulin pump may cause inappropriate insulin delivery. Your insulin pump must be programmed with your personal settings prior to starting insulin pump therapy. Do not operate your insulin pump without knowing your personal settings. If you are not sure about your personal settings or you are less experienced, have your doctor or healthcare team check them. Work with your doctor or healthcare team to

determine the timing as well as the amount and type of bolus you need to deliver. Make sure you know your personal carbohydrate insulin ratio and your correction bolus ratio.

- ➔ **When your insulin pump is in STOP**, it does not deliver any insulin. Put your insulin pump into **RUN** to continue the insulin delivery.
  
- ➔ After inserting a new battery, **always ensure that the time and date of your insulin pump are set correctly**. Incorrect programming of the time and date may cause incorrect insulin delivery.
  
- ➔ **If your insulin pump “times out” before you finish programming a bolus, no bolus is delivered**. Check the bolus history screens on your insulin pump for the bolus amount delivered with time and date and program a new bolus if necessary.
  
- ➔ **If you do not set the correct battery type** on your insulin pump, the alert A2: BATTERY LOW may not occur in time to provide you sufficient warning to replace your battery.
  
- ➔ **If you flip the display orientation of your insulin pump by 180°, the  and  keys will also reverse their functions**. In comparison to the display orientation, the upper key will be  and the lower key will be . This change applies for all  and  functions including turning on the backlight. The  and  keys will not change in their function, but will remain the same, regardless of your screen orientation.
  
- ➔ **If data transfer between your PC and your insulin pump is disrupted in any way**, the configuration may be incomplete and an error E12: DATA INTERRUPTED occurs. The data transfer must be completed successfully before you can put your insulin pump into **RUN**.



## 2.1.2 Warnings and cautions regarding hygiene

- ➔ Your insulin pump was designed to use an Accu-Chek 3.15 ml plastic cartridge. **The cartridge is a sterile product intended for single use only.** Sterility is guaranteed for unopened packaging up to the indicated expiration date. **Do not use sterile products if the package is damaged.**
  
- ➔ **Do not reuse single use materials.** Reuse of single use materials might lead to a malfunction of your insulin pump and to incorrect insulin delivery and/or an infection.
  
- ➔ **Always handle the items with clean hands.** Avoid any contact of the infusion set and especially the connecting parts of your insulin pump with antiseptics, antibiotic creams, soaps, perfumes, deodorants, body lotions or any other cosmetics. They may contaminate these parts.
  
- ➔ **Avoid deliberate contact with water.** In case of deliberate contact with water disconnect and take off your insulin pump.
  
- ➔ **Check daily that your insulin pump and its sterile products are not chipped, cracked or damaged in any way and that the battery cover and the adapter are correctly closed.** In the presence of chips and cracks, water, dust, insulin, or other foreign substances may enter your insulin pump and lead to malfunction.
  
- ➔ **Only use sterile products and accessories that are designed for use with your insulin pump.** Other sterile products and accessories have not been tested for compatibility with your insulin pump and may therefore endanger your health if used.

➔ **Always have extra sterile products and accessories with you.**

This allows you to exchange components when needed. Materials designed for single use (such as cartridges and infusion sets) may not be reused due to the increased risk of infections, malfunctions and/or incorrect insulin delivery.



➔ **Always remove the cartridge and adapter and place your insulin pump in STOP while cleaning.** Avoid pressing the keys of your insulin pump during cleaning, as this may accidentally change your settings. Do not use alcohol, solvents, strong detergents, bleaching agents, scouring pads or sharp instruments (knife, key, screw driver, etc.) for cleaning as they may damage your insulin pump.

### 2.1.3 Warnings and cautions on proper use

➔ **Training and use of your insulin pump requires the support of an experienced doctor or healthcare team.** Regular visits with your doctor or healthcare team are absolutely essential during insulin pump therapy. Change your personal settings only after consulting your doctor or healthcare team. Always follow the instructions given by your doctor or healthcare team.

➔ **Do not perform any servicing or repair on your insulin pump by yourself.** Do not use any lubricant on the insulin pump mechanism. For additional questions contact your local Customer Care Area.

➔ **Never press any of the keys with a sharp or pointed item such as the tip of your fingernail.** This could damage the casing of your insulin pump. Use the pad of your finger tip instead to prevent puncturing the keys.

- 
- **Check the amount of insulin remaining in the cartridge at least once a day.** Before going to bed, make sure that the cartridge contains enough insulin to last through the night.
  
  - **Contact your doctor or healthcare team for an alternate therapy plan** when you interrupt the insulin pump therapy for a longer period of time.
  
  - **The bolus increment that can be programmed into your insulin pump determines the “Quick” Standard bolus amount using the  and  keys of your insulin pump.** To ensure correct insulin delivery be sure the bolus increment is set appropriate for your therapy.
  
  - **Do not use the Accu-Chek Spirit insulin pump if you are not able to read the display, hear the beeps or you are not able to feel the vibrations.** The display, beeps and vibrations indicate critical alert information and must be addressed immediately by the insulin pump user. Failure to see, hear or feel the alerts from the insulin pump could result in serious harm or death of the user.
  
  - **If insulin administration is interrupted for any reason** (e.g. stopping the pump by the user, technical problem of the insulin pump, leakage in the cartridge, occlusion of the infusion set tube or of the infusion set needle, infusion set needle has slipped out of the infusion site), you must be prepared to replace the missing insulin immediately.
  
  - **Always carry spare sterile products and accessories (infusion set, insulin cartridge, batteries) as well as an insulin pen/syringe and insulin with you.** Without insulin, diabetic ketoacidosis may develop and this may require in-patient hospital treatment.

## 2.1.4 Warnings and cautions for secure use

- ➔ **The interruption of insulin delivery (e.g. due to a leak, occlusion or loss of insulin potency) or the malfunction of your insulin pump can result in a rapid rise of your blood glucose level.** Although your insulin pump has an internal security check system, it cannot alert you if your infusion set is leaking or if the insulin you are using has lost its potency.
  
- ➔ **Do not use any adapter from another insulin pump** (e.g. H-TRON / H-TRONplus adapters) for your Accu-Chek Spirit.
  
- ➔ **If an error E4: OCCLUSION occurs, immediately check your blood glucose level,** because insulin delivery has been interrupted. If your blood glucose level is high, take appropriate actions according to your doctor's or healthcare team's instructions.
  
- ➔ **Properly tighten the infusion set to the adapter to avoid leakage.** Turn the infusion set luer-lock connector only as far as to the stop. Do not turn it any further and do not use any auxiliary tools as this may crack the infusion set luer-lock connector and result in leakage.
  
- ➔ Your insulin pump cannot detect infusion set leakage. **You must inspect all parts of your infusion set at least every three hours during the day and before you go to sleep.** Should you detect any loss of insulin and all parts are correctly tightened, immediately replace the leaking component. Immediately check your blood glucose level, because insulin delivery has been interrupted.
  
- ➔ **Air bubbles in the cartridge and infusion set may cause the infusion of air instead of insulin.** Your body does not receive the required amount of insulin. An error E4: OCCLUSION may be delayed. Remove these bubbles while filling the cartridge and priming the infusion set and without having the infusion set connected to your body.

- ➔ Cold insulin may release air when warmed. **Use only insulin at room temperature when filling the cartridge and priming the infusion set.** Inspect the cartridge and the infusion set for air bubbles at least every three hours during the day and before you go to sleep. Remove any air bubbles and, if necessary, change system components.
  
- ➔ **Only use your insulin pump for insulin delivery if it is fully and correctly equipped** (including cartridge, adapter and infusion set) and programmed with your personal settings approved by your doctor or healthcare team.
  
- ➔ **Do not carry out the CHANGE THE CARTRIDGE function if the cartridge compartment is not completely dry.** During the rewinding of the piston rod liquid might enter your insulin pump and could lead to malfunction.
  
- ➔ **Programming your insulin pump for a relatively low basal rate (0.1 U/h) can result in a delayed error E4: OCCLUSION.** Roche Diagnostics recommends to use plastic cartridges when a low basal rate is required for your therapy.
  
- ➔ **Never prime an infusion set that is connected to your body. You risk uncontrolled insulin delivery into your body.** With disconnectable infusion sets, make sure you have disconnected the tubing from your site prior to priming. Always follow the instructions provided with your infusion set.
  
- ➔ **To prevent free flow of insulin move the piston rod to the proper position and twist the adapter until the cartridge plunger sits flush with the end plate of the piston rod.** Free flow of insulin from the cartridge or infusion set can occur when the cartridge plunger and the piston rod are not properly connected and you position your insulin pump at a higher level than the infusion site.

- ➔ **If the piston rod does not return completely to its base after the cartridge has been removed**, immediately contact your local Customer Care Area for support.
- ➔ **After the occurrence of an error, your insulin pump will be in STOP** and the insulin delivery is interrupted. In order to maintain insulin delivery, you must act immediately according to the instructions given for each error code. Put your insulin pump into **RUN** to continue the insulin delivery.
- ➔ **Your insulin pump and its sterile products and accessories include small component pieces that could pose a choking hazard to children.** Keep sterile products and accessories out of reach of children.

### 2.1.5 Warnings and cautions about environmental hazards

- ➔ **Avoid electromagnetic fields** of radar or antenna installations, high-voltage sources, X-Ray sources, MRI, CAT scan or other sources of electrical current. Do not use your insulin pump in such areas. Electromagnetic fields may cause malfunction to your insulin pump. Always stop and remove your insulin pump prior to entering these areas. In other cases, the insulin delivery may be stopped immediately and an error E7: ELECTRONIC ERROR occurs. See the section 10.2.7 “Error E7: ELECTRONIC ERROR” for further information.
- ➔ **Do not use your insulin pump in hyperbaric chambers and in hazardous areas of any classification** (such as areas where explosive or flammable gases or vapours could exist) as this might interfere with insulin delivery and/or lead to harmful situations.

- ➔ **Your insulin pump is designed to work in normal barometric conditions from 70 to 106 kPa (700 to 1060 mbar).** Do not exceed 3000 meters (10000 feet) above sea level. It has not been tested for use in hazardous areas of any classification. Always stop and remove your insulin pump prior to entering these areas.
  
- ➔ **Your insulin pump has not been tested with other electronic medical devices.** Therefore, do not use the Accu-Chek Spirit insulin pump with another electronic medical device unless advised by your doctor or healthcare team.
  
- ➔ **Temperatures over +40°C (+104°F) and below +5°C (+41°F) may damage the insulin,** the electronics of your insulin pump and may cause a malfunction of the battery. Do not place your insulin pump in direct sunlight. Overheating of the insulin and your insulin pump must be avoided.
  
- ➔ **Protect your insulin pump from direct exposure to cold wind.** In cold weather wear your insulin pump under your clothes or directly on your body. Refer to the instructions for use of the insulin you are using for information on the admissible temperature range for the insulin.

### 2.1.6 Warnings and cautions regarding safe battery use

- ➔ **Always have extra AA batteries available with you and pay attention to the following:**
  - Always keep batteries in their original package until ready to use.
  - Prevent early discharge by avoiding contact between batteries themselves and between batteries and any other metal objects (e.g. coins, keys).
  - To operate your insulin pump, the operating temperature of the battery must be +5°C to +40°C (+41°F to +104°F).
  - Do not use a battery that has been dropped.

- ➔ **Change the battery in a dry environment.** Tighten or loosen the battery cover with the Accu-Chek Spirit battery key (sharp objects could damage the insulin pump). Do not over-tighten. Over-tightening may cause damage to the battery cover and to the insulin pump casing. The battery cover is correctly inserted and tightened when the battery cover is even with the insulin pump casing.
- ➔ **Do not use the insulin pump if the battery cover is not correctly inserted.** See the section 5.1 “Insert and change the battery” for information on changing the battery and battery cover.
- ➔ **Using AA batteries other than those supplied or recommended by Roche Diagnostics may significantly reduce battery usage time and void the warranty.** They may leak and corrode the battery contacts within your insulin pump.
- ➔ **Always make sure that you remove the battery if your insulin pump is not in use** for a longer time period in order to preserve the battery lifetime. Old or used batteries may cause a failure at the insulin pump start-up.
- ➔ **Change the battery only in a dry environment** to prevent water from entering the casing, and make sure that the seal is not worn out or missing and that the battery is properly inserted.

### 2.1.7 Warnings and cautions about wearing the insulin pump

- ➔ **Avoid any contact with objects that may damage or accidentally press the keys of your insulin pump** (e.g. keychain, garment buttons, pocket knives, coins).
- ➔ **If you cannot hear the beeps your insulin pump needs to be returned since** you may not be notified of any changes of your insulin pump in a suitable amount of time. Contact your local Customer Care Area.



- ➔ **If you cannot feel the vibration alert, your insulin pump needs to be returned** since you may not be notified of any changes of your insulin pump in a suitable amount of time. Contact your local Customer Care Area.
- ➔ **If you discover incomplete characters, numbers or symbols on the display, your insulin pump needs to be returned** since you may not be notified of any changes of your insulin pump in a suitable amount of time. Contact your local Customer Care Area.
- ➔ **Check the display of your insulin pump at least every three hours during the day**, before you go to sleep, and especially if for any reason you might be unable to hear the beeps or to feel the vibrations. This is the only way you will be notified of any changes of your insulin pump in a suitable amount of time.
- ➔ **If the buttons on your pump are not functioning properly or you are unable to identify the buttons**, disconnect your insulin pump and contact your local Customer Care Area.

**If you have any additional questions about these warnings and cautions please feel free to contact your local Customer Care Area.**

## 2.2 Guarantee

- ➔ **Any changes or modification to the devices not expressly approved by Roche Diagnostics could void your operating warranty for the Accu-Chek Spirit insulin pump.**



## 3 The Accu-Chek Spirit insulin pump

### Caution

Check your insulin pump and its sterile products and accessories for chips and cracks at least once per day, especially if they have been dropped. Do not use them if chipped or cracked. In the presence of chips and cracks, water, dust, insulin, or other foreign substances may enter your insulin pump and lead to malfunction.

### 3.1 Display



Your insulin pump has a full graphic LCD (Liquid Crystal Display) which displays important current and historical information. **Check the display of your insulin pump at least every three hours during the day, before you go to sleep, and especially if for any reason you might be unable to hear the beeps or to feel the vibrations.** This is the only way you will be notified of any changes of your insulin pump in a suitable amount of time.

For your convenience you can flip the display orientation by 180° and adjust the display contrast if desired. See the sections 7.7.6 “Display orientation” and 8.5.7 “Display contrast” for further information.

### Caution

If you discover incomplete characters, numbers or symbols on the display, put your insulin pump into **STOP** and perform a system check by simply removing the battery for a few seconds. Upon reinsertion of the battery the system check will occur. If the incomplete characters, numbers or symbols remain, your insulin pump needs to be returned since you may not receive the correct information for operating your pump. Put your insulin pump into **STOP** and contact your local Customer Care Area immediately.

## 3.2 Backlight

The backlight helps you use your insulin pump or read the information in poorly lighted areas. You can turn on the backlight from the **RUN** screen, from the **STOP** screen or when browsing the menus, by pressing . The backlight turns on automatically when an alert or error message appears. The backlight turns off automatically after 20 seconds  if no further keys are pressed.

**Note** In areas of poor lighting, turn the backlight on before you start programming.

## 3.3 Keys and key combinations








Your insulin pump's four keys are all you need to program the beeps, vibration alerts and messages which appear on the display.

### Caution

Never press any of the keys with a sharp or pointed item such as the tip of your fingernail. This could damage the casing of your insulin pump. Use the pad of your finger tip instead to prevent puncturing the keys.



The use of your insulin pump requires the use of four keys and three key combinations. Each key or key combination has a specific function:

Key	Name	Function
	<b>Menu</b>	<b>move</b> through menus, function and information screens
	<b>Check</b>	<b>select</b> a menu  <b>save</b> changes and <b>exit</b> the function and information screen  <b>view</b> the quick info screen
	<b>Up</b>	<b>increase</b> or <b>decrease</b> a setting
	<b>Down</b>	<b>move forward</b> or <b>backward</b> in the information screens  <b>program</b> a Standard bolus  <b>cancel</b> a Standard bolus  <b>turn off</b> the <b>STOP</b> -Warning (press and hold  or  for 3 seconds until you hear a melody)
	<b>Up</b>	<b>turn on</b> the backlight

**A plus symbol (+) connecting 2 key functions means that both keys must be pressed simultaneously:**

Key Combination	Name	Function
+	<b>Menu and Up</b>	<b>Exit</b> a menu or function and information screens. <b>Loop backward</b> in the menu structure.
+	<b>Menu and Down</b>	Unlock the KeyLock (press and hold simultaneously until you hear three beeps in <b>RUN</b> or one beep when in <b>STOP</b> ).
+	<b>Up and Down</b>	<b>Copy</b> an hourly basal rate to the next hour(s).
		Press either key when two key symbols are unconnected to set a desired value.
()		You can program some sequence starting either with the  or  key. For the recommended sequences press the key without parenthesis, here the  key. For the alternative sequence press the key in the parenthesis, here the () key.

You will hear a beep each time a key is pressed, unless the beep tone has been set to 0.






**Warning**

If the buttons on your pump are not functioning properly or you are unable to identify the buttons, disconnect your insulin pump and contact your local Customer Care Area.



### Symbol Timeout ☺

For your safety and convenience during programming, your insulin pump automatically returns to the **RUN** or **STOP** screen if no key is pressed within 20 seconds ☺. For your safety, the changes you made before your insulin pump “times out” will not be saved. Within this user guide the symbol ☺ represents the timeout feature.







## 3.4 KeyLock

Symbol	Beep	Display	Comment
	A sequence of beeps (you will hear 3 beeps in <b>RUN</b> and 1 beep in <b>STOP</b> ).	 means that all 4 keys are locked (except  for turning on the backlight).	If no KeyLock symbol appears in the <b>RUN</b> or <b>STOP</b> screen, the KeyLock function is turned off.
		 means that all 4 keys are unlocked.	

The KeyLock function allows you to lock all four keys on your insulin pump and serves as an additional safety measure against unintentional activation of functions (e.g. while sleeping, during sports). To make use of the KeyLock function, it must be turned on in the SETUP MENU STANDARD.





If the KeyLock function is turned on, it needs to be unlocked or turned off before any programming of your insulin pump can be performed. All instructions in this user guide are given with the assumption that the KeyLock is NOT turned on. To unlock the KeyLock, press  +  simultaneously (you will hear three beeps in **RUN**, or one beep in **STOP**) and hold for 3 seconds until you hear again three beeps (in **RUN**) or one beep (in **STOP**).

### 3.5 Scrolling

In order to enter larger or smaller values, the  or  key must normally be pressed many times with small increments per press. However, in most menus, you can hold the  or  key (scrolling) until the desired number is displayed. This value can be corrected by a single press of  or  if needed.

If activated, your insulin pump will beep once when you start scrolling.

---

Key	Beep	Display	Comment
	A short initial beep	Pressing the key	By pressing the  or
	identifies the	results in a change	the  key you can
	scrolling function.	of display information.	change a value. Each
			press of the key will
			increase or decrease the
			value by one increment.

---



## 3.6 STATUS screen

Beeps, vibrations and messages on the display inform you of the status of your insulin pump. Your insulin pump beeps when a key is pressed, regardless of how the alarm signals are set. You can turn off either the beeps or the vibrations as alarm signals, but not both at the same time. The beep volume can be adjusted to different levels.

See the sections 7.7.2 “Adjust the beep volume”, 7.7.3 “Alarm signals” and 16.2 “Beeps and melodies” for further information.

**Note** If the beep volume is turned off, your insulin pump will not beep when a key is pressed.

### Caution


Check the display of your insulin pump at least every three hours during the day, before you go to sleep, and especially if for any reason you might be unable to hear the beeps or to feel the vibrations. This is the only way you will be notified of changes in your insulin pump in a suitable amount of time.

## 3.7 Alarm signals

Your insulin pump communicates the alerts and errors using beeps and vibrations. You can turn off either the beeps or the vibrations, but not both at the same time for the initial alarm signal.

### 3.8 RUN and STOP modes


During normal use, your insulin pump is in **RUN** and insulin is delivered (basal rate) at the rate you have programmed. The **RUN** screen is considered the “home” screen and will be displayed whenever the pump is in **RUN**. After any programming or review of information, the insulin pump will return to the **RUN** screen. On this screen, you will see the current time, the current basal rate profile and the current hourly basal rate in units per hour. In addition, possible reminder information (e.g. cartridge low) and, information on special functions (e.g. a running Extended bolus or KeyLock locked or unlocked) may appear on the **RUN** screen.



See the sections 7.3 “Start insulin delivery”, 7.7.3 “Alarm signals” and 16.4 “Symbols” for further information on the definition of these symbols.

Some functions must be programmed with your insulin pump in **STOP**, during which insulin will not be delivered. Your insulin pump must be in **STOP** to change the cartridge. Always disconnect your pump from the infusion site and place the pump in **STOP** when you prime the infusion set and/or transfer data.

The time and date, the **STOP** symbol, possible reminder information (e.g. battery low) and special functions (e.g. KeyLock locked or unlocked) appear on the **STOP** screen. See the section 7.4 “Stop insulin delivery” for further information.



### 3.9 STOP-Warning

The **STOP**-Warning is turned on each time your pump is changed from **RUN** into **STOP** or when a battery is inserted. This reminder is to inform you that insulin delivery is currently interrupted. If the **STOP**-Warning is turned on, it occurs every minute with a long beep and a vibration. Turning off the **STOP**-Warning stops this warning function.

#### To turn off the STOP-Warning

Action	Result	
Press and hold ▲ or ▼ (for 3 seconds).	You hear a melody. The <b>STOP</b> -Warning is turned off.	If you press ▲ or ▼ too briefly, the <b>STOP</b> -Warning is not suppressed and occurs again.
		The volume of the <b>STOP</b> -Warning is independent of the programmed beep volume. It always occurs at the maximum level.

The **STOP**-Warning will be reactivated the next time your pump is changed from **RUN** into **STOP** or with the next battery change.

The operation of your insulin pump is logical and simple. For your convenience a “map” of these functions is found in the back of this user guide.

---

## 4 Accessories, disposables and software

Accu-Chek sterile products, accessories and software complete your Accu-Chek Spirit insulin pump system. They are specially designed for safe and convenient insulin pump therapy.

### **Caution**

Only use sterile products and accessories that are designed for use with your insulin pump. The proper functionality of your insulin pump can only be guaranteed in conjunction with Accu-Chek sterile products and Accu-Chek Spirit accessories. All Accu-Chek sterile products and accessories designed for the Accu-Chek Spirit insulin pump have been thoroughly tested and approved for use with your insulin pump. Other sterile products and accessories have not been tested for compatibility with your insulin pump and may therefore endanger your health if used.

Your Accu-Chek Spirit insulin pump, sterile products, accessories and software should be used for the first time in the presence of your doctor or healthcare team. Regular medical examinations are required. Always follow the instructions given by your doctor or healthcare team and the instructions for use of the sterile products and accessories you are using.

**Caution**

Always have extra sterile products and accessories with you. This allows you to exchange components when needed. Materials designed for single use (such as cartridges and infusion sets) should not be reused due to the increased risk of infections, malfunctions and/or incorrect insulin delivery.

Consult our sterile products brochures and accessories catalogue or your local Customer Care Area for more information.

**Warning**

Your insulin pump and its sterile products and accessories include small component pieces that could pose a choking hazard to children. Keep sterile products and accessories out of the reach of children.

## **4.1 Sterile products**

Accu-Chek sterile products are a critical part of your Accu-Chek Spirit insulin pump system and therapy. Sterile products are designed for single use only to enable proper function of your Accu-Chek Spirit insulin pump and to avoid infections.


### 4.1.1 Cartridge

Your insulin pump was designed to use an Accu-Chek 3.15 ml plastic cartridge. The cartridge is a sterile product intended for single use only. Sterility is guaranteed for unopened packaging up to the indicated expiration date. Do not use sterile products if the package is damaged.

#### Warning


Do not reuse single use materials. Reuse of single use materials might lead to a malfunction of your insulin pump and to incorrect insulin delivery and/or an infection. Always handle the items with clean hands.

Use specified cartridges only. The use of other cartridges can endanger your health and may void the warranty.



- ▶ Look at the sections 5.2.1 “Prepare the cartridge”, 5.3.1 “Connect the cartridge, adapter and infusion set” and 5.7 “Change the cartridge and infusion set” for further information.

## 4.1.2 Infusion sets



The Accu-Chek infusion sets connect your insulin pump to your body and are a crucial part of your insulin pump therapy. The insulin is delivered from the cartridge in your insulin pump through the infusion set into your subcutaneous tissue. The cannula or needle of the infusion set is usually placed in the abdomen.

Infusion sets can be:

- ▶ disconnectable
- ▶ non-disconnectable

Use only infusion sets with a luer-lock connector. All current Accu-Chek infusion sets have luer-lock connectors and are the preferred use for your Accu-Chek Spirit insulin pump. All Accu-Chek infusion sets are PVC free. They are made of materials which are skin-friendly and neutral to insulin. The diameter of the Accu-Chek infusion set tubing has been reduced in order to achieve faster and more economical priming and to minimize waste of insulin.




Roche Diagnostics offers a wide range of infusion sets:

Common lengths of the tubing	Filling volume with U100 insulin
30 cm / 12"	7 units
60 cm / 24"	10 units
80 cm / 31"	13 units
110 cm / 43"	18 units

See the sections 5.3.1 “Connect the cartridge, adapter and infusion set”, 5.4 “Prime the infusion set”, 5.5 “Prepare the infusion site”, 5.6 “Change the infusion set” and 4 “Accessories, disposables and software” for further information.

## 4.2 Adapter



The adapter physically connects the cartridge and the infusion set. It has two seals and forms an efficient seal against water ingress for the cartridge compartment of your insulin pump. The two small vents on the adapter allow air pressure to equalize. The adapter only works correctly when the vents are not plugged or dirty and when the seals are not worn out or missing. The adapter should be replaced with at least every 10<sup>th</sup> cartridge change.

See the sections 5.3.1 “Connect the cartridge, adapter and infusion set” and 5.8 “Change the adapter” for further information.

## 4.3 Battery



Your insulin pump has been designed for use with AA ALKALINE batteries with a minimum capacity of 2500 mAh. Do not use lithium, carbon zinc or nickel cadmium (NiCd) batteries. If you prefer rechargeable batteries, Roche Diagnostics recommends NiMH batteries with a minimum capacity of 1500 mAh. If you use rechargeable batteries, always use a battery charger recommended by the battery manufacturer.

The battery life is approximately four weeks for alkaline batteries and one week for rechargeable batteries if used in a typical usage pattern (50 U/day; operating temperature  $22^{\circ}\text{C} \pm 3^{\circ}\text{C}$  [ $72^{\circ}\text{F} \pm 6^{\circ}\text{F}$ ]).

**Note** If you change from using regular batteries to rechargeable batteries or vice versa, you must change the battery type in the SETUP MENU STANDARD.

See the section 7.7.5 “Battery type” for further information.


### Caution

Using AA batteries other than those supplied or recommended by Roche Diagnostics may significantly reduce battery usage time. Batteries other than those recommended could leak and corrode the battery contacts within your insulin pump. Therefore, the use of batteries not supplied or recommended by Roche Diagnostics may void the warranty.

**Note** Many types of batteries on the market are not designed to provide adequate power for your insulin pump. To ensure that the battery lasts as long as possible, be sure to use alkaline batteries with a minimum capacity of 2500 mAh or NiMH (rechargeable) batteries with a minimum capacity of 1500 mAh.

Alkaline batteries supplied by Roche Diagnostics are the correct batteries for providing a maximum battery life. The operating temperature of the battery must be between +5°C and +40°C (+41°F and +104°F). The lifetime of the battery is affected by usage of your insulin pump, your personal settings, delivery rates, temperature and other factors.

## 4.4 Battery cover



The battery cover forms an efficient seal for the battery compartment of your insulin pump. The battery cover should be replaced with at least every 4<sup>th</sup> battery change.

- Change the battery in a dry environment to prevent water from entering the casing.
- Make sure that the seal is not worn out or missing and that the battery is properly inserted into your insulin pump.
- Tighten or loosen the battery cover only with the Accu-Chek Spirit battery key (using knives, screwdrivers or other sharp objects may harm your insulin pump).
- Do not over-tighten as this may cause damage to the battery cover and to the casing. The battery cover is correctly inserted and tightened when the battery cover is even with the casing.

See the section 5.1 “Insert and change the battery” for further information.

**Caution**

Always make sure that you remove the battery if your insulin pump is not in use for a longer time period in order to preserve the battery lifetime. To prevent water from entering the casing, change the battery only in a dry environment and make sure that the seal of the battery cover is not worn out or missing and that the battery is properly inserted.

## 4.5 Battery key



Use the battery key to tighten and to loosen the battery cover. You may also use it to loosen the luer-lock connection between your infusion set and adapter if it cannot be loosened by hand. The battery key has a notch that fits the luer-lock connectors of the Accu-Chek FlexLink and the Accu-Chek TenderLink infusion sets. Never use the battery key to attach or tighten an infusion set.

## 4.6 Emergency kit

It is recommended that you take a small amount of supplies with you in case of an emergency. This allows you to exchange components when needed. Such an emergency kit may contain:

- A new Accu-Chek infusion set
- A new AA Alkaline battery
- A new 3.15 ml plastic cartridge
- A pen or syringe for alternate therapy
- A vial of insulin
- Treatment for low blood glucose levels
- Blood glucose monitoring supplies (such as an Accu-Chek blood glucose monitor)
- A disinfectant for skin
- A new site dressing
- The battery key

The above list is only an example of what this emergency kit should contain. Ask your doctor or healthcare team how to compose your personal emergency kit.

## 4.7 Software

The Accu-Chek Insulin Pump Configuration Software is an optional insulin pump programming tool for the Accu-Chek Spirit insulin pump\*. This software allows for fast and easy change of parameters and settings directly from your Microsoft Windows compatible PC.

**Accu-Chek Insulin Pump Configuration Software** is designed to allow insulin pump users to conveniently program new settings to their insulin pump using a Microsoft Windows compatible PC, to transfer settings and information to and from their insulin pump.

**Accu-Chek Insulin Pump Configuration Software Pro for healthcare professionals** is designed to allow healthcare professionals to program the insulin pump and to effectively manage patients' insulin pump data. Using Accu-Chek Insulin Pump Configuration Software Pro for healthcare professionals allows healthcare professionals to program additional settings.

Contact your local Customer Care Area for further information or refer to the Accu-Chek Insulin Pump Configuration Software user manual.

---

\* Available in selected countries depending on registration approval.

## 5 Prepare the pump for use

This section will show you how to get your pump up and running.

### 5.1 Insert and change the battery

Your insulin pump turns on as soon as you insert the battery. When you remove the battery, your insulin pump keeps track of the time and date for approximately 1 hour. Your insulin pump settings (such as the hourly basal rates, bolus increment and active user menu) and the event memory (bolus and alarm history, history of daily insulin totals and temporary basal rates) are saved, regardless of battery voltage and the time your insulin pump has been without a battery.

#### **Keep the time/date settings correct**

Always ensure that the time and date of your insulin pump are set correctly. Incorrect programming of the time and date may cause incorrect insulin delivery. If you or your healthcare team store and analyse your therapy data electronically and the time and date of the devices used are not set identically, the gathered data might not be meaningful.

**Caution**

To prevent water from entering the casing, change the battery only in a dry environment and make sure that the seal is not worn out or missing and that the battery is properly inserted. Tighten or loosen the battery cover with the Accu-Chek Spirit battery key (use of knives, screwdrivers or other sharp objects may harm your insulin pump). Do not over-tighten as this may cause damage to the battery cover and to the casing. The battery cover is correctly inserted and tightened when the battery cover is level with the insulin pump casing.


Your insulin pump is powered with a single 1.5 volt AA battery. See chapter 4.3 “Battery” for more detailed information on recommended batteries.



Make sure your insulin pump is in **STOP** and the KeyLock is unlocked or turned off.


Remove or disconnect your infusion set from the infusion site to ensure that no risk of insulin delivery may occur.






Remove the battery cover from your insulin pump. Use the battery key to turn the battery cover counter-clockwise as indicated by the arrow on the battery cover.

Make sure the opening of the battery compartment and the seal are clean and undamaged.



Put the battery cover on the positive end of the battery.

Insert the battery with cover, with the flat (negative) end first, into the battery compartment.



Carefully push the battery with cover in and turn the battery cover clockwise. The battery cover is correctly inserted and tightened when the battery cover is even with the casing. Do not over-tighten.

Your insulin pump will now perform its start-up procedure.

A melody signals the end of the start-up procedure, and your insulin pump returns to the **STOP** screen.

**Note** If your insulin pump was in **RUN** before you removed the battery, an error E8: POWER INTERRUPT occurs.

Press twice  to confirm and to turn off the error.

Check the time and date in the **STOP** screen and correct, if necessary.

Put your insulin pump into **RUN**, if necessary.

**Caution**

When you remove the battery your insulin pump remembers the settings you made and saved. If your Accu-Chek Spirit insulin pump remains without a battery for more than one hour, check the time and date upon reinsertion of a battery.

Using AA batteries other than those supplied or recommended by Roche Diagnostics may significantly reduce the battery usage time. Batteries other than those recommended may leak and corrode the battery contacts within your insulin pump. Therefore, the use of batteries not supplied or recommended by Roche Diagnostics may void the warranty. Do not insert old or used batteries into your insulin pump as this may cause a failure at the insulin pump start-up.

## 5.2 Start-Up procedure

Your insulin pump performs a start-up procedure (internal tests) when you insert a battery and when you change the cartridge. When the battery is changed the entire start-up procedure will be performed. When only changing the cartridge, the start-up procedure will begin at the SELF-TEST step.

### Warning

Do not interrupt the start-up procedure by pressing keys or by any handling on the insulin pump. Interruption of the start-up procedure may lead to malfunctions of your insulin pump.

If your insulin pump does not beep and vibrate as described, or if individual numbers, letters, symbols or lines are shown incompletely or not at all, contact your local Customer Care Area for support.

**Note** When your insulin pump is in use, its safety system monitors its functioning continuously. If it detects a deviation from the normal state, an alert (warning instruction) or an error (error message) occurs.

SW VERSION  
V2.XX


The software version (SW VERSION) of your insulin pump appears on the display.

**Accu-Chek**

The Accu-Chek logo appears.

SELF-TEST

SELF-TEST appears.

BEEP TEST  


Your insulin pump beeps and BEEP TEST appears.

Listen for the beeps.

**Caution**

If you cannot hear the beeps your insulin pump needs to be returned since you may not be notified of any changes of your insulin pump in a suitable amount of time. Contact your local Customer Care Area.

VIBRATION TEST



Your insulin pump vibrates and VIBRATION TEST appears.

Check the vibrations.

**Caution**

If you cannot feel the vibrations your insulin pump needs to be returned since you may not be notified of any changes of your insulin pump in a suitable amount of time. Contact your local Customer Care Area.



The display turns white.  
Check that the display is empty.  
The backlight should turn on during this step.



The display shows a pattern.  
Check that the pattern is uniform.



The display turns black.  
Check that the display is completely black.

### Caution

If you discover incomplete characters, numbers or symbols on the display, your insulin pump needs to be returned since you may not be notified of any changes of your insulin pump in a suitable amount of time. Contact your local Customer Care Area.

Accu-Chek Spirit

Additionally, you will see an information screen. This data can be modified by you or your physician with the Accu-Chek Insulin Pump Configuration Software. Please see the Accu-Chek Insulin Pump Configuration Software user manual for additional information.

**Note** If your insulin pump detects an alert or error during the start-up procedure, it occurs at this point.

See the section 10 “Alerts and errors” for further information.

Your insulin pump goes into **STOP**.



The **STOP**-Warning occurs every minute with a long beep and a vibration to remind you that no insulin is delivered in **STOP**.

Press and hold  or  until you hear a melody to turn the **STOP**-Warning off.

**Note** If the voltage of the inserted battery is too low, the start-up procedure cannot be performed completely. Your insulin pump beeps a five tone and keeps restarting the start-up procedure until you remove the battery or the battery is depleted.

## 5.2.1 Prepare the cartridge


### **Warning**

Your insulin pump was designed to use an Accu-Chek 3.15 ml plastic cartridge. The cartridge is a sterile product intended for single use only. Sterility is guaranteed for unopened packaging up to the indicated expiration date.


- ▶ Do not use sterile products if the package is damaged.
- ▶ Do not reuse single use materials. Reuse of single use materials might lead to a malfunction of your insulin pump and to incorrect insulin delivery and/or an infection.
- ▶ Always handle the items with clean hands.
- ▶ Use specified cartridges only. The use of other cartridges can endanger your health and may void the warranty.

A cartridge set consists of the following items:

1. 3.15 ml plastic cartridge with protective cap, plunger rod and filling needle.



## 2. Instructions for use




### 5.2.2 Fill the cartridge

A filling aid by Roche Diagnostics can be used to help you fill an empty 3.15 ml plastic cartridge from your vial of insulin. For detailed instructions for using a filling aid consult the instructions included with that item.


To fill a cartridge, have the following materials ready:

- ▶ New 3.15 ml plastic cartridge with protective cap, plunger rod and filling needle
- ▶ insulin at room temperature to minimize formation of air bubbles.



Wash your hands.


Set the insulin vial on a hard and flat surface, like a table.



Clean the rubber membrane of the vial with an antiseptic wipe.

Remove the cartridge, the cartridge protective cap and the filling needle with its protective cap from their packaging.






Place the filling needle with protective cap onto the cartridge tip and make sure the connection is tight.

Move the plunger rod twice back and forth in the cartridge to distribute the lubricant.  
Pull the plunger back to fill the cartridge with air.

Remove the protective cap from the needle.  
Ensure that you do not touch anything with the needle.




Push the tip of the filling needle into the centre of the insulin vial's rubber membrane.

Push the plunger rod down so that all air from the cartridge enters the insulin vial.

While holding steady pressure against the plunger rod with your thumb, invert the insulin vial so that the filling needle and cartridge are pointing upwards into the insulin vial.

Make sure that the tip of the filling needle stays in the insulin.

Slowly release the pressure against the plunger rod and allow the insulin to flow into the cartridge.




Do not pull or push the plunger rod while it is moving on its own (this would cause the formation of air bubbles).

Slowly pull the plunger rod straight down to fill the cartridge.

Remove any air bubbles by tapping on the cartridge to dislodge them. Push the air bubbles back into the insulin vial with the plunger rod.

The cartridge is completely filled when there are no air bubbles in the cartridge and the plunger is at the bottom of the cartridge.

Remove the filling needle from the insulin vial.



Place the protective cap onto the filling needle.

Remove the plunger rod by turning it clockwise out of the plunger.

Remove the filling needle with its protective cap from the cartridge by turning counter-clockwise.

Push the cartridge protective cap tightly onto the cartridge tip until it clicks.

The cartridge is now ready for use.

**Note** Do not pull or push the plunger rod during removal.

### Caution

Air bubbles in the cartridge and infusion set may cause the infusion of air instead of insulin. Your body does not receive the required amount of insulin. An error E4: OCCLUSION may be delayed. Remove these bubbles while filling the cartridge and priming the infusion set and without having the infusion set connected to your body.

Cold insulin may release air when warmed. Use only insulin at room temperature when filling the cartridge and priming the infusion set. Inspect the cartridge and the infusion set for air bubbles at least every three hours during the day and before you go to sleep. Remove any air bubbles and, if necessary, change system components.

## 5.3 Insert the cartridge, adapter and infusion set

**Note** Prior to changing the cartridge, adapter and infusion set, care should be taken to ensure that your insulin pump is in good, working order. See the section 14.1 “System check” for the inspection checklist.

Free flow of insulin from the cartridge or infusion set can occur when the cartridge plunger and the piston rod are not properly connected and you position your insulin pump at a higher level than the infusion site. Prevent free flow by correctly inserting the cartridge: move the piston rod to the proper position and twist the adapter until the cartridge plunger sits flush with the end plate of the piston rod.

When inserting an approved cartridge and adapter into the cartridge compartment twist the adapter clockwise until the adapter sits flush with the cartridge compartment. In this position the cartridge plunger must also sit flush against the end plate of the piston rod. This flush fit provides an additional safety measure to prevent the free flow of insulin and minimizes the occlusion volume.

### 5.3.1 Connect the cartridge, adapter and infusion set

Have the following materials ready:

- ▶ A filled Accu-Chek 3.15 ml plastic cartridge with a luer connection
- ▶ An Accu-Chek Spirit adapter
- ▶ A new Accu-Chek infusion set


### Warning

The cartridge and infusion set are **sterile products intended for single use only**. Sterility is guaranteed for unopened packaging up to the indicated expiration date. **Do not use sterile products if the packaging is damaged.**

Do not reuse single use materials. Reuse of single use materials might lead to a malfunction of your insulin pump and to incorrect insulin delivery and/or an infection. Always handle the items with clean hands. Avoid any contact of the infusion set and especially the connecting parts of your insulin pump with antiseptics, antibiotic creams, soaps, perfumes, deodorants, body lotions or any other cosmetics. They may contaminate these parts.

Remove the cartridge protective cap.

Make sure that you do not touch the cartridge tip.



Push the adapter onto the cartridge tip as far as to the stop (a).

Carefully prepare a new Accu-Chek infusion set for use.

Hold the adapter and twist the infusion set by hand clockwise into the adapter (b).

Do not overtighten.

**Caution**

To avoid leakage, properly tighten the infusion set to the adapter. Turn the luer-lock connector only as far as to the stop. Do not turn it any further and do not use any auxiliary tools as this may crack the infusion set luer-lock connector and result in leakage.


**Note** Your insulin pump cannot detect infusion set leakage. You must inspect all parts of the infusion set at least every three hours during the day and before you go to sleep. Should you detect any loss of insulin and all parts are correctly tightened, immediately replace the leaking component. Because insulin delivery has been interrupted, check your blood glucose level and take appropriate actions according to your doctor's or healthcare team's instructions.

### 5.3.2 Insert the cartridge

**Caution**

Do not carry out the CHANGE THE CARTRIDGE function if the cartridge compartment is not completely dry. During the rewinding of the piston rod liquid might enter your insulin pump and could lead to malfunction.

Should the piston rod not return completely, contact your local Customer Care Area for support.



Have the following materials ready:


- ▶ Your Accu-Chek Spirit insulin pump
- ▶ A filled 3.15 ml plastic cartridge with a connected Accu-Chek Spirit adapter and a new Accu-Chek infusion set.



Make sure your insulin pump is in **STOP**.



Press to move to the CHANGE THE CARTRIDGE menu.




Select with .

Remove the current cartridge, if necessary.




Press and hold (for 3 seconds) until you hear a melody.


Your insulin pump starts to return the piston rod.






When the piston rod has returned completely, the CARTRIDGE VOLUME screen displays the maximum cartridge volume (315 U).



Hold your insulin pump upright. Hold the new filled cartridge with connected adapter and infusion set pointing upwards and hold them parallel and close to the cartridge compartment. Make sure that the bottom edge of the coloured part of the adapter is level with the top of the cartridge compartment.



Press or scroll  or  to move the plunger rod forward until the end plate of the piston rod is level with the bottom of the cartridge plunger. The cartridge volume decreases accordingly on the display.





Position your insulin pump upright. Insert the cartridge with attached adapter and infusion set into the cartridge compartment. Twist (do not push) the adapter clockwise until it sits flush with the cartridge compartment. It is not necessary to apply pressure as the twisting will properly position the cartridge. The cartridge has been correctly inserted when the end plate of the piston rod sits flush against the cartridge plunger. Do not overtighten the adapter as this may make removal of the cartridge difficult.

If the end plate of the piston rod does not sit flush against the cartridge plunger, remove the cartridge from the cartridge compartment.


Hold the cartridge and adapter next to the cartridge compartment.



Press or scroll  or  to move the piston rod forward until the base plate of the piston rod and the cartridge plunger are level. The displayed cartridge volume decreases accordingly. Re-insert the cartridge into the cartridge compartment.

- Continue with these instructions when the cartridge is correctly inserted (the end plate of the piston rod sits flush against the cartridge plunger), or
- Repeat this step until the cartridge is correctly inserted (the end plate of the piston rod sits flush against the cartridge plunger).

SELF-TEST

When the cartridge has been inserted correctly press . Your insulin pump performs a self test.

See the section 5.2 “Start-Up procedure” for further information.

Check that the infusion set luer-lock connector is still correctly connected to the adapter. Tighten the infusion set by hand clockwise until it sits tightly in the adapter.

PRIME THE  
INFUSION SET




When the start-up procedure is finished, the PRIME THE INFUSION SET menu is displayed. See the section 5.4 “Prime the infusion set” for further information. Start from the point in those instructions where the insulin pump screen matches the menu to the left.

## 5.4 Prime the infusion set

### Warning


Never prime an infusion set that is connected to your body. You risk uncontrolled insulin delivery into your body. With disconnectable infusion sets, make sure you have disconnected the tubing from your site prior to priming. Always follow the instructions for use of the infusion set you are using.

Your insulin pump automatically stops priming after 25 units  of insulin. If insulin has not emerged from the infusion set needle after priming, repeat the prime function. When insulin emerges from the tip of the needle, press any of your insulin pump's keys to stop priming. The amount of insulin used for the priming is not added to the history of daily insulin totals.

During the priming process your insulin pump should be positioned upright, with the adapter pointing upwards so as to maximize the removal of air bubbles from the cartridge and the infusion set.

### Warning

Air bubbles in the cartridge and infusion set may cause the infusion of air instead of insulin. Your body will not receive the required amount of insulin. An error E4: OCCLUSION may be delayed. Inspect the cartridge and the infusion set for air bubbles at least every three hours during the day and before you go to sleep. Remove these bubbles while filling the cartridge and priming the infusion set and without having the infusion set connected to your body. Ensure that the connections between the infusion set and the adapter as well as between the adapter and the pump are properly tightened.




Make sure your insulin pump is in **STOP**.



Press to move to the PRIME THE INFUSION SET menu.

Select with .



START PRIMING and the priming amount of 25 units of insulin is shown on the display.


Press and hold until you hear a melody and the piston rod moves forward. The priming process starts and a count-up is shown on the display. Your insulin pump automatically stops priming after 25 units of insulin.



When the priming process is finished, your insulin pump returns to **STOP**.

The infusion set is properly primed when no air bubbles are visible in the tubing anymore and a bubble-free flow of insulin emerges from the tip of the needle.

## Stop priming



The priming process can be stopped at any time by pressing any key. Your insulin pump returns to the **STOP** screen.

## 5.5 Prepare the infusion site

Always follow the instructions given by your doctor or healthcare team and the instructions for use of the infusion set you are using.

### Infusion site selection

Your doctor or healthcare team will help you select your infusion site areas according to a planned infusion site rotation pattern. You should avoid your waistline, bones, recent infusion sites, bruises and sores. Make sure to choose an infusion site at least 2.5 cm (one inch) away from the navel and previous infusion sites.

### Infusion site preparation


Proper infusion site preparation is essential to reduce the risks of infection. Your doctor or healthcare team will give you infusion site preparation guidelines.

Have the following materials ready:

- ▶ Your Accu-Chek Spirit insulin pump with filled cartridge, connected Accu-Chek Spirit adapter and a new primed infusion set
- ▶ A disinfectant for skin




Wash your hands thoroughly.



Desinfect the infusion site and let it dry completely.



Make sure your insulin pump is in **STOP**.



Insert your Accu-Chek infusion set according to the corresponding instructions for use.

**Note** If you use a soft cannula-type infusion set, you must bolus to fill the air space in the cannula after removing the introducer needle and connecting the tubing according to the instructions for the infusion set you are using. Failure to do so could result in missed insulin dosing.

See the section 7.5 “Program a bolus” for further information. Inspect your infusion site at least once or twice each day for irritation or infection. Signs of infection may include, but are not limited to, pain, lumps, redness, heat or seepage.

If you see redness or swelling, change your infusion set and infusion site immediately and contact your doctor or healthcare team.

You should change your infusion site according to the instructions for use of the infusion set you are using and your doctor’s or healthcare team’s recommendations.

Make sure your insulin pump is correctly equipped (including cartridge, adapter and infusion set) and programmed with your personal settings approved by your doctor or healthcare team.

After you have carried out all steps required for the selection and preparation of the infusion site, your insulin pump is ready for operation. Program your insulin pump with your personal settings prior to starting insulin pump therapy. See the section 7 “STANDARD user menu” for further information.

**Warning**

Incorrect programming of your insulin pump may cause inappropriate insulin delivery. Your insulin pump must be programmed with your personal settings prior to starting insulin pump therapy. Do not operate your insulin pump without knowing your personal settings.

## 5.6 Change the infusion set

When changing the infusion set, please keep the following points clearly in mind:


- ▶ Infusion sets are sterile products intended for single use only. Sterility is guaranteed for unopened packaging up to the indicated expiration date. Do not use sterile products if the package is damaged.  
Do not reuse single use materials. Reuse of single use materials might lead to a malfunction of your insulin pump and to incorrect insulin delivery and/or an infection. Always handle the items with clean hands.
- ▶ Properly tighten the infusion set to the adapter to avoid leakage. Turn the infusion set luer-lock connector only as far as to the stop. Do not turn it any further and do not use any auxiliary tools as this may crack the infusion set luer-lock connector and result in leakage.
- ▶ Your insulin pump cannot detect infusion set leakage. You must inspect all parts of your infusion set at least every three hours during the day and before you go to sleep. Should you detect any loss of insulin and all parts are correctly tightened, immediately replace the leaking component. Immediately check your blood glucose level, because insulin delivery has been interrupted. Check your blood glucose level and take appropriate actions according to your doctor's or healthcare team's instructions.
- ▶ Never change the cartridge or prime with an infusion set inserted into your body. You risk uncontrolled insulin delivery into your body. With disconnectable infusion sets, make sure you have disconnected the tubing from your site prior to changing or priming.

Have the following materials ready:

- ▶ Your Accu-Chek Spirit insulin pump
- ▶ A new Accu-Chek infusion set
- ▶ A disinfectant for skin



Make sure your insulin pump is in **STOP**.



Remove your infusion set from your infusion site.

Remove the infusion set from the adapter and dispose of it properly.

**Note** The battery key has a notch that fits the luer-lock connectors of the Accu-Chek FlexLink and the Accu-Chek TenderLink infusion sets. You may use the battery key to loosen the luer-lock connector between your infusion set and adapter if it cannot be loosened by hand.

Carefully prepare a new Accu-Chek infusion set for use.

Turn the Accu-Chek infusion set by hand clockwise into the adapter. Properly tighten the luer-lock connector by hand.

Prime the infusion set.

When finished, put your insulin pump into **RUN** (see section 7.3 “Start insulin delivery”).




## 5.7 Change the cartridge and infusion set

Have the following materials ready:

- ▶ Your Accu-Chek Spirit insulin pump
- ▶ A filled Accu-Chek 3.15 ml plastic cartridge
- ▶ A new Accu-Chek infusion set
- ▶ A new adapter (must be changed at least with every 10<sup>th</sup> cartridge)
- ▶ A disinfectant for skin




Make sure your insulin pump is in **STOP**.



Remove your infusion set from your infusion site.

Remove the infusion set from the adapter and dispose of it properly.

**Note** The battery key has a notch that fits the luer-lock connectors of the Accu-Chek FlexLink and the Accu-Chek TenderLink infusion sets. You may use the battery key to loosen the luer-lock connector between your infusion set and adapter if it cannot be loosened by hand.



- a. Hold your Accu-Chek Spirit with the adapter pointing downwards, to prevent the remaining insulin from flowing into the cartridge compartment.
- b. Remove the adapter and cartridge from your insulin pump by unscrewing the adapter from the insulin pump casing.

**Note** Ensure that the cartridge plunger has been completely unscrewed from the piston rod before you pull the cartridge out of the cartridge compartment. When unscrewing the adapter, the cartridge plunger also needs to turn.

- c. Remove the cartridge from the adapter and dispose of it properly.
- d. Hold the adapter up to the light to check for signs of wear or dirt (especially on the seals inside and outside of the adapter).
- e. Clean the adapter with water and dry it, if necessary. If it shows signs of wear or dirt, replace it immediately.
- f. Set up your insulin pump with a new cartridge and a new infusion set.

See the sections 5.2.1 “Prepare the cartridge”, 5.3 “Insert the cartridge, adapter and infusion set”, 5.4 “Prime the infusion set” and 7.3 “Start insulin delivery” for further information.

**Note** Check the amount of insulin remaining in the cartridge at least once a day. Before going to sleep, make sure that the cartridge contains enough insulin to last through the night.

## 5.8 Change the adapter

It is recommended to replace the adapter at least after every 10<sup>th</sup> insulin cartridge. Replace an used adapter with a new one while changing a cartridge. See the section 5.3.1 “Connect the cartridge, adapter and infusion set” for further information.

## 5.9 How to wear your pump

### **Caution**

When wearing your insulin pump make sure to prevent any contact with objects that may damage or accidentally press the keys of your insulin pump (e.g. keychain, garment buttons, pocket knives, coins).

Attach your insulin pump safely to your body or clothes to prevent damage. Specially designed carrying systems allow you to carry your insulin pump on or underneath your clothes. All specified accessories have been tested and approved for use with your Accu-Chek Spirit insulin pump. The proper functionality of your insulin pump can only be guaranteed in conjunction with Accu-Chek Spirit accessories.

Consult our sterile products brochures and accessories catalogue for more information or contact your local Customer Care Area.

---

## 6 Select a user menu (STANDARD, ADVANCED or CUSTOM)


### 6.1 The three levels of the Accu-Chek Spirit

The operation of your insulin pump is divided into three levels.

<b>1<sup>st</sup> level</b>	<b>2<sup>nd</sup> level</b>	<b>3<sup>rd</sup> level</b>
<b>RUN screen and STOP screen</b>	<b>Menus</b>	<b>Function screens and information screens</b>
This level is the home screen for your insulin pump's operation. You can get to all other functions from this level.	This level consists of all available menus. Move through the menus to get to a specific function or information.	This level consists of all functions.

When no keys are pressed, your pump automatically returns to the 1<sup>st</sup> level.

## 6.2 Navigate through the menus





1. Press ☰ to move to the menus (2<sup>nd</sup> level).
2. When you reach the desired menu, press ✓ to enter (3<sup>rd</sup> level).
3. Press ↑ or ↓ to increase or decrease the selected value.
4. Press ☰ again to move to the next function screen you want to change.
5. When finished, press ✓ to confirm the changes you have made.

Your insulin pump saves the changes and returns to the **RUN** or **STOP** screen.

**Note** Flow Chart information printed in black is described in the instructions.  
 Gray printed information illustrates additional options.  
 – Three dots (...) indicates alternative ways to exit a menu.

Your insulin pump is equipped with advanced features to make programming even easier.



## 6.2.1 Scrolling



You can hold  or  pressed to enter larger or smaller values until the desired number is displayed. This value can be corrected by single key presses up or down if necessary.

## 6.2.2 Looping

The menus and screens are set up in a “loop” so that you automatically return to the first function or information screen within a menu once you have reached the last choice within a menu.

## 6.2.3 Backward loop

Pressing  +  simultaneously allows you to move backwards in the menu structure or to return to menus you just skipped.

**Note** Within a function and information screen (3<sup>rd</sup> level) pressing  +  simultaneously allows you to move up to the corresponding menu (2<sup>nd</sup> level) without saving the current changes.

## 6.2.4 Exit options



At any function screen you have three options to exit:

if you want to **confirm** and **save the changes**

– press . Your insulin pump returns to the **RUN** screen.

if you want to **undo the changes**

– wait for your insulin pump to return to the **RUN** screen (timeout), or

– press  +  simultaneously (exit feature) to exit the current function screen.

## 6.2.5 Select a user menu

A user menu comprises a fixed or an individually composed selection of the Accu-Chek Spirit menus. Your insulin pump offers a choice of 3 different user menus:

**Note** In the flap at the end of this user guide you will find diagrams for the STANDARD and ADVANCED user menus.

- **STANDARD user menu**

offers all the functions you need for a successful insulin pump therapy and the ability to change between the user menus. We recommend that patients new to pump therapy first use this feature. Later, as your comfort with pump therapy grows, you may wish to choose the additional features offered in the ADVANCED user menu.

- **ADVANCED user menu**

offers the full range of the Accu-Chek Spirit functions. The ADVANCED user menu includes all functions from the STANDARD user menu, plus a wide range of additional functions for the more experienced user.

- **CUSTOM user menu**

because the Accu-Chek Spirit has so many features, you may want to selectively choose those features that you see when scrolling through the menus. The CUSTOM user menu can be adjusted to your individual needs by you (using Accu-Chek Insulin Pump Configuration Software) or by your doctor (using Accu-Chek Insulin Pump Configuration Software Pro for healthcare professionals). Adjusting this feature allows you to display or hide menus. Menus that are fundamental for insulin pump therapy will always be seen and the function settings within a hidden menu remain as set (e.g. on or off) when menus are hidden.



For your personalized CUSTOM user menu you can turn the following menus on or off by using Accu-Chek Insulin Pump Configuration Software:


- EXTENDED BOLUS
- MULTIWAVE BOLUS
- TEMPORARY BASAL RATE (TBR)
- CHANGE BASAL RATE PROFILE
- PROGRAM BASAL RATE PROFILE 2, 3, 4 and/or 5
- ALARM CLOCK
- SETUP MENU STANDARD
- SETUP MENU ADVANCED

Using Accu-Chek Insulin Pump Configuration Software Pro for healthcare professionals, your healthcare professionals can also turn on or off the following menus and functions:

- PROGRAM BASAL RATE PROFILE 1
- SELECT USER MENU

See the Accu-Chek Insulin Pump Configuration Software Pro for healthcare professionals user manual for further information on customizing your user menus.

**Note** When you change your user menu, your current basal rate profile will not appear if it is not activated. Make sure that your desired basal rate profile number(s) is/are activated with Accu-Chek Insulin Pump Configuration Software Pro for healthcare professionals, or select the previous or the ADVANCED user menu.



Press to move to the SELECT USER MENU.


Select with .

The current user menu is displayed.

**Note** If the SELECT USER MENU is not available, contact your doctor or healthcare team.




Press or to select the desired user menu.



Confirm with .


The newly selected basal rate profile and its daily insulin total appear.



or

if the previously turned on basal rate profile is not available from the newly selected user menu, dashes appear instead.

Press or to select an available basal rate profile, if necessary.

BASAL RATE TOTAL  
 2 24.0 U

Check the (new) basal rate profile and its daily insulin total.

Save and exit with .

The new user menu and basal rate profile are active immediately.

**Caution**

A running temporary basal rate increase or decrease will remain active, even if you changed the basal rate profile.

---

## 7 STANDARD user menu

This chapter will familiarize you with the basic functions of your Accu-Chek Spirit insulin pump that are necessary for a successful therapy.

### 7.1 Set time and date

You must take care to set the time and date correctly, because the basal rates and all history information are stored on the basis of time and date.

**When travelling across time zones make sure that the time and date are set correctly.**

The format can be either European or American:

American	European
12-hour clock month/day/year	24-hour clock day.month.year


See the sections 8.5.1 “Time format” and 8.5.2 “Date format” for further information.

**Warning**

Incorrect setting of time and date can lead to incorrect insulin delivery. Make sure that the time and date of your insulin pump is programmed correctly to ensure the correct insulin delivery and data memory. Do not make therapy decisions based upon a single result in the insulin pump's memory.


If you, your doctor or healthcare team review your therapy data electronically, it is essential that time and date of your insulin pump, your blood glucose measurement system (such as an Accu-Chek Blood Glucose Meter), PC and other used devices, are set identically. If they are not, the gathered data might not be meaningful. Check regularly that the time and date of your insulin pump, your blood glucose measurement system, PC and other used devices are identical.

## Set time and date




Press to move to the SET TIME AND DATE menu.


Select with .



Press or to set the **hour**.



Press to move to the minutes.  
Press or to set the **minutes**.



Press to move to the year.  
Press or to set the **year**.



Press to move to the month.  
Press or to set the **month**.



Press to move to the day.  
Press or to set the **day**.  
Save and exit with .

**Note** At any screen you have three options to exit:

if you want to **confirm** and **save the changes**  
– press . Your insulin pump returns to the **RUN** screen.

if you want to **undo the changes**  
– wait for your insulin pump to return to the **RUN** screen (timeout), or  
– press + simultaneously (exit feature) to exit the current function screen.

## Review the time and date

The **time** is displayed on the **RUN** screen and **STOP** screen.  
The **date** (and time) can be reviewed on the **STOP** screen.



## 7.2 Basal rate

Your insulin pump delivers insulin every three minutes, 20 equal doses each hour, 24-hours a day. This flow of insulin, measured in international units per hour (U/h), is called the basal rate and is calculated to meet your basic insulin needs.

An Accu-Chek Spirit basal rate profile consists of up to 24 different hourly basal rates. Each hourly rate may be changed independently. The sum of all 24 hourly basal rates in one basal rate profile is called the (daily) basal rate total.

Your basal rate profile should be reviewed:

- after any programming change performed with your insulin pump or the Accu-Chek Insulin Pump Configuration Software,
- after a new battery has been inserted and
- after confirming an error E7: ELECTRONIC ERROR.



Insulin delivery is not interrupted during basal rate programming.

**Note** The insulin amounts and other values used in this user guide are examples only. Your personal values may differ.


### **Caution**

Programming your insulin pump for a relatively low basal rate (0.1 U/h) can result in a delayed error E4: OCCLUSION. Roche Diagnostics recommends to use plastic cartridges when a low basal rate is required for your therapy.


## 7.2.1 Program a basal rate profile




Press to move to the PROGRAM BASAL RATE PROFILE 1 menu.




Select with .  
The daily BASAL RATE TOTAL is displayed.




Press to move to the first hour.  
The first hour, which always begins at midnight (00:00–01:00 or 12:00 AM–01:00 AM if the American time format is used) is shown.



Press or or scroll to set the hourly basal rate given by your doctor or healthcare team.



Press to move to the next hour. Continue using and the and keys to set the hourly basal rates for the remaining hours. In this way you can program your individual basal rate profile hour-by-hour. Continue until all 24 hours are programmed.



Confirm with .

The new daily BASAL RATE TOTAL is displayed. Check the new daily BASAL RATE TOTAL.

Save and exit with .

If the currently active basal rate profile is the same as the newly programmed basal rate, it is active immediately.

**Note** At any screen you have three options to exit:

if you want to **confirm** and **save the changes**

- press . Check the new daily basal rate total. Press again. Your insulin pump returns to the **RUN** screen.

if you want to **undo the changes**

- wait for your insulin pump to return to the **RUN** screen (timeout), or
- press + simultaneously (exit feature) to exit the current function screen.

## 7.2.2 Copy an hourly basal rate


An hourly basal rate can be copied for 1 or more hours.




Press to move to the PROGRAM BASAL RATE PROFILE 1 menu.



Select with . The daily BASAL RATE TOTAL is displayed.




Press to move to the hourly basal rate which you want to copy.



Press + simultaneously. A special beep confirms that you are pressing the keys correctly.

The same hourly basal rate is copied to the next hour or hours, if you press + repeatedly.

Press and hold + (scroll) to copy the same hourly basal rate for multiple hours.



Confirm with . Check the new daily BASAL RATE TOTAL.

Save and exit with . If the currently active basal rate profile is the same as the newly programmed basal rate, it is active immediately.


### 7.2.3 Basal rate and bolus

- During normal use, your insulin pump is in **RUN** and insulin is continuously delivered (basal rate).
- A bolus can be programmed when needed. Unless you are currently changing settings, giving a bolus or reviewing information, the **RUN** screen is displayed.

**Note** Since only short-acting insulin or fast-acting insulin analogues are used for insulin pump therapy, there is only a small insulin reserve in the body. If insulin administration is interrupted for any reason (e.g. stopping the insulin pump by the user, technical problem with the insulin pump, leakage in the cartridge, occlusion of the infusion set tube or of the infusion set needle, infusion set needle has slipped out of the infusion site), you must be prepared to replace the missing insulin immediately.

Always carry spare sterile products and accessories (infusion set, insulin cartridge, batteries) as well as an insulin pen/syringe and insulin with you. Without insulin, diabetic ketoacidosis may develop and may require in-patient hospital treatment.

## 7.3 Start insulin delivery



Insulin delivery starts at the moment you put your insulin pump into **RUN**.



Start from the **STOP** screen.



Press to move to the START YOUR ACCU-CHEK menu.



Confirm with . Your insulin pump shows the **RUN** screen.

Insulin delivery begins within the next 3 minutes at the hourly basal rate shown on the display.

**Note** Check the amount of insulin remaining in the cartridge at least once a day.

Before going to sleep, make sure that:

- the cartridge contains enough insulin to last through the night.
- the time and date are set correctly.

## 7.4 Stop insulin delivery


For your safety there are functions and actions that require that your insulin pump is in **STOP**. These functions and actions include:

- changing the cartridge
- connecting and disconnecting the adapter and/or infusion set
- priming the infusion set and
- data transfer from the insulin pump to the PC and vice versa using Accu-Chek Insulin Pump Configuration Software.


### Caution


When your insulin pump is in **STOP**, it does not deliver any insulin. Put your insulin pump into **RUN** to continue the insulin delivery.



Press  to move to the STOP YOUR ACCU-CHEK menu.



Confirm with . Your insulin pump returns to the **STOP** screen and insulin delivery stops.

Press  from the **RUN** screen or **STOP** screen (quick info screen) to check for symbols of activated functions and features (like the beep and vibration alarm signals).

**Note** When an Extended bolus, a MultiWave bolus and/or a temporary basal rate (TBR) is/are active and you put your insulin pump into **STOP**, an alert A6: TBR CANCELLED and/or an alert A8: BOLUS CANCELLED occur(s).



Press  twice for each appearing alert to confirm and turn off.

See the sections 10.1.6 “Alert A6: TBR CANCELLED” and 10.1.8 “Alert A8: BOLUS CANCELLED” for further information.

### **STOP-Warning**

The **STOP**-Warning is turned on each time your insulin pump is put from **RUN** into **STOP** or when a battery is inserted. To inform you that the insulin delivery is interrupted, your insulin pump reminds you every minute with a long beep and a vibration. Turning off the **STOP**-Warning suppresses this warning function.

### **Turn off the STOP-Warning**

Press and hold  or  until you hear a melody. The **STOP**-Warning is turned off.

**Note** The volume of the **STOP**-Warning is independent of the programmed beep volume. It always occurs at the maximum level.

The **STOP**-Warning will be reactivated the next time your pump is changed from **RUN** to **STOP** or with the next battery change.



## 7.5 Program a bolus

The bolus type and amount is determined by your doctor or healthcare team's guidelines, your blood glucose level, your food intake, your health condition and your activity level. Discuss the timing, amount and type of bolus you need to deliver with your doctor or healthcare team. The time, date and amount of the last 30 boluses can be reviewed in the bolus history. See the section 7.8.1 "Review the bolus history" for further information.

Your insulin pump allows for the delivery of three types of boluses:

- Standard bolus** (fast delivery)
- Extended bolus** (delivery over a programmed period of time)
- MultiWave bolus** (combines fast delivery with delivery over a programmed period of time)

For instructions on how to program an Extended bolus and a MultiWave bolus please refer to section 8 "ADVANCED user menu".

### Warning



- ▶ Incorrect programming of your insulin pump may cause inappropriate insulin delivery.
- ▶ Your insulin pump must be programmed with your personal settings prior to starting insulin pump therapy.
- ▶ Do not operate your insulin pump without knowing your personal settings.
- ▶ If you are not sure about your personal settings or you are less experienced, have your doctor or healthcare team check your personal settings.
- ▶ Work with your doctor or healthcare team to determine the timing, amount and type of bolus you need to deliver.
- ▶ Make sure you know your personal carbohydrate insulin ratio and your correction bolus ratio.

## 7.5.1 Standard bolus

Your insulin pump offers you two possibilities to program a Standard bolus:





---

**“Quick” Standard bolus**


Using the  and  keys of your insulin pump, with bolus increments determined by the user.




---

**“Scroll” Standard bolus**

Menu-guided using the  and  keys with the  and  keys for programming of the amount.

---

The bolus amount per delivery is limited to 25 units  of insulin.


Programming a Standard bolus includes a short start delay of 5 seconds before the actual bolus delivery. If necessary, this allows you to cancel the bolus before the bolus delivery has started by the press of the  or  key. Alert A8: BOLUS CANCELLED will occur. Press  twice to confirm and turn off the alert. See the section 10.1.8 “Alert A8: BOLUS CANCELLED” for further information.

## 7.5.2 Program a “Quick” Standard bolus

The following programming sequence involves only the  $\uparrow$  and  $\downarrow$  keys of your insulin pump. The sequence is “reversible”, meaning you can start with either key. The recommended sequence is always described with the first key appearing (without parenthesis), e.g.  $\uparrow$ . For the alternative sequence press the key in the parenthesis, e.g.  $(\downarrow)$ .



**Note** Because your insulin pump beeps and vibrates to safely guide you through the programming, you may program the “Quick” bolus with one finger. Once the bolus is programmed, you can count the beeps and vibrations of your insulin pump to verify the correct programming.

The bolus increment for the “Quick” Standard bolus is initially set to 0.5 units  $\text{⌚}$ . It can be changed in the **SETUP MENU ADVANCED** or by using Accu-Chek Insulin Pump Configuration Software.





STANDARD BOLUS  
 $\perp$  0.0u






From the **RUN** screen press and hold  $\downarrow$  ( $\uparrow$ ) until you hear a beep sequence and feel a vibration. This activates the “Quick” STANDARD BOLUS function.

**Note** We recommend starting with the  key in areas of good lighting and, starting with the  key in areas of poor lighting as this will also turn on the backlight.





Press  () repeatedly until the desired bolus amount is reached.

Each press of the  () key adds one bolus increment to the bolus amount. Your insulin pump simultaneously beeps and vibrates once for each programmed bolus increment.

5 seconds  after the last press of  () , your insulin pump confirms the total bolus amount with one simultaneous beep and vibration for each bolus increment programmed. For 5 seconds  (bolus delivery start delay), the Standard bolus symbol () blinks.

After this bolus delivery start delay, your insulin pump beeps three times and starts to deliver the total bolus amount programmed. The countdown of the remaining bolus amount appears on the display.

### **Warning**

The bolus increment that can be programmed into your insulin pump determines the “Quick” Standard bolus amount using the  and  keys of your insulin pump. To ensure correct insulin delivery be sure the bolus increment is set appropriate for your therapy.

### 7.5.3 Cancel a “Quick” Standard bolus



- During programming (the bolus amount blinks):
  - Press the key you started programming with (⏏ [⏏]) to reset the bolus amount to 0.0 units. You hear a melody.  
If no new bolus is programmed, your insulin pump returns to the **RUN** screen and beeps three times after 5 seconds 🗘.  
No bolus is delivered.
  
- During confirmation (when the beeps and vibrations occur), or during the start delay (⏏ blinks):
  - Press ⏏ or ⏏. You hear a melody. Your insulin pump returns to the **RUN** screen. An alert A8: BOLUS CANCELLED occurs.  
Press ⏏ twice to confirm and turn off the alert. See the section 10.1.8 “Alert A8: BOLUS CANCELLED” for further information.  
No bolus is delivered.
  
- During bolus delivery (countdown of bolus amount):
  - Press and hold the ⏏ or ⏏ key for 3 seconds until you hear a melody.  
An alert A8: BOLUS CANCELLED occurs. Press ⏏ twice to confirm and turn off the alert. See the section 10.1.8 “Alert A8: BOLUS CANCELLED” for further information.

#### **Bolus delivery is interrupted.**


The actual bolus amount delivered prior to the cancellation can be reviewed in the bolus history. Please refer to the section 7.8.1 “Review the bolus history” for further information.


Ensure that the cancellation was intended and program a new Standard bolus, if necessary.


### 7.5.4 Program a menu-guided “Scroll” Standard bolus


The bolus increments for the “Scroll” Standard bolus are fixed to 0.1 unit. You can program this bolus in the STANDARD BOLUS menu by holding the  or  key (scrolling) until the desired bolus amount is displayed.







Press  to move to the STANDARD BOLUS menu.



Select with . The BOLUS AMOUNT screen appears.



Press or scroll  to increase or  to decrease the bolus amount.





Press  to confirm the bolus amount. During 5 seconds  (bolus delivery start delay), the Standard bolus symbol (⊥) blinks.

After this bolus delivery start delay, your insulin pump beeps three times and starts to deliver the total bolus amount programmed. The countdown of the remaining bolus amount appears on the display.



## 7.5.5 Cancel a menu-guided “Scroll” Standard bolus

– During programming (the bolus amount blinks):

You have three options for cancelling a “Scroll” Standard bolus during programming:

- If you do not press any key for 20 seconds , your insulin pump returns to the **RUN** screen (timeout).
- Exit to the STANDARD BOLUS menu by pressing  +  simultaneously.
- Set the bolus amount to 0.0 units. Exit with . No bolus is delivered.




– During the start delay ( $\perp$  blinks):

- Press  or . You hear a melody. Your insulin pump returns to the **RUN** screen.

An alert A8: BOLUS CANCELLED occurs. Press  twice to confirm and turn off the alert. No bolus is delivered.

Please refer to section 10.1.8 “Alert A8: BOLUS CANCELLED” for further information.

– During bolus delivery (countdown of bolus amount):

- Press and hold the  or  key for 3 seconds until you hear a melody. An alert A8: BOLUS CANCELLED occurs. Press  twice to confirm and turn off the alert. Bolus delivery is interrupted. Please refer to section 10.1.8 “Alert A8: BOLUS CANCELLED” for further information.

### **Bolus delivery is interrupted.**

The actual bolus amount delivered prior to the cancellation can be reviewed in the bolus history. Please refer to the section 7.8.1 “Review the bolus history” for further information.

Ensure that the cancellation was intended and program a new Standard bolus, if necessary.



## 7.6 Temporary basal rate (TBR)

Your insulin pump allows you to temporarily increase or decrease your basal rate. This function is helpful in matching changing insulin needs due to increased or decreased activity level, illness or stress.

### Warning

Incorrect programming of your insulin pump may cause inappropriate insulin delivery. Your insulin pump must be programmed with your personal settings prior to starting insulin pump therapy. Do not operate your insulin pump without knowing your personal settings. If you are not sure about your personal settings or you are less experienced, have your doctor or healthcare team check your personal settings.

Typically, a basal rate is set to 100%, but you can

- ▶ increase it to 200% 
- ▶ decrease it to 0% for up to 24 hours 

When you increase or decrease your basal rate, each hourly basal rate for the duration you program is effected.


The characteristic shape of your basal rate profile remains if the change is larger than 0%. With a temporary basal rate of 0% no insulin (except boluses) is delivered during the programmed duration.

Discuss programming a temporary basal rate with your doctor or healthcare team.

**Note** A running temporary basal rate increase or decrease remains active, even if you change the basal rate profile.



## 7.6.1 Program a temporary basal rate



TEMPORARY  
BASAL RATE (TBR)  
d b %

Make sure your insulin pump is in **RUN**.

Press to move to the TEMPORARY BASAL RATE (TBR) menu.

TBR PERCENTAGE  
d b 100%

Select with . The TBR PERCENTAGE screen appears.

**Note** If no temporary basal rate is currently active, the percentage is set to 100%. If a temporary basal rate is active, its duration and percentage appear on the **RUN** screen.

TBR PERCENTAGE  
d b 120%  
→ 04:00

Press to increase or to decrease the basal rate.

**Note** When you increase (or decrease) the temporary basal rate the duration of your previous temporary basal rate increase (or decrease) appears automatically, or when using your insulin pump for the first time, a default value appears.



Press to move to the TBR DURATION.

Press or to correct and set the duration.

**Note** Press to loop between programming the basal rate percentage and duration, if necessary.



Save and exit with .

The programmed temporary basal rate is immediately activated.

**Note** At any screen you have three options to exit:

if you want to **confirm** and **save the changes**

– press . Your insulin pump returns to the **RUN** screen.

if you want to **undo the changes**

– wait for your insulin pump to return to the **RUN** screen (timeout), or


– press + simultaneously (exit feature) to exit the current function screen.

– set percentage to 100% and press .

**Note** During the TBR delivery, the increase or decrease in percentage, the increased or decreased hourly basal rate (hourly basal rate plus/minus TBR percentage) and the remaining time are displayed in the **RUN** screen.

An arrow pointing up symbolizes a temporary basal rate increase, an arrow pointing down a temporary basal rate decrease.

**Note** If an Extended bolus or a MultiWave bolus is active at the same time, the remaining time and amount of the bolus and the amount of the increased or decreased hourly basal rate are displayed in the **RUN** screen.

At the end of a temporary basal rate, an alert A7: TBR OVER occurs .





Press  twice to confirm and turn off the alert.

See the section 10.1.7 “Alert A7: TBR OVER” for further information.

## 7.6.2 Cancel a temporary basal rate

– During programming:

You have three options for cancelling a temporary basal rate during programming:

- If you do not press any key for 20 seconds , your insulin pump returns to the **RUN** screen (timeout).
- Exit to the TEMPORARY BASAL RATE (TBR) menu by pressing  +  simultaneously.
- Set the TBR percentage to 100%. Save and exit with .


**No temporary basal rate is active!**


– During delivery:


You have two options for cancelling a temporary basal rate during delivery:

- Set the TBR percentage to 100%:





Press  to move to the TEMPORARY BASAL RATE (TBR) menu.



Select with . The TBR PERCENTAGE screen appears.



Press  or  to return the basal rate to 100%.

Save and exit with .

- Put your insulin pump into **STOP**.

**The temporary basal rate delivery is interrupted.**



An alert A6: TBR CANCELLED occurs. Press  twice to confirm and turn off the alert.

See the section 10.1.6 “Alert A6: TBR CANCELLED” for further information.

**Caution**

When your insulin pump is in **STOP**, it does not deliver any insulin. Put your insulin pump into **RUN** to continue the insulin delivery.

**Note** If an Extended bolus or a MultiWave bolus is active at the same time, it is also cancelled when your insulin pump is put into **STOP**. An alert A8: BOLUS CANCELLED and an alert A6: TBR CANCELLED occur.

Press  twice to confirm and turn off the first alert. Then the second alert appears on the display. Press  twice to confirm and turn off the second alert. Both alerts are recorded in the alarm history.

See the sections 10.1.6 “Alert A6: TBR CANCELLED” and 10.1.8 “Alert A8: BOLUS CANCELLED” for further information.

Ensure that the cancellation was intended and program a new temporary basal rate (and/or an Extended or MultiWave bolus), if necessary.

## 7.7 Setup menu standard

### Warning

Incorrect programming of your insulin pump may cause inappropriate insulin delivery. Your insulin pump must be programmed with your personal settings prior to starting insulin pump therapy. Do not operate your insulin pump without knowing your personal settings. If you are not sure about your personal settings or you are less experienced, have your doctor or healthcare team check your personal settings.



For all SETUP MENU STANDARD settings the Accu-Chek Spirit can save or exit using the same simple method. For all of the features in this section you may do the following when you are finished with a section:

if you want to **confirm** and **save the changes**

– press . Your insulin pump returns to the **RUN** screen.

if you want to **undo the changes**

– wait for your insulin pump to return to the **RUN** screen (timeout), or


– press  +  simultaneously (exit feature) to exit the current function screen.

### 7.7.1 KeyLock function

The KeyLock function allows you to lock the four keys of your insulin pump and serves as an additional safety measure against unintentional activation of functions (e.g. while sleeping, during contact sports).

To make use of the KeyLock function, it must be turned on.

See the section 3.3 “Keys and key combinations” for further information.



### Turn the KeyLock on or off

Make sure the KeyLock is unlocked (press and hold + simultaneously for 3 seconds until it is unlocked) or turned off.



Press to move to the SETUP MENU STANDARD.



Select with . The current KeyLock status (ON or OFF) appears.





Press or to turn the KeyLock ON or OFF. Save and exit with .


If no further key is pressed for 10 seconds after your insulin pump has returned to the **RUN** screen or **STOP** screen, all four keys are automatically locked if the KeyLock is turned on.

**Note** The only function which can still be activated is the backlight (press .

## Unlock the KeyLock



From the **RUN** screen or **STOP** screen, press + simultaneously (you hear three beeps in **RUN**, respectively one beep in **STOP**) and hold for 3 seconds until you hear again three beeps (in **RUN**), respectively one beep (in **STOP**) to confirm that the KeyLock has been unlocked.



If no further key is pressed for 10 seconds after your insulin pump has returned to the **RUN** screen or **STOP** screen, all four keys are automatically locked.

### 7.7.2 Adjust the beep volume

Your insulin pump beeps when a key is pressed or when an alert or error occurs. The volume of these beeps can be changed.










Press to move to the SETUP MENU STANDARD.

Select with .



Press to move to the BEEP VOLUME screen. The active beep volume is displayed.

Press  or  to adjust the beep volume.  
You can choose between 5 volume levels:




-  no beeps (beeps are turned off)
-  very low
-  low
-  normal
-  maximum

Save and exit with .


**Note** The volume of the **STOP-Warning** is independent of the programmed beep volume. It always occurs at the maximum level. When the beeps are turned off and an alert or error occurs, the beeps will become activated after a short period of time to ensure that you are aware of the alarm or error. These beeps will steadily increase in volume until the maximum level occurs if the alert or error is not cancelled

### 7.7.3 Alarm signals

The Accu-Chek Spirit allows you to choose how you want your pump to alert you when an alarm occurs. You have three options:


-  beep only
-  vibrate only
-  beeps and vibrate



Press  to move to the SETUP MENU STANDARD.

Select with .





Press to move to the ALARM SIGNALS screen. The active alarm signals are displayed. Press or to set the alarm signals as desired.

Save and exit with .

## 7.7.4 Automatic off

### Caution

Discuss the use of the automatic off function with your doctor or healthcare team.

Automatic off is a safety feature which stops all insulin delivery by triggering an error E3: AUTOMATIC OFF if no keys are pressed within a programmed time period in **RUN**.

The automatic off can either be set to OFF or programmed up to 24 hours in 1 hour intervals in the SETUP MENU STANDARD.




Press to move to the SETUP MENU STANDARD.

Select with .



Press to move to the AUTOMATIC OFF screen. The active status is displayed.



Press or scroll to increase or to decrease the duration for the automatic off in 1 hour intervals. If you would like to set the automatic off function to OFF, press or scroll until 0 hours and OFF appear on the display.

Save and exit with .

## 7.7.5 Battery type


### Warning

If you do not set the correct battery type on your insulin pump, the alert A2: BATTERY LOW may not occur in time to provide you sufficient warning to replace your battery.

Your insulin pump has the option of using 1.5 volt AA Alkaline batteries or rechargeable AA NiMH batteries. If you change from one battery type to the other, you must change the battery type accordingly in the SETUP MENU STANDARD.


See the section 4.3 “Battery” for further information on recommended batteries.





Press  to move to the SETUP MENU STANDARD.

Select with .




Press  to move to the BATTERY TYPE screen. The current battery type is displayed.




Press  or  to select the desired battery type (ALKALINE or NiMH [rechargeable] battery).


Save and exit with .

## 7.7.6 Display orientation



Your insulin pump allows you to flip the display orientation by 180° so that you can view it more conveniently depending on how you are wearing the insulin pump.






Press  to move to the SETUP MENU STANDARD.

Select with .











Press  to move to the ORIENTATION screen.



Press  or  to select the desired display orientation.

Save and exit with .

### Caution

If you flip the display orientation of your insulin pump by 180°, the  and  keys will also reverse their functions. In comparison to the display orientation, the upper key will be the  and the lower key will be . This change applies for all  and  functions including turning on the backlight. The  and  keys will not change in their function, but will remain the same, regardless of your screen orientation.

## 7.8 Review data memory (INFORMATION)


The memory of your insulin pump stores all events (i.e. alerts and errors, programming operations and insulin delivery records) for up to 4500 events. This data memory corresponds typically to about the last 90 days of usage and is accessible via Accu-Chek information management products.

The following data can be reviewed directly on your insulin pump:

- Bolus history last 30 boluses
- Alarm history last 30 alerts and errors
- History of daily insulin totals last 30 daily totals of insulin delivered
- Temporary basal rate history last 30 increases and decreases
- Pump Timer remaining time in days until the pump timer expires
- Quick info screen remaining cartridge content

### 7.8.1 Review the bolus history


The bolus history enables review of up to the last 30 bolus deliveries starting with the most recent entry in chronological order.




```
BOLUS HISTORY
┌ 0.7U      01/08
⊙ 19:40   ☒ 29.11
```

Each bolus history screen displays the

- bolus type (└ Standard bolus, ☐ Extended bolus or ☒ MultiWave bolus),
- bolus duration (➔ only displayed for an Extended or a MultiWave bolus),
- time (⊙),
- date (☒) and
- entry number (01 is the most recent) of total entries (e.g. 01/30) of a delivered bolus.




Press to move to the INFORMATION menu.



Select with .

The BOLUS HISTORY screen appears. The most recent bolus information is displayed.



Press or to review the entries.

Exit with .

## 7.8.2 Review the alarm history


The alarm history enables you to review the last 30 alerts or errors starting with the most recent entry in chronological order.

Each alarm history screen displays the


- number and type of the alert or error (e.g.  $\Delta$  A6),
- alert or error (e.g. TBR CANCELLED),
- time (🕒),
- date (📅) and
- entry number (01 is the most recent) of total entries (e.g. 01/30) of an alert or error.



Press  $\square$  to move to the INFORMATION menu.  
Select with  $\checkmark$ .



Press  $\square$  to move to the ALARM HISTORY screen.



Press  $\uparrow$  or  $\downarrow$  to review the entries.

Exit with  $\checkmark$ .

## 7.8.3 Review the daily insulin totals history


This history enables you to review the last 30 daily totals of insulin delivered (from midnight to midnight; basal rate plus bolus deliveries) starting with the most recent entry in chronological order.

Each daily insulin total screen displays the following history


- daily total of insulin delivered ( $\Sigma$ ),
- date (📅) and
- entry number (01 is the most recent) of total entries (e.g. 01/30) of a daily insulin total.



Press to move to the INFORMATION menu.  
Select with .



Press to move to the DAILY TOTALS screen.



Press or scroll or to review the entries.

Exit with .

## 7.8.4 Review the temporary basal rate history


The temporary basal rate history enables you to review the last 30 TBR increases or decreases starting with the most recent entry in chronological order.

Each TBR history screen displays the


- TBR increase () or decrease () in percentage,
- TBR duration (),
- time when TBR was over (),
- date () and
- entry number (01 is the most recent) of total entries (e.g. 01/30) of a delivered temporary basal rate.



Press to move to the INFORMATION menu.  
Select with .



Press to move to the TBR HISTORY screen.



Press or scroll or to review the entries.

Exit with .

## 7.8.5 Review the time remaining

Your insulin pump has been designed to be a highly reliable system that will provide worry-free performance for a long period of time. To ensure optimal performance, a pump timer will limit the operating time of the pump. This timer counts the number of days of operation remaining in the Accu-Chek Spirit insulin pump.


Before the pump timer expires, an alert occurs to remind you of the upcoming end of operating time so that you can make the necessary arrangements. The pump timer symbol (⌚) appears in the **RUN** screen or **STOP** screen as a reminder. When the timer expires, an error occurs and your insulin pump goes into **STOP**. It can no longer be put into **RUN**. See the sections 10.1.5 “Alert A5: PUMP TIMER” and 10.2.5 “Error E5: END OF OPERATION” for further information.


The time remaining in days until the pump timer expires can be reviewed in the INFORMATION menu.



Press  to move to the INFORMATION menu.

Select with .





Press  to move to the PUMP TIMER screen. The time remaining in days until the pump timer expires is shown.

Exit with .





## 7.8.6 Quick info screen



With the press of the  key, your insulin pump will display the remaining cartridge content in units and the insulin concentration. Symbols for activated functions and features (like the beep and vibration alarm signals) are displayed.

### To get to the quick info screen

- ▶ Press  from the **RUN** screen or **STOP** screen.
- ▶ After 8 seconds , your insulin pump automatically returns to the **RUN** screen or **STOP** screen.

---

## 8 ADVANCED user menu

This chapter introduces you to the advanced functions of your insulin pump. The ADVANCED user menu activates all of the therapy options you have with your pump.

The Accu-Chek Spirit insulin pump allows you to choose the menus that you wish to view. This allows you to adjust your pump's features to better fit your level of diabetes management.

To turn on all functions:

- ▶ choose the SELECT USER MENU function,
- ▶ select the ADVANCED user menu.

To use the functions individually:



- ▶ choose the SELECT USER MENU function,
- ▶ select the CUSTOM user menu and turn the individual features on or off as desired using Accu-Chek Insulin Pump Configuration Software.

As was done with the STANDARD user menu settings on the Accu-Chek Spirit, you can save or exit your programming using the same three simple methods. For all of the features in this section you may do the following when you are finished:

if you want to **confirm** and **save the changes**

- press . Your insulin pump returns to the **RUN** screen.

if you want to **undo the changes**

- wait for your insulin pump to return to the **RUN** screen (timeout), or
- press  +  simultaneously (exit feature) to exit the current function screen.

**Warning**

Incorrect programming of your insulin pump may cause inappropriate insulin delivery. Your insulin pump must be programmed with your personal settings prior to starting insulin pump therapy. Do not operate your insulin pump without knowing your personal settings. If you are not sure about your personal settings or you are less experienced, have your doctor or healthcare team check your personal settings.

## 8.1 Extended bolus


The Extended bolus is only available in the ADVANCED user menu and may be turned on or off in the CUSTOM user menu. The Extended bolus function allows you to program a bolus to be delivered over a period of time. This function can be helpful during long meals, during dinners or receptions, or having meals which are digested slowly. The use of the Extended bolus may also be appropriate with people who have gastroparesis (delayed digestion).

The duration of bolus delivery can be programmed in 15 minute ⌚ intervals up to 12 hours ⌚ and begins immediately after confirmation.


**Note** You can add a Standard bolus to an ongoing Extended bolus. If you cancel the Standard bolus, the Extended bolus remains running.

During the delivery of an Extended bolus, the programming of another Extended bolus or a MultiWave bolus is blocked. If you want to program a different Extended bolus, put your insulin pump into **STOP** to cancel the current bolus and program a new one.


## 8.1.1 Program an Extended bolus



Press to move to the EXTENDED BOLUS menu.




Select with . The BOLUS AMOUNT screen appears.



Press to increase or to decrease the bolus amount.

The bolus duration of the last Extended bolus delivered appears, or when using your insulin pump for the first time, a default value appears.


Press to move to the BOLUS DURATION.




Press to increase or decrease the bolus duration.

**Note** Press to loop between programming the bolus amount and bolus duration, if necessary.


Check the bolus amount and duration programmed on the display.

Press  to confirm the bolus amount and duration.

You hear a melody and the bolus delivery begins within the next 3 minutes.



**Note** During the entire bolus delivery, the remaining time and amount of the Extended bolus and the current hourly basal rate are displayed in the **RUN** screen.







**Note** If a temporary basal rate is active at the same time, the remaining time and amount of the bolus and the amount of the increased or decreased hourly basal rate are displayed in the **RUN** screen.

### 8.1.2 Cancel an Extended bolus


– During programming (the bolus amount or bolus duration blinks):

You have three options for cancelling an Extended bolus during programming:

- If you do not press any key for 20 seconds , your insulin pump returns to the **RUN** screen (timeout).
- Exit to the EXTENDED BOLUS menu by pressing  +  simultaneously.
- Set the bolus amount to 0.0 units and press .

**No bolus is delivered.**

– During bolus delivery:

- If the Extended bolus delivery has begun, it can be cancelled by putting your insulin pump into **STOP**. This cancels the bolus delivery and an alert A8: BOLUS CANCELLED occurs. Press  twice to confirm and turn off the alert. See the section 10.1.8 “Alert A8: BOLUS CANCELLED” for further information.  
**Bolus delivery is interrupted.**



The bolus amount delivered up to the cancellation can be reviewed in the BOLUS HISTORY screen.

See the section 7.8.1 “Review the bolus history” for further information.

### Caution

When your insulin pump is in **STOP**, it does not deliver any insulin. Put your insulin pump into **RUN** to continue the insulin delivery.

**Note** If a temporary basal rate is active at the same time, it is also cancelled when your insulin pump is put into **STOP**. An alert A8: BOLUS CANCELLED and an alert A6: TBR CANCELLED occur.

Press  twice to confirm and turn off the first alert. Then the second alert appears on the display. Press  twice to confirm and turn off the second alert. Both alerts are recorded in the alarm history.

See the sections 10.1.6 “Alert A6: TBR CANCELLED” and 10.1.8 “Alert A8: BOLUS CANCELLED” for further information.

Ensure that the cancellation was intended and program a new Extended bolus (and/or temporary basal rate), if necessary.

## 8.2 MultiWave bolus

The MultiWave bolus is only available in the ADVANCED user menu and depending on your personal settings in the CUSTOM user menu. It is designed to better simulate a body's insulin delivery. It combines an immediate bolus delivery followed by an Extended bolus delivery. This function can be helpful when having meals with both rapid and slowly absorbed carbohydrates (e.g. pizza).

The duration of bolus delivery can be programmed in 15 minute ⌚ intervals up to 12 hours ⌚ and begins immediately after confirmation.



**Note** You can add a Standard bolus to an ongoing MultiWave bolus.

If you cancel the Standard bolus, the MultiWave bolus remains running.

During the delivery of a MultiWave bolus, the programming of another MultiWave bolus or an Extended bolus is blocked. If you need a different MultiWave bolus, put your insulin pump into **STOP** to cancel the current bolus and program a new one.




## 8.2.1 Program a MultiWave bolus



Press **➤** to move to the MULTIWAVE BOLUS menu.




Select with **⬅**. The BOLUS AMOUNT screen appears.





Press **⬆** to increase or **⬇** to decrease the total bolus amount.


With the first press of the **⬆** key, the bolus duration of the last MultiWave bolus appears, or when using your insulin pump for the first time, a default value appears.

In addition, as you increase or decrease the total bolus amount, the suggested immediate bolus amount increases or decreases proportionally.





Press  to move to the IMMEDIATE BOLUS amount.







Press  or  to correct and set the immediate bolus amount.


Press  to move to the BOLUS DURATION.



Press  or  to correct and set the Extended bolus duration in increments of 15 minutes  up to 12 hours .

**Note** Press  to loop between programming the total bolus amount, immediate bolus amount and Extended bolus duration, if necessary.

- Check the total and immediate bolus amounts and the Extended bolus duration programmed on the display. The set values can be viewed in the total bolus amount, immediate bolus and bolus duration screen.
- Press  to confirm both bolus amounts and the bolus duration.
- For 5 seconds  (bolus delivery start delay), the MultiWave bolus symbol () blinks.
- Your insulin pump beeps three times and the immediate bolus delivery begins. The countdown of the remaining bolus continues to appear on the display.








**Note** Once the immediate bolus is delivered, the remaining time, bolus amount and the current hourly basal rate are displayed in the **RUN** screen.




If a temporary basal rate is active at the same time, the remaining time and amount of the bolus and the amount of the increased or decreased hourly basal rate are displayed in the **RUN** screen.

## 8.2.2 Cancel a MultiWave bolus



– During programming (the bolus amount or bolus duration blinks):  
You have three options for cancelling a MultiWave bolus during programming:

- If you do not press any key for 20 seconds , your insulin pump returns to the **RUN** screen (timeout).
  - Exit to the MULTIWAVE BOLUS menu by pressing  +  simultaneously.
  - Set the total bolus amount to 0.0 units and press .
- No bolus is delivered.**


– During the start delay ( blinks):

- Press  or . Your Accu-Chek Spirit beeps and vibrates. Your insulin pump returns to the **RUN** screen. An alert A8: BOLUS CANCELLED occurs. Press  twice to confirm and turn off the alert. See the section 10.1.8 “Alert A8: BOLUS CANCELLED” for further information.
- No bolus is delivered.**

– During bolus delivery:

- The **immediate delivery** can be cancelled by pressing and holding  or  for 3 seconds until you hear a melody. This cancels the whole bolus (immediate and extended bolus delivery). An alert A8: BOLUS CANCELLED occurs.

**Note** If the alert A8: BOLUS CANCELLED does not occur, the bolus delivery was not cancelled. Cancel the extended delivery (see below).

- The **extended delivery** can be cancelled by putting your insulin pump into **STOP**. This cancels the Extended bolus. An alert A8: BOLUS CANCELLED occurs. Press  twice to confirm and turn off the alert.

See the section 10.1.8 “Alert A8: BOLUS CANCELLED” for further information.

**Bolus delivery is interrupted.**



The actual bolus amount delivered prior to the cancellation can be reviewed in the bolus history.

See the section 7.8.1 “Review the bolus history” for further information.

**Caution**

When your insulin pump is in **STOP**, it does not deliver any insulin. Put your insulin pump into **RUN** to continue the insulin delivery.


**Note** If a temporary basal rate is active at the same time, it is also cancelled when your insulin pump is put into **STOP**. An alert A8: BOLUS CANCELLED and an alert A6: TBR CANCELLED occur.

Press  twice to confirm and turn off the first alert. Then the second alert appears on the display. Press  twice to confirm and turn off the second alert. Both alerts are recorded in the alarm history.

See the sections 10.1.6 “Alert A6: TBR CANCELLED” and 10.1.8 “Alert A8: BOLUS CANCELLED” for further information.

Ensure that the cancellation was intended and program a new MultiWave bolus (and/or temporary basal rate), if necessary.

## 8.3 Basal rate profiles

Your insulin pump offers the option of delivering up to five  different basal rate profiles in order to easily meet changing insulin needs (for example Monday-Friday versus Exercise Day versus Sleeping Late Day).


Discuss programming additional basal rate profiles with your doctor or healthcare team.

Consult your doctor or healthcare team before changing basal rate profiles as this may change how you utilize your insulin pump.

If you and your doctor or healthcare team decide that additional basal rate profiles are not necessary you may hide basal rate profiles 2 to 5 in the CUSTOM user menu using Accu-Chek Insulin Pump Configuration Software Pro for healthcare professionals (see the Accu-Chek Insulin Pump Configuration Software user manual for additional information).

### 8.3.1 Set an additional basal rate profile




1. Make sure the PROGRAM BASAL RATE PROFILE menu to be set is turned on.
2. Press  to move to the desired PROGRAM BASAL RATE PROFILE menu.

**Note** Only basal rate profiles that are accessible in your current user menu are available.

3. The hourly basal rates for each additional basal rate profile are set in the same way as for basal rate profile 1. See the section 7.2.1 “Program a basal rate profile” for further information.


### 8.3.2 Select a basal rate profile



Make sure the CHANGE BASAL RATE PROFILE menu is turned on.




Press to move to the CHANGE BASAL RATE PROFILE menu.



Select with .  
The current basal rate profile is displayed.

**Note** Only basal rate profiles that are accessible in your current user menu are available.



Press or to select the desired basal rate profile. The newly selected basal rate profile and its basal rate total appear.

Save and exit with .  
The new basal rate profile is active immediately.

**Note** When you change the user menu your current basal rate profile might not appear. Make sure that your desired basal rate profile numbers are activated with Accu-Chek Insulin Pump Configuration Software.

Should your desired basal rate profile number(s) not appear,

- select the ADVANCED user menu, or
- add or have your doctor or healthcare team add your desired basal rate profile number(s) to the CUSTOM user menu (see Accu-Chek Insulin Pump Configuration Software user manual), or
- program the desired basal rate profile to an available basal rate profile number.

A running temporary basal rate increase or decrease remains active, even if you change the basal rate profile on your pump.

## 8.4 Alarm clock


The alarm clock can be set for a single alarm and for multiple alarms. Multiple alarms are repeated every day. These alarms may be utilized to remind you of blood glucose testing or other important personal events.

### Set single or multiple alarms

Make sure the ALARM CLOCK menu is turned on.



Press  to move to the ALARM CLOCK menu.



Select with .

The current alarm clock status (OFF, ONCE [single alarm] or EVERY DAY) and alarm time appear.



Press or to set the alarm clock to OFF (), to a single alarm () or to every day ().



Press to move to the SET ALARM HOUR screen.

Press or to set the hour.



Press to move to the SET ALARM MINUTE screen.

Press or to set the minute.

Save and exit with .

## To turn off the alarm clock

When the alarm clock goes off, an alert A4: ALARM CLOCK occurs. Press twice to confirm and turn off the alert.

See the section 10.1.4 “Alert A4: ALARM CLOCK” for further information.

## 8.5 Setup menu advanced

### 8.5.1 Time format

The time can be set in the following formats:


European: 24-hour clock (00:00–23:59), e.g. 13:39

American: 12-hour clock (AM/PM), e.g. 1:39 PM




Make sure the SETUP MENU ADVANCED is turned on.





Press  to move to the SETUP MENU ADVANCED.

Select with .



The TIME FORMAT screen appears.

Press  or  to select the desired time format.

Save and exit with .

## 8.5.2 Date format


The date can be set in the following formats:

European: dd.mm.yy, e.g. 26.04.06

American: mm/dd/yy, e.g. 04/26/06

Make sure the SETUP MENU ADVANCED is turned on.





Press  to move to the SETUP MENU ADVANCED.

Select with .





Press  to move to the DATE FORMAT screen.




Press  or  to select the desired date format.

Save and exit with .

### 8.5.3 Bolus increment



#### Warning


The bolus increment that can be programmed into your insulin pump determines the “Quick” Standard bolus amount using the  and  keys of your insulin pump. To ensure correct insulin delivery be sure the bolus increment is set appropriate for your therapy.


To address patients that may need large or small amounts of insulin, the bolus increment for the Accu-Chek Spirit is adjustable. The “Quick” Standard bolus is initially set to 0.5 units  per key press but can be changed in the SETUP MENU ADVANCED or by using Accu-Chek Insulin Pump Configuration Software to 0.1, 0.2, 0.5, 1.0 or 2.0 units. Units for the “Scroll” Standard bolus, which may be adjusted by simply holding down the  or  keys, will move at a constant 0.1 units. See the section 7.5.1 “Standard bolus” for further information.



Make sure the SETUP MENU ADVANCED is turned on.



Press  to move to the SETUP MENU ADVANCED.  
Select with .




Press  to move to the BOLUS INCREMENT screen.

Press  or  to select the desired bolus increment.

Save and exit with .


## 8.5.4 Prime quantity


The prime quantity necessary to completely fill an infusion set depends on the length of the infusion set tubing. The shorter your infusion set tubing, the less insulin will be required to prime the infusion set. The default setting for the prime quantity is 25  units of insulin.

**Note** The amount of insulin used for the priming is not added to the history of daily insulin totals.


You can press any of your insulin pump's keys to stop priming.


Make sure the SETUP MENU ADVANCED is turned on.





Press  to move to the SETUP MENU ADVANCED.

Select with .



Press  to move to the PRIME QUANTITY screen.

Press  or  to select the desired prime quantity.


Save and exit with .

## 8.5.5 Lock basal rate profiles


Your insulin pump allows you to change your hourly basal rates as needed. This feature provides additional protection against accidentally changing a basal rate during normal operation. When this function is activated, your hourly basal rates cannot be changed.


Make sure the SETUP MENU ADVANCED is turned on.







Press  to move to the SETUP MENU ADVANCED.









Select with .



Press  to move to the BASAL RATE LOCK screen.

Press  or  to turn the basal rate lock ON () or OFF (.

Save and exit with .


**Note** If the basal rate lock is set to ON () , the programming of the basal rate profiles 1, 2, 3, 4 and 5 is blocked. To check if the basal rate lock is ON, move with  to a PROGRAM BASAL RATE menu and select with . The  symbol in the BASAL RATE TOTAL screen indicates the lock of the basal rate programming (The  symbol appears in the BASAL RATE TOTAL screen only). Using  the hourly basal rates can be viewed but they can not be changed by using  and .


## 8.5.6 Language

Your insulin pump is programmed in several languages. To use the language of your choice, do the following:


Make sure the SETUP MENU ADVANCED is turned on.



Press  to move to SETUP MENU ADVANCED.

Select with .



Press  to move to the LANGUAGE screen.


Press  or  to select the desired language.


Save and exit with .


## 8.5.7 Display contrast

Your insulin pump allows you to adjust the display contrast.


Make sure the SETUP MENU ADVANCED is turned on.



Press  to move to the SETUP MENU ADVANCED.

Select with .



Press  to move to the DISPLAY CONTRAST screen.

Press  or  to select the desired contrast.


Save and exit with .

## 9 Data transfer

### Caution


If data transfer between your PC and your insulin pump is disrupted in any way, the configuration may be incomplete and an error E12: DATA INTERRUPTED occurs. The data transfer must be completed successfully before you can put your insulin pump into **RUN**. See section 10.2.11 “Error E12: DATA INTERRUPTED” for further information.

The built-in infrared interface in the bottom cover of your insulin pump allows for wireless data transfer from your insulin pump to a PC. Pump programming information may be exchanged between Accu-Chek Insulin Pump Configuration Software on your computer and your Accu-Chek Spirit insulin pump using this data transfer. Infrared communication with your PC is established via an infrared adapter, which is connected to your PC's RS232 serial interface.



Refer to the user manuals of the corresponding Accu-Chek information management products for detailed infrared adapter information and for further instructions on the data transfer feature.

## Setup for data transfer




Prepare your PC and Accu-Chek Insulin Pump Configuration Software for data transfer (see Accu-Chek Insulin Pump Configuration Software user manual).




Make sure your insulin pump is in **STOP**.



Press  to move to the COMMUNICATION menu.




Select with .



The DATA TRANSFER screen appears and a melody signals that your insulin pump is ready to receive data from or deliver data to your PC. Follow the instructions in the Accu-Chek Insulin Pump Configuration Software user manual. When the data transfer is finished, press  to return to **STOP**.



Check the configuration on your insulin pump to make sure that all parameters are written correctly.

Set up your insulin pump with a cartridge, adapter and a new infusion set and put it into **RUN**, if necessary.

**Note** Your insulin pump automatically returns to **STOP** 15 minutes  after the last data transfer respectively when the **DATA TRANSFER** screen has been selected.

You may press  +  simultaneously to exit while data is being transferred. The data transfer is interrupted and an error E12: DATA INTERRUPTED occurs.

See section 10.2.11 “Error E12: DATA INTERRUPTED” for further information.

---

## 10 Alerts and errors

When your insulin pump is in use, a two microprocessor safety system continuously monitors and controls the insulin pump. Over 9 million safety checks are performed daily to secure all relevant functions. If your insulin pump detects a deviation of its normal state, an alert (warning instruction) or an error (error message) will occur. Beeps, vibrations and messages on the display inform you of the status or malfunction of your insulin pump. **Check the display of your insulin pump at least every three hours during the day, before you go to sleep, and especially if for any reason you might be unable to hear the beeps or to feel the vibrations.** This is the only way you will be notified of any changes of your insulin pump in a suitable amount of time.

In general, alerts and errors are signalled by both beeps and vibrations. If you have turned off either the beeps or vibrations, only the turned on alarm signal goes off when the alert or error occurs for the first time. If you do not confirm the alert or error within 60 seconds, it occurs again with both beeps and vibrations. When the beeps are turned off or programmed at a lower than maximum volume level, the beep volume increases within 10 seconds up to the maximum volume level. The beeps and vibrations will continue at this level until confirmed.

**Note** The volume of the **STOP-Warning** is independent of the programmed beep volume. It always occurs at the maximum level.


The KeyLock is automatically unlocked each time an alert or error occurs.

## To confirm an alert or error









Press  to turn off the beeps and vibrations.

The alert or error code remains on the display.

Confirm that you have seen and understood the alert or error code and messages by pressing  again.

Take the appropriate action, if necessary.

If more than one alert and/or error occur at the same time, press  twice for each alert or error to turn off and confirm.

When you confirm A1: CARTRIDGE LOW , E1: CARTRIDGE EMPTY , A2: BATTERY LOW , E2: BATTERY DEPLETED , A5: PUMP TIMER  or E5: END OF OPERATION, the symbol for that alert and error stays in the **RUN** screen or **STOP** screen as a reminder. For the other alerts and errors, the alarm information is cleared from the display and remains only in the alarm history and the event memory.



After the occurrence of an alert or error, it is important to make sure that your insulin pump is returned to **RUN**, if desired.

## 10.1 Alerts

### Warning

After the occurrence of an alert, your insulin pump may be in **STOP** and the insulin delivery may be interrupted. In order to maintain insulin delivery, you must act immediately according to the instructions given for each error code or put your insulin pump into **RUN** to continue the insulin delivery.

### 10.1.1 Alert A1: CARTRIDGE LOW

If the remaining content in the insulin cartridge reaches 20 units  of insulin, an alert A1: CARTRIDGE LOW occurs. Your insulin pump displays that the cartridge is almost empty. Be prepared to replace the cartridge as quickly as possible. The cartridge low symbol () appears in the **RUN** or **STOP** screens as a reminder until you change the cartridge.





Press  twice to turn off and confirm the alert.

Change the cartridge before it is empty.

Put your insulin pump into **RUN**, if necessary.

### 10.1.2 Alert A2: BATTERY LOW

Your insulin pump checks the battery voltage before each insulin delivery (at least every 3 minutes). If it falls below a certain voltage value, the alert A2: BATTERY LOW occurs. Be prepared for a battery change as quickly as possible. The battery low symbol () appears in the **RUN** screen or **STOP** screen as a reminder until you insert a new battery.




Press  twice to turn off and confirm the alert.


Replace the battery within the next hours.

Put your insulin pump into **RUN**, if necessary.

### 10.1.3 Alert A3: REVIEW TIME AND DATE

If the battery is removed for more than 1 hour, it may be necessary to enter the time and date again. Other Accu-Chek Spirit settings and the event memory are not affected and are preserved, regardless of battery voltage and time your insulin pump has been without a battery.



Press  twice to turn off and confirm the alert.


Review the time and date and correct them, if necessary.

Put your insulin pump into **RUN**, if necessary.

### 10.1.4 Alert A4: ALARM CLOCK

An alert A4: Alarm clock will occur at the time that you have programmed into the pump.



Press  twice to turn off and confirm the alert.

Put your insulin pump into **RUN**, if necessary.


**Note** If you have set the alarm occurrence to **EVERY DAY**, the same alert occurs again every 24 hours. Turn the alarm clock function to **OFF**, if you do not want additional alarms at that same time.

### 10.1.5 Alert A5: PUMP TIMER

To ensure that insulin delivery is optimized, your pump features a pump timer.

To inspect the pre-set lifetime of your Accu-Chek Spirit pump before you start using it, please refer to the “Accu-Chek Spirit Pump Protocol” included in the starter kit. During use, the remaining lifetime in days until the expiration date of your insulin pump can be reviewed in the **INFORMATION** menu (see “Review the time remaining”).


To ensure you are aware of this timing feature, the alert A5: PUMP TIMER will remind you 60 days prior to the end of operation to check your pump’s programmed remaining operation time.




When the PUMP TIMER has reached zero, the pump will change into **STOP** and will no longer operate.

An **E5: END OF OPERATION** will appear on the display.




Press  twice to turn off and confirm the alert.

Put your insulin pump into **RUN**.



Check the INFORMATION screen on your pump for the PUMP TIMER display. The amount shown is the remaining number of programmed days in your insulin pump.

After the alert A5: PUMP TIMER occurs, the Accu-Chek Spirit insulin pump will display the  symbol to remind you to check the timer. Contact your local HelpDesk to discuss your options for further use of the Accu-Chek Spirit insulin pump.

**Note** The alert A5: PUMP TIMER and the error E5: END OF OPERATION are based on running time (period with pump in **RUN**) and not calendar years. Your actual warranty is based upon the date the pump was first purchased.

### 10.1.6 Alert A6: TBR CANCELLED (Temporary Basal Rate cancelled)

An alert A6: TBR CANCELLED occurs when a temporary increase or decrease of the basal rate has been cancelled. The original basal rate (100%) will resume automatically.



Press  twice to turn off and confirm the alert.

Put your insulin pump into **RUN**, if necessary.

Ensure that the cancellation was intended and, if necessary, program a new temporary basal rate.



### 10.1.7 Alert A7: TBR OVER (Temporary Basal Rate over)

An alert A7: TBR OVER occurs when a temporary basal rate has ended. Once the alert has been confirmed, the original basal rate (100%) will resume automatically.



Press  twice to turn off and confirm the alert.

Decide whether a further temporary basal rate change is appropriate and program, if necessary.

**Note** If desired, the alert A7: TBR OVER can be turned off using Accu-Chek Insulin Pump Configuration Software. If turned off, the alert does not occur when a temporary basal rate ends.

### 10.1.8 Alert A8: BOLUS CANCELLED

If a bolus is cancelled during the start delay or after delivery has begun, an alert A8: BOLUS CANCELLED occurs.



Press  twice to turn off and confirm the alert.

Put your insulin pump into **RUN**, if necessary.

Ensure that the cancellation was intended and, if necessary, program a new bolus.

**Note** The actual bolus amount delivered prior to the cancellation can be reviewed in the bolus history. See the section 7.8.1 “Review the bolus history” for further information.

## 10.2 Errors


### Warning

After the occurrence of an error, your insulin pump will be in **STOP** and the insulin delivery is interrupted. In order to maintain insulin delivery, you must act immediately according to the instructions given for each error code.

### 10.2.1 Error E1: CARTRIDGE EMPTY

If the cartridge is empty, an error E1: CARTRIDGE EMPTY occurs. Change the cartridge immediately.



Press  twice to turn off and confirm the error.

Disconnect or remove the infusion set from your infusion site.

Change the cartridge and prime the new infusion set.

Reconnect the infusion set (tubing) at your infusion site or insert the new infusion set.


When finished, put your insulin pump into **RUN**.


## 10.2.2 Error E2: BATTERY DEPLETED

### Warning

After inserting a new battery, always ensure that the time and date of your insulin pump are set correctly. Incorrect programming of the time and date may cause incorrect insulin delivery. If you or your healthcare team store and analyse your therapy data electronically and the time and date of the devices used are not set identically, the gathered data might not be meaningful.

If the battery is depleted, an error E2: BATTERY DEPLETED occurs. Replace it immediately.



Press  twice to turn off and confirm the error.

Replace the battery.

When finished, put your insulin pump into **RUN**.

**Note** Your insulin pump stores the settings you made and saved before the battery was depleted. If your insulin pump remains without power for more than one hour, check the time and date.

## 10.2.3 Error E3: AUTOMATIC OFF

This safety function interrupts insulin delivery if no keys are pressed within a certain period of time.

Refer to section 7.7.4 “Automatic off” for more information.



Press  twice to turn off and confirm the error.

Put your insulin pump into **RUN**.

## 10.2.4 Error E4: OCCLUSION

### Warning

If an error E4: OCCLUSION occurs, immediately check your blood glucose level, because insulin delivery has been interrupted. If your blood glucose level is high, take appropriate actions according to your doctor's or healthcare team's instructions.

An error E4: OCCLUSION occurs if insulin (a maximum of 3.5 units of insulin when using plastic cartridges or a maximum of 6 units of insulin when using Aventis Insuman Infusat glass cartridges\*) is failed to be delivered. This may be caused by:


- a blocked infusion set or
- a reused cartridge or
- a dirty or damaged piston rod.

**This error is dependent upon several factors including your current hourly basal rate and boluses.** For example, a blocked infusion set – with a basal rate of 2.0 U/h (no additional boluses) would have an error E4: OCCLUSION in approximately three hours or less using Accu-Chek 3.15 ml plastic cartridges, depending on the type and severity of the blockage.

\* Insuman® Infusat produced by Aventis Pharma Deutschland GmbH, a member and a brand of the sanofi-aventis group.

## With a disconnectable infusion set do the following:



Press  twice to turn off and confirm the error. Disconnect the infusion set tubing from the infu-sion site. Activate the PRIMING function.

If no error E4: OCCLUSION occurs, the **occlusion is within the infusion set cannula**. Change the infusion set cannula and reconnect the infusion set tubing. Put you insulin pump into **RUN**.

If however still an error E4: OCCLUSION occurs, remove the infusion set tubing from the insulin pump. Activate the PRIMING function without an attached infusion set. Hold your insulin pump with the adapter pointing downwards to avoid that insulin flows through the adapter into the cartridge compartment.

If no error E4: OCCLUSION occurs, the **occlusion was within the infusion set tubing**. Change the infusion set tubing. Prime the new infusion set tubing and connect it to your infusion site. Put you insulin pump into **RUN**.


If however still an error E4: OCCLUSION occurs, the cartridge can be the cause of the occlusion. Remove the cartridge and activate the PRIMING function (without cartridge and without attached adapter and infusion set).

If no error E4: OCCLUSION occurs, the **occlusion was within the cartridge**. Activate the CHANGE THE CARTRIDGE function. Insert a new filled cartridge with attached adapter and a new infusion set tubing into the insulin pump. Prime the new infusion set tubing and connect it to your infusion site. Put you insulin pump into **RUN**.

If all these corrective actions do not help and the error E4: OCCLUSION continues to occur, contact your doctor or healthcare team for an alternate therapy plan. Contact your local Customer Care Area for further assistance.

## With a non-disconnectable infusion set do the following:



Press  twice to turn off and confirm the error. Remove your infusion set from the infusion site and afterwards from your insulin pump. Activate the PRIMING function without an attached infusion set. Hold your insulin pump with the adapter pointing downwards to avoid that insulin flows through the adapter into the cartridge compartment.

If no error E4: OCCLUSION occurs, the **occlusion was within the infusion set**. Change the infusion set. Prime the new infusion set and insert it in a new infusion site. Put you insulin pump into **RUN**.

If however still an error E4: OCCLUSION occurs, the cartridge can be the cause of the occlusion. Remove the cartridge and activate the PRIMING function (without cartridge and without attached adapter and infusion set).

If no error E4: OCCLUSION occurs, the **occlusion was within the cartridge**. Activate the CHANGE THE CARTRIDGE function. Insert a new filled cartridge with attached adapter and a new infusion set into the insulin pump. Prime the new infusion set and insert it into a new infusion site. Put you insulin pump into **RUN**.


If all these corrective actions do not help and the error E4: OCCLUSION continues to occur, contact your doctor or healthcare team for an alternate therapy plan. Contact your local Customer Care Area for further assistance.

## 10.2.5 Error E5: END OF OPERATION

The error E5: END OF OPERATION indicates that your pump timer has reached zero. The pump will no longer operate. Before this error occurs an alert A5: PUMP TIMER will occur to ensure that you are able to replace your pump within sufficient time.



When the PUMP TIMER has reached zero, it will change into **STOP** and will no longer operate. An E5: END OF OPERATION will appear on the display.

Press  twice to turn off and confirm the error.

To ensure that insulin delivery is optimized, your pump features a pump timer.

To inspect the pre-set lifetime of your Accu-Chek Spirit pump before you start using it, please refer to the “Accu-Chek Spirit Pump Protocol” included in the starter kit. During use, the remaining lifetime in days until the expiration date of your insulin pump can be reviewed in the INFORMATION menu (see “Review the time remaining”).



Contact your local HelpDesk to discuss your options for further use of the Accu-Chek Spirit insulin pump.

**Note** The alert A5: PUMP TIMER and the error E5: END OF OPERATION are based on running time (period with pump in **RUN**) and not calendar years. Your actual warranty is based upon the date the pump was first purchased.

## 10.2.6 Error E6: MECHANICAL ERROR



With every insulin delivery (at least every 3 minutes) and each time you put your insulin pump into **RUN**, the safety system of your insulin pump reviews all mechanical parameters. If a mechanical error is discovered during this process, the insulin delivery is stopped and an error E6: MECHANICAL ERROR occurs.



Press  twice to turn off and confirm the error. Press  to get to the Quick Info Screen. Note down the remaining cartridge content.

Remove or disconnect the infusion set from your infusion site.

Remove the cartridge, adapter, infusion set (tubing) and battery from your insulin pump. Insert a new battery. Rewind the piston rod.

Start the CHANGE THE CARTRIDGE function. When the cartridge volume is displayed, press or scroll  or  to move the plunger rod forward until it reaches the remaining cartridge content, which was noted down.


Set up your insulin pump with the cartridge and a new infusion set. Prime the infusion set. Reconnect the infusion set (tubing) at your infusion site.

When finished put your insulin pump into **RUN**.

Should the error E6: MECHANICAL ERROR occur again in spite of these actions, contact your doctor or healthcare team for an alternate therapy plan. Contact your local Customer Care Area for further assistance.



## 10.2.7 Error E7: ELECTRONIC ERROR

The safety system of your insulin pump continuously monitors the functioning of your insulin pump. If an electronic error is discovered, the insulin delivery is stopped and an error E7: ELECTRONIC ERROR occurs. You **cannot** turn off and confirm an error E7: ELECTRONIC ERROR by pressing , remove the battery instead.



Remove or disconnect the infusion set from your infusion site.

Insert a new battery into the pump.

Prime the infusion set.

Reconnect the infusion set (tubing) at your infusion site or insert the new infusion set.

When finished, put your insulin pump into **RUN**.


If an error E7: ELECTRONIC ERROR occurs during programming (i.e. basal rate, time and date) check your settings.

If the error E7: ELECTRONIC ERROR occurs again, contact your local Customer Care Area for assistance. Contact your doctor or healthcare team for an alternate therapy plan, if necessary.

## 10.2.8 Error E8: POWER INTERRUPT

A power interrupt may occur when you change the battery without putting your insulin pump into **STOP** first, or when your insulin pump was dropped.



Press  twice to turn off and confirm the error.

Your insulin pump goes into **STOP**.

Check the time and date and correct, if necessary.

If a bolus and/or a temporary basal rate was interrupted by the error, review the bolus and/or the temporary basal rate history for the delivered amount and duration.

When finished, put your insulin pump into **RUN**, if necessary.


If necessary, program a new bolus and/or temporary basal rate.

If an error E8: POWER INTERRUPT occurs during programming, check your settings.

## 10.2.9 Error E10: CARTRIDGE ERROR

If the CHANGE THE CARTRIDGE function was not correctly carried out, an error E10: CARTRIDGE ERROR will occur.



Press  twice to turn off and confirm the error.

Remove or disconnect the infusion set from your infusion site.

Remove the cartridge. Start the CHANGE THE CARTRIDGE menu. Reinsert the cartridge.

Prime the infusion set.

Reconnect the infusion set (tubing) at your infusion site or insert the new infusion set.

When finished, put your insulin pump into **RUN**.

If an error E10: CARTRIDGE ERROR occurs during the winding of the piston rod, insert a new battery into the pump and repeat the steps above.


## 10.2.10 Error E11: SET NOT PRIMED

### Warning

Never prime an infusion set that is connected to your body. You risk uncontrolled insulin delivery into your body. With disconnectable infusion sets, make sure you have disconnected the tubing from your site prior to priming. Always follow the instructions provided with your infusion set.

If you changed the cartridge and infusion set but did not prime the infusion set, an error E11: SET NOT PRIMED occurs when you try to put your insulin pump into **RUN**.



Press  twice to turn off and confirm the error.

Remove or disconnect the infusion set from your infusion site.

Prime the infusion set.

Reconnect the infusion set (tubing) at your infusion site or insert the new infusion set.

When finished, put your insulin pump into **RUN**.

## 10.2.11 Error E12: DATA INTERRUPTED


### Caution

If data transfer between your PC and your insulin pump is disrupted in any way, the configuration may be incomplete and an error E12: DATA INTERRUPTED occurs. The data transfer must be completed successfully before you can put your insulin pump into **RUN**.

When data transfer between your insulin pump and a PC is interrupted an error E12: DATA INTERRUPTED occurs. Restart communication between your insulin pump and the PC and make sure that the data transfer is completed correctly. Before you can put your insulin pump into **RUN**, the data transfer must be successfully concluded.

Refer to the Accu-Chek Insulin Pump Configuration Software user manual for further instructions on using the data transfer feature.



Press  twice to turn off and confirm the error.  
Restart the data transfer.

### 10.2.12 Error E13: LANGUAGE ERROR

When an error E13: LANGUAGE ERROR occurs check the language and reset if necessary.



Press  twice to turn off and confirm the error.

Refer to the ADVANCED SETUP MENU for instructions on how to change the language.

When finished, put your insulin pump into **RUN**, if necessary.


# 11 Troubleshooting

This section describes some potential problems that may occur while using your Accu-Chek Spirit insulin pump.

See sections 10 “Alerts and errors”, 12 “Life with your Accu-Chek Spirit insulin pump” and 14 “Care of your Accu-Chek Spirit insulin pump” for further information.

**Note** The possible solutions given are only suggestions. Do not use these procedures without approval of your doctor or healthcare team. Always follow your doctor’s or healthcare team’s recommendations, and contact them if you have any concerns with your insulin pump therapy.

Be aware that the interruption of insulin delivery (e.g. due to a leak, occlusion or loss of insulin potency) or the malfunction of your insulin pump can result in a rapid rise of your blood glucose level. Although your insulin pump has an internal security check system, it cannot alert you if your infusion set is leaking or if the insulin you are using has lost its potency.



## Problems with your Accu-Chek Spirit

<b>Problem</b>	<b>Recommended Actions</b>
Individual characters, numbers or symbols on the display are shown incompletely or not at all.	Contact your local Customer Care Area.
Your insulin pump does not beep or vibrate.	Check that the beep volume and the alarm signals are turned on. Check and change the battery. If your insulin pump still does not beep or vibrate, contact your local Customer Care Area.
Self-test incomplete (after an error E7: ELECTRONIC ERROR).	Remove the battery, insert it again after a few seconds. If the self-test remains incomplete, contact your local Customer Care Area.
Incorrect insulin delivery due to incorrect time setting.	The correct hourly basal rate cannot be delivered if the time is not set correctly. Check and reset the time as necessary. Do the same check when changing the battery.
Incorrect insulin delivery due to incorrect basal rate programming.	Always review any changes made in basal rate programming. Confirm that you are using your correct basal rate profile.
Incorrect insulin delivery due to incorrect insulin concentration used.	Your insulin pump has been developed exclusively for the continuous subcutaneous delivery of U100 short-acting insulin or U100 fast-acting insulin analogue. Do not use your insulin pump for the delivery of medications other than U100 short-acting insulin or U100 fast-acting insulin analogue.



<b>Problem</b>	<b>Recommended Actions</b>
Air bubbles in the cartridge and/or infusion set	<p>In the cartridge: Disconnect and re-prime your infusion set. Should the air bubbles remain, replace the cartridge.</p> <p>In your infusion set: Disconnect and prime your infusion set. Make sure all connections are tight and secure.</p>
Blood in infusion set tubing	Change your infusion set and infusion site.
Empty cartridge	Replace the cartridge. Cartridges are <b>for single use only</b> .
Disconnected or dislodged infusion set	<p>Test the connections for their security. Inspect connections for leaking insulin. Check your blood glucose level. Change your infusion set and infusion site.</p>
Leak in system	Check for leaking insulin at all connections and on the skin. Change your infusion set, infusion site, cartridge and adapter. Check your blood glucose level.
Occlusion in system	<p>Disconnect and re-prime your infusion set. If the occlusion remains, carry out the PRIME THE INFUSION SET function without the cartridge inserted. See the section “Error E4: OCCLUSION”.</p>
Failure to complete prime	<p>Disconnect and re-prime the infusion set. Always prime until a bubble-free insulin flow emerges from the tip of the needle. Make sure the tubing is free of air.</p>
Poor absorption of insulin	<p>Choose infusion sites without scar tissue, bruises, or tissue build-up.</p> <p>Do not use an infusion set or infusion site longer than recommended by your doctor or healthcare team.</p>

<b>Problem</b>	<b>Recommended Actions</b>
Infusion site becomes sore, red or swollen	Change your infusion set and infusion site immediately. Use proper insertion technique and always follow the planned infusion site rotation pattern and advice given by your doctor or healthcare team.

---

### **Therapy-related issues**

Many additional factors, such as alcohol consumption, medications other than insulin, ineffective or expired insulin, decreased activity, illness and stress can also potentially affect blood glucose levels. Contact your doctor or healthcare team for recommendations on how to manage these, or any other therapy related issues.

## 12 Life with your Accu-Chek Spirit insulin pump

### 12.1 Recommendations for daily use

#### Caution

Training and use of your insulin pump require the support of an experienced doctor or healthcare team. Regular visits with your doctor or healthcare team are absolutely essential during insulin pump therapy. Only change your personal settings after consulting your doctor or healthcare team. Always follow the instructions given by your doctor or healthcare team.

Successful insulin pump therapy requires frequent self-monitoring of blood glucose levels. It is recommended to check your blood glucose at least four times a day, or as directed by your doctor or healthcare team.

#### 12.1.1 Short interruption of insulin pump therapy

Consult your doctor or healthcare team about when and for how long insulin pump therapy can be interrupted.

1. Put your insulin pump into **STOP**.
2. Remove or disconnect your infusion set from your infusion site.  
When using a disconnectable infusion set, disconnect at your infusion site and use the protective covers provided for the infusion set.

Measure your blood glucose level regularly during any insulin interruption. Use a syringe or insulin pen to inject insulin according to the instructions of your doctor or healthcare team.

### 12.1.2 Continuation of insulin pump therapy

1. Attach a new Accu-Chek infusion set and prime it prior to insertion. When using a disconnectable infusion set, reconnect the tubing.
2. Follow the instructions for the infusion set you are using, and put your insulin pump into **RUN**.

Measure your blood glucose level within 2 hours in order to check the functionality of the system and to ensure insulin delivery.

### 12.1.3 Interruption of insulin pump therapy for a longer period of time

Before using a different Accu-Chek Spirit insulin pump, always check your personal settings in order to avoid incorrect insulin delivery. Make sure that the time and date are set correctly.

#### **Caution**

Contact your doctor or healthcare team for an alternate therapy plan when you interrupt the insulin pump therapy for a longer period of time.

There may be times when you may interrupt your insulin pump therapy for more than 1 day.

#### **To interrupt the use of your insulin pump**

1. Put your insulin pump into **STOP**.
2. Remove the cartridge, adapter and infusion set.
3. Reinsert the adapter.
4. Remove the battery and reinsert the battery cover.
5. Store your insulin pump properly. (See the section 14.3 “Storing your Accu-Chek Spirit insulin pump” for further information.)

## 12.2 Weather conditions

Wear your insulin pump under your clothes or directly on your body in cold and rainy weather. In case of doubt regarding these conditions, check the operating conditions of your insulin pump and contact your local Customer Care Area.

### **Caution**

Do not place your insulin pump in direct sunlight. Overheating of the insulin and your insulin pump must be avoided. Protect your insulin pump from direct exposure to cold wind. Temperatures over 40°C (104°F) and below 5°C (41°F) may damage the insulin and the electronics of your insulin pump and may cause a malfunction of the battery.

See the instructions for use of the insulin you are using for information on the admissible temperature range for the insulin.

See the section 15.1 “General technical data” for further information of operating conditions.

## 12.3 Your insulin pump and water

### Caution

Avoid deliberate contact with water. Check daily that your insulin pump is not chipped, cracked or damaged in any way and that the battery cover and the adapter are correctly closed. In the presence of chips and cracks, water, dust, insulin, or other foreign substances may enter your insulin pump and lead to malfunction. In case of deliberate contact with water disconnect and take off your insulin pump.

### 12.3.1 Daily situations

Your insulin pump is protected against unintentional water contact such as splashes or rain. For deliberate contact with water like bathing (e.g. tubs, jacuzzis or whirlpools), taking showers, swimming or any other water activities disconnect and take off your insulin pump. Avoid deliberate exposure with high humidity (e.g. saunas, etc.) as this may also harm your insulin pump. For your convenience we recommend to use disconnectable Accu-Chek infusion sets.

Consult your doctor or healthcare team for the length of time insulin pump therapy can safely be interrupted.

### 12.3.2 Accidental contact with water

You and your insulin pump do not have to worry about accidental temporary water contact. Examples of accidental temporary water contact include:

- rain, snow;
- water splashes while biking, jogging, hiking and similar activities;
- accidental immersions into the sink or bath tub.

#### Caution

You must inspect your insulin pump immediately after an accidental immersion in water.

### 12.3.3 What to do after water contact

Put your insulin pump into **STOP**. Disconnect your insulin pump prior to inspection. Use a soft cloth to dry the outside casing of your insulin pump. Check the battery compartment and the cartridge compartment for any water ingress. Should the battery compartment or cartridge compartment get wet, turn your insulin pump upside down to let the water run out and let it dry. Do not use warm air (e.g. a hair dryer) to dry as this could damage the internal electronics of your insulin pump. Do not insert the battery or the cartridge before these compartments are completely dry. See the section 14.2.1 “Cleaning your Accu-Chek Spirit insulin pump” for further information.

### 12.3.4 Other liquids

You and your insulin pump do not have to worry about:

- sweat
- saliva

You must check your insulin pump immediately after an accidental contact with other liquids such as:

- cleaning solutions
- alcohol
- beverages

See the section 14.2.1 “Cleaning your Accu-Chek Spirit insulin pump” for further information.

#### **Caution**

Avoid any contact of your insulin pump, the infusion set and the connecting parts of your insulin pump with health and beauty items (e.g. antiseptics, antibiotic creams, soaps, perfumes, deodorants, body lotions or any other cosmetics).



## 13 On the go with your Accu-Chek Spirit insulin pump

### 13.1 Electromagnetic fields and hazardous areas

Avoid electromagnetic fields of radar or antenna installations, high-voltage sources, X-Ray sources, MRI, CAT scan or other sources of electrical current. Do not use your insulin pump in such areas. Electromagnetic fields may cause malfunction to your insulin pump. Always stop and remove your insulin pump prior to entering these areas. In other cases, the insulin delivery may be stopped immediately and an error E7: ELECTRONIC ERROR occurs.

For further information regarding electromagnetic immunity, please refer to the section 15 “Technical data”.

#### **Warning**

Your insulin pump has not been tested with other electronic medical devices. Therefore do not use the Accu-Chek Spirit insulin pump with another electronic medical device unless advised by your doctor or healthcare team.

Your insulin pump has been tested and found to comply with the regulations concerning unintentional electromagnetic interference. Security systems for transit checks in airports as well as anti-theft monitoring devices (e.g. department stores, etc.) should not affect the functionality of your insulin pump.

Because of the high variety of devices with electromagnetic radiation such as cellular phones, their influence on your insulin pump and vice versa cannot be ruled out completely. It is advisable to keep your insulin pump at a distance of at least 10 cm (4 inches) from such devices while the other device is being operated.

Do not use your insulin pump in hyperbaric chambers and in hazardous areas of any classification (such as areas where explosive or flammable gases or vapours could exist) as this might interfere with insulin delivery and/or lead to harmful situations.

Your insulin pump is designed to work in normal barometric conditions from 70 to 106 kPa (700 to 1060 mbar). Do not exceed 3000 meters (10000 feet) above sea level.

Your pump has not been tested for use in hazardous areas of any classification. Always stop and remove your insulin pump prior to entering these areas.

If you have additional questions, contact your local Customer Care Area.

## 13.2 Sports

Exercise is an important aspect in the management of diabetes. Of course, you can perform sporting activities with the Accu-Chek Spirit insulin pump. However, prior to starting these types of activities be sure to protect your insulin pump. Do not wear your insulin pump during sports with body contact (e.g. boxing, football, hockey or rugby) as this type of activity could damage your insulin pump.

Consult our sterile products brochures and accessories catalogue for additional information regarding carrying systems or contact your local Customer Care Area for further information.

## 13.3 Travelling

Ask your doctor or healthcare team about any special precautions needed prior to travelling. Make sure you take extra insulin pump and blood-glucose testing supplies. It is highly recommended to make sure you know where you can obtain additional supplies while travelling.

When travelling across time zones make sure that the time and date of your insulin pump are set correctly. For further information please refer to the section 7.1 “Set time and date”.

## 14 Care of your Accu-Chek Spirit insulin pump

### Caution


Do not perform any servicing or repair on your insulin pump by yourself. For any questions contact your local Customer Care Area.

Roche Diagnostics as the manufacturer guarantees the functions of your insulin pump according to its specifications, on the condition that all servicing of your insulin pump is done and/or authorized by Roche Diagnostics.

### 14.1 System check

Only a well maintained system guarantees an accurate insulin delivery.

**Check the display of your insulin pump at least every three hours during the day, before you go to sleep, and especially if for any reason you might be unable to hear the beeps or to feel the vibrations.** Use only sterile products and accessories designed for use with your insulin pump. Replace and discard these items according to your doctor's or healthcare team's recommendations and according to the specifications in the corresponding instructions for use.



Check your insulin pump daily. Make sure that:

- The pump casing, the display and the cartridge are free of chips or cracks and the display is free of incomplete or abnormal letters or symbols.
- You visually inspect the cartridge. Make sure that the actual amount of insulin in the cartridge equals the displayed cartridge amount. This may be checked on the quick info screen.
- **You inspect all parts of your infusion set at least every three hours during the day and before you go to sleep.** Should you detect any loss of insulin immediately replace the leaking component.
- The battery cover is correctly tightened (even with the pump casing).
- The adapter is inserted properly and correctly tightened.
- Your infusion set is primed, free of air bubbles and tight in the adapter.
- Your infusion set is inserted according to the corresponding instructions for use.
- Your infusion site is secure, comfortable and free of irritation or infection.
- The cartridge is free of air bubbles.
- Your insulin pump is in **RUN**.
- The basal rates are set correctly according to your doctor's or healthcare team's recommendations.
- The time and date are correctly set.
- The temporary basal rate changes are set according to your doctor's or healthcare team's recommendations.
- The beeps and/or vibrations are set as desired.
- You have your personal emergency kit with you.

## 14.2 Maintenance and cleaning

### 14.2.1 Cleaning your Accu-Chek Spirit insulin pump

The ideal time to clean your insulin pump is during the change of a cartridge. Only use a soft, dry cloth for cleaning. If the insulin pump compartments are excessively dirty, consult your local Customer Care Area for further instructions.

#### Caution

Always remove the cartridge and adapter and place your insulin pump in **STOP** while cleaning. Avoid pressing the keys of your insulin pump during cleaning, as this may accidentally change your settings. Do not use alcohol, solvents, strong detergents, bleaching agents, scouring pads or sharp instruments for cleaning as they may damage your insulin pump.

### 14.2.2 Battery information

#### Caution

To prevent water from entering the pump casing, change the battery only in a dry environment and make sure that the seal of the battery cover is not worn out or missing and that the battery is properly inserted.

Always have an extra battery available and pay attention to the following:

- Use only 1.5 volt AA Alkaline batteries with a minimum capacity of 2500 mAh supplied and/or recommended by Roche Diagnostics or NiMH rechargeable batteries with a minimum capacity of 1500 mAh.

- **There are many types of AA batteries that are not designed to provide adequate power for your insulin pump. To ensure that the battery lasts as long as possible, be sure to use alkaline or NiMH (rechargeable) batteries in your Accu-Chek Spirit insulin pump.** Alkaline batteries supplied by Roche Diagnostics are the correct batteries for providing a maximum battery life.
- To operate your insulin pump, the operating temperature of the battery must be +5°C to +40°C (+41°F to +104°F).
- **Tighten or loosen the battery cover with the Accu-Chek Spirit battery key** (use of knives, screwdrivers or other sharp objects may harm your insulin pump). Do not over-tighten as this may cause damage to the battery cover and to the pump casing. The battery cover is correctly inserted and tightened when the battery cover is even with the pump casing.

The battery life is approximately four weeks for alkaline batteries and one week for rechargeable batteries if used in a typical usage pattern (50 U/day using U100 insulin; operating temperature 22°C ±3°C [72°F ±6°F]).

### 14.3 Storing your Accu-Chek Spirit insulin pump

#### Caution

When your insulin pump will not be in use for a long period of time, it should be properly stored, to prevent a malfunction.

To store your insulin pump:

- Remove the battery to preserve the battery lifetime
- Remove the cartridge
- Insert the battery cover and the adapter
- Store your insulin pump in its shipping case

## Storage conditions

Temperature	+5 to +45°C (+41 to +113°F)
Air humidity	5 to 85% relative air humidity
Barometric pressure	70 to 106 kPa (700 to 1060 mbar)

### 14.4 When your insulin pump is dropped

Dropping may harm your insulin pump and damage the seal against water. Make sure that your insulin pump is not dropped. Use the proper Accu-Chek Spirit carrying systems for the various daily situations to prevent dropping.

If your insulin pump is dropped:

- Check that all connections regarding your infusion set are still tight and reconnect, if necessary
- Check your insulin pump, its sterile products and accessories for chips and cracks  
AND
- Change the cartridge, if necessary.

#### Caution

Check your insulin pump and its sterile products and accessories for chips and cracks at least once per day, especially if they have been dropped. Do not use them if chipped or cracked. In the presence of chips and cracks, water, dust, insulin, or other foreign substances may enter your insulin pump and lead to malfunction.

In case of doubt, contact your local Customer Care Area.

## 14.5 Repair

### Caution

Do not perform any servicing or repair on your insulin pump by yourself. Do not use any lubricant on the pump mechanism. For additional questions contact your local Customer Care Area.

When experiencing problems with your insulin pump, contact your local Customer Care Area. The local Customer Care Area can address your concerns via the telephone or internet. Your insulin pump may need to be returned to Roche Diagnostics when approved by the local Customer Care Area. Return of your insulin pump would occur if an alert or error could not be cleared by following the procedures described in the section 10 “Alerts and errors”.

Pack your insulin pump with the used cartridge, battery, battery cover, adapter and infusion set safely to avoid damage during transportation. Place your insulin pump back in the shipping case and put the case in a package that will not rip or tear. It is best to send your insulin pump with a carrier that tracks the shipment.

Remember to include a note with

- a description of the reason for returning your insulin pump
  - your name and address
  - your daytime phone
  - the serial number of your insulin pump
- AND
- any RMA (Return Material Authorization) number received from the local Customer Care Area.

By returning your insulin pump to Roche Diagnostics, you allow Roche Diagnostics to undertake any testing (including destructive testing) on your insulin pump necessary for a proper inspection.



## 14.6 Disposal

### **Accu-Chek Spirit insulin pump**

If necessary, return your insulin pump to your local Customer Care Area for professional disposal.

### **Battery**

Dispose of depleted batteries at your nearest battery disposal center.

### **Cartridge, infusion set, adapter and additional accessories**

Dispose of these items in your normal garbage. To prevent injury to others, replace the needle protective cap or insert the needle of the infusion set into a disposable object, or put it into a hard shell container.

---

## 15 Technical data

### 15.1 General technical data

Maximum dimensions (without adapter)	Approx. 81 × 55 × 20 mm (3.2 × 2.2 × 0.8 inches)	
Weight	<ul style="list-style-type: none"> <li>• an empty insulin pump: approx. 80 g (2.8 ounces)</li> <li>• insulin pump inclusive battery, full plastic cartridge and infusion set: approx. 110 g (4.0 ounces)</li> </ul>	
Pump casing	Shock- and scratch-resistant plastic, resistant to pharmaceuticals, all edges rounded	
Temperature ranges*	During operation	+5 to +40°C (+41 to +104°F)
	Storage in its shipping case	+5 to +45°C (+41 to +113°F)
Air humidity	During operation	20 to 90% rel.
	Storage in its shipping case	5 to 85% rel.
Barometric pressure	During operation or storage in its shipping case	70 to 106 kPa (700 to 1060 mbar) (do not exceed 3000 meters [10000 feet] above sea level)

\* See the instructions for use of the insulin you are using for information on the acceptable temperature range for storage and handling.

Power supply	One 1.5 volt AA Alkaline battery or one rechargeable NiMH AA battery. Alkaline batteries should have a minimum capacity of 2500 mAh and NiMH AA rechargeable batteries a minimum capacity of 1500 mAh. Use only a battery charger officially recommended by the battery manufacturer.
Lifetime of battery	If used in a typical usage pattern (50 U/day using U100 insulin; operating temperature $22^{\circ}\text{C} \pm 3^{\circ}\text{C}$ [ $72^{\circ}\text{F} \pm 6^{\circ}\text{F}$ ]) the battery life is approximately four weeks for alkaline batteries and one week for rechargeable batteries.
Data storage time	The time and date is safely stored in the memory for about 1 hour after the battery has been removed. Your insulin pump's settings (the hourly basal rates, remaining cartridge content, bolus increments and active user menu) and the event memory (bolus history, history of daily insulin totals, temporary basal rate history, alarm history) are saved, regardless of battery voltage and time your insulin pump has been without a battery.
Delivery	1/20 of the current hourly basal rate in 3 minute intervals.
Basal rate	Min. 0.1 U/h, max. 25 U/h. There are 24 hourly basal rates adjustable in 0.1 unit increments.

---

Bolus	The max. bolus amount per delivery is 25 insulin units. The bolus amount for the “Quick” Standard bolus is adjustable in increments of 0.1, 0.2, 0.5, 1.0 and 2.0 units. For the “Scroll” Standard bolus, the Extended bolus and the MultiWave bolus, the amount is adjustable in fixed increments of 0.1 units. The duration of the Extended bolus and the MultiWave bolus is adjustable in intervals of 15 minutes (15 minutes up to 12 hours).
Temporary basal rate	Adjustable in 10% increments, 0–90% for decreases, 110–200% for increases. The duration is adjustable in 15 minute intervals, up to a maximum of 24 hours. The last programmed duration is given by default for the next temporary basal rate change programmed.
Maximum time for an error E4: OCCLUSION**	<p>Plastic cartridges:</p> <ul style="list-style-type: none"> <li>– at med. basal rate 1.0 U/h:     ≤ 5 h</li> <li>– at min. basal rate 0.1 U/h:     ≤ 50 h</li> </ul> <p>The typical time for an error was 3.5 hours.</p> <p>Aventis Insuman Infusat glass cartridges*:</p> <ul style="list-style-type: none"> <li>– at med. basal rate 1.0 U/h:     ≤ 10 h</li> <li>– at min. basal rate 0.1 U/h:     ≤ 100 h</li> </ul>

\* Insuman® Infusat produced by Aventis Pharma Deutschland GmbH, a member and a brand of the sanofi-aventis group.

\*\* Determined using measurement methods according to IEC 60601-2-24-1998.

## Technical data

---

---

Maximum volume before error E4: OCCLUSION**	Plastic cartridges: $\leq 3.5$ U The typical volume for an error was 2.3 U.
--	---

Aventis Insuman Infusat glass cartridges*:	$\leq 6.0$ U
---	--------------

---

Maximum pressure	400 kPa (4.0 bar), for plastic and Aventis Insuman Infusat glass cartridges*.
------------------	---

---

Flow (delivery rate)	During infusion set priming and bolus 0.2 U/sec.
----------------------	---

---

Maximum quantity delivered at a single fault condition	$\leq 1.0$ U
---	--------------

---

Cartridge	Accu-Chek 3.15 ml plastic cartridges with a luer connection and by Roche Diagnostics approved third party products. Aventis Insuman Infusat cartridges* have been evaluated and approved for use with the Accu-Chek Spirit insulin pump.
-----------	--

---

Infusion sets	Accu-Chek infusion sets with a luer-lock connector. Accu-Chek Infusion sets are PCV and latex free.
---------------	--

---

Data transfer	Infrared interface
---------------	--------------------

---

\* Insuman<sup>®</sup> Infusat produced by Aventis Pharma Deutschland GmbH, a member and a brand of the sanofi-aventis group.

\*\* Determined using measurement methods according to IEC 60601-2-24-1998.

---

Safety system	<p>Alert system, beeps, information on the display, vibrations, dual microprocessors.</p> <p>Your insulin pump is controlled by two microprocessors. The safety concept is based on one processor (supervisor processor) supervising the other (main processor). Whenever a defect or fault occurs in the main processor, it is identified by the supervisor processor. The motor is immediately switched off and an error E7: ELECTRONIC ERROR occurs. Conversely, the main processor can also determine at any time whether the supervisor is working correctly. The motor also constitutes an important safety component, as the combination of the main and supervisor processor and brushless motor provides the best possible reliability and accuracy in insulin delivery.</p>
IPX8	<p>Protected against the effects of temporary immersion in water under standardized conditions (up to 60 minutes and 2.5 meters [8 feet])</p>

---

## 15.2 Technical standards on electromagnetic emissions

The standard concerning electromagnetic compatibility of medical equipment (IEC 60601-1-2) requires the concrete specification of the corresponding levels which refer to specified electromagnetic interferences.

---

### Guidance and manufacturer's declaration – electromagnetic emissions

---

Emissions test RF emissions CISPR 11	Compliance Group 1	Electromagnetic environment - guidance The Accu-Chek Spirit uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Accu-Chek Spirit is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	–
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	–

---



## 15.3 Technical standards on electromagnetic immunity

### Guidance and manufacturer's declaration – electromagnetic immunity

The Accu-Chek Spirit is intended for use in the electromagnetic environment specified below. The customer or the user of the Accu-Chek Spirit should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 8 kV contact ± 15 kV air (IEC 60601-2-24)	Avoid any contact with synthetic materials. A relative humidity of 10% will allow ESD to be more likely.
Power frequency (50/60Hz) magnetic field  IEC 61000-4-8	3 A/m 400 A/m	400 A/m  (IEC 60601-2-24)	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

### Guidance and manufacturer's declaration – electromagnetic immunity

The Accu-Chek Spirit is intended for use in the electromagnetic environment specified below. The customer or the user of the Accu-Chek Spirit should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	Not applicable	

Radiated RF	3 V/m	10 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the Accu-Chek Spirit, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
IEC 61000-4-3	80 MHz to 2.5 GHz	80 MHz to 2.5 GHz (IEC 60601-2-24)	

Recommended separation distance:

$$d = 0.4 \cdot \sqrt{P} \quad \begin{array}{l} 80 \text{ MHz} \\ \text{to } 800 \text{ MHz} \end{array}$$

$$d = 0.7 \cdot \sqrt{P} \quad \begin{array}{l} 800 \text{ MHz} \\ \text{to } 2.5 \text{ GHz} \end{array}$$

where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and  $d$  is the recommended separation distance in meters (m).

Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey (a), should be less than the compliance level in each frequency range (b).

Interference may occur in the vicinity of equipment marked with the following symbol:



**Note 1** At 80 MHz and 800 MHz, the higher frequency range applies.

**Note 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM an FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To asses the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Accu-Chek Spirit is used exceeds the applicable RF compliance level above, the Accu-Chek Spirit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Accu-Chek Spirit.


b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

<b>V<sub>1</sub> in V:</b>	10.00	<b>E<sub>1</sub> in V/m:</b>	10.00
Rated maximum output power of transmitter W	150 kHz to 80 MHz out of ISM (Industrial, Scientific, Medical Band)	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	<b>Separation distance according to frequency of transmitter m</b>		
0.01	0.04	0.04	0.07
0.1	0.11	0.11	0.22
1	0.35	0.35	0.70
10	1.11	1.11	2.21
100	3.50	3.50	7.00

## Accuracy of flow rate for both Accu-Chek 3.15 ml plastic cartridges and Aventis glass cartridges\*\*\*

### Trumpet curve plotted from data after the end of the stabilization period

The trumpet curve shows the accuracy of the delivery rate in relation to the observation period.




The maximum deviation of the delivered amount (overall mean percentage flow error) for U100 insulin is  $\leq \pm 5\%^*$ .

\*\*\* Insuman® Infusat produced by Aventis Pharma Deutschland GmbH, a member and a brand of the sanofi-aventis group.

## Start-up graph over the stabilization period

The start-up graph shows changes in the flow rate over the stabilization time.



## Influence of height on delivery accuracy

The maximum deviation of the delivered amount (overall mean percentage flow error) for U100 insulin is  $\leq \pm 10\%$ \* in case your Accu-Chek Spirit insulin pump is located  $\leq 1$  m (3.3 feet) above/below the infusion site.

## Bolus

For U100 insulin the maximum deviation of a maximum bolus is  $\leq \pm 5\%$  and the maximum deviation of a minimum bolus is  $\leq \pm 30\%$ \*\*.

\* The measurements were made at a medium basal rate of 1.0 U/h according to IEC 60601-2-24:1998 with a Disetronic Classic infusion set PC 16/110 at room temperature.

\*\* The measurements were made according to IEC 60601-2-24:1998 with a Disetronic Classic infusion set PC 16/110 at room temperature.

## 15.4 Configuration parameters

Your insulin pump is provided by Roche Diagnostics with a standard configuration (factory settings). A list containing the important factory settings is added to the original insulin pump package when your insulin pump is first shipped. You may need to adapt the factory settings to meet your individual needs. Set the configuration parameters directly on your insulin pump and/or by using the Accu-Chek Insulin Pump Configuration Software. Always consult your doctor or healthcare team before changing parameters.

The Accu-Chek Spirit allows you to customize your user menus. If a menu does not appear, see the section 6 “Select a user menu (STANDARD, ADVANCED or CUSTOM)” and the Accu-Chek Insulin Pump Configuration Software user manual for further information on addressing all menus on your insulin pump.

The following list provides you with the full range of the adjustable configuration parameters on your insulin pump. Additionally, included are examples for typical range and parameter settings.

	Typical range available on the pump and standard parameter settings*	Maximum parameter range adjustable with the Accu-Chek Insulin Pump Configuration Software
<b>Basal rate</b>		
– hourly basal rate + max. TBR (hourly basal rate combined with temporary basal rate, max. increase)	20.0 IU/h	max. amount 0 – 62.50 IU/h
– hourly basal rates increment	0.1 IU	
– hourly basal rates	10.0 IU/h	max. amount 0 – 25.0 IU/h
– temporary basal rate TBR increase or decrease in increments of	0 – 200%	0 – 250%
– TBR duration	10% 15 min – 24 h	10% 15 min – 24 h
– TBR duration increments	15 min	15 min, 30 min, 1 h
– basal rate profiles available	1 – 5	1 – 5
– basal rate lock	on or off <b>presetting: off</b>	on or off
<b>Bolus</b>		
– max. bolus amount per bolus delivery	25.0 IU	max. amount 0 – 25.0 IU
– bolus amount increments for “Quick” Standard bolus	0.1, 0.2, 0.5, 1.0 or 2.0 IU <b>presetting: 0.5 IU</b>	0.1, 0.2, 0.5, 1.0 or 2.0 IU
– bolus duration (Extended bolus, MultiWave bolus) in increments of	15 min – 12 h 15 min	15 min – 24 h 15 min, 30 min, 1 h
Prime volume	0 – 30.0 IU <b>presetting: 25 IU</b>	max. amount 0 – 30.0 IU

\* Note: these factory settings may vary in some countries and depend on the changes you, your doctor or healthcare team made using your insulin pump or pump programming software for healthcare professionals.

	<b>Typical range available on the pump and standard parameter settings*</b>
Display orientation	standard or inverted <b>presetting: standard</b>
Display contrast	7 steps <b>presetting: medium</b>
User menus	3 (STANDARD, ADVANCED or CUSTOM) <b>presetting: ADVANCED</b>
Volume (for beeps)	off, 1 – 4 steps <b>presetting: medium</b>
Alarm signals turned on	beeps, vibrations or both <b>presetting: beeps and vibration</b>
Battery type	AA ALKALINE or AA NiMH rechargeable battery <b>presetting: AA ALKALINE battery</b>
KeyLock	turned on or off <b>presetting: off</b>
Automatic off	turned on (1 – 24 h) or off (0 h) <b>presetting: off</b>
Alarm clock	turned on or off when on: once or every day's occurrence <b>presetting: off</b>
Time format	European (24 h) or American (12 h, am/pm) <b>presetting: European</b>
Date format	European (dd.mm.yy) or American (mm/dd/yy) <b>presetting: European</b>

\* Note: these factory settings may vary in some countries and depend on the changes you, your doctor or healthcare team made using your insulin pump or pump programming software for healthcare professionals.



# 16 Annexes

## 16.1 Abbreviations

approx.	Approximately
BR	Basal rate
BRP	Basal rate profile
DDM	Diabetes Data Management
h	Hour(s)
IEC	International Electrotechnical Commission
incl.	Inclusive
IR	Infrared
I.U., IU	International Units in context with biological effectiveness of a certain insulin amount
kPa	Kilo pascals
LCD	Liquid Crystal Display
LED	Light Emitting Diode for infrared communication with e.g. a PC
NiMH	Nickel metal hydride (used with rechargeable batteries)
rel.	Relative
sec.	Second(s)
TBR	Temporary Basal Rate
U	See I.U.
U/h	Amount of International Units of insulin delivered per hour

---

U100	This is the insulin concentration. Each millilitre of liquid contains 100 International Units of insulin.
------	---

---

dd.mm.yy	European date format: day.month.year
----------	--------------------------------------





















---

mm/dd/yy	American date format: month/day/year
----------	--------------------------------------

---

## 16.2 Beeps and melodies

You can turn off the beeps or the vibrations, but not both at the same time. Please refer to sections 7.7.2 “Adjust the beep volume” and 7.7.3 “Alarm signals” for further information.




	Each press of 
	Each press of 
	Each press of  , respectively confirmation of each programmed bolus increment when programming a “Quick” Standard bolus
	Each press of  , respectively confirm the <b>STOP</b> Warning
	Temporary basal rate active
	Unlock the KeyLock in <b>RUN</b> , enter the <b>RUN</b> or <b>INFO</b> screen, start a temporary basal rate
	Lock the KeyLock in <b>RUN</b>
	Unlock the KeyLock in <b>STOP</b> , enter <b>STOP</b> screen.
	Lock the KeyLock in <b>STOP</b>
	Exit a screen
	Cancel a running bolus or stop priming
	Start-up procedure successfully concluded
	Copy an hourly basal rate
	Alert and error
	Maximum amount reached
	Minimum amount reached

## 16.3 Vibrations

Your insulin pump communicates alerts and errors using beeps and vibrations. You can turn off either the beeps or the vibrations, but not both at the same time. When you program a “Quick” Standard bolus, your Accu-Chek Spirit confirms all programming with vibrations. You cannot turn off these vibrations.

### Your Accu-Chek insulin pump vibrates when

---
















 (Short vibration)	Your Accu-Chek Spirit performs a start-up procedure
 (medium vibration)	Your Accu-Chek Spirit confirms the programmed bolus amount of the “Quick” Standard bolus.
 (Long vibration)	You <b>reset</b> the “Quick” Standard <b>bolus amount</b> to zero, You <b>cancel</b> the “Quick” Standard <b>bolus</b> , or When the <b>STOP</b> -Warning occurs.

---

See the sections 7.5.1 “Standard bolus” and 7.7.3 “Alarm signals” for further information.

## 16.4 Symbols

### 16.4.1 General symbols

		Read user guide
		Sterilized using ethylene oxide
		Sterilized using radiation
		Sterilized using antiseptics
		Year of manufacture
		Batch number
		Expiration date
<b>REF</b>		Item number
<b>SN</b>		Serial number
		Admissible temperature range during operation
		Admissible temperature range
		Admissible humidity range
		Admissible air pressure range
		Fragile – handle with care
		Protect against moisture
		Protect from heat and sunlight



For single use only

---



Does not contain pyrogen

---



Does not contain PVC

---



Recycling

---



Do not throw away

---



Do not use if the package is damaged

---



Flammable

---



See

---



Electronic device of type BF according the standard IEC 60601-1. Protection against electrical shock.

---

**IPX8**

Symbol for protection against the effect of temporary immersion of water (up to 60 minutes and 2.5 meters [8 feet]), according to IEC 60529.

---



Marking of conformity according to the European Medical Device Directive MDD 93/42/EEC with the number of the notified body.

---

**Rx only**

Federal law (USA) restricts this device to sale by or on the order of a physician.



















---























Manufacturer

---

## 16.4.2 Display symbols

	Start your Accu-Chek Spirit
	<b>STOP</b> screen and Stop insulin delivery
	Time, and Set the time and date menu
	Date
	Locked keys when the KeyLock function is turned on
	Unlocked keys when the KeyLock function is turned on
	Change the cartridge menu, and cartridge content
	Cartridge low warning
	Empty cartridge
	Low battery
	Empty battery
	Regular AA battery
	Rechargeable AA battery
U/h	Amount of International Units of insulin delivered per hour
	Basal rate profile
	Temporary basal rate
%	TBR percentage
	Daily basal rate total
	Basal rate programming unlocked
	Basal rate programming locked















---

	Standard bolus and Bolus increment programming
	Extended bolus and extended delivery of the MultiWave bolus
	MultiWave bolus
	Immediate bolus delivery of MultiWave bolus
	Decreased temporary basal rate (0 – 90%)
	Increased temporary basal rate (110 – 200%  )
	Remaining duration of an ongoing Extended bolus, MultiWave bolus, TBR or automatic off
	Daily total of insulin delivered (basal rate plus boluses)
	Prime the infusion set (menu)
	Alert or error occurring
	Alert
	Error
	Information menu
	Pump Timer
	Select the user menus STANDARD, ADVANCED or CUSTOM
	SETUP MENU STANDARD menu
	SETUP MENU ADVANCED menu
	Beeps turned on
	Vibrations turned on

---



---

	Beeps and vibrations turned on
	Set beep volume
	Screen orientation
<b>EU</b>	European date format
<b>US</b>	American date format
	Alarm clock turned on
	Alarm clock turned off
	Single alarm when set an alarm clock
	Every day repeat alarm when set an alarm clock
	PC communication
	Menu key
	Check key
	Up key
	Down key
	Timeout of a menu
	This setting may have been changed (by using your pump programming software).

---

## 16.5 Sterile products and accessories

### 16.5.1 Sterile products

#### Cartridges

Name	Comments
Accu-Chek 3.15 ml Plastic Cartridge	Plastic cartridges are single-use items only. Do not reuse cartridges. Refer to the instructions for use of the insulin you are using for information on the acceptable temperature range for storage and handling.
Aventis Insuman Infusat Pre-filled	
3.15 ml Glass Cartridges* (U100)	

#### Accu-Chek infusion sets

Name	Comments
Accu-Chek RapidLink	All Accu-Chek infusion sets are available in various tubing and needle lengths. Talk with your doctor and/or healthcare team to find the infusion set that fits you best.  Roche Diagnostics recommends that infusion sets are changed every two to three days, or at the recommendation of your doctor or healthcare team.
Accu-Chek TenderLink	
Accu-Chek FlexLink	
Accu-Chek Rapid-D Link	

\* Insuman® Infusat produced by Aventis Pharma Deutschland GmbH, a member and a brand of the sanofi-aventis group.

## 16.5.2 Accessories

Name	Comments
Adapter	Replace your adapter with every 10 <sup>th</sup> cartridge.
Batteries One 1.5 volt AA Alkaline battery with a minimum capacity of 2500 mAh.	If used in a typical usage pattern (50 U/day using U100 insulin; operating temperature 22°C ±3°C [72°F ±6°F]) the battery life is approximately four weeks.
Rechargeable Batteries AA NiMH rechargeable batteries with a minimum capacity of 1500 mAh.	If used in a typical usage pattern (50 U/day using U100 insulin; operating temperature 22°C ±3°C [72°F ±6°F]) the rechargeable battery life is approximately one week.  Rechargeable batteries are not provided by Roche Diagnostics. Use only a battery charger officially recommended by the battery manufacturer.
Battery Cover	Replace your battery cover with every 4 <sup>th</sup> battery.
Carrying Systems	Roche Diagnostics offers a wide range of carrying system solutions to best fit your lifestyle. If used in a normal usage pattern, the carrying systems have a life of approximately one year.
Accu-Chek Insulin Pump Configuration Software	Two versions of the pump configuration software are available; a consumer version and a healthcare professional version.
Travel Size Guide	This is a convenient pocket size reference that you may take with you to assist with your insulin pump questions.

---

## 17 Glossary

### **Accu-Chek Insulin Pump Configuration Software**

Roche Diagnostic's insulin pump programming tool. This software allows the fast and easy change of parameters and settings directly from your Microsoft Windows compatible PC.

### **Adapter**

The adapter physically connects the cartridge and infusion set. It has two seals and forms an efficient seal for the cartridge compartment of your insulin pump. The two small vents on the adapter allow air pressure to equalize.

### **Basal rate**

The amount of insulin delivered per hour that is needed to cover your basic insulin needs. In insulin pump therapy, your basal rate is determined together with your doctor or healthcare team and can be adjusted to meet your individual physiological needs throughout the day. Your basal rate is delivered by your insulin pump according to the curve of your personal basal rate profile or profiles.

### **Basal rate profile**

Your insulin pump offers the option of delivering up to 5 different basal rate profiles in order to easily meet changing insulin needs (for example during the week versus on the weekend). A basal rate profile consists of 24 programmed hourly basal rates.

### **Basal rate total**

The sum of all 24 hourly basal rates in one basal rate profile is called the (daily) basal rate total.

## **Bolus**

The amount of insulin delivered (in addition to the basal rate) to cover the intake of food and correct high blood glucose levels. The bolus amount is determined by your doctor or healthcare team's guidelines, your blood glucose level, your food intake and your activity level.

## **Carrying system**

A wide variety of carrying systems made of different materials are specially designed to suit your individual needs while wearing your insulin pump.

## **Cartridge**

The insulin reservoir of your insulin pump. It holds 3.15 ml (315 IU) of short acting insulin or fast acting insulin analogue.

## **Cartridge compartment**

The opening in your insulin pump for the cartridge.

## **Daily insulin total**

The total amount of insulin delivered (basal rate plus boluses) delivered in a 24 hour day, beginning at midnight. This amount does not include any insulin needed for the priming of infusion sets.

## **Diabetes Data Management (DDM)**

Diabetes Data Management is the recording of the therapy relevant data gained from your insulin delivery system and your blood glucose measurement system (such as an Accu-Chek blood glucose monitor) in order to analyse and illustrate this data on a PC or other communication devices.

## **Factory settings**

Your insulin pump is provided by Roche Diagnostics with a standard configuration (factory settings). These settings can be adapted to your individual needs directly on your insulin pump or by using the Accu-Chek Insulin Pump Configuration Software.

**Filling aid**

The filling aid can be used to help you fill an empty Accu-Chek 3.15 ml plastic cartridge from your vial of insulin.

**Free flow**

Free flow of insulin from the infusion may occur when two conditions are met. First, the cartridge with a connected infusion set and the piston rod of your insulin pump are not correctly connected and second, you position your insulin pump at a higher level than the infusion site.

**Hourly basal rate**

An hourly basal rate is the amount of insulin delivered by your insulin pump in 3 minute intervals for a given hour.

**Infusion set**

Infusion sets connect your insulin pump to your body. The insulin is delivered from the cartridge through the infusion set tubing and cannula or needle into your subcutaneous tissue.

Infusion sets can be disconnectable and non-disconnectable.

**Infusion site**

The place where the infusion set cannula or needle is inserted into your subcutaneous tissue to deliver insulin.

**Insulin**

A hormone that helps cells transform glucose into energy. Insulin is produced in the beta cells of the pancreas (also called islets of Langerhans).

**Insulin, fast-acting (insulin analogue)**

A type of insulin created using recombinant DNA technology. The onset time of insulin analogue is 5 to 15 minutes.

**Insulin, short-acting (regular insulin)**

A type of insulin created using recombinant DNA technology. The onset time of regular insulin is 30 to 45 minutes.

### **IPX8 according to IEC 60529**

Protection against the effects of temporary immersion in water. Ingress of water in quantities causing harmful effects is not possible when the device is temporarily immersed in water under standardized conditions (maximal 60 minutes a day and 2.5 meters [8 feet]).

### **Luer-lock connector**

A standardized fitting at the end of the infusion set, the cartridge and the adapter which allows them to be locked together without leaking.

### **Personal settings**

Your insulin pump must be programmed with your personal settings prior to starting insulin pump therapy. Personal settings include the basal rate profile(s), the correct time and date and all other variable values which can be tailored to your individual needs on your insulin pump or using Accu-Chek Insulin Pump Configuration Software.

### **RUN**

During normal use, your insulin pump is in **RUN** and insulin is continuously delivered. Boluses, temporary basal rates and almost all functions can be programmed from **RUN**.

### **RUN screen**

The **RUN** screen is the starting point for all functions which can be accessed while your insulin pump delivers insulin. Your insulin pump displays the **RUN** screen during normal use when no programming is currently being performed. The time, current hourly basal rate, basal rate profile and all functions currently activated appear on the **RUN** screen.

### **Settings**

Settings are individually programmable values and parameters that affect the way your insulin pump works.

### **Scrolling**

Scrolling allows a fast and easy programming of larger values.



**STOP**

When your insulin pump is in **STOP**, it does not deliver any insulin. The insulin delivery is only stopped if an error occurs or if your insulin pump is put into a necessary **STOP**, (e.g. when changing the cartridge, adapter or infusion set or for data transfer). Functions such as Extended bolus or temporary basal rate are interrupted by putting your insulin pump into **STOP**.

**STOP screen**

The **STOP** screen is the starting point for all functions that require that your insulin pump does not deliver insulin.

**Temporary basal rate**

Temporary increase or decrease of your basal rate profile in percentages (from 0 – 200%) to match changing insulin needs due to increased or decreased activity level, illness or stress.

**Timeout**

For your safety and convenience, your insulin pump automatically returns to the **RUN** or **STOP** screen if no key is pressed within a certain time limit. Any changes made are not saved.

**U100**

This is the insulin concentration. Each millilitre of liquid contains 100 International Units of insulin. Your insulin pump has been developed exclusively for delivery of U100 short-acting insulin or fast-acting insulin analogue.



## Index


- Accessories . . . . . 35
- Accu-Chek Insulin Pump
- Configuration Software . . . . 44, 78
- Adapter. . . . . 39
- Advanced setup menu. . . . . 134
- Advanced user menu . . . . 78, 121
- Alarm . . . . . 31, 110
- Alarm clock. . . . . 133, 149
  - turning off . . . . . 134, 149
- Alarm history . . . . . 116
- Alarm signals . . . . . 110
- Alcohol consumption. . . . . 168
- Alerts and errors . . . . . 145
  - confirming . . . . . 146
- Automatic off . . . . . 153
  
- Backlight . . . . . 26
- Backward loop . . . . . 77
- Basal rate . . . . . 87
  - and Bolus . . . . . 91
- Basal rate profile. . . . . 88
  - reviewing. . . . . 87
- Battery . . . . . 112, 179
  - battery depleted. . . . . 153
  - battery low . . . . . 148
  - type. . . . . 112, 179
- Beeps . . . . . 109
  - adjusting the volume . . . . 109
  - beeps and vibrations . . . . 145
  
- Blood glucose . . . . . 168
  - checking . . . . . 169
- Bolus . . . . . 95
  - cancel a “Quick” Standard . . 99
  - cancel a MultiWave . . . . . 129
  - cancel an Extended . . . . . 124
  - cancelled. . . . . 151
  - Extended . . . . . 95, 122
  - increment . . . . . 136
  - MultiWave . . . . . 95, 126
  - program a “Quick”  
Standard . . . . . 95, 97
  - program a MultiWave . . . . 127
  - program an Extended. . . . 123
- Bolus history. . . . . 114
  
- Cartridge. . . . . 37, 52
  - air bubbles in the. . . . . 57
  - cartridge empty . . . . . 152
  - cartridge low . . . . . 147
  - change the cartridge . . 71, 161
  - filling the . . . . . 54
  - inserting . . . . . 58, 60
  - plunger rod . . . . . 62
- Check key. . . . . 27
- Cleaning . . . . . 179
  - see also “Maintenance”
- Connector. . . . . 38
- Connect the cartridge, adapter  
and infusion set . . . . . 58
- Copy an hourly basal rate . . . . 90

- Daily insulin totals history . . . . . 116
- Daily totals history . . . . . 116
- Data interrupted . . . . . 163
- Data memory . . . . . 114
- Data transfer. . . . . 141
  - interruption of . . . . . 141
  - setting up . . . . . 142
- Date . . . . . 114, 116, 117
- Date and time setting; see “Time and date” . . . . . 83
- Date and time settings . . . . . 45
- Date format . . . . . 135
- Display . . . . . 25
  - contrast. . . . . 140
  - orientation . . . . . 14, 25, 113
- Disposal . . . . . 183
- Down key . . . . . 27
  
- Electromagnetic
  - emissions . . . . . 190
  - fields . . . . . 175
  - immunity . . . . . 191
- Electronic error . . . . . 159
- Emergency kit. . . . . 43
- Entry number . . . . . 114, 116, 117
- Errors . . . . . 152
- Exit options. . . . . 77
  
  
  
  
  
  
  
  
  
  
- Graphic LCD . . . . . 25
- Guarantee. . . . . 23
  
  
  
  
  
  
  
  
  
  
- Hygiene . . . . . 67
  
  
  
  
  
  
  
  
  
  
- Infection . . . . . 66
- Infusion set. . . . . 38
  - changing . . . . . 69
  - disconnectable. . . . . 38
  - instructions for use . . . . . 67
  - non-disconnectable . . . . . 38
  - not primed. . . . . 162
  - priming . . . . . 64
- Infusion site . . . . . 66
  - preparation . . . . . 66
  - selection . . . . . 66
- Insulin. . . . . 11, 116, 213
  - expired . . . . . 168
- Insulin delivery . . . . . 92, 93
  - start. . . . . 92
  - stop . . . . . 93
- Insulin pump. . . . . 25, 169, 175, 177
- Insulin pump therapy
  - continuation. . . . . 170
  - interruption . . . . . 169
  
  
  
  
  
  
  
  
  
  
- KeyLock; see also “Keys” . . . . . 29
  - making use of . . . . . 29
  - symbols . . . . . 29
  - turn on/off . . . . . 108
  - unlocking. . . . . 109
- Keys . . . . . 26
  - and key combinations . . . . . 26
  
  
  
  
  
  
  
  
  
  
- Language . . . . . 139
- Language error . . . . . 164
- Liquids . . . . . 174
- Looping . . . . . 77
- Luer-lock . . . . . 38, 42, 214

- 
- Maintenance . . . . . 179
- Malfunction . . . . . 37  
    reusing single use materials
- Mechanical error . . . . . 158
- Memory . . . . . 114
- Menu . . . . . 76  
    3 different . . . . . 78  
    navigating through . . . . . 76  
    navigation levels . . . . . 76  
    personalized . . . . . 79
- Menu key . . . . . 27
- Microsoft Windows; see “Software”
- MultiWave . . . . . 95, 126  
    cancelling . . . . . 129
- Occlusion . . . . . 154
- Orientation; see “Display” . . . . . 113
- Power interrupt . . . . . 160
- Prepare the pump for use . . . . . 45
- Prime quantity . . . . . 137
- Priming; see “Infusion set” . . . . . 64  
    stop . . . . . 66
- Products . . . . . 35
- Profile . . . . . 87, 131
- Programming . . . . . 88, 103  
    additional basal rates . . . . . 131
- Pump . . . . . 25  
    stopping the . . . . . 93  
    storing conditions . . . . . 181  
    storing . . . . . 180  
    wearing . . . . . 73
- Pump programming software . . . . . 44
- Pump therapy . . . . . 9, 11
- Pump timer . . . . . 114, 118, 149
- Quick info screen . . . . . 119
- Repair . . . . . 182
- Review time and date . . . . . 148
- RUN . . . . . 32
- Safety . . . . . 29
- Screen . . . . . 49
- Scrolling . . . . . 77
- Selecting profile . . . . . 132
- Self-test . . . . . 49
- Setup menu advanced . . . . . 134
- Setup menu standard . . . . . 107
- Software . . . . . 44
- Standard user menu . . . . . 83
- Start-up procedure . . . . . 48
- STATUS screen . . . . . 31
- Sterile product . . . . . 35
- STOP . . . . . 32, 93
- Stop-Warning . . . . . 33, 94  
    turning off . . . . . 94
- Supplies . . . . . 35
- System check . . . . . 177
-

TBR; see “Temporary basal rate” . . . . .	102
Temporary basal rate . . . . .	102
cancelled . . . . .	150
decrease . . . . .	102
duration . . . . .	104, 117
history . . . . .	117
increase . . . . .	102
interruption of . . . . .	106
over . . . . .	151
Time . . . . .	114, 116, 117, 118
Time and date	
american . . . . .	83
european . . . . .	83
review . . . . .	86
setting . . . . .	83
Time format . . . . .	134
Time remaining . . . . .	118
Troubleshooting . . . . .	165
Tubing . . . . .	39
Up key . . . . .	27
User menu	
advanced . . . . .	78, 121
custom . . . . .	78, 121
standard . . . . .	78, 83
Water . . . . .	172
accidental contact . . . . .	173
Wireless communication . . . . .	141

## Your Insulin pump STANDARD MENU




### Standard bolus

(A) **“Quick” Standard bolus** with direct access from the RUN screen using the  $\leftarrow$  and  $\rightarrow$  keys of your insulin pump.

(B) **“Scroll” Standard bolus** menu-guided using the STANDARD BOLUS menu and the scroll function of the  $\uparrow$  and  $\downarrow$  keys for programming the bolus amount.

For a complete review of your insulin pump features as well as the warnings and precautions associated with its use, please refer to your Accu-Chek Spirit Insulin pump user’s guide.

## Your Insulin pump ADVANCED MENU



---

## 19 Alert overview

In all cases, turn off and confirm the alert prior to addressing the issue.

No	Alert	What you have to do	Page
A1	CARTRIDGE LOW	Change cartridge before it is completely empty.	147
A2	BATTERY LOW	Replace the battery as soon as possible.	148
A3	REVIEW TIME AND DATE	Set the time and date.	148
A4	ALARM CLOCK	Put your Accu-Chek Spirit into <b>RUN</b> if necessary.	149
A5	PUMP TIMER	In your INFORMATION menu check for the insulin pump operation remaining under the menu PUMP TIMER. Arrange for replacement of pump within the remaining pump operation time (60 days).	149
A6	TBR CANCELLED	Put your Accu-Chek Spirit into <b>RUN</b> if necessary. Ensure the cancellation was intended and program a new temporary basal rate if necessary.	150
A7	TBR OVER	Decide whether a further temporary basal rate change is appropriate and program one if necessary.	151
A8	BOLUS CANCELLED	Put your Accu-Chek Spirit into <b>RUN</b> if necessary. Ensure that the cancellation was intended and program a new bolus if necessary.	151

No	Error	What you have to do	Page
E1	CARTRIDGE EMPTY	Change the cartridge.	152
E2	BATTERY DEPLETED	Change the battery.	153
E3	AUTOMATIC OFF	Put your Accu-Chek Spirit in <b>RUN</b> if necessary.	153
E4	OCCLUSION	See the appropriate section in this user guide for detailed information.	154
E5	END OF OPERATION	Arrange for immediate replacement of your insulin pump. Contact your healthcare professional for alternative insulin therapy options.	157
E6	MECHANICAL ERROR	See the appropriate section in this user guide for detailed information.	158
E7	ELECTRONIC ERROR	See the appropriate section in this user guide for detailed information.	159
E8	POWER INTERRUPT	See the appropriate section in this user guide for detailed information.	160
E10	CARTRIDGE ERROR	See the appropriate section in this user guide for detailed information.	161
E11	SET NOT PRIMED	Prime your infusion set.	162
E12	DATA INTERRUPTED	Restart the data transfer.	163
E13	LANGUAGE ERROR	See the appropriate section in this user guide for detailed information.	164

---





ACCU-CHEK, ACCU-CHEK SPIRIT,  
ACCU-CHEK TENDERLINK, ACCU-CHEK FLEXLINK,  
ACCU-CHEK RAPIDLINK and DISETRONIC  
are trademarks of Roche.



Disetronic Medical Systems AG  
3401 Burgdorf, Switzerland  
[www.accu-chek.com](http://www.accu-chek.com)

CE 0123

**Rx only**  
Swiss made

04540344001/V06/10.08