

# IntelliVue MP2 Patient Monitor

## Philips M8102A Technical data sheet

The IntelliVue MP2 portable patient monitor is compact in size, ergonomic, and modular in design. It provides an easy-to-use touchscreen user interface, is highly customizable and shares a technological platform with the Philips IntelliVue MP5-MP90 patient monitors.

The IntelliVue series offers a complete monitoring solution that is flexible and modular, designed to suit a broad spectrum of monitoring needs.

### Measurement Features

- Compact, rugged, lightweight monitor with built in measurements.
- ECG monitoring using any combination of three to 10 electrodes.

- 12-lead ECG monitoring with five electrodes using the EASI method or with 10 electrodes using the conventional method.
- Multi-lead arrhythmia and ST segment analysis at the bedside on all available leads.
- Mainstream or Sidestream CO<sub>2</sub>
- Choice of FAST SpO<sub>2</sub> or Nellcor™ OxiMax™ SpO<sub>2</sub><sup>1</sup>
- Invasive Pressure and Temperature measurement
- The monitor can operate using battery power for up to 3 hours with basic monitoring configuration to let you safely and easily monitor patients during

<sup>1</sup> The following are trademarks of Covidien AG and/or its affiliates: Nellcor™, Durasensor™, Dura-Y™, Oxiband™, OxiCliq™, OxiMax™, MAXFAST™.

# PHILIPS

in-hospital transfer. AC power is provided by an external power supply.

- Telemetry devices can be connected via short range radio to monitor telemetry data (ECG/SpO<sub>2</sub>) on the MP2 screen (Telemetry as a parameter (TAAP)).
- IntelliVue Cableless Measurement Devices can be connected via short range radio to monitor data from the IntelliVue CL SpO<sub>2</sub> Pod or IntelliVue CL NBP Pod on the MP2 screen. The Cableless Measurement Devices can also be controlled from an assigned MP2 via short range radio.

### Usability Features

- Touchscreen and hardkeys as input device.
- Intuitive user interface.
- Simple menu hierarchy gives fast access to all basic monitoring tasks.
- Patient data management with tabular and graphic trends.
- Settings “Profiles” for rapid case turnover.
- Patented automatic alarm limits help clinicians provide care more efficiently.
- 3.5" TFT flat panel display with QVGA (320 x 240) resolution, wide viewing angle, large numerics, permanently visible alarm limits, and up to three real-time waves.
- Capable of functioning in a wireless infrastructure (IIT)

### Intended Use

The monitor is intended to be used for monitoring and recording of, and to generate alarms for, multiple physiological parameters of adults, pediatrics, and neonates in a hospital environment and during patient transport inside and outside of hospitals. The MP2 when used with the TRx4841A/TRx4851A IntelliVue Telemetry System Transceiver is intended for use in a hospital environment and during patient transport inside the hospital environment.

The monitor is intended for use by health care professionals.

The monitor is only for use on one patient at a time. It is not intended for home use. Not a therapeutic device.

Rx only: U.S. Federal Law restricts this device to sale by or on the order of a physician.

ST segment monitoring is intended for use with adult patients only and is not clinically validated for use with neonatal and pediatric patients.

The ECG measurement is intended to be used for diagnostic recording of rhythm and detailed morphology of complex cardiac complexes (according to AAMI EC 11).

### Hospital Environment:

The monitor is suitable for use in all medically used rooms which fulfill the requirements regarding electrical installation according to IEC60364-7-710 “Requirements for special installations or locations - Medical locations”, or corresponding local regulations.

### Upgradability

The MP2 monitor allows new capabilities to be added in the future as your monitoring requirements evolve. This upgradability gives the security of knowing that the monitors can be enhanced and updated as practices and technologies advance, and it protects long-term investments.

### Main Components

#### Monitor

The monitor has a color LCD TFT display with a wide viewing angle, providing high resolution waveform and data presentation.

The display, processing unit and measurements are integrated into one device. An external power supply provides power for the monitor.

#### User Interface

The user interface is designed for fast and intuitive operation. The color graphical user interface ensures that clinicians quickly feel at ease using the monitor.

Configurable SmartKeys with intuitive icons allow monitoring tasks to be performed quickly and easily, directly on the monitor screen.

Waves and numerics are color-coded.

The monitor displays up to three measurement waves simultaneously. For 12-lead ECG monitoring it can display 12 real-time ECG waves, with a rhythm strip and all ST values.

Flexible screen layout allows optimal use of the available display space, for example, waves can be overlapped or wave size can adjust dynamically depending on the number of waves configured for the space.

The MP2 monitor is supplied with a resistive touchscreen.

#### Simulated Keyboard

If alpha or numeric data entry is required, for example to enter patient demographics, an on-screen keyboard will automatically appear on the screen.

#### Mounting

The mounting options available enable flexible, space saving placement of the monitors for an ergonomic work space. The monitor is shipped with a low cost mounting plate if not specified otherwise.

### Application Features

#### Critical and Cardiac Care Features

- The monitor performs multi-lead **arrhythmia detection** analysis on the patient's ECG waveform at the bedside. It analyzes for ventricular arrhythmias, calculates heart rate, and generates alarms, including asystole, bradycardia, and ventricular fibrillation.

- Up to 12 leads of **ST segment analysis** can be performed on adult patients at the bedside, measuring ST segment elevation and depression and generating alarms and events. The user can trend ST changes, set high and low alarm limits, and set both ST and isoelectric measurement points. Using ST Snippets, one-second wave segments can be compared with a baseline segment for each measured ST lead.
- optional **ST Map** application shows ST changes over time in two multi-axis spider diagrams.
- **QT/QTc interval monitoring** provides the measured QT interval, the calculated heart-rate corrected QTc value and a  $\Delta$ QTc value, which tracks variation in the QT interval in relation to a baseline value.
- optional **12-lead ECG** data can be measured, using either the EASI placement method with five standard electrodes or conventional electrode placement with 10 electrodes.<sup>1</sup>  
12 realtime ECG waveforms can be displayed simultaneously. Diagnostic 12-lead ECG can be captured, reviewed and stored on the patient monitor before it is sent to the Information Center. Local printout is available, in harmonized layout.
- High performance pulse oximetry technologies perform accurately even in cases with low perfusion.
- Choice of sidestream or mainstream **CO<sub>2</sub> monitoring** for high quality measurements with intubated and non-intubated patients.
- **Telemetry devices** (TRx4841A/TRx4851A TRx/TRx+ IntelliVue Transceiver) can be connected via short range radio to the MP2 to monitor telemetry data (ECG/SpO<sub>2</sub>) on the MP2 screen.

#### Ease of Use

- **Screen layouts** are easily adjustable, allowing flexible display of measurement information.
- Temperature, height, and weight can be configured either in metric or imperial **units**. Pressure measurements can be displayed in kPa or mmHg. Gases can be displayed in kPa, mmHg.

#### Trends

- The **trend database** stores patient data from up to 16 measurement numerics. The measurement information can be sampled every 12 seconds, one minute, or five minutes, and stored for a period ranging from four to 48 hours.
  - Each NBP measurement generates a column in the Vital Signs trend table. The values for the other measurements are added to provide a complete vital signs set for the NBP measurement time.
- **Horizon Trends** show the deviation from a stored baseline.

#### Transport Features

- The monitor's portable design means it can be used for in and out-of-hospital transport: a basic monitor weighs 1.5 kg.
- The monitor can operate using battery power for up to 3 hours, to let you safely and easily monitor patients during procedures or in-hospital transfer.
- Specially-designed mounting solutions let you quickly disconnect the monitor for transport and reconnect to the mount after transport.
- The Universal Admit, Discharge and Transfer (ADT) feature means that all ADT information is shared between the networked monitor and the Information Center. Information need only be entered once.

#### Patient Data Documentation

- An extensive range of **Patient Reports** can be printed:
  - 12-lead ECG Reports
  - Alarm Limit Reports
  - Vital Signs
  - Graphic Trends
  - Realtime Wave Reports
 Report templates can be defined in advance, enabling print-outs tailored to each hospital's specific requirements to be started quickly. Reports can be printed on centrally-connected printers or via the IntelliVue PC Printing Solution, and they can be initiated manually or automatically at user-defined intervals.
- The IntelliVue PC Printing Solution allows printing of reports, waveform captures and trends from the MP2 to a standard off-the-shelf printer or to an electronic file.

#### Alarms

The alarm system can be configured to present either the traditional HP/Agilent/Philips alarm sounds or sounds compliant with the IEC 60601-1-8 Standard.

Alarm limits are permanently visible on the main screen. The Alarm Limits page provides a graphic depiction of alarm limits in relation to the currently monitored measurement values and lets you adjust alarm limits. It also lets you preview wide and narrow automatic alarm limits before you apply them.

When an alarm limit is exceeded, it is signalled by the monitor in the following ways:

- an alarm tone sounds, graded according to severity
- an alarm message is shown on the screen, color-coded according to severity
- the numeric of the alarming measurement flashes on the screen
- alarm lamps flash for red and yellow alarms and are illuminated for technical INOPs

A "SmartAlarm Delay" algorithm helps to reduce the number of pulse oximetry nuisance alarms.<sup>2</sup>

<sup>1</sup> EASI-derived 12-lead ECGs and their measurements are approximations to conventional 12-lead ECGs. As the 12-lead ECG derived with EASI is not exactly identical to the 12-lead conventional ECG obtained from an electrocardiograph, it should not be used for diagnostic purposes.

If the monitor is connected via a network to a central monitoring station, alarming is simultaneous at the monitor and at the Information Center.

Alarms are graded and prioritized according to severity:

- **Red Alarms\*\*\*** identify a potentially life threatening situation for a patient.
- **Yellow Alarms\*\*** indicate conditions violating preset vital signs limits.
- **Technical Alarms (INOPS)** are triggered by signal quality problems, equipment malfunction or equipment disconnect. The Silence/Pause Alarms function (equivalent to Silence/Suspend with previous monitor generations) allows the user to switch off alarm tones with one touch.

All alarms can be paused for a period of one, two, three, five, or 10 minutes or turned off indefinitely.

Alarm strip recordings are available on a centrally-connected recorder or via the IntelliVue PC Printing Solution.

Patented automatic alarm limits automatically adapt the alarm limits to the patient's currently measured vital signs within a safe margin defined individually for each patient.

Visual and/or audible latching and non-latching alarm handling is available.

## Profiles

Profiles are predefined configuration settings for Screens, measurement settings, and monitor settings. Each Profile can be designed for a specific application area and patient category, for example OR adult, or ICU neonatal. Profiles enable a quick reaction to patient and care location changes: activating a Profile with a particular patient category (Adult, Pediatric or Neonatal) automatically applies suitable alarm and safety limits and saves time usually spent carrying out a complete set-up procedure.

Profiles can be created directly on the monitor or remotely on a personal computer and transferred to the monitor using the IntelliVue Support Tool. A selection of Profiles for common monitoring situations is provided with the monitor. These profiles can be changed, added to, renamed, or deleted.

## Optional Networking Capabilities

The monitor can operate as part of a wired or wireless hospital network system, using the Philips IntelliVue Clinical Network interface. This includes:

- DHCP protocol support (as an alternative to BootP in certain network designs)
- 802.1x basic support on wireless networks

- WMM on wireless networks
- QoS Tagging

## Service Features

- The Support Tool helps technical personnel to
  - carry out configuration, upgrades and troubleshooting via the network, or on an individual monitor
  - share configuration settings between monitors
  - back up the monitor settings.
- A password-protected Service Mode ensures that only trained staff can access service tests and tasks.
- The Configuration Mode is password-protected and allows trained users to customize the monitor configuration.

## Device Connections

The monitor can be connected to:

- an Information Center (for example M3150B)
- a PC
- MMS Extensions (M3012A, M3014A, M3015A/B, )<sup>1</sup>

## Network Interface

The network interface provides the system with networking capability via a wired or wireless network connection.

## Wireless Network

The monitor can function within a wireless infrastructure based on an IEEE 802.11 a/g network in the 2.4 GHz / 5 GHz bands (ISM). Additionally, the monitor can function within a telemetry infrastructure compatible with the Philips Cellular Telemetry System (CTS) in the WMTS and ISM bands. Additional components are required to complete the system. Please refer to the M3185A IntelliVue Clinical Network Technical Data Sheet for further information.

A short range radio interface for an IEEE 802.15.4 network in the 2.4 GHz (ISM) band is also available. This allows a telemetry device with a short range radio adapter or IntelliVue Cableless Measurement Devices to be assigned to the monitor.

# Monitor Specifications

## Safety Specifications

The monitor complies with the Medical Device Directive 93/42/EEC (CE0366) and with IEC 60601-1:1988 + A1:1991 + A2:1995; EN60601-1:1990 + A1:1993 + A2:1995; UL 60601-1:2003; CAN/CSA

<sup>2</sup> The "SmartAlarm Delay" is not available in the U.S.A. and territories relying on FDA Market clearance.

The Smart Alarm Delay functionality is currently not available in China or in clinical environments under SFDA control.

<sup>1</sup> The MMS Extensions will only function when they are connected to the Philips Battery Extension, or the monitor is connected to external power.

C22.2#601.1-M90 + Suppl. No 1-94 + Am.2; JIS T 0601-1:1999; IEC 60601-1-1:2000; EN 60601-1-1:2001.

All applied parts are Type CF unless otherwise specified. They are protected against damage from defibrillation and electrosurgery. The possibility of hazards arising from software errors was minimized in compliance with ISO 14971:2000, EN60601-1-4:1996 + A1:1999 and IEC 60601-1-4:1996 + A1:1999.

The monitor complies with the EMC standards IEC 60601-1-2:2001; EN 60601-1-2:2001

This ISM device complies with Canadian ICES-001. Cet appareil ISM est conforme à la norme NMB-001 du Canada.

The MP2 patient monitor with measurements and interfaces other than those listed below cannot be used for patient transport outside of the hospital environment.

The MP2 patient monitor with the following measurements and interfaces:

- ECG/Respiration, NBP, SpO<sub>2</sub>, Pressure, Temperature, CO<sub>2</sub><sup>1</sup>
- LAN, Battery

can be used in a transport environment such as road ambulance, airplane or helicopter. For this purpose, the monitor fulfills the following additional mechanical, EMC and environmental requirements:

- **Shock Tests** according to IEC TR 60721-4-7, Class 7M3. Test procedure according to IEC/EN 60068-2-27 (peak acceleration up to 100g).
- **Random Vibration** according to IEC TR 60721-4-7, Class 7M3. Test procedure according to IEC/EN 60068-2-64 (RMS acceleration 5g).
- **Sinusoidal Vibration** according to IEC TR 60721-4-7, Class 7M3. Test procedure according to IEC/EN 60068-2-6 (acceleration up to amplitude 2g).
- **Bump Test** according to IEC/EN60068-2-29 (peak acceleration 15g, 1000 bumps).
- **Free Fall Test** according to EN1789 (covers also IEC TR 60721-4-7 and Class 7M3). Test procedure according to EN 60068-2-32 (height 0.75 m).
- Specification for degrees of protection provided by enclosures according to **IEC/EN 60529: IP 32**
- **EN 1789 +A1:2003** Medical vehicles and their equipment - Road ambulances (chapter 6 - Medical Devices).
- **EN13718-1** Air, water and difficult terrain ambulances. Medical devices interface requirements for the continuity of patient care. For Ambulances, Patient transport equipment, Emergency vehicles, Ambulance services, Rough-terrain vehicles, Water transport, Air transport, Medical equipment, Medical instruments, Interfaces, Performance.

- **Radiated susceptibility 20 V/m** according to EN ISO 9919 (SpO<sub>2</sub>) and EN ISO 21647 (CO<sub>2</sub>).
  - **Altitude Range** from -500 to 3000 m operating and -500 to 4600 m storage and transportation.
  - Extended radiated susceptibility tests
- The MP2 patient monitor with its out-of-hospital parameter set provides a general immunity level of 20 V/m with only few restrictions. Details are as listed below:
- GSM 900: Immunity at 900 MHz (uplink mobile phone), 20 V/m, duty cycle 1:8
  - GSM 1800: Immunity at 1800 MHz (uplink mobile phone), 20 V/m, duty cycle 1:8.
  - DECT: Immunity at 1800 MHz (digital cordless phone), 20 V/m, duty cycle 1:24
  - AM: 1 kHz Immunity from 80 MHz to 2.5 GHz (any radio communication unit, broadcasting and TV transmitter), 20 V/m, modulation factor 80%. (ECG: 20 V/m except 0.8-1.2 GHz where it is 10 V/m)
- **Operating ambient temperature** testing over the range from 0°C to 40°C (32°F to 104°F).
  - **Operating ambient humidity** testing up to 95% RH at 40°C (104°F), non condensing.

#### US Army Airworthiness Certification Granted

- U.S. Army Airworthiness Certification and Evaluation (ACE) program of U.S. Army Aeromedical Research Laboratory (USAARL)
- Tests performed in accordance with the following standards:
  - MIL-STD-461E Electromagnetic interference characteristics requirements and limits.
  - MIL-STD-810F Department of Defense test method standard for environmental engineering considerations and laboratory tests.
  - MIL-STD-1472F Human Engineering.
  - ANSI/AAMI HE48-1993 HF Engineering guidelines & preferred practices for the design of medical devices.
  - ANSI/AAMI ES1-1993 Safe current limits for electromedical apparatus.

<sup>1</sup> The [MP2 Patient Monitor](#) in combination with the M2741A LoFlo CO<sub>2</sub> Sensor meets CISPR 11, Group 1, Class A emission limits. Thus it is suitable only for use in all establishments other than domestic and those connected to a low power network which supplies buildings used for domestic purposes.

## Physical Specifications

a sufficient for flight altitudes up to 12,000 m with pressurized cabins

| Product                                     | Max Weight | W x H x D            |
|---|------------|----------------------|
| M8102A                                      | <1.25 kg   | < 188 x 99 x 86 mm   |
| IntelliVue MP2 (without handle and options) | (2.8 lb)   | (7.4 x 3.9 x 3.4 in) |

## Environmental Specifications

| Item   | Condition                      | Range  |
|--|--------------------------------|--|
| <b>Temperature Range</b>                           | Operating                      | 0°C to 40°C<br>(32°F to 104°F)   |
|  | Storage<br>(incl. Transport)   | -20°C to 60°C<br>(-4°F to 140°F)   |
| <b>Temperature Range when charging the battery</b> | Operating                      | 0°C to 35°C<br>(32°F to 95°F)  |
| <b>Humidity Range</b>                              | Operating                      | 15% to 95% Relative Humidity (RH) (non condensing)   |
|  | Storage and Transport          | 5% to 95% Relative Humidity (RH)   |
| <b>Altitude Range</b>                              | Operating                      | -500 m to 3000 m<br>(10000 ft)   |
|  | Storage and Transport          | -500 m to 4600 m<br>(15000 ft) <sup>a</sup>  |
| <b>Ingress Protection</b>                          | Monitor                        | IP32 (protected against the ingress of solid foreign objects 2.5 mm in diameter or larger, and the ingress of water when the water is dripping vertically and the monitor is tilted up to 15°).  |
|  | External Power Supply (M8023A) | IP31(protected against the ingress of solid foreign objects 2.5 mm in diameter or larger, and the ingress of water when the water is dripping vertically) when rested on its rubber feet on a flat, level surface.<br><br>IP32 when mounted with the connectors facing downwards |

## Performance Specifications

| Monitor Performance Specifications |  |  |
|------------------------------------|--|--|
| <b>Power Specifications</b>        | Power consumption  | <40 W average, <65 W peak  |
|                                    | Line Voltage   | 100 to 240 V ~   |
|                                    | Current  | 1.3 to 0.7 A   |
|                                    | Frequency  | 50/60 Hz   |
| <b>Battery Specifications</b>      | Operating Time (with new, fully charged battery at 25°C)   | Basic monitor configuration: 3 hours   |
|                                    | Charge Time  | When MP2 is off: 2 h<br>When MP2 is in use and connected to the external power: 12 h approx. |
| <b>Indicators</b>                  | Alarms Off   | red LED  |
|                                    | Alarms   | red/yellow/cyan LED  |
|                                    | On/Standby/Error   | green/red LED  |
|                                    | AC Power   | green LED  |
|                                    | Battery  | yellow (charging)/red blinking (empty) LED   |
| <b>Sounds</b>                      | External Power   | green LED  |
|                                    | Audible feedback for user input. Prompt tone. QRS tones, or SpO <sub>2</sub> modulation tone. Four different alarm sounds. |  |

### Trends:

12, or 16 numerics @ 12 sec, 1 minute, 5 minute resolution.  
Multiple choices of number of numerics, resolution and duration depending on trend option and application area.

|                        |   |   |
|------------------------|---|---|
| <b>Alarm Signal</b>    | System delay  | less than 3 seconds                                   |
|                        | Pause duration  | 1,2,3 minutes or infinite, depending on configuration |
|                        | Extended alarm pause  | 5 or 10 minutes                                       |
| <b>Review Alarms</b>   | Information: all alarms / inops, main alarms on/off, alarms acknowledged and time of occurrence   |   |
|                        | capacity  | 500 items   |
| <b>Real Time Clock</b> | Range: from: January 1, 1997, 00:00 to: December 31, 2080, 23:59  |   |
|                        | Accuracy: < 4 seconds per day (typically)<br>Hold Time: infinite if powered by host monitor or external power supply; otherwise at least 48 hours |   |

## Monitor Performance Specifications

|                        |  |
|------------------------|--|
| <b>Buffered Memory</b> | Contents: Active settings, trends, patient data, realtime reports, review alarms     |
|                        | Hold Time: infinite if powered by external power supply; otherwise at least 48 hours |

Restart time: After power interruption, an ECG wave will be shown on the display after 30 seconds maximum.

## M8023A External Power Supply Performance Specifications

| M8023A External Power Supply Performance Specifications |              |                |
|---|--------------|----------------|
| <b>Power Specifications</b>                             | Power        | < 12 W average |
|   | Consumption  | < 30 W peak    |
|   | Line Voltage | 100 to 240 V ~ |
|   | Current      | 0.7 to 0.4 A   |
| <b>Indicators</b>                                       | Frequency    | 50/60 Hz ~     |
|   | AC Power     | green LED      |
|   |              |                |

## Interface Specifications

| MP2 (M8102A) Interface Specifications    |  |   |
|--|--|---|
| <b>Measurement Link (MSL)</b>            | Connectors                               | Female ODU (Proprietary)  |
|  | Power                                    | 30 V to 60 V input  |
|  | Power Sync                               | RS-422 compliant input<br>78.125 kHz (typical)                                  |
|  | LAN signals                              | IEEE 802.3 10-Base-T complaint  |
|  | Serial signals                           | RS-422 compliant  |
|  | Local signals                            | Provided for connecting MMS extensions  |
| <b>ECG Sync Pulse Output<sup>a</sup></b> | Cable Detection                          | Yes   |
|  | Marker In                                | No  |
|  | Wave Output                              | No  |
|  | Connector                                | Binder Series 709/719   |
|  | Output Levels                            | Output low <0.8V @ I = -4 mA<br>Output high >2.4 V @ I = 4 mA                   |
|  | Isolation                                | None  |
|  | Pulse Width                              | 100 +/- 10 ms (high)  |
|  | Delay from R-wave peak to start of pulse | 20 ms maximum per AAMI EC13   |
|  | Minimum required R-wave amplitude        | 0.5 V   |
|  |  |   |
| <b>802.11 Bedside Adapter</b>            | Wireless Technology                      | IEEE 802.11 a/b/g   |
|  | Frequency Band                           | 2.4 GHz and 5 GHz ISM   |
|  | Technology                               | compatible with Philips Cellular Telemetry System (CTS) cellular infrastructure |
| <b>Internal WMTS Adapter (US only)</b>   | Frequency Band                           | WMTS, 1395-1400 MHz and 1427-1432 MHz   |
|  | Technology                               | compatible with Philips Cellular Telemetry System (CTS) cellular infrastructure |
| <b>Internal ISM Adapter</b>              | Frequency Band                           | 2.4 GHz ISM   |
|  | Technology                               | compatible with Philips Cellular Telemetry System (CTS) cellular infrastructure |

## MP2 (M8102A) Interface Specifications

|  |                          |                             |
|--|--------------------------|-----------------------------|
| <b>Short Range Radio Interface<sup>b</sup></b> | Type                     | Dual internal SRR Interface |
|  | Technology               | IEEE 802.15.4               |
|  | Frequency                | 2.4 GHz ISM                 |
|  | Band                     | (2,400 - 2,483 GHz)         |
|  | Modulation               | DSSS (O-QPSK)               |
|  | Technique                |                             |
|  | Effective Radiated Power | max. 0 dBm (1 mW)           |

<sup>a</sup> ECG Sync Pulse not available if ECG is sourced from the telemetry device

<sup>b</sup> The short range radio interface is compatible with the IntelliVue Cableless Measurements and the following telemetry devices: TRx4841A/TRx4851A IntelliVue Telemetry System Transceiver.

## M8023A External Power Supply Interface Specifications

|                               |                |  |
|-------------------------------|----------------|--|
| <b>Measurement Link (MSL)</b> | Connectors     | Male ODU (Proprietary)                       |
|                               | Power          | 48 V output                                  |
|                               | Power Sync.    | RS-422 compliant output 78.125 kHz (typical) |
|                               | LAN signals    | IEEE 802.3 10-Base-T compliant               |
|                               | Serial signals | RS-422 compliant output                      |
|                               | Local signals  | Not connected                                |

## Display Specifications

|                                |                   |                             |
|--------------------------------|-------------------|-----------------------------|
| <b>Integrated QVGA Display</b> | Sweep Speeds      | 6.25, 12.5, 25 and 50 mm/s; |
|                                | Resolution        | 320 x 240                   |
|                                | Refresh frequency | 60 Hz                       |
|                                | Useful screen     | 72 x 54 mm (2.8 x 2.1 in)   |
|                                | Pixel size        | 0.22 x 0.22 mm              |



| MP2 (M8102A) Compatible Devices                               |                |  |
|---|----------------|--|
| IntelliVue Instrument Telemetry Wireless Network (USA only)   |                |  |
| <b>Internal WMTS Adapter</b>                                  | Technology     | compatible with Philips Cellular Telemetry System (CTS), cellular infrastructure |
|   | Frequency Band | WMTS, 1395-1400 MHz and 1427-1432 MHz  |
| IntelliVue Instrument Telemetry Wireless Network (except USA) |                |  |
| <b>Internal ISM Adapter</b>                                   | Technology     | compatible with Philips Cellular Telemetry System (CTS), cellular infrastructure |
|   | Frequency Band | 2.4 GHz ISM  |

| M4607A Battery Specifications   |  |
|---------------------------------|--|
| Physical Specifications         |  |
| W x D x H                       | 66 mm (2.36 in) x 80 mm (3.15 in) x 20 mm (0.79 in)  |
| Weight                          | 160 g ±5%  |
| Performance Specifications      |  |
| Nominal Voltage                 | 10.8 Volt  |
| Rated Capacity at discharge C/5 | 1000 mAh (typical)   |
| Environmental Specifications    |  |
| Temperature Range               | Discharge 0°C to 60°C (32°F to 140°F)<br>Charge 0°C to 60°C (32°F to 140°F)<br>Storage and Transportation: -20°C to 65°C (-4°F to 149°F) |
| Humidity Range                  | Operating: 15% to 95% Relative Humidity (RH)<br>Storage and Transportation: 5% to 95% Relative Humidity (RH)                             |
| Battery Type                    | Lithium Ion Mangan, 10.8 V, 1000 mAh,  |
| Safety                          | complies with UL 2054 (UL Recognized)  |
| Communication Standard          | complies with the SMBus specification v1.1   |

| M4605A Battery Specifications   |  |
|---------------------------------|--|
| Physical Specifications         |  |
| W x D x H                       | 149 mm (5.866 in) x 89 mm (3.504 in) x 19.8 mm (0.78 in)   |
| Weight                          | 490 g (1.08 lb)  |
| Performance Specifications      |  |
| Nominal Voltage                 | 10.8 Volt  |
| Rated Capacity at discharge C/5 | 6000 mAh (typical)   |
| Continuous Discharge Capability | 6.5 A  |
| Environmental Specifications    |  |
| Temperature Range               | Discharge 0°C to 50°C (32°F to 122°F)<br>Charge 0°C to 50°C (32°F to 122°F)<br>Storage and Transportation: -20°C to 65°C (-4°F to 149°F) |
| Humidity Range                  | Operating: 15% to 95% Relative Humidity (RH)<br>Storage and Transportation: 5% to 95% Relative Humidity (RH)                             |
| Battery Type                    | Smart Battery 10.8 V, 6000 mAh, Lithium Ion  |
| Safety                          | complies with UL 2054 (UL Recognized)  |
| Communication Standard          | complies with the SMBus specification v1.1   |

## Measurement Specifications

### ECG/Arrhythmia/ST/QT

Complies with IEC 60601-2-25:1993 + A1:1999 /EN60601-2-25:1995 + A1:1999, IEC 60601-2-27:2005/EN60601-2-27:2006, IEC 60601-2-51:2003 /EN 60601-2-51:2003 and AAMI EC11/EC13:1991/2002.

| ECG/Arrhythmia/ST Performance Specifications |                                |   |
|--|--------------------------------|---|
| <b>Cardiotach</b>                            | Range                          | Adult/pedi:<br>15 to 300 bpm<br>Neo range:<br>15 to 350 bpm       |
|  | Accuracy                       | ±1% of range  |
|  | Resolution                     | 1 bpm   |
|  | Sensitivity                    | ≥200 $\mu V_{peak}$   |
| <b>PVC Rate</b>                              | Range                          | 0 to 300 bpm  |
|  | Resolution                     | 1 bpm   |
| <b>ST Numeric</b>                            | Range                          | -20 to +20 mm   |
|  | Accuracy                       | ±0.5 mm or 15%,<br>whichever is greater                           |
|  | Resolution                     | 0.1 mm  |
| <b>QT Numeric</b>                            | Range                          | 200 to 800 ms   |
|  | Accuracy                       | ±30 ms  |
|  | Resolution                     | 8 ms  |
| <b>QTc Numeric</b>                           | Range                          | 200 to 800 ms   |
|  | Resolution                     | 1 ms  |
| <b>ΔQTc Numeric</b>                          | Range                          | -600 to +600 ms   |
|  | Resolution                     | 1 ms  |
| <b>Q-T-HR Numeric</b>                        | Range - adult                  | 15 to 300 bpm   |
|  | Range - pediatric and neonatal | 15 to 350 bpm   |
| <b>Sinus and SV Rhythm Ranges</b>            | Brady                          | Adult: 15 to 60 bpm<br>Pedi: 15 to 80 bpm<br>Neo: 15 to 90 bpm    |
|  | Normal                         | Adult: 60 to 100 bpm<br>Pedi: 80 to 160 bpm<br>Neo: 90 to 180 bpm |
|  | Tachy                          | Adult: >100 bpm<br>Pedi: >160 bpm<br>Neo: >180 bpm                |

| ECG/Arrhythmia/ST Performance Specifications   |                          |  |
|--|--------------------------|--|
| <b>Bandwidth</b>   | Diagnostic Mode          | Adult/neo/pedi: 0.05 to 150 Hz   |
|  | Extended Monitoring Mode | Neo/pedi: 0.5 to 150 Hz  |
|  | Monitoring Mode          | Adult: 0.5 to 40 Hz<br>Neo/pedi: 0.5 to 55 Hz  |
|  | Filter Mode              | Adult/neo/pedi: 0.5 to 20 Hz   |
| <b>Bandwidth when ECG is transmitted from a telemetry device via short range radio</b> | Diagnostic Mode          | Adult/neo/pedi: 0.05 to 40 Hz  |
|  | Extended Monitoring Mode | Neo/pedi: 0.5 to 40 Hz   |
|  | Monitoring Mode          | Adult: 0.5 to 40 Hz<br>Neo/pedi: 0.5 to 40 Hz  |
|  | Filter Mode              | Adult/neo/pedi: 0.5 to 20 Hz   |
| <b>Differential Input Impedance</b>  |                          | >2 M $\Omega$ RA-LL leads (Resp)<br>>5 M $\Omega$ at all other leads (at 10 Hz including patient cable)                              |
| <b>Common Mode Rejection Ratio</b>   |                          | Diagnostic mode:<br>>86 dB (with a 51 k $\Omega$ /47 nF imbalance).<br>Filter mode: >106 dB (with a 51 k $\Omega$ /47 nF imbalance). |
| <b>Electrode Offset Potential Tolerance</b>  |                          | ±500 mV  |
| <b>Auxiliary Current (Leads off Detection)</b>   |                          | Active electrode:<br><100 nA<br>Reference electrode:<br><900 nA  |
| <b>Input Signal Range</b>  |                          | ±5 mV  |

| ECG/<br>Arrhythmia/<br>ST Alarm<br>Specifications | Range  | Adjustment  |
|---|--|---|
| <b>HR</b>   | 15 to 300 bpm<br>maximum delay:<br>10 seconds<br>according to AAMI<br>EC 13-1992<br>standard | Adult:1 bpm steps (15 to<br>40 bpm)<br>5 bpm steps (40 to<br>300 bpm)<br>Pedi/Neo:1 bpm steps<br>(15 to 50 bpm)<br>5 bpm steps (50 to<br>300 bpm) |
| <b>Extreme<br/>Tachy</b>                          | Difference to high<br>limit 0 to 50 bpm<br>Clamping at 150 to<br>300 bpm                     | 5 bpm steps<br>5 bpm steps  |
| <b>Extreme<br/>Brady</b>                          | Difference to low<br>limit 0 to 50 bpm<br>Clamping at 15 to<br>100 bpm                       | 5 bpm steps<br>5 bpm steps  |
| <b>Run PVCs</b>                                   | 2 PVCs   | Not adjustable by user  |
| <b>PVCs Rate</b>                                  | 1 to 99 PVCs/<br>minute  | 1 PVC   |
| <b>Vent Tach HR</b>                               | 20 to 300 bpm  | 5 bpm   |
| <b>Vent Tach<br/>Run</b>                          | 3 to 99 PVCs/<br>minute  | 1 PVC   |
| <b>Vent Rhythm<br/>Run</b>                        | 2 to 99 PVCs/<br>minute  | 1 PVC   |
| <b>SVT HR</b>                                     | 120 to 300 bpm   | 5 bpm   |
| <b>SVT Run</b>                                    | 3 to 99 SV beats   | 1 SV beat   |
| <b>ST High</b>                                    | -19.8 to +20 mm  | 0.2 mm  |
| <b>ST Low</b>                                     | -20 to +19.8 mm  | 0.2 mm  |
| <b>QTc High</b>                                   | 200 ms to 800 ms   | 10 ms steps   |
| <b>ΔQTc High</b>                                  | 30 ms to 200 ms  | 10 ms steps   |

| ECG/Arrhythmia/ST Supplemental Information as<br>required by AAMI EC11/13 |  |
|---|--|
| <b>Respiration Excitation<br/>Waveform</b>                                | Sinusoidal signal, 260 $\mu$ A, 40.5 kHz           |
| <b>Noise Suppression</b>  | RL drive gain 44 dB max., max.<br>voltage 1.8 Vrms |

| ECG/Arrhythmia/ST Supplemental Information as<br>required by AAMI EC11/13 |                      |   |
|---|----------------------|---|
| <b>Time to<br/>Alarm for<br/>Tachy-<br/>cardia</b>                        | Vent                 | Gain 0.5, Range 6.5 to 8.4 seconds,<br>Average 7.2 seconds  |
|   | Tachycardia          |   |
|   | 1 mV <sub>pp</sub> , | Gain 1.0 Range 6.1 to 6.9 seconds,<br>Average 6.5 seconds   |
|   | 206 bpm              |   |
|   |                      | Gain 2.0, Range 5.9 to 6.7 seconds,<br>Average 6.3 seconds  |
|   |                      |   |
|   | Vent                 | Gain 0.5, Range 5.4 to 6.2 seconds,<br>Average 5.8 seconds  |
|   | Tachycardia          |   |
|   | 2 mV <sub>pp</sub> , | Gain 1.0, Range 5.7 to 6.5 seconds,<br>Average 6.1 seconds  |
|   | 195 bpm              |   |
|   |                      | Gain 2.0, Range 5.3 to 6.1 seconds,<br>Average 5.7 seconds  |
|   |                      |   |
| <b>Tall T-Wave Rejection<br/>Capability</b>                               |                      | Exceeds ANSI/AAMI EC 13 Sect.<br>3.1.2.1(c)<br>minimum recommended 1.2 mV T-<br>Wave amplitude  |
| <b>Heart Rate Averaging<br/>Method</b>                                    |                      | Three different methods are used:<br>Normally, heart rate is computed<br>by averaging the 12 most recent<br>RR intervals.<br>For runs of PVCs, up to 8 RR<br>intervals are averaged to compute<br>the HR.<br>If each of 3 consecutive RR<br>intervals is greater than 1200 ms<br>(that is, rate less than 50 bpm),<br>then the 4 most recent RR<br>intervals are averaged to compute<br>the HR. |
| <b>Response Time of Heart<br/>Rate Meter to Change in<br/>Heart Rate</b>  |                      | HR change from 80 to 120 bpm:<br>Range: [6.4 to 7.2 seconds]<br>Average: 6.8 seconds<br>HR change from 80 to 40 bpm:<br>Range: [5.6 to 6.4 sec] Average:<br>6.0 seconds   |
| <b>Heart Rate Meter<br/>Accuracy and Response<br/>to Irregular Rhythm</b> |                      | Ventricular bigeminy: 80 bpm<br>Slow alternating ventricular<br>bigeminy: 60 bpm<br>Rapid alternating ventricular<br>bigeminy: 120 bpm<br>Bidirectional systoles: 90 bpm  |
| <b>Accuracy of Input Signal<br/>Reproduction</b>                          |                      | Methods A and D were used to<br>establish overall system error and<br>frequency response.   |

## Respiration

| Respiration Performance Specifications |            |   |
|--|------------|---|
| <b>Respiration Rate</b>                | Range      | Adult/pedi: 0 to 120 rpm<br>Neo: 0 to 170 rpm                   |
|  | Accuracy   | at 0 to 120 rpm $\pm 1$ rpm<br>at 120 to 170 rpm<br>$\pm 2$ rpm |
|  | Resolution | 1 rpm   |
| <b>Bandwidth</b>                       |            | 0.3 to 2.5 Hz (-6 dB)   |
| <b>Noise</b>                           |            | Less than 25 m $\Omega$ (rms)<br>referred to the input          |

| Respiration Alarm Specifications | Range   | Adjustment  | Delay   |
|----------------------------------|---|---|---|
| <b>High</b>                      | Adult/pedi: 10 to 100 rpm<br>Neo: 30 to 150 rpm | under 20 rpm: 1 rpm steps<br>over 20 rpm: 5 rpm steps | max. 14 seconds   |
| <b>Low</b>                       | Adult/pedi: 0 to 95 rpm<br>Neo: 0 to 145 rpm    | under 20 rpm: 1 rpm steps<br>over 20 rpm: 5 rpm steps | for limits from 0 to 20 rpm: max. 4 seconds<br>for limits above 20 rpm: max. 14 seconds |
| <b>Apnea Alarm</b>               | 10 to 40 seconds                                | 5 second steps  |   |

### Philips FAST SpO<sub>2</sub>

Complies with EN ISO 9919:2005 (except alarm system; alarm system complies with IEC 60601-2-49:2001).

Measurement Validation: The SpO<sub>2</sub> accuracy has been validated in human studies against arterial blood sample reference measured with a CO-oximeter. Pulse oximeter measurements are statistically distributed, only about two-thirds of the measurements can be expected to fall within the specified accuracy compared to CO-oximeter measurements. Display Update Period: Typical: 2 seconds, Maximum: 30 seconds. Max. with NBP INOP suppression on:

60 seconds. For SpO<sub>2</sub> specifications of a connected telemetry device, see the specifications document of the telemetry device.

| SpO <sub>2</sub> Performance Specifications |          |   |
|---|----------|---|
| <b>SpO<sub>2</sub>*</b>                     | Range    | 0 to 100%   |
|   | Accuracy | <b>Philips Reusable Sensors:</b><br>M1191A, M1191AL, M1191B, M1191BL, M1192A: 2% (70% to 100%)<br>M1193A, M1194A, M1195A, M1196A: 3% (70% to 100%)<br><b>Philips Reusable Sensors with M1943A(L):</b><br>M1191T, M1192T, M1193T (Adult), M1196T: 3% (70% to 100%)<br>M1193T (Neonate): 4% (70% to 100%)<br><b>Philips Disposable Sensors with M1943A(L):</b><br>M1132A, M1133A, M1134A (adult/infant): 2%<br>M1131A, M1133A, M1134A (neonate), M1901B, M1902B, M1903B, M1904B: 3% (70% to 100%)<br><b>Nellcor™ Sensors with M1943A(L):</b><br>MAXA, MAXAL, MAXP, MAXI, MAXN, D-25, D-20, I-20, N-25, OxiCliq A, P, I, N: 3% (70% to 100%) |

| SpO <sub>2</sub> Performance Specifications |            |   |
|---|------------|---|
| SpO <sub>2</sub> *                          | Accuracy   | <b>Masimo Reusable Sensors® with LNOP MP12 or LNC MP10:</b><br>LNOP DCI, LNOP DCIP, LNOP YI, LNCS DCI, LNCS DCIP:<br>2% (70% to 100%)<br>LNOP TC-I, LNCS TC-I:<br>3.5% (70% to 100%)<br><b>Masimo Disposable Sensors® with LNOP MP12 or LNC MP10:</b><br>LNOP Adt, LNOP Adtx, LNOP Pdt, LNOP Pdtx, LNOP Inf-L, LNCS Adtx, LNCS Pdtx, LNCS Inf-L: 2% (70% to 100%)<br>LNOP Neo-L, LNOP NeoPt-L, LNCS Neo-L, LNCS NeoPt-L: 3% (70% to 100%) |
|   | Resolution | 1%  |
|   | Pulse      |   |
|   | Range      | 30 to 300 bpm   |
| Pulse                                       | Accuracy   | ±2% or 1 bpm, whichever is greater  |
|   | Resolution | 1 bpm   |
| Sensors                                     |            | Wavelength range: 500 to 1000 nm<br>Emitted Light Energy: ≤15 mW<br>Information about the wavelength range can be especially useful to clinicians (for instance, when photodynamic therapy is performed)  |
| Pulse Oximeter Calibration Range            |            | 70 - 100%   |

\*The specified accuracy is the root-mean-square (RMS) difference between the measured values and the reference values

| SpO <sub>2</sub> Alarm Specifications | Range  | Adjustment  | Delay                           |
|---------------------------------------|--|---|---------------------------------|
| SpO <sub>2</sub>                      | Adult: 50% to 100%<br>Pedi/Neo: 30 to 100%                         | 1% steps  | (0, 1, 2, 3,... 30) + 4 seconds |
| Desat                                 | Adult: 50% to Low alarm limit<br>Pedi/Neo: 30% to Low alarm limit  | 1% steps  |                                 |
| Pulse                                 | 30 to 300 bpm  | Adult:<br>1 bpm steps (30 to 40 bpm)<br>5 bpm steps (40 to 300 bpm)<br>Pedi/Neo:<br>1 bpm steps (30 to 50 bpm)<br>5 bpm steps (50 to 300 bpm) | max. 14 seconds                 |
| Tachycardia                           | Difference to high limit 0 to 50 bpm<br>Clamping at 150 to 300 bpm | 5 bpm steps<br>5 bpm steps  | max. 14 seconds                 |
| Bradycardia                           | Difference to low limit 0 to 50 bpm<br>Clamping at 30 to 100 bpm   | 5 bpm steps<br>5 bpm steps  | max. 14 seconds                 |

**Nellcor OxiMax SpO<sub>2</sub> Specifications  
(M8102A #SP4)**

| Pulse Oximetry Performance Specifications |  |
|---|--|
| <b>SpO<sub>2</sub></b>                    |  |
| <b>Measurement Range</b>                  | 1 to 100%  |
| <b>Resolution</b>                         | 1%   |
| <b>Accuracy</b>                           | see table below  |
| <b>Low Perfusion Accuracy<sup>a</sup></b> | 2% (70 - 100%)   |
| <b>Pulse</b>                              |  |
| <b>Measurement Range</b>                  | 25 to 300 bpm  |
| <b>Resolution</b>                         | 1 bpm  |
| <b>Accuracy</b>                           | +/- 3 bpm (25 to 250 bpm)  |
| <b>Low Perfusion Accuracy<sup>a</sup></b> | +/- 3 bpm (25 to 250 bpm)  |
| <b>Sensors</b>                            |  |
| <b>Sensors</b>                            | Wavelength range: 500 to 1000 nm<br>Emitted Light Energy: ≤15 mW<br>Information about the wavelength range can be especially useful to clinicians (for instance, when photodynamic therapy is performed) |
| <b>Numeric Update Rate</b>                |  |
| <b>Numeric Update Rate</b>                | typical 1 second, max ≤ 60 sec   |

a Specification applies to Monitor performance. Reading accuracy in the presence of low perfusion (detected IR pulse modulation amplitude 0.03% - 1.5%) was validated using signals supplied by a patient simulator. SpO<sub>2</sub> and pulse rate values were varied across the monitoring range over a range of weak signal conditions and compared to the known true saturation and pulse rate of the input signals.

The specified accuracy is the root-mean-square (RMS) difference between the measured values and the reference values.

| SpO <sub>2</sub> Accuracy Table |                                    |      |                                       |
|---------------------------------|------------------------------------|------|---------------------------------------|
|                                 | SaO <sub>2</sub> Range:<br>70-100% |      | SaO <sub>2</sub><br>Range:<br>60%-80% |
| <b>MAXFAST</b>                  | 2%                                 | NA   | 3%                                    |
| <b>MAXR<sup>b</sup></b>         | 3.5%                               | NA   | NA                                    |
| <b>SC-A</b>                     | 2%                                 | NA   | NA                                    |
| <b>SC-PR-I<sup>c</sup></b>      | NA                                 | 2%   | NA                                    |
| <b>SCNEO-I<sup>c</sup></b>      | NA                                 | 2%   | NA                                    |
| <b>OxiCliq A</b>                | 2.5%                               | NA   | NA                                    |
| <b>OxiCliq P</b>                | 2.5%                               | NA   | NA                                    |
| <b>OxiCliq N<sup>d</sup></b>    | 2.5%                               | 3.5% | NA                                    |
| <b>OxiCliq I</b>                | 2.5%                               | NA   | NA                                    |
| <b>D-YS<sup>d</sup></b>         | 3%                                 | 4%   | NA                                    |
| <b>D-YS &amp; D-YSE</b>         | 3.5%                               | NA   | NA                                    |
| <b>D-YSPD</b>                   | 3.5%                               | NA   | NA                                    |
| <b>DS-100A-1</b>                | 3%                                 | NA   | NA                                    |
| <b>OXI-A/N<sup>d</sup></b>      | 3%                                 | 4%   | NA                                    |
| <b>OXI-P/I</b>                  | 3%                                 | NA   | NA                                    |

a M1901B/MAXN:

Clinical functionality has been demonstrated on a population of hospitalized neonate patients. The observed SpO<sub>2</sub> accuracy was 2.5% in a study of 42 patients with ages of 1 to 23 days, weight from 750 to 4100 grams, and 63 observations made spanning a range of 85 to 99% SaO<sub>2</sub> while monitored with Nellcor OxiMax N-595 pulse oximeters.

b The accuracy specification has been determined between saturations of 80%–100%.

c SoftCare SC-PR-I, SCNEO-I:

Clinical functionality has been demonstrated on a population of hospitalized neonate and infant patients. The observed SpO<sub>2</sub> accuracy was 3.0% in a study of 57 patients with ages of 24 to 40 weeks, weight from 710 to 5,000 grams, and 185 observations made spanning a range of 63 to 100% SaO<sub>2</sub> while monitored with Nellcor OxiMax N-595 pulse oximeters.

d Neonatal accuracy: When sensors are used on neonatal subjects as recommended, the specified accuracy range is increased by ± 1 digit, as compared to adult usage, to account for the theoretical effect on oximeter measurements of fetal hemoglobin in neonatal blood. For example, OxiCliq N accuracy on neonates is ± 3.5 digits, rather than ± 2.5.

| SpO <sub>2</sub> Accuracy Table |                                    |                |                                       |
|---------------------------------|------------------------------------|----------------|---------------------------------------|
|                                 | SaO <sub>2</sub> Range:<br>70-100% |                | SaO <sub>2</sub><br>Range:<br>60%-80% |
| <b>Sensor</b>                   | <b>Adult/Infant</b>                | <b>Neonate</b> | <b>Adult</b>                          |
| <b>M1901B<sup>a</sup></b>       | Identical to OxiMax MAXN           |                |                                       |
| <b>M1902B</b>                   | Identical to OxiMax MAXI           |                |                                       |
| <b>M1903B</b>                   | Identical to OxiMax MAXP           |                |                                       |
| <b>M1904B</b>                   | Identical to OxiMax MAXA           |                |                                       |
| <b>MAXA, MAXAL</b>              | 2%                                 | NA             | 3%                                    |
| <b>MAXN<sup>a</sup></b>         | 2%                                 | 2%             | 3%                                    |
| <b>MAXP</b>                     | 2%                                 | NA             | 3%                                    |
| <b>MAXI</b>                     | 2%                                 | NA             | 3%                                    |

## NBP

Complies with IEC 60601-2-30:1999/EN60601-2-30:2000.

| NBP Performance Specifications         |            |  |
|--|------------|--|
| <b>Measure-<br/>ment Ranges</b>        | Systolic   | Adult: 30 to 270 mmHg (4 to 36 kPa)<br>Pedi: 30 to 180 mmHg (4 to 24 kPa)<br>Neo: 30 to 130 mmHg (4 to 17 kPa)   |
|  | Diastolic  | Adult: 10 to 245 mmHg (1.5 to 32 kPa)<br>Pedi: 10 to 150 mmHg (1.5 to 20 kPa)<br>Neo: 10 to 100 mmHg (1.5 to 13 kPa)   |
|  | Mean       | Adult: 20 to 255 mmHg (2.5 to 34 kPa)<br>Pedi: 20 to 160 mmHg (2.5 to 21 kPa)<br>Neo: 20 to 120 mmHg (2.5 to 16 kPa)   |
|  | Pulse Rate | Adult: 40 to 300<br>Pedi: 40 to 300<br>Neo: 40 to 300  |
| <b>Accuracy</b>                        |            | Max. Std. Deviation: 8 mmHg (1.1 kPa)<br>Max. Mean Error: $\pm 5$ mmHg ( $\pm 0.7$ kPa)  |
| <b>Pulse Rate Measurement Accuracy</b> |            | 40 to 100 bpm: $\pm 5$ bpm<br>101 to 200 bpm: $\pm 5\%$ of reading<br>201 to 300 bpm: $\pm 10\%$ of reading (average over NBP measurement cycle)                               |
| <b>Heart Rate Range</b>                |            | 40 to 300 bpm  |
| <b>Measurement Time</b>                |            | Typical at HR > 60 bpm<br>Auto/manual: 30 seconds (adult)<br>25 seconds (neonatal)<br>Stat: 20 seconds<br>Maximum time: 180 seconds (adult/pediatric)<br>90 seconds (neonates) |
| <b>Cuff Inflation Time</b>             |            | Typical for normal adult cuff: Less than 10 seconds<br>Typical for neonatal cuff: Less than 2 seconds  |

| NBP Performance Specifications         |                 |   |
|--|-----------------|---|
| <b>Initial Cuff Inflation Pressure</b> |                 | Adult: 165 $\pm 15$ mmHg<br>Pedi: 130 $\pm 15$ mmHg<br>Neo: 100 $\pm 15$ mmHg |
| <b>Auto Mode Repetition Times</b>      |                 | 1, 2, 2.5, 3, 5, 10, 15, 20, 30, 45, 60 or 120 minutes                        |
| <b>STAT Mode Cycle Time</b>            |                 | 5 minutes   |
| <b>Venipuncture Mode Inflation</b>     |                 |   |
| <b>Inflation Pressure</b>              | Adult           | 20 to 120 mmHg (3 to 16 kPa)  |
|  | Pediatric       | 20 to 80 mmHg (3 to 11 kPa)   |
|  | Neonatal        | 20 to 50 mmHg (3 to 7 kPa)  |
| <b>Automatic deflation after</b>       | Adult/pediatric | 170 seconds   |
|  | Neonatal        | 85 seconds  |

**Measurement Validation:** In adult and pediatric mode, the blood pressure measurements determined with this device comply with the American National Standard for Electronic or Automated Sphygmomanometers (ANSI/AAMI SP10 - 1992) in relation to mean error and standard deviation, when compared to intra-arterial or auscultatory measurements (depending on the configuration) in a representative patient population. For the auscultatory reference the 5th Korotkoff sound was used to determine the diastolic pressure. In neonatal mode, the blood pressure measurements determined with this device comply with the American National Standard for Electronic or Automated Sphygmomanometers (ANSI/AAMI SP10 - 1992 and AAMI/ANSI SP10A -1996) in relation to mean error and standard

deviation, when compared to intra-arterial measurements in a representative patient population.

| NBP Alarm Specifications | Range                                    | Adjustment                         |
|--------------------------|--|------------------------------------|
| <b>Systolic</b>          | Adult: 30 to 27 mmHg<br>(4 to 36 kPa)    | 10 to 30 mmHg:<br>2 mmHg (0.5 kPa) |
|                          | Pedi: 30 to 180 mmHg<br>(4 to 24 kPa)    | > 30 mmHg:<br>5 mmHg (1 kPa)       |
|                          | Neo: 30 to 130 mmHg<br>(4 to 17 kPa)     |                                    |
| <b>Diastolic</b>         | Adult: 10 to 245 mmHg<br>(1.5 to 32 kPa) |                                    |
|                          | Pedi: 10 to 150 mmHg<br>(1.5 to 20 kPa)  |                                    |
|                          | Neo: 10 to 100 mmHg<br>(1.5 to 13 kPa)   |                                    |
| <b>Mean</b>              | Adult: 20 to 255 mmHg<br>(2.5 to 34 kPa) |                                    |
|                          | Pedi: 20 to 160 mmHg<br>(2.5 to 21 kPa)  |                                    |
|                          | Neo: 20 to 120 mmHg<br>(2.5 to 16 kPa)   |                                    |

| NBP Overpressure Settings |                                |                     |
|---------------------------|--------------------------------|---------------------|
| <b>Adult</b>              | > 300 mmHg (40 kPa)<br>> 2 sec | not user adjustable |
| <b>Pedi</b>               | > 300 mmHg (40 kPa)<br>> 2 sec |                     |
| <b>Neo</b>                | > 150 mmHg (20 kPa)<br>> 2 sec |                     |

### Invasive Pressure and Pulse

Complies with IEC 60601-2-34:2000/EN60601-2-34:2000.

| Invasive Pressure Performance Specifications |            |   |
|--|------------|---|
| <b>Measurement Range</b>                     |            | -40 to 360 mmHg                                       |
| <b>Pulse Rate</b>                            | Range      | 25 to 350 bpm   |
|  | Accuracy   | ±1% Full Range  |
|  | Resolution | 1 bpm   |
| <b>Input Sensitivity</b>                     |            | Sensitivity: 5 $\mu$ V/V/mmHg<br>(37.5 $\mu$ V/V/kPa) |
|  |            | Adjustment range: ±10%                                |

| Invasive Pressure Performance Specifications |  |  |
|--|--|--|
| <b>Transducer</b>                            |  | Load Impedance: 200 to 2000 $\Omega$ (resistive)<br>Output Impedance: $\leq$ 3000 $\Omega$ (resistive) |
|  | <b>Frequency Response</b>              | dc to 12.5 Hz or 40 Hz   |
| <b>Zero Adjustment</b>                       | Range                                  | ±200 mmHg (±26 kPa)  |
|  | Accuracy                               | ±1 mmHg (±0.1 kPa)   |
|  | Drift                                  | Less than 0.1 mmHg/°C<br>(0.013 kPa/°C)  |
| <b>Gain Accuracy</b>                         | Accuracy                               | ±1%  |
|  | Drift                                  | Less than 0.05%/°C   |
|  | Non linearity and Hysteresis           | Error of $\leq$ 0.4% FS (@CAL 200 mmHg)  |
| <b>Overall Accuracy</b>                      | (including transducer)                 | ± 4% of reading or ± 4 mmHg<br>(± 0.5 kPa), whichever is greater                                       |
|  | <b>Volume displacement of CPJ840J6</b> | 0.1 mm <sup>3</sup> /100 mmHg  |

| Invasive Pressure Alarm Specifications | Range                                 | Adjustment  | Delay              |
|--|---------------------------------------|---|--------------------|
| <b>Pressure</b>                        | -40 to 360 mmHg<br>(-5.0 to 48 kPa)   | -40 to 30 mmHg<br>2 mmHg<br>(0.5 kPa)<br>> 30 mmHg<br>5 mmHg<br>(1 kPa) | max.<br>12 seconds |
| <b>Extreme High</b>                    | Difference to high limit 0 to 25 mmHg | 5 mmHg steps<br>(0.5 kPa)   |                    |
|  | Clamping at -40 to 360 mmHg           | 5 mmHg steps<br>(1.0 kPa)   |                    |
| <b>Extreme Low</b>                     | Difference to low limit 0 to 25 mmHg  | 5 mmHg steps<br>(0.5 kPa)   |                    |
|  | Clamping at 40 to 360 mmHg            | 5 mmHg steps<br>(1.0 kPa)   |                    |



| Invasive Pressure Alarm Specifications | Range  | Adjustment  | Delay              |
|--|--|---|--------------------|
| <b>Pulse</b>                           | 25 to 300 bpm  | Adult:<br>1 bpm steps<br>(25 to 40 bpm)<br>5 bpm steps<br>(40 to 300 bpm)<br>Pedi/Neo:<br>1 bpm steps<br>(25 to 50 bpm)<br>5 bpm steps<br>(50 to 300 bpm) |                    |
| <b>Tachycardia</b>                     | Difference to high limit 0 to 50 bpm<br>Clamping at 150 to 300 bpm | 5 bpm steps   | max.<br>14 seconds |
| <b>Bradycardia</b>                     | Difference to low limit 0 to 50 bpm<br>Clamping at 25 to 100 bpm   | 5 bpm steps   | max.<br>14 seconds |

## Temp

Complies with EN 12470-4:2000

| Temp Performance Specifications |            |   |
|---------------------------------|------------|---|
| <b>Temp</b>                     | Range      | -1 °C to 45 °C (30 °F to 113 °F)  |
|                                 | Resolution | 0.1 °C (0.2 °F)   |
|                                 | Accuracy   | ±0.1 °C (±0.2 °F)   |
| <b>Average Time Constant</b>    |            | Less than 10 seconds  |
| <b>Alarms</b>                   | Range      | -1 °C to 45 °C (30 °F to 113 °F)  |
|                                 | Adjustment | -1 to 35 °C (30 to 95 °F):<br>0.5 °C (1.0 °F) steps<br>35 to 45 °C (95 to 113 °F):<br>0.1 °C (0.2 °F) steps |

| Temp Alarm Specifications   | Range                            | Adjustment   |
|-----------------------------|----------------------------------|--|
| <b>Temp High/Low Alarms</b> | -1 °C to 45 °C (30 °F to 113 °F) | -1 °C to 35 °C (30 °F to 95 °F),<br>0.5 °C (1.0 °F) steps<br>35 °C to 45 °C (95 °F to 113 °F), 0.1 °C (0.2 °F) steps |

## CO<sub>2</sub>

The CO<sub>2</sub> measurement in the monitor, M3014A and M3015A/B complies with EN ISO 21647:2004 + Cor.1:2005 (except alarm system; alarm system complies with IEC 60601-2-49:2001).

| M3015A/B Microstream CO <sub>2</sub> Performance Specifications |            |  |
|---|------------|--|
| <b>CO<sub>2</sub></b>   | Range      | 0 to 98 mmHg (0 to 13 kPa), or 13% CO <sub>2</sub> , whichever is lower  |
|   | Accuracy   | Up to 5 minutes during warm-up:<br>±4 mmHg or 12%, whichever is greater<br>After 5 minutes warm-up:<br>0 to 40 mmHg (0 to 5.3 kPa):<br>±2.2 mmHg (±0.3 kPa)<br>Above 40 mmHg (5.3 kPa): ±(5.5% + (0.08%/mmHg above 40 mmHg)) of reading<br>These specifications are valid for 21% O <sub>2</sub> and N <sub>2</sub> balance, up to 35 °C ambient temperature, up to 60 rpm in adult mode and 100 rpm in neonatal mode. Outside of these conditions the accuracy reaches at a minimum ±4 mmHg or ±12% of the reading, whichever is greater. |
|   | Resolution | Numeric: 1.0 mmHg (0.1 kPa)<br>Wave: 0.1 mmHg (0.01 kPa)   |
|   | Stability  | Included in Accuracy specifications  |
|   | awRR       |  |
|   | Range      | 0 to 150 rpm   |
|   | Accuracy   | 0 to 40 rpm: ±1 rpm<br>41 to 70 rpm: ±2 rpm<br>71 to 100 rpm: ±3 rpm<br>>100 rpm: ±5% of reading   |
| <b>Warm-up Time</b>   |            | 5 minutes for full accuracy specification  |

### M3015A/B Microstream CO<sub>2</sub> Performance Specifications

|                                   |  |
|-----------------------------------|--|
| <b>Rise Time</b>                  | 190 ms for neonatal mode<br>(measured with FilterLine H for neonatal)<br>240 ms for adult mode<br>(measured with FilterLine H for adult) |
| <b>Sample Flow Rate</b>           | 50 + 15/-7.5 ml/minute   |
| <b>Gas Sampling Delay Time</b>    | Typical: 2.3 seconds<br>Maximum: 3 seconds   |
| <b>Sound Pressure</b>             | Acoustic noise: <45 dBA  |
| <b>Total System Response Time</b> | The total system response time is the sum of the delay time and the rise time.   |

### M3014A Mainstream CO<sub>2</sub> Performance Specifications

|                       |  |   |
|-----------------------|--|---|
| <b>CO<sub>2</sub></b> | <b>Range</b>   | 0 to 150 mmHg (0 to 20.0 kPa)   |
|                       | <b>Accuracy</b>  | after 2 minutes warmup:<br>For values between 0 and 40 mmHg: ±2.0 mmHg (±0.29 kPa)<br>For values from 41 to 70 mmHg: ±5% of reading<br>For values from 71 to 100 mmHg: ±8% of reading<br>The specifications are valid for standard gas mixtures, balance air, fully hydrated at 35°C, P <sub>abs</sub> = 760 mmHg, flow rate = 2 l/min. |
|                       | <b>Resolution</b>  | Numeric: 1.0 mmHg (0.1 kPa)<br>Wave: 0.1 mmHg (0.01 kPa)  |
|                       | <b>Stability:</b>  |   |
|                       | Short term drift<br>Long term drift  | ±0.8 mmHg over four hours<br>Accuracy specification will be maintained over a 120 hour period   |
| <b>awRR</b>           | <b>Range</b>   | 2 to 150 rpm  |
|                       | <b>Accuracy</b>  | ±1 rpm  |
| <b>Warm-up Time</b>   | 2 minutes with CO <sub>2</sub> transducer attached for full accuracy specification |   |
| <b>Response Time</b>  | Less than 60 ms (with adult or infant reusable or disposable adapter)              |   |

### M3014A Sidestream CO<sub>2</sub> Performance Specifications

|                                   |  |  |
|-----------------------------------|--|--|
| <b>CO<sub>2</sub></b>             | <b>Range</b>   | 0 to 150 mmHg (0 to 20.0 kPa)  |
|                                   | <b>Accuracy</b>  | after 2 minutes warmup:<br>For values between 0 and 40 mmHg: ±2.0 mmHg (±0.29 kPa)<br>For values from 41 to 70 mmHg: ±5% of reading<br>For values from 71 to 100 mmHg: ±8% of reading<br>For values from 101 to 150 mmHg: ±10% of reading<br>At respiration rates above 80 rpm, all ranges are ±12% of actual. The specifications are valid for gas mixtures of CO <sub>2</sub> , balance N <sub>2</sub> , dry gas at 760 mmHg within specified operating temperature range. |
|                                   | <b>Resolution</b>  | Numeric: 1.0 mmHg (0.1 kPa)<br>Wave: 0.1 mmHg (0.01 kPa)   |
|                                   | <b>Stability:</b>  |  |
|                                   | Short term drift<br>Long term drift  | ±0.8 mmHg over four hours<br>Accuracy specification will be maintained over a 120 hour period  |
| <b>awRR</b>                       | <b>Range</b>   | 2 to 150 rpm   |
|                                   | <b>Accuracy</b>  | ±1 rpm   |
| <b>Warm-up Time</b>               | 2 minutes with CO <sub>2</sub> sensor attached for full accuracy specification |  |
| <b>Sample Flow Rate</b>           | 50 ±10 ml/minute   |  |
| <b>Total System Response Time</b> | 3 seconds  |  |
| <b>Operating Temperature</b>      | 0°C to 40°C (32°F to 104°F)  |  |

### M3014A Mainstream and Sidestream CO<sub>2</sub> Humidity Correction Factor

Either BTPS or STPD can be selected as the humidity correction factor for the CO<sub>2</sub> readings. The formula for the correction calculation is:

$$P_{STPD} = P_{BTPS} \cdot \frac{P_{abs}}{P_{abs} - P_{H2O}}$$

Where  $p$  = partial pressure,  $P_{abs}$  = absolute pressure, and  $P_{H_2O}$  = 42 mmHg @35°C and 100% RH.

|  |  |
|--|--|
|  |  |
|  |  |
|  |  |
|  |  |
|  |  |
|  |  |
|  |  |
|  |  |
|  |  |

#### Mainstream CO<sub>2</sub> Humidity Correction Factor

Either BTPS or STPD can be selected as the humidity correction factor for the Mainstream CO<sub>2</sub> readings. The formula for the correction calculation is:

$$P_{STPD} = P_{BTPS} \cdot \frac{P_{abs}}{P_{abs} - P_{H_2O}}$$

Where  $p$  = partial pressure,  $P_{abs}$  = absolute pressure, and  $P_{H_2O}$  = 47 mmHg @37°C and 100% RH.

| CO <sub>2</sub> Alarm Specifications | Range                         | Adjustment                | Delay   |
|--------------------------------------|-------------------------------|---------------------------|---|
| <b>etCO<sub>2</sub> High</b>         | 20 to 95 mmHg (2 to 13 kPa)   | 1 mmHg (0.1 kPa)          | M3014A: less than 14 seconds                                    |
| <b>etCO<sub>2</sub> Low</b>          | 10 to 90 mmHg (1 to 12 kPa)   |                           | M3015A/B: less than 18 seconds.                                 |
| <b>imCO<sub>2</sub> High</b>         | 2 to 20 mmHg (0.3 to 3.0 kPa) | steps of 1 mmHg (0.1 kPa) | M3014A: less than 14 seconds<br>M3015A/B: less than 18 seconds. |

| CO <sub>2</sub> Alarm Specifications | Range   | Adjustment  | Delay   |
|--------------------------------------|---|---|---|
| <b>awRR High</b>                     | Adult/pedi: 10 to 100 rpm<br>Neo: 30 to 150 rpm | under 20 rpm: 1 rpm steps<br>over 20 rpm: 5 rpm steps | M3014A: less than 14 seconds<br>M3015A/B: less than 18 seconds.   |
| <b>awRR Low</b>                      | Adult/pedi: 0 to 95 rpm<br>Neo: 0 to 145 rpm    |   | M3015A/B: settings <20 rpm: less than 8 seconds<br>>20 rpm: less than 18 seconds<br>M3014A settings <20 rpm: less than 4 seconds<br>>20 rpm: less than 14 seconds |
| <b>Apnea delay</b>                   | 10 to 40 seconds                                | 5 second steps  | set apnea delay time + 4 seconds (M3014A) or 8 seconds (M3015A/B)   |

## Ordering Information

Ordering information for the M8102A patient monitor is given here.

| Parameters   | M8102A |
|--|--------|
| <b>Order one Bxx option</b>                        |        |
| ECG, Resp, NBP, SpO <sub>2</sub>                   | B20    |
| ECG, Resp, NBP, SpO <sub>2</sub> , Press/Temp      | B22    |
| ECG, Resp, NBP, SpO <sub>2</sub> , CO <sub>2</sub> | B23    |

## Application OptionsXDS Connectivity

| Application Options                    | M8102A |
|--|--------|
| Full Arrhythmia Capability             | C01    |
| 12-Lead ECG Application (conventional) | C12    |
| ST Map                                 | C13    |
| Full Networking                        | C15    |

| Options                  | M8102A           |
|--------------------------|------------------|
| 4-Wave XDS Connectivity  | X04              |
| 6-Wave XDS Connectivity  | X06 <sup>a</sup> |
| XDS Remote Control       | X20              |
| XDS Clinical Workstation | X30              |

<sup>a</sup> BU approval required

## SpO<sub>2</sub> Technology Choice

| Options   | M8102A   |
|---|----------|
| Philips FAST SpO <sub>2</sub>   | Standard |
| Substitute Philips FAST SpO <sub>2</sub> with Nellcor OxiMax SpO <sub>2</sub> | SP4      |

## Hardware Options

| Hardware Add-Ons           | M8102A |
|----------------------------|--------|
| Anti-slip pad              | E18    |
| Carrying strap             | E19    |
| MMS Mount                  | E20    |
| Protective cover           | E23    |
| Add 1X Lithium-Ion battery | E24    |
| Add 2X Lithium-Ion battery | E26    |
| SN3 ECG Sync Cable         | SN3    |

## Interface Options

| Interfaces                        | M8102A           |
|-----------------------------------|------------------|
| IntelliVue 802.11 Bedside Adapter | J35 <sup>a</sup> |
| Instrument Telemetry 1.4 GHz      | J45 <sup>a</sup> |
| Instrument Telemetry 2.4 GHz      | J47 <sup>a</sup> |
| Short Range Radio Interface       | J46 <sup>a</sup> |

<sup>a</sup> May not be available in all geographies

## Upgrade Options M8102AU

| Options                              | MP2<br>M8102AU   |
|--------------------------------------|------------------|
| <b>Application</b>                   |                  |
| <u>Full Arrhythmia capability</u>    | <u>C01</u>       |
| conv. 12-lead ECG                    | C12              |
| <u>Full Networking</u>               | <u>C15</u>       |
| <b>XDS external display solution</b> |                  |
| 4-wave XDS connectivity              | X04              |
| XDS Remote Control                   | X20              |
| XDS Clinical Workstation             | X30              |
| <b>Interfaces</b>                    |                  |
| 802.11 Wireless Interface            | J35 <sup>a</sup> |
| Instrument Telemetry 1.4 GHz         | J45 <sup>a</sup> |
| Short Range Radio                    | J46 <sup>a</sup> |
| Instrument Telemetry 2.4 GHz         | J47 <sup>a</sup> |
| <b>Software Upgrade</b>              |                  |
| Current Software Revision            | SUJ <sup>9</sup> |

<sup>a</sup> May not be available in all geographies

## Sensors and Disposables

| Accessory   | M8102A |
|---|--------|
| 3-lead Accessories Bundle ICU-AAMI                          | G06    |
| Tyco low cost cable   |        |
| 3-lead Accessories Bundle ICU-IEC                           | G07    |
| Tyco low cost cable   |        |
| 5-lead Accessories Bundle ICU-AAMI                          | G08    |
| Tyco low cost cable   |        |
| 5-lead Accessories Bundle ICU-IEC                           | G09    |
| Tyco low cost cable   |        |
| 5-lead Accessories Bundle ICU-AAMI                          | H06    |
| 5-lead Accessories Bundle ICU-IEC                           | H07    |
| 5-lead Accessories Bundle OR-AAMI                           | H08    |
| 5-lead Accessories Bundle OR-IEC                            | H09    |
| Accessories Bundle Neonatal-AAMI                            | H14    |
| Accessories Bundle Neonatal-IEC                             | H15    |
| 3-lead Accessories Bundle ICU-AAMI                          | H16    |
| 3-lead Accessories Bundle ICU-IEC                           | H17    |
| 3-lead Accessories Bundle OR-AAMI                           | H18    |
| 3-lead Accessories Bundle OR-IEC                            | H19    |
| CO <sub>2</sub> Mainstream Sensor                           | N01    |
| Reusable Adult Airway Adapter (msCO <sub>2</sub> )          | N02    |
| Reusable Infant Airway Adapter (msCO <sub>2</sub> )         | N03    |
| Single Use Adult Airway Adapter (msCO <sub>2</sub> )        | N04    |
| Single Use Infant Airway Adapter (msCO <sub>2</sub> )       | N05    |
| CO <sub>2</sub> Sidestream Sensor                           | N11    |
| Non-intubated Adult Airway Adapter (ssCO <sub>2</sub> )     | N12    |
| Non-intubated pediatric Airway Adapter (ssCO <sub>2</sub> ) | N13    |
| Intubated Adult Airway Adapter (ssCO <sub>2</sub> )         | N14    |
| Intubated Pediatric Airway Adapter (ssCO <sub>2</sub> )     | N15    |

## Related Products

M3086A Support Tool

## Mounting Information

The Intellivue MP2 Roll Stand Mounting Kit (Order No. 989803153021) is compatible with the table top mount and the standard mounting plate. For information on other mounting hardware, contact your local Philips sales representative. For GCX mounting hardware information, see [www.gcx.com/philips](http://www.gcx.com/philips).

## Documentation

All documentation is available in .pdf format on documentation CD-ROM. Additionally, a printed copy of the Instructions for Use and Quick Guide ships with each monitor.

- Instructions for Use (printed)
- Quick Guide (printed)
- Installation and Service Guide
- Configuration Guide
- Documentation CD-ROM
- Training Guide (printed)
- Computer Based Training (optional)

## Carry Case & NVG Display Filter

| Accessory                      | Part No.     |
|--------------------------------|--------------|
| MP2 Carry Case Std Red         | 989803163331 |
| MP2 Carry Case Std Blk         | 989803163341 |
| MP2 Carry Case Mini Red        | 989803163351 |
| MP2 Carry Case Mini Blk        | 989803163361 |
| MP2 Carry Case Large Blk       | 989803163371 |
| MP2 Carry Case Replacement Kit | 989803163631 |
| MP2 NVG Display Filter         | 989803163381 |

## Battery Extension

The Philips Battery Extension provides additional battery runtime of up to 6 hours for intra-hospital patient transport and concurrent CO<sub>2</sub> measurement with additional invasive blood pressure and temperature measurement.

| Accessory                 | Part No. |
|---------------------------|----------|
| Philips Battery Extension | 865297   |

## ECG Accessories



This symbol indicates that the cables and accessories are designed to have special protection against electric shocks (particularly regarding allowable leakage currents), and are defibrillator proof.

## Trunk Cables

|          | 3-Electrode Cable Set | 5-Electrode Cable Set | 6-Electrode Cable Set | 10-Electrode Cable set<br>(5+5) | 10-Electrode Cable set<br>(6+4) |
|----------|-----------------------|-----------------------|-----------------------|---------------------------------|---------------------------------|
| Part No. | M1669A                | M1668A                | M1667A                | M1663A                          | M1665A                          |
| Length   | 2.7 m                 | 2.7 m                 | 2.7 m                 | 2.0 m                           | 2.7 m                           |

## 3-Electrode Cable Sets

| Description           | Length | AAMI Part No. | IEC Part No. |
|-----------------------|--------|---------------|--------------|
| OR Grabber shielded   | 1.0 m  | M1675A        | M1678A       |
| ICU Grabber shielded  | 1.0 m  | M1671A        | M1672A       |
| ICU snap shielded     | 1.0 m  | M1673A        | M1674A       |
| ICU Clip non-shielded | 0.45 m | M1622A        | —            |
| ICU Clip non-shielded | 0.7 m  | M1624A        | M1626A       |

## 5-Electrode Cable Sets

| Description               | Length      | AAMI Part No. | IEC Part No. |
|---------------------------|-------------|---------------|--------------|
| OR Grabber shielded       | 1.0 m/1.6 m | M1973A        | M1974A       |
| ICU Grabber shielded      | 1.0 m/1.6 m | M1968A        | M1971A       |
| ICU Snap shielded         | 1.0 m/1.6 m | M1644A        | M1645A       |
| ICU Miniclip non-shielded | 0.7 m/1.3 m | M1647A        | M1648A       |

## 6-Electrode Cable Sets

| Description | Length      | AAMI Part No. | IEC Part No. |
|-------------|-------------|---------------|--------------|
| OR Grabber  | 1.0 m/1.6 m | M1684A        | M1685A       |
| ICU Grabber | 1.0 m/1.6 m | M1680A        | M1681A       |
| ICU Snap    | 1.0 m/1.6 m | M1682A        | M1683A       |

## 10-Electrode (5+5) Cable Sets

| Description                               | Length | AAMI Part No. | IEC Part No. |
|---|--------|---------------|--------------|
| ICU Grabber, chest, shielded              | 1.0 m  | M1976A        | M1978A       |
| ICU Snap, chest, shielded                 | 1.0 m  | M1602A        | M1604A       |
| OR Grabber, chest, shielded               | 1.0 m  | M1979A        | M1984A       |
| For Limb Leads see 5-electrode cable sets |        |               |              |

## 10-Electrode (6+4) Cable Sets

| Description                               | Length | AAMI Part No. | IEC Part No. |
|---|--------|---------------|--------------|
| ICU Grabber, chest, shielded              | 1.0 m  | M1532A        | M1533A       |
| ICU Snap, chest, shielded                 | 1.0 m  | M1537A        | M1538A       |
| OR Grabber, chest, shielded               | 1.0 m  | M1557A        | M1558A       |
| For Limb Leads see 6-electrode cable sets |        |               |              |

## One-piece Cables

| Description         | Length | AAMI Part No. | IEC Part No. |
|---------------------|--------|---------------|--------------|
| 3-lead Grabber, ICU | 1.0 m  | 989803143181  | 989803143171 |
| 5-lead Grabber, ICU | 1.0 m  | 989803143201  | 989803143191 |

### Radio-translucent Cables







Pack of five single wires, radio-translucent, 0.9m, M1649A





## Set Combiners and Organizers

| Set combiners and organizers                                       |             | Part No.         |
|--|-------------|------------------|
| Set combiner   | 3-electrode | M1501A           |
|  | 5-electrode | M1502A           |
| Set organizer for shielded leadsets - grabber and snap             | 3-electrode | M1503A           |
|  | 4-electrode | M1664A           |
|  | 5-electrode | M1504A           |
| Set organizer for non-shielded lead sets - miniclip                | 3-electrode | M1636A           |
|  | 5-electrode | M1638A           |
| Bedsheet clip  |             | M1509A           |
| Replacement red cover for trunk cable (for 5-electrode cable sets) |             | 98980814886<br>1 |





## Philips FAST SpO<sub>2</sub> Accessories

### Philips Reusable Sensors

| Part Number | Description   | Connector Type   |
|-------------|---|--|
| M1191A/B    | Adult Sensor (2 m cable)  | Philips 8-pin<br>  |
| M1191AL/BL  | Adult Sensor (3 m cable)  |  |
| M1191T      | Adult Sensor (requires M1943A (1.1 m) or M1943AL (3 m) adapter cable)                 | Generic D-Sub<br> |
| M1192A      | Small Adult/Pediatric sensor (1.5 m cable)  | Philips 8-pin<br> |
| M1192T      | Small Adult Pediatric sensor (requires M1943A (1.1 m) or M1943AL (3 m) adapter cable) | Generic D-Sub<br> |
| M1193A      | Neonatal Hand/Foot Sensor (1.5 m cable)   | Philips 8-pin<br> |
| M1193T      | Neonatal Sensor (requires M1943A (1.1 m) or M1943AL (3 m) adapter cable)              | Generic D-Sub<br> |

| Part Number | Description  | Connector Type   |
|-------------|--|--|
| M1194A      | Adult/Pediatric Clip Sensor (ear) (1.5 m cable)                            | Philips 8-pin<br> |
| M1195A      | Infant Sensor (1.5 m cable)  | Philips 8-pin<br> |
| M1196A      | Adult Clip Sensor (3 m cable)  | Philips 8-pin<br> |
| M1196T      | Adult Clip Sensor (requires M1943A (1.1 m) or M1943AL (3 m) adapter cable) | Generic D-Sub<br> |

### Philips Disposable Sensors

| Part Number | Description   | Connector Type   |
|-------------|---|--|
| M1131A      | Adult/Pediatric Sensor (requires M1943A (1.1 m) or M1943AL (3 m) adapter cable)                     | Generic D-Sub<br> |
| M1132A      | Infant Sensor (requires M1943A (1.1 m) or M1943AL (3 m) adapter cable)                              | Generic D-Sub<br> |
| M1133A      | Adult/Infant/Neonatal Sensor (requires M1943A (1.1 m) or M1943AL (3 m) adapter cable)               | Generic D-Sub<br> |
| M1134A      | Adhesive-free Neonatal/Infant/Adult Sensor (requires M1943A (1.1 m) or M1943AL (3 m) adapter cable) | Generic D-Sub<br> |



## Nellcor Accessories

### NELLCOR™ Disposable Sensors<sup>1</sup>:

Purchase Nellcor sensors directly from Nellcor.

| Product Number                 | Description                            | Philips Part Number |
|--------------------------------|--|---------------------|
| OxiMax MAXA <sup>a</sup>       | Adult SpO <sub>2</sub> Sensor          | M1904B <sup>b</sup> |
| OxiMax MAXAL <sup>a</sup>      | Adult XL SpO <sub>2</sub> Sensor       | n/a                 |
| OxiMax MAXP <sup>a</sup>       | Pediatric SpO <sub>2</sub> Sensor      | M1903B <sup>b</sup> |
| OxiMax MAXI <sup>a</sup>       | Infant SpO <sub>2</sub> Sensor         | M1902B <sup>b</sup> |
| OxiMax MAXN <sup>a</sup>       | Neonatal-Adult SpO <sub>2</sub> Sensor | M1901B <sup>b</sup> |
| Oxisensor II D-25 <sup>a</sup> | Adult Sensor                           | n/a                 |
| Oxisensor II D-20 <sup>a</sup> | Pediatric Sensor                       | n/a                 |
| Oxisensor II I-20 <sup>a</sup> | Infant Sensor                          | n/a                 |
| Oxisensor II N-25 <sup>a</sup> | Neonatal Sensor                        | n/a                 |
| OxiCliq A <sup>c</sup>         | Adult SpO <sub>2</sub> Sensor          | n/a                 |
| OxiCliq P <sup>c</sup>         | Pediatric SpO <sub>2</sub> Sensor      | n/a                 |
| OxiCliq I <sup>c</sup>         | Infant SpO <sub>2</sub> Sensor         | n/a                 |
| OxiCliq N <sup>c</sup>         | Neonatal-Adult SpO <sub>2</sub> Sensor | n/a                 |

a Requires M1943 A(L) adapter cable

b not available from Philips in the U.S.A.

c Requires M1943 A(L) and Nellcor OC3 adapter cables

## Masimo Accessories

Adapter cables are available from Philips and also from Masimo.

Sensors are available directly from Masimo.

### MASIMO LNOP<sup>®2</sup> Reusable Sensors:

| Product Number | Description                      |
|----------------|----------------------------------|
| DCI            | Reusable Finger Sensor           |
| DCIP           | Reusable Pediatric Finger Sensor |
| YI             | Reusable Multi-Site Sensor       |
| TC-I           | Reusable Ear Sensor              |

### MASIMO LNCS<sup>®1</sup> Reusable Sensors:

| Product Number | Description         |
|----------------|---------------------|
| LNCS DCI       | Adult Sensor        |
| LNCS DCIP      | Pediatric Sensor    |
| LNCS-TCI       | Reusable Ear Sensor |

### MASIMO LNOP<sup>®</sup> Disposable Adhesive Sensors:

| Product Number | Description                            |
|----------------|--|
| Adt            | Adult Adhesive Sensor                  |
| Adtx           | Adult Adhesive Sensor                  |
| Pdt            | Pediatric Adhesive Sensor              |
| Pdtx           | Pediatric Adhesive Sensor              |
| INF-L          | Infant Adhesive Sensor                 |
| Neo-L          | Neonate Adhesive Sensor                |
| NeoPt-L        | Sensitive Skin Neonate Adhesive Sensor |

### MASIMO LNCS<sup>®</sup> Disposable Adhesive Sensors:

| Product Number | Description                                  |
|----------------|--|
| Adtx           | Adult Adhesive Sensor                        |
| Pdtx           | Pediatric Adhesive Sensor                    |
| INF-L          | Infant Toe Sensor                            |
| Neo-L          | Neo Foot Sensor or Adult Finger Sensor       |
| NeoPt-L        | Neo Pre-Term Sensitive Skin Adhesive Sensors |



The Philips M8102A with Philips FAST SpO<sub>2</sub> technology uses Masimo certified pulse oximetry for reduced noise and low perfusion performance with Masimo Sensors under the Masimo NR&LP protocol available from Masimo.

<sup>1</sup> Nellcor™, Durasensor™, Dura-Y™, Oxiband™, OxiCliq™, OxiMax™, MAXFAST™, are trademarks of Covidien AG and/or its affiliates

<sup>2</sup> LNOP and LNCS are federally registered trademarks of Masimo Corporation

## Extension/Adapter Cables:

| Part Number              | Description   |
|--------------------------|---|
| M1941A                   | Extension Cable (2m) (8-pin to 8-pin)   |
| M1943A                   | Adapter Cable (1.1 m) for Philips and Nellcor disposable sensors (8-pin to 9-pin D-Sub) |
| M1943AL                  | Adapter Cable (3m) for Philips and Nellcor disposable sensors (8-pin to 9-pin D-Sub)    |
| Nellcor OC3              | Adapter cable for OxiCliq Sensors (available from Nellcor only)                         |
| LNOP MP12 (451261000761) | LNOP MP Series Patient Cable (3.6 m) Adapter Cable for Masimo LNOP Sensors              |
| LNC MP10 (989803148221)  | LNCS MP Series Patient Cable (3.0 m) Adapter Cable for Masimo LNCS Sensors              |

## Nellcor OxiMax SpO<sub>2</sub> Accessories (for M8102A #SP4)

| Product Number                      | Description <sup>a/</sup><br>Preferred Application Site | Comments   |
|-------------------------------------|---|--|
| Philips Disposable Sensors          |   |  |
| M1904B <sup>b</sup>                 | Adult Sensor  | Requires M1943NL adapter cable (included with #A02).   |
| M1903B <sup>b</sup>                 | Pediatric Sensor  |  |
| M1902B <sup>b</sup>                 | Infant Sensor   |  |
| M1901B <sup>b</sup>                 | Neonatal/Adult Sensor                                   |  |
| NELLCOR OxiMax Sensors <sup>c</sup> |   |  |
| MAXA                                | Adult SpO <sub>2</sub> Sensor                           | Requires M1943NL adapter cable.  |
| MAXAL                               | Adult XL SpO <sub>2</sub> sensor                        |  |
| MAXP                                | Pediatric SpO <sub>2</sub> Sensor                       |  |
| MAXI                                | Infant SpO <sub>2</sub> Sensor                          |  |
| MAXN                                | Neonatal-Adult Sensor                                   |  |
| OxiCliq A                           | Adult   | Requires M1943NL adapter cable. <b>Additionally</b> the Nellcor OC3 adapter cable is needed. |
| OxiCliq P                           | Pediatric   |  |
| OxiCliq I                           | Infant  |  |
| OxiCliq N                           | Neonatal  |  |

| Product Number           | Description <sup>a/</sup><br><i>Preferred Application Site</i> | Comments                        |
|--------------------------|--|---------------------------------|
| MAXR                     | Adult SpO <sub>2</sub> Nasal Sensor                            | Requires M1943NL adapter cable. |
| MAXFAST                  | Forehead SpO <sub>2</sub> Sensor                               |                                 |
| Oxiband OXI-A/N          | Adult-Neonatal SpO <sub>2</sub> Sensor with Wraps              |                                 |
| Oxiband OXI-P/I          | Pediatric-Infant SpO <sub>2</sub> Sensor with Wraps            |                                 |
| SoftCare SC-A            | Adult SpO <sub>2</sub> Sensor                                  |                                 |
| SoftCare SCNEO-I         | Neonatal SpO <sub>2</sub> Sensor                               |                                 |
| SC-PR-I                  | Preemie SpO <sub>2</sub> Sensor                                |                                 |
| Durasensor DS100A-1      | Adult SpO <sub>2</sub> Sensor                                  |                                 |
| Dura-Y D-YS              | SpO <sub>2</sub> Sensor  |                                 |
| Extension/Adapter Cables |  |                                 |
| M1943NL                  | Adapter Cable (3 m)  | Available from Nellcor only.    |
| Nellcor OC3              | Adapter Cable for Nellcor OxiCliqu sensors                     |                                 |

a For application site, please refer also to the Instructions for Use provided with the sensors.

b Philips disposable sensors M1901B, M1902B, M1903B and M1904B are not available in the USA.

c Can only be ordered from Nellcor.

## Non Invasive Blood Pressure Accessories



These cuffs and tubings are designed to have special protection against electric shocks (particularly regarding allowable leakage currents), and are defibrillator proof.

| <b>Multi-Patient Comfort Cuffs and Disposable Cuffs</b> |                 |               |
|---|-----------------|---------------|
| Patient Category  | Disposable cuff | Reusable cuff |
| Adult (Thigh)   | M1879A          | M1576A        |
| Large Adult   | M1878A          | M1575A        |
| Adult   | M1877A          | M1574A        |
| Small Adult   | M1876A          | M1573A        |
| Pediatric   | M1875A          | M1572A        |
| Infant  | M1874A          | M1571A        |

Tubing: Use M1598B or M1599B

| Reusable Cuff Kits                     | Part No. |
|--|----------|
| Infant, pediatric, small adult, adult  | M1577A   |
| Small adult, adult, large adult, thigh | M1578A   |

| Reusable Cuff Kits  | Part No. |
|---|----------|
| Infant, pediatric, small adult, adult, large adult, thigh | M1579A   |

| Adult/Pediatric Antimicrobial Coated Reusable cuffs |                    |                       |                      |
|---|--------------------|-----------------------|----------------------|
| Cuff Size (color)                                   | Circumference (cm) | Bladder Width         | Single-Hose Part No. |
| Infant (orange)                                     | 9.0 - 14.8         | 5.4 cm<br>2.1 inches  | M4552A               |
| Pediatric (green)                                   | 13.8 - 21.5        | 8.0 cm<br>3.1 inches  | M4553A               |
| Small Adult (royal blue)                            | 20.5 - 28.5        | 10.6 cm<br>4.2 inches | M4554A               |
| Adult (navy blue)                                   | 27.5 - 36.5        | 13.5 cm<br>5.3 inches | M4555A               |
| Adult X-long (navy blue)                            | 27.5 - 36.5        | 13.5 cm<br>5.3 inches | M4556A               |
| Large Adult (burgundy)                              | 35.5 - 46.0        | 17.0 cm<br>6.7 inches | M4557A               |
| Large Adult X-long (burgundy)                       | 35.5 - 46.0        | 17.0 cm<br>6.7 inches | M4558A               |
| Thigh (grey)  | 45 - 56.5          | 21.0 cm<br>8.3 inches | M4559A               |

Tubing: Use M1598B or M1599B

| Adult/Pediatric Soft Single Patient Single-Hose Disposable Cuffs |                    |               |                          |
|--|--------------------|---------------|--------------------------|
| Patient Category   | Limb Circumference | Bladder Width | Disposable cuff Part No. |
| Adult (Thigh)  | 45.0 - 56.5 cm     | 20.4 cm       | M4579A                   |
| Large Adult X-long   | 35.5 - 46.0 cm     | 16.4 cm       | M4578A                   |
| Large Adult  | 35.5 - 46.0 cm     | 16.4 cm       | M4577A                   |
| Adult X-long   | 27.5 - 36.5 cm     | 16.4 cm       | M4576A                   |
| Adult  | 27.5 - 36.5 cm     | 13.1 cm       | M4575A                   |
| Small Adult  | 20.5 - 28.5 cm     | 10.4 cm       | M4574A                   |
| Pediatric  | 15.0 - 21.5 cm     | 8.0 cm        | M4573A                   |
| Infant   | 9.0 - 15.0 cm      | 5.6 cm        | M4572A                   |

Tubing: Use M1598B or M1599B

| Neonatal/Infant Cuffs (Disposable, non-sterile) |                    |               |          |
|---|--------------------|---------------|----------|
| Cuffs   | Limb Circumference | Bladder Width | Part No. |
| Size 1  | 3.1 to 5.7 cm      | 2.2 cm        | M1866A   |
| Size 2  | 4.3 to 8.0 cm      | 2.8 cm        | M1868A   |
| Size 3  | 5.8 to 10.9 cm     | 3.9 cm        | M1870A   |
| Size 4  | 7.1 to 13.1 cm     | 4.7 cm        | M1872A   |

Tubing: Use M1596B or M1597B

| Cuff Tubing |             |        |
|-------------|-------------|--------|
| Adult       | 1.5 m /4.9' | M1598B |
|             | 3.0 m/9.8'  | M1599B |
| Neonatal    | 1.5 m /4.9' | M1596B |
|             | 3.0 m/9.8'  | M1597B |

## Temperature Accessories

| Temperature Probes                            | Part No. |
|---|----------|
| Reusable                                      |          |
| General purpose probe                         | 21075A   |
| Small flexible vinyl probe (Infant/Pediatric) | 21076A   |
| Attachable surface probe                      | 21078A   |
| Disposable                                    |          |
| General purpose probe                         | M1837A   |
| Skin probe                                    | 21091A   |
| Esophageal/Stethoscope Probe (12 French)      | 21093A   |
| Esophageal/Stethoscope Probe (18 French)      | 21094A   |
| Esophageal/Stethoscope Probe (24 French)      | 21095A   |
| Foley Catheter Probe (12 French)              | M2255A   |
| Foley Catheter Probe (16 French)              | 21096A   |
| Foley Catheter Probe (18 French)              | 21097A   |
| Adapter cable 1.5 m/4.9'                      | 21082B   |
| Adapter cable 3.0 m/9.8'                      | 21082A   |

## PRESS Accessories



These transducers and accessories are designed to have special protection against electric shocks (particularly regarding allowable leakage currents), and are defibrillator proof.

| Pressure Transducers and Accessories                            | Part No. |
|---|----------|
| <b>Reusable</b>   |          |
| Reusable pressure transducer<br>5 $\mu$ V/V/mmHg sensitivity    | CPJ840J6 |
| Sterile disposable pressure domes for<br>CPJ840J6 (pack of 50)  | CPJ84022 |
| Transducer holder for CPJ840J6<br>(pack of 4)                   | CPJ84046 |
| IV pole mount for CPJ840J6                                      | CPJ84447 |
| <b>Disposable (EU/EFTA only. Not available in USA)</b>          |          |
| Single channel disposable sensor kit (20)                       | M1567A   |
| Dual channel disposable sensor kit (20)                         | M1568A   |
| Transducer holder for M1567/8A                                  | M2271A   |
| IV pole mount for M1567/8A                                      | M2272C   |
| Adapter cable for disposable sensor kit,<br>3.0 m, for M1567/8A | M1634A   |



**Philips Healthcare is part of  
Royal Philips Electronics**

**How to reach us**

[www.philips.com/healthcare](http://www.philips.com/healthcare)

[healthcare@philips.com](mailto:healthcare@philips.com)

fax: +31 40 27 64 887

Asia

+852 2821 5888

Europe, Middle East, Africa

+49 7031 463 2254

Latin America

+55 11 2125 0744

North America

+1 425 487 7000

800 285 5585 (toll free, US only)



M8102A complies with the requirements of  
the Council Directive 93/42/EEC of 14 June  
1993 (Medical Device Directive).

Please visit [www.philips.com/](http://www.philips.com/)



© 2012 Koninklijke Philips Electronics N.V.  
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Published in The Netherlands.  
4522 962 81041 \* FEB 2012