Owner's Manual
for model: **BD-3000S**

SAMSUNG AMERICA

Made in Taiwan

105 Challenger Road, 6th Floor
Ridgefield Park, NJ 07660

[www. SAMUNGHHealthyLiving.com](http://www.SAMUNGHHealthyLiving.com)

DELUXE

Deluxe Inflate Blood Pressure Monitor with Heart Sense®

- Large Three Row Digital Display
- 3 Zones of 30 Memory - 90 Memory Total
 - Automatic Power Off
 - Time and Date Indications
 - Pulse Display
- AC Capability and Adapter

Important Warnings/Medical Disclaimer

Before using your blood pressure monitor, you should read and understand all instructions and follow all warnings.

Information in this manual is provided for informational purposes only. This manual and product are not meant to be a substitute for the advice provided by your own physician or other medical professional. You should not use the information contained herein or this product for diagnosing or treating a health problem, disease, or prescribing any medication. If you have or suspect that you have a medical problem, promptly contact your healthcare provider.

In This Manual

Healthy Living®	3-4
Commonly Asked Questions About Blood Pressure	5-8
Intended Use	8
Measurement Method	9
Main Unit and Accessories	10
Special Features	11
Accuracy	12
Special Conditions	12
Digital Display Symbols	13
Battery Installation	13
AC Capability and Adapter	14
Selecting Time and Date	14
Important Reminders	15-16
Applying the Cuff	17
Taking a Blood Pressure Measurement	18-19
3 Zones of 30 Memory	20
Taking Care of Your Blood Pressure Monitor	21
Troubleshooting	22-23
Five Year Limited Warranty	24-25
Specifications	26
Blood Pressure Diary	27-30

Healthy Living®

This blood pressure monitor is the first step to living a healthier life. But monitoring your blood pressure at home is just the beginning, it's also important to:

EXERCISE REGULARLY - Start out by taking a quick walk each day. Before you know it, you'll feel energized. That's because along with being an easy way to shed a few pounds, walking also helps bones and heart to become stronger.

EAT RIGHT - Remember, you are what you eat! So when you reach for a snack, try grabbing a piece of fruit or a vegetable. When cooking meals at home try to use fats and oils sparingly. Sometimes a healthy diet is the best medicine of all! It's also a good idea to monitor your sodium intake since about one-third of individuals who have high blood pressure are also sodium sensitive. Consult your doctor regarding over-the-counter supplements and herbs such as onion and garlic. Both of these plants have been historically linked to producing mild anti-hypertension effects.

BE HAPPY! - A positive attitude is a big part of a healthy body. Try not to get too stressed out over things. Daily relaxation techniques such as gentle yoga or even quiet time by yourself can do a lot of good for the mind and body. Stay focused on the good things in life. Begin by congratulating yourself for that short walk you just took or the smart food choice you made. It's your first step to a healthier life!

KEEP AN EYE ON YOUR WEIGHT! - Surprisingly, losing even a small amount of weight can often times help control high blood pressure. Talk with your doctor to find a "healthy" weight for you and strive to come within a range of it.

Remember to consult your physician prior to beginning exercise, vitamin/mineral, dietary changes or programs and blood pressure monitoring.

Commonly Asked Questions About Blood Pressure


Why Should I Monitor My Blood Pressure At Home?

By monitoring at home, you can manage your health between doctor visits. Home monitoring also provides you with a record of blood pressure measurements over time to give to your doctor - a valuable tool in monitoring hypertension.

Home monitoring may also help alleviate *White Coat Hypertension*. This term refers to the medical phenomena which causes many individuals' blood pressure to be higher when taken in a medical office setting. That's because many times the stress of being in the doctor's office may cause blood pressure to rise as much as 20-50 points higher than normal. When monitoring at home, you may get a truer indication of your blood pressure in a normal, relaxed state.

What Is Blood Pressure?

The commonly used term "blood pressure" refers to the force or *pressure* that is exerted on the body's arteries as blood flows through them. Each time the heart beats or *contracts*, it produces pressure in the arteries. This blood pressure is what moves blood through the body, supplying oxygen and nutrients to every organ.

Commonly Asked Questions About Blood Pressure

When the heart contracts it sends blood out into the body. This pressure is referred to as **systolic**.

The pressure that is created as the heart pulls blood back through is called **diastolic**.


Blood pressure is recorded as systolic over diastolic, as shown here.

What Is Considered Normal Blood Pressure?

The Seventh Report of the Joint National Committee on Prevention, Detection, Evaluation, and Treatment of High Blood Pressure from the National Heart, Lung, and Blood Institute has set forth the following guidelines regarding high blood pressure. Guidelines provided are for information purposes only. Please consult with your physician for proper diagnosis.

Commonly Asked Questions About Blood Pressure

Why Does My Blood Pressure Fluctuate?

Because blood pressure changes with every beat of the heart it is in constant fluctuation throughout any given day. In addition to these natural changes there are other factors that may affect blood pressure such as:

- Temperature (too hot or cold)
- Humidity
- Restlessness
- Lack of sleep
- Physical exercise
- Constipation
- Stress
- Fear
- Anger
- Anxiety
- Food and beverage consumption
- Individual physical conditions

In addition to these factors you should talk with your physician to see what may be affecting you.

How Do I Know I Have High Blood Pressure?

It's nearly impossible to know if your blood pressure is elevated without having it checked. People have the common misconception that one can "feel" their blood pressure rise. It is because of

Commonly Asked Questions About Blood Pressure

this mystery that hypertension has been named “the silent killer.” Don’t let it sneak up on you! Start monitoring your blood pressure every day with your new Healthy Living[®] monitor.

Why Is It Important To Control My Blood Pressure?

Left uncontrolled, hypertension can place an incredible amount of stress on your heart and arteries, both of which become forced to work harder to keep blood flowing in a healthy manner throughout the body. Over time, this stress may result in health problems including over-enlargement of the heart, loss of elasticity of the arteries, heart disease, stroke and even death. Making some simple lifestyle changes now may help you to stay healthy in the future.

Intended Use


Healthy Living[®] blood pressure monitors are convenient for home use. They are noninvasive, meaning no part of the monitor enters your body, and provide systolic, diastolic and pulse measurements. Healthy Living[®] monitors are recommended for use by people over the age of 18 and are not suitable for clinical use.

Measurement Method


Healthy Living[®] monitors measure blood pressure and pulse by using what’s called an *oscillometric* method, meaning they measure the fluctuations in pressure. This is done by first wrapping a fitted cuff around the upper arm. The monitor then automatically fills the cuff with air creating pressure around the arteries inside of the arm. Mounted inside of the cuff is a special gauge that is capable of sensing the small *oscillations* (fluctuations) in pressure. These fluctuations are produced as a result of the pressure the cuff has placed around the arm and are, in a basic sense, the arteries contracting (getting larger and smaller) with each heart beat. The monitor then measures how high and how quickly these contractions occur and converts that information into a digital value. A special valve, called an *electrode*, deflates the monitor automatically when the measurement is complete.

**Main Unit and
Accessories**

Main Unit

Accessories

4 AA Alkaline Batteries
Storage Pouch

Special Features

Heart Sense®

The Heart Sense® Inflation System senses your blood pressure before the measurement is complete - the monitor then adjusts the air pressure amount according to your individual level. The results are quick, comfortable measurements each time you check your blood pressure.

Noise Control System

A unique low-noise, rolling motor pump allows our monitors to offer quiet inflation.

Large Three Row Display

A large LCD with three rows displaying time, date, and your blood pressure - systolic over diastolic - makes measuring your blood pressure easy to read.

Memory Feature

Your Healthy Living® Deluxe Inflate monitor features a 3 party memory capability or 3 "zones" - each zone stores 30 measurements each for a total capacity of 90 measurements. This unique feature allows 3 users to use and store their blood pressure measurements separately all on one monitor! This capability is also useful for one individual to track their blood pressure throughout the day - morning, afternoon, and night - for quick and easy comparisons.

Time and Date Indications

This monitor automatically keeps track of the time, month and day.


Accuracy

Healthy Living® monitors have been clinically tested against a scientific device called a *sphygmomanometer*, considered the gold standard in blood pressure measurement. All Healthy Living® monitors have performed equivalent to measurements taken with this scientific device and are within the accuracy limits prescribed by the American National Standard for Electronic or Automated Sphygmomanometers.

Special Conditions

All Healthy Living® blood pressure monitors are safe and accurate to use with arrhythmias such as atrial or ventricular premature beats and atrial fibrillation.

Digital Display Symbols


Battery Installation

When the display shows  , it's time to replace the batteries. This monitor uses 4 AA batteries.

Insert the batteries into the monitor making certain to align the polarities correctly.

Close the battery cover securely. It will snap into place when it is fully closed.


* It's a good idea to remove the batteries if you know you will not be using your monitor for an extended amount of time.

AC Capability and Adapter

This Healthy Living® blood pressure monitor is equipped with AC capability and AC adapter should you choose to use electric power to operate your unit instead of batteries.

AC Adapter Output: 6V DC 1A

*It is suggested that you remove the batteries from the compartment if you are going to use AC adapter for an extended amount of time.

Selecting Time and Date

Your Healthy Living® monitor will keep track of the time, month, and day unless the batteries have run down or are removed from the monitor.


To set the month simply press the "MODE" button. The month will begin to flash. Press the "SET" button until the correct month appears on the display, then press the "MODE" button to store, follow this same process to set the day and time.

Important Reminders

Follow your doctor's advice on how and when to take your blood pressure. Otherwise, here are some important hints:

Before Taking a Measurement

- ✓ Sit quietly for about 10 minutes. This will allow your body to return to its normal, resting state.
- ✓ Take your blood pressure in a comfortable environment - not too cold, not too hot-as temperature may affect your measurement.
- ✓ If you've just consumed beverages containing caffeine such as coffee, non-herbal tea, or cola wait about 30-45 minutes. Also, don't smoke just before taking a measurement.


Cuff Plug

- ✓ Plug the cuff tube into the monitor before taking a measurement.

Important Reminders


While Taking a Measurement

- ✓ Sit still and quietly while measuring. Talking or moving may elevate measurements.
- ✓ For consistency, it's a good idea to measure your blood pressure using the same arm and at about the same time each day. If possible, we suggest using your left arm.
- ✓ When taking multiple measurements right after each other, make sure you wait at least 10 minutes in between. Waiting allows your blood vessels to return to their normal state.
- ✓ Sit with your legs uncrossed and your feet flat on the floor. Do not touch the cuff or monitor at any time during the measurement. Relax.

Applying the Cuff

- Before you put the cuff on your arm make certain that you are not wearing a thick coat or sweater.
- If you are wearing a long sleeved shirt, roll up the sleeve, but make sure that it is not constricting circulation to your arm.
- The cuff should be positioned on your bicep and the red arterial marker should be in the center just above your elbow. Instructions are also printed directly on the cuff for reference.
- Wrap the cuff around your upper arm (usually the left) so it fits securely but is not too tight. You should be able to place 2 fingers in between the cuff and your arm. The velcro tab should fall within the "OK" range printed on the cuff for proper use.
- Place your arm on a table with your palm facing upwards. The inflation tube should line up with your arm. Do not allow the tube to kink or wrap around anything; it should fall down the center of your arm.

Taking a Blood Pressure Measurement

- Once the cuff is properly applied and you are seated with your feet flat on the floor you are ready to take a measurement.

- Press the "SET" button and you will see a number appear: 1, 2, or 3. These are each of the 3 memory zones/users explained on page 11. From here on, the 3 party memory capability will be referred to as "memory zones" or "zones". Press the "SET" button to toggle between each of three memory zones and to select the zone in which you would like to store your measurement. Then press the "START/STOP" button to begin the measurement process. The monitor will automatically inflate to the level that is right for you.

*If you select zone No. 1, your blood pressure will be recorded as one of the 30 readings available for that zone. The same is true for zones 2 and 3. When taking measurements, be sure and select the correct zone number so that the measurement is registered to the appropriate zone.

- You will feel the cuff get tighter around your arm - this is normal. It should not inflate to the point of discomfort. If this occurs, the cuff is secured too tightly. Press "START/STOP" and rewrap the cuff leaving more room between your arm and the cuff.

Taking a Blood Pressure Measurement

- Numbers will appear on the digital display. This is the normal process the monitor goes through to determine your blood pressure. Don't be alarmed if you notice the numbers quickly change, or go higher than you believe your blood pressure to be. 
- When your measurement is complete the monitor will deflate automatically, releasing the air pressure in the cuff. 
- The entire measurement process should take less than a minute. Once measurement is complete, your blood pressure and pulse measurements will be displayed. The monitor will automatically turn off after approximately one minute. 

NOTE: Blood pressure constantly changes during the course of a day. Occasionally high or low readings do not necessarily mean your monitor is malfunctioning or that you are ill. It is important to speak with your doctor to learn what your normal blood pressure is and if you experience related abnormal measurements.

3 Zones of 30 Memory

When you take a measurement, the monitor automatically stores it in memory space number one. As you continue to take measurements, the monitor will push the previous measurement back one space. In other words, the value in space one will always be the most recent blood pressure measurement taken.

Your monitor features a 3 party memory capability or 3 "zones" - each zone stores 30 measurements each for a total capacity of 90 measurements. This unique feature allows 3 users to use and store their blood pressure measurements separately all on one monitor. This feature is also useful for one user to track their blood pressure throughout the various periods of the day: morning, afternoon, and night.

Press the "SET" button to toggle between memory zone No.1, 2, and 3. Once you select a zone that you want to view the measurements in, press the "MEMORY" button once to view the most recent measurement along with the time and date that it was taken. Press the "MEMORY" button again in succession and you will see all of the measurements registered to that zone. You can then press the "SET" button to switch to the next zone to view measurements. The monitor will automatically power off after about one minute of inactivity.

***When taking measurements, be sure and select the correct zone number using the SET button so that the measurement is registered to the appropriate zone.**

Taking Care of Your Blood Pressure Monitor

It's important to treat your monitor with care. Here are some basic tips:

- ✓ Keep the monitor out of direct sunlight, high humidity, and extreme temperature changes. 
- ✓ Try not to drop it on the floor.
- ✓ Do not take the monitor or the cuff apart. None of the monitor's parts are fixable at home. If your monitor needs repairs see details in warranty section.
- ✓ Keep the monitor free of dust and other materials. Avoid getting any part of the monitor or cuff extremely wet. If the monitor needs cleaning it's best to wipe it with a cloth slightly dampened with water. Healthy Living® recommends cleaning your monitor about once a month.
- ✓ Do not harshly fold or bend the cuff as this may damage it on the inside. 
- ✓ Do not inflate the monitor's cuff if it is not wrapped around your arm.
- ✓ Always remove the cuff tube from the monitor with care.
- ✓ For your safety, do not wrap the cuff around any other body part other than your arm.

Troubleshooting

Error/Symbol	What Does It Mean?	How To Fix It
ERROR	Monitor may have been inflated without the cuff around the arm <i>OR</i> Cuff may not be positioned correctly on the arm.	Read page 17 (Applying the Cuff).
	The cuff tube may be folded <i>OR</i> You may have moved while the cuff was inflating.	Read pages 15 & 16 (Important Reminders).
	You may have moved or spoken during the measurement process.	Remember that it is important not to move and speak as it may elevate measurements.
	You may have been startled by a noise or vibration while measuring.	Make sure you are in a quiet and comfortable place when taking a measurement.

Troubleshooting

Error/Symbol	What Does It Mean?	How To Fix It
ERROR	You have been ignoring the flashing replace battery symbol ().	Replace with four new "AA" batteries. Make sure the polarities are positioned correctly in the battery compartment before continuing to measure your blood pressure.
	The batteries have run down.	Replace with 4 AA batteries quickly so you may continue to measure your blood pressure without any interruptions.
Blood pressure value seems too high or too low	Cuff may not be aligned correctly	Read page 17 (Applying the Cuff).
Blood pressure value seems too high or too low	Circuitry locked.	Remove then reinsert the batteries and then proceed to take measurement again.
Unable to take a measurement after receiving an error symbol.	The cuff may not be positioned correctly.	Read page 17 (Applying the Cuff).

Five Year Limited Warranty

Samsung America, Inc. will, at its option, repair or replace your monitor, including the cuff, with the same or comparable model free of charge (except you must pay for shipping charges set forth below) for a period of (5) five years from the date of original purchase in the event of a defect in materials or workmanship occurring during normal use. This product is for home use only and is not intended for commercial or clinical use. This warranty only extends to the original retail purchaser. Purchase receipt or other proof of date of original purchase is required before warranty performance. The warranty on any replacement product provided under this warranty shall be for the unexpired portion of the warranty period applicable to the original product.

This warranty extends solely to failures due to defects in materials or workmanship occurring during normal use. It does not cover normal wear of the product. This warranty is void if: the product housing has been removed, if the product's label, logo or serial number have been removed or if the product fails to function properly as a result of accident, misuse, abuse, neglect, mishandling, misapplication, defective batteries, faulty installation, setup, adjustments, improper maintenance, alteration, maladjustment of controls, modification, power surges, commercial or clinical use of the product, service by anyone other than an authorized service center or acts beyond the manufacturer's or distributor's control.

The warranty and remedy provided herein are exclusive and in lieu of all other express and implied warranties and unless stated herein, any statements or representations made by any other person or entity are void. NEITHER THE MANUFACTURER, DISTRIBUTOR NOR THEIR AFFILIATES SHALL BE LIABLE FOR ANY LOSS, INCONVENIENCE, INJURY, OR DAMAGE INCLUDING DIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES RESULTING FROM THE USE, IMPROPER USE, FAILURE TO HEED ANY WARNING, DISCLAIMER

Five Year Limited Warranty

OR INSTRUCTION OR INABILITY TO USE THIS PRODUCT, WHETHER ARISING OUT OF ANY BREACH OF THIS WARRANTY OR ANY OTHER LEGAL THEORY. ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING THE WARRANTIES OF MERCHANTABILITY, ARE LIMITED TO THE APPLICABLE WARRANTY PERIOD SET FORTH ABOVE. Manufacturer and Distributor assume no responsibility for injuries, damages or penalties incurred or suffered resulting from the use of this product in an improper manner or manner or location other than for that for which it is intended or from failure to heed any warning, disclaimer or instruction.

Some states do not allow the exclusion or limitation of incidental or consequential damages, or limitations on how long an implied warranty lasts, so the above exclusions or limitations may not apply to you. This warranty gives your specific legal rights and you may also have other rights which vary from state to state.

Should your product prove defective during this warranty period, please return the product in the original carton prepaid to:

SAMSUNG AMERICA, INC.
c/o Healthy Living Division
105 Challenger Road, 6th Floor
Ridgefield Park, NJ 07660
Attn: Customer Service

Include \$11.95 for return shipping & insurance. (CHECK OR MONEY ORDER ONLY PLEASE, NON-REFUNDABLE). Sorry, we cannot deliver to P.O. Boxes. Please allow 6-8 weeks for delivery.

Include with your monitor a letter stating: your name, street address, daytime phone number, detailed description of the problem, a copy of your dated cash register receipt or other proof of purchase.

Specifications

Model	BD-3000S
Monitor	Deluxe inflate
Measurement Method	Oscillometric
Display	Liquid crystal digital display
Measurement Range	Pressure: 0-260 mm Hg Pulse: 40-199 beats/Minute
Accuracy	Pressure: within ± 3mm Hg Pulse: within ± 5%
Memory	Stores up to 30 measurements for each zone
Automatic Power Off	About 1 minute after measurement is complete
Battery Type	4 AA batteries (Alkaline)
Battery Life	About 6 months when used once a day, less when using a large cuff
Storage Temperature	Temperature: -4°F to 140°F Relative Humidity: less than 95%
Operating Temperature	Temperature: 50°F to 104°F Relative Humidity: less than 95%
Dimensions (approx.)	Monitor: 4 1/4" (l) x 2 1/2" (w) x 5 3/4" (h) Cuff: fits arms 9"-13" around
Weight (approx.)	14 oz (including batteries)

For purposes of improvement, specifications are subject to change without notice.

Blood Pressure Diary

Date	Time	Systolic/Diastolic	Pulse
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	

Blood Pressure Diary

Date Time Systolic/Diastolic Pulse

Date	Time	Systolic/Diastolic	Pulse
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	

Blood Pressure Diary

Date Time Systolic/Diastolic Pulse

Date	Time	Systolic/Diastolic	Pulse
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	

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