Owner's Manual
for model:
BW-325S

SAMSUNG AMERICA

Made in Taiwan

105 Challenger Road, 6th Floor
Ridgefield Park, NJ 07660

[www. SAMUNGHHealthyLiving.com](http://www.SAMUNGHHealthyLiving.com)

WRIST

Wrist Blood Pressure Monitor with Heart Sense®

- Large digital display panel
- Measures pulse
- Memory stores 48 measurements
- Automatic power off

Important Warnings/Medical Disclaimer

Before using your blood pressure monitor, you should read and understand all instructions and follow all warnings.

Information in this manual is provided for informational purposes and this manual and product are not meant to be a substitute for the advice provided by your own physician or other medical professional. You should not use the information contained herein or this product for diagnosing or treating a health problem, disease, or prescribing any medication. If you have or suspect that you have a medical problem, promptly contact your healthcare provider.

In This Manual

Healthy Living® _____	3-4
Commonly Asked Questions About Blood Pressure _____	5-8
Intended Use of Healthy Living® Blood Pressure Monitors _____	8
Measurement Method _____	9
Main Unit _____	10
Special Features _____	11
Accuracy _____	12
Special Conditions _____	12
Digital Display Symbols _____	13
Battery Installation _____	13
Selecting Time and Date _____	14
Important Reminders _____	15-16
Applying the Cuff _____	17
Taking a Blood Pressure Measurement _____	18-19
Forty-Eight Measurement Memory _____	19
Taking Care of Your Blood Pressure Monitor _____	20
Troubleshooting _____	21-22
Five Year Limited Warranty _____	23-24
Specifications _____	25
Blood Pressure Diary _____	26-30

Healthy Living®

This blood pressure monitor is the first step to living a healthier life. But monitoring your blood pressure at home is just the beginning, it's also important to:

STAY ACTIVE Start out by taking just a short walk each day. Over time you will begin to feel more energy. Walking not only helps you lose weight, it also builds stronger bones and a healthy heart!

EAT HEALTHY Adding more fruits and vegetables to your diet is a great way to start. Studies have shown that a diet high in these foods, as well as whole grains, fish and low-fat dairy products may be the best medication of all. It's also important to keep a check on your sodium intake as about one-third of people who have high blood pressure are also sensitive to this mineral. You may want to explore the potential benefits of incorporating certain botanicals into your diet such as onion and garlic, which have long been linked to having mild anti-hypertension effects. As always, consult your physician before beginning any dietary changes. Remember, just like you've always heard, you are what you eat!

Healthy Living®

STAY POSITIVE A healthy mind goes hand in hand with a healthy body. Daily relaxation techniques such as meditation and gentle yoga may be just what you need to alleviate unwanted stress from your life. Keep your mind focused on the positive things in life. Start out by giving yourself credit for those short walks or that apple you just ate. It's your first step.

WATCH YOUR WEIGHT! Sometimes, just a small amount of weight loss is all it takes to help control high blood pressure. If you've already started eating healthy and staying active, your weight just might take care of itself!

Remember to consult your physician prior to beginning an exercise program, vitamin/mineral, dietary changes or programs and blood pressure monitoring.

Commonly Asked Questions About Blood Pressure


Why Should I Monitor My Blood Pressure At Home?

By monitoring at home, you can manage your health in between doctor visits. Home monitoring also provides you with a record of blood pressure measurements over time to give to your doctor - a valuable tool in hypertensives' monitoring.


Home monitoring may also help alleviate *White Coat Hypertension*. This term refers to the medical phenomena which causes many individuals' blood pressure to be higher when taken in a medical office setting. That's because many times the stress of being in the doctor's office may cause blood pressure to rise as much as 20-50 points higher than normal. When monitoring at home, you may get a truer indication of your blood pressure in a normal, relaxed state.

What Is Blood Pressure?

The commonly used term "blood pressure" refers to the force or *pressure* that is exerted on the body's arteries as blood flows through them. Each time the heart beats or *contracts*, it produces pressure in the arteries. This blood pressure is what moves blood through the body, supplying oxygen and nutrients to every organ. (*con't*)


Commonly Asked Questions About Blood Pressure

When the heart contracts it sends blood out into the body. This pressure is referred to as **systolic**.

The pressure that is created as the heart pulls blood back through it is called **diastolic**.


Blood pressure is recorded as systolic over diastolic, as shown here.

What Is Considered Normal Blood Pressure?

The Seventh Report of the Joint National Committee on Prevention, Detection, Evaluation, and Treatment of High Blood Pressure from the National Heart, Lung, and Blood Institute has set forth the following guidelines regarding high blood pressure. Guidelines provided are for information purposes only. Please consult with your physician for proper diagnosis.

Commonly Asked Questions About Blood Pressure

Why Does My Blood Pressure Fluctuate?

Because blood pressure changes with every beat of the heart it is in constant fluctuation throughout any given day. In addition to these natural changes there are other factors that may affect blood pressure such as:

- Temperature (too hot or cold)
- Humidity
- Restlessness
- Lack of sleep
- Physical exercise
- Constipation
- Stress
- Fear
- Anger
- Anxiety
- Food and beverage consumption
- Individual physical conditions

In addition to these factors you should talk with your physician to see what may be affecting you.

How Do I Know I Have High Blood Pressure?

It's nearly impossible to know if your blood pressure is elevated without having it checked. People have the common misconception that one can "feel" their blood pressure rise. It is because of

Commonly Asked Questions About Blood Pressure

this mystery that hypertension has been named “the silent killer.” Don’t let it sneak up on you! Start monitoring your blood pressure every day with your new Healthy Living[®] monitor.

Why Is It Important To Control My Blood Pressure?

Left uncontrolled, hypertension may place an incredible amount of stress on your heart and arteries, both of which become forced to work harder to keep blood flowing in a healthy manner throughout the body. Over time, this stress may result in health problems including over enlargement of the heart, loss of elasticity of the arteries, heart disease, stroke and even death. Making some simple lifestyle changes now may help you to stay healthy in the future.

Intended Use of Healthy Living[®] Blood Pressure Monitors


Healthy Living[®] blood pressure monitors offer convenient, home use. They are noninvasive, meaning no part of the monitor enters your body, and provide systolic, diastolic and pulse measurements. Healthy Living[®] monitors are recommended for use by people over the age of 18 and are not suitable for clinical use.

Measurement Method


Healthy Living[®] monitors measure blood pressure and pulse by using what’s called an *oscillometric* method, meaning they measure the fluctuations in pressure. This is done by first wrapping a fitted cuff around the upper arm (or wrist, depending on the model). The monitor then automatically fills the cuff with air creating pressure around the arteries inside of the arm (or wrist). Mounted inside of the cuff is a special gauge that is capable of sensing the small *oscillations* (fluctuations) in pressure. These fluctuations are produced as a result of the pressure the cuff has placed around the arm and are, in a basic sense, the arteries contracting (getting larger and smaller) with each heart beat. The monitor then measures how high and how quickly these contractions occur and converts that information into a digital value. A special valve, called an *electrode*, deflates the monitor automatically when the measurement is complete.

Main Unit

Main Unit

Accessories

2 AAA Batteries Alkaline Storage Case

Special Features

Heart Sense®

The Heart Sense® Inflation system senses your blood pressure before the measurement is complete - the monitor then adjusts the air pressure amount according to your individual level. The results are quick, comfortable measurements each time you check your blood pressure.

Noise Control System

A unique and Healthy Living® exclusive low-noise, rolling motor pump allows our monitors to offer quiet inflation.

Memory Feature

Healthy Living® special memory feature allows the wrist monitor to store up to 48 measurements along with the time and date for quick and easy comparisons.

Time and Date

This monitor automatically keeps track of the time, month and day.


Accuracy

Healthy Living® monitors have been clinically tested against a scientific device called a *sphygmomanometer*, considered the gold standard in blood pressure measurement. All Healthy Living® monitors have performed equivalent to measurements taken with this scientific device and are within the accuracy limits prescribed by the American National Standard for Electronic or Automated Sphygmomanometers.

Special Conditions

All Healthy Living® blood pressure monitors are safe and accurate to use with arrhythmias such as atrial or ventricular premature beats and atrial fibrillation.

Digital Display Symbols


Battery Installation

When the display shows , it's time to replace the batteries. This monitor uses 2 AAA batteries.

Insert the batteries into the monitor making certain to align the polarities correctly.

Close the battery cover securely. It will snap into place when it is fully closed.

* If you know you will not be using your monitor for an extended amount of time it's a good idea to remove the batteries.

Selecting the Time and Date

Your Healthy Living® monitor will keep track of the time, month, and day unless the batteries have run down or are removed from the monitor.

To set the time simply press the "MODE" button. The month will begin to flash. Press the "SET" button until the correct month appears on the display, then press the "MODE" button to store. Follow this same process to set the day and time.

Important Reminders

Follow your doctor's advice on how and when to take your blood pressure. Here are some important hints:

Before Taking A Measurement

- ✓ Sit quietly for about 10 minutes. This will allow your body to return to its normal, resting state.
- ✓ Take your blood pressure in a comfortable environment - not too cold, not too hot as temperature may affect your measurement.
- ✓ If you've just consumed beverages containing caffeine such as coffee, non-herbal tea, or cola wait about 30-45 minutes. Also, don't smoke just before taking a measurement.

While Taking A Measurement

- ✓ Sit still and quietly while measuring. Talking or moving may elevate measurements.
- ✓ For consistency, it's a good idea to measure your blood pressure using the same wrist and at about the same time each day. If possible, we suggest using your left wrist.

Important Reminders

While Taking A Measurement, con't


- ✓ When taking multiple measurements right after each other, make sure you wait at least 10 minutes in between. Waiting allows your blood vessels to return to their normal state.
- ✓ Sit with your legs uncrossed and your feet flat on the floor. Do not touch the cuff or monitor at any time during the measurement. Relax.


NOTE: Blood pressure constantly changes during the course of a day. Occasionally, high or low readings do not necessarily mean your monitor is malfunctioning or that you are ill. It is important to speak with your doctor to learn what your normal blood pressure is and if you experience related abnormal measurements.

Applying The Cuff

- Before you put the cuff on your wrist make certain that you are not wearing a thick coat or sweater.
- If you are wearing a long sleeved shirt, roll up the sleeve, but make sure that it is not constricting circulation to your arm.
- Wrap the cuff around your left wrist making sure that it is secure but is not too tight.


- Place your arm on a table with your palm facing upwards; the monitor should be facing you. It is important to make sure that the wrist cuff is at the same level as your heart. You can rest your wrist on the closed storage case provided with your monitor.


Taking a Blood Pressure Measurement

- Once the cuff is properly applied and you are seated with your feet flat on the floor you are ready to take a measurement.
- To take a measurement, press "START/STOP" button. The monitor will automatically inflate to 195mmHg. If the monitor senses you require a higher inflation level, it will automatically increase the air pressure to the level that's right for you.


- You will feel the cuff get tighter around your wrist - this is normal. It should not inflate to the point of discomfort. If this occurs, the cuff is secured too tightly. Press the "START/STOP" button to stop inflation and rewrap the cuff leaving more room between your wrist and the cuff.
- Numbers will appear on the digital display. This is the normal process the monitor goes through to determine your blood pressure. Don't be alarmed if you notice the numbers quickly change, or go higher than you believe your blood pressure to be.

Taking a Blood Pressure Measurement

- When your measurement is complete, the monitor will deflate automatically, releasing the air pressure in the cuff.
- The entire measurement process should take less than a minute. Once measurement is complete your blood pressure and pulse measurements will be displayed.

Forty-Eight Measurement Memory

When you take a measurement the monitor automatically stores it in memory space number one. As you continue to take measurements, the monitor will push the previous measurement back one space. In other words, the value in space one will always be the most recent blood pressure measurement taken. Your monitor will store up to 48 measurements. To review measurement in the memory, depress the "MEMORY" button. Continue pressing "MEMORY" to review all 48 of your last blood pressure measurements.

Taking Care of Your Blood Pressure Monitor

It's important to treat your monitor with care. Here are some common sense tips:

- ✓ Keep the monitor out of direct sunlight, high humidity, and extreme temperature changes.   
- ✓ Try not to drop it on the floor.
- ✓ Do not take the monitor or the cuff apart. None of the monitor's parts are fixable at home. If your monitor needs repairs see details in warranty section.
- ✓ Keep the monitor free of dust and other materials. Avoid getting any part of the monitor or cuff extremely wet. If the monitor needs cleaning it's best to wipe it with a cloth slightly dampened with water. Healthy Living® recommends cleaning your monitor about once a month.
- ✓ Do not harshly fold or bend the cuff as this may damage it on the inside.
- ✓ Do not inflate the monitor's cuff if it is not wrapped around your wrist.
- ✓ For your safety, do not wrap the cuff around any other body part other than your wrist.

Troubleshooting

Error/Symbol	What Does It Mean?	How To Fix It
ERROR	Monitor may have been inflated without the cuff around the wrist <i>Or</i> Cuff may not be positioned correctly on the wrist.	Read page 17 (Applying The Cuff).
	You may have moved while the cuff was inflating, it is now attempting to reinflate. Check the cuff to make sure it's positioned correctly.	Read pages 15 & 16 (Important Reminders).
	You may have moved or spoken during the measurement process.	Remember that it is important to not move or speak while taking a measurement.
	You may have been startled by a noise or vibration while measuring.	Make sure you are in a quiet and comfortable place when taking a measurement.

Troubleshooting

Error/Symbol	What Does It Mean?	How To Fix It
	You have been ignoring the flashing replace battery symbol ().	Replace with 2 new "AAA" batteries. Make sure the polarities are positioned correctly in the battery compartment before continuing to measure your blood pressure.
	The batteries have run down.	Replace with 2 "AA" batteries quickly so you may continue to measure your blood pressure without any interruptions.
Blood pressure value seems too high or too low	Cuff may be positioned too high or too low on the wrist.	Read page 17 (Applying The Cuff).
The monitor keeps reinflating	You may have moved during inflation <i>Or</i> If you have a very weak pulse the monitor may continue to reinflate until it senses your pulse.	Relax for a few minutes then take a measurement.
Unable to take a measurement after receiving an error symbol.	The cuff may not be positioned correctly.	Read page 17 (Applying The Cuff).

Five Year Limited Warranty

Samsung America, Inc. will, at its option, repair or replace your monitor, including the cuff, with the same or comparable model free of charge (except you must pay for shipping charges set forth below) for a period of (5) five years from the date of original purchase in the event of a defect in materials or workmanship occurring during normal use. This product is for home use only and is not intended for commercial or clinical use. This warranty only extends to the original retail purchaser. Purchase receipt or other proof of date of original purchase is required before warranty performance. The warranty on any replacement product provided under this warranty shall be for the unexpired portion of the warranty period applicable to the original product.

This warranty extends solely to failures due to defects in materials or workmanship occurring during normal use. It does not cover normal wear of the product. This warranty is void if: the product housing has been removed, if the product's label, logo or serial number have been removed or if the product fails to function properly as a result of accident, misuse, abuse, neglect, mishandling, misapplication, defective batteries, faulty installation, setup, adjustments, improper maintenance, alteration, maladjustment of controls, modification, power surges, commercial or clinical use of the product, service by anyone other than an authorized service center or acts beyond the manufacturer's or distributor's control.

The warranty and remedy provided herein are exclusive and in lieu of all other express and implied warranties and unless stated herein, any statements or representations made by any other person or entity are void. NEITHER THE MANUFACTURER, DISTRIBUTOR NOR THEIR AFFILIATES SHALL BE LIABLE FOR ANY LOSS, INCONVENIENCE, INJURY, OR DAMAGE INCLUDING DIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES RESULTING FROM THE USE, IMPROPER USE, FAILURE TO HEED ANY WARNING, DISCLAIMER

Five Year Limited Warranty

OR INSTRUCTION OR INABILITY TO USE THIS PRODUCT, WHETHER ARISING OUT OF ANY BREACH OF THIS WARRANTY OR ANY OTHER LEGAL THEORY. ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING THE WARRANTIES OF MERCHANTABILITY, ARE LIMITED TO THE APPLICABLE WARRANTY PERIOD SET FORTH ABOVE. Manufacturer and Distributor assume no responsibility for injuries, damages or penalties incurred or suffered resulting from the use of this product in an improper manner or manner or location other than for that for which it is intended or from failure to heed any warning, disclaimer or instruction.

Some states do not allow the exclusion or limitation of incidental or consequential damages, or limitations on how long an implied warranty lasts, so the above exclusions or limitations may not apply to you. This warranty gives you specific legal rights and you may also have other rights which vary from state to state. This warranty is only valid in the United States.

Should your product prove defective during this warranty period, please return the product in the original carton prepaid to:

SAMSUNG AMERICA, INC.
c/o Healthy Living Division
105 Challenger Road, 6th Floor
Ridgefield Park, NJ 07660
Attn: Customer Service

Include \$11.95 for return shipping & insurance. (CHECK OR MONEY ORDER ONLY PLEASE, NON-REFUNDABLE). Sorry, we cannot deliver to P.O. Boxes or outside the U.S. Please allow 6-8 weeks for delivery.

Include with your monitor a letter stating: your name, street address, daytime phone number, detailed description of the problem, a copy of your dated cash register receipt or other proof of purchase.

Specifications

Model	BW-325S
Monitor Type	Wrist cuff
Measurement Method	Oscillometric
Display	Liquid crystal digital display
Measurement Range	Pressure: 0-260 mm Hg Pulse: 40-199 beats/Minute
Accuracy	Pressure: within \pm 3mm Hg Pulse: within \pm 5%
Memory	Stores up to 48 measurements
Automatic Power Off	About 2 minutes after measurement is complete
Battery Type	2 AAA batteries (Alkaline)
Battery Life	About 6 months when used once a day
Storage Temperature	Temperature: -4°F to 140°F Relative Humidity: less than 95%
Operating Temperature	Temperature: 50°F to 104°F Relative Humidity: less than 95%
Dimensions (approx.)	Monitor: 3/4" (L) x 2 1/2" (W) x 1 3/8" (H) Cuff: fits wrists 5" - 8" around
Weight (approx.)	5 oz. (including batteries)

For purposes of improvement, specifications are subject to change without notice.

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