



Cytokine Secretion Assays

Miltenyi Biotec GmbH
Friedrich-Ebert-Str. 68
51429 Bergisch Gladbach, Germany
Phone: +49 2204 83060
Fax: +49 2204 85197
macs@miltenyibiotec.de

Miltenyi Biotec Inc.
2303 Lindbergh Street,
Auburn CA 95602, USA
Phone: 800 FOR MACS, +1 530 888 8871
Fax: +1 530 888 8925
macs@miltenyibiotec.com

Miltenyi Biotec Pty. Ltd. (Australia)
Phone: +61 2 8877 7400
Fax: +61 2 9889 5044
macs@miltenyibiotec.com.au

Miltenyi Biotec Shanghai Office
Phone: +86 21 62351005
Fax: +86 21 62350953
miltenyibiotec@china.com

Miltenyi Biotec (France)
Phone: +33 1 56 98 16 16
Fax: +33 1 56 98 16 17
macs@miltenyibiotec.fr

For further information refer to our website www.miltenyibiotec.com

For technical questions, please contact your local distributor or our
Technical Support Team in Germany:

e-mail: macsTec@miltenyibiotec.de
phone: +49-2204-830 6 830.

Miltenyi Biotec S.r.l. (Italy)
Phone: +39 51 64 60 411
Fax: +39 51 64 60 499
macs@miltenyibiotec.it

Miltenyi Biotec K.K. (Japan)
Phone: +81 3 56 46 8910
Fax: +81 3 56 46 89 11
macs@miltenyibiotec.jp

**Miltenyi Biotec Asia Pacific Pte. Ltd.
(Singapore)**
Phone: +65 6238 8183
Fax: +65 6238 0302
macs@miltenyibiotec.com.sg

Miltenyi Biotec S.L. (Spain)
Phone: +34 91 512 12 90
Fax: +34 91 512 12 91
macs@miltenyibiotec.es

Miltenyi Biotec Ltd. (UK)
Phone: +44 1483 799 800
Fax: +44 1483 799 811
macs@miltenyibiotec.co.uk

Miltenyi Biotec

This MACS® product is for *in vitro* research use only and not for diagnostic or therapeutic procedures.




Cytokine Secretion Assays

IFN- γ Secretion Assay Cell Enrichment and Detection Kit (PE) human

For 50 tests with 10^7 cells

Order No. 130-054-201



Miltenyi Biotec



Index

1. Description

Index


1. Description
 - 1.1 Principle of the IFN- γ Secretion Assay
 - 1.2 Background and product applications
 - 1.3 Reagent and instrument requirements
2. Protocol overview
3. Experimental set-up
 - 3.1 Controls
 - 3.2 Kinetics of restimulation and proposed time schedule
 - 3.3 Counterstaining of cytokine secreting cells
 - 3.4 Two color cytokine analysis
 - 3.5 Combination with peptide-MHC tetramer staining
 - 3.6 Detection without prior enrichment
4. Protocol for the IFN- γ Secretion Assay
 - 4.1 Cell preparation
 - 4.2 (Antigen-specific) In vitro stimulation
 - 4.3 Cytokine Secretion Assay
 - 4.4 Magnetic labeling
 - 4.5 Magnetic separation
5. Detection and analysis of IFN- γ secreting antigen-specific T cells
6. References
7. Appendix
 - A: Flask and dish sizes for stimulation
 - B: Detection and enrichment of cytokine secreting cells from whole blood

1. Description

Components	<p>1 mL IFN-γ Catch Reagent: anti-IFN-γ monoclonal antibody (mouse IgG1) conjugated to cell surface specific monoclonal antibody (mouse IgG2a).</p> <p>1 mL IFN-γ Detection Antibody: anti-IFN-γ monoclonal antibody (mouse IgG1) conjugated to PE (R-phycoerythrin).</p> <p>1 mL Anti-PE MicroBeads: colloidal superparamagnetic MicroBeads conjugated to monoclonal mouse anti-PE antibody (mouse IgG1).</p>
Size	For 50 tests with 10^7 cells
Product format	All components are supplied as a suspension containing stabilizer and 0.05% sodium azide.
Storage	Store protected from light at 4–8 °C. Do not freeze. The expiration dates are indicated on the vial labels.

Miltenyi Biotec





1.1 Principle of the IFN- γ Secretion Assay

Antigen-specific T cells are analyzed and isolated using the IFN- γ Secretion Assay starting from whole blood, PBMC or other leukocyte containing single-cell preparations. The cells are restimulated for a short period of time with specific peptide, protein or other antigen preparations.

Subsequently, an IFN- γ -specific **Catch Reagent** is attached to the cell surface of all leukocytes. The cells are then incubated for a short time at 37 °C to allow cytokine secretion. The secreted IFN- γ binds to the IFN- γ Catch Reagent on the positive, secreting cells. These cells are subsequently

labeled with a second IFN- γ -specific antibody, the **IFN- γ Detection Antibody** conjugated to R-phycoerythrin (PE) for sensitive detection by flow cytometry.

The IFN- γ -secreting cells can now be magnetically labeled with **Anti-PE MicroBeads** and enriched over a MACS® Column which is placed in the magnetic field of a MACS Separator. The magnetically labeled cells are retained in the MACS Column while the unlabeled cells run through. After the column has been removed from the magnetic field, the magnetically retained cells can be eluted as positively selected cell fraction, enriched for cytokine secreting cells. The cells can now be used for cell culture or analysis. Since viable cells are analyzed, non-specific background can be minimized by dead cell exclusion. This provides highest sensitivity of analysis.

1.2 Background and product applications

The IFN- γ Secretion Assay - Cell Enrichment and Detection Kit is designed for the detection, isolation and analysis of viable IFN- γ secreting leukocytes. It is specially developed for the **detection and isolation of antigen-specific T cells**. After restimulation with specific antigen *in vitro* secretion of IFN- γ is induced. IFN- γ is predominantly secreted by activated CD4⁺ and CD8⁺ memory and effector T cells and by NK cells upon activation.

Quantitative analysis of antigen-specific T cell populations can provide important information on the natural course of immune responses. MACS enrichment of the antigen-specific T cells increases the sensitivity of analysis, allowing detection of frequencies as low as one in a million cells.

Miltenyi Biotec

4

140-000-274.07

140-000-274.07

5



The MACS enrichment also enables further functional characterization of the antigen-specific cells and downstream experiments, as well as the expansion of antigen-specific cells allowing research on potential future immunotherapeutic applications.

Examples of applications

- Detection and enrichment of viable IFN- γ -secreting leukocytes.
- Detection and enrichment of IFN- γ -secreting, antigen-specific T cells for enumeration and phenotypic analysis as well as for expansion and functional characterization.
- Monitoring and analysis of antigen-specific T cell immunity, e.g. in infection, autoimmunity, cancer, allergy or alloreactivity.
- Isolation and expansion of antigen-specific T cells for research in immunotherapy.
- Enrichment and analysis of IFN- γ secreting cells for determination of functional antigens in disease and for T cell receptor (TCR) epitope mapping.
- Analysis or cloning of TCR repertoire of antigen-specific T cells.

1.3 Reagent and instrument requirements

- Buffer (degassed): Prepare a solution containing PBS (phosphate buffered saline) pH 7.2, 0.5% BSA and 2 mM EDTA by diluting MACS BSA Stock Solution (# 130-091-376) 1:20 with autoMACS™ Rinsing Solution (# 130-091-222). Keep buffer cold (4–8 °C).
- Culture medium, e.g. RPMI 1640 (# 130-091-440), containing 5% human serum, like autologous or AB serum (do not use BSA or FCS because of non-specific stimulation!).
- Propidium iodide (PI) or 7-AAD for flow cytometric exclusion of dead cells. For cell fixation and flow cytometric exclusion of dead cells, the Fixation and Dead Cell Discrimination Kit (# 130-091-163) is recommended.
- (Optional) Staining reagents such as CD4-FITC (# 130-080-501) or CD8-FITC (# 130-080-601) and CD14-PerCP™.

Miltenyi Biotec

6

140-000-274.07

140-000-274.07

7




- MACS Columns and MACS Separators:

Column	max. number of labeled cells	max. number of total cells	Separator
MS	10^7	2×10^8	MiniMACS, OctoMACS; with Column Adapter: VarioMACS, SuperMACS
LS	10^8	2×10^9	MidiMACS; with Column Adapter: VarioMACS, SuperMACS
autoMACS	2×10^8	4×10^9	autoMACS

▲ **Note:** Column adapters are required to insert certain columns into VarioMACS™ Separator or SuperMACS™ Separator. For details, see MACS Separator data sheets.

- Refrigerated centrifuge (4–8 °C).
- Rotation device for tubes: MACSmix™ tube rotator (# 130-090-753).
- (Optional) Pre-Separation Filters (# 130-041-407) to remove cell clumps.

2. Protocol overview



Miltenyi Biotec



3. Experimental set-up

3.1 Controls

Negative control

For accurate detection of IFN-γ-secreting antigen-specific cells, a negative control sample should always be included. This will provide information about IFN-γ secretion unrelated to the specific antigen-stimulation, but e.g. due to ongoing in vivo immune responses. The control sample should be treated exactly the same as the antigen-stimulated sample except for the addition of antigen, or by using a control antigen.

Positive control

When setting up a new experiment, it is recommended to include a positive control. As a positive control, a sample stimulated with the superantigen Staphylococcal Enterotoxin B (Sigma) 1 µg/mL for 3–16 hours, may be included in the experiment.

▲ **Note:** Mitogens like PHA or PMA/Ionomycin are not recommended for stimulation of a positive control, as the resulting high frequencies of IFN-γ secreting cells do not allow conclusions on the performance (e.g. sensitivity) of the IFN-γ Secretion Assay.

3.2 Kinetics of restimulation and proposed time schedule

Peptides

Upon stimulation with peptide, the cells can be analyzed for IFN-γ secretion 3–6 hours later.

It is possible to prepare the cells first and take them into culture overnight, but without adding the antigen (see 4.2 step 2.). Peptide is then added the next morning for 3 hours of stimulation, directly followed by the IFN-γ Secretion Assay.

Proteins

Upon stimulation with protein, the cells can be analyzed for IFN-γ secretion 6–16 hours later.

It is possible to start the stimulation of the cells late in the afternoon, and to perform the IFN-γ Secretion Assay the following morning.

Costimulation

The addition of costimulatory agents like CD28 or CD49d antibody may enhance the response to the antigen. If costimulatory agents are added to the antigen sample, they also have to be included in the control sample.

3.3 Counterstaining of cytokine secreting cells

The IFN-γ secreting cells are stained with PE-conjugated IFN-γ Detection Antibodies. To identify cells of interest, counterstaining for T cells with e.g. CD4-FITC (# 130-080-501) or CD8-FITC (# 130-080-601) is important.

Miltenyi Biotec



▲ Do **not** use tandem conjugates of phycoerythrin, like Cy-Chrome® (PharMingen), PE-Cy5 (Serotec), ECD, PC5 (Coulter-Immunotech) etc., they may also be recognized by the Anti-PE MicroBeads.

▲ Upon activation of T cells, TCR and some associated molecules, like CD3, might be down-regulated.

▲ The samples should be stained with propidium iodide (PI) or 7-AAD prior to acquisition, to exclude dead cells from analysis. This will reduce non-specific background staining and increase sensitivity.

▲ For optimal sensitivity, we recommend labeling of undesired non-T cells such as monocytes with antibodies conjugated to PerCP™, e.g. CD14-PerCP™. These cells can then be excluded together with PI stained dead cells by gating.

3.4 Two color cytokine analysis

IFN-γ-secreting-cells can be analyzed simultaneously for IL-2 or IL-10 production by two color cytokine analysis combining the IFN-γ Secretion Assay with the IL-2 Secretion Assay - Detection Kit (APC) (# 130-090-763), or the IL-10 Secretion Assay - Detection Kit (APC) (# 130-090-761). Detailed protocols are included in the data sheets of the Cytokine Secretion Assay - Detection Kits (APC) and are available from our website www.miltenyibiotec.com.

Miltenyi Biotec

12

140-000-274.07

3.5 Combination with peptide-MHC tetramer staining

IFN-γ-secreting cells can be analyzed simultaneously for peptide-MHC tetramers by combining the IFN-γ Secretion Assay (PE) with APC-conjugated peptide-MHC tetramers. For combination with PE-conjugated peptide-MHC tetramers the IFN-γ Secretion Assay - Detection Kit (APC) (# 130-090-762) and the IFN-γ Secretion Assay - Detection Kit (FITC) (# 130-090-433) are available. Detailed recommendations for the experimental setup and the procedure are included in the data sheets of the Cytokine Secretion Assay - Detection Kits (APC) and are available from our website www.miltenyibiotec.com.

3.6 Detection without prior enrichment

(Optional) If the sample contains more than 0.01–0.1% of IFN-γ-secreting cells, you may be able to analyze IFN-γ-secreting cells without prior enrichment (see also: IFN-γ Secretion Assay - Detection Kit (PE), # 130-054-202). The assay can also be performed directly starting from whole blood. A detailed protocol is included in the data sheet of the IFN-γ Secretion Assay - Detection Kit (PE) and is available from our website www.miltenyibiotec.com.



13

4. Protocol for the IFN-γ Secretion Assay

4.1 Cell preparation

For the detection and isolation of cytokine secreting cells, best results are achieved by starting the assay with fresh PBMC, or other leukocyte containing single-cell preparations from tissues or cell lines. Alternatively, frozen cell preparations can be used.

▲ **Note:** PBMC may be stored over night. The cells should be resuspended and incubated in culture medium as described in 4.2 step 2., but without addition of antigen. The antigen is then added to the culture on the next day.

▲ **Note:** Remove platelets after density gradient separation. Resuspend cell pellet, fill tube with buffer and mix. Centrifuge at 200xg for 10–15 minutes at 20°C. Carefully remove supernatant.

Special protocols for whole blood: You can start the IFN-γ Secretion Assay directly from whole blood. For details on the procedure, see 7. Appendix B: Detection and enrichment of cytokine secreting cells from human whole blood. This special protocol is also available from our website www.miltenyibiotec.com.

4.2 (Antigen-specific) *In vitro* stimulation

▲ Always include a **negative control** in the experiment. A positive control may also be included (see 3.1).

▲ Do **not** use media containing any **non-human** proteins, like BSA or FCS because of non-specific stimulation.

Miltenyi Biotec

14

140-000-274.07



Protocol for *in vitro* stimulation

1. Wash cells by adding medium, centrifuge at 300xg for 10 minutes.
2. Resuspend cells in culture medium, containing 5% human serum, adjust to 10⁷ cells/mL and 5×10⁶ cells/cm² (see 7. Appendix A: Flask and dish sizes for stimulation).
3. Add antigen or control reagent:
 peptide: 3–6 hours at 37 °C, 5–7% CO₂, e.g. 1–10 µg/mL
 protein: 6–16 hours at 37 °C, 5–7% CO₂, e.g. 10 µg/mL
 SEB: 3–16 hours at 37 °C, 5–7% CO₂, e.g. 1 µg/mL
 For comparison of different experiments, the stimulation time should always be the same (see 3.2).
4. Collect cells carefully by using a cell scraper, or by pipetting up and down when working with smaller volumes. Rinse the dish with cold buffer. Check microscopically for any remaining cells, if necessary, rinse the dish again.

4.3 Cytokine Secretion Assay

General considerations

▲ The assay is optimized for cell samples containing < 5% of total IFN-γ-secreting cells. If ≥ 5% of IFN-γ-secreting cells are expected, it is necessary to dilute the cells further during the cytokine secretion period, and therefore a larger test tube will be needed



15

(see table below). The dilution prevents non-specific staining of cells not secreting IFN- γ during this period.

- ▲ For each test with 10^7 total cells, prepare:

100 mL of cold buffer (4–8 °C)

100 μ L of cold medium (4–8 °C)

10 mL (or 100 mL; see table below) of warm medium (37 °C).

▲ Work fast, keep the cells cold, use pre-cooled solutions which will prevent capping of antibodies on the cell surface and a non-specific cell labeling (exception: warm medium during secretion period).

▲ Volumes shown below are for 10^7 total cells. When working with fewer than 10^7 cells, use the same volumes as indicated. When working with higher cell numbers, scale up all reagent volumes and total volumes, accordingly (e.g. for 2×10^7 total cells, use twice the volume of all indicated reagent volumes and total volumes).

▲ Do not remove supernatant by decanting. This will lead to cell loss and incorrect incubation volumes. Pipette off or aspirate supernatant.

▲ Dead cells may bind non-specifically to MACS® MicroBeads or antibodies. Therefore, when working with cell preparations containing large amounts of dead cells, they should be removed before starting the IFN- γ Secretion Assay, e.g. by density gradient centrifugation or by using the Dead Cell Removal Kit (# 130-090-101).

Miltenyi Biotec

16

140-000-274.07



Labeling cells with IFN- γ Catch Reagent

1. Use 10^7 total cells in a 15 mL closable tube per sample.

▲ **Note:** For larger cell numbers, scale up all volumes accordingly. For fewer than 10^7 cells, use same volumes.

2. Wash cells by adding 10 mL of cold buffer, centrifuge at $300 \times g$ for 10 minutes at 4–8 °C, pipette off supernatant completely.

▲ **Note:** Do not remove supernatant by decanting. This will lead to cell loss and incorrect incubation volumes.

3. Resuspend cell pellet in 80 μ L of cold medium per 10^7 total cells.

4. Add 20 μ L of IFN- γ Catch Reagent per 10^7 total cells, mix well and incubate for 5 minutes on ice.



IFN- γ secretion period

1. Add warm (37 °C) medium to dilute the cells according to the following table:

Expected number of IFN- γ secreting cells	Dilution	Amount of medium to add per 10^7 total cells
< 5%	10^6 cells/mL	10 mL
\geq 5%	$\leq 10^5$ cells/mL	100 mL

▲ **Note:** For frequencies of cytokine secreting cells \gg 20% the cells need to be further diluted, e.g. by a factor of 5.



17

2. Incubate cells in closed tube for 45 minutes at 37 °C under slow continuous rotation using the MACSmix™ tube rotator (# 130-090-753), or turn tube every 5 minutes to resuspend settled cells.

▲ **Note:** During this step it is crucial to prevent contact of cells to avoid cross contamination with cytokines.



Labeling cells with IFN- γ Detection Antibody

1. Put the tube **on ice**.
2. Wash the cells by filling up the tube with **cold buffer**, and centrifuge at $300 \times g$ for 10 minutes at 4–8 °C. Pipette off supernatant completely.
▲ **Note:** If the volume of the cell suspension was higher than the volume of added buffer, repeat wash step.
3. Resuspend cell pellet in 80 μ L of **cold buffer** per 10^7 total cells.
4. Add 20 μ L of **IFN- γ Detection Antibody (PE)** per 10^7 total cells.
5. (Optional) Add additional staining reagents, e.g. 10 μ L of CD4-FITC (# 130-080-501) or 10 μ L of CD8-FITC (# 130-080-601) and CD14-PerCP™.
6. Mix well and incubate for 10 minutes **on ice**.
7. Wash cells by adding 10 mL of **cold buffer**, centrifuge at $300 \times g$ for 10 minutes at 4–8 °C, pipette off supernatant.

Miltenyi Biotec

18

140-000-274.07

4.4 Magnetic labeling



Magnetic labeling with Anti-PE MicroBeads

1. Resuspend cell pellet in 80 μ L of cold buffer per 10^7 total cells.
2. Add 20 μ L of **Anti-PE MicroBeads** per 10^7 total cells, mix well and incubate for 15 minutes at 4–8 °C.
▲ **Note:** Incubate in refrigerator at 4–8 °C, do not work on ice during this step.
3. Wash cells by adding 10 mL of **cold buffer**, centrifuge at $300 \times g$ for 10 minutes at 4–8 °C. Pipette off supernatant.
4. Resuspend cell pellet in 500 μ L of **cold buffer**. For higher cell numbers than 5×10^7 use a dilution of 10^8 cells/mL.
5. (Optional) Take an aliquot for flow cytometric analysis and cell count of the fraction before enrichment.
6. Proceed to magnetic separation (see 4.5).

4.5 Magnetic separation



Magnetic separation using MS or LS Columns

▲ Choose an appropriate MACS® Column and MACS Separator according to the number of total cells (see table in 1.3).



19

▲ When enriching antigen-specific T cells, **always perform two consecutive column runs** to achieve best results.

1. Prepare two columns per sample by rinsing with cold buffer:
MS: 500 μ L LS Column: 3 mL
and discard effluent.
2. Place the first column into the magnetic field of a MACS Separator (use column adapter with VarioMACS or SuperMACS Separator).
3. (Optional) Pass the cells through Pre-Separation Filters (# 130-041-407) to remove clumps.
4. Apply cell suspension onto the column.
5. Collect unlabeled cells which pass through and wash with appropriate amount of cold buffer. Perform washing steps by adding buffer successively once the column reservoir is empty.
MS: 3 \times 500 μ L LS: 3 \times 3 mL
Collect total effluent. This is the unlabeled cell fraction.
6. Remove the first column from separator, place the second column into the separator, and put the first column on top of the second one.
7. Pipette appropriate amount of cold buffer onto the first column. Immediately flush out the fraction with the magnetically labeled cells by firmly applying the plunger, supplied with the column. directly onto the second column.
MS: 1 mL LS: 5 mL

8. Collect unlabeled cells that pass through and wash with appropriate amount of cold buffer. Perform washing steps by adding buffer successively once the column reservoir is empty.
MS: 3 \times 500 μ L LS: 3 \times 3 mL
9. Remove the second column from separator, place the column on a suitable collection tube.
10. Pipette appropriate amount of cold buffer onto the column. Immediately flush out the fraction with the magnetically labeled cells by firmly applying the plunger, supplied with the column.
MS: 500 μ L LS: 5 mL
▲ **Note:** For subsequent cell culture, the cells can also be eluted with medium. If part of the cells are analyzed by flow cytometry, the medium should **not contain** phenol red.
11. Proceed to analysis (see section 5.), cell culture or other subsequent experiment.



Magnetic separation using the autoMACS™ Separator

▲ Refer to the autoMACS™ User Manual for instructions on how to use the autoMACS Separator.

1. Prepare and prime autoMACS Separator.
2. (Optional) Pass cells through Pre-Separation Filters (# 130-041-407) to remove clumps.

3. Place tube containing magnetically labeled cells in autoMACS Separator. Choose separation program "Posseld". Collect the separated fractions from outlet port "pos2".
4. Proceed to analysis (see section 5.), cell culture or other subsequent experiment.

5. Detection and analysis of IFN- γ -secreting T cells

▲ Add propidium iodide (PI) or 7-AAD to a final concentration of 0.5 μ g/mL **just prior** to acquisition for exclusion of dead cells from flow cytometric analysis. Incubating with PI for longer periods will affect the viability of the cells.




Do not fix the cells when using PI or 7-AAD.

▲ For optimized sensitivity, an appropriate number of viable cells has to be acquired from the antigen stimulated sample as well as from the control sample.

- Acquire 2 \times 10⁵ viable cells from the fraction before enrichment (see 4.4 step 5.).
- For **enumeration** of low frequent IFN- γ -secreting cells, acquire all of the positive fraction. For **preparative purposes**, acquire an aliquot of the positive fraction to determine the performance of the cell enrichment.


To illustrate the analysis, we describe the detection of IFN- γ -secreting T cells using the IFN- γ Secretion Assay. The detailed description, including how to set gates, should serve as a model for the analysis of your own sample.

1. 10⁷ human PBMC of a CMV⁺ donor have been restimulated for 16 hours with and without CMV lysate (5 μ g/mL; Biowhittaker).
2. The IFN- γ Secretion Assay was performed on the stimulated and the unstimulated sample.
3. **Counterstaining of T cells** was performed using CD4-FITC.
4. **Monocytes** were stained with CD14-PerCP[™].
5. **Dead cells** were stained with propidium iodide (PI), which was added just prior to flow cytometric analysis to a final concentration of 0.5 μ g/mL.
6. 200,000 viable cells of the fractions before enrichment and the complete enriched fractions were acquired by flow cytometry, from the stimulated and the unstimulated samples.
7. A **lymphocyte gate** based on forward and side scatter (FSC/SSC) properties was activated prior to further gating to exclude monocytes and debris (see A.).
8. Dead cells and monocytes were excluded according to PI- and CD14-PerCP[™]-staining in a fluorescence 2 (PE) versus fluorescence 3 plot (PerCP) (see B.).
- The **dead cell exclusion** is crucial for the analysis of rare antigen-specific T cells, as dead cells may bind non-specifically to antibodies or MicroBeads. This could lead to false positive events.
- The sensitivity of detection is further enhanced by exclusion of undesired non-T cells, like monocytes which may cause non-specific background staining.
9. Analysis of secreted IFN- γ (PE) versus CD4-FITC staining by viable lymphocytes is displayed (see C.).

A. Lymphocyte gate in the forward versus side scatter plot**B. Dead cell and monocyte exclusion in FL-2 versus FL-3****C. Antigen-specific CD4⁺ T cells stained for secreted IFN-γ****Sample stimulated with CMV lysate**

$$\% \text{ IFN-}\gamma^+ \text{ cells among CD4}^+ = \frac{\# \text{ of IFN-}\gamma^+ \text{CD4}^+ \text{ cells in the analyzed sample}}{\# \text{ of total CD4}^+ \text{ cells in the analyzed sample}} \times 100$$

$$\% \text{ IFN-}\gamma^+ \text{ cells among CD4}^+ = \frac{\text{abs. \# of IFN-}\gamma^+ \text{CD4}^+ \text{ cells in the enriched fraction}}{\text{abs. \# of total CD4}^+ \text{ cells before enrichment}} \times 100$$

Unstimulated control sample**Miltenyi Biotec**

6. References

6. References

- Manz, R; Assenmacher, M; Pflüger, E; Miltenyi, S; Radbruch, A (1995) Analysis and Sorting of Live cells According to Secreted Molecules Relocated to a Cell-Surface Affinity Matrix. Proc. Natl.Acad.Sci. USA 92: 1921–1925. [139]
- Assenmacher, M; Löhning, M; Scheffold, A; Manz, RA; Schmitz, J; Radbruch, A. (1998) Sequential production of IL-2, IFN-γ and IL-10 by individual staphylococcal enterotoxin B-activated T helper lymphocytes. Eur. J. Immunol. 28: 1534–1543. [483]
- Brosterhus, H; Brings, S; Leyendeckers, H; Manz, RA; Miltenyi, S; Radbruch, A; Assenmacher, M; Schmitz, J (1999) Enrichment and detection of live antigen-specific CD4⁺ and CD8⁺ T cells based on cytokine secretion. Eur. J. Immunol. 29: 4053–4059. [573]
- Oelke, M; Moehrl, U; Chen, JL; Behringer, D; Cerundolo, V; Lindemann, A; Mackensen, A (2000) Generation and purification of CD8⁺ Melan-A-Specific Cytotoxic T Lymphocytes for Adoptive Transfer in Tumor Immunotherapy. Clin. Cancer Res. 6: 1997–2005. [663]
- Oelke, M; Kurokawa, T; Hentrich, I; Behringer, D; Cerundolo, V; Lindemann, A; Mackensen, A (2000) Functional Characterization of CD8⁺ Antigen-Specific Cytotoxic T Lymphocytes after Enrichment Based on Cytokine Secretion: Comparison with the MHC-Tetramer Technology. Scand. J. Immunol. 52: 544–549. [970]
- Bickham, K; Münz, C; Tsang, ML; Larsson, M; Fonteneau, J-F; Bhardwaj, N; Steinmann, R (2001) EBNA1-specific CD4⁺ T cells in healthy carriers of Epstein-Barr virus are primarily Th1 in function. J. Clin. Invest. 107: 121–130. [1035]
- Pittet, MJ; Zippelius, A; Speiser, DE; Assenmacher, M; Guillaume, P; Valmori, D; Lienard, D; Lejeune, F; Cerottini, JC; Romero, P (2001) Ex vivo IFN-γ secretion by circulating CD8 T lymphocytes: Implications of a novel approach for T cell monitoring in infectious malignant diseases. J. Immunol. 166: 7634–7640. [1037]
- Becker, C; Pohla, H; Frankenberger, F; Schüller, T; Assenmacher, M; Schendel, DJ; Blankenstein, T (2001) Adoptive tumor therapy with T lymphocytes enriched through an IFN-γ capture assay. Nature Medicine 7, 10: 1159–1162. [1207]

For further information visit our website www.miltenyibiotec.com.

Miltenyi Biotec

7. Appendix

7. Appendix:**A: Flask and dish sizes for stimulation**

For *in vitro* stimulation (see 4.2 step 2.) the cells should be resuspended in culture medium, containing 5% of human serum, at 10⁷ cells/mL and 5×10⁶ cells/cm². Both the dilution and the cell density are important to assure optimum stimulation.

The following table lists culture plate, dish and flask sizes suitable for different cell numbers. It also indicates the appropriate amount of medium to add.

total cell number	medium volume to add	culture plate	well diameter
0.15 × 10 ⁷	0.15 mL	96 well	0.64 cm
0.5 × 10 ⁷	0.5 mL	48 well	1.13 cm
1 × 10 ⁷	1 mL	24 well	1.6 cm
2 × 10 ⁷	2 mL	12 well	2.26 cm
5 × 10 ⁷	5 mL	6 well	3.5 cm
total cell number	medium volume to add	culture dish	dish diameter
4.5 × 10 ⁷	4.5 mL	small	3.5 cm
10 × 10 ⁷	10 mL	medium	6 cm
25 × 10 ⁷	25 mL	large	10 cm
50 × 10 ⁷	50 mL	extra large	15 cm
total cell number	medium volume to add	culture flask	growth area
12 × 10 ⁷	12 mL	50 mL	25 cm ²
40 × 10 ⁷	40 mL	250 mL	75 cm ²
80 × 10 ⁷	80 mL	720 mL	162 cm ²
120 × 10 ⁷	120 mL	900 mL	225 cm ²

B: Detection and enrichment of cytokine secreting cells from whole blood

B1. Reagent and instrument requirements

B2. Protocol

- B 2.1 (Antigen-specific) *in vitro* stimulation
- B 2.2 Cytokine Secretion Assay
- B 2.3 Magnetic labeling
- B 2.4 Magnetic separation

The following special protocol can be used in combination with one of the Cytokine Secretion Assay - Cell Enrichment and Detection Kits for human cells.

B 1. Reagent and instrument requirements

- Cytokine Secretion Assay Kit, for example:

IFN- γ Secretion Assay - Cell Enrichment and Detection Kit (PE)	(# 130-054-201)
IL-2 Secretion Assay - Cell Enrichment and Detection Kit (PE)	(# 130-090-488)
IL-4 Secretion Assay - Cell Enrichment and Detection Kit (PE)	(# 130-054-101)
IL-10 Secretion Assay - Cell Enrichment and Detection Kit (PE)	(# 130-090-435)

- **Anticoagulant:** sodium heparin
- **Buffer** (degassed): Prepare a solution containing PBS (phosphate buffered saline) pH 7.2, 0.5% BSA and 2 mM EDTA by diluting MACS BSA Stock Solution (# 130-091-376) 1:20 with autoMACS™ Rinsing Solution (# 130-091-222). Keep buffer cold (4–8 °C).
- **Culture medium**, e.g. RPMI 1640 (# 130-091-440) containing 20% of human serum, like autologous serum or AB serum (**do not use** BSA or FCS because of non-specific stimulation).
- **Erythrocyte lysing solution (1x):**
 - prepare freshly from 10 \times stock solution.
 - **10 \times stock solution:** 41.4 g NH₄Cl (1.55 M), 5 g KHCO₃ (100 mM), 1 mL 0.5 M EDTA (1 mM), adjust pH to 7.3, fill up to 500 mL with dd H₂O.
- ▲ **Note:** Do not use FACS Lysing solution™.
- (Optional) **Staining reagents:** CD4-FITC (# 130-080-501) or CD8-FITC (# 130-080-601) and CD14-PerCP™.
- ▲ **Note:** Do **not** use tandem conjugates of phycoerythrin, like Cy-Chrome® (PharMingen), PE-Cy5 (Serotec), ECD, PC5 (Coulter-Immunotech) etc., they may also be recognized by the Anti-PE MicroBeads.
- ▲ **Note:** Upon activation of T cells, TCR and some associated molecules, like CD3, might be down-regulated.
- ▲ **Note:** For optimal sensitivity, we recommend labeling of undesired non-T cells such as monocytes with antibodies conjugated to PerCP™, e.g. CD14-PerCP™. These cells can then be excluded together with PI stained dead cells by gating.

Miltenyi Biotec



- **Propidium iodide** (PI) or 7-AAD to exclude dead cells from analysis.
- MACS Columns and MACS Separators:

Column	max. number of labeled cells	max. number of total cells	Separator
MS	10 ⁷	2 \times 10 ⁸	MiniMACS, OctoMACS; with Column Adapter: VarioMACS, SuperMACS
autoMACS	2 \times 10 ⁸	4 \times 10 ⁹	autoMACS

▲ **Note:** Column adapters are required to insert certain columns into VarioMACS™ Separator or SuperMACS™ Separator. For details, see MACS Separator data sheets.

- (Optional) Rotation device for tubes: MACSmix™ tube rotator (# 130-090-753)
- (Optional) Pre-Separation Filters (# 130-041-407) to remove cell clumps.

B 2. Protocol

B 2.1 (Antigen-specific) *in vitro* stimulation

▲ The peripheral blood should not be older than 20 hours and should be supplemented with anticoagulant **sodium heparin**. **Do not use EDTA, or ACD.** Lymphocyte activation and secretion of cytokines requires calcium, and is consequently inhibited by chelating anticoagulants.

▲ **Note:** Whole blood may be stored overnight at **room temperature**.

▲ Always include a **negative control** sample in the experiment. A **positive control** with e.g. Staphylococcal Enterotoxin B (SEB) may be included in the experiment (see also detailed protocol provided with the Cytokine Secretion Assay Kits).

▲ **Do not use** media containing any **non-human** proteins, like BSA or FCS because of non-specific stimulation.



Protocol for *in vitro* stimulation

1. Start with **5 mL of fresh, sodium heparinized, human blood** (containing about 10⁷ lymphocytes) in a 50 mL conical polypropylene tube.

Miltenyi Biotec



2. Add the antigen or, as a positive control, 1 µg/mL SEB for 3–16 hours at 37 °C, 5–7% CO₂ (for details on the kinetics of cytokine secretion and on concentrations of antigen to add, refer to Cytokine Secretion Assay data sheet, 3.1–3.2).
3. A negative control sample, treated exactly the same as the antigen-stimulated sample but without addition of antigen, should always be included in the experiment.
4. (Optional) Costimulatory agents like CD28 and CD49d antibodies may be added.

B 2.2 Cytokine Secretion Assay

▲ This protocol is optimized for cell samples containing < 5% of total cytokine secreting cells. If ≥ 5% of cytokine secreting cells are expected, it is necessary to dilute the cells further during the cytokine secretion period, and therefore a larger test tube will be needed. The dilution avoids non-specific staining of cells not secreting cytokines during this period.

▲ For each sample with 5 mL whole blood prepare:

- 100 mL of **cold buffer** (4–8 °C)
- 200 µL of **cold medium** (4–8 °C)
- 7 mL of **warm medium** (37 °C)
- 45 mL of **erythrocyte lysing solution** (room temperature).

▲ Work fast, keep the cells cold, use pre-cooled solutions which will prevent capping of antibodies on the cell surface and a non-specific cell labeling (exception: warm medium during secretion period and room temperature during lysing step).

▲ Do not remove supernatant by decanting. This will lead to cell loss and incorrect incubation volumes. Pipette off or aspirate supernatant.

▲ Dead cells may bind non-specifically to MACS[®] MicroBeads or antibodies. Therefore, when working with cell preparations containing large amounts of dead cells, they should be removed before starting the Cytokine Secretion Assay, e.g. by density gradient centrifugation or by using the Dead Cell Removal Kit (# 130-090-101).

▲ Higher temperatures and longer incubation times for staining should be avoided. This will lead to non-specific cell labeling.



Lysis of erythrocytes

1. After stimulation add 45 mL of erythrocyte lysing solution to 5 mL whole blood sample.
2. Mix gently and incubate for 10 minutes at **room temperature**. Rotate tube continuously using the MACSmix[™] tube rotator (# 130-090-753), or turn tube several times during incubation.
3. Centrifuge cells at 300×g for 10 minutes at **room temperature**, remove supernatant **completely**.

Miltenyi Biotec



Labeling cells with Cytokine Catch Reagent

1. Resuspend cell pellet in 15 mL of **cold buffer**, and transfer into a 15 mL conical propylene tube.
2. Centrifuge at 300×g for 10 minutes at 4–8 °C. Pipette off supernatant completely.
3. Resuspend pellet in 160 µL of **cold medium**.
4. Add 40 µL of **Cytokine Catch Reagent**, mix well and incubate for 5 minutes **on ice**.



Cytokine secretion period

1. Add 7 mL of **warm medium** (37 °C) to dilute the cells.

▲ **Note:** For frequencies of cytokine secreting cells ≥ 5% the cells need to be further diluted, e.g. by a factor of 5.
2. Incubate cells in a closed tube for 45 minutes at 37 °C under slow continuous rotation using the MACSmix tube rotator, or turn tube every 5 minutes to resuspend settled cells.

▲ **Note:** During this step it is crucial to prevent contact of cells to avoid cross contamination with cytokines.



Labeling cells with Cytokine Detection Antibody

1. Put the tube **on ice**.
2. Wash cells by adding 8 mL of **cold buffer**, centrifuge at 300×g for 10 minutes at 4–8 °C. Pipette off supernatant completely.
3. Resuspend cell pellet in 160 µL of **cold buffer**.
4. Add 40 µL of **Cytokine Detection Antibody (PE)**.
5. (Optional) Add additional staining reagents, e.g. 20 µL of CD4-FITC (# 130-080-501) or CD8-FITC (# 130-080-601) and CD14-PerCP[™].
6. Mix well and incubate for 10 minutes **on ice**.
7. Wash cells by adding 10 mL of **cold buffer**, centrifuge at 300×g for 10 minutes at 4–8 °C. Pipette off supernatant completely.

B 2.3 Magnetic labeling



Magnetic labeling with Anti-PE MicroBeads

1. Resuspend cell pellet in 160 µL of **cold buffer**.
2. Add 40 µL of **Anti-PE MicroBeads**, mix well and incubate for 15 minutes at 4–8 °C.

▲ **Note:** Incubate in refrigerator at 4–8 °C; do not work on ice during this step.

Miltenyi Biotec



3. Wash cells by adding 10 mL of **cold buffer**, centrifuge at 300×g for 10 minutes at 4–8 °C. Pipette off supernatant completely.
4. Resuspend cell pellet in 500 µL of **cold buffer**.
5. (Optional) Take an aliquot for flow cytometric analysis and cell count of the fraction before enrichment.
6. Proceed to magnetic separation.

B 2.4 Magnetic separation



Magnetic separation using MS Columns

▲ When enriching antigen-specific T cells, **always perform two consecutive MS Columns** to achieve best results.

1. Prepare **two MS Columns** per sample by rinsing with 500 µL **cold buffer**, discard effluent.
2. Place first column into the magnetic field of a MACS® Separator (use column adapter with VarioMACS or SuperMACS Separator).
3. (Optional) Pass cells through Pre-Separation Filters (# 130-041-407) to remove clumps.
4. Apply cell suspension onto the column.

5. Collect unlabeled cells which pass through and wash with 3×500 µL of cold buffer. Perform washing steps by adding buffer successively once the column reservoir is empty. Collect total effluent. This is the unlabeled cell fraction.
6. Remove first column from separator, place second column into the separator, and put the first column on top of the second one.
7. Pipette 1 mL of cold buffer on top of the first column. Immediately flush out the fraction with the magnetically labeled cells by firmly applying the plunger, supplied with the column, directly onto the second column.
8. Collect unlabeled cells that pass through and wash with 3×500 µL of cold buffer. Perform washing steps by adding buffer successively once the column reservoir is empty.
9. Remove second column from separator, place column on a suitable collection tube.
10. Pipette 500 µL of cold buffer on top of the column. Immediately flush out the fraction with the magnetically labeled cells by firmly applying the plunger, supplied with the column.

▲ **Note:** For subsequent cell culture, the cells can also be eluted with medium. If part of the cells are analysed by flow cytometry, the medium should **not contain** phenol red.
11. Proceed to flow cytometric analysis (see detailed protocol), cell culture or other subsequent experiment.

Miltenyi Biotec



Magnetic separation using the autoMACS™ Separator

▲ Refer to the autoMACS™ User Manual for instructions on how to use the autoMACS Separator.

1. Prepare and prime autoMACS Separator.
2. (Optional) Pass cells through Pre-Separation Filters (# 130-041-407) to remove clumps.
3. Place tube containing magnetically labeled cells in autoMACS Separator. Choose separation program “Posseld”. Collect the separated fractions from outlet port “pos2”.
4. Proceed to flow cytometric analysis (see detailed protocol), cell culture or other subsequent experiment.

Warnings

Reagents contain sodium azide. Under acidic conditions sodium azide yields hydrazoic acid, which is extremely toxic. Azide compounds should be diluted with running water before discarding. These precautions are recommended to avoid deposits in plumbing where explosive conditions may develop.

Warranty

The products sold hereunder are warranted only to be free from defects in workmanship and material at the time of delivery to the customer. Miltenyi Biotec GmbH makes no warranty or representation, either expressed or implied, with respect to the fitness of a product for a particular purpose. There are no warranties, expressed or implied, which extend beyond the technical specifications of the products. Miltenyi Biotec GmbH's liability is limited to either replacement of the products or refund of the purchase price. Miltenyi Biotec GmbH is not liable for any property damage, personal injury or economic loss caused by the product.

Cy-Chrome® is a trademark of PharMingen.

Peridin Chlorophyll Protein (PerCP™) is a trademark of Becton Dickinson.

MACS® is a registered trademark of Miltenyi Biotec GmbH.

Miltenyi Biotec

