
Revisions

SO 0191-5


Rev from	Rev to	ECO #
F	G	3528-05

Notes:


- BD Cat. Number 442003, 442288
- Blank (Sheet) Size : Length: 24" Width: 27.5"
 Number of Pages: 30 Number of Sheets: 1
 Page Size: Length 8" Width 5.5" Final Folded Size: 5.5" x 8"
- Style (see illustrations below): # 7


#1


#2


#3


#4

#5

#6

#7

- See Specification Control Number N/A for Material Information
- Ink Colors: Printed two sides ☒ Yes ☐ No
 No. of Colors: 1 PMS# Black
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

VS Controlled by BD Caribe, Ltd.

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: PP124JAA		Category and Description Package Insert, BACTEC Myco/F Lytic	Sheet: 1 of 31	A
			Scale: N/A	

BACTEC™ Myco/F Lytic Culture Vials

Supplemented Middlebrook 7H9 and Brain Heart Infusion Broth

FOR USE WITH BACTEC 9000MB

English:	pages	1 – 4	Español:	páginas	15 – 18
Français :	pages	4 – 7	Dansk:	side	18 – 21
Deutsch:	Seiten	8 – 11	Português:	páginas	21 – 25
Italiano:	pagine	11 – 14	Svenska:	sidan	25 – 28

PP124JAA(G)
2005/05

See symbol glossary at end of insert. / Viz popis symbolů na konci příbalového letáku. / Se symbolglossaret i slutningen af indlægssedlen. / Zie lijst met symbolen aan het einde van de bijsluiter. / Vaadake sümbolite seletust infolehe lõpus. / Katso pakkauselosteen lopussa olevaa kuvamerkkien sanastoa. / Voir le glossaire des symboles à la fin de la notice. / Siehe Symbol-Erklärungen am Ende der Packungsbeilage. / Δείτε το γλωσσάριο των συμβόλων στο τέλος του ένθετου. / A jelmagyarazát a használati utasítás végén találhatja. / Vedere il glossario dei simboli alla fine del foglio illustrativo. / Ζητουμε informacijo lapelo pabalgoje pateikiamą simbolių glosarijų. / Se i symbolforklaringen på slutten av produktvedlegget. / Zobacz objaśnienie symboli na końcu ulotki. / Consulte o glossário de símbolos no fim do folheto informativo. / Pozri slovník symbolov na konci letáka. / Consulte el glosario de símbolos al final del prospecto. / Se symbolförteckningen vid slutet av bipacksedeln.

Pokyny vám poskytne miestni zástupce spoločnosti BD. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Instrukcie ziskate u miestneho zástupcu spoločnosti BD.

NAME

BACTEC™ culture vials type Myco/F Lytic (a modified Middlebrook 7H9 broth) support the growth and detection of mycobacteria.

INTENDED USE

BACTEC Myco/F Lytic when used with the BACTEC 9000MB instrument is a non-selective culture medium for the qualitative culture and recovery of mycobacteria from blood specimens.

SUMMARY AND EXPLANATION OF THE TEST

Since the mid-1980s and spread of the AIDS epidemic, the incidence of septicemia caused by opportunistic mycobacteria has risen. *Mycobacterium tuberculosis* (MTB) and mycobacteria other than tuberculosis (MOTT), especially *Mycobacterium avium* complex (MAC), have become resurgent. From 1985 to 1992, the number of MTB cases reported increased 18%. Between 1981 and 1987, AIDS case surveillances indicated that 5.5% of the patients with AIDS had disseminated nontuberculous mycobacterial infections, e.g., MAC. By 1990, the increased cases of disseminated nontuberculous mycobacterial infections had resulted in a cumulative incidence of 7.6%.

The Centers for Disease Control and Prevention (CDC) have recommended that every effort must be made for laboratories to use the most rapid methods available for diagnostic mycobacteria testing. These recommendations include the use of a liquid medium for mycobacterial culture.^{1,2,3}

The BACTEC 9000MB System is designed for the rapid detection of mycobacteria in clinical specimens. BACTEC Myco/F Lytic Culture medium is a Middlebrook 7H9 and Brain Heart Infusion broth formulation for the recovery of mycobacteria from blood specimens. Specific modifications were made to enhance the growth and recovery of mycobacteria. These modifications include ferric ammonium citrate to provide an iron source for specific strains of mycobacteria, the addition of saponin as a blood lysing agent and the addition of specific proteins and sugars to provide nutritional supplements. Each vial contains a sensor which can detect decreases in oxygen concentration in the vial resulting from microorganism metabolism and growth. The sensor is monitored by the BACTEC 9000MB System for increasing fluorescence which is proportional to the decrease in oxygen. A positive determination indicates the presumptive presence of viable microorganisms in the vial.

PRINCIPLES OF THE PROCEDURE

The BACTEC Myco/F Lytic culture vial is designed for the rapid detection of mycobacteria in blood. Blood specimens are inoculated into the BACTEC Myco/F Lytic vial either with a syringe or direct draw with a needle and tubing. The vial is placed into the BACTEC 9000MB System and is continuously incubated at 37°C, with agitation once every ten minutes for maximum recovery. Each vial contains a sensor which can detect decreases in oxygen concentration in the vial resulting from microorganism metabolism and growth. The sensor is monitored by the BACTEC 9000MB System every ten minutes. Analysis of the rate of oxygen decrease as measured by increasing fluorescence enables the BACTEC fluorescent series instrument to determine if the vial is instrument positive. A positive determination indicates the presumptive presence of viable microorganisms in the vial. Detection is limited to microorganisms that will grow in the medium at 37°C. The medium is not selective and will support the growth of other aerobic organisms including yeast, fungi and bacteria which may interfere, if present, with the recovery of mycobacteria. Culture vials which remain negative for 42 days and which show no visible sign of positivity are removed from the instrument and sterilized prior to discarding.

REAGENTS

Each BACTEC Myco/F Lytic culture vial contains the following active ingredients prior to processing:

List of Ingredients

Processed Water	40 mL qs
7H9 Middlebrook Broth Base without phosphate salts	0.12% w/v
Brain Heart Infusion	0.5% w/v
Casein Hydrolysate	0.10% w/v
Supplement H	0.10% w/v
Inositol	0.05% w/v
Glycerol	0.10% w/v
Sodium Polyanetholsulfonate	0.025% w/v
Tween 80	0.0025% w/v
Pyridoxal HCl	0.0001% w/v
Ferric Ammonium Citrate	0.006% w/v
Potassium Phosphate	0.024% w/v
Saponin	0.24% w/v
Antifoam	0.01% w/v

This BACTEC media is dispensed with added CO₂ and O₂.

Composition may have been adjusted to meet specific performance requirements.

BACTEC Myco/F Lytic medium requires no supplement addition. Each 40 mL vial of BACTEC Myco/F Lytic is ready for use when received. The appearance of the media upon receipt shall be clear and light amber in color.

WARNINGS

Precautions: *in vitro* diagnostic.

This Product Contains Dry Natural Rubber.

POTENTIAL INFECTIOUS TEST SPECIMEN. Observe "Universal Precautions"^{4,5} and institutional guidelines when handling and disposing of infectious materials.

BACTEC Myco/F Lytic vials will accept more than the recommended maximum of 5 mL of specimen volume, monitoring of fill volume should be conducted.

Biosafety Level 2 practice, containment equipment and facilities are recommended for preparing acid-fast stains and for culturing clinical specimens. For activities involving the propagation and manipulation of *Mycobacterium tuberculosis* or *Mycobacterium bovis* grown in culture, Biosafety Level 3 practice, containment equipment and facilities are recommended.^{5,6,7}

Prior to use, each vial should be examined for evidence of contamination such as cloudiness, bulging or depressed septum, or leakage. **DO NOT USE** any vial showing evidence of contamination, leakage or damage. Vial contamination may not be readily apparent. A contaminated vial could contain positive pressure. If a contaminated vial is used for direct draw, gas or contaminated culture media could be refluxed into the patient's vein. On rare occasions, the glass bottle neck may be cracked and the neck may break during removal of the flip-off cap or in handling. Also, on rare occasions, a vial may not be sealed sufficiently. In both cases the contents of the vials may leak or spill, especially if the vial is inverted.

To minimize the potential of leakage during inoculation by syringe of specimen into culture vials, use syringes with **LUER-LOK™** brand tips. A one-handed inoculation technique and a suitable vial holder should be employed to prevent accidental needle stick injury.

Before discarding, sterilize all inoculated **BACTEC Myco/F Lytic** vials by autoclaving.

Before sampling **positive culture vials for subculturing or staining, etc.:** It is necessary to release gas which often builds up due to microbial metabolism. Sampling must be performed in a biological safety cabinet, and appropriate protective clothing, including gloves and masks, should be worn. See Procedure Section for more information on subculturing.

LEAKING OR BROKEN VIALS

CAUTION: Because an inoculated leaking or broken vial may produce an aerosol of mycobacteria, including *M. tuberculosis* or other bacteria, appropriate handling should be observed.

If an inoculated vial is found to be leaking or is accidentally broken during collection or transport, use the established procedure in your facility for dealing with mycobacterial spills. As a minimum, "Universal Precautions" should be employed. Vials should be discarded in an appropriate manner.

In the rare instance where a vial is found to have leaked contents into the instrument proper, or if a vial is accidentally broken, turn off the instrument immediately. Vacate the affected area. Contact your facility's Safety or Infection Control Officer(s). Determine the necessity of turning off or modifying the settings of the air handling units serving the affected area. Do not return to the area until any potential aerosols have settled or have been removed by appropriate ventilation. Becton Dickinson and Company should be notified by calling 1-800-544-7434 in the USA or the appropriate Becton Dickinson representative in your area. Guidelines for proper handling of accidental mycobacterial contamination due to breakage of culture tubes or broth suspensions have been issued by the CDC.^{5,6,7}

STORAGE INSTRUCTIONS

Store at 2° – 25°C in a dry location out of direct sunlight.

SPECIMEN COLLECTION

NOTE: It is recommended that this procedure be reviewed with appropriate personnel prior to use of medium to insure proper specimen collection techniques as described in this section.

The specimen must be collected using sterile technique to reduce the chance of contamination. The range of blood volume which can be cultured is 1.0 mL to 5.0 mL. It is recommended that the specimen be inoculated at bedside. Most commonly, a syringe with a **LUER-LOK™** brand tip is used to draw the specimen. If appropriate, a **VACUTAINER™** brand Needle Holder and a **VACUTAINER™** brand Blood Collection Set, **VACUTAINER™ SAFETY-LOK™** Blood collection Set or other tubing "butterfly" set may be used. If using a needle and tubing (direct draw), carefully observe the direction of the blood flow when starting sample collection. Prior to inoculation, the medium fill volume should be noted on the label with a pen or marker to indicate the starting point of specimen collection. The vacuum in the bottle will usually exceed 5 mL, so the user should monitor the volume collected by means of the 5 mL graduation marks on the vial label. When the desired 1 – 5 mL of blood has been drawn, the flow should be stopped by crimping the tubing and removing the needle from the **BACTEC** vial. The **BACTEC** vial should be transported as quickly as possible to the laboratory and placed in the **BACTEC** instrument. A yellow-top **VACUTAINER** Brand Tube containing SPS may also be used to collect the blood sample from the patient. The tube should be transported to the laboratory as quickly as possible for transfer into the **BACTEC** culture vial.

PROCEDURE

Materials Provided: **BACTEC Myco/F Lytic Culture Vials.**

Materials Required But Not Provided: Biological Safety Cabinet, Autoclave, venting unit, mycobacterial disinfectant, 70% isopropyl alcohol, Quality Control Organisms (*Mycobacterium intracellulare*, ATCC 13950; *Mycobacterium kansasii*, ATCC 12478; and *Mycobacterium fortuitum*, ATCC 6841), microscope and materials for staining slides and subculturing vials.

CAUTION: **BACTEC Myco/F Lytic vials must be used with instrument software version 3.6 or higher.**

Inoculation of **BACTEC Myco/F Lytic Culture Vials**

1. Remove the flip-off cap from the **BACTEC** vial top and inspect the vial for cracks, leaks, contamination, excessive cloudiness, and bulging or indented septum. **DO NOT USE** if any defect is noted.
2. Label culture vial with specimen identification and mark medium fill graduation line on vial label.
3. Before inoculating, swab the septum with alcohol. Aseptically inject with a syringe or draw directly with the aid of the graduation lines on the vial label 1 – 5 mL of specimen per vial (see the section on Limitations of the Procedure). **Inoculated vials should be placed into the BACTEC 9000MB instrument as soon as possible** for incubation and monitoring.
4. Vials entered into the instruments will be automatically tested for the duration of the testing protocol. Positive vials will be identified by the **BACTEC** Fluorescent System (see the **BACTEC** User's Manual, MA-0092). The sensor inside the vial may not appear visibly different in positive or negative vials; however, the **BACTEC** Fluorescent System can determine a difference in sensor fluorescence.
5. Positive vials should be subcultured and an appropriate smear prepared. All positive vials should be handled using BSL III practices and containment facilities.

Processing an instrument-positive vial

- a) Remove the vial from the instrument and transport to an area using BSL III practices and containment facilities.
- b) Invert vial to mix contents.
- c) In biological safety cabinet, vent the vial to equilibrate vial pressure with atmosphere.
- d) Remove aliquot from vial (approx. 0.1 mL) for stain preparations (AFB and Gram).
- e) Inspect smear and report preliminary results only after smear evaluation.

If at the end of six weeks incubation an instrument negative vial appears to be positive (i.e., bulging septum, or very darkened blood), it should be subcultured, AFB stained and treated as a presumptive positive, providing the stain result is positive.

Subculturing of Vial: Subculturing should be performed in a biological safety cabinet, and appropriate clothing, including gloves and masks, should be worn. Prior to subculturing, place the vial in an upright position, and place an alcohol wipe over the septum. To release any positive pressure in the vial which could be caused by growth of possible contaminants, insert a sterile 25-gauge (or smaller) needle equipped with an appropriate filter or pledget through the alcohol wipe and septum. The needle should be removed after any pressure is released and before sampling the vial for subculture. The insertion and withdrawal of the needle should be done in a straight-line motion, avoiding any side-to-side motions which could permanently damage the septum. **Do not re-cap the needle. Discard needles and syringes in a puncture-resistant biohazard container.**

QUALITY CONTROL

Quality Control Certificates are provided with each carton of media.

It is recommended that each new shipment or lot of **BACTEC Myco/F Lytic** media be tested with the ATCC control organisms identified in the chart below as a positive control, and an uninoculated vial as a negative control.

Organism	Range of Time-to-detection (days)
<i>Mycobacterium intracellulare</i> , ATCC 13950	8 to 16
<i>Mycobacterium kansasii</i> , ATCC 12478	3 to 13
<i>Mycobacterium fortuitum</i> , ATCC 6841	1 to 3

The positive vials should be inoculated using a 1:100 dilution of a McFarland #1 suspension grown on solid medium. Inoculate the vial with 0.1 mL of the diluted culture. The vials and an uninoculated control vial should be scanned into the instrument and tested. The inoculated vial should be detected as positive by the instrument within the test protocol. The negative control should remain negative. If expected results for Quality Control are not obtained, do not use the medium and contact Becton Dickinson Technical Services (in the U.S. only: 1-800-638-8663) or your local Becton Dickinson Representative for further assistance.

For information on quality control for the **BACTEC** System, refer to the User's Manual (MA- 0092).

Reporting of RESULTS

An instrument positive vial may be confirmed by acid-fast smear or Gram stain. A positive result indicates the presumptive presence of viable microorganisms in the vial.

If AFB smear positive, subculture to solid media and report as: instrument-positive, AFB smear positive, ID pending.

If microorganisms other than acid-fast bacilli are present, subculture to solid media and report as: instrument positive, AFB smear negative, ID pending.

If no microorganisms are present on the smears, subculture to solid media, re-enter the vial into the instrument as an ongoing negative vial within 5 hours of removal and allow to complete test protocol. No reportable result.

Perform subcultures from the **BACTEC Myco/F Lytic** vial for identification and susceptibility testing.

LIMITATIONS OF THE PROCEDURE

Detection of mycobacterial species in blood specimens is dependent on the number of organisms present in the specimen, specimen collection methods, and patient factors such as presence of symptoms and prior treatment.

Mycobacteria may vary in acid-fastness depending on strain, age of culture and other variables.

Care must be taken to prevent contamination of the sample during collection and inoculation into the **BACTEC** vial. A contaminated vial will give a positive instrument reading, but will not indicate a relevant clinical result. Such a determination must be made by the user, based on such factors as, stain results, type of organism recovered, occurrence of the same organism in multiple cultures, patient history, etc.

BACTEC Myco/F Lytic vials are not selective and will support the growth of other aerobic organisms besides mycobacteria. Positive vials may contain one or more species of mycobacteria and/or other non-mycobacterial species. If present, fast growing organisms may mask the detection of slower growing mycobacteria. Subculture and additional procedures are required. The consistency of microscopic morphology in **BACTEC Myco/F Lytic** has not been established.


Optimum recovery of isolates will be achieved by adding 1 – 5 mL of blood to each vial. Use of lower or higher volumes may adversely affect recovery, detection times and/or specificity. False positivity most likely will increase when the blood volume is above 5 mL.

Blood may contain antimicrobials or other inhibitors which may slow or prevent the growth of microorganisms.

BACTEC Myco/F Lytic vials are incubated at 37°C potentially precluding the recovery of mycobacteria requiring other incubation temperatures (e.g., *M. marinum*, *M. ulcerans*, *M. haemophilum*). Recovery of such organisms requires additional culture methods. The following isolates were detected as positive in the **BACTEC 9000MB** instrument using **BACTEC Myco/F Lytic** medium during internal studies and/or clinical trials: *M. tuberculosis*, *M. kansasii*, *M. fortuitum*, *M. avium*, *M. intracellulare*, *M. bovis*, *M. terrae*, *M. simiae*, *M. goodii*, *M. celatum*, *M. abscessus*, *M. malmoense*. During internal studies, *M. xenopi* and *M. szulgai* exhibited unsatisfactory recovery with **BACTEC Myco/F Lytic** culture medium.

EXPECTED RESULTS

Frequency distribution of recovery times (TTD) for clinical trial blood specimens positive with the **BACTEC Myco/F Lytic** Culture medium is illustrated in the following figure.

PERFORMANCE CHARACTERISTICS

The BACTEC Myco/F Lytic culture medium was evaluated with the BACTEC 9000MB instrument at two clinical sites considered large tertiary care teaching hospitals in geographically diverse areas. The site populations included patients suspected of a mycobacterial infection, immunocompromised patients and transplant patients. The BACTEC Myco/F Lytic culture medium was compared to the BACTEC 13A culture medium for the recovery and detection of mycobacteria from blood specimens. A total of 284 compliant blood specimens were tested during the study. The total number of pathogenic mycobacterial isolates recovered in the study was 39 (See TABLE 1). Of these positives, five (13%) were recovered in the BACTEC Myco/F Lytic culture medium only and two (5%) were recovered by BACTEC 13A culture medium only. A total of 28 BACTEC Myco/F Lytic vials were over filled with specimen (between 6 to 20 mL) during the clinical evaluation and were not included in this study since they were above the maximum fill volume (not compliant). Of these 28 BACTEC Myco/F Lytic vials, 16 (57%) were identified as false positive.

Of the 284 blood specimens tested in the clinical study, one BACTEC Myco/F Lytic vial (0.4%) was determined to be false positive (instrument-positive, smear and/or subculture-negative). Of the 38 instrument positive Myco/F Lytic vials, 1 (2.6%) was determined to be false positive. The false negative rate (instrument-negative, smear and/or subculture-positive) was determined to be 0% based on terminal subcultures of $\geq 50\%$ of negative vials. The contamination rate during this evaluation was determined to be 0.9%.

TABLE 1: SUMMARY OF MYCO/F LYTIC CULTURE MEDIUM ISOLATE RECOVERY DURING CLINICAL TRIALS

Organism	Total Isolates	Myco/F Lytic Medium Only	13A Medium Only	Both
All Pathogenic Mycobacteria:				
<i>Mycobacterium avium</i>	30	3	1	26
<i>Mycobacterium tuberculosis</i>	6	0	0	6
<i>Mycobacterium kansasii</i>	3	2	1	0
Total	39	5	2	32

REFERENCES

- Horsburg Jr., C.R. 1991. *Mycobacterium Avium* Complex Infection in the acquired immunodeficiency syndrome. New England Journal of Medicine 324:1332-1338.
- Tenover, F.C., et al, 1993. The resurgence of Tuberculosis: Is Your Laboratory Ready? Journal of Clinical Microbiology 31:767-770
- Isolation and Identification of *Mycobacterium tuberculosis*: A Guide for the Level II Laboratory. U.S. Department of Health and Human Services, Atlanta, GA, 1981.
- Kent, P.T., and G.P. Kubica. 1985. Public Health Mycobacteriology: A Guide for the Level III Laboratory. U.S. Department of Health and Human Services, Centers for Disease Control and Prevention, Atlanta, GA.
- Recommendations for preventing transmission of Human Immunodeficiency Virus and Hepatitis B Virus to patients during exposure-prone invasive procedures. MMWR 1991, Vol. 40, No. RR-8.
- Bloodborne pathogens. Code of Federal Regulations, Title 29, Part 1910.1030, Federal Register 1991, 56:64175-64182.
- Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services Public Health Service/ Centers for Disease Control and Prevention, Atlanta, GA, May, 1993.

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ATCC is a registered trademark of the American Type Culture Collection.
BD, BD Logo, BACTEC, VACUTAINER, SAFETY-LOK and LUER-LOK are trademarks of Becton, Dickinson and Company. © 2005 BD.
BENEX Limited