# RAPIDEC® staph

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# RAPIDEC® staph

Identification of the main staphylococci

IVD

## SUMMARY AND EXPLANATION

RAPIDEC staph is a standardized micromethod for the 2-hour identification of the main staphylococci isolated from human specimens (*Staphylococcus aureus*, *Staphylococcus epidermidis*, *Staphylococcus saprophyticus*) using colonies. RAPIDEC staph can also be used for the identification of *S. aureus* directly with blood culture broths.

## PRINCIPLE

The RAPIDEC staph strip incorporates 4 biochemical tests to identify 3 clinically significant staphylococci : *S. aureus*, *S. epidermidis* and  $\beta$ -galactosidase positive *S. saprophyticus*. *S. intermedius* (animal origin), *S. xylosus* and *Staphylococcus* spp. are presumptively identified. The biochemical tests are read and interpreted according to the sequence of the reactions.

The "aurease" (AUR) test is a spontaneous fluorescent reaction, specific for *S. aureus*, which is detected with the use of a UV lamp (365 nm). The alkaline phosphatase (PAL) test is a spontaneous reaction expressed by a color change. The  $\beta$ -galactosidase ( $\beta$ -GAL) test requires the addition of a reagent, FVB. Both tests have been widely used to identify *Staphylococcus* species. The conventional catalase (CAT) test may also be performed to confirm the identification of staphylococci.

Results are obtained after 2 hours of incubation at 37°C.

## CONTENT OF THE KIT (Kit for 30 tests) :

- 10 strips (3 tests each)
- 10 lids
- 1 package insert

## COMPOSITION OF THE STRIP

The composition of the RAPIDEC staph strip is given below in the list of tests :

CUP.	TESTS	SUBSTRATES	QTY (mg/cup.)	REACTIONS
C		Turbidity control (4 McFarland)		
S		Preparation of the bacterial suspension		
0		Negative control for fluorescence		
1	AUR	Prothrombin (human origin)*	0.000625 UI	Detection of «AUREase», enzyme specific for <i>Staphylococcus aureus</i>
		Enzymatic substrate	0.0096	
2	PAL	4-nitrophenyl-phosphate-2CHA	0.084	Alkaline Phosphatase
3	$\beta$ -GAL	2-naphthyl- $\beta$ D-galactopyranoside	0.092	$\beta$ -GALacto-pyranosidase

\* This product has been tested and shown to be negative for HBs antigen, antibodies to HIV1, HIV2 and HCV. However, since no existing test method can totally guarantee their absence, this product must be treated as potentially infectious. Therefore, usual safety procedures should be observed when handling.

- The quantities indicated may be adjusted depending on the titer of the raw materials used.
- Certain cupules contain products of animal origin, notably peptones.

## REAGENTS AND MATERIAL REQUIRED BUT NOT PROVIDED

### Reagents :

- API® Suspension Medium, 2 ml (Ref. 70 700)
- FVB reagent (Ref. 70 500)
- Hydrogen peroxide (3 %)

### Material

- Wooden or glass applicator sticks
- Pipettes or PSIpettes, or precision pipette
- Ampule rack
- Ampule protector
- General microbiology laboratory equipment including Ultra violet lamp (365 nm)

## WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- For professional use only.
- This kit contains products of human origin. No known analysis method can totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions; refer to the Laboratory biosafety manual - 1993, 2nd edition WHO Geneva.
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled observing the usual safety precautions (do not ingest or inhale).
- All specimens, microbial cultures and inoculated products should be considered infectious and handled appropriately. Aseptic technique and usual precautions for handling the bacterial group studied should be observed throughout this procedure. Refer to "CLSI® M29-A, Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline - Current revision". For additional handling precautions, refer to "Biosafety in Microbiological and Biomedical Laboratories - CDC/NIH - Latest edition", or to the regulations currently in use in each country.
- Do not use reagents past the expiration date.
- Before use, check that the packaging of the various components is intact.
- Do not use strips which have been damaged : cupules deformed, desiccant sachet open, etc.

- The performance data presented were obtained using the procedure indicated in this package insert. Any change or modification in the procedure may affect the results.
- Interpretation of the test results should be made taking into consideration the patient history, the source of the specimen, colonial and microscopic morphology of the strain and, if necessary, the results of any other tests performed.

## STORAGE CONDITIONS

The strips should be stored at 2-8°C until the expiration date indicated on the packaging.

## SPECIMENS (COLLECTION AND PREPARATION)

RAPIDEC staph is not for use directly with clinical or other specimens.

The microorganisms to be identified must first be isolated on a suitable culture medium according to standard microbiological techniques.

## INSTRUCTIONS FOR USE

### Preparation of the strip

- Remove the strip from its packaging.
- Write the specimen reference numbers on the strip next to the letters A, B or C.
- Place the lid on the strip.  
If the strip is not completely used, cut off the number of tests required using a pair of scissors. Replace the remainder of the strip in the refrigerator in its packaging and use it within 2 days.

### Preparation of the inoculum

- Open an ampule of API® Suspension Medium (2 ml) as indicated in the paragraph "Warnings and Precautions" of the package insert for this product, or use any tube containing 2 ml of distilled water without additives : dispense 250 µl into each of the cupules C and S. The level of the water must be 1 mm below the top edge of the cupule.
- Gently homogenize the contents of cupule C with an applicator stick. Proceed immediately to the next step.
- Using colonies :  
With the end of a new stick, pick up several colonies of the same morphology and homogenize them in cupule S until a turbidity equivalent to that in cupule C (4 McFarland) is obtained, i.e. about 2 to 5 colonies.  
It is recommended to use young cultures (18-24 hours old).
- Using blood culture broth (detection of *S. aureus* only) :
  - Transfer 2-3 ml of blood culture broth (excluding red blood cells) which has formed a deposit into a cone-shaped tube.
  - Add 2-3 ml of distilled water.
  - Centrifuge for 5-10 minutes (3000 r.p.m.)
  - The centrifugation pellet is used to prepare the bacterial suspension in cupule S in order to obtain a turbidity equivalent to that in cupule C (4 McFarland).

### NOTES :

- Do not use saline medium for rehydration of cupules C and S (it would prevent the appearance of turbidity in cupule C).
- Any spatter of liquid from cupule C into cupules S, 0, 1, 2 or 3 will kill the bacteria and prevent the enzymatic reactions.
- It is easier to compare turbidities by placing the strip on a black background.

### Inoculation of the strip

- Transfer 50 µl of the bacterial suspension (S) into the corresponding cupules 0, 1, 2 and 3.
- Put the excess suspension back into cupule S.
- Withdraw the liquid from cupule C to avoid contamination of the other cupules.
- Place a lid on the strip.
- Incubate the RAPIDEC staph strip for 2 - 2 ¼ hours at 36°C ± 2°C.

## READING AND INTERPRETATION

### Reading the strip

Refer to the Reading Table.

#### • Fluorescent reaction :

Place the strip under an Ultra Violet lamp (365 nm) and read the AUR test.

Only *S. aureus* is AUR (+) and not the other coagulase

(+) veterinary species : *S. intermedius* and *S. hyicus*.

To interpret the AUR test, compare the fluorescence in cupules 1 and 0 (negative control for fluorescence). The AUR test is positive when it is more fluorescent than the control.

#### • If the AUR test is (-) :

- Read the PAL test.

- Reveal the β-GALactosidase test (β-GAL) :

Add 1 drop of FVB reagent to cupule 3. Read the reaction after 2-3 minutes.

#### • Confirmatory test : CATALase test (CAT)

Add 1 drop of hydrogen peroxide (3 %) to cupule 2 (after reading the PAL reaction). Read the reaction after 1-2 minutes.

Cupule 2 (CAT) negative : no bubbles  
positive : **bubbles**.

Staphylococci are catalase (+).

### NOTES :

- The incubation time can exceed 2 - 2 ¼ hours, but should not exceed 4 hours. After this length of time, it is difficult to interpret the reactions and the results could be false.
- Reading is made easier by placing the strip on a white background.
- If RAPIDEC staph is used directly with blood culture broth, **only read the AUR reaction**.

### Interpretation

Refer to the Reading Table.

### Continuation of the analysis

After reading, the bacterial suspension which remains in cupule S can be used for complete identification (API strip), an antimicrobial susceptibility test (ATB™ strip) or a subculture.

## QUALITY CONTROL

The media, strips and reagents are systematically controlled at various stages of their manufacture. For those users who wish to perform their own quality control tests with the strip, it is preferable to use the strain **1. *Staphylococcus aureus* ATCC® 25923** or else the following strain :

**2. *Staphylococcus saprophyticus*** ATCC 43867

ATCC : American Type Culture Collection, 10801 University Boulevard, Manassas, VA 20110-2209, USA.

	Fluorescent reaction	Spontaneous reaction	Reaction requiring a reagent
	Cupule 1	Cupule 2	Cupule 3
	AUR (UV)	PAL	β-GAL (FVB)
<b>1.</b>	<b>+</b>	<b>+</b>	<b>—</b>
<b>2.</b>	<b>—</b>	<b>—</b>	<b>+</b>

It is the responsibility of the user to perform Quality Control in accordance with any local applicable regulations.

## LIMITATIONS OF THE METHOD

- RAPIDEC staph is designed for the identification of the main staphylococci isolated from **human** specimens. **It should only be used for staphylococci** (check morphology, catalase reaction, anaerobic glucose fermentation and susceptibility to lysostaphin).
- The following isolation media can be used : blood agar or Baird-Parker. Media containing lactose (BCP, MacConkey, CLED...) should be used with caution as the β-GAL test for *S. epidermidis*, although negative, can give a slightly pink or red coloration. Take this into account when reading : only a dark pink or red coloration should be considered a positive reaction when these media are used.
- When blood culture broths are used directly, only the "aurease" reaction should be read and interpreted to determine whether *S. aureus* is present or not. The presence of blood products may affect the PAL and β-GAL reactions and therefore their results should not be taken into consideration. Other staphylococcal species should be identified after culture and complete identification.
- Some strains belonging to the veterinary species *S. intermedius* can give a weakly positive aurease reaction due to the presence of coagulase activity.
- Only pure cultures of a single organism should be used.

## RANGE OF EXPECTED RESULTS

Consult the Reading Table at the end of this package insert for the range of expected results for the various biochemical reactions.

## PERFORMANCE

- Results obtained during an evaluation conducted in Europe using colonies (197 strains) :

Species	Number of profiles	Taxon		
		Correct	<i>Staph</i> sp	Incorrect
<i>S. aureus</i>	81	80	0	1
<i>S. epidermidis</i>	85	63	22	0
<i>S. saprophyticus</i>	7	6	1	0
<i>S. xylosus</i>	1	1	0	0
Other staphylococci of human origin *	23	N/A	22	1

\* *S. haemolyticus*, *S. hominis*, *S. warneri*, *S. cohnii*, *S. capitis*

- Direct identification of *S. aureus* using blood culture broths has been evaluated in various studies (see Literature References 1, 10 and 14).

## WASTE DISPOSAL

Dispose of used or unused reagents as well as any other contaminated disposable materials following procedures for infectious or potentially infectious products.

It is the responsibility of each laboratory to handle waste and effluents produced according to their type and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

## WARRANTY

bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.

### READING TABLE

• Using colonies :

Fluorescence control	Fluorescent reaction	Spontaneous reaction	Reaction requiring a reagent	Results
Cupule 0	Cupule 1	Cupule 2	Cupule 3	
	<b>AUR (UV)</b>	<b>PAL</b>	<b>β-GAL (FVB)</b>	
	FLUORESCENT	YELLOW or colorless	██████████	<i>S. aureus</i>
	colorless	YELLOW	PINK / RED	<i>S. intermedius</i> * <i>S. xylosus</i>
			colorless / yellow	<i>S. epidermidis</i>
	colorless	colorless	PINK / RED	<i>S. saprophyticus</i> **
			colorless / yellow	<i>Staphylococcus</i> spp.

(\*) veterinary species.

(\*\*) *S. simulans*, very rare in clinical laboratories, gives the same profile. In case of doubt, inoculate an API® identification strip.

If the combination of colors obtained is not listed in the table, it may be due to a mixture of organisms or to the presence of a species other than *Staphylococcus*. Check the purity of the strain and inoculate an API identification strip.


• Using blood culture broth (detection of *S. aureus* only) :

RAPIDEC staph used directly with blood culture broth : only read the AUR reaction (cupule 1) :

- AUR (UV) negative : no fluorescence ➡ absence of *S. aureus*.
- AUR (UV) positive : fluorescence ➡ presence of *S. aureus*.

PROCEDURE	p. I
LITERATURE REFERENCES	p. II
INDEX OF SYMBOLS	p. III

CLSI is a used, pending and/or registered trademark belonging to Clinical Laboratory and Standards Institute, Inc.  
ATCC is a used, pending and/or registered trademark belonging to American Type Culture Collection.

 **bioMérieux® SA**  
au capital de 12 029 370 €  
673 620 399 RCS LYON  
69280 Marcy-l'Etoile / France  
Tel. 33 (0)4 78 87 20 00  
Fax 33 (0)4 78 87 20 90  
<http://www.biomerieux.com>

**bioMérieux, Inc**  
Box 15969,  
Durham, NC 27704-0969 / USA  
Tel. (1) 919 620 20 00  
Fax (1) 919 620 22 11  
Printed in France

bioMérieux, the blue logo, RAPIDEC, API and ATB are used, pending and/or registered trademarks belonging to bioMérieux SA or one of its subsidiaries.

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



1. CHAPIN K., MUSGNUG M., Evaluation of Three Rapid Methods for the direct identification of *Staphylococcus* from Positive Blood Cultures, (2003) *J. Clin. Microbiol.* 41, 4324-4327.
2. DOUCET M., MONGET D., GUICHERD M. RAPIDEC staph : A New Miniaturized System for the 2-Hour Detection of the Main Staphylococci. (1987) 4-6 Novembre FLORENCE, Fifth International Symposium on Rapid Methods and Automation in Microbiology and Immunology.
3. GEARY C., STEVENS M. A rapid test to detect the most clinically significant *Staphylococcus* species. (1991) *Medical Laboratory Sciences*, 48, 99-105.
4. KLOOS W.E., SCHLEIFER K.H. Simplified Scheme for Routine Identification of Human *Staphylococcus* Species. (1975) *J. Clin. Microbiol.* 1, 82-88.
5. KRIEG N.R., HOLT J.G. Bergey's Manual of Determinative Bacteriology. Eighth Edition. (1984) Williams and Wilkins Co., Baltimore, Md.
6. LE MINOR L., VERON M. Bactériologie Médicale. 2ème édition. (1989) Flammarion Médecine-Sciences. Paris.
7. LOTTENBERG R., CHRISTENSEN U., JACKSON C.M., COLEMAN P.L. Assay of Coagulation Proteases Using Peptide Chromogenic and Fluorogenic Substrates. (1981) *Methods in Enzymology* 80, 341-361
8. MacFADDIN J.F. Biochemical Tests for Identification of Medical Bacteria. Second edition. (1980) Williams and Wilkins, Baltimore, Md.
9. MATHIEU D., PICARD V. Comparative Evaluation of Five Agglutination Techniques and a New Miniaturized System for Rapid Identification of Methicillin-Resistant Strains of *Staphylococcus aureus*. (1991) *Zbl. Bakt.*, 276, 46-53.
10. MITCHELL C.J., GEARY C., STEVENS M. Detection of *Staphylococcus aureus* in Blood Cultures : Evaluation of a Two-Hour Method. (1991) *Medical Laboratory Sciences*, 48, 106-109.
11. MONGET D., DOUCET M., DESMONCEAUX M., GUICHERD M. 2-Hour Detection of the main Staphylococci with the RAPIDEC staph System. (1988) 8-13 May MIAMI, American Society of Microbiology.
12. MURRAY P.R., BARON E.J., PFALLER M.A., TENOVER F.C., YOLKEN R.H. Manual of Clinical Microbiology. 7th Edition. (1999) Amer. Soc. Microbiol., Washington, D.C.
13. PIEMONTE Y., ORSONI A., BOUVET A. et al. Evaluation of RAPIDEC Staph (API) Versus Conventional Staphylococcal Identification Methods. (1989) NICE, 4th European Congress of Clinical Microbiology, ABSTRACT 504/PP23.
14. SPEERS D.J., OLMA T.R., GILBERT G.L. Evaluation of Four Methods for Rapid Identification of *Staphylococcus aureus* from Blood Cultures. (1998) *J. Clin. Microbiol.* 36, 1032-1034.

**TABLE DES SYMBOLES / INDEX OF SYMBOLS / SYMBOLE /  
CUADRO DE SIMBOLOS / TABELLA DEI SIMBOLI / QUADRO DOS SÍMBOLOS / SYMBOLER /  
SYMBOLFORTEGNELSE / TABELA SYMBOLI**

Symbole / Symbol Símbolo / Simbolo	Signification / Meaning / Bedeutung Significado / Significato Betydelse / Betydning / Znaczenie
<b>REF</b>	Référence du catalogue Catalogue number (GB) / Catalog number (US) Bestellnummer / Número de catálogo / Numero di catalogo Referência de catálogo Katalognummer / Katalognummer / Numer katalogowy
<b>IVD</b>	Dispositif médical de diagnostic in vitro In Vitro Diagnostic Medical Device In Vitro Diagnostikum Producto sanitario para diagnóstico in vitro Dispositivo medico-diagnóstico in vitro Dispositivo médico para diagnóstico in vitro Medicintekniska produkter för in vitro diagnostik Medicinsk udstyr til in vitro-diagnostik Wyrób do diagnostyki In Vitro
	Fabricant / Manufacturer / Hersteller / Fabricante Fabbricante / Tillverkare / Producent
	Limites de température / Temperature limitation Temperaturbegrenzung / Límite de temperatura Limiti di temperatura / Limites de temperatura / Temperaturbegränsning Temperaturbegrænsning / Przestrzegać zakresu temperatury
	Utiliser jusque / Use by / Verwendbar bis Fecha de caducidad / Utilizzare entro / Prazo de validade Använd före / Holdbar til / Użyć przed
<b>LOT</b>	Code du lot / Batch code Chargenbezeichnung / Código de lote Codice del lotto / Código do lote / Lot nummer / Lotnummer / Kod partii
	Consulter les instructions d'utilisation Consult Instructions for Use Gebrauchsanweisung beachten Consulte las instrucciones de uso Consultare le istruzioni per l'uso Consulte as instruções de utilização Se handhavandebeskrivningen / Se brugsanvisning Sprawdź w instrukcji obsługi
	Contenu suffisant pour "n" tests Contains sufficient for <n> tests Inhalt ausreichend für <n> Prüfungen Contenido suficiente para <n> ensayos Contenuto sufficiente per "n" saggi Conteúdo suficiente para "n" ensaios Räcker till "n" antal tester Indeholder tilstrækkeligt til "n" test Wystarczy na wykonanie <n> testów