# api<sup>®</sup> 50 CHL Medium

IVD

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*Lactobacillus* and related genera

## SUMMARY AND EXPLANATION

API 50 CHL Medium, intended for the identification of the genus *Lactobacillus* and related genera, is a ready-to-use medium which allows the fermentation of the 49 carbohydrates on the API 50 CH strip to be studied.

## PRINCIPLE

A suspension is made in the medium with the microorganism to be tested and each tube of the strip is then inoculated with the suspension. During incubation, the carbohydrates are fermented to acids which produce a decrease in the pH, detected by the change in color of the indicator. The results make up the biochemical profile which is used by the identification software to identify the strain.

## CONTENT OF THE KIT (Kit for 10 tests)

- 10 ampules of API 50 CHL Medium
- 1 package insert

## COMPOSITION OF THE MEDIUM

<b>API 50 CHL Medium</b> 10 ml	Polypeptone (bovine/porcine origin)	10 g
	Yeast extract	5 g
	Tween 80	1 ml
	Dipotassium phosphate	2 g
	Sodium acetate	5 g
	Diammonium citrate	2 g
	Magnesium sulfate	0.20 g
	Manganese sulfate	0.05 g
	Bromcresol purple	0.17 g
	Demineralized water	1000 ml
	pH : 6.7-7.1	

The quantities indicated may be adjusted depending on the titer of the raw materials used.

## REAGENTS AND MATERIAL REQUIRED BUT NOT PROVIDED

### Reagents


- API 50 CH strips (Ref. 50 300)
- McFarland Standard (Ref. 70 900), point 2 on the scale or DENSIMAT (Ref. 99 234) or ATB™ Densitometer
- **apiweb™** identification software (Ref. 40 011) (consult bioMérieux)
- Mineral oil (Ref. 70 100)
- MRS agar (Ref. 42 602)
- API Suspension Medium, 2 ml (Ref. 70 700) and 5 ml (Ref. 20 150)

### Material

- Pipettes or PSlpettes
- Ampule rack
- Small and large ampule protectors
- Swabs
- General microbiology laboratory equipment


## WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use and microbiological control.**
- **For professional use only.**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled observing the usual safety precautions (do not ingest or inhale).
- All specimens, microbial cultures and inoculated products should be considered infectious and handled appropriately. Aseptic technique and usual precautions for handling the bacterial group studied should be observed throughout this procedure. Refer to "CLSI<sup>®</sup> M29-A, *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline - Current revision*". For additional handling precautions, refer to "Biosafety in Microbiological and Biomedical Laboratories - CDC/NIH - Latest edition", or to the regulations currently in use in each country.
- Do not use media past the expiry date.
- Before use, check that the ampules are intact.
- Allow media to come to room temperature before use.
- Open ampules carefully as follows :
  - Place the ampule in the ampule protector.
  - Hold the protected ampule in one hand in a vertical position (white plastic cap uppermost).
  - Press the cap down as far as possible.

\* Model 1 :

- Cover the flattened part of the cap with the upper part of the thumb.
- Apply thumb pressure to the base of the flattened part of the cap to snap off the top of the ampule.

\* Model 2 :

- Position the thumb tip on the striated part of the cap and press forward to snap off the top of the ampule.

- Take the ampule out of the ampule protector and put the protector aside for subsequent use.
- Carefully remove the cap.
- The performance data presented were obtained using the procedure indicated in this package insert. Any change or modification in the procedure may affect the results.
- Interpretation of the test results should be made taking into consideration the patient history, the source of the specimen, colonial and microscopic morphology of the strain and, if necessary, the results of any other tests performed, particularly the antimicrobial susceptibility patterns.

**STORAGE CONDITIONS**

The media should be stored at 2-8°C until the expiry date indicated on the packaging.

**SPECIMENS (COLLECTION AND PREPARATION)**

API 50 CHL Medium is not for use directly with clinical or other specimens.

The microorganisms to be identified must first be isolated on a suitable culture medium according to standard microbiological techniques.

**INSTRUCTIONS FOR USE**

**Selection of the colonies**

- Check the purity of the strain.
- Culture it on MRS agar, and incubate anaerobically for 24 hours at 30 or 37°C. The temperature of incubation depends on the origin of the strain.
- Check that it belongs to the group of lactic bacteria : Gram (+), catalase (-), non-sporeforming, anaerobic (facultative or occasionally obligate) bacteria, which grow on MRS agar.
- If freeze-dried or frozen strains are used, subculture twice in MRS broth before isolation on MRS agar.

**Preparation of the strip**

See the API 50 CH package insert.

**Preparation of the inoculum**

- If the DENSIMAT or ATB™ Densitometer is used :
  - Open an ampule of API 50 CHL Medium as indicated in the paragraph "Warnings and Precautions".
  - Pick up several identical colonies.
  - Prepare a suspension with a turbidity equivalent to 2 McFarland in the ampule of API 50 CHL Medium. This suspension must be used immediately after preparation.
- If the DENSIMAT or ATB Densitometer is not used :
  - Open an ampule of API Suspension Medium (2 ml) as indicated in the paragraph "Warnings and Precautions" or use any tube containing sterile distilled water without additives.
  - Pick up all the bacteria from the culture using a swab.
  - Prepare a heavy suspension (S) in the ampule.

- Open an ampule of API Suspension Medium (5 ml) as indicated in the paragraph "Warnings and Precautions".
- Prepare a suspension with a turbidity equivalent to 2 McFarland by transferring a certain number of drops of suspension S into the ampule : record this number of drops (n).
- Open an ampule of API 50 CHL Medium as indicated in the paragraph "Warnings and Precautions" and inoculate by transferring twice the number of drops of suspension S (i.e. 2n) into the ampule. This suspension must be used immediately after preparation.

- Homogenize.

**Inoculation of the strip**

- Fill the tubes (not the cupules) with the inoculated API 50 CHL Medium, and cover all the tests with mineral oil.
- Incubate aerobically at 29°C ± 2°C or 36°C ± 2°C, for 48 hours (± 6 hours).

**READING AND INTERPRETATION**

**Reading the strip**

(see the API 50 CH package insert)

- Read after 48 hours of incubation.
- A positive test corresponds to acidification revealed by the bromcresol purple indicator contained in the medium changing to YELLOW. For the esculin test (tube no. 25), a change in color from purple to BLACK is observed.
- Record the results on the result sheet.

**Interpretation**

The biochemical profile obtained for the strain can be identified using the **apiweb™** identification software with database (V5.1).

**NOTE :**

The biochemical profile may also be used with other results for a taxonomic study.

**QUALITY CONTROL**

The media and strips are systematically controlled at various stages of their manufacture. For those users who wish to perform their own quality control tests with the strip, it is preferable to use the strain **1. *Lactobacillus plantarum* ATCC® 14917** or the following strain :

- 2. *Lactobacillus paracasei* ssp *paracasei* NCFB 206 or ATCC BAA-52

ATCC : American Type Culture Collection, 10801 University Boulevard, Manassas, VA 20110-2209, USA.

NCFB : National Collection for Food Bacteria (=NCDO), Institute of Food Research, Reading Laboratory, Earley Gate, Reading RG6 6BZ, ENGLAND

	0	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49		
1.	24	-	-	-	-	+	+	-	-	-	-	+	+	+	+	-	-	-	-	+	-	-	-	+	+	+	+	+	+	+	-	+	+	+	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
	48	-	-	-	-	+	+	-	-	-	-	+	+	+	+	-	-	-	-	+	+	+	+	+	+	+	+	+	V	+	+	+	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
2.	24	-	-	-	-	+	-	-	-	-	+	+	+	+	-	-	V	-	+	V	-	+	+	-	V	+	V	V	+	-	-	+	+	+	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	48	-	-	-	-	+	-	-	-	-	+	+	+	+	-	-	+	+	+	+	+	+	+	+	+	+	+	+	+	+	-	-	+	+	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-

It is the responsibility of the user to perform Quality Control in accordance with any local applicable regulations.

**LIMITATIONS OF THE METHOD**

- The API 50 CHL system is intended uniquely for the identification of those species included in the database (see Identification Table at the end of this package insert). It cannot be used to identify any other microorganisms or to exclude their presence.
- Only pure cultures of a single organism should be used.

**RANGE OF EXPECTED RESULTS**

Consult the Identification Table at the end of this package insert for the range of expected results for the various biochemical reactions.

**PERFORMANCE**

944 collection strains and strains of various origins belonging to species included in the database were tested :

- 81.36% of the strains were correctly identified (with or without supplementary tests).
- 6.99% of the strains were not identified.
- 11.65% of the strains were misidentified.

**WASTE DISPOSAL**

Dispose of used or unused reagents as well as any other contaminated disposable materials following procedures for infectious or potentially infectious products.

It is the responsibility of each laboratory to handle waste and effluents produced according to their type and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

**WARRANTY**

bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.

PROCEDURE	p. I
IDENTIFICATION TABLE	p. II
LITERATURE REFERENCES	p. III
INDEX OF SYMBOLS	p. IV

CLSI is a used, pending and/or registered trademark belonging to Clinical Laboratory and Standards Institute, Inc.  
ATCC is a used, pending and/or registered trademark belonging to American Type Culture Collection.

**bioMérieux® SA**

au capital de 12 029 370 €  
673 620 399 RCS LYON  
69280 Marcy-l'Etoile / France  
Tel. 33 (0)4 78 87 20 00  
Fax 33 (0)4 78 87 20 90  
<http://www.biomerieux.com>

**bioMérieux, Inc**

Box 15969,  
Durham, NC 27704-0969 / USA  
Tel. (1) 919 620 20 00  
Fax (1) 919 620 22 11

Printed in France

bioMérieux, the blue logo, API, ATB and **apiweb** are used, pending and/or registered trademarks belonging to bioMérieux SA or one of its subsidiaries.

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