ID 32 STAPH

IVD

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ID 32 STAPH

IVD

Identification system for staphylococci

SUMMARY AND EXPLANATION

ID 32 STAPH is a standardized system for the identification of the genera *Staphylococcus*, *Micrococcus* and related genera, *Rothia* and *Aerococcus*, which uses 26 miniaturized biochemical tests, as well as a specific database. The complete list of organisms that it is possible to identify with this system can be found in the Identification Table at the end of this package insert. Reading and interpretation are carried out automatically or manually.

PRINCIPLE

The ID 32 STAPH strip consists of 32 cupules, 26 of which are used as test cupules and contain dehydrated test substrates.

After 24 hours of incubation, the reactions are read either using the ATB™ Expression™ or *mini API*® instruments, or visually.

Identification is obtained using the identification software.

CONTENT OF THE KIT (Kit for 25 tests) :

- 25 ID 32 STAPH strips
- 25 incubation lids
- 1 package insert

COMPOSITION OF THE STRIP

The composition of the ID 32 STAPH strip is given in the Reading Table of this package insert.

REAGENTS AND MATERIAL REQUIRED BUT NOT PROVIDED

Reagents / Instrumentation

- API Suspension Medium, 2 ml (Ref. 70 700) or 3 ml (Ref. 70 640) if the ATB Inoculator is used
- Reagents : VP A + VP B (Ref. 70 572)
NIT 1 + NIT 2 (Ref. 70 442)
FB (Ref. 70 562)
- Mineral oil (Ref. 70 100)
- DENSIMAT (Ref. 99 234) or ATB Densitometer or McFarland Standard (Ref. 70 900), point 0.5 on the scale
- ATB Electronic Pipette (consult bioMérieux) or ATB Inoculator and Tips (Ref. 15 710)
- ATB Expression or *mini API*, or *apiweb*™ identification software (Ref. 40 011) (consult bioMérieux)

Material

- Ampule rack
- Ampule protector
- Air-tight box
- General microbiology laboratory equipment

WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use and microbiological control.**
- **For professional use only.**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled observing the usual safety precautions (do not ingest or inhale).
- All specimens, microbial cultures and inoculated products should be considered infectious and handled appropriately. Aseptic technique and usual precautions for handling the bacterial group studied should be observed throughout this procedure. Refer to "CLSI/NCCLS M29-A, *Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids, and Tissue; Approved Guideline* - Current revision". For additional handling precautions, refer to "Biosafety in Microbiological and Biomedical Laboratories - CDC/NIH - Latest edition", or to the regulations currently in use in each country.
- Do not use reagents past the expiry date.
- Before use, check that the packaging of the various components is intact.
- Do not use strips which have been damaged : cupules deformed, desiccant sachet open, etc.
- The performance data presented were obtained using the procedure indicated in this package insert. Any change or modification in the procedure may affect the results.
- Interpretation of the test results should be made taking into consideration the patient history, the source of the specimen, colonial and microscopic morphology of the strain and, if necessary, the results of any other tests performed, particularly the antimicrobial susceptibility patterns.

STORAGE CONDITIONS

The strips should be stored at 2-8°C until the expiry date indicated on the packaging.

SPECIMENS (COLLECTION AND PREPARATION)

ID 32 STAPH is not for use directly with clinical or other specimens.

The microorganisms to be identified must first be isolated on a suitable culture medium according to standard microbiological techniques.

INSTRUCTIONS FOR USE

Selection of the colonies

- Harvest the colonies obtained on blood agar, Mannitol Salt agar or Baird Parker.

NOTE: MacConkey may also be used. This culture medium is nevertheless not very suitable for the growth of staphylococci and micrococci : the growth of certain species may be inhibited, particularly if the medium contains crystal violet.

Preparation of the strip

- Remove the strip from its packaging.
- Discard the desiccant.
- Place the lid on the strip.
- Record the strain reference on the elongated flap of the strip. (Do not record the reference on the lid as it may be misplaced during the procedure.)

Preparation of the inoculum

- Open an ampule of API® Suspension Medium, 2 ml (3 ml if the ATB™ Inoculator is used) as indicated in the paragraph "Warnings and Precautions" of the API Suspension Medium package insert, or use any tube containing sterile distilled water without additives.
- Pick up several identical colonies. It is recommended to use young cultures.
- Prepare a suspension with a turbidity equivalent to 0.5 McFarland : measure with the ATB Densitometer, or the DENSIMAT, or compare with a turbidity control (McFarland Standard). This suspension must be used immediately after preparation.

NOTE : If the strip is to be read AUTOMATICALLY, the ATB Densitometer or the DENSIMAT must be used to adjust the turbidity of the bacterial suspension.

Inoculation of the strip

- **AUTOMATIC inoculation** :
 - Place the strip, the inoculated ampule of API Suspension Medium and a Tip on the ATB Inoculator tray.
 - The inoculator will automatically homogenize the contents of the ampule and fill the cupules (55 µl / cupule).
- **MANUAL inoculation** :
 - Homogenize the ampule of inoculated API Suspension Medium and dispense 55 µl of the suspension into each cupule of the strip using the ATB Electronic Pipette.
- Cover the tests URE, ADH and ODC with **2 drops** of mineral oil (cupules 1.0, 1.1 and 1.2).
- Place the lid on the strip.
- Incubate at 36°C ± 2°C for 24 hours (± 2 hours) in **aerobic** conditions.

NOTE : Some ventilated incubators may completely dehydrate the medium in the cupules. In this case, place the strip in an air-tight box with a receptacle containing a small volume of water. A humid atmosphere is thus created, which prevents the tests from drying out.

READING AND INTERPRETATION

Reading the strip

Reveal all the reactions in row 0 by adding 1 drop of the following reagents:

- NIT test (cupule 0.0) : NIT 1 and NIT 2 reagents.
- VP test (cupule 0.1) : VP A and VP B reagents.
- Tests βGAL to PyrA (cupules 0.2 to 0.5) : FB reagent.

Read after 5 minutes (do not exceed 10 minutes) :

- **AUTOMATIC reading** using ATB Expression™ or **mini API** :
 - check that the middle part of the strip is clean so that the reader can recognize the strip code.
 - The reader records the color for each cupule and transmits the information to the computer.
- **VISUAL reading** :
 - refer to the Reading Table. Record the results on the result sheet.


NOTE : In certain cases, the VP reaction becomes even stronger after the 10 minute period indicated ; this test can be read again after 12 minutes.

Interpretation

Identification is obtained using the database (V2.1) :

- **AFTER AUTOMATIC READING** :
 - the results transmitted to the computer are interpreted by the ATB Expression or **mini API** identification software.
 - Check that the name printed on the strip corresponds to the strip name displayed by the software.
- **AFTER VISUAL READING** :
 - the reactions obtained are coded into a **numerical profile** :
 - On the result sheet, the tests are separated into groups of 3 and a number 1, 2 or 4 is indicated for each. The values corresponding to positive reactions are then added together within each group.

Identification is obtained using the **apiweb™** identification software by manually entering the 9-digit numerical profile : the 4 digits from the upper row (1.0-1.B), followed by the 4 digits from the lower row (0.0-0.B), and completed by the 9th digit for the tests RIB and CEL (1.C and 1.D).

2661 2461 0 Staphylococcus haemolyticus

QUALITY CONTROL

The strips are systematically controlled at various stages of their manufacture. For those users who wish to perform their own quality control tests with the strip, it is preferable to use the strain **1. Staphylococcus gallinarum ATCC® 700401** or else one of the following strains :

2. *Staphylococcus lugdunensis*

ATCC 49576

3. *Staphylococcus aureus*

ATCC 29213

ATCC : American Type Culture Collection, 10801 University Boulevard, Manassas, VA 20110-2209, USA.

	URE	ADH	ODC	ESC	GLU	FRU	MNE	MAL	LAC	TRE	MAN	RAF	RIB	CEL	NIT	VP	βGAL	ArgA	PAL	PyrA	NOVO	SAC	NAG	TUR	ARA	βGUR
1.	+	-	-	+	+	+	+	+	+	+	+	+	-	+	+	-	+	-	+	V	+	+	+	+	+	+
2.	V	-	+	-	+	+	+	+	+	+	-	-	-	-	+	+	-	-	-	+	-	+	+	-	-	-
3.	V	+	-	-	+	+	+	+	+	+	+	-	-	-	+	+	-	-	+	V	-	+	+	+	-	-

Profiles obtained after culture of the strains on Columbia agar + 5 % sheep blood and automatic reading of the results.

It is the responsibility of the user to perform Quality Control in accordance with any local applicable regulations.

RECOMMENDATIONS

To obtain the best results with the ID 32 STAPH strip, it is important to scrupulously respect the following points of the procedure :

- Use the isolation media recommended in this package insert (see paragraph "Selection of the colonies").
- Precisely adjust the inoculum to 0.5 McFarland. The ATB Densitometer or the DENSIMAT must be used if the strip is to be read and interpreted by ATB™ Expression™ or *mini API*®.
- Dispense exactly 55 µl per cupule with the ATB Electronic Pipette or the ATB Inoculator (essential if the strip is to be read and interpreted by ATB Expression or *mini API*).
- Respect the incubation time and the reading time.
- If, after automatic reading of the strip, the result of the VP test is against the identification proposed and is negative, read the VP test again visually : if the reaction is positive or weakly positive, manually enter the profile with a positive result for the VP reaction.
- The reagents should be of good quality : check the expiry date and storage conditions and use within one month of opening the ampules.

LIMITATIONS OF THE METHOD

- The ID 32 STAPH system is intended uniquely for the identification of the species included in the database (see Identification Table at the end of this package insert). It cannot be used to identify any other microorganisms or to exclude their presence.
- Only pure cultures of a single organism should be used.

RANGE OF EXPECTED RESULTS

Consult the Identification Table at the end of this package insert for the range of expected results for the various biochemical reactions.

PERFORMANCE

2312 collection strains and strains of various origins belonging to species included in the database were tested :

- 95.81 % of the strains were correctly identified (with or without supplementary tests).
- 2.98 % of the strains were not identified.
- 1.21 % of the strains were misidentified.

WASTE DISPOSAL

Dispose of used or unused reagents as well as any other contaminated disposable materials following procedures for infectious or potentially infectious products.

It is the responsibility of each laboratory to handle waste and effluents produced according to their type and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

WARRANTY

bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.

READING TABLE

CUPULE	TEST	ACTIVE INGREDIENTS	QTY (mg/cup.)	REACTIONS/ENZYMES	RESULT	
					NEGATIVE	POSITIVE
1.0	URE	urea	1.12	UREase	yellow	orange red-violet
1.1	ADH	L-arginine	0.76	Arginine DiHydrolase	yellow	orange-red
1.2	ODC	L-ornithine		Ornithine DeCarboxylase		
1.3	ESC	esculin ferric citrate	0.224 0.032	Hydrolysis (ESCulin)	colorless- pale grey	brown-black
1.4	GLU	D-glucose	0.56	fermentation (GLUcose)	red red-orange	yellow yellow-orange
1.5	FRU	D-fructose	0.56	fermentation (FRUctose)		
1.6	MNE	D-mannose	0.56	fermentation (MaNnosE)		
1.7	MAL	D-maltose	0.56	fermentation (MALtose)		
1.8	LAC	D-lactose (bovine origin)	0.56	fermentation (LACtose)		
1.9	TRE	D-trehalose	0.56	fermentation (TREhalose)		
1.A	MAN	D-mannitol	0.56	fermentation (MANnitol)		
1.B	RAF	D-raffinose	0.56	fermentation (RAFfinose)		
1.C	RIB	D-ribose	0.56	fermentation (RIBose)		
1.D	CEL	D-cellobiose	0.56	fermentation (CELLobiose)		
1.E 1.F	Empty cupules					
0.0	NIT	potassium nitrate	0.054	reduction (NITrates)	<u>NIT 1 + NIT 2 / 5 min < 10 min</u> colorless pink-purple	
0.1	VP	sodium pyruvate	0.475	acetoin production (Voges Proskauer)	<u>VP A + VP B / 10 min < 12 min</u> colorless pink-red	
0.2	βGAL	2-naphthyl-βD-galactopyranoside	0.0364	β GALactosidase	<u>FB / 5 min < 10 min (βGAL → PyrA)</u> colorless purple pale purple pale orange	
0.3	ArgA	L-arginine β-naphthylamide	0.0172	Arginine Arylamidase	colorless pale orange	orange
0.4	PAL	2-naphthyl phosphate	0.0123	Alkaline Phosphatase	colorless pale purple pale orange	purple
0.5	PyrA	pyroglutamic acid- β-naphthylamide	0.0128	Pyrrolidonyl Arylamidase	colorless pale orange	orange
0.6	NOVO	novobiocin	0.0018	resistance (NOVObiocin)	red red-orange	yellow yellow-orange
0.7	SAC	D-saccharose (sucrose)	0.56	fermentation (SACcharose)		
0.8	NAG	N-acetyl-glucosamine	0.56	fermentation (N-Acetyl- Glucosamine)		
0.9	TUR	D-turanose	0.56	fermentation (TURanose)		
0.A	ARA	L-arabinose	0.56	fermentation (ARAbinose)		
0.B	βGUR	4-nitrophenyl-βD-glucuronide	0.0158	β GlucURonidase	colorless	yellow
0.C 0.D 0.E 0.F	Empty cupules					

- The quantities indicated may be adjusted depending on the titer of the raw materials used.
- Certain cupules contain products of animal origin, notably peptones.

PROCEDURE	p. I
IDENTIFICATION TABLE	p. II
LITERATURE REFERENCES	p. III
INDEX OF SYMBOLS	p. IV
RESULT SHEET	p. V

ATCC is a used, pending and/or registered trademark belonging to American Type Culture Collection.

bioMérieux® SA
 au capital de 12 029 370 €
 673 620 399 RCS LYON
 69280 Marcy-l'Etoile / France
 Tel. 33 (0)4 78 87 20 00
 Fax 33 (0)4 78 87 20 90
<http://www.biomerieux.com>

bioMérieux, Inc
 Box 15969,
 Durham, NC 27704-0969 / USA
 Tel. (1) 919 620 20 00
 Fax (1) 919 620 22 11

Printed in France

bioMérieux, the blue logo, API, ATB, Expression and **apiweb** are used, pending and/or registered trademarks belonging to bioMérieux SA or one of its subsidiaries.

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