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Identification system for *Enterobacteriaceae* and other non-fastidious Gram-negative rods

SUMMARY AND EXPLANATION

API 10 S is a standardized identification system for *Enterobacteriaceae* and other non-fastidious Gram-negative rods, which uses 11 miniaturized biochemical tests and a database. The complete list of those organisms that it is possible to identify with this system is given in the Identification Table at the end of this package insert.

PRINCIPLE

The API 10 S strip consists of 10 microtubes containing dehydrated substrates. These tests are inoculated with a bacterial suspension which reconstitutes the media. During incubation, metabolism produces color changes that are either spontaneous or revealed by the addition of reagents.

The reactions are read according to the Reading Table and identification is obtained by consulting the list of profiles in this package insert or using the identification software.

CONTENT OF THE KIT (Kit for 50 tests)

- 50 API 10 S strips
- 50 incubation boxes
- 50 result sheets
- 1 clip seal
- 1 package insert

COMPOSITION OF THE STRIP

The composition of the API 10 S strip is given in the Reading Table of this package insert.

REAGENTS AND MATERIAL REQUIRED BUT NOT PROVIDED

Reagents

- API NaCl 0.85 % Medium, 5 ml (Ref. 20 230) or API Suspension Medium, 5 ml (Ref. 20 150)
- Reagents : TDA (Ref. 70 402)
JAMES (Ref. 70 542)
NIT 1 + NIT 2 (Ref. 70 442)
- Oxidase (Ref. 55 635*)
* reference not sold in certain countries : use an equivalent reagent.
- Mineral oil (Ref. 70 100)
- **apiweb**[™] Identification software (Ref. 40 011)
(consult bioMérieux)

Material

- Pipettes or PSlpettes
- Ampule protector
- Ampule rack
- General microbiology laboratory equipment

WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use and microbiological control.**
- **For professional use only.**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled observing the usual safety precautions (do not ingest or inhale).
- All specimens, microbial cultures and inoculated products should be considered infectious and handled appropriately. Aseptic technique and usual precautions for handling the bacterial group studied should be observed throughout this procedure. Refer to "CLSI/NCCLS M29-A, *Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids, and Tissue; Approved Guideline* - Current revision". For additional handling precautions, refer to "Biosafety in Microbiological and Biomedical Laboratories- CDC/NIH - Latest edition", or to the regulations currently in use in each country.
- Do not use reagents past the expiry date.
- Before use, check that the packaging of the various components is intact.
- Do not use strips which have been damaged : cupules deformed, desiccant sachet open, etc.
- The performance data presented were obtained using the procedure indicated in this package insert. Any change or modification in the procedure may affect the results.
- Interpretation of the test results should be made taking into consideration the patient history, the source of the specimen, colonial and microscopic morphology of the strain and, if necessary, the results of any other tests performed, particularly the antimicrobial susceptibility patterns.

STORAGE CONDITIONS

The API 10 S strips are supplied in an aluminum pouch with desiccant sachets.

Once opened (*), the pouch should be re-sealed using the clip seal (included in the kit) to preserve the remaining strips with the desiccant sachets : place the open end of the pouch along the seal and carefully clamp between the two parts. The strips may then be kept for up to **10 months after the pouch has been opened**, at 2-8°C (or until the expiry date indicated on the packaging, if this comes before).

(* *Recommended method for opening the pouches* : cut open the pouch just below the seal while holding the pouch upright, in order to avoid damaging the desiccant sachets.

SPECIMENS (COLLECTION AND PREPARATION)

API 10 S is not for use directly with clinical or other specimens.

The microorganisms to be identified must first be isolated on a suitable culture medium according to standard microbiological techniques.

INSTRUCTIONS FOR USE

Oxidase test

The oxidase test must be performed according to the manufacturer's instructions for use. The result should be recorded on the result sheet as it is an integral part of the final profile (11th identification test).

Preparation of the strip

- Prepare an incubation box (tray and lid) and distribute about 3 ml of distilled water or demineralized water [or any water without additives or chemicals which may release gases (e.g. Cl₂, CO₂, etc.)] into the honey-combed wells of the tray to create a humid atmosphere.
- Record the strain reference on the elongated flap of the tray. (Do not record the reference on the lid as it may be misplaced during the procedure.)
- Remove the strip from its packaging.
- Place the strip in the incubation box.

NOTE : API 10 S should only be used with *Enterobacteriaceae* and/or non-fastidious Gram-negative rods. Fastidious organisms having demanding nutritional requirements and requiring appropriate handling precautions (i.e. *Brucella* and *Francisella*) are not included in the API 10 S database. Alternative procedures must be used to exclude or confirm their presence.

Preparation of the inoculum

- Open an ampule of API NaCl 0.85 % Medium (5 ml) or an ampule of API Suspension Medium (5 ml) as indicated in the paragraph "Warnings and Precautions" of the package insert for these products, or use any tube containing 5 ml of sterile saline or sterile distilled water, without additives.
- Using a pipette or PSipette, remove a single well-isolated colony from an isolation plate. It is recommended to use young cultures (18-24 hours old).
- Carefully emulsify to achieve a homogeneous bacterial suspension.
This suspension must be used immediately after preparation.

NOTE : most *Vibrio* species are halophilous. If a *Vibrio* is suspected, suspend the bacteria in API NaCl 0.85 % Medium.

Inoculation of the strip

- With the same pipette, distribute the bacterial suspension into the tubes of the strip (to avoid the formation of bubbles at the base of the tubes, tilt the strip slightly forwards and place the tip of the pipette or PSipette against the side of the cupule):
 - For the [CIT] test, fill both tube and cupule ,
 - For the other tests, fill only the tubes (and not the cupules),
 - for the tests LDC, ODC, H₂S and URE, create anaerobiosis by overlaying with mineral oil.
- Close the incubation box.
- Incubate at 36°C ± 2°C for 18-24 hours.

READING AND INTERPRETATION

Reading the strip


- After the incubation period, read the strip by referring to the Reading Table.
- Record all spontaneous reactions on the result sheet.
- Reveal the tests which require the addition of reagents :
 - TDA Test : add 1 drop of TDA reagent. A **reddish brown** color indicates a **positive** reaction to be recorded on the result sheet.
 - IND Test : add 1 drop of JAMES reagent. A **pink** color developed in the whole cupule indicates a **positive** reaction to be recorded on the result sheet.
 - NO₂ Test : add 1 drop each of NIT 1 and NIT 2 reagents to the GLU tube. Wait 2 to 5 minutes. A **red** color indicates a **positive** reaction to be recorded on the result sheet.

NOTE : The nitrate reduction and indole production tests must be performed last since these reactions release gaseous products which interfere with the interpretation of other tests on the strip. The plastic incubation lid should not be replaced after the addition of these reagents.

Interpretation

Identification is obtained with the **numerical profile**.

- Determination of the numerical profile :
On the result sheet, the tests are separated into groups of three and a value 1, 2 or 4 is indicated for each. By adding together the values corresponding to positive reactions within each group, a 4-digit numerical profile is obtained for the 10 tests of the API 10 S strip. (The oxidase reaction constitutes the 11th test and reduction of nitrates to nitrites (NO₂) the 12th test).
- Identification
This is performed using the database (V3.1)
 - * with the numerical profile :
 - Look up the numerical profile in the list of profiles in this package insert.
 - * with the **apiweb™** identification software :
 - Enter the 4-digit numerical profile manually via the keyboard.

7 504 *Klebsiella pneumoniae* ssp *pneumoniae*

Further tests may be proposed in case of low discrimination. Refer to the identification software.

QUALITY CONTROL

The media, strips, and reagents are systematically quality controlled at various stages of their manufacture. For those users who wish to perform their own quality control tests with the strip, it is preferable to use the strain

1. *Escherichia coli* ATCC® 25922 or else one of the following strains :

- | | | | |
|--------------------------------|------------|--|------------|
| 2. <i>Enterobacter cloacae</i> | ATCC 13047 | 4. <i>Stenotrophomonas maltophilia</i> | ATCC 51331 |
| 3. <i>Proteus mirabilis</i> | ATCC 35659 | | |

ATCC : American Type Culture Collection, 10801 University Boulevard, Manassas, VA 20110-2209, USA.

	ONPG	GLU	ARA	LDC	ODC	CIT	H ₂ S	URE	TDA	IND	OX	NO ₂
1.	+	+	+	+	+	-	-	-	-	+	-	+
2.	+	+	+	-	V	+	-	-	-	-	-	+
3.	-	+	-	-	+	V	+	+	+	-	-	+
4.	+	-	-	V	-	V	-	-	-	-	-	-

Profiles obtained after 18-24 hours of incubation, after culture of the strains on Trypticase Soy agar + sheep blood

It is the responsibility of the user to perform Quality Control in accordance with any local applicable regulations.

LIMITATIONS OF THE METHOD

- The API 10 S system is designed uniquely for the identification of *Enterobacteriaceae* and non-fastidious, Gram-negative rods included in the database (see Identification Table at the end of this package insert) which may be encountered in clinical specimens.
- Supplementary tests are sometimes necessary to differentiate between two species. To choose between these two species, it is necessary to take into consideration the patient history, the source of the specimen, colonial and microscopic morphology of the strain and, if necessary, the results of any other tests performed, as well as the identification percentages (% ID) and typicity (T) index (e.g. : for a clinical specimen, in case of poor discrimination between *Escherichia coli* and *Serratia odorifera*, *Escherichia coli* will be indicated, particularly if the % ID and T index values are high and in favor of this species, since *E. coli* are frequently encountered in clinical specimens whereas *S. odorifera* are very rare). Some identifications may be extended by use of the API 20 E strip which provides 10 extra tests compared to the API 10 S strip. In this case, the microorganisms can be identified using the API 20 E Analytical Profile Index.
- Only pure cultures of a single organism should be used.

RANGE OF EXPECTED RESULTS

Consult the Identification Table at the end of this package insert for the range of expected results for the various biochemical reactions.

PERFORMANCE

- *Enterobacteriaceae*
4775 collection strains and strains of various origins belonging to species included in the database were tested :
 - 94.01 % of the strains were correctly identified (with or without supplementary tests).
 - 2.03 % of the strains were not identified.
 - 3.96 % of the strains were misidentified.
- Other non-fastidious Gram-negative rods :
1286 collection strains and strains of various origins belonging to species included in the database were tested :
 - 95.72 % of the strains were correctly identified (with or without supplementary tests).
 - 1.48 % of the strains were not identified.
 - 2.80 % of the strains were misidentified.

WASTE DISPOSAL

Dispose of used or unused reagents as well as any other contaminated disposable materials following procedures for infectious or potentially infectious products.

It is the responsibility of each laboratory to handle waste and effluents produced according to their type and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

WARRANTY

bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.

READING TABLE

TESTS	ACTIVE INGREDIENTS	QTY (mg/cup.)	REACTIONS/ENZYMES	RESULTS	
				NEGATIVE	POSITIVE
ONPG	2-nitrophenyl-βD-galactopyranoside	0.223	β-galactosidase (Ortho-NitroPhenyl-βD-Galactopyranosidase)	colorless	yellow (1)
GLU	D-glucose	1.9	fermentation / oxidation (GLUcose) (3)	blue / blue-green	yellow / yellow-grey
ARA	L-arabinose	1.9	fermentation / oxidation (ARAbinose) (3)	blue / blue-green	yellow
<u>LDC</u>	L-lysine	1.9	Lysine DeCarboxylase	yellow	red / orange
<u>ODC</u>	L-ornithine	1.9	Ornithine DeCarboxylase	yellow	red / orange
<u>CIT</u>	trisodium citrate	0.756	CITrate utilization	pale green / yellow	blue-green / blue (2)
<u>H₂S</u>	sodium thiosulfate	0.075	H ₂ S production	colorless / greyish	black deposit / thin line
<u>URE</u>	urea	0.76	UREase	yellow	red / orange
TDA	L-tryptophane	0.38	Tryptophane DeAminase	TDA / immediate yellow reddish brown	
IND	L-tryptophane	0.19	INDole production	JAMES / immediate colorless pink pale green / yellow	
OX	(see oxidase test package insert)	-	cytochrome-OXidase	(see oxidase test package insert)	
NO ₂	(GLU tube)	-	NO ₂ production	NIT 1 + NIT 2 / 2-5 min yellow red	

(1) A very pale yellow should also be considered positive.

(2) Reading made in the cupule (aerobic).

(3) Fermentation begins in the lower portion of the tubes, oxidation begins in the cupule.

- The quantities indicated may be adjusted depending on the titer of the raw materials used.
- Certain cupules contain products of animal origin, notably peptones.

IDENTIFICATION TABLE	p. I
LIST OF NUMERICAL PROFILES	p. II
LITERATURE REFERENCES	p. IV
INDEX OF SYMBOLS	p. V

ATCC is a used, pending and/or registered trademark belonging to American Type Culture Collection.

bioMérieux® SA
 au capital de 12 029 370 €
 673 620 399 RCS LYON
 69280 Marcy-l'Etoile / France
 Tel. 33 (0)4 78 87 20 00
 Fax 33 (0)4 78 87 20 90
<http://www.biomerieux.com>

bioMérieux, Inc
 Box 15969,
 Durham, NC 27704-0969 / USA
 Tel. (1) 919 620 20 00
 Fax (1) 919 620 22 11

Printed in France

bioMérieux, the blue logo, API and **apiweb** are used, pending and/or registered trademarks belonging to bioMérieux SA or one of its subsidiaries.

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1. MEHTAR S., AFSHAR S.A.