# CRM

Cardiac Rhythm Management

External Devices

Technical Manual

Manual técnico

Manual técnico

## Reliatty

Measuring device for analysis of lead functions during the implantation of pacemakers and ICDs

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**BIOTRONIK**  
excellence for life

# Table of Contents

<b>Safety When Using the Device</b> .....	<b>5</b>
Intended Medical Use .....	6
Required Expertise .....	8
Contraindications .....	9
Possible Complications .....	9
Installation – Setup Requirements .....	10
Maintenance, Service, Disposal .....	16
General Safety Instructions .....	19
<b>Handling of the Device</b> .....	<b>25</b>
Device Overview .....	26
Transport and Setup .....	29
Power Supply .....	30
Switching On and Off .....	35
Connections and Cables .....	36
Keys, Displays and Signals .....	45
<b>Using the Software</b> .....	<b>49</b>
Overview of the Software .....	50
Handling the Operating Elements .....	52
Measuring P/R Wave Amplitudes .....	60
Measuring Pacing Threshold and Lead Impedance ..	62
Measuring Conduction Times .....	64
Measuring the Wenckebach Point .....	67
Printing/Exporting a Report .....	70
Delivering Burst Pacing .....	74
Delivering Maximum Pacing .....	77
Delivering the Safe Program (Safe VVI Pacing) .....	79
Viewing or Changing System Settings .....	81
Viewing or Changing Preferred Settings .....	83

<b>Appendix</b> .....	<b>87</b>
Technical Data .....	88
Fixed parameter values .....	92
Programmable parameter values .....	95
Scope of Delivery and Accessories .....	97
Messages .....	99
Electromagnetic Compatibility in compliance with IEC 60601-1-2:2007 .....	101
Country-Related Information .....	105
Symbols on Label .....	106
Index .....	108

## 1

## Safety When Using the Device

**What's in this chapter?**

This chapter contains the following topics:

<b>Topic</b>	<b>Page</b>
Intended Medical Use	6
Required Expertise	8
Contraindications	9
Possible Complications	9
Installation – Setup Requirements	10
Maintenance, Service, Disposal	16
General Safety Instructions	19

## Intended Medical Use

---

**Note:** Federal (USA) law restricts this device to sale by, or on the order of, a physician (or properly licensed practitioner).

### Indication

This device is indicated for use in pacing lead system analysis during the implantation of pacemakers and defibrillators.

**Note:** During implantation of cardiac pacemaker and ICD systems the device does not measure shock impedance and defibrillation threshold (DFT).

### System description

This device is a portable medical diagnostic device intended to be used during the implantation of pacemakers and defibrillators to evaluate the placement and integrity of pacing leads and to determine the appropriate pacing parameters for the implanted device.

The user controls the operations of the device through its LCD touch-screen, control-knob and control-buttons. The user monitors device measurements, heart activity and device status through its LCD display and LED indicators.

The device provides 3-chamber pacing and sensing. It can display up to 3 channels of real-time IEGMs and markers. The device provides standard PSA functionality such as lead impedance measurements, intrinsic amplitude, timing measurements, and pacing thresholds. It has a report feature that provides a summary of device measurements along with optional IEGM traces. It has a freeze/timing feature that allows the user to snapshot IEGM traces and contrast and analyze IEGM morphologies.

The device paces and senses the patient's heart through patient cables that connect the output of the device to the implanted pacing leads.

The device can operate from line-power or AA alkaline batteries. When operating from batteries, the device

continuously displays the battery status. Visual and audio low-battery indicators advise the user when the remaining battery capacity is getting short.

The device has a VGA output port that can be connected to an external CRT/LCD display. It also has a USB connector that can be used to export measurement reports to a USB flash memory stick or print those reports to a Bluetooth printer via an attached USB Bluetooth adapter.

### Intended use

This device is intended to be used during the implantation of pacemakers and defibrillators to evaluate the placement and integrity of pacing leads and to determine the appropriate pacing parameters for the implanted device.

Only trained medical personnel may use the device. While the device is in use, the medical staff must continuously monitor the patient with the aid of a surface ECG monitor and always have cardiac emergency equipment (e.g., external pacemaker, defibrillator) in an operational status available for immediate life support. Consider additional pre-emptive measures in patients where loss of pacing could cause life-threatening danger.

In support of this intended use, the device provides the following measurement/diagnostic capabilities:

- For the sensing of intrinsic events of the heart:
  - P/R wave amplitudes and slew rate
  - Rates (PP, RR interval)
  - Intrinsic AV delay (PR interval)
  - Conduction times
  - Wenckebach point
  - VV delay
- For the pacing of the heart:
  - Pacing threshold in up to 3 chambers
  - Lead impedances
  - Burst pacing

**Note:** The device may not be used as a life-support system. During the implantation, the device is suitable for temporary external pacing in up to 3 chambers while the medical staff must continuously monitor the patient with the aid of a surface ECG monitor.

## Required Expertise

---

**Required knowledge**

Users must be thoroughly familiar with electrotherapy of the heart. Only qualified medical specialists having the special knowledge required for the proper use of the device are permitted to use it.

Users must have normal or corrected to normal eyesight and hearing.

**Medical product operator's ordinance**

---

National regulations may require that the user, manufacturer or manufacturer representative perform and document safety checks of the device during installation. They may also require that the manufacturer or its representative provide training to users on the proper use of the device and its accessories.

- If you do not know the national regulations in your country, please contact your local BIOTRONIK representative.
-

## Contraindications

---

The following applications are contraindicated:

- With AV conduction disorders:
    - Atrial single-chamber pacing
  - With competing intrinsic rhythms:
    - Asynchronous modes
  - With chronic atrial tachycardia as well as chronic atrial fibrillation or flutter:
    - Modes with atrial control (DDD, VDD)
  - With poor tolerance of high ventricular rates (e.g., with angina pectoris):
    - Tracking modes (i.e., atrial control modes) and propensity for atrial tachycardia
  - Use as an external pacemaker outside of the implantation procedure.
- 

## Possible Complications

---

Depending on the patient's condition and depending on the scope and type of pacing program, the following possible complications associated with the use of pacing system analyzers are reported in medical references: life-threatening atrial and ventricular arrhythmia, bradycardia, tachycardia, and asystole.

---

## Installation – Setup Requirements

### Shipping and storage

- If the packaging in which the device is shipped is damaged, do not use the device and contact BIOTRONIK immediately.
- To protect the device from damage, it should be placed in its inner carrying case when stored, carried or shipped.

#### CAUTION

##### Functional impairment due to external damage.

Mechanical impact, for example dropping the unit – unpackaged, from a height of over 5 cm (2 inches) – can permanently impair the function of the system.

- Do not use the device if there is apparent damage.
- If damage has occurred, contact BIOTRONIK to test, and if necessary, repair the device.

**Note:** Remove the batteries from the device if it will not be used for a prolonged period.

### Ambient conditions during shipping and storage

The device may only be stored in an environment that meets the following ambient conditions:

Parameter	Value
Temperature	0 ... 50°C 32 ... 122°F
Relative humidity	30 ... 75 % (non-condensing)
Atmospheric pressure	50 ... 106 kPa 7.2 ... 15.4 lbf/in. <sup>2</sup>

**Ambient conditions during operation**

The device may only be operated under the following ambient conditions:

Parameter	Value
Temperature	10 ... 40°C 50 ... 104°F
Relative humidity	30 ... 75 % (non-condensing)
Atmospheric pressure	70 ... 106 kPa 10.1 ... 15.4 lbf/in. <sup>2</sup>

**Installation site**

The device may only be operated in an environment that meets the following conditions:

- The ambient temperature, relative humidity and atmospheric pressure are within the specified operating conditions.
- Sterile operating conditions can be maintained.
- There are no explosive gases in the vicinity of the device.
- There is always a cardiac emergency equipment (e.g., external pacemaker, defibrillator) in an operational status available for immediate life support. Consider additional pre-emptive measures in patients where loss of pacing could cause life-threatening danger.
- The patient is monitored with a surface ECG monitor.

The device should stand on a level, dry surface. It should be placed so that it cannot slide, especially when cables are connected.

---

### Applying language foil

Fig. 1: The labeling language can be changed by applying a language foil.

- To apply the language foil, remove its adhesive backing, position the foil as shown in the above figure, and press it firmly onto the front panel.

### Power supply

The device can be operated with the following power sources:

- External power supply unit as a primary source – for continuous operation.
- Replaceable batteries as a secondary source – for protection in case of a power failure and for operation when line-power is not readily available. Refer to section "Power Supply" for battery longevity information.

#### **WARNING**

##### **Danger from loss of power**

Operating the device with depleted batteries, unapproved battery types or no batteries can endanger the patient if line-power is temporarily interrupted.

- Do not use the device with depleted batteries, unapproved battery types, or when the battery magazines are not fully populated with batteries.

**Cable and lead connections**

- Inspect cables prior to use. Replace cables if they are worn or damaged.
- Arrange cables to avoid entanglement with equipment or medical staff.
- Cable connections and plugs must be cleaned as defined in section "Maintenance, Service, Disposal". Soiled contacts can lead to signal distortions and false diagnoses.
- All electrical contacts must be dry.
- Make sure that cables are securely connected. When disconnecting cables from the device, clasp the locking connector and pull. Do not pull on the cable.
- All lead connections are swap-safe and encoded at the lead connectors. Do not force cables to connect together or force lead connectors into the connector ports.

**Initial operation**

- Make sure that the device is free of visible damage or contamination and all safety related markings are legible.

**CAUTION****Danger from unauthorized modification of the equipment**

Modification of the device or its accessories may place the health and safety of the operator and/or patient at risk.

- Do not modify the device or its accessories.

**Default settings – start parameters**

When turned on, the device's operational settings are set to the factory settings (see section "Fixed Parameter Values").

**Note:** Press the [**Pref. Set.**] key on the front panel to change the device's operation settings to user-defined settings (see section "Viewing or Changing Preferred Settings").

**Possible electro-magnetic interferences**

When used with its approved accessories and configured and operated in accordance with the instructions in this document, this device meets the electromagnetic compatibility requirements of IEC 60601-1-2.

However, strong electromagnetic interference can occur in close proximity to other electrical equipment. Therefore the device should not be placed adjacent to or stacked on other electrical equipment.

Electromagnetic interference may cause interference with the device, including:

- Resetting of the device.
- Noise on the IEGM traces, spurious sense events. If electromagnetic interference occurs, try the following, as appropriate:
  - Switching off the source of interference.
  - Moving the source of interference away from the device.
  - Disconnecting any electrical connections between the device and the source of interference.
  - Switching the device off and back on, again.
  - If the interference continues, contact your local BIOTRONIK representative immediately.

**Note:**

- Portable and mobile RF communication equipment can affect the operation of the device.
- Use of accessories that are not approved by BIOTRONIK can increase the device's electromagnetic emissions and/or electromagnetic susceptibility.

**Note:** Electromagnetic interference may cause noise on the IEGM traces and spurious sense events and as a consequence may influence the timing of the pacing delivered to the patient.

Lower sensitivity settings increase the susceptibility of the device to interfere with electromagnetic fields of other devices or equipment. A sensitivity setting at or above 1 mV is strongly recommended, if clinically suitable.

**Note:** Prolonged exposure to electromagnetic interference may cause the device to pace asynchronously.

---

## Maintenance, Service, Disposal

When maintaining, cleaning or disposing of the device, the following warnings and cautions must be observed.

### WARNING

#### Fatal injury if exposed to fluids

Before cleaning and disinfecting device surfaces, disconnect the external power supply.

### WARNING

#### Danger of explosion if exposed to cleaning and disinfecting agents

Before operating the device, let cleaning and disinfection agents evaporate.

### CAUTION

#### Damage by cleaning agents

Do not use strong, abrasive cleaning agents or other organic solvents such as ether or benzene. These agents corrode the surface of the device.

### Cleaning and disinfecting

- Use soft, lint-free cloth.
- Housing: Clean with a damp cloth and mild soap solution or a mixture of 70 % from isopropanol and 30 % from water. Disinfect with alcohol or aldehyde-based agents such as Aerodesin 2000, Fugaten spray, Lysoformin 2000, or Aldasan 2000.
- Visually inspect the connectors: Make sure that the ports for all connections and the plugs for the cables are free of debris.
- Patient cables and adapters: Clean with a mixture of 70 % from isopropanol and 30 % from water. Disinfect by exposing to water with a 2 % concentration of Lysoformin for 15 minutes.

**Not to be sterilized**

- The device cannot be sterilized.
  - The following patient cables and adapters are intended for single use only and cannot be sterilized:
    - PK-155
- 

**Sterilization**

- The following patient cables and adapters can be sterilized with steam for 20 minutes at 121°C (250°F) and 1.1 bar with 20 reesterilization cycles permitted:
    - PK-67-L
    - PK-67-S
    - PA-1-B
    - PA-2
    - PA-4
  - The following patient cables and adapters can be sterilized with steam for 5 minutes at 134°C (273°F) and 2.2 bar, or for 20 minutes at 121°C (250°F) and 1.1 bar with 50 reesterilization cycles permitted:
    - PK-141
- 

**Checkup prior to use**

- Perform a checkup on the device and the approved accessories before each use:
    - Inspect the housing for mechanical damage, dents, loose parts, cracks or tampering.
    - Inspect the legibility of the labeling.
    - Inspect cables for damage to the insulation, connectors or electrical contacts.
    - Replace damaged or worn cables.
    - Turn the device on.
      - Verify that all indicator lights flash briefly.
      - Wait 15 seconds to allow an internal self-test to complete.
      - Verify that no error messages are displayed.
  - If the device has mechanical damage or does not pass its internal self-test, do not use the device and contact your local BIOTRONIK representative immediately.
-

**Periodic inspection**

- A complete inspection of the device needs to be performed:
  - in suspicion of malfunction;
  - every two years.

Excessive electrical power, e.g. combined use of HF surgery equipment or a defibrillator on a patient during use of the pacing system analyzer, may damage the device. Check the operation of the device after such usage and initiate this inspection if any signs of malfunction exist, especially if failure messages during startup selftest occur on the display.
- This inspection needs to be performed according to manufacturer’s instructions. These instructions can be provided on demand and state all required test steps and devices related to the check.
- If you have any questions, please contact your local BIOTRONIK representative.

**Disposal**

- 
- This device contains material that must be disposed of in an environmentally safe manner. The European Directive 2002/96/EC regarding electric and electronic equipment (WEEE) applies.
  - The symbol on the label – a crossed out garbage can – indicates that the device must be disposed of in accordance with the WEEE directive. The black bar underneath the garbage can symbol indicates that the device was sold after the national implementation of the WEEE directive was implemented in your country.
  - Return devices that are no longer used to BIOTRONIK.
  - Dispose of the batteries in an environmentally safe manner.
-

## General Safety Instructions

---

<b>Technical manual</b>	The device may be used only in accordance with this technical manual.
<b>Risks of improper handling</b>	<p>Disregarding the safety warnings can endanger the patient, the user, or others, as well as the equipment.</p> <div style="border: 1px solid black; padding: 5px;"><p><b>Note:</b> Failure to observe the safety precautions voids all damage claims and manufacturer liability.</p></div> <p>The following dangers may arise in the event of improper use:</p> <ul style="list-style-type: none"><li>● Failure of important device functions.</li><li>● Personal endangerment due to electrical effects.</li></ul>
<b>No modification of equipment</b>	Do not modify the device or its accessories. Modifications may put the health and safety of the operator and/or patient at risk.
<b>Physician supervision</b>	The device may only be operated under the constant supervision of a physician. During a procedure, the patient must be continuously monitored by medical personnel with the aid of a surface ECG monitor.
<b>Emergency equipment</b>	<p>During the procedure, always have cardiac emergency resuscitation equipment (e.g., external pacemaker, defibrillator) in an operational status available for immediate life support.</p> <p>Consider additional pre-emptive measures in patients where loss of pacing could cause life-threatening danger.</p>
<b>External pacing</b>	The device may not be used as a life-support system. During the duration of the implantation, the device is suitable for temporary external pacing while the patient is being continuously monitored by medical personnel.

---

<b>Touching clips of cables and leads</b>	Do not touch the metal clips on the patient cable or the pacing lead. The device is in electrical contact with the patient's heart and blood via the implanted leads. Touching the metal clips on the patient cable or the pacing lead may expose the patient's heart to dangerous electrical currents.
<b>Connecting cables to leads</b>	Verify that RV patient cable is connected to the RV lead before selecting cross-chamber LV pacing polarities. Connecting the patient cable to the wrong lead may result in ineffective sensing and pacing behavior and loss of pacing support.
<b>Patient cables</b>	Patient cables should not be connected to the device before it has reached the ready-for-operation status (Ready-for-operation status, p. 35). After the device is "ready-for-operation", securely attach all patient cable connections.
<b>Wet cables</b>	Do not use wet cables.
<b>Cables and wounds</b>	Do not allow the cable to come into contact with an open wound.
<b>Disposal of Cables</b>	Dispose of non-contaminated cables according to the European Directive 2002/96/EC regarding electric and electronic equipment (WEEE).
	<b>Note:</b> Cables which have had contact to blood must be disposed of as contaminated medical waste.
<b>Unused cable connections</b>	Attach unused cable connections close to the patient.
<b>Protective sleeves of cables</b>	Before connecting cables ensure correct position of protective sleeves.
<b>Liquids</b>	Avoid spilling liquids on the device or its accessories. While the device is designed for limited protection against the ingress of particulate matter, it is not protected against the ingress of fluids.

<b>Electrostatic potentials</b>	Avoid delivering dangerous electrostatic shocks to the patient or equipment. Before handling the device, the patient cable or the corresponding leads, the electrostatic potential between physician or medical technicians and the patient must be equalized, for instance by touching the patient at a point as far as possible from the leads.
<b>Leakage currents</b>	<p>Avoid leakage currents between all connected devices, the patient cable, and the patient if line-powered devices are used in the vicinity of the patient. Such leakage currents may trigger lethal arrhythmias.</p> <p>Potential equalization cables, if present, must be attached to all connected components.</p> <p>National and international regulations concerning the use of electromedical devices also apply to patient cables.</p>
<b>Self-test</b>	Be aware that when the device is turned on, it performs an internal self-test for approximately 15 seconds. During this time the device's pacing outputs are inactive and the device is unavailable for use.
<b>Pacing threshold test</b>	Consider the health of the patient prior to performing a pacing threshold test. A loss of capture, asystole and pacing during vulnerable periods can occur.
<b>Termination of pacing</b>	Do not abruptly terminate pacing. The sudden termination of pacing can lead to extended periods of asystole in some patients. Gradually decrease the pacing rate until the patient's intrinsic rate is detected.
<b>External ECG device</b>	During the implantation the medical staff must continuously monitor the patient with the aid of a surface ECG monitor.
<b>Pacing mode selection</b>	Select a pacing mode that is consistent with the patient cable connections to the leads. Loss of pacing may be a consequence of inconsistent selection and may cause life-threatening danger to the patient.

<b>Wenckebach test</b>	Because atrial single chamber pacing is contraindicated for use in patients with no AV conduction, the Wenckebach test may not be performed on such patients.
<b>Burst pacing</b>	Burst pacing can induce or accelerate dangerous arrhythmias. Always have cardiac emergency equipment (e.g., external pacemaker, defibrillator) in an operational status available for immediate life support when using this feature. Consider additional pre-emptive measures in patients where loss of pacing could cause life-threatening danger.
<b>High-frequency surgery</b>	<p>The device is equipped with protective circuitry to prevent damage when used with HF surgery. Although this circuitry has been tested to and exceeds standard requirements, its efficiency is limited and depends on the strength, waveform and conduction path of the induced current.</p> <p>In addition use of HF surgery may induce dangerous currents into the patient cables which may be conducted into patient's heart.</p> <p>Therefore:</p> <ul style="list-style-type: none"><li>• Disconnect the patient cables from the device when performing HF surgery procedures.</li><li>• Check the operations of the device if HF surgery has been used.</li></ul>
<b>Defibrillation</b>	While the device is designed and tested to be defibrillator safe, whenever possible, disconnect the device from the patient when defibrillating. If the device is attached during defibrillation, check its operations afterwards.
<b>Basic advice</b>	Only the manufacturer may perform corrective maintenance, enhancements or modifications to the device.

<b>Connecting accessories</b>	Do not connect a USB flash memory stick, a Bluetooth adapter or an external VGA monitor while the device is operating on low battery power (i.e., when the battery gauge is displayed red). Connecting one of these devices may cause a sudden increase in power consumption and cause the device to prematurely shut off.
<b>Sensitivity</b>	<p>Electromagnetic interference may cause noise on the IEGM traces and spurious sense events and as a consequence may influence the timing of the pacing delivered to the patient.</p> <p>Lower sensitivity settings increase the susceptibility of the device to interfere with electromagnetic fields of other devices or equipment. A sensitivity setting at or above 1 mV is strongly recommended, if clinically suitable.</p>
<b>Replacement parts and accessories</b>	Ensure safety compliance by using only original replacement parts and accessories authorized by BIOTRONIK. Using any other parts voids the liability for subsequent events, the guarantee, and the warranty.
<b>Defects</b>	Do not put defective or damaged devices into operation.

## 2

# Handling of the Device

### What's in this chapter?

This chapter contains the following topics:

<b>Topic</b>	<b>Page</b>
Device Overview	26
Transport and Setup	29
Power Supply	30
Switching On and Off	35
Connections and Cables	36
Keys, Displays and Signals	45

## Device Overview

### Operating elements on the front panel

Fig. 2: Operating elements on the front panel

Designation	Description
Display, indicators, keys and control knob	See section "Keys, Displays and Signals"

### Operating elements on right

Fig. 3: Operating elements on the right of the device

### Operating elements on left and back

Designation	Description
	On/off switch
	Right-side battery compartment (secondary power source)

Fig. 4: Operating elements on the left and back of the device

Designation	Description
	Left-side battery compartment (secondary power source)
	USB port used to: <ul style="list-style-type: none"> <li>– export reports to a USB flash memory stick</li> <li>– print reports via Bluetooth adapter to an external printer</li> </ul>
	VGA port used to output the image on the display to an external monitor
 7.5 V	Power supply connector (primary power source)
 A / RV	Patient cable port used to connect the atrium/right-ventricle patient cable
 LV	Patient cable port used to connect the left-ventricle patient cable
Support stand	Support stand used to adjust the viewing angle of the device

Designation	Description
Markings on bottom side:	
	CE mark
	Device must be disposed of in accordance with the WEEE directive

**External Power supply of the device**

Fig. 5: External power supply of the device

Designation	Description
Indicator (green)	Power indicator of the external power supply (primary energy source)
Markings on bottom side:	
	Manufacturer identification
	Manufacturer and distributor
	CE mark
	Follow the instructions for use!
	Safety class II
	ENEC certification by VDE
	Underwriters Laboratories certification for Recognized Component (USA and Canada)

## Transport and Setup

### Transporting the device

- The device's carrying case is designed to protect the device from damage during transport. Transporting the device outside its carrying case may result in damage to the device.

### Setting up the device

#### **WARNING**

##### **Danger from contamination**

The device cannot be sterilized.

- Do not allow the device to enter a sterile area.
- The device should stand on a level, dry surface. It should be placed so that it cannot slide, especially when cables are connected.
- Use the support stand on the back of the device to change the device's viewing angle. Under some lighting conditions, this may improve operation and the visibility of the display.

## Power Supply

### Primary power source

The primary power source is the external power supply, a medical AC–DC power source:

- Type: FRIWO FW 7555M/08,
- input rating: 100 ... 240 V AC, 50/60 Hz, 400 mA,
- output rating: 7.5 V DC, 1.7 A.

### Secondary power source

The secondary power source are 8 standard batteries:

- Duracell MN1500, type AA/LR6 or Energizer E91, type AA/LR6.
- Stored in 2 battery magazines, 4 batteries per magazine.
- Replaceable
  - with the battery magazine, and
  - without special tools.
- With a service time of at least 12 hours – under the following Standard Operating Conditions:
  - 2-chamber pacing,
  - a pacing rate of 70 ppm,
  - a pacing amplitude of 5.0 V,
  - a pacing width of 0.5 ms,
  - at 500 Ohm pacing impedance,
  - with long-life illumination setting,
  - without Bluetooth printer operation,
  - without USB flash memory stick,
  - without active VGA port,
  - at a battery temperature of 20 ... 25°C (68 ... 77°F).

### Supply concept

When the primary power source is unavailable the device automatically switches between the primary and secondary power sources without loss of function.

### Connecting the power supply unit

**Note:** Only operate the device with the designated external power supply. Do not use a different external power supply.

#### WARNING

##### Danger from electrical leakage currents

Electrical currents can harm the patient.

- Do not touch the patient and the device's electrical contacts at the same time.

Fig. 6: Connecting the external power supply to the device

- Plug the DC voltage plug of the external power supply into the DC socket of the device.
- Plug the AC voltage plug of the external power supply into a suitable AC power outlet (for specifications see section "Technical Data").
- Check that the power indicator of the external power supply is on.

Fig. 7: External power supply with power indicator

### Disconnecting the AC power supply

To disconnect the device from AC power, pull the AC voltage plug from the AC power outlet.

**Note:** The AC voltage plug must be easily accessible. Do not block the connection between AC plug and outlet.

### Determining the battery status

When operating on battery power, the device provides an always visible battery level indicator which indicates the remaining battery service time under the Standard Operating Conditions.

Fig. 8: Battery level indicator on the display

- If the batteries have more than 4 hours of remaining service time, the battery level indicator is green.
- The battery level indicator turns yellow when the remaining service time is between 2 and 4 hours.
- The battery level indicator turns red when the remaining service time is only 30 minutes. When the battery level indicator is red, the device emits an audible beep every 20 seconds.
- The battery level indicator displays a question mark  when the remaining service time is being calculated.
- A warning displays when the device is switched on if:
  - The remaining service time is less than 30 minutes.
  - The batteries are completely depleted or there are no batteries in the device.

**Note:** Replace the depleted batteries with new batteries and acknowledge this warning before continuing to use the device.

#### **WARNING**

##### **Danger from loss of power**

When the battery level indicator is red, the device has less than 30 minutes of remaining battery service time.

- For additional patient safety connect the external power supply at this point.
- When operating on battery power, do not attempt to replace the batteries in either battery magazine when the battery level indicator is red.
- Do not plug in a USB device or attach an external VGA monitor when the battery level indicator is red.

**Note:** Prolonged attachment of a USB device or an external VGA monitor will decrease the service time of the batteries.

When operating on battery power and the power in one of the battery-magazines nears exhaustion, the device will visually designate that magazine as needing replacement.

Fig. 9: Display indicates: Right battery magazine is depleted and should be replaced.

- When operating on battery power and the battery level indicator is NOT red, the batteries in the designated battery magazine may be replaced with fresh batteries.

Fig. 10: Display indicates: Both battery magazines are depleted, do not remove any magazine.

- When operating on battery power and the battery level indicator is red, do not replace the batteries, unless you first connect the power supply.

## Exchanging the batteries

The battery compartments are located on the right and left sides of the device.

Each battery compartment contains a battery magazine with four batteries (for specifications see section "Technical Data").

- When operating on battery power and replacing batteries, replace the batteries in the depleted magazine first.
- Only use the approved batteries (see section "Secondary power source"). Do not use rechargeable batteries.
- Only use new batteries. Always replace all four batteries in a battery magazine at the same time.
- Remove the batteries if the device will not be used for a prolonged period of time.
- Use and store the batteries at room temperature.

Fig. 11: Battery compartments with battery magazines

- To open a battery compartment, lift the locking mechanism and pull the battery magazine out.
  - Remove the old batteries from the battery magazine and dispose of them properly.
  - Insert the new batteries with the terminals (+/-) aligned as shown on the battery magazine.
  - Push the battery magazine into the battery compartment until the locking mechanism latches.
  - Turn on the device and check that the battery level indicator displays sufficient battery capacity.
-

## Switching On and Off

### Prerequisites

Check the following requirements before switching on:

- No patient cables should be connected to the device before it has reached the ready-for-operation status (Ready-for-operation status, p. 35).
- All batteries must be inserted and ready for use – even on line-power – to avoid loss of function in the event of line-power interruptions.
- Do not operate the device with the external power supply if the device is missing batteries or the batteries do not have at least 30 minutes of remaining service time (see section "Determining the battery status").

### WARNING

#### Danger from loss of pacing support

After the device is switched on, the pacing functions are switched off for approximately 15 seconds while a self-test is conducted, so that no pacing is possible in this time.

- Keep a separate pacing device ready for emergencies.

### Switching on the device

The on/off button  is located on the right side of the device.

- Press the on/off button in to switch on the device.

### Start screen and self-test

The start screen appears shortly after the device is switched on. A self-test then runs for approximately 15 seconds.

### Ready-for-operation status

After the self-test, the device is ready for use.

- Continuously monitor all displays and signals for plausibility. If a display is not plausible, look for the cause. If necessary, switch the device off and then on again.

### Switching off the device

- Press the on/off button again to switch off the device.

## Connections and Cables

**Note:** Do not force the plugs into the ports. When disconnecting the plugs, do not pull on the cable.

### Connecting an external monitor

The image on the display can also be sent to an external monitor.

**Note:** When connecting an external monitor, the user is creating a medical electrical system and is responsible that this system complies with the requirements of IEC 60601-1.

- When using an external monitor, measure the leakage currents.

#### WARNING

##### Danger from electrical leakage currents

Electrical currents can be dangerous to the patient.

- Never touch the patient and the device's electrical contacts at the same time.

#### WARNING

##### Danger from electrostatic charges

The lead system is in electrical contact with the patient's heart and blood. Touching the metal clips on the patient cable or the pacing lead may induce dangerous electrical currents in the patient's heart.

- Do not touch the metal clips on the patient cable or the pacing lead.

**Note:** It is possible to connect or disconnect the external monitor while the device is on.

Fig. 12: Connecting an external monitor

- Connect the VGA plug of the monitor to the VGA port .

When the device detects a connection to a VGA monitor, it will enable its external VGA port. The device may not detect some VGA monitors, so it may be necessary to manually enable the external VGA port.

To manually enable or disable the VGA port, proceed as follows:

- Press the **[Pref. set.]** button to go to the Preferred Settings screen.
- Press the **[VGA]** button repeatedly to activate one of the following settings:
  - **[OFF]**
  - **[ON]**
  - **[auto detect.]** (factory setting)

#### Preparing for printing or exporting reports

The USB port of the device can be used as interface for printing or exporting reports to USB/Bluetooth devices (see section "Printing or Exporting a Report").

- Using a Bluetooth adapter, reports can be sent to an external printer.
- Using a USB flash memory stick, reports can be saved as PDF files.

#### Connecting a Bluetooth adapter

Only connect Bluetooth adapters

- Complying with the Microsoft Bluetooth Stack standard
- With no electrical connection with the line power or other conducting parts

#### WARNING

##### Danger from electrical leakage currents

Electrical currents can be dangerous to the patient.

- Never touch the patient and the device's electrical contacts at the same time.

**Note:** It is possible to connect or disconnect a Bluetooth adapter while the device is on.

Fig. 13: Connecting a USB Bluetooth adapter for printing reports

- Connect a Bluetooth adapter to the USB connector port .

**Connecting a USB flash memory stick**

Only connect USB flash memory sticks with the following features :

- No electrical connection with the line power or other conducting parts
- No customized functions like activated Boot manager, Password protection or Hardware encryption

**WARNING**

**Danger from electrical leakage currents**

Electrical currents can be dangerous to the patient.

- Never touch the patient and the device's electrical contacts at the same time.

**Note:** It is possible to connect or disconnect a USB flash memory stick while the device is on.

Fig. 14: Connecting a USB flash memory stick for exporting reports

- Connect a USB flash memory stick to the USB connector port .

## Patient cable overview

The following patient cables and adapters can be used for the device.

**Note:** Only the patient cables and adapters listed here may be used.

**Note:** EMC conformity of these cables is warranted for use with this device.

### Patient cables and adapters

Patient cable	Adapter	Patient connection
PK-141, sterile, can be resterilized	–	Insulated alligator clips
PK-67-L, sterile, can be resterilized	PK-155, sterile, single use	Insulated alligator clips
PK-67-S, sterile, can be resterilized	PK-155, sterile, single use	Insulated alligator clips
PK-67-L, sterile, can be resterilized	PA-1-B, sterile, can be resterilized	2 mm pin connector
PK-67-L, sterile, can be resterilized	PA-2, sterile, can be resterilized	IS-1 connector
PK-67-L, sterile, can be resterilized	PA-4, sterile, can be resterilized	Insulated alligator clips
PK-67-S, sterile, can be resterilized	PA-4, sterile, can be resterilized	Insulated alligator clips


Fig. 15: Patient cable PK-141 (2.8 m)

Fig. 16: Patient cable PK-67-L (2.6 m)

Fig. 17: Patient cable PK-67-S (0.8 m)

Fig. 18: Patient adapter PK-155 (2.0 m)

Fig. 19: Patient adapter PA-1-B

Fig. 20: Patient adapter PA-2


Fig. 21: Patient adapter PA-4

### Connecting the patient cables

**Note:** Sterility can not be guaranteed if the sterile packaging was damaged or improperly handled or stored.

If the sterile package is wet, punctured, opened or otherwise damaged, do not use the cable.

### WARNING

#### Danger from loss of function

Moisture from wet cable can impair cable function and endanger the patient.

- Do not use wet cables.

**WARNING****Danger from loss of function**

Incorrect positioning of the protective sleeves over the cable clip(s) can cause unintended electrical connections that can impair cable function and endanger the patient.

- Before connecting cables ensure correct position of protective sleeves.

**WARNING****Danger from allergic reaction**

The contact of the cable and open wounds may result in allergic reaction of the patient.

- Do not allow the cable to come into contact with an open wound.

**WARNING****Danger from electrical currents**

Unused cable connections can induce electrical currents into the patient's heart.

- Attach unused cable connections close to the patient.

**Note:** Before connecting the patient cable to the leads, ensure that the leads are securely implanted in the patient's heart.

- Patient cables may only be used by healthcare professionals that are qualified for intracardial examinations and therapy.
- Precautionary measures must be maintained while conducting an intracardial examination. Use in suitable rooms featuring X-ray facility and cardiac emergency equipment (e.g., external pacemaker, defibrillator) in an operational status available for immediate life support. Consider additional pre-emptive measures in patients where loss of pacing could cause life-threatening danger.

**WARNING****Danger from defibrillation**

While the device is designed and tested to be defibrillator safe, the patient can be endangered and the device can be damaged.

- Whenever possible disconnect the device from the patient when defibrillating.
- If the device is connected to the patient during defibrillation, check its operations afterwards.

**WARNING****Danger from HF surgery**

The device is equipped with protective circuitry to prevent damage when used with HF surgery. Although this circuitry has been tested to and exceeds standard requirements, its efficiency is limited and depends on the strength, waveform and conduction path of the induced current.

In addition use of HF surgery may induce dangerous currents into the patient cables which may be conducted into patient's heart.

Therefore:

- Disconnect the patient cables from the device when performing HF surgery procedures.
- Check the operations of the device if HF surgery has been used.

**WARNING****Danger from electrostatic charges**

Electrical currents can be dangerous to the patient.

- Discharge any electrical static charge on your person before touching the patient, the patient cables or the device.

- Please observe the technical manual of the patient cables and adapters that are used.
- If patient cables and adapters are used, first connect the patient cables with the adapters.
- Only hold alligator clips by the protective sleeve.

- No patient cables should be connected to the device before it has reached the ready-for-operation status (Ready-for-operation status, p. 35). Current parameter settings should be verified, prior to connecting the patient cables to the device.

**Note:** It is possible to connect or disconnect the patient cables while the device is on.

Fig. 22: Connecting a patient cable (here: A/RV) to the device

For the right side of the heart (A/RV):

- Insert the patient cable's locking connector into port  A/RV. Ensure that the plug has locked into the port.
- Connect the patient connections of the patient cable or adapter to the A/RV leads of the patient.

For the left side of the heart (LV):

- Insert the patient cable's locking connector into port  LV. Ensure that the plug has locked into the port.
- Connect the patient connections of the patient cable or adapter to the LV leads of the patient.

### Removing the patient cables

- 
- Disconnect the patient cable or adapter from the pacing lead by opening the clips and removing the pacing lead.
  - Disconnect the patient cable from the device by pulling on the cable's locking connector.
-

**Disposal of cables**

**Note:** Dispose of used cables in the following manner:

- Sterile cables: as contaminated medical waste.
  - Unsterile cables: according to the European Directive 2002/96/EC regarding electric and electronic equipment (WEEE).
-

## Keys, Displays and Signals

### Keys on the front panel

The keys on the front panel have specific functions and are labeled accordingly.

Fig. 23: Keys on the front panel

Designation	Description
	Key used to initiate the Safe Program
Save:	Keys used to save the measurements of the selected channel in the report:
- <b>[Sensing]</b>	- P/R measurements
- <b>[Threshold]</b>	- Pacing amplitude
- <b>[Impedance]</b>	- Lead impedance
- <b>[Timing/ Freeze]</b>	Key used to open the IEGM analysis window, the Timing/Freeze screen
<b>[Max Stim]</b>	Key used to initiate pacing at maximum amplitude and pulse width in the selected channel
<b>[Pref. Set.]</b>	Key used to collectively change the value of parameters to a set of user defined values

- To activate a function, briefly press the corresponding key.
- Exception: To pace with maximum amplitude and pulse width, press and hold the **[Max Stim]** key.

**Indicators on the front panel**

Indicator lights for atrial, right-ventricular, and left-ventricular senses and paces are located above the display.

Fig. 24: Indicators on the front panel

Designation	Description
<b>[As]</b> (green)	Atrial sense indicator light
<b>[Ap]</b> (yellow)	Atrial pace indicator light
<b>[RVs]</b> (green)	Right-ventricular sense indicator light
<b>[RVp]</b> (yellow)	Right-ventricular pace indicator light
<b>[LVs]</b> (green)	Left-ventricular sense indicator light
<b>[LVp]</b> (yellow)	Left-ventricular pace indicator light

**Display**

The display is a touch-sensitive screen (touch screen), which is operated with the touch of a finger.

**Note:** Operating the touch screen with sharp objects (e.g., pens, surgical instruments) may damage the display and will reduce the lifetime of the display.

The following is displayed on the screen:

- Parameters and measurements
- IEGM and marker channel
- Buttons

**Buttons on the display**

The display buttons are labeled and have specific functions. The buttons can be used to:

- navigate between windows;
- turn on/off a function;
- select from one of several possible values by pressing repeatedly;
- select a value, which can then be modified by turning the control knob on the front panel.

**Control knob on the front panel**

Use the control knob to modify the selected parameter highlighted on the display. The new value is effective immediately.

- Before turning the control knob, ensure that the desired parameter is highlighted in the display.

**Acoustic signal**

- The device emits a high pitched tone when pressing a key or button has a functional effect.
- The device emits a low pitched tone when pressing a key or button has no functional effect.

**Setting the illumination brightness**

The illumination brightness of the display is preset to a basic value (setting `Normal`) to provide an optimal compromise between light intensity of the display and service time for battery operation.

**Note:** Changing the display illumination setting to **[Bright]** will decrease the battery service life.

To change the illumination brightness, proceed as follows:

- Press the **[Pref. set.]** button to go to the `Preferred Settings` screen.
- Press the **[Illumination]** button repeatedly to activate one of the following settings:
  - **[Normal]** (factory setting)
  - **[Bright]**
  - **[Long life]**

## 3 Using the Software

### What's in this chapter?

This chapter contains the following topics:

Topic	Page
Overview of the Software	50
Handling the Operating Elements	52
Measuring P/R Wave Amplitudes	60
Measuring Pacing Threshold and Lead Impedance	62
Measuring Conduction Times	64
Measuring the Wenckebach Point	67
Printing/Exporting a Report	70
Delivering Burst Pacing	74
Delivering Maximum Pacing	77
Delivering the Safe Program (Safe VVI Pacing)	79
Viewing or Changing System Settings	81
Viewing or Changing Preferred Settings	83

## Overview of the Software

### Description

The device provides the following measurement/diagnostic capabilities:

- For the sensing of intrinsic events of the heart:
  - P/R wave amplitudes and slew rate,
  - Rates (PP, RR interval),
  - Intrinsic AV delay (PR interval),
  - Conduction times,
  - Wenckebach point,
  - VV delay.
- For the pacing of the heart:
  - Pacing thresholds in up to 3 chambers,
  - Lead impedances,
  - Burst pacing.

**Note:** Before each implantation the battery status should be checked.

**Note:** After switching on the device, a self-test is initiated, which detects and displays problems with the device. Depending on the display, please inform your representative immediately.

The device always starts with the same parameters. These start parameters are predefined at the factory and cannot be modified.

- You can save an individual parameter set in the Preferred settings (Pref. set.) screen. You can activate these settings by using the **[Pref. Set.]** key on the front panel.

The software application has the following main screens: Tests, Report, Burst, Preferred settings (Pref. set.), and Timing/Freeze.



- From the Tests screen, you can test for extra-cardiac stimulation and measure P/R wave amplitudes, lead impedances, and pacing thresholds.
  - From the Timing/Freeze screen, you can perform detailed analyzes of IEGM traces.
- From the Report screen, you can view, print and export reports containing measurement values.

- From the `Burst` screen, you can deliver high-rate pacing in a selected channel.
  - From the `Pref. set.` screen, you can change system settings and configure your preferred device settings.
-

## Handling the Operating Elements

**General properties** The display has a touch-sensitive surface (touch screen). To avoid damaging the display, do not touch the surface with sharp objects.

In the various screens of the software application, the following general operating principles apply.

Example	Description
	Control button: Function can be selected.
	Control button: Function is selected.
	Channel Control button: Channel is OFF.
	Channel Control button: Channel is ON.
	Channel Control button: Channel is ON and selected.
	ON/OFF Control button: Function is OFF.
	ON/OFF Control button: Function is ON.
	Parameter Control button: Parameter is not active.
	Parameter Control button: Parameter is active and can be selected.

Example	Description
	Parameter Control button: Parameter is selected, active, and can be modified with the control knob on the front panel; the new value becomes immediately effective.
	Parameter Control button: Parameter is selected, but not active, and can be modified with the control knob on the front panel.

### Operating elements in all screens

The following figures and explanations describe the operating elements (buttons, icons and data fields) of the software application.


Fig. 25: Operating elements in all screens

Designation	Description
Display: active mode	Displays the active mode.
Display: current rates [A/RV/LV]	Displays the sensed rates in bpm for all connected leads. Refractory senses are used in the calculation of the rate. The displays are updated for each event. Rates lower than 30 bpm are not displayed.

Designation	Description
Display: power source	Displays the type of power supply and the state of the batteries:  : Line power operation;  : Battery status;  : Battery status check (shortly after turning on the device or after changing the battery magazine).
Display: left battery magazine depleted	Indicates that the batteries in the left battery magazine are depleted: 
Display: right battery magazine depleted	Indicates that the batteries in the right battery magazine are depleted: 
Display: battery magazines must not be replaced	Indicates that the batteries in both battery magazines must not be replaced, unless the power supply is connected:  
<b>[Tests]</b> button	Switches to the Tests screen.
<b>[Report]</b> button	Switches to the Report screen.
<b>[Burst]</b> button	Switches to the Burst screen (high-rate pacing).
<b>[Pref. set.]</b> button	Switches to the Preferred settings screen (system settings and preferred settings).

### Operating elements on the left in the Tests screen


Fig. 26: Operating elements on the left in the Tests screen

Designation	Description
<b>[Modes]</b> button	Opens a list of the commonly used single-chamber and dual-chamber modes.
Channel ON OFF buttons <b>[A]/[RV]/[LV]</b>	<p>Activate or deactivate the corresponding channel:</p> <ul style="list-style-type: none"> <li>● : activated;</li> <li>○ : deactivated.</li> </ul> <p><b>Note:</b> In order to activate or deactivate a channel, the channel has to be selected (Channel Control button highlighted) first.</p> <p>Sensing and pacing is disabled in deactivated channels, and the IEGM traces and markers are not displayed. The mode is changed accordingly.</p>
<b>[Pace]</b> buttons	<p>Enable or disable pacing for the corresponding channel:</p> <ul style="list-style-type: none"> <li>● : enabled;</li> <li>○ : disabled.</li> </ul> <p>The mode is changed accordingly.</p>
<b>[Sense]</b> buttons	<p>Enable or disable sensing for the corresponding channel:</p> <ul style="list-style-type: none"> <li>● : enabled;</li> <li>○ : disabled.</li> </ul> <p>The mode is changed accordingly.</p>
<b>[Tracking]</b> button	Enables or disables tracking. Switches the mode from DDD to DDI, e.g.

### WARNING

#### Danger from abruptly terminating pacing

Abruptly terminating pacing may result in extended periods of asystole in some patients.

- Gradually decrease the pacing rate until the patient's intrinsic rate is detected for a controlled transition from pacing to intrinsic action.

### Operating elements on the right in the Tests screen


Fig. 27: Operating elements on the right in the Tests screen

Designation	Description
<b>[Rate]</b> button	Displays the pacing rate. To change, press the button and set the value with the control knob on the front panel. The new value is immediately effective.
Channel parameter buttons <b>[A]/[RV]/[LV]</b>	Display the pacing parameters for the respective channel. To pace in a channel, the channel must be active and pacing enabled. To select a channel, press the button. The button can display the following states: <ul style="list-style-type: none"> <li>– Dark button with white text, other text black: Parameter is active and can be modified with the control knob on the front panel; the new value is immediately effective.</li> <li>– Dark button with white text, other text light gray: Parameter is not active, can however be modified with the control knob on the front panel; the new value does not take effect until the parameter becomes active.</li> </ul>
<b>[Amplitude]</b> buttons	Display the pace amplitude in the respective channel. To change, press the button and set the value with the control knob on the front panel. The new value is immediately effective.
<b>[Pulse width]</b> buttons	Display the pulse width in the respective channel. To change, press the button and set the value with the control knob on the front panel. The new value is immediately effective.

Designation	Description
<b>[More]</b> button	Opens a window with additional parameter buttons: <ul style="list-style-type: none"> <li>- <b>[AV delay]</b>;</li> <li>- <b>[RV/LV delay]</b>;</li> <li>- <b>[Sensitivity]</b> A/RV/LV;</li> <li>- <b>[Polarity LV]</b> (BIPL, CRBP, IVBP, RRPB);</li> <li>- <b>[Back]</b>: back to the channel parameter buttons.</li> </ul>

### WARNING

#### Danger from loss of pacing support

Connecting the patient cable to the wrong lead may result in ineffective sensing and pacing behaviour and loss of pacing support.

- Verify that the RV patient cable is connected to the RV lead before selecting cross-chamber LV pacing.

#### Operating elements in the IEGM window


Fig. 28: Operating elements in the IEGM window

Designation	Description
Marker channels (A/RV/LV)	Display a real-time sequence of event markers. Markers are displayed for all activated channels (see channel ON/OFF buttons): <ul style="list-style-type: none"> <li>- Display P: pace marker;</li> <li>- Display S: sense marker;</li> <li>- Display RS: refractory sense marker.</li> </ul>

Designation	Description
IEGM traces (A/RV/LV)	Display a real-time plot of the IEGM signals. Signals are displayed for all activated channels (see channel ON/OFF buttons). The IEGM traces in the IEGM window are wide range filtered.
Display: Measured values	Displays P/R wave measurements and pace amplitudes. The P/R wave amplitudes (baseline-to-peak) are displayed in mV at the respective events. The pacing amplitudes are displayed in V at the respective events.
 Buttons for amplitude gain	Control the gain of the respective IEGM trace. Gains: 0.5, 1.0, 2.0, 5.0, 10.0 mm/mV. Cycle through the available gains by pressing the button repeatedly.
Display: Lead impedance	Displays the lead impedance measurement from the last pace, in the respective channel. For each pacing pulse the lead impedance is determined and displayed in Ohm for the respective paced channel: <ul style="list-style-type: none"> <li>- Display 123<math>\Omega</math>: determined impedance (example);</li> <li>- Display &lt;100<math>\Omega</math>: lead not correctly connected, low impedance lead;</li> <li>- Display &gt;3k <math>\Omega</math>: no lead or lead not correctly connected, high-impedance lead;</li> <li>- Display -- <math>\Omega</math>: no pacing.</li> </ul> When pacing is turned off in a channel or a non-pacing mode is activated, the lead impedance measurement will be cleared within 2 seconds.

### Operating elements in the Timing/Freeze screen


Fig. 29: Operating elements in the Timing/Freeze screen

Designation	Description
 Buttons for amplitude gain	Control the gain of the respective IEGM trace. Gains: 0.5, 1.0, 2.0, 5.0, 10.0 mm/mV. Cycle through the available gains by pressing the button repeatedly.
<b>[Speed] [+]/[-]</b> buttons	Control the time scale. Time scale: 25, 50, 100 mm/s.
 Slider	Moves the displayed timing interval of the IEGM window to the right or left.
Calipers	Measure the time interval between the calipers in ms and ppm. A caliper can be repositioned by dragging it, or by first selecting it and then using the control knob on the front panel to move it left or right.
	Move the respective caliper (e.g., left caliper or right caliper) to the left or right.
	Move both calipers to the left or right by one time interval.
<b>[Grid]</b> button	Plots a grid (1000 ms by 5 mm) for use in measuring amplitudes and intervals.
<b>[Save]</b> button	Saves the currently displayed, 5-second snapshot of the IEGM window to the report. In the Report screen press the <b>[Details]</b> button to view the measured values and the snapshot of the IEGM window.
<b>[Close]</b> button	Closes the Timing/Freeze screen and returns to the Tests screen.

## Measuring P/R Wave Amplitudes

### Description

P/R wave amplitude measurements (e.g., intrinsic amplitudes) are used in assessing the ability of the implanted pacemaker or ICD to sense the intrinsic activity of the heart at the current lead position within the heart.

If the intrinsic rate of the patient is above the set pacing rate or pacing is disabled, the heart's intrinsic activity can be measured.

Measurements of P/R wave amplitudes are made from the filtered IEGM signal. The P wave or R wave amplitude is, together with the rate, displayed in the IEGM window.

An intrinsic amplitude measurement is made on each sense event. Intrinsic amplitude measurements, for a channel, are displayed continuously in the IEGM window, next to the IEGM trace for that channel.

### Performing the measurement

- Switch to the **Tests** screen by pressing the **[Tests]** button.
- Deactivate all channels to which no leads are connected by pressing the channel **ON|OFF** buttons **[A]/[RV]/[LV]**.

To measure intrinsic amplitudes in a channel, proceed as follows:


Fig. 30: P/R wave amplitudes: Tests screen

- Use the **[Pace]** button to disable pacing in the channel, or select a non-pacing mode. Alternatively, use the **[Rate]** button to set the pacing rate below the patient's intrinsic rate. Measured intrinsic amplitudes are displayed next to each sense event of the channel's IEGM trace.
  - Save the measured values and a 5-second snapshot of the IEGM window in the report by making sure that the channel is selected and then pressing the **[Sensing]** key on the front panel. To perform measurements for additional channels, repeat the previous steps for each channel.
- 

#### Terminating the measurement

- To view, print or export the report, go to the **Report** screen using the **[Report]** button.
-

## Measuring Pacing Threshold and Lead Impedance

### Description

Pacing threshold measurements are used in assessing the ability of the implanted pacemaker or ICD to effectively and efficiently pace the heart at the current lead position within the heart.

### WARNING

#### Danger from loss of capture

Pacing threshold testing implies loss of capture. At loss of capture, asystole and pacing during vulnerable periods can occur.

- Consider the health of the patient prior to performing a pacing threshold test.

Lead impedance measurements are used in assessing the integrity of the lead and the connection of the lead's electrodes with the myocardium.

A lead impedance measurement is made on each pace. The last lead impedance measurement, in a channel, is displayed on the left-side of the IEGM window next to the IEGM trace for that channel.

### Performing the measurement

- Switch to the **Tests** screen by pressing the **[Tests]** button.

To measure the pacing threshold and lead impedance in a channel, proceed as follows:


Fig. 31: Pacing threshold: Tests screen

- Activate the desired channel by pressing the channel ON/OFF button **[A]**, **[RV]** or **[LV]**. Enable pacing by pressing the **[Pace]** button.
- Press the **[Rate]** button and increase the pacing rate with the control knob on the front panel until it is higher than the intrinsic rate of the patient. In the IEGM window, pace amplitudes are displayed next to each pace event of the channel's IEGM trace.
- Press the **[Amplitude]** button and reduce the pacing amplitude with the control knob on the front panel until capture is lost. Increase the pacing amplitude again until the pacing just starts to be effective.
- Save the current pacing parameters (pulse amplitude and pulse width) and a 5-second snapshot of the IEGM window in the report by making sure that the channel is selected and then pressing the **[Threshold]** key on the front panel.
- Reset the pacing amplitude to a safe value.
- Save the current pace impedance measurement and a 5-second snapshot of the IEGM window in the report by making sure that the channel is selected and then pressing the **[Impedance]** key on the front panel.

To perform measurements for additional channels, repeat the previous steps for each channel.

#### Terminating the measurement

- Press the **[Rate]** button and change the pacing rate to the original value with the control knob on the front panel. Or press the **[Pref. Set.]** key on the front panel to return to your preferred settings.
- To view, print or export the report, go to the Report screen using the **[Report]** button.

## Measuring Conduction Times

### Description

Conduction measurements are used in determining appropriate pacing parameters for the implanted pacemaker or ICD.

The retrograde conduction is measured by pacing in the ventricles and measuring the time interval between the ventricular pace and any subsequently induced P wave.

### Performing the measurement

- Switch to the **Tests** screen by pressing the **[Tests]** button.

To measure the retrograde conduction in a channel, proceed as follows:

- Use the **[Mode]** button to set the mode to VDI. Use the **[Rate]** button to set the pacing rate above the patient's intrinsic rate.


Fig. 32: Conduction times: Tests screen

- To document the result of the conduction test and to freeze the IEGM, press the **[Timing/Freeze]** key on the front panel. The **Timing/Freeze** screen is displayed.


Fig. 33: Conduction times: Timing/Freeze screen

The screen shows the frozen IEGM with markers, from which – depending on the current mode – VA, AV and VV conduction times can be evaluated.

### Operating the Timing/Freeze screen

The controls on the Timing/Freeze screen can be used to isolate and analyze portions of the frozen IEGM:

- To change the time scale of the IEGM, press the **[Speed] [+]/[-]** buttons.
- To move the viewing window of the IEGM left or right, move the slider.
- To change the amplitude gain of the IEGM trace for one channel, press the corresponding **[Gain]** button repeatedly.
- To measure the time between a ventricular pace and the ensuing atrial sense, use the calipers. Reposition a caliper by dragging it, or by first selecting it and then using the control knob on the front panel. The time between calipers is displayed next to the calipers.
- To superimpose a measurement grid on the IEGM traces, press the **[Grid]** button.
- To save the evaluated IEGM window to the report, press the **[Save]** button.

**Terminating the measurement**

- To end the analysis and return to the Tests screen, press the **[Close]** button.
  - Change the mode and the pacing rate back to their original values. Or press the **[Pref. Set.]** key on the front panel to return to your preferred settings.
  - To view, print or export the report, go to the Report screen using the **[Report]** button. In the Report screen press the **[Details]** button to view the measured values and the snapshot of the IEGM window.
-

## Measuring the Wenckebach Point

### Description

Wenckebach point (WKB) measurements are used in determining appropriate pacing parameters for the implanted pacemaker or ICD.

The Wenckebach point is measured in AAI mode with simultaneous observation of the ventricular IEGM. The pacing rate is increased until Wenckebach behavior is displayed in the IEGM.

### WARNING

#### Danger from loss of pacing support

Single chamber atrial modes are contraindicated for patients with impaired AV conduction.

- If the patient has impaired AV conduction, a Wenckebach test must not be performed.

### Performing the measurement

- Switch to the **Tests** screen by pressing the **[Tests]** button.

To determine the Wenckebach point, proceed as follows:

- Use the **[Mode]** button to set the mode to AAI. Use the **[Rate]** button to set the pacing rate above the patient's intrinsic rate.


Fig. 34: Wenckebach procedure: Tests screen

- Increase the pacing rate while observing the ventricular IEGM. Look for the point at which Wenckebach behavior occurs.


Fig. 35: Wenckebach point: Tests screen

- To document the result of the Wenckebach test and to freeze the IEGM, press the **[Timing/Freeze]** key on the front panel. The *Timing/Freeze* screen is displayed.


Fig. 36: Wenckebach point: Timing/Freeze screen

The screen shows the frozen IEGM, where the PP interval, displayed next to the atrial markers in the marker channel, can be used to calculate the Wenckebach rate.

### Operating the Timing/Freeze screen

The controls on the *Timing/Freeze* screen can be used to isolate and analyze portions of the frozen IEGM.

- See section "Measuring Conduction Times" for instructions on operating the *Timing/Freeze* screen.

**Terminating the measurement**

- To end the analysis and return to the `Tests` screen, press the **[Close]** button.
  - Change the mode and the pacing rate back to their original values. Or press the **[Pref. Set.]** key on the front panel to return to your preferred settings.
  - To view, print or export the report, go to the `Report` screen using the **[Report]** button. In the `Report` screen press the **[Details]** button to view the measured values and the snapshot of the IEGM window.
-

## Printing/Exporting a Report

---

### Description

The following information can be saved in a report:

- For every channel (A/RV/LV)
  - Measured P/R wave amplitudes,
  - Pacing parameters at the pacing threshold,
  - Lead impedance measurements.
- The IEGM window of the Timing/Freeze screen.

Up to 15 reports can be stored in the report list. Each report is displayed in the report list with date and time or a name.

**Note:**

- Once 15 reports have been stored and a new one is added, the oldest report is automatically deleted.
- The report list remains stored when the device is switched off.

You can print or export reports.

---

### Preparing for printing or exporting reports

To prepare the necessary resources for printing or exporting reports, proceed as follows:

- If you want to print reports to an external Bluetooth printer, be sure a Bluetooth adapter is connected to the USB connector port (see section "Connections and Cables"). Be sure the designated printer is on.
  - If you want to export reports as a PDF file, be sure a USB flash memory stick is connected to the USB connector port (see section "Connections and Cables").
-

## Summary report

To print, export, view or manage reports, proceed as follows:

- Go to the Report screen by pressing the **[Report]** button. The current summary report is displayed with the saved measurements.

ADI	A 89 bpm	RV 89 bpm	LV 89 bpm
Report	A RV LV		
John Johnson			
Software version: M2_0_1			
Serial number: 123456789			
<b>Sensing</b>			
Min. amplitude	8.7	15.8	15.8 mV
Mean amplitude	8.8	16.0	16.0 mV
Max. amplitude	8.9	16.1	16.1 mV
Slew rate	0.9	1.0	1.0 V/s
Measuring time	12:07:11 PM	12:08:18 PM	12:08:15 PM
<b>Threshold</b>			
Pulse amplitude	5.0	5.0	5.0 V
Pulse width	0.5	0.5	0.5 ms
Impedance	421	503	384 Ω
Pacing polarity	---	---	---
Pacing current	11	9	13 mA
Measuring time	12:07:40 PM	12:07:29 PM	12:07:33 PM
<b>Impedance</b>			
Impedance	422	503	385 Ω
Pulse amplitude	5.0	5.0	5.0 V
Pulse width	0.5	0.5	0.5 ms
Measuring time	12:07:41 PM	12:07:30 PM	12:07:34 PM

Fig. 37: Summary report in the Report screen

- To view the list of saved reports, press the **[Report list]** button.
- To print or export the currently displayed report, press the **[Print/Export]** button. If the currently displayed report is empty, the button is inactive.
  - Follow the instructions in the displayed window to use the designated USB device. If no USB device information is displayed, disconnect and connect the USB device again.
  - For printing, follow the instructions in the displayed window to use the designated printer or search for new printers.
  - If an error message occurs, follow the troubleshooting instructions (see section "System Messages").
- To view a detailed report with both measured values and 5 second snapshots of the IEGM window, press the **[Details]** button. A message window indicates briefly that the detailed report is being prepared.

## Detailed report


Fig. 38: Detailed report in the Report screen

- To view the summary report, press the **[Details]** button.
- To print or export the currently displayed report, press the **[Print/Export]** button. If the currently displayed report is empty, the button is inactive.
- To view the list of saved reports, press the **[Report list]** button.

## Report list

Fig. 39: Report list in the Report screen

- To view a summary report, select the desired entry from the list and press the **[View Report]** button.
- To rename a report, select the desired entry from the list and press the **[Rename Report]** button. Change the name of the report using the on-screen keyboard (see figure below). Close the on-screen keyboard by pressing the return button.

- To add a new report – for saving new measurements and IEGMs – press the **[New Report]** button. Enter the name of the report using the on-screen keyboard (see figure below). Close the on-screen keyboard by pressing the return button.

**Note:** When a new report is added, all subsequent measurements are saved to that report – until another report is added.

Fig. 40: On-screen keyboard for the Report screen

---

## Delivering Burst Pacing

### Description

Burst pacing can be used to terminate certain types of arrhythmias. During burst pacing, the currently selected channel is stimulated with a fixed pulse amplitude and an adjustable pacing rate. Burst pacing is possible on all three channels, but on only one channel at a time.

**Note:** During burst pacing the keys **[Save...]**, **[Max Stim]** and **[Pref. Set.]** on the front panel are disabled.

### Settings

The burst pacing is performed with the settings shown in the table below.

#### Burst pacing settings

Parameters	Setting
Pacing rate, preset	100 ppm (600 ms)
Pacing rate, adjustable range	80 ... (10) ... 1000 ppm (750 ... 60 ms)
Pulse amplitude, fixed setting	7.5 V
Pulse width, fixed setting	1.0 ms
Duration	max. 30 s

### Backup program

While burst pacing in the atrium, the device will provide ventricular backup pacing if the pacing mode is a ventricular pacing mode when starting the burst. If the pacing mode is not a ventricular pacing mode (e.g., AAI), ventricular pacing support will not be provided during the burst.

If sensing is switched off in the ventricle, e.g. in the D00 mode, the backup pacing is performed asynchronously.

## Delivering burst pacing

To deliver burst paces in a channel, proceed as follows:

- Press the **[Burst]** button. A confirmation window is displayed.

Fig. 41: Burst pacing: Confirmation window

### WARNING

#### Danger from induction

Always have cardiac emergency equipment (e.g., external pacemaker, defibrillator) in an operational status available for immediate life support. Consider additional pre-emptive measures in patients where loss of pacing could cause life-threatening danger. Only activate the burst pacing when all necessary provisions have been made.

- Press the **[OK]** button to continue, or press the **[Cancel]** button to return to the previous screen. In the **Burst** screen the burst rate buttons **[Rate]/[Interval]** for the channels A, RV and LV are displayed with the current settings.

Fig. 42: Burst pacing: Selecting a channel

- Select the channel in which to deliver the burst paces by pressing the corresponding burst rate button **[A]**, **[RV]** or **[LV]**. The burst deliver button **[Burst A]**, **[Burst RV]** or **[Burst LV]** for that channel is displayed.


Fig. 43: Burst pacing: Delivering burst paces (here: in channel A)

- Change the burst rate **[Rate]/[Interval]** for that channel using the control knob on the front panel.
- Deliver the burst paces by pressing and holding the burst deliver button **[Burst A]**, **[Burst RV]** or **[Burst LV]**.

**Note:** The maximum duration of a burst is limited to 30 seconds.

#### Terminating burst pacing

- Terminate burst pacing by releasing the burst deliver button.

## Delivering Maximum Pacing

### Description

Maximum pacing is used to assess whether the extra-cardiac muscle can be inadvertently stimulated by paces with the lead in its current position within the heart. During maximum pacing, the currently selected channel is paced with maximum amplitude and pulse width.

### Settings

The maximum pacing is performed with the settings shown in the table below.

#### Maximum pacing settings

Parameters	Setting
Pulse amplitude, fixed setting	10.0 V
Pulse width, fixed setting	2.0 ms

### Delivering maximum pacing

To deliver maximum pacing in a channel, proceed as follows:

- If the channel is not already active, activate the channel by pressing the channel ON/OFF button **[A]**, **[RV]** or **[LV]**.
- If pacing is not enabled in the channel, enable pacing by pressing the **[Pace]** button.
- Select the channel in which to deliver maximum pacing.
- Deliver maximum pacing by pressing the **[Max Stim]** key on the front panel and holding it. While maximum pacing in a channel, the pace amplitude and pulse width are displayed on the channel's parameter buttons.

Fig. 44: Maximum pacing: Max Stim key on the front panel


Fig. 45: Maximum pacing: Activating the pacing (here: in channel RV)

### Terminating maximum pacing

- Terminate maximum pacing by releasing the **[Max Stim]** key on the front panel.

## Delivering the Safe Program (Safe VVI Pacing)

### Description

The Safe program function sets critical pacing parameters to preset, safe values.

**Note:** You can initiate the Safe program in every operating mode of the device. The Safe program immediately interrupts all other functions of the device.

### Settings

The Safe program settings are shown in the table below.

#### Safe program settings

Parameters	Setting
Mode	VVI RV
Pacing rate, preset	60 ppm (1000 ms)
Pulse amplitude, preset	7.5 V
Pulse width, preset	1.0 ms
Sensitivity RV, preset	2.5 mV
Refractory period RV, fixed setting	250 ms

### Delivering the Safe program

To deliver the Safe program, proceed as follows:

Fig. 46: Safe program: Safe program key

- Press the Safe program key  on the front panel. The Safe program is active immediately, without additional confirmation.

## Viewing or Changing System Settings

### Description

System settings are general settings of the device such as date, time, language, brightness, etc.

- A change in the system settings can be made without affecting the current operation of the device.
- All system settings are saved permanently in the device, even after it is turned off.

### Viewing/changing system settings

To change a system setting proceed as follows:

- Go to the Preferred settings screen by pressing the **[Pref. set.]** button.


Fig. 48: System settings in the Preferred settings screen

- To change a setting, select the value and change it using the control knob on the front panel. The new value is immediately effective.

### Date, time, date format, time format

On the report of the saved measured values the corresponding date and time are denoted.

You can set the following values:

- **[Date]**
  - Day, month, year
- **[Time]**
  - Hour, minutes
  - Time format
    - 12h (e.g. 01 : 30 pm)
    - 24h (e.g. 13 : 30)

- **[Date format]**
    - mm/dd/yyyy (American)
    - dd/mm/yyyy (European)(dd = day, mm = month, yyyy = year)
- 

**Language**

Buttons and screen messages can be displayed in different languages.

You can set the following values:

- **[Language]**
    - English
    - Español
    - Français
    - Italiano
    - Deutsch
- 

**Illumination**

The brightness of the display is preset to provide an optimal compromise between light intensity of the battery service time.

**Note:** The service time for battery operation is reduced when the display is made brighter.

You can set the following values:

- **[Illumination]**
    - Long life
    - Normal
    - Bright
- 

**Resetting system settings**

You can reset all system settings to the factory default values.

**Note:** All preferred settings will also be reset (see section "Viewing or Changing Preferred Settings").

To reset all settings to their factory values proceed as follows:

- Press the **[Reset pref. settings]** button. The system settings are changed to the factory settings. The current operational settings are not affected by this change.
-

## Viewing or Changing Preferred Settings

### Description

Preferred settings are a set of user configurable values that can be quickly activated (e.g., recalled) with a press of a single key – the **[Pref. Set.]** key on the front panel.

**Note:** These settings can be changed in the Preferred settings screen. For these settings to take effect, you have to activate them using the **[Pref. Set.]** key on the front panel.

- A change in the preferred settings can be made without affecting the current operation of the device.
- All preferred settings are saved permanently in the device, even after it is turned off.

### Viewing/changing preferred settings

To change a preferred setting proceed as follows:

- Go to the Preferred settings screen by pressing the **[Pref. set.]** button.
- Use the scroll bar to scroll through the available settings.


Fig. 49: Preferred settings in the Preferred settings screen (upper section)


Fig. 50: Preferred settings in the Preferred settings screen (lower section)

- To change an ON/OFF value, tap on the corresponding field: ● = ON, ○ = OFF. The new value is saved.
- To change a parameter, select the value and change it using the control knob on the front panel. The new value is saved.

### Adjustable values

You can set the following values:

### Preferred settings values

Parameters	Range of values	Factory value
A Lead	ON / OFF	ON
A Pacing (Ap)	ON / OFF	OFF
A Sensing (As)	ON / OFF	ON
A Pulse width (APW)	0.1 ... (0.1) ... 2 ms	0.5 ms <sup>a)</sup>
A Amplitude	0.1 ... (0.1) ... 10 V	5 V <sup>a)</sup>
A Sensitivity	0.2 ... (0.1) ... 20 mV	1.0 mV
RV Lead	ON / OFF	ON
RV Pacing (RVp)	ON / OFF	OFF
RV Sensing (RVs)	ON / OFF	ON
RV Pulse width (RVPW)	0.1 ... (0.1) ... 2 ms	0.5 ms <sup>b)</sup>
RV Amplitude	0.1 ... (0.1) ... 10 V	5 V <sup>b)</sup>
RV Sensitivity	0.5 ... (0.1) ... 20 mV	2.5 mV

Parameters	Range of values	Factory value
LV Lead	ON / OFF	ON
LV Pacing (LVp)	ON / OFF	OFF
LV Sensing (LVs)	ON / OFF	ON
LV Pulse width (LVPW)	0.1 ... (0.1) ... 2 ms	0.5 ms <sup>c)</sup>
LV Amplitude	0.1 ... (0.1) ... 10 V	5 V <sup>c)</sup>
LV Sensitivity	0.5 ... (0.1) ... 20 mV	2.5 mV
Tracking	ON / OFF	OFF
Rate (Basic interval)	30 ... (1) ... 100 ... (2) ... 180 ppm (2000 ... 333 ms)	90 ppm (667 ms)
AV delay (AVD)	0 ... (5) ... 300 ms	120 ms <sup>b)</sup>
RV/LV delay (VVD)	-100 ... (5) ... +100 ms	+5 ms <sup>c)</sup>
Polarity LV (LV lead polarity)	- LV bipolar: BIPL - Common ring bipolar: CRBP - LV Inverse bipolar: IVBP - Ring-ring bipolar: RRBP	- LV bipolar: BIPL <sup>c)</sup>
IEGM Filter	ON / OFF	OFF <sup>d)</sup>
	<p>a) <b>Not active, as Ap is OFF</b></p> <p>b) <b>Not active, as RVp is OFF</b></p> <p>c) <b>Not active, as LVp is OFF</b></p> <p>d) <b>In certain environments, switching the IEGM Filters on may result in a reduction of EMI interferences or noise on the IEGM traces. However, activated IEGM Filters may attenuate peak IEGM potentials from cardiac activity, thus influencing the morphology of the IEGM signal. The default setting for the IEGM Filters is off upon start of the device.</b></p>	

**Resetting  
preferred settings**

You can reset all preferred settings to the factory default values.

**Note:** All system settings will also be reset (see section "Viewing or Changing System Settings").

To reset all settings to their factory values proceed as follows:

- Press the **[Reset pref. settings]** button. The preferred settings are changed to the factory settings. The current operational settings are not affected by this change – unless you activate the preferred settings using the **[Pref. Set.]** key on the front panel.
-

## 4

**Appendix****What's in this chapter?**

This chapter contains the following topics:

<b>Topic</b>	<b>Page</b>
Technical Data	88
Fixed parameter values	92
Programmable parameter values	95
Scope of Delivery and Accessories	97
Messages	99
Electromagnetic Compatibility in compliance with IEC 60601-1-2:2007	101
Country-Related Information	105
Symbols on Label	106
Index	108

## Technical Data

### General characteristics

Category	Specification
Designation	Reliaty
Type	PSA – Pacing System Analyzer (Three-Chamber Threshold Analyzer)
Classification	AIMD, according to EC Directive 90/385/EEC

### Physical properties

#### Dimensions

Parameter	Specification
Length x Width x Height	220 x 180 x 60 mm
Weight (incl. 2 filled battery magazines)	1.2 kg
Housing material	ABS flame resistant

#### Display / Touch screen

Parameter	Specification
Size / diagonal	7 in.
Resolution	800 x 480 pixel
Backlight	LED
Screen brightness	350 cd/m <sup>2</sup>
Touchscreen	Resistive
Pace/sense indicators	6 LEDs

#### Input / Output

Parameter	Specification
USB port	<ul style="list-style-type: none"> <li>• USB 2.0 standard (12 MBit/s)</li> <li>• Bluetooth adapter for Windows-Stack or</li> <li>• USB-Memory-Stick</li> </ul>
Bluetooth interface	Support for external HP PCL5 compatible printer

Parameter	Specification
VGA interface	DE-15 DDC2 only for medical monitor
External power supply	Socket for barrel connector 5.5 mm / 2.5 mm
Patient connectors	<ul style="list-style-type: none"> <li>• 2 Redel connectors PKN06</li> <li>• Connection of patient cables according to the accessory list</li> </ul>

### Power supply

#### External power supply

Parameter	Specification
Designation / type	FRIWO FW 7555M/08
Voltage range	100 ... 240 V (AC)
Frequency range	50 ... 60 Hz
Rated input current	400 mA
Rated output voltage	7.5 V (DC)
Rated output current	1.7 A
Safety class	II
Protection class	IP 41
Connector type	Angled coaxial connector
Outer connector diameter	5.5 mm
Inner connector diameter	2.5 mm
Plus pole	Inner
Cable length	2.5 m (model for USA, Canada, China, Taiwan, Japan) 2.0 m (model for Europe and other regions)

## Batteries

Parameter	Specification	Remark
Battery magazines	2	Can be replaced independently
Batteries per magazine	4	
Battery type	Mignon (AA LR6)	
Operation time (with 2 fully loaded magazines = 8 batteries)	> 12 hours	<ul style="list-style-type: none"> <li>• With Duracell MN1500 AA LR6 or Energizer E91 AA LR6</li> <li>• Two chamber stimulation at 70 ppm, 5 V, 0,5 ms, 500 Ohm</li> <li>• No external devices are connected at the VGA and the USB connector</li> </ul>
Remaining operation time at EOS warning	<ul style="list-style-type: none"> <li>• Stage 1: 3 hours</li> <li>• Stage 2: 30 minutes</li> </ul>	<ul style="list-style-type: none"> <li>• Visual display (battery symbol)</li> <li>• Visual display, and warning signal every 20 seconds</li> </ul>

## Environmental conditions

## Shipping and Storage

Parameter	Specification
Temperature	0 ... 50°C 32 ... 122°F
Relative humidity	30 ... 75 % (non-condensing)
Atmospheric pressure	50 ... 106 kPa 7.2 ... 15.4 lbf/in. <sup>2</sup>

## Operation

Parameter	Specification
Temperature	10 ... 40°C 50 ... 104°F
Relative humidity	30 ... 75 % (non-condensing)
Atmospheric pressure	70 ... 106 kPa 10.1 ... 15.4 lbf/in. <sup>2</sup>

## Safety equipment

Category	Specification
Equipment class	CF (cardiac floating)
Defibrillation protection	Yes
Safety class (with external power supply)	II
Protection class – PSA	IP 30
Protection class – external power supply	IP 41
Resistance against impact	According to IEC 60601-1 and UL 60601-1

## Fixed parameter values

### Factory settings

Parameter	Value
Mode	VVI
Active channels	<ul style="list-style-type: none"> <li>• A: ON</li> <li>• RV: ON</li> <li>• LV: ON</li> </ul>
Pace	<ul style="list-style-type: none"> <li>• A: OFF</li> <li>• RV: OFF</li> <li>• LV: OFF</li> </ul>
Sense	<ul style="list-style-type: none"> <li>• A: ON</li> <li>• RV: ON</li> <li>• LV: ON</li> </ul>
Tracking	OFF
AV delay	120 ms
VV delay (RV > LV)	5 ms
Upper tracking rate	425 ms (+/- 2 ms)
Basic rate	666 ms 90 ppm
Atrium: Pulse amplitude	5 V
Atrium: Pulse width	0.5 ms
Atrium: Sensitivity	1.0 mV
Atrium: Refractory period (TARP)	425 ms (+/- 2 ms)
Atrium: Refractory period (PVARP)	425 ms (+/- 2 ms)
RV: Pulse amplitude	5 V
RV: Pulse width	0.5 ms
RV: Sensitivity	2.5 mV
RV: Refractory period	250 ms (+/- 2 ms)
LV: Pulse amplitude	5 V
LV: Pulse width	0.5 ms
LV: Sensitivity	2.5 mV

Parameter	Value
LV: Refractory period	250 ms (+/- 2 ms)
LV: Pacing polarity	BIPL
Date format	mm/dd/yy
Time format	12 h
Language	English
Screen brightness	Normal

**Maximum pacing**

Parameter	Value
Pulse amplitude	10.0 V
Pulse width	2.0 ms

**Safe program**

Parameter	Value
Mode	VVI RV
Active channels	<ul style="list-style-type: none"> <li>● A: ON</li> <li>● RV: ON</li> <li>● LV: ON</li> </ul>
Pace	<ul style="list-style-type: none"> <li>● A: OFF</li> <li>● RV: ON</li> <li>● LV: OFF</li> </ul>
Sense	<ul style="list-style-type: none"> <li>● A: OFF</li> <li>● RV: ON</li> <li>● LV: OFF</li> </ul>
Tracking	ON
Basic rate	1000 ms 60 ppm
Pulse amplitude RV	7.5 V
Pulse width RV	1.0 ms
Sensitivity RV	2.5 mV
Refractory time RV	250 ms

**Intracardiac  
measurements**

<b>Parameter</b>	<b>Value / Range</b>
Measured intervals	Accuracy: +/- 2 ms
R-wave amplitudes	0.5 ... (0.1) ... 30 mV (+/- 10 % or +/- 0.2 mV, whichever is greater)
P-wave amplitudes	0.2 ... (0.1) ... 30 mV (+/- 10 % or +/- 0.2 mV, whichever is greater)
Impedance	100 ... (1) ... 3000 Ohm (+/- 15 % at pace ampli- tude 0.2 ... 1.9 V, +/- 10 % at pace amplitude 2.0 ... 10 V)
Slew rate	0.5 ... (0.1) ... 3.0 V/s (+/- 30 % or +/- 0.3 V/s, whichever is greater)
IEGM capture window (in Timing/Freeze screen)	5 s
Number of stored reports (in Report List)	up to 15 (remain stored when the device is switched OFF)

## Programmable parameter values

### Parameters and value ranges

Parameter	Value range
Mode	DDD, DDI, VDD, VDI, VVI, V00, AAI, A00, OAO, OVO, ODO
Active channels	<ul style="list-style-type: none"> <li>• A: ON, OFF</li> <li>• RV: ON, OFF</li> <li>• LV: ON, OFF</li> </ul>
Pace	<ul style="list-style-type: none"> <li>• A: ON, OFF</li> <li>• RV: ON, OFF</li> <li>• LV: ON, OFF</li> </ul>
Sense	<ul style="list-style-type: none"> <li>• A: ON, OFF</li> <li>• RV: ON, OFF</li> <li>• LV: ON, OFF</li> </ul>
Tracking	ON, OFF
AV delay	0 ... (5) ... 300 ms (+/- 2 ms)
VV delay (RV > LV)	-100 ... (5) ... 100 ms (+/- 2 ms)
Basic rate (A, RV, LV)	30 ... 180 ppm, with 30 ... (1) ... 100 ppm and 100 ... (2) ... 180 ppm (+/- 1 ppm)
Pulse amplitude (A, RV, LV)	0.1 ... (0.1) ... 10 V (+/- 10 % or +/- 0.1 V, whichever is greater, at 500 Ohm)
Pulse width (A, RV, LV)	0.1 ... (0.1) ... 2.0 ms (+/- 10 % or +/- 0.02 ms, whichever is greater)
Sensitivity, atrium	0.2 ... (0.1) ... 20 mV (+/- 10 % or +/- 0.2 mV, whichever is greater)
Sensitivity, ventricle	0.5 ... (0.1) ... 20 mV (+/- 10 % or +/- 0.2 mV, whichever is greater)

<b>Parameter</b>	<b>Value range</b>
LV: Pacing polarity	BIPL CRBP IVBP RRBP
Burst frequency	80 ... (10) ... 1000 ppm (+/- 1 ppm, +/- 2 ms)
Burst time maximum	30 s (+/- 1 s)
High rate protection	200 ppm (+/- 15 ppm)
Blanking after pace – cross-chamber	25 ms (+/- 2 ms)
Blanking after pace – same chamber	125 ms (+/- 2 ms)
Autoshort after pace (A, RV, LV)	15 ms (+/- 2 ms)
Measurement time sense (A, RV, LV)	60 ms (+/- 2 ms)
Noise interval	70 ms (+/- 2 ms)

<b>Parameter</b>	<b>Value / Range</b>
Sense calibration (Test waveforms)	Cenelec

## Scope of Delivery and Accessories

### Scope of delivery

Name	Quantity	Order no.	Details
Reliatty	1	365383	Device (when ordered separately: without battery magazines)
Bag	1	370109	-
Power supply brick	1	367905	FRIWO FW 7555M/08
Power cord	1	351665 for: NK-29 / 2.5 m	
		345475 for: NK-30 / 2.0 m UK	
		342415 for: NK-31 / 2.0 m EU	
		367906 for: NK-32 / 2.0 m AUS	
		379321 for: NK-34 / 2.0 m BR	
Battery magazines	3	365382	2 inside the device
Batteries	8	370736	Duracell MN1500 AA LR6
Patient cables			
- USA	2	128085	PK-67-S
- All other countries	2	353181	PK-141
Adhesive foil (labels) in different languages	4	DE: 370728 ES: 370730 FR: 370731 IT: 370732	
Technical manual	1	Language specific order no.	
Quick reference guide	1	Language specific order no.	

**Available  
accessories**

<b>Name</b>	<b>Order no.</b>	<b>Details</b>
Patient cable PK-67-L	123672	Applicable with PK-155
Patient cable PK-67-S	128085	Applicable with PA-1-B, PA-2, PA-4
Patient cable PK-141	353181	-
Patient adapter PA-1-B	123751	-
Patient adapter PA-2	123157	-
Patient adapter PA-4	123090	-
Patient adapter PK-155	337358	-
NK-29 / 2.5 m US	380080	Power cord

## Messages

### System Messages Meanings and measures

Message text	Meaning	Measures
Error while creating the *.pdf - file!	PDF file is corrupt.	
Storage medium out of space!	Insufficient remaining memory on storage device.	<ul style="list-style-type: none"> <li>• Delete unused data on storage device.</li> <li>• Use storage device with greater capacity.</li> </ul>
Storage device not ready!	Files cannot be exported to storage device.	Ensure the following: <ul style="list-style-type: none"> <li>• Storage device is compatible.</li> <li>• Storage device is properly connected.</li> </ul>
Storage device write-protected!	Files cannot be exported to storage device.	Ensure that the storage device does not have write protection enabled.
Bluetooth device not ready!	Data cannot be transferred to Bluetooth adapter.	Ensure the following: <ul style="list-style-type: none"> <li>• Bluetooth adapter is compatible.</li> <li>• Bluetooth adapter is properly connected.</li> </ul>
Printer communication error!	Data cannot be transferred to Bluetooth printer.	Ensure the following: <ul style="list-style-type: none"> <li>• Printer is Bluetooth enabled, powered on, and within close proximity.</li> </ul>
Printer not ready!	Data cannot be transferred to Bluetooth printer.	Ensure the following: <ul style="list-style-type: none"> <li>• Printer is Bluetooth enabled, powered on, and within close proximity.</li> </ul>

Message text	Meaning	Measures
No storage/bluetooth device connected!	Data cannot be transferred to storage device or Bluetooth adapter.	Ensure the following: <ul style="list-style-type: none"><li>• Storage device or Bluetooth adapter is compatible.</li><li>• Storage device or Bluetooth adapter is properly connected.</li></ul>

**Critical Error Messages**

Triggered by a built-in selftest routine, the device can also show a critical error message. The message consists of a condition code number and a proposed measure. In case the device repeatedly shows such error message, note the condition number and contact your local BIOTRONIK representative.

---

## Electromagnetic Compatibility in compliance with IEC 60601-1-2:2007

- As the user, you must ensure that the device is operated in a suitable electromagnetic environment.
- The following guidelines may not be applicable in all cases. The propagation of electromagnetic values is, for example, affected by absorption and reflection by structures, objects and people. This data is for your personal information.

Devices with the warning sign "Beware of non-ionizing radiation" must not be operated in the environment of the device due to potential interferences.

### Electromagnetic Emissions (Table 1)

Measuring the emitted interference	Compliance	Guidelines for the electromagnetic environment
High-frequency interference according to CISPR 11	Group 1	The device uses HF energy exclusively for its own function. Therefore, the HF interference is very low and not likely to cause any interference in nearby electronic equipment.
High-frequency interference according to CISPR 11	Class B	The device is suitable for use in all establishments. This includes residences and facilities directly connected to the public power supply network that supplies buildings used for domestic purposes.
Interference of harmonics according to IEC 61000-3-2	Class A	
Emitted interference of voltage fluctuations according to IEC 61000-3-3	Complies	

**Recommended safety distances (Table 6)**

- Safety distances help prevent interference if you maintain a minimum distance between transmitters such as mobile RF telecommunication devices and the device. The necessary distance depends on the respective power output of the transmitter.

**Note:** At 80 MHz and at 800 MHz, the higher frequency range applies.

Transmission frequency	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
Maximum output power of the transmitter [W]	Safety distance [m]		
0.01	0.12	0.12	0.24
0.1	0.37	0.37	0.74
1	1.17	1.17	2.34
10	3.70	3.70	7.40
100	11.7	11.7	23.4

- For transmitters whose maximum output power is not indicated in the table, the recommended safety distance *d* can be calculated in meters using an equation that is suitable for the respective transmission frequency range. *P* is the maximum output power of the transmitter in watts[W] according to the specification of the transmitter’s manufacturer.

Transmission frequency	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
Equation	$d = 1.17 \sqrt{P}$	$d = 1.17 \sqrt{P}$	$d = 2.34 \sqrt{P}$

**Electromagnetic Interference Resistance (Tables 2 and 4)**

- When the measured field strength exceeds the specified compliance level at the operating location of the device, observe the device in order to determine whether it is functioning properly.
- If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device. In the frequency range of 150 kHz to 80 MHz, ensure that field strengths are lower than 3 V/m.

**Note:**  $U_T$  is the mains alternating voltage before applying the test levels.

Test of resistance to interference	Test level according to IEC 60601-1-2	Compliance	Guidelines for the electromagnetic environment
Electrostatic discharge (ESD) according to IEC 61000-4-2	$\pm 6$ kV contact discharge $\pm 8$ kV air discharge	Same as test level	<ul style="list-style-type: none"> <li>Operate the devices on floors made of wood, concrete, or ceramic tile. If the floor is covered with synthetic material, the relative humidity must be at least 30 %.</li> </ul>
Fast transient electric interferences (bursts) according to IEC 61000-4-4	$\pm 2$ kV for power supply lines $\pm 1$ kV for input and output lines	Same as test level	<ul style="list-style-type: none"> <li>Ensure that the power supply quality is that of a typical commercial and/or hospital environment.</li> </ul>
Surges according to IEC 61000-4-5	$\pm 1$ kV push-pull voltage $\pm 2$ kV for common-mode voltage		
Voltage drops, brief interruptions and fluctuations in the supply voltage according to IEC 61000-4-1	$< 5\%$ $U_T$ for 1/2 cycle $> 95\%$ drop $40\%$ $U_T$ for 5 cycles $60\%$ drop $70\%$ $U_T$ for 25 cycles $30\%$ drop $< 5\%$ $U_T$ for 5 s $> 95\%$ drop	Same as test level	<ul style="list-style-type: none"> <li>Ensure that the power supply quality is that of a typical commercial and/or hospital environment.</li> <li>If you require continued operation during power supply interruptions, connect the device to an uninterruptible power supply or use a battery for operation.</li> </ul>

Test of resistance to interference	Test level according to IEC 60601-1-2	Compliance	Guidelines for the electromagnetic environment
Magnetic field at the supply frequencies (50/60 Hz) according to IEC 61000-4-8	3 A/m	Same as test level	<ul style="list-style-type: none"> <li>• Ensure that the magnetic field strengths are at levels characteristic of a location in a typical commercial and/or hospital environment.</li> </ul>

Test of resistance to interference	Test level according to IEC 60601-1-2	Compliance	Guidelines for the electromagnetic environment
Conducted HF interferences according to IEC 61000-4-6	3 V <sub>eff</sub>	3 V	<ul style="list-style-type: none"> <li>• Maintain safety distance of mobile radio equipment to the device; see table 6.</li> </ul>
Radiated HF interferences according to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	<ul style="list-style-type: none"> <li>• The field strength of stationary transmitting devices must be measured on site and must be lower than the compliance level at all frequencies: consider conducting a study of the site.</li> <li>• The field strength must be lower than 3 V/m over the frequency range of 150 kHz to 80 MHz.</li> </ul>

## Country-Related Information

---

### UL certification

Reliaty has been certified by Underwriters Laboratories Inc. with respect to electrical shock, fire and mechanical hazards only in accordance with:

- UL 60601-1
- CAN/CSAC22.2 No. 601.1
- IEC 60601-1: 1988 + A1 + A2
- IEC 60601-2-31: 1994 + A1

UL-certified devices are identified as follows:

## Symbols on Label

### Device label

The label icons symbolize the following:

	Pacing system analyzer
	Catalogue number BIOTRONIK
	Serial number of device
	Date of manufacture
	Storage and transport temperature limitation
	Storage and transport atmospheric pressure limitation
	Storage and transport humidity limitation
	CE conformity marking
	Caution: Federal (U.S.A.) law restricts this device to sale by, or on the order of, a physician.
	Contents
	Consult instructions for use
	Device contains materials that must be correctly disposed of in accordance with environmental protection regulations. European Directive 2002/96/EC regarding waste electrical and electronic equipment (WEEE) applies. Return devices that are no longer used to BIOTRONIK.

	Battery magazine
	External power supply (power supply brick plus power cord)
	Patient with implanted lead
<b>Distributor:</b>	Distributor
	Manufacturer

## Index

---

### A

- Accessories, 23, 97
  - connecting, 23
- Ambient conditions, 10, 11
- Applying language foil, 12

### B

- Battery status, 31
- Bluetooth adapter, connecting, 37
- Burst Pacing, 74
  - delivering, 75
  - settings, 74
  - terminating, 76

### C

- Cable and lead connections, 13
- Cables, 36
  - connecting, 20
  - disposal, 20, 44
  - touching clips, 20
- Characteristics, general, 88
- Cleaning, 16
- Complications, possible, 9
- Conduction Times
  - performing the measurement, 64
  - terminating the measurement, 66
  - Timing/Freeze screen, 65
- Connections, 36
- Contraindications, 9
- Control knob, 47
- Critical Error Messages, 100

### D

- Default settings, 13
- Device
  - overview, 26
  - switching on/off, 35
- Disinfecting, 16
- Display, 45, 46
  - buttons, 46
  - illumination, 47
- Disposal, 16, 18

**E**

- Electromagnetic Compatibility, 101
- Electromagnetic Emissions, 101
- Electromagnetic Interference Resistance, 102
- Electromagnetic interferences, 14
- Electrostatic potentials, 21
- Emergency equipment, 19
- Environmental conditions, 90
- Expertise, required, 8
- Exporting, 70
- External monitor, connecting, 36

**F**

- Fixed parameter values, 92
  - factory settings, 92
  - intracardiac measurements, 94
  - maximum pacing, 93
  - safe programm, 93

**I**

- Indication, 6
- Indicators, front, 46
- Initial operation, 13
- Inspections, 17, 18
- Installation, 10
- Installation site, 11
- Intended use, 6, 7
- Intracardiac measurements
  - fixed parameter values, 94

**K**

- Keys, front, 45

**L**

- Leads
  - connecting, 20
  - touching clips, 20

**M**

- Maintenance, 16
- Maximum Pacing, 77
  - delivering, 77
  - fixed parameter values, 93
  - settings, 77
  - terminating, 78
- Measuring
  - Conduction Times, 64
  - P/R Wave Amplitudes, 60
  - Pacing Threshold and Lead Impedance, 62
  - Wenckebach Point, 67
- Medical product operator's ordinance, 8
- Messages
  - Critical Error Messages, 100
  - System Messages, 99

**O**

- Operating elements
  - all screens, 53
  - back, 27
  - front, 26
  - handling, 52
  - IEGM window, 57
  - left, 27
  - left in the Test screen, 54
  - right, 26
  - right in the Test screen, 56
  - Timing/Freeze screen, 59

**P**

- P/R Wave Amplitudes, 60
  - performing the measurement, 60
  - terminating the measurement, 61
- Pacing Threshold and Lead Impedance
  - performing the measurement, 62
  - terminating the measurement, 63
- Patient cables
  - connecting, 40
  - overview, 39
- Physical properties, 88

- Power supply, 12, 30
  - batteries, 30
  - battery status, 31
  - connecting, 30
  - external, 28, 30
  - technical data, 89
- Preferred Settings, 83
  - adjustable values, 84
  - changing, 83
  - resetting, 86
  - viewing, 83
- Printing, 70
- Programmable parameter values, ranges, 95
- Properties, 88

## R

- Recommended safety distances, 102
- Replacement parts, 23
- Report
  - detailed, 72
  - exporting, 37
  - list, 72
  - printing, 37
  - printing/exporting, 70
  - summary, 71

## S

- Safe Programm, 79
  - delivering, 79
  - fixed parameter values, 93
  - settings, 79
  - stopping, 80
- Safety equipment, 91
- Safety Instructions, general, 19
- Scope of Delivery, 97
- Self-test, 21, 35
- Service, 16
- Setup, 29
- Setup Requirements, 10
- Shipping, 10
- Signals, 45
  - acoustic, 47
- Software
  - overview, 50

Start parameters, 13  
Start screen, 35  
Sterilization, 17  
Storage, 10  
Supervision, physician, 19  
Symbols on Label, 106  
System Messages, 99  
System Settings, 81  
    changing, 81  
    date, 81  
    date format, 81  
    illumination, 82  
    language, 82  
    resetting, 82  
    time, 81  
    time format, 81  
    viewing, 81  
System, description, 6

## **T**

Technical Data, 88  
Technical safety check, 17, 18  
Timing/Freeze, operating, 65  
Transport, 29

## **U**

USB flash memory stick, connecting, 38

## **W**

Wenckebach Point  
    performing the measurement, 67  
    terminating the measurement, 69  
Timing/Freeze screen, 68

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All rights reserved. Specifications are subject to modification, revision and improvement.