

# User manual

## **Check-Direct ESBL for BD MAX<sup>™</sup>**

Real time PCR kit for the detection of ESBL-producing Enterobacteriaceae

Version 1.0 Date of issue: 10.07.2014



## EU: **C € IVD**

U.S.: For Research Use Only Not for use in diagnostic procedures

## Contents

Intended use	2
Introduction and principle of the method	2
Kit contents (for 24 reactions)	2
Materials required but not supplied with the kit	2
Storage, handling, and stability	3
Good laboratory practices	3
Sample preparation procedures	4
BD MAX operation	4
Results Interpretation	5
Frequently asked questions (FAQ) & Troubleshooting	6
Limitations	7
Key to symbols used	7
Technical assistance	7
Appendix 1: creating the Check-Direct ESBL test program	8
Appendix 2: performance characteristics	9



## Intended use

Check-Direct ESBL for BD MAX is a qualitative *in vitro* diagnostic test for the rapid detection of Extended Spectrum Beta-Lactamase (ESBL) genes in *Enterobacteriaceae*. The assay is performed on the BD MAX<sup>TM</sup> system using bacterial cell suspensions from clinical specimens of individuals at risk of colonization with ESBL-producing *Enterobacteriaceae*. Check-Direct ESBL detects the presence of the ESBL gene families CTX-M1, CTX-M2, CTX-M9 and SHV-ESBL, presently the primary cause of extended spectrum  $\beta$ -lactam resistance in *Enterobacteriaceae*. The assay uses the BD MAX system for extraction of DNA and subsequent real-time PCR employing the reagents provided combined with universal reagents and disposables for the BD MAX system. Check-Direct ESBL for BD MAX can be used as an aid to identify, prevent and control ESBL-producing *Enterobacteriaceae* that colonize patients in healthcare settings. Check-Direct ESBL for BD MAX is not intended to diagnose infections with ESBL-producing *Enterobacteriaceae* nor to guide or monitor treatment for these infections. Parallel cultures are necessary to recover organisms for epidemiological typing, susceptibility testing and for further confirmatory identification.

## Introduction and principle of the method

The worldwide emergence and dissemination of extended spectrum  $\beta$ -lactam resistance in *Enterobacteriaceae* is a serious threat to public health. These organisms are associated with high mortality rates and have the potential to spread widely. The most common cause of extended spectrum  $\beta$ -lactam resistance in *Enterobacteriaceae* is the production of Extended Spectrum Beta-Lactamases, i.e. ESBL's. Most Enterobacteriaceae producing ESBL's have resistance to all  $\beta$ -lactam antibiotics except carbapenems. Presently, the vast majority of clinically relevant ESBL's are expressed from one of the following plasmid-encoded gene families: CTX-M1, CTX-M2, CTX-M9 and SHV-ESBL.

Check-Direct ESBL is a multiplex real-time PCR assay for detection of the CTX-M1, CTX-M2, CTX-M9 and SHV-ESBL gene families. The assay is based on specific recognition and amplification of target sequences by PCR, and the simultaneous detection of the accumulation of PCR amplification products by fluorescent DNA probes. The CTX-M1, CTX-M2, CTX-M9 and SHV-ESBL gene families have many gene variants and Check-Direct ESBL has been designed to reliably detect all variants. Check-Direct ESBL for BD MAX employs five different fluorescent probes and enables detection and discrimination of the 4 ESBL gene families and the control target SPC, that monitors DNA extraction and PCR amplification.

## Kit contents (for 24 reactions)

Components (Mat. No.)	Description	Storage conditions
ESBL reagent tubes (9-0065)	24 sealed tubes (red seal)	+ 4°C, store in the dark
ESBL positive control (9-0064)	1 tube (red cap)	+ 4°C
User Manual (9-0109)	Leaflet – download from website	Not critical

## Materials required but not supplied with the kit

Supplies	Equipment
<ul> <li>BD MAX ExK<sup>™</sup> DNA-1 Extraction Kit (Ref:442818)</li> <li>BD MAX DNA MMK Master Mix (Ref: 442848)</li> <li>BD MAX PCR Cartridges (Ref: 437519)</li> <li>Disposable laboratory (powder-free) gloves/Lab coat</li> <li>Pipettes &amp; disposable (filter-) tips for volumes of 10 to 1000 µl</li> <li>Milli-Q water or aqua bidest</li> </ul>	<ul> <li>Real-time PCR instrument: BD MAX System, software version 2.96A</li> <li>Densitometer suitable for bacterial suspensions</li> <li>Vortex mixer</li> </ul>



## Storage, handling and stability

The Check-Direct ESBL kit is shipped at ambient temperature and should be stored at + 4°C upon receipt. Please visually inspect the product upon initial opening to ensure that its contents are intact. Do not use this product if the packaging is damaged upon arrival and do not use reagents if their protective pouches are open or broken upon arrival. Do not use reagents if desiccant is not present or broken inside and do not remove desiccant from protective pouches. Store all opened reagents at + 4°C until expiration date. Store in the dark. Close protective pouches promptly with the zip seal after each use. Remove any access air in the pouches prior to sealing. Please contact the Check-Points office at *support@check-points.com* if you have any further questions.

## **Good laboratory practices**

## **Recommendations for best results**

The quality of the results depends on strict compliance with the following good laboratory practices, especially concerning PCR practices.

- The test must be performed by adequately trained personnel.
- Do not use reagents after their expiration date
- Follow recommendations for storage and handling to preserve the quality of the kit's reagents.
- Protect reagents from light to avoid photo-bleaching of the dyes.
- Periodically, verify the accuracy and precision of pipettes, as well as correct functioning and calibration of the instruments.

### **Prevention of contaminations**

**Use separate rooms**: a sample preparation room and a PCR room with the BD MAX system. Never transfer items from the PCR room to the sample preparation room.

#### To keep laboratory free of PCR product contamination:

- Use pipettes with hydrophobic filter tips.
- Make sure to always use a new pipette tip when adding solutions, test samples, and controls to a reaction tube to avoid contamination.
- Follow proper pipette-dispensing techniques to prevent aerosols.
- Wear clean disposable gloves and clean lab coats for the different steps of the test.
- Change gloves whenever you suspect that they are contaminated.
- Keep the tubes of all kit components and samples closed as much as possible.
- Clean the lab benches and all equipment regularly with a 0,5% sodium hypochlorite solution.

## Please read the full protocol before starting the test

## Sample preparation procedures

#### Test preparation for bacteria from culture

1. Inoculate nutrient agar plates with the clinical samples or the bacterial strains to be tested and incubate overnight at 37°C. Typical growth media include blood agar, MacConkey agar and Tryptic Soy agar.

Check-Direct ESBL

- 2. Prepare a bacterial cell suspension of McFarland 0.5 ( $\approx 1 \times 10^8$  CFU/ml) from the bacterial cells grown on the agar plate. Cell suspension buffer (*e.g.* PBS or 10mM Tris.HCl pH8.0) or Milli-Q water or aqua bidest may be used.
- 3. Pipette 500 of Milli-Q water or aqua bidest and 10  $\mu$ L of the bacterial cell suspension ( $\approx$ 1 x 10<sup>6</sup> CFU) into one DNA Sample Buffer Tube SB-1. (supplied by BD with the DNA extraction kit, refer to *Materials required but not supplied with the kit*).
- 4. Close the Sample Buffer Tube with a septum cap and vortex 10 second at low speed.
- 5. Transfer the Sample Buffer Tube with the bacterial cell suspensions to be analyzed to the PCR room.

#### **Preparation of control reactions**

To validate the run, perform positive and negative control reactions for each Check-Direct ESBL PCR run. The positive control is supplied with the kit.

Positive control:

Pipette 10  $\mu$ L of the positive control and 500  $\mu$ L of Milli-Q water or aqua bidest into one Sample Buffer Tube. Vortex for 10 seconds.

• Negative control:

Pipette 500 µL of Milli-Q water or aqua bidest into one Sample Buffer Tube. Vortex for 10 seconds.

## **BD MAX operation**

#### 1. Multiplex real-time PCR setup

Table 1 presents the multiplex real-time PCR setup with the targets detected in each detector channel of the BD MAX System.

Table 1: Multiplex qPCR setup

Detector	475/520	530/565	585/630	630/665	680/715
Channel	1	2	3	4	5
Target	CTX-M1	CTX-M2	SHV-ESBL	СТХ-М9	SPC*

\*SPC: Sample Processing Control

When the test is performed for the first time create the PCR test program "Check-Direct ESBL" as described in Appendix 1.



Figure 1: DNA Unitized Reagent Strip setup.



#### 2. BD MAX Rack set-up

- 2.1. Load the BD MAX system racks with the number of DNA Unitized Reagents Strips necessary for the number of samples to test. Gently tap each strip to make sure all liquids are at the bottom of their container.
- 2.2. Prepare Unitized Reagents Strips:
- 2.2.a. Snap a DNA extraction BD Exk-1 Reagent tube (white seal) into position 1 of the DNA Strip, see Figure 1.
- 2.2.b. Snap a DNA MMK Master Mix tube (green/yellow seal) into position 2 of the DNA Strip, see Figure 1.
- 2.2.c. Snap a ESBL reagent tube (red seal) into position 3 of the DNA Strip, see Figure 1.

#### 3. BD MAX instrument set-up

- 3.1 Open the Run screen of the BD MAX System software v2.96A.
- 3.2 In the Assay menu select Check-Direct ESBL.
- 3.3 Enter the Sample Buffer Tube barcode using the barcode scanner (you can also enter the barcode manually). Start with position 1 of rack A.
- 3.4 Place each of the Sample Buffer Tubes in their corresponding position in the BD MAX racks (with septum cap).
- 3.5 Enter the specimen or patient identification information into the work list. Check that each specimen or patient information correspond to its specific Sample Buffer Tubes in the Rack.
- 3.6 Load the Rack(s) into the BD MAX System. (Rack A is positioned on the left side of the instrument and Rack B on the right side).
- 3.7 Load the BD MAX PCR cartridge(s).
- 3.8 Close the instrument door and select Start Run.

## **Results Interpretation**

**Important points before starting**: For a detailed description on how to analyze data, refer to *BD MAX™ System User's manual*.

Always visually inspect the amplification plot for each sample tested versus C<sub>T</sub> values obtained with the software.

#### **1. Reported results**

The BD MAX<sup>m</sup> software reports C<sub>T</sub> values and amplification curves for each detector channel of each specimen tested in the following way:

- $C_T$  value of **0** indicates that there was no  $C_T$  value calculated by the software. Amplification curve of the sample showing a "0"  $C_T$  value must be checked manually.
- $C_{T}$  value of -1 indicates that no valid amplification process has occurred. Check that there is no amplification curve for the sample with a  $C_{T}$  value of -1 on the graphical results.
- Any other  $C_{T}$  value should be interpreted in correlation with the amplification curve (**PCR Analysis** tab) and according to the interpretation method outlined in Tables 2 and 3.

#### 2. Interpretation

#### 2.1 Run validation

Verify that the real-time PCR run is valid before data interpretation of the results. Check that there is no report of BD MAX System failure. Check the positive and negative control amplification curves. Table 2 shows criteria for a valid real-time Check-Direct ESBL run on the BD MAX<sup>TM</sup> System. If the  $C_T$  values of the controls are not as expected refer to FAQ and Troubleshooting **"3"**.

 Table 2: Criteria for a valid run with Check-Direct ESBL test.

Sample Type*	Ст 475/520 СТХ-М1	С <sub>т</sub> 530/565 СТХ-М2	С <sub>т</sub> 585/630 SHV-ESBL	Ст 630/665 СТХ-М9	С <sub>т</sub> 680/715 SPC
Positive controls	33 ±3	30 ±3	30 ±3	30 ±3	29 ±3
Negative sample	-1	-1	-1	-1	29 ±3

#### 2.2 Results interpretation

If the run has been validated, interpret results as positive, negative or invalid with the  $C_T$  values obtained for the samples following the guidelines summarized in Table 3. Invalid runs should be retested.



Ct-values of bacterial cells will be in a specific Ct-window for each target because of the well-defined amount of cells used as input material for the test. However, Ct-values will differ between individual strains. Ct-values out of this Ct-window are regarded as invalid and should be retested. A higher Ct-value above the Ct-window suggest contamination of the sample or a strain that is not pure. A lower Ct-value suggests an aberrant amplification plot. Please check the amplification plot in such a case to confirm this.

С <sub>т</sub> 475/520 СТХ-М1	Ст 530/565 СТХ-М2	C <sub>T</sub> 585/630 SHV-ESBL	Ст 630/665 СТХ-М9	C <sub>τ</sub> 680/715 SPC	Interpretation
≤ 28	≤ 28	≤ 30	≤ 29	29 ±3	Positive
-1	-1	-1	-1	29 ±3	Negative
> 28	> 28	> 30	> 29	29 ±3	Invalid
-1 or YES	-1 or YES	-1 or YES	-1 or YES	<26 or >32	Invalid
-1 or YES	-1 or YES	-1 or YES	-1 or YES	-1	Invalid

#### Table 3: Data interpretation guidelines for bacterial cells

## Frequently asked questions (FAQ) & Troubleshooting

Refer to "the troubleshooting" section of the BD MAX<sup>™</sup> System User's Manual for additional information

- **1.** Real-time results show no C<sub>T</sub> values or interpretation indicates that the sample is invalid. Possible causes and troubleshooting:
  - The PCR reaction has been inhibited by exogenous or endogenous substances. Please repeat sample testing. When still inhibited a lower amount of input sample may improve the results.
  - The DNA extraction failed since the SPC was not detected.
  - The BD DNA MMK may have expired.
  - An error in liquid handling has occurred: check unitized reagent strips and PCR cartridge to determine where liquid handling problem has occurred (example: air bubble in the cartridge) and re-run the sample. If the problem persists, contact your local BD representative.

#### 2. Troubleshooting for invalid results.

For Invalid results: Repeat test with the original specimen by preparing a new Sample Buffer Tube. Alternatively, test newly collected specimen or use a lower amount of specimen.

- **3.** Real-time results show no C<sub>T</sub> values for the positive control or interpretation indicating that sample is invalid? Possible causes and troubleshooting:
  - The positive control solution was not added.
  - The BD DNA MMK may have expired.
  - Air bubbles have occurred in the PCR reaction chamber of the positive control.
- 4. Real-time results show very low fluorescent signals in all samples and detector channels including the SPC signal.

Possible causes and troubleshooting:

- The ESBL reagent tubes containing the fluorescent probes and primers may be degraded. Please check expiration date and make sure that the ESBL tubes have been stored correctly.
- The BD MAX<sup>™</sup> System can be responsible for these results. Please refer to BD MAX<sup>™</sup> User's manual or contact your BD local representative.

#### 5. The BD MAX<sup>™</sup> System states an error or failure.

Refer to the BD MAX<sup>™</sup> instrument user manual or contact your BD local representative.

#### 6. Duplicate samples tested with Check-Direct ESBL test do not yield identical results.

 $C_{T}$  values of identical samples may vary between individual reactions. Large variations, > 2  $C_{T}$  values, suggest pipetting errors or other differences between the duplicate samples.



## Limitations

Check-Direct ESBL uses a range of specific DNA markers to detect the presence of the CTX-M1, CTX-M2, CTX-M9 and SHV-ESBL genes, which currently represent the clinically most prevalent ESBL's. The test detects almost all clinically relevant variants of CTX-M1, CTX-M2, CTX-M9 and SHV-ESBL. It should be noted that some more rare variants of CTX-M1, CTX-M2, CTX-M9 and SHV-ESBL. It should be noted that some more rare variants of CTX-M1, CTX-M2, CTX-M9 and SHV-ESBL are not detected and that other ESBL gene families like e.g. TEM, GES, VEB and PER are also not detected. Bacterial cell suspensions should be used as input material for the test.

The quality of the input DNA is an important factor for obtaining reliable results with Check-Direct ESBL. For cell suspensions the correct cell densities are an important factor to obtain reliable results and the procedure described in this manual must be strictly followed. The assay has been tested extensively with DNA purified from gram-negative bacteria, such as *Escherichia, Salmonella, Klebsiella, Enterobacter, Citrobacter* and *Pseudomonas,* with excellent results. However, it may never be excluded that other Gram-negative bacteria or certain strains of the above species will yield poor results. Check-Direct ESBL cannot and does not make any representation or warranty that it is capable of correctly detecting the ESBL genes in all gram-negative species, subspecies or types or in all clinical samples. Results may need to be confirmed by additional methodologies in specific cases (e.g. for regulatory samples). Due to the high variability of bacterial genomes it is possible that certain subtypes might not be detected. The test reflects the state of knowledge of Check-Points Health B.V.

The presence of multiple bacterial species in a sample may hamper the interpretation of the test. As with other diagnostic assays, the results of this test may only be interpreted in combination with additional laboratory and clinical data available to the responsible person. Use of this assay is limited to appropriately qualified personnel, well-trained in performing DNA-based molecular detection methods.

Symbol	Definition			
CONTROL ESBL	ESBL control			
IVD	For In Vitro Diagnostic Use			
REF	Catalog number			
LOT	Batch code			
IFU	IFU number			
8	Use before YYYY-MM			
(II)	Consult instructions for use			
(A)	Manufacturer			
X	Temperature limitation			
$\overline{\mathbb{V}}$	Contains sufficient for < n > tests			

## Key to symbols used

## **Technical assistance**

### support@check-points.com +31 317 453 908

Despite the utmost care in the development and preparation of the protocol Check-Points cannot take any responsibility for errors, omissions and/or future changes herein.

Literature Citation: When describing a procedure for publication using this product, please refer to it as the Check-Direct ESBL.

#### Notice to Purchaser:

This product is sold under license from PHRI Properties and may be used under PHRI Properties patent rights only for human in vitro diagnostics, food testing, veterinary testing, or research.

Dye & quencher compounds in this product are sold under license from Biosearch Technologies, Inc. and protected by U.S. and world-wide patents either issued or in application. The license grant covers human *in vitro* diagnostic (IVD) applications **Trademarks** 

BD, BD MAX are trademarks Becton Dickinson GmbH

**Check-Points Health BV** Binnenhaven 5 6709 PD Wageningen The Netherlands Tel: +31 317 453 908 Fax: +31 317 210 147 info@check-points.com www.check-points.com







## Appendix 1: Creating the Check-Direct ESBL test program

**Important points before starting:** Refer to BD MAX System User's Manual for detailed instructions on how to operate the BD MAX System and software version 2.96A.

1. Create a new Test, select **Create test**, and enter the following parameters:

- <u>Test Name</u>: type Check-Direct ESBL
- <u>Extraction Type</u>: Select *Exk DNA-1 (Plasma/Serum)*
- Master Mix Format: choose Type 1: BD MMK or MMK(SPC) and Dried Primers & Probes
- Channel detector Settings: set Gain and Threshold with parameters presented in Table A
- GardRail: select Default
- <u>Test details</u>: enter the PCR profile, see Table B
- Spectral Cross Talk tab: enter parameters presented in Table C

2 Select Save Test.

Table A: Gain parameters.

Detector	Gain	Threshold
475/520	80	100
530/565	80	100
585/630	30	200
630/665	80	100
680/715	30	150

Table B: Real-time protocol parameters.

Step Name	Profile Type	Cycles	Time (s)	Temp(°C)	Detect
Denaturation	Hold	1	600	98	NO
Amplification & Detection	2 - temperature	45	15	98	NO
		45	62	60	YES

Table C: Spectral cross-talk parameters.

	False Receiving Channel					
		475/520	530/565	585/630	630/665	680/715
	475/520		1.6	0.0	0.0	0.0
Excitation	530/565	1.2		1.2	0.0	0.0
Channel	585/630	0.0	0.0		10.7	0.0
Channel	630/665	0.0	0.0	6.3		5.6
	680/715	0.0	0.0	0.0	7.8	



## **Appendix 2: Performance Characteristics**

### Limit of Detection with control DNA targets

The analytical limit of detection (LoD) of Check-Direct CPE was determined using linearized control plasmids containing target sequences for CTX-M1, CTX-M3 & M15, CTX-M10, CTX-M2, CTX-M9, SHV-2 and SHV-5 & -12. Of each target DNA 6 concentrations were tested, i.e.  $10^5$ ,  $10^4$ ,  $10^3$ ,  $10^2$ , 10, 1 DNA copies/µl. Samples were tested on the BD MAX system in the PCR only mode. First Check-Direct ESBL Snap tubes were rehydrated with 12.5µl of a solution containing Tris EDTA buffer, MMK PCR diluent and SPC target DNA. Next 12.5 µl of target DNA solution was added. Subsequently, the resulting 25µl of each reaction mix this mix was transferred to the MMK tube to rehydrate the master mix. Finally, 12µl of each solution was manually transferred into the individuals well of a BD MAX PCR cartridge. The cartridge was then placed into the BD MAX instrument and the PCR only protocol was started using the Check-Direct ESBL test program.

The assays with  $10^5$ ,  $10^4$ ,  $10^3$ ,  $10^2$  DNA copies/µl were always positive. The assays with 10 and 1 DNA copies/µl reached the true LoD of the test: results are presented in the table below.

Target	Input DNA copies	DNA copies/PCR	Success rate
CTXM-1	12.5	2	4 out of 4
CTXM-2	12.5	2	3 out of 4
CTXM-3/15	12.5	2	2 out of 2
CTXM-10	125	20	2 out of 2
CTXM-9	12.5	2	2 out of 2
SHV-2	12.5	2	2 out of 2
SHV-5/12	12.5	2	2 out of 2

#### In silico Specificity

The specificity of the Check-Direct ESBL real-time diagnostic test is ensured by the selection of the correct primers and probes, as well as the selection of stringent reaction conditions. Primers and Probes sequences were designed to specifically identify the gene variants listed in the table below. A 100% sequence match with the primers and probes by *in silico* analysis was assumed to warrant reliable detection of each of the depicted variants. Single mismatches with the primers and probes exist in some variants, of which we expected that detection would not be compromised. This was confirmed by testing such variants in comparison with variants which were 100% homologous.

Primers and Probes sequences were tested for potential homologies with genes from other organisms using all gene sequences present in the international gene bank on April 1<sup>st</sup>, 2014. (GenBank<sup>®</sup>, NIH genetic sequence database). using sequence comparison analysis. No cross homology was found with other organisms for the selected primers and probes.

ESBL gene	Variants detected
CTX-M1 group	1,3,10-12,15,22,23,28-30, 32-34,36,37,42,52-55,57,58,60-62,66,68,69,71,72,88,96,101,107-109,114,116,117,132 133,136
CTX-M2 group	2,4-7,20,31,35,43,44,56,59,76,77,95,124,131
CTX-M9 group	9,13,14,16-19,21,24,27,38,45-51,55,57,64,65,67,81,93,102,104,105,106,110,111,112,113,121,122,123,126,134
SHV-ESBL	2,2a,3, 4,5,7,9,10,12,15,20-23,30,34,39, 45,46,55,64,66,86,90,105,106,123,124,128,129,134,141,152,153,152,160, 163-165