QUICK START GUIDE

Prepare the Reservoir:



Replace

Replace the reservoir onto the base, PRESSING DOWN FIRMLY to ensure full seating.



Prior to Very First Use of Unit, Prime Unit: Set pressure control slide switch on handle to high, point the handle and tip into sink. Turn unit ON until water flows from the handle, turn unit OFF.





2 Select and Insert Tip:





Press tip firmly into handle. When tip is properly inserted there should be no gap between control knob and tip. To remove tip, simultaneously hold control knob down and lift tip up.

4 Adjust Pressure and Pause Control:

For first time users, adjust pressure control to lowest position. Gradually increase pressure as you prefer.

You can temporarily stop water flow by pressing the pause button on the handle.



continued on back

Using your Waterpik® Water Flosser:

Always place tip in mouth before turning unit on.

Lean over sink and close lins enough to prevent splashing. while still allowing water to flow from mouth into the sink.









Technique Tips for Best Results:

Aim the tip just above the gum line at a 90 degree angle.



Starting with the back teeth, follow the gumline and pause briefly between teeth.

Continue until all areas around and between teeth have been cleaned.

REGISTER YOUR PRODUCT ONLINE NOW!

> Receive up to \$44 in savings.

www.waterpik.com/register