



04520

Please read the following instructions and warnings carefully.  
Keep this instruction manual for future reference.

**For Technical Support Call 1-800-268-6237 or E-Mail at [customerservice@summerinfant.com](mailto:customerservice@summerinfant.com)**

For use with children from birth and up.



## INSTRUCTION MANUAL

# DELUXE BABY'S ESSENTIALS KIT

Thank you for buying the Deluxe Baby's Essentials Kit from Summer Infant. This comprehensive kit contains essential health and hygiene items for your baby. Convenient carry bag keeps these items organized and easily accessible.

## **WARNING:**

**DO NOT TIE MEDI-DISPENSER AROUND CHILD'S NECK AS IT PRESENTS A STRANGULATION DANGER.**

## **CAUTION:**

The components of this product, except for the teether, are not toys. Do not allow children to play with them. Cool teether only in a domestic refrigerator; do not place in the freezer compartment.

---

## BRUSH & COMB:

---

Soft brush bristles are gentle on baby's scalp. Wide-toothed comb makes untangling hair easier. Handle design with soft grip texture is easy for toddlers to hold and use.

**To use:**

- Wash brush & comb before using.
- Untangle hair from the roots out. Keep a hand on the roots while brushing or combing to avoid pulling on child's scalp.
- Encourage your toddler to try brushing his or her own hair with your assistance (approximately 3 years old).

**To clean:** Wash with mild soap and warm water. Rinse thoroughly; air dry.

---

## NAIL CLIPPERS:

---

This standard size nail clipper's handle rotates 360 degrees so you can safely trim nails from any angle. Soft over-molded handle is comfortable and helps you maintain a steady grip.

**Note:** Do not clip too close to baby's skin. We recommend clipping baby's nails while he or she is sleeping.

**To clean:** Sterilize with isopropyl rubbing alcohol and rinse in warm water. Wipe dry.

---

## MEDICINE SYRINGE:

---

Easily dispense and measure medicine.

---

## MEDICINE SYRINGE:

---

**To use:** Clean before first use. Push the inner half of the syringe all the way in and place the tip of the syringe into the liquid medicine. Pull the inner half out until the medicine reaches the correct marked dosage. Make sure your child is elevated to avoid choking and then place one finger inside their mouth and gently pull the cheek outward. Place the syringe toward the back of the cheek (not toward the back of the mouth to avoid gagging) and push down on the inner plunger to dispense the liquid.

**To clean:** Remove the inner plunger and wash both halves in warm, soapy water. Air dry. Top rack dishwasher safe. Do not place in boiling water.

---

## 6 EMERY BOARDS:

---

- Gently glide across baby's fingernails to smooth them.

---

# NASAL ASPIRATOR:

---

**To use:**

- Hold baby upright.
- While squeezing the aspirator bulb, gently place the tip of the aspirator into baby's nostril and tilt it to the side. Never point the tip directly in the center of the nostril. Be careful not to over-insert the aspirator.
- Slowly release the bulb allowing the suction to draw the mucus out of baby's nostril.
- Repeat with other nostril.
- Repeat as necessary until baby's experiences relief from congestion.
- To clean, put the aspirator under soapy, warm water. Squeeze and release the bulb several times to clean the inside of the bulb. Rinse with warm water. Squeeze out all excess water from the bulb.

---

# PACIFIER MEDICINE DISPENSER:

---

The Medicine Dispenser provides parents with a convenient way to dispense medicine to their child. The unique design reduces the likelihood of spillage and helps ensure proper dosage. The soft, contoured applicator is comfortable for baby.

**To use:** Boil in water for 2 minutes before first use. Twist the clear cap counterclockwise and pivot open. Pour in the desired amount of liquid using the dosage marks. Holding the cap upright, pivot the pacifier shield down on the cap, and then twist clockwise to close. Place the Summer medicine dispenser into baby's mouth. Works best when the child is in a reclined position.

**To clean:** Immediately after use, open the cap and wash with hot soapy water. Rinse thoroughly. Dishwasher safe.

---

# TEETHER:

---

- Water-filled teether soothes baby's gums while teething.

---

# FOREHEAD THERMOMETER & CASE:

---

This temperature indicator is reusable and comes with a convenient carrying case. The parent simply holds it in place on the child's forehead and gets a reading in less than 15 seconds.

---

# 5 ALCOHOL SWABS:

---

- Use to clean and disinfect.

---

## EMERGENCY CARD:

---

- Keep important phone numbers and medical information handy for quick reference.
  - Keep in a safe place. Familiarize family members, baby-sitters and care givers with the location of the card and the information on it.
- 

## CARRY BAG:

---

- Keep items organized with the convenient storage case.

**To clean:** Wash with mild soap and warm water. Rinse thoroughly; air dry.

---

## DIGITAL THERMOMETER:

---

For complete instructions for use and care, please see separate instruction manual enclosed.

---

## FINGERBRUSH & CASE:

---

This unique item fits compactly into its own sanitary case. Whether a parent is massaging sore teething gums on the go, or maintaining a good oral brushing hygiene after eating out, this wonderful little product fills the need.

- Clean infant finger toothbrush before first use and after each use.
  - Place on parent's fingertip.
  - Gently clean and massage gums for emerging teeth.
- 

## GUM MASSAGER:

---

- Clean gum stimulator before first use and after each use.
  - As your child's first tooth emerges: massage gums to ease tooth emergence and relieve discomfort.
- 

## 6 EMERY BOARDS:

---

- Gently glide across baby's fingernails to smooth them.
- 

## 5 ALCOHOL SWABS:

---

- Use to clean and disinfect.

We at Summer Infant stand behind all of our products. If you are not completely satisfied or have any questions, please contact our Customer Service Department at 1-800-268-6237 or e-mail us at [customerservice@summerinfant.com](mailto:customerservice@summerinfant.com).

### **Summer Infant, Inc.**

6 Blackstone Valley Place, Suite 206  
Lincoln, RI 02865 USA 1-800-263-6237

[www.summerinfant.com](http://www.summerinfant.com)

© 2005 Summer Infant, Inc.

Please retain information for  
future reference.

### **Summer Infant Europe, LTD.**

20 Sparrows Herne  
Bushey, Hertfordshire  
WD23 1FX UK 44 (0)20 8420 4429  
Colors and styles may vary.

**MADE IN CHINA**

6/05