

User Manual

Body Scale

GBS-1269BT



- Thank you very much for selecting the IC Sciences WeighWise™ Body Scale GBS-1269BT.
- Please read the user manual carefully and thoroughly to ensure the safe usage of this product, and keep the manual available for further reference.



Download and install My BP Coach Pro
from iTunes Store or Google play



Distributed by:
Informed Clinical Sciences Corporation
1 Joy Street, Boston, USA MA 02108

WWW.ICSCIENCES.COM

SCIENCES

Table of Contents

Symbol Meaning
Symbol Meaning 2

Safety Information
Your Scale and Its Environment 3
Efficient Use of Your Scale 3

Overview
General Instructions 4
Device Components 4
LCD Display 5

Initial Start-Up
Insert the Batteries 6
Install App and Pair-up 7

Start Measuring
Select Measurement Unit 8
Daily Measurement 8

Data Management
Data Transmission 9

Troubleshooting
Error Prompt 10
When Measuring 11
When Data Transmitting 11







Specifications 12

Maintenance 13

Warranty 13

♥ Symbol Meaning

The warning signs and symbols are essential to ensure your correct and safe use of this product and protect you and others from injury. Review the warning signs and symbols below:

	Symbol for "THE OPERATION GUIDE MUST BE READ"		Symbol for "MANUFACTURER"
	Data transmitting		Symbol for "COMPLIES WITH EU REQUIREMENTS"
SN	Symbol for "SERIAL NUMBER"		Symbol for "ENVIRONMENT PROTECTION – Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice"
	Symbol for "DIRECT CURRENT"		

♥ Your Scale and Its Environment

To ensure your safety as well as the service life of your scale, please avoid using the scale under the following circumstances:

- Slippery floor such as tile floor
- Weighing immediately after bath or with wet hands
- Near a cell phone or microwave oven

Avoid storage in the following locations:

- Near water
- Where the device may be exposed to extreme temperatures, humidity, moisture, direct sunlight, dust, or salt air
- Where there is risk of shock or dropping
- Where you store chemicals or corrosive gases
- Within reach of infants or children

♥ Efficient Use of Your Scale


To ensure accuracy, please follow the instructions below when you weigh.

- Place the scale on a flat, hard surface. Soft surfaces, such as carpet, will affect the performance of the scale.
- Weigh at least two hours after waking or eating.
- Avoid weighing immediately after strenuous exercise, sauna or bath, drinking, and eating.
- Always weigh at the same time and on the same scale located on the same flat, hard surface.

♥ General Instructions

IC Sciences WeighWise Body Scale GBS-1269BT applies SENSE ON technology. It's not necessary to turn on or turn off the scale, just stand on the scale for measurement.


♥ Device Components






♥ List

- 1. Body Scale GBS-1269BT
- 2. Three AAA-size Batteries (1.5V each)
- 3. User Manual

♥ LCD Display



kg	Kilogram		Data transmitting
lb	Pound		Data pending to transmit to wireless device
	Low Battery		

♥ Insert the Batteries

- Open the battery door in the back of the scale.
- Insert the batteries (3 x 1.5V AAA) into the battery compartment according to the polarity indications marked inside the compartment.


* The digits "8888" will be shown on the LCD.



- Close the battery door and wait until the digits "0.0lb" are shown on the LCD.



CAUTION

- When the symbol  appears, the device will power off in four seconds. Replace batteries with a new set. Please replace all three batteries at the same time. Do NOT mix the old batteries with the new one.
- Worn batteries are hazardous waste. Do NOT dispose of them in household garbage. Please refer to local ordinances and recycling instructions in your area regarding disposal of the worn batteries.

♥ Download and install My BP Coach app from iTunes app store or Google Play


♥ Install App and Pair-Up


With advanced Bluetooth 4.0 technology, your device is capable to receive your personal health data.

Simply install the My BP Coach app and pair up your scale with your mobile device.

Note: It is only available for the device which supports Bluetooth 4.0 technology.

- Search and install the app in your smartphone.
- Turn on Bluetooth and App. Make sure both are ON when pair-up is proceeding.
- Press and hold "UNIT" button in the back of the scale to start pair-up.


Symbol  and symbol  will be shown on the LCD alternatively, indicating pair-up is proceeding.



If successful, symbol  will be shown on the LCD.




If the device fails to connect, symbol "E1" will be shown on the LCD.



♥ Select Measurement Unit

With batteries correctly installed, press “UNIT” button in the back of the scale to select measurement unit. The default measurement unit is “lb”. You may press “UNIT” button to choose between pound and kilogram.



♥ Daily Measurement

• With original SENSE ON patent technology, the scale will automatically switch on as you step on the platform barefooted.




• Stand still and maintain full contact with the scale until the LCD displays the weight and stops blinking.



• When your scale is successfully paired with your smartphone and the Bluetooth is ON, it will process data transmission automatically.


(Please refer to [Data Transmission](#) for more details.)






♥ Data Transmission

• With the scale successfully paired up with your smartphone, the measurement data will be automatically transmitted to your mobile via Bluetooth.







• The symbol  will disappear after successful data transmission, and you may check your personal health data stored in your smartphone.

• If the data transmission fails, the symbol  will remain. The pending measurement data will be temporarily kept in the scale and transmitted to your smartphone when next measurement is complete.

	Data transmitting
	Data ready to transmit to the App: -If successful, the symbol disappears; -If the scale fails to transmit the data, the symbol remains.

Note: Whether the data transmission is successful or not, the device will shut down automatically after 5 seconds.

♥ Error Prompt

Error	Description	Solution
	Overload.	Stop using this scale for measurement.
	Low Battery. The device will power off in four seconds.	Replace all three batteries at the same time. Please purchase the authorized batteries for replacement.
	Failure of pairing up your scale with your smartphone. Or failure of activation.	Please check below items: -Bluetooth is ON. -App is ON. -Both devices are within the transmission distance of Bluetooth.
	Data doesn't upload. (The notebook icon will be displayed.)	Upload the data.

♥ When Measuring ...

Problem	Root Cause	Solution
No display on LCD when the device powers on.	Batteries not yet installed.	Install the batteries. (Please refer to Insert the Batteries)
	Worn batteries.	Replace all three batteries at the same time. Please purchase the authorized batteries for replacement.
The device powers off automatically.	Low battery.	Replace all three batteries at the same time. Please purchase the authorized batteries for replacement.

♥ When Transmitting Data ...

Problem	Root Cause	Solution
Data transmission failed.	Bluetooth is OFF.	Turn ON the Bluetooth via the settings.
	App is OFF.	Press the icon to turn ON your app.
	Out of range of Bluetooth transmission.	Place your smartphone closer to the scale.

♥ Specifications

Product Name	Body Scale (GBS-1269BT)
Dimension	Scale: 11 x 9.6 x 0.9 in (280x245x23mm) (Approximately)
Net Weight	Approximately 3.2 lb (1.46kg) (Excluding the dry cells)
Display	Blue LCD with White Backlight V.A.: 81x41.6mm
Measurement Unit	Kilogram / Pound
Measurement Range	2.5kg to 150kg / 5.6lb to 330lb
Division	0.1kg / 0.2lb
Accuracy	2.5-100kg: ±0.2kg; 100-150kg: ±0.3kg; 5.5-220lb: ±0.4lb; 220-330lb: ±0.6lb;
Working Environment	Temperature: 5°C to 40 °C Humidity: 15% RH to 93%RH 41° F to 104° F
Storage Environment	Temperature: -20 °C to 60 °C Humidity: ≤93% RH -4° F to 140° F
Power Source	4.5V (Three AAA-size Batteries)
Auto-ON	SENSE ON technology
Auto-OFF	About 15 seconds after the weight data is locked
Accessories	1. Three AAA-size Batteries 2. User Manual
Mode of Operation	Continuous Operation

About the Accuracy of This Product

This product passes strict inspection before delivery and therefore its accuracy is guaranteed by the manufacturer. Please refer to the above table for the descriptions on accuracy.

WARNING: No modification of this equipment is allowed.

♥ Maintenance

Follow these tips to maintain your scale:

- Use a dry soft cloth to wipe the dust.
- Use a wet soft cloth, dipped into water and wrung out, to wipe dirt. Then use a dry soft cloth to dry up the device.
- DON'T wash the device with water or immerse it in water.
- DON'T use propellant, abrasive or other chemicals to wipe dirt to avoid discoloration or malfunction.
- DON'T disassemble this device. If you have any problems, please contact Transtek. (Please refer to Warranty for contact information)

♥ Warranty

- Informed Clinical Sciences, Corp warrants its products free of defects in materials and workmanship in normal use for a period of TWO years from the date of retail purchase.
- This warranty does NOT cover damages caused by misuse or abuse, including but not limited to:
 - Failure caused by unauthorized repairs or modifications;
 - Damage caused by shock or drop during transportation;
 - Failure caused by improper operation inconsistent with the instructions stated in this user manual;
 - Malfunction or damage from failure to provide the recommended maintenance;
 - Damage caused by improper use of power supply.
- Should this device require maintenance (or replacement at our option) under warranty, please deliver the original package to Informed Clinical Sciences Corporation prepaid. Please return the store receipt (with the retail purchase date) and a note with reasons to return on it as well.

Distributed by: Informed Clinical Sciences Corporation
1 Joy Street
Boston, MA USA 02108

Website: <http://www.icsciences.com>