

Infra-Heat Lumbar Cushion & Extender



User's Manual



Infrared Heating Element



On/Off Switch

Adapter Input



*AC (house) and DC (car)
adapters included*

Infra-Heat Lumbar Cushion & Extender by Wagan Tech

Thank you for purchasing the Infra-Heat Lumbar Cushion & Extender by Wagan Tech. With normal care and proper treatment it will provide years of reliable service. Please read all operating instructions carefully before use.

WARNING—READ ALL INSTRUCTIONS BEFORE USING

- Always unplug the appliance from the electrical outlet or vehicle immediately after using and before cleaning.
- **DO NOT** reach for an appliance that has fallen into water. Unplug it immediately.
- **DO NOT** use while bathing or in the shower.
- **DO NOT** place or store appliance where it can fall or be pulled into a tub or sink.
- **NEVER** use pins or other metallic fasteners with this appliance.
- Keep dry—**DO NOT** operate in a wet or moist environment.
- An appliance should never be left unattended when plugged in. Unplug from outlet when not in use.
- Close supervision is necessary when this appliance used by, on or near children or disabled persons.
- Use this appliance only for its intended use as described in this manual.
- **NEVER** operate this appliance if it has a damaged cord or plug.
- Keep cord away from heated surfaces.
- **NEVER** use while sleeping or fall asleep while using the product.
- **NEVER** drop or insert any object into any opening.
- **DO NOT** operate where aerosol (spray) products are being used or where oxygen is being administered.

- **DO NOT** operate under a blanket or pillow. Excessive heating can occur and cause fire, electrocution or injury to persons.
- This appliance is designed for personal, non-professional use only.
- **DO NOT** use on an infant, child, disabled person, sleeping or unconscious person.
- **DO NOT** use on insensitive skin or on a person with poor blood circulation.
- **DO NOT** crush.
- If you have any concerns regarding your health, consult your doctor before using this product.
- Individuals with pacemakers should consult a doctor before use.
- Never leave the appliance unattended, especially if children are present.
- **DO NOT** use this product for more than 20 minutes at a time.
- Extensive use could lead to the product's excessive heating and shorten the product's life.
- If excessive heating should occur, discontinue use and allow the unit to cool before operating.
- **NEVER** use this product directly on swollen or inflamed areas or skin lesions.
- **DO NOT** use this product as a substitute for medical attention.
- This product should **NEVER** be used by any individual suffering from any physical ailment that would limit the user's capacity to operate the controls or who has sensory deficiencies in the lower half of their body.

Infra-Heat Lumbar Cushion & Extender by Wagan Tech

OPERATION

For Vehicle Use:

- Remove the supplied 12V adapter from the package.
- Plug in the output side of the 12V cable into the Infra-Heat Lumbar Cushion.
- Plug the 12V tip into the 12V socket in your vehicle.
NOTE: Be sure to unplug the unit when the vehicle is off to avoid possibly depleting your vehicle's battery.

For Home Use:

- Remove the supplied 110V AC adapter from the package.
- Plug in the output side of the 12V cable into the Infra-Heat Lumbar Cushion.
- Plug the AC adapter into your wall outlet.

To turn on, press the button to the "I" position.

To turn off, press the button to the "O" position

USES

Lumbar Extender:

Place the cushion with the flat side of the unit on the floor. Place your body on top of the cushion with your lower back directly on the lumbar cushion. Relax in this position for 20 minutes to allow the lumbar cushion to penetrate and rejuvenate your muscles.



Foot Relief:

Place the cushion with the flat side of the unit on the floor. While seated in a chair, gently place your feet on top of the Infra-heat Lumbar Cushion to enjoy heat and relaxation after a hard day of work.



Infra-Heat Lumbar Cushion & Extender by Wagan Tech

Lower Back Support:

Place the flat side of Infra-Heat Lumbar Cushion onto your chair/couch. Seat yourself on the chair/couch and adjust the angle of the cushion to achieve the most comfortable position for you.



MAINTENANCE

To Clean:

- Before cleaning, remove the power cord from the outlet and allow unit to cool.
- Wipe with a soft cloth. **DO NOT** use abrasive cleaners.
- Never submerge the unit in any liquid. Keep away from all solvents and harsh detergents.

Storage:

Place massager in a safe, dry and cool place. Avoid contact with sharp edges or pointed objects which might cut or puncture the surface. To avoid breakage, **DO NOT** wrap the power cord around the unit.

WAGAN Corp. Limited Warranty

The WAGAN Corp. warranty is limited to products sold only in the United States.

All Wagan Tech products must be registered within (30) days of purchase to activate its warranty. To register your product, please visit <http://tinyurl.com/waganwarranty>. Be sure to keep the original receipt as it will be required when returning a product under the warranty.

Warranty Duration: This product is warranted to the original purchaser for a period of one (1) year from the original purchase date, to be free of defects in material and workmanship. WAGAN Corporation disclaims any liability for consequential damages. In no event will WAGAN Corporation be responsible for any amount of damages beyond the amount paid for the product at retail.

Warranty Performance: During the above one (1) year warranty period, a product with a defect will be replaced with a comparable model when the product is returned to WAGAN Corporation with an original store receipt. The replacement product will be in warranty for the balance of the original one (1) year warranty period.

To return a defective item, please contact WAGAN Corporation at (800) 231-5806 to obtain a Returned Merchandise Authorization number (RMA#), and return instructions. Each item returned will require a separate RMA#. After you have received the RMA# and the return instructions from WAGAN Corporation, please follow the instructions and send the item with PREPAID SHIPPING, along with all of the required documentation, a complete explanation of the problem, your name, address and daytime phone number. WAGAN Corporation will, at its option, replace or repair the defective part.

A Returned Merchandise Authorization number (RMA#) is REQUIRED when sending in any defective item. WAGAN Corporation is not responsible for any item(s) returned without an official Returned Merchandise Authorization number. The item(s) must be returned with prepaid shipping. WAGAN Corporation is not responsible for any shipping charges incurred in returning the item(s) back to the company for repair or replacement. This warranty is void if the product has been damaged by accident, in shipment, unreasonable use, misuse, neglect, improper service, commercial use, repairs by unauthorized personnel or other causes not arising out of defects in materials or workmanship. This warranty does not extend to any units which have been used in violation of written instructions furnished.

Warranty Disclaimers: This warranty is in lieu of all warranties expressed or implied and no representative or person is authorized to assume any other liability in connection with the sale of our products. There shall be no claims for defects or failure of performance or product failure under any theory of tort, contract or commercial law including, but not limited to negligence, gross negligence, strict liability, breach of warranty and breach of contract.

©2008 WAGAN Corp. REV2012



www.wagan.com

customerservice@wagan.com

US & Canada Toll Free: 1.800.231.5806

Tel: +1.510.471.9221

31088 San Clemente Street

Hayward, CA 94544

U.S.A.