

trulyclear

LIGHT THERAPY

Welcome to the **trulyclear 4 in 1 Skin Renewal System**
Ultrasound/Phototherapy/Ionotherapy

for your skin



User Manual

Congratulations on your purchase of a new **trulyclear 5 in 1 Skin Renewal System**

Please read this manual carefully before using your new

trulyclear 5 in 1 Skin Renewal System

Recommended for home and personal use only

Not for treatment of medical conditions

If you have medical conditions please see a medical professional

specifications

Weight: 150 g. (electrical unit) 165 g. (rechargeable unit)
Size: Height 2 inches/ Length 7.3 inches / Width 2 inches
Temperature : 40 - 100 degrees Farenheit

Power Supply:

Input-Power Cord Unit: AC100-240V 50/60HZ; Output: DC15V, 0.4A
Input-Rechargeable Unit: AC100-240V 50/60HZ; Output: DC5.1V 1.0A
Li-ion battery volume: 1150 mAh
Charging Time: We recommend approximately three (3) hours for rechargeable units
Electrical power supply classification: Class II-double insulation

Product Information

Model/Name 5 in 1 Skin Care Renewal System
Manufacturer: accelebrite, subsidiary of milliondollarsmile

LED Information

Red 625± 3nm
Green 525± 3nm
Blue 465± 3nm
Yellow: 592± 3nm

Warranty Covers Unit for One (1) Year

Congratulations on the purchase of your new trulyclear 5 in 1 skin renewal system. Please read the information included in this user manual for complete instructions on getting the most from your system. Keep the manual handy for future reference.

Thank you for ordering our new, revolutionary trulyclear 5 in 1 skin renewal system, an advanced Photon Galvanic Rejuvenation System that includes a skin renewal device designed for home use that is safe, effective and easy to use. Prepare to experience real gain with no pain as you employ our revolutionary system.

Our trulyclear 5 in 1 skin renewal system will deliver phenomenal results with minimum effort on your part. The only scientifically proven home device that actually produces dramatic results, the trulyclear 5 in 1 skin renewal system stimulates your body's own healing powers to heal damaged skin, reduce the appearance of fine lines and wrinkles, firm, smooth and restore a youthful radiance.

Powered by cutting edge, patented BIOWAVE technology, the trulyclear 5 in 1 skin renewal system focuses four low power light beams deep into the dermis, much like a laser lens would. The BIOWAVE stimulates dermal activity, tightening collagen fibers and spurring new collagen production. Skin becomes firmer and more supple, lines and wrinkles are reduced and the effects of the aging process are virtually reversed, for every type and color of skin.

An effective alternative to surgery, botox and chemical peels
A proven, relaxing treatment in the privacy of your own home
Scientific, clinical proof of a high level of success for long-lasting effects

Find a new you.....today!!

important safety notice

- Do not open or disassemble the **trulyclear 5 in 1 skin renewal system** yourself.
- Keep the **trulyclear 5 in 1 skin renewal system** out of reach of children.
- Do not expose **trulyclear 5 in 1 skin renewal system** to open flame
- Do not use the **trulyclear 5 in 1 skin renewal system** in or near water and store in a dry location.
- Dropping the **trulyclear 5 in 1 skin renewal system** device may cause irreparable damage.
- Unplug docking base when not in use

Warning: If you are currently under medical care for a skin condition, consult your doctor before using the trulyclear 5 in 1 skin renewal system. The trulyclear 5 in 1 skin renewal system is designed for personal use only by a single user. Do not share with others.

trulyclear 5 in 1 skin renewal system: Core Technology and Recommendations for Use
trulyclear 5 in 1 skin renewal system effectively focuses BLOWAVE heat deep beneath the skin to support and harness the body's own healing powers.

trulyclear 5 in 1 skin renewal system increases dermal activity and supports healing for a noticeable improvement in skin tone, firming, rejuvenation and the appearance of fine lines and wrinkles. Skin becomes smoother and healthier in a very short time.

- For optimum effect, it is recommended to use maximum effect, trulyclear 5 in 1 skin renewal system three (3) times per week for six (6) to eight (8) weeks.
- For long-lasting effects, repeat the trulyclear 5 in 1 skin renewal system regimen twice per week.
- Safe and effective for all skin types and skin color

Note: It is important for users to read this user manual carefully before use of the trulyclear 5 in 1 skin renewal system.

trulyclear 5 in 1 skin renewal system

List of Parts

Figure 1: trulyclear 5 in 1 skin renewal system device

Figure 2: Adaptor (rechargeable 2-1) (AC Power 2-2)

Figure 3: Manual and Warranty Card

Figure 4: Eye Masks

Figure 5: Gift box

Figure 6: Cotton pad

Figure 7: Cotton holder

High quality stainless steel design of 3M HZ Ultrasonic ION (Galvanic) microwave massage

Dynamic lights: Red LED: 625+/-3nm

Blue LED: 465+/-3 nm

Green LED: 525+/-3nm

On/Off Button

function LCD display

trulyclear 5 in 1 skin renewal system works to resolve your skin issues:

The trulyclear 5 in 1 skin renewal system accomplishes a wide array of functions, including:

- Heal and reduce blemishes and blackheads
- Reduces swelling, edema and puffiness around the eyes
- Helps to heal and clear acne and other skin conditions
- Reduces the appearance of fine lines, wrinkles, and scarring
- renews skin, firms, smoothes, increases skin's capacity to absorb nutrients and moisture

Functions of ultrasonic therapy

The ultrasonic function of the trulyclear 5 in 1 skin renewal system massages skin and increases cell efficiency, stimulating the healing power of the cytomembrane, increasing the permeability of the skin while raising the overall temperature of deep layers of skin to promote absorption of nutrients, accelerate metabolism and reinforce the regenerative processes in tissue. Besides its regenerative abilities, ultrasonic also works to decompose waste in the surrounding fat as well as aiding in the repair of lymph and capillary vessels.

Galvanic (ION⁺ ION⁻)

The Galvanic function of the trulyclear 5 in 1 skin renewal system is effective in establishing a healing electric field between the skin surface and the germinal layer of the dermis, removing dirt from pores and increasing the healing of skin conditions such as acne and rosacea, while moisturizing, firming and reducing the appearance of fine lines and wrinkles. This function is also highly effective in replenishing deep skin healing processes, increasing absorption of essential nutrients and moisture at a level much higher than oral treatments can achieve.

Additionally, this positive electrode therapy is known to serve in increasing blood flow to the skin, strengthening skin tissue, shrinking pores and reducing redness. Concurrently, negative electrode therapy generates alkali action, stimulates nerves, promotes blood flow, firms, moisturizes and reduces the aging process in skin.

trulyclear 5 in 1 skin renewal system LED Light Functions for Skin Rejuvenation

Bio Red LED Light: Promotes the reduction of blemishes, softens skin and reduces the appearance of fine lines and wrinkles, firms and increases the production of collagen

Bio Blue LED Light: Effective in diminishing inflammation, eliminating bacteria that causes acne and healing acne and other scarring

Bio Gree LED Light: Neutralizes, restores balance and calms reddened, irritated skin, easing the effects of skin stress and trauma, while reducing swelling

trulyclear 5 in 1 skin renewal system's Soothing Microwave Massage

Applying a safe and pleasant vibrating microwave massage at a frequency of 12,000 rpm to the deep dermis, trulyclear 5 in 1 skin renewal system effectively cleanses the skin, softening the residue of cosmetics in pores to facilitate cleaning, breaks up oils, reduces the incidence of blackheads and clogged pores and helps to facilitate the discharge of skin cleaning oils. trulyclear 5 in 1 skin renewal system's microwave massage is a relaxing, rejuvenating deep skin massage that reduces swelling, renews and refreshes, erases fine lines and wrinkles, replenishes collagen, firms, smoothes and heals redness, acne, rosacea and chronic skin conditions.

instructions for use of your accelebrite 5 in 1 skin renewal system

1. Make sure the trulyclear 5 in 1 skin renewal system unit is plugged in to an electrical socket for electrical power units. Make sure the trulyclear 5 in 1 skin renewal system portable rechargeable unit is fully connected to the charger and charged (recommended at least 3 hours) before using. You may charge your unit by connecting it to the charger device and connecting the charger to an electrical wall socket, or you may connect it to your computer using a USB cord and USB port. An indicator light on the unit is displayed when charging and turns off when unit is fully charged.
2. Remove the protective film from the head of the trulyclear 5 in 1 skin renewal system prior to using.
3. Before using on face, test an area of skin on your wrist for sensitivity, following the directions for use. If swelling, reddening, itching, pain or other unusual reactions occur, discontinue use. If the condition persists, call your doctor.
4. Remove the protective film from the head of the trulyclear 5 in 1 skin renewal system unit.
5. Turn the power on for the trulyclear 5 in 1 skin renewal system using the power indicator button on the device handle. For rechargeable units, make sure the unit is fully charged. For electrical units, make sure the unit's cord is connected to an electrical outlet.
6. Select the mode preferred using the MODE indicator button on the device handle.
7. Adjust the level by using the SKIN indicator button on the device handle. Levels of operation are CONTINUOUS level and PULSATING level. The timer should be set for one (1) to fifteen (15) minutes.
8. Confirm your settings using the CONFIRM indicator button on the device handle. During your treatment, the CONFIRM button may be used to pause the process, during which you can adjust the mode and levels before resuming the treatment.
9. Apply four (4) to five (5) grams of the trulyclear 5 in 1 skin renewal system essence formula to the tip of the renewal system unit. (Is this correct? In the manual they show a graphic of a substance labeled "Essence" being applied to the tip of the probe, but they're calling it a "cosmetic." I noticed it's not listed in the list of parts included so perhaps its just the user's face

cream.)

10. Begin the treatment by placing the ultrasonic probe of the device against your facial skin, with the complete surface of the head in contact with the skin (do not allow the unit head to tip so that only part of it is in contact with skin. Make sure it is fully in contact with skin.)
11. Move the head across the targeted area of your skin and enjoy the benefits of a relaxing, rejuvenating facial using state of the art ultrasound technology for amazing results.

LCD display settings

Battery charge indicator for rechargeable units

Ultrasonic indicator for five (5) levels of 3M Hz ultrasonic therapy

Ion indicator for Ion and skin cleansing function, and nutrition absorption

Relax indicator for micro-vibration massage, including continuous and pulse modes

Photon indicator for Bio-wave pure red light, blue light and green light

Time indicator for time set, from one (1) to fifteen (15) minutes for safe and efficient operation

Using the trulyclear 5 in 1 skin renewal system

For Eye swelling, puffiness and edema (once weekly treatments are recommended)

1. Clean your skin thoroughly, particularly around the eyes, using a clean cloth and gently massaging with mild soap and water, rinse well and pat dry with a soft cloth.
2. On your trulyclear 5 in 1 skin renewal system device, select the ULTRASOUND therapy and choose a weak intensity level. Indicate the red, blue or green light mode and set the time at ten (10) minutes. During the photon therapy, it is recommended to wear the included eye mask to protect the eyes. (Looking directly at the LED light application may cause discomfort, dizziness and possible injury to eyes.)
3. For a relaxing eye massage, select the RELAX therapy indicator. Choose either continuous or intermittent massage and set time for five (5) minutes.

For once a week care of sensitive skin around the eyes

1. Select the ultrasonic function on your trulyclear 5 in 1 skin renewal system and set functions the same as your daily care regimen settings
2. Choose green and blue light therapy, alternating the two during the session. It is recommended to set the timer for ten (10) minutes for each light application.

Results of trulyclear 5 in 1 skin renewal system treatments

- Some skin redness and minor swelling is expected immediately after the treatment and is a natural part of the healing process. Both will disappear within a few moments.
- Following trulyclear 5 in 1 skin renewal system therapy sessions, users will notice results that

include visible improvement in the elasticity and texture of the skin, which will be smoother and more supple.

- The appearance of fine lines and wrinkles will be significantly reduced.
- In a clinical study, proven results showed a 30 percent reduction in fine lines, a 16 percent reduction in wrinkling within three (3) to four (4) weeks. In addition, overall skin texture was greatly improved. (Note - These results differ in some people, and depend on the individual skin characteristics of each.)

Care and Maintenance of the trulyclear 5 in 1 skin renewal system

- Make sure the power is OFF before cleaning
- Gently wipe the trulyclear 5 in 1 skin renewal system unit with a damp cloth after each use, taking care to remove any residue of gel.
- Store your trulyclear 5 in 1 skin renewal system in its protective case in a cool, dry place

troubleshooting

If trulyclear 5 in 1 skin renewal system does not work properly, please try these common maintenance tasks before calling our repair number at (??)

For Electrical Models:

- Make sure the power cord is properly connected at the unit and at the wall socket.
- Make sure the power button is lit and the power to the unit is on.
- Make sure the proper function is selected according to the instructions in this manual.
- If the power cord is damaged, contact us at (??) for a replacement.
- If you have not selected any operation to begin within 30 seconds, the unit will automatically shut down to protect the circuitry.

For Rechargeable Models:

- Make sure the device is completely charged to full capacity.
- Make sure the device is turned on and the LCD display light is on.
- If the display is not lit, the built-in lithium battery may be low. Try charging the unit to full capacity (about 3 hours).

- If the battery indicator light is flashing, it indicates that the built-in lithium battery is low. Charge the unit to full capacity (about 3 hours).
- If you have not selected any operation to begin within 30 seconds, the unit will automatically shut down to protect the circuitry.