



**OMEGA**

*Redefining Health & Luxury*



**OM-510 User Guide**

Thank you for choosing OM-510 massage chair. Please read this user manual carefully before using the massage chair. Also, please pay careful attention to the section regarding Safety Precautions.

Please keep this user manual handy for future reference.

## Table of Contents

---

- +Safety Precautions
- ◆ **Component Names and Functions**
- +Characteristics and Effectiveness
- +Massage Chair Setup
- +Troubleshooting Guide
- +Daily Maintenance
- +Product Specifications

## Safety Precautions

### 1.Environment

- Do NOT use the massage chair in a damp environment, or one that is full of dust as seepage or mechanical failure may occur.
- Do NOT use the massage chair in an environment above 104" F.
- Do NOT place the massage chair near the stove, burner, or in direct sunlight.
- Please place the massage chair on a smooth and flat surface on the ground.
- Do NOT use the massage chair in wet areas such as the bathroom and kitchen.

#### **Note**

- Keep the massage chair in a cool place. Do NOT increase the room temperature above normal room temperature. If the massage chair is moved to an environment with a significant change in temperature, let it sit for at least an hour before operating it.

### 2.Users with any of the conditions, including but not limited to the following, MUST consult their physicians before using the massage chair:

- People with heart diseases,
- People with malignant tumors,
- Women who are pregnant or menstruating,
- People with spinal abnormalities or curvature,
- People who are being treated or have just completed treatment of any injuries or health problems,
- People with osteoporosis, and
- People who require rest.

### 3.Notes to users

- Users should operate the massage chair according to the user manual
- Do NOT use other treatment instruments or electronics concurrently while operating the massage chair as this may result in injuries.
- Children and pets should NOT play on or around the massage chair.
- Do NOT place heavy articles on the armrests, footrest, or chair-back of the massage chair.
- Do NOT use the massage chair when body is wet.

## Safety Precautions

- Do NOT operate the massage chair when your hands are wet.  
Do NOT massage one spot for long periods of time.

### **Attention**

- Keep your body upright during a massage session.
- To lessen the intensity beyond its preset limits, place a towel, blanket, or pillow between the chair-back and yourself.
- Do NOT use the massage chair as a bed.
- The recommended time for each massage session should not exceed fifteen (15) minutes, and five (5) minutes for each spot on the body.

### **4. Stop** and discontinue using the massage chair if...

- The massage chair is wet and cannot be operated. Switch all power on the massage chair OFF, and unplug its power cord from the wall outlet. Contact customer support.  
You feel uncomfortable during the massage.
- The massage chair is not operating. Switch all power on the massage chair OFF, and unplug its power cord from the wall outlet. Contact customer support.
- It has an electrical failure. Switch all power on the massage chair OFF, and unplug its power cord from the wall outlet. Contact customer support.
- There is a severe thunderstorm. Switch all power on the massage chair OFF, and unplug its power cord from the wall outlet.

### **Note:**

- Occasionally, check the fabric between the chair-back upholstery and the chair-back for wear and tear.

### **Warning**

- Do NOT use the massage chair if the fabric between the **chair-back** upholstery and the chair-back is torn or ripped.

## Safety Precautions

---

### **5.Repair** and maintenance

- Do NOT disassemble and disconnect the chair-back from its base.
- Do NOT attempt to repair the chair yourself. Consult customer support for instructions.

### **6.Electrical** plug

- Do NOT touch the electrical plug or wires when your hands are wet.
- Always hold the electrical wire by its plug when attempting to disconnect it from the wall outlet.
- Do NOT pinch the electrical wire under the weight of the massage chair.
- Do NOT use the massage chair if the electrical wires are damaged. Consult customer support.
- Do NOT use the massage chair if the electrical plug is loose.

## Component Names and Functions


**Front of massage chair**



**Back of massage chair**



**Inside of seat cushion and footrest**



- 1.Headrest pillow
- 2.Chair-back upholstery
- 3.Massage rollers and vibrator(internal)
- 4.Armrest(right)
- 5.Armrest(left)
- 6.Seat cushion
- 7.Footrest
- 8.Chair frame
- 9.The side table
- 10.Remote control

- 11.Remote control holder
- 12.Main power switch
- 13.Foot bracket
- 14.Power cord
- 15.Chair-back locking lever
- 16.Chair-back net cover
- 17.Buttocks massage airbags(3 pieces)
- 18.Thighs massage airbags(1 piece)
- 19.Calf massage airbags(6 pieces)
- 20.Foot massage airbags(8 pieces)

## Component Names and Functions

### Remote control




## Component Names and Functions

1. **Main power ON/OFF** - turn the massage chair ON/OFF.
2. **Timer** - switch from massage session of **10, 20, or 30** minutes.
3. **Back massage ON/OFF** - turn just the back massage ON/OFF, while the rest of the chair continues to massage.
4. **Back massage mode** - choose from automatic (Auto) or manual (Hand), overall (All), partial (Recycle), or fixed (Location) massage.
5. **Technique (Tech)** - choose from kneading (KD), finger press (FP), knocking (KK), or kneading and knocking (SYNC) style of massage.
6. **Strength (back)** - choose back massage intensity from low, medium, or high.
7. **Shake** - turn ON/OFF the back vibration massage.
8. **Width** - adjust the width of the rollers from narrow, regular, **or wide**.
9. **Rollers up arrow** - in manual mode, press and hold down to move the massage rollers up.
10. **Rollers down arrow** - in manual mode, press and hold down to move the massage rollers down.
11. **Seat & footrest massage mode** - choose from automatic (Auto) or manual (Hand) seat and footrest massage.
12. **Strength (seat & footrest)** - choose seat & footrest air massage intensities from low, medium, or high.
13. **Seat & footrest massage ON/OFF** - turn just the seat & footrest air massage ON/OFF, while the rest of the chair continues to massage.
14. **Air adjustment** - choose air massage in the Feet and Calf, Thigh, or Buttocks. Note: Waist is not a selectable option.
15. **Footrest down arrow** - press and hold down to lower the footrest.
16. **Footrest up arrow** - press and hold down to raise the footrest.
17. **Chair-back up arrow** - press and hold down to raise the chair-back.
18. **Chair-back down arrow** - press and hold down to lower the chair-back.



## Component Names and Functions



1. A 2-digit display of the number of minutes remaining in the massage session.

2. **Massage** rollers widths display.

3. **Massage** rollers strengths (intensity) display.

4. **Massage** rollers techniques display.


5. **Seat & footrest** air massage display.

6. **Back** massage display.

7. **Chair-back** and footrest movements display.

### Note

- Do NOT stand, sit, or place heavy objects on the remote control.
- Do NOT drop or throw the remote control to the ground. Please place the remote control in its holder when not in use.
- Do NOT aggressively pull or tug on the remote control wire.
- Do NOT operate the remote control with wet hands.
- Do NOT use any sharp objects to operate the remote control.



1. Air pressure massage indicators.

2. **Massage** rollers position indicators.

## Characteristics and effectiveness

---

### Characteristics:

- Experience the wave-like roller massage on the upper body while the lower body relaxes to the mild and continuous air-pressure massage.
- Choose from the automatic program or a series of manual functions, designed by researchers to simulate professional massages and techniques.
- Four (4) massage rollers to provide your upper back with the best massages.
- Eighteen (18) special-designed air pockets to provide your lower body with a soothing massage experience, which is great for relieving lower back and hip pains.
- LCD remote control lights up the screen to display all the features.
- The built-in back vibrator adds another dimension to the massage.
- Standard time of twenty (20) minutes for automatic massage.
- Effortlessly reclines the chair-back to a maximum of 170 degrees.
- Footrest can be raised simultaneously with the chair-back, or independently.
- Detachable chair-back upholstery makes it easy to clean or replace.
- Remote control holder for your remote control so it is always within reach.

### Effectiveness:

Daily thorough massage may achieve the following effects:

- Dissipate fatigue
- Increase blood circulation
- Reduce muscle fatigue
- Relieve muscle pain
- Relieve muscle strain
- Reduce nerve pain

## Massage Chair Setup and Instructions


### 1. How to bring the chair-back to its upright position.

- The chair-back is folded forward during packaging. Please lift the chair-back to its upright position until a "click" sound is heard, indicating the chair-back is locked in place.



### 2. How to place and position the massage chair.

- Make sure there is enough room for the chair-back and footrest to move. Allocate at least 50 inches behind the chair-back and 30 inches in front of the footrest to allow the chair-back and footrest to move up and down.



Note:

To fold the chair-back forward, pull back the locking mechanism beneath the chair-back and push the chair-back forward.

### ⚠ Caution

- Please place a carpet underneath the massage chair to avoid damage to the floorboard.
- Please be sure the immediate front and back areas of the massage chair are free from any barriers.
- Do NOT drop the massage chair during a move as it may damage the floorboard or carpet.
- Please make sure there is at least a distance of three (3) feet between the massage chair and any audio/video electronics to avoid signal interference.

## Massage Chair Setup and Instructions

### 3. Moving the massage chair.

- Make sure the chair-back is in the upright position and locked. Power down the massage chair and unplug its power cord from the wall outlet. Lift the massage chair by the footrest and tilt it on its casters as shown in the figure. The massage chair can now be rolled to a different location.



### **Caution**

- Do NOT attempt to relocate the massage chair while it is occupied.

### 4. Turn on the power.

- Insert the power cord into the wall outlet
- Power on the chair from the right side of the chair-back

### **Caution**

- Please make sure the power cord is not pinched underneath the massage chair.
- Please make sure the power switch on the chair is turned to OFF before plugging in its power cord.
- If the massage rollers are located in the middle of the chair-back, they will return to the top when the massage chair is first turned on.

### 5. How to operate the massage chair.

- Make sure the footrest is lowered.
- The massage rollers should be located near the top part of the chair, the "storage position" .
- If the massage rollers are not at the storage position, turn the massage chair OFF immediately and return the massage rollers to its storage position. Thereafter, sit on the center of the massage chair slowly.

### 6. After usage.

- Turn OFF the main power switch on the massage chair when it is not in use to avoid accidents and damage to the massage chair.

## Massage Chair Setup and Instructions

### 7. **Reclining** the chair-back.

The power chair-back and footrest adjust buttons are located on the Lower front side of the remote control.



### **Reclining/Raising** the chair-back

- Press and hold down the chair-back down/up button to lower or raise the chair-back, respectively

### **Raising/Lowering** the footrest

- Press and hold down the footrest up/down button to raise or lower the footrest, respectively.

#### **Caution**

- The chair-back and footrest can be adjusted and controlled without having to switch the remote control on.
- Keep hands and feet away from the footrest areas while the footrest is raised or lowered.
- Make sure the front and back areas of the massage chair are free from any barriers.

#### **Warning**

- Do NOT sit or stand on the footrest while it is raised to avoid malfunction and damage to the footrest, or injuries to the user.

## Troubleshooting Guide

Situation	What to do
<ul style="list-style-type: none"> <li>The massage rollers stopped functioning</li> </ul>	<ul style="list-style-type: none"> <li>Check the power cord and ensure it is connected properly.</li> <li>Check the power switch of the massage chair and ensure it is switched on.</li> <li>Make sure the massage rollers are in the storage position (at the top).</li> <li>Press the main power button on the remote control and start the massage again.</li> </ul>
<p>The chair-back is moving in the opposite direction of the button I pressed.</p>	<ul style="list-style-type: none"> <li>Make sure you do not have both up/down buttons pressed.</li> <li>Release the buttons and try again.</li> </ul>
<ul style="list-style-type: none"> <li>The air massage in the footrest stops when I raise or lower the footrest.</li> </ul>	<ul style="list-style-type: none"> <li>This is normal, and a safety feature when the footrest is in motion.</li> <li>Air massage should resume after adjustments of the footrest completes.</li> </ul>
<p>Noises:</p> <ul style="list-style-type: none"> <li>Sound made by the frictions between the cloth and the rollers when back massage is working.</li> <li>Sound made by the back massage when TAPPING function is working.</li> <li>Sound made by the air pump during air massage.</li> <li>Sound made when function changes from one to the next.</li> <li>Sound made by the vibration motor.</li> </ul>	<ul style="list-style-type: none"> <li>Do not worry, as these noises are common to the normal working conditions of the massage chair.</li> </ul>
<ul style="list-style-type: none"> <li>The left and right rollers are not vertically aligned and at the same level</li> </ul>	<ul style="list-style-type: none"> <li>Do not worry, as this is common to the normal working conditions of the massage chair.</li> </ul>
<ul style="list-style-type: none"> <li>The remote control has no control over the massage chair.</li> </ul>	<ul style="list-style-type: none"> <li>Check the power cord and ensure it is connected properly.</li> <li>Check the power switch of the massage chair and ensure it is switched on.</li> <li>Make sure the massage rollers are in the storage position (at the top).</li> <li>Start the massage again.</li> </ul>
<ul style="list-style-type: none"> <li>Massage rollers return to the top after massage ends.</li> </ul>	<ul style="list-style-type: none"> <li>Do not worry, as this is common to the normal working conditions of the massage chair.</li> </ul>

### 1.Storage method

- The massage chair should be kept in a cool and dry place, and away from dust
- If the massage chair **is not going** to be used for a long period of time, please cover it to **avoid** dust.
- Do NOT place the massage chair in a high temperature environment, or in direct sunlight as this may cause the color of the upholstery to fade.


### 2.Maintenance of massage chair

- Turn the massage chair OFF and unplug its power cord from the wall outlet before cleaning it.
- Use only leather-friendly detergent to clean the upholstery of the massage chair. Quickly wipe clean with a dry cloth.
- Use soft dry cloth for daily cleaning.
- Use only a dry cloth to wipe and clean the remote control.

### 3.Cleaning the massage chair-back upholstery

Please use an appropriate detergent, that is leather-friendly, and a soft cloth to clean the chair-back upholstery when necessary.

- Do NOT use an **electric iron** on the surface of the massage chair.
- Do NOT attempt to wash the chair-back upholstery in a washing machine.



### Caution

- Please lift up the chair-back upholstery every so often to inspect for any wear and tear in the chair-back cloth section.
- Please refrain from using the massage chair if there are tears and rips on the chair-back cloth section, or the chair-back upholstery.



## Product specifications

Product name/model number	OM-510
Electric voltage	120V/50-60Hz, 220V/50-60Hz
Consumption	120W-200W
Frequency	50-60Hz
Massage time range	10-30 minutes
Default auto-timer	20 minutes

### Upper-body massage specifications

• Massage path	24 inches		
• Massage width	Narrow	Middle	Wide
— Top rollers	3"	6"	8"
— Bottom rollers	4"	4.5"	<b>6.5"</b>
• Kneading intensities	• 1-3 levels; 16-40 rpm		
• Tapping intensities	• 1-3 levels; 307-610 rpm		
• Vibration intensities	• Automatic		
• Manual massage options	• Upper-body options available • Lower-body options available		

## Product specifications

### Lower-body massage specifications

• Seat massage intensities	• 3 levels; low, medium, high
• Foot massage intensities	• 3 levels: low, medium, high

Chair-back recline angle	120-170 degrees
Massage chair dimension	Upright - 33" x 45" x 47" (W x L x H) Fully reclined 33" x 71" x 33" (W x L x H)
Net weight	159 lbs.
Gross weight	185 lbs.

Note: Actual products may differ in colors, designs, features, and specifications. All product colors, designs, features, and specifications are subject to change without notice.



**OMEGA**

*Redefining Health & Luxury*

OMEGA MASSAGE CHAIR, INC.

Tel: 1-800-659-3650

[www.omegamassage.com](http://www.omegamassage.com)

Email Address: [info@omegamassage.com](mailto:info@omegamassage.com)