HoMEDICS®

HoMEDICS®

The beauty of living well.™

Instruction Manual

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
HoMedics
Service Center Dept. 168
3000 Pontiac Trail
Commerce Township,
MI 48390

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HealthStation™

KnowledgeBook™ Technology measures Body Fat, Body Water, Muscle Mass & Bone Mass

Know exactly what you're made of.™

Getting Started

Battery Installation

1. Uses one 9-Volt alkaline battery (not included).
2. Remove Battery Compartment Cover on unit bottom and attach battery onto the 9-Volt battery snap.
3. The unit needs to be initialized before first use or after battery replacement.
Press lightly on platform. The display will show "0000", then will automatically turn off. Your HealthStation is now ready for use.
4. Replace battery when "LO" is displayed.
5. Remove battery if product is not used for a long time.

Weight Only

The weight only function can be used anytime by just stepping on the HealthStation with or without shoes.

Body Fat, Body Water, Muscle Mass & Bone Mass Functions

To use the Body Composition readings, you must first enter your personal information.

When entering your information, you must move through each step **quickly**. You have 15 seconds to enter data before the unit shuts off to conserve battery power.

Enter Personal Data into Memory

- 1) Press *MODE*
- 2) Set *Height* – Press the ▲ or ▼ button to find your height
When your *Height* is correct, press *MODE*
- 3) Set *Age* – Press the ▲ or ▼ button to find your age
When your *Age* is correct, press *MODE*
- 4) Set *Gender* – Press the ▲ or ▼ button to select MALE or FEMALE
When *Gender* is correct, press *MODE*
- 5) Set *Normal* or *Athlete* – Press the ▲ or ▼ button to select your level of activity
After you choose, press *MODE*
You will hear one long *Beep*
- 6) Set Memory Location – Press the ▲ or ▼ button to select your Memory Location (1-4) After you choose, press *MODE* to store your information into memory
You will hear one long *Beep*
NOTE: If a black box (■) is shown below a Memory Location number before setting a memory location, that Memory Location is assigned to another person. Select a different Memory Location for yourself.

Taking a Body Composition Reading

Before You Begin:

Be sure the scanner is on a *firm, flat surface*

You must have *bare feet* to take this reading

Wipe your feet with a damp cloth, leaving them *slightly moist*, to get the most accurate and consistent reading

NOTE: You *must wait* several hours before taking a reading when:
drinking coffee or alcohol, taking diuretic medications or exercising.

Recall Personal Data from Memory:

- 1) While the power is off, press the ▲ or ▼ button to toggle between Memory Locations.
The ▲ button recalls locations 1 and 3, The ▼ button recalls location 2 and 4.
Once you have selected the correct memory location, press *MODE* to confirm. A "beep" sounds.
- 2) Step on the unit with bare feet. Position your feet as shown.
- 3) Continue to stand still while the unit measures your weight.

The scanner will take your Body Composition reading.
You will hear one long *Beep*

Your weight, Body Fat %, Body Water %, Muscle Mass and Bone Mass are displayed on the screen.

The Status Bar at the bottom of the screen will show your Body Fat and Body Water Health Status. (See Reference Charts on pages 5 - 6)

To Clear Memory:

1. While power is off press either the ▲ or ▼ button to toggle between memory locations.
2. If the selected memory has a (■) under the memory number, press both the ▲ and ▼ buttons at the same time to clear the content of the selected memory location.

Why Should I Monitor My Body Composition?

Typically, total body weight is used to determine if a person is at a healthy weight. A change in weight alone does not tell you whether you have a gain/loss of muscle or fat. In weight management, you want to maintain muscle mass while body fat is lost.

How Does This Body Composition Scanner Work?

This unit uses BIA (Bioelectrical Impedance Analysis) to measure Body Fat, Body Water, Muscle Mass and Bone Mass. BIA sends a low intensity signal through the body. The signal is very low and causes no bodily harm. The signal flows with a different degree of difficulty depending upon the individual's body fat level. This difficulty is called electrical impedance.

Who Can Use This Body Composition Scanner?

This HealthStation is designed for individuals aged 10 – 99 years.

WARNINGS: DO NOT use this product if you have a pacemaker or other internal medical device. This product is not intended for use by children under 10 years of age or pregnant women. The measurements provided are not substitutes for physician assessments of bone mineral density or osteoporosis. Always consult your physician to determine what is right for you.

When Should I Use This Body Composition Scanner?

For maximum accuracy and repeatability, it is recommended that the unit's functions be used at approximately the same time of day under the same circumstances, such as first thing in the morning before you step in the shower. It is also a good practice to avoid swings in the body's hydration level prior to taking the measurement. By establishing the pattern of your body's fluctuation over a period of one month, you can determine your average level of body fat, body water, muscle mass and bone mass. Women will have more fluctuation in water content than men and need to monitor their reading over a longer time period to establish a personal average. **You must determine this as a reference point before measuring progress in a health program.**

You must wait several hours before taking a reading when: drinking coffee or alcohol, taking diuretic medications or exercising. These affect your level of hydration and the accuracy of this scanner.

This instrument is a personal monitor and should be used in a consistent manner for the most accurate reading. While readings may vary between this scanner and other methods of measurement (hydro-densitometry or hand held calipers), your changes in body fat, body water, muscle mass and bone mass shown by this instrument will be reflected accurately.

What Is Body Fat Percentage?

Every body carries two types of body fat: essential fat and storage fat.

Essential fat is needed for normal, everyday physiological functions within the human body. Storage fat is found just below the skin's surface. Though some storage fat is necessary for normal body functions, excessive amounts can be detrimental to your health and overall quality of life. Use this chart as a guide:

Note: Ideal body fat content is not the same for all people. Age, sex and heredity are variables in this measurement. Always consult your physician to determine what is most ideal for you.

BODY FAT % RANGES FOR MEN †

Rating	Age				
	20-29	30-39	40-49	50-59	60+
Low	<13	<14	<16	<17	<18
Normal	14-20	15-21	17-23	18-24	19-25
Moderate	21-23	22-24	24-26	25-27	26-28
High	>23	>24	>26	>27	>28

BODY FAT % RANGES FOR WOMEN †

Rating	Age				
	20-29	30-39	40-49	50-59	60+
Low	<19	<20	<21	<22	<23
Normal	20-28	21-29	22-30	23-31	24-32
Moderate	29-31	30-32	31-33	32-33	33-35
High	>31	>32	>33	>34	>35

† Source: University of Illinois Department of Food Science and Human Nutrition. Body Fat Percentage Calculator. www.ag.uiuc.edu/~food-lab/ai/bfc.html

Your goal in weight loss should be to increase lean muscle mass and to decrease body fat for a higher burning metabolism and improved body tone.

Why Is It Important to Monitor Your Total Body Water (TBW%) Percentage?

Monitoring the level of body water can be a useful tool for one’s health maintenance.

Water is an essential component of the body and is one of the body’s overall health indicators. Water makes up approximately 50-70% of the body’s weight and is present proportionally more in lean tissue as compared to fat tissue.

- For example:
- 73% of lean muscle tissue is water
 - 80% of blood is water
 - 25% of fat is water
 - 22% of bone tissue is water

Additionally, water is also a medium for biochemical reactions that regulate body functions.

The estimated TBW% may vary according to your hydration level. For better accuracy, avoid fluctuations in hydration level prior to taking a measurement.

The accuracy of the unit in estimating TBW% will also decrease with individuals suffering from diseases or illnesses that tend to accumulate water in the body. Please remember that the unit represents only a good approximation of your total body water.

The normal TBW% of an individual varies according to body fat percentage and gender. The table below may be used as a guide:

	Rango de % BF	Rango normal de TBW%
Hombres	4 a 14%	70 a 63%
	15 a 21%	63 a 57%
	22 a 24%	57 a 55%
	25% y m+s	55 a 37%
Mujeres	4 a 20%	70 a 58%
	21 a 29%	58 a 52%
	30 a 32%	52 a 49%
	33% y m+s	49 a 37 %

†† Source: Derived from Wang & Deurenberg: "Hydration of fat-free body mass". American Journal Clin Nutr 1999; 69 833-841.

Why Should I Monitor My Skeletal Muscle Mass?

According to the American College of Sports Medicine (ACSM), lean muscle mass may decrease by nearly 50 percent between the ages of 20 and 90. If you don't do anything to replace that loss, you're losing muscle and increasing fat.

It is also important to know your muscle mass during weight reduction. At rest, the body burns approximately 50 additional calories for each pound of muscle gained. Some advantages of gaining muscle mass include:

- Reversing the decline in strength, bone mineral and muscle mass that accompanies age
- Maintenance of flexible joints
- Guide weight reduction when combined with a healthy diet

This body composition scanner tells you the total weight of skeletal muscle in your body. The normal weight of skeletal muscle mass for men is approximately 40% of total body weight. For women, normal muscle mass is approximately 30% of total body weight (International Commission on Radiological Protection, 1975). To find your muscle mass %, follow this formula:

$$\frac{\text{* Muscle Mass (in pounds)}}{\text{Total Body Weight}} = \text{Muscle Mass \%}$$

Use the chart below to see if your muscle mass is low, normal or high.

MUSCLE MASS % RANGES †

		Masa Muscular (en libras)							
		Normal Mujeres				Normal Hombres			
Peso (en libras)	330	66	83	99	116	132	149	165	
	320	64	80	96	112	128	144	160	
	310	62	78	93	109	124	140	155	
	300	60	75	90	105	120	135	150	
	290	58	73	87	102	116	131	145	
	280	56	70	84	98	112	126	140	
	270	54	68	81	95	108	122	135	
	260	52	65	78	91	104	117	130	
	250	50	63	75	88	100	113	125	
	240	48	60	72	84	96	108	120	
	230	46	58	69	81	92	104	115	
	220	44	55	66	77	88	99	110	
	210	42	53	63	74	84	95	105	
	200	40	50	60	70	80	90	100	
	190	38	48	57	67	76	86	95	
	180	36	45	54	63	72	81	90	
	170	34	43	51	60	68	77	85	
	160	32	40	48	56	64	72	80	
	150	30	38	45	53	60	68	75	
	140	28	35	42	49	56	63	70	
	130	26	33	39	46	52	59	65	
	120	24	30	36	42	48	54	60	
	110	22	28	33	39	44	50	55	
	100	20	25	30	35	40	45	50	
		20%	25%	30%	35%	40%	45%	50%	
		(% de Masa Muscular							

† Source: HoMedics, 2006. Figures based on the following formula: Weight of individual's skeletal muscle (in pounds)/total body weight (in pounds).

How to Use Muscle Mass Chart

- 1) Find your weight (approximate) on the left side of the chart.
- 2) Follow that row across until you find your Muscle Mass in pounds, as shown on the scale.
- 3) Lastly, follow that column down to Muscle Mass %.

Remember, 30% is normal for a female and 40% is normal for a male.

* You may also use the formula on page 7 to better approximate your Muscle Mass %.

Why Should I Monitor Bone Mineral Content?

Your bones contain blood vessels, nerve cells and living bone cells. These are all held together by a framework of hard, non-living material containing calcium and phosphorous.

This body composition scanner estimates the weight of the non-living bone mineral content. This is important information for monitoring and maintaining healthy bones through exercise and a calcium-rich diet.

A healthy level of bone mineral content differs according to age and sex. For example, a 26-year-old FEMALE should have approximately 6 pounds of bone mineral content, where a 26-year-old MALE is considered to have healthy bone mineral content between 6 and 8 pounds. (Rico et al. 1993). For children, 10-16 years old, healthy bone mineral content is between 1.76 pounds and 4.85 pounds (Faulkner et al. 1993)

WARNINGS: The measurements provided are not substitutes for physician assessments of bone mineral density or osteoporosis. Always consult your physician to determine what is right for you.

What is the definition of an Athlete?

An athlete is considered a person who does 10 hours or more per week of aerobic activity and has a resting heart rate of 60 beats per minute. These individuals should select Athlete Mode for the most accurate results. The Athlete Mode is available only for adults 18 years of age or older.

Trouble Shooting

You must have bare feet when using body composition measurements. To get the most accurate and consistent reading, wipe your feet with a damp cloth, leaving them slightly moist before stepping on the HealthStation.

Be sure your feet are centered on the electrodes.

The condition of the skin on the bottom of your feet can affect the reading. The natural effects of aging or activity can make this skin hard. Take a reading with clean, slightly moist feet for best accuracy.

1. Display **ERR0**: Initialization error. The user is advised to take the measurement again for an accurate result.
2. Display **ERR1**: Wrong usage or body not stable. Repeat measurement again.
3. Display **ERR2**: Overload Warning. Remove the weight immediately; otherwise permanent damage to the unit will occur.
4. Display **ERR3**: Contact error. Impedance cannot be measured. Make sure that you are standing still and maintain maximum contact between your feet and metal contacts. Clean bottom of feet with a damp cloth, leave slightly moist, and repeat measurement again.

5. Display **ERR4**: Body fat percentage is beyond the technical limit. Clean bottom of feet with a damp cloth, leave slightly moist, and repeat measurement again.
6. Display **ERR5**: Total body water percentage is out of range. Moistening your feet may help to improve the contact.
7. Display **L0**: Replace the battery when "L0" is displayed.

Care and Maintenance

1. Do not disassemble the HealthStation other than replacing the battery; it contains no user serviceable parts. Damage to the unit may occur as a result of improper handling.
2. Remove battery when the unit is not used for a prolonged period of time.
3. Clean the unit after use with a lightly dampened cloth. Do not use solvents or immerse the unit in water.
4. Avoid excessive impact or vibration to the unit, such as dropping it on the floor.


The product is intended for home/consumer use only; it is not intended for professional use in hospitals or medical facilities.

Note: If you are having any problems operating this HealthStation, please call 800-466-3342, between the hours of 8:30 a.m. to 5:00 p.m. EST or contact Customer Service by e-mail: cservice@homedics.com

Specifications

- Uses BIA (Bioelectrical Impedance Analysis) to determine body composition
- High precision STRAIN GAUGE technology to determine weight
- Weight Capacity: 350 Lb, 160 Kg or 25 St
- Weight Graduation: 0.2 Lb, 100 g, 1/4 Lb (Dashes shown before weight reading in Stone mode)
- Body Fat Graduation: 0.1% (In the range from 4%–50%)
- Body Water Graduation: 0.1% (In the range from 37%–70%)
- Muscle Mass Graduation: 0.2Lb
- Bone Mass Graduation: 0.2Lb
- Height Range: 3' 3.4"–7' 2" (100.0–220.0 cm)
- Age Range: 10–99 years
- Measures weight, body fat, total body water, muscle mass and bone mass simultaneously in about 10 seconds
- Athlete or Normal mode selection
- Weight-only feature
- Auto on/auto off function for convenience and conserving power
- Data can be stored for up to 4 individuals
- User selectable units of measure: (Lb/in) pounds, (Kg/cm) kilograms, or (St/in) stone-pounds
- Low-battery indicator
- Operates on one 9-Volt alkaline battery (not included)

LCD Display

HOMEDICS®

HOMEDICS®
The beauty of living well.™

MANUAL DE INSTRUCCIONES

LIMITED LIFETIME WARRANTY

HoMedics sells its products with the intent that they are free of defects in manufacture and workmanship from the date of original purchase, for as long as it is owned by the original retail purchaser, except as noted below. HoMedics warrants that its products will be free of defects in material and workmanship under normal use and service.

To obtain warranty service on your HoMedics product, mail the product and your dated sales receipt (as proof of purchase), postpaid, to the following address: HoMedics Service Center, Dept. 168, 3000 Pontiac Trail, Commerce Township, MI 48390. No COD's will be accepted.

HoMedics does not authorize anyone, including, but not limited to, Retailers, the subsequent consumer purchaser of the product from a Retailer or remote purchasers, to obligate HoMedics in any way beyond the terms set forth herein.

This warranty does not cover damage caused by misuse or abuse; accident; the attachment of any unauthorized accessory; alteration to the product; improper installation; unauthorized repairs or modifications; improper use of electrical/power supply; loss of power; dropped product; malfunction or damage of an operating part from failure to provide manufacturer's recommended maintenance; transportation damage; theft; neglect; vandalism; or environmental conditions; loss of use during the period the product is at a repair facility or otherwise awaiting parts or repair; or any other conditions whatsoever that are beyond the control of HoMedics.

This warranty is effective only if the product is purchased and operated in the United States of America. A product that requires modifications or adoption to enable it to operate in any other country than the country for which it was designed, manufactured, approved and/or authorized, or repair of products damaged by these modifications is not covered under this warranty.

**For service or repair -
do not return this unit
to the retailer.
Contact HoMedics
Consumer Relations at
1-800-HOMEDIC
(1-800-466-3342).**

Mail To:

**HoMedics
Service Center Dept. 168
3000 Pontiac Trail
Commerce Township, MI
48390**

**e-mail:
cservice@homedics.com**

THE WARRANTY PROVIDED HEREIN SHALL BE THE SOLE AND EXCLUSIVE WARRANTY. THERE SHALL BE NO OTHER WARRANTIES EXPRESS OR IMPLIED INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS OR ANY OTHER OBLIGATION ON THE PART OF THE COMPANY WITH RESPECT TO PRODUCTS COVERED BY THIS WARRANTY. HOMEDICS SHALL HAVE NO LIABILITY FOR ANY INCIDENTAL, CONSEQUENTIAL OR SPECIAL DAMAGES. IN NO EVENT SHALL THIS WARRANTY REQUIRE MORE THAN THE REPAIR OR REPLACEMENT OF ANY PART OR PARTS WHICH ARE FOUND TO BE DEFECTIVE WITHIN THE EFFECTIVE PERIOD OF THE WARRANTY. NO REFUNDS WILL BE GIVEN. IF REPLACEMENT PARTS FOR DEFECTIVE MATERIALS ARE NOT AVAILABLE, HOMEDICS RESERVES THE RIGHT TO MAKE PRODUCT SUBSTITUTIONS IN LIEU OF REPAIR OR REPLACEMENT.

This warranty does not extend to the purchase of opened, used, repackaged and/or resealed products, including but not limited to sale of such products on Internet auction sites and/or sales of such products by surplus or bulk resellers. Any and all warranties or guarantees shall immediately cease and terminate as to any products or parts thereof which are altered, or modified, without the prior express and written consent of HoMedics.

This warranty provides you with specific legal rights. You may have additional rights which may vary from state to state. Because of individual state regulations, some of the above limitations and exclusions may not apply to you.

For more information regarding our product line in the USA, please visit: www.homedics.com

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