Precautions

1. Clean exposed parts with a soft, slightly damp cloth. To remove stains, use a mild soap. Never use detergents, excess water, treated cloths, harsh cleaning agents, or sprays. Do not attempt to lubricate or open scale casing as this will void your warranty.
2. Although your scale is designed to be maintenance-free, it contains sensitive electronic parts. Avoid rough treatment. Do not drop, kick, or jump on it. Treat it with care to ensure the best performance.
3. This scale is designed to read weights between 25 and 400 lbs / 12 and 180 kgs / 2 and 28 st 8 lbs. To prevent damage, do not place more than 400 lbs / 180 kgs / 28 st 8 lbs. on it.
4. Do not store scale where you store cleaning chemicals. The vapors from some household products may affect the electronic components of your scale. Do not store scale on its side.
5. This scale is an extremely sensitive weighing device. In order to prevent run down of the battery, do not store anything on the scale.
6. Do not dispose of batteries in fire. Batteries may explode or leak. Remove the batteries if the scale will not be used for a long period of time.

Five (5) Year Limited Warranty

This scale is warranted against defects in materials of workmanship (excluding battery) for five (5) years for the original purchaser from date of retail purchase. It does not cover damages or wear resulting from accident, misuse, abuse, commercial use, or unauthorized adjustment and/or repair. Do not return to retailer. Should this scale require service (or replacement at our option) while under warranty, please pack the item in the original packaging and return it prepaid, along with store receipt showing date of purchase and a note explaining reason for return to:

TAYLOR Precision Products
2220 Entrada Del Sol, Suite A
Las Cruces, New Mexico
88001 USA

There are no express warranties except as listed above. This warranty gives you specific legal rights, and you may have other rights which vary from state to state.

For additional product information, or warranty information in Canada or elsewhere outside the USA, please contact us through www.taylorusa.com.

©2010 Taylor Precision Products and its affiliated companies, all rights reserved. Taylor® and Leading the Way in Accuracy® are registered trademarks of Taylor Precision Products and its affiliated companies. All rights reserved.

Not legal for trade
Made to our exact specifications in China.
www.taylorusa.com

CP7407 8.10

TAYLOR®

7407

Lithium Electronic Scale

Instruction Manual

Thank you for purchasing a Taylor® Lithium Electronic scale. All of our precision instruments for weight measurement are engineered and designed to meet the highest quality standards ...to assure you years of uncompromising accuracy and consistently dependable, convenient performance. Please read this instruction manual carefully before use. Keep these instructions handy for future reference.

General Operation

1. There is a strip under the Lithium battery cover to prevent battery drain in shipping. Please remove the strip before continuing with these instructions. Also, some models have a static cling label on the lens to prevent scratching. Please remove before use.
2. Your scale has been set at the factory to read weight in pounds. You may select kilogram or stone (1 stone = 14 pounds) units of measurement by moving the slide switch on the bottom of the scale.
3. Place scale on flat, hard surface. Carpeted or uneven floors may affect accuracy.
4. Step on the scale (disregarding the reading) and step off. This process insures the scale is set and is now ready for use. Do this whenever the scale is moved. **At all other times step straight onto the scale.**
5. Step onto the scale platform and remain still. The display will show **"8888"** running from left to right, then display the weight.
6. If you step off of the scale, the display will continue to show your weight for 3 seconds, then **"0.0"** then will turn off automatically. If you continue to stand on the scale, the display will show your weight for 4 seconds and then the scale will turn off automatically. **NOTE:** If you've weighed yourself and a **"C"** is displayed after the zero, weigh yourself again for a more accurate reading. Your weight cannot be locked-in if you do not stand still. If the scale displays **"Err"**, you must step-off and weigh yourself again.

Battery Information

Your Taylor® scale has lithium batteries installed to provide extra long energy life for your scale. Under heavy usage (Health Club, Wellness Clinic, etc.) the batteries may eventually need to be replaced.

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