# tempTalk™

## Voice Readout Digital Thermometer OPERATING INSTRUCTIONS

Temperature Range: 89.6°F – 109.2°F 32.0°C – 42.9°C

### 1. To TURN ON the thermometer:

Press the **ON/OFF** button for **1 Second**.

After releasing the ON/OFF button, the LCD screen will display 37.0°C (98.6°F), then “LO°C (LO°F)” and you will hear it say “Please measure.”

### 2. To disable or activate the VOICE READOUT function:

While unit is off press the **ON/OFF** button for **4 Seconds**. After 3 seconds at one setting normal mode resumes.

### 3. To change voice readout LANGUAGE:

While unit is off press the **ON/OFF** button for **8 Seconds**.

Then press the ON/OFF button again to cycle to the preferred language readout to one of the following options:

Code	Language	Code	Language
US	English	ES	Spanish
FR	French	RU	Russian
IT	Italian	DE	German

After 3 seconds at one setting normal mode resumes.

### 4. To change readouts to CELSIUS (C) or FAHRENHEIT (F):

While unit is off press the **ON/OFF** button for **12 Seconds**.

Then press the ON/OFF button to cycle to the preferred mode of measurement. After 3 seconds at one setting normal mode resumes.

### 5. To take an ORAL TEMPERATURE reading:

(3 months and older) Place the probe tip **under the tongue** as near as possible to the area **furthest back**.

When you hear beeping the measurement is complete and the thermometer will readout both, visibly on the display screen, and audibly (unless disabled).

### 6. To take an UNDERARM TEMPERATURE reading:

Dry the underarm with a towel and place the probe tip **high in the armpit** so the thermometer is perpendicular to the body, and position the arm across the chest so the tip is fully covered by bare skin.

When you hear beeping the measurement is complete and the thermometer will readout both, visibly on the display

screen, and audibly (unless disabled).

### 7. To take a RECTAL TEMPERATURE reading:

(newborn and older) Coat the probe tip with petroleum jelly. With your free hand, spread baby's buttocks apart. Gently insert the probe tip about half an inch to an inch (1.3 to 2.5 centimeters) into the rectum, or until the silver tip disappears. When you hear beeping the measurement is complete and the thermometer will readout both, visibly on the display screen, and audibly (unless disabled).

### 8. To LISTEN AGAIN to the temperature readout:

Press the **Speaker** symbol button (🔊) and your thermometer will recall and announce the last measurement taken.

### 9. To TURN OFF the thermometer:

Press the **ON/OFF** button for **less than 1 Second**. The unit will also turn off automatically after five minutes from the last use.

### Battery Replacement

1. Replace the batteries in your thermometer when the low battery indicator appears on the LCD display.
2. Remove the battery cover from the back of the unit and remove all three batteries.
3. Load three new batteries (size AG10) so they are lined up with the + sides to the right (nearest the latch indentation) inside and replace the cover.

### Precautions

- The thermometer probe tip should be disinfected with alcohol or washed with soap and warm water; the LCD display screen should only be cleaned with a towel as that part is not waterproof; **Do Not** place upper part of thermometer into water or wipe it with an overly moist towel.

### Tips

- Always place the thermometer tip in the correct position to consistently measure temperature accurately.
- Oral and rectal readings are more precise than underarm measurements.
- Be sure to replace the batteries after 18 months of use (or after approximately 350 uses) for continued accurate readings.

## Mobi Technologies, Inc.

5913 Blackwelder Street  
Culver City, CA 90232

[www.getmobi.com](http://www.getmobi.com)

Mobi and TempTalk are trademarks of Mobi Technologies, Inc.

© 2011 Mobi Technologies, Inc. All rights reserved.

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