


OMRON




Digital thermometer Model Eco Temp II (Water-resistant)

- Digitales Fieberthermometer
- Thermomètre digital
- Termometro digitale
- Termómetro digital
- Termómetro digital
- Digitale termometer
- Digitaalinen kuumemittari
- Digital termometer
- Digital feberthermometeret
- Digitale koortsthermometer
- Ψηφιακό θερμόμετρο
- Dijital klinik termometresi
- Цифровой термометр



A Good Sense of Health

IM-MC-241W-E-03-08/08
1664195-7C



Digital Thermometer English

3. Use of the temperature

- A. Sensor Tip
B. Display Window
C. ON/OFF Button
D. Storage Case

Introduction

Thank you for purchasing an OMRON Eco Temp II Digital Thermometer. The thermometer you purchased offers a safe, accurate and quick temperature reading. You can measure your temperature either in the anus (rectal) or in the mouth (orally) or in the armpit (axillary).

Thoroughly read this instruction manual before using your OMRON Eco Temp II Digital Thermometer. Contact your doctor if you have any questions regarding specific temperatures.

- (37.0°C/98.6°F) is not displayed during the measurement preparation period as a test function. (Fig. 9)
- Note: When a temperature other than 37.0°C/98.6°F is displayed, please refer to the "Troubleshooting" section for details to correct the problem.

2. Applying the thermometer to the measurement site

2. Apply the thermometer to the measurement site as described in "2. Correct measurement".
- Note: Either temperature or "L" is displayed while taking a measurement. (Fig. 10)

3. When the reading is ready, the thermometer will emit a (beep) sound three times after approx. 60 to 120 seconds (depending on place of measurement and environmental conditions). Remove the thermometer from the measurement site and read the result. (Fig. 11)

Important Safety Notes

To assure the correct use of the product basic safety measures should always be followed including the precautions listed below.

⚠ Warning

Please note that the temperature reading can be affected by various physical circumstances. Please consult a doctor if the thermometer measures abnormal temperatures over a prolonged period of time.

- Please keep still during measuring.
- Keep the thermometer and its small parts out of reach of children. Do not allow children to take their temperatures unattended.
- Do not attempt measurements when the thermometer is wet.
- Do not attempt rectal measurements on persons with rectal disorders.

⚠ Caution

Do handle the thermometer with care and do not use the thermometer for anything else than the intended use on human beings. Do not expose the thermometer to shocks.

- Changing and disposal of the battery has to be handled with reference to this manual and has to take place in accordance with the regulations of the country. Do not attempt to disassemble or repair the thermometer. After each use, disinfect the thermometer especially in case the device is used by more than one person. Please refer to "7. Cleaning and Disinfecting".

⚠ Attention

The tip of the thermometer contains Nickel.

8. Specifications

Digital thermometer

- Product type:** Eco Temp II
Product Name: MC-241W
Power Supply: 1.5V DC, 1 LR41 Alkaline-Magnesium Button Battery
Power Consumption: 0.1mW
Battery Life: Approx. 1100 readings (If used for 10 minutes per day).
Sensing Unit: Thermistor
Measurement Method: Actual measurement (non-predictive)
Temperature Display: 4-digit, +/-(°F) display in 0.1 degree increments
Battery Replacement Indicator (Fig. 10): Replace the battery when the battery replacement indicator appears when the thermometer is switched on.

- Measurement Range:** 32.0 to 42.0°C (89.6 to 107.6°F)
Usage Environment Temp and Humidity: Surrounding temperature +10 to +40°C (+50 to +104°F), Relative Humidity 30-85%RH
Storage Environment Temp and Humidity: Surrounding temperature -20 to +60°C (-4 to +140°F), Relative Humidity 10-95%RH


- Weight:** Approx. 11g (with battery installed)
External Dimensions: 18.4mm(w) x 130.0mm(l) x 11.3mm(d)
Package Content: Test Battery (Alkaline-Magnesium Button Battery LR41), Storage Case, Instruction Sheet.

- Notes:**
• This OMRON product is produced under the strict quality system of OMRON Healthcare Co., Ltd., Japan.
• The specifications may be changed without prior notice.

• This unit is water-resistant.
• This thermometer does not contain any serviceable part. OMRON recommends that the accuracy is verified every two years by OMRON's representative.

- Warranty:**
This product is guaranteed by OMRON for 3 years after date of purchase. Any damage caused by improper handling shall not be covered by warranty.

Batteries and packaging are also excluded from the warranty. All other damage claims excluded. A warranty claim must be submitted with the purchase receipt.



= Type B

Read the Instructions in this manual carefully.

This device fulfills the provisions of the EC directive 93/42/EEC (Medical Device Directive) and the European Standard EN12470:2000, Clinical thermometers - Part 3: Performance of compact electrical thermometers (non-predictive and predictive) with maximum device.

The following may result in incorrect measurements:
• Heavy armpit perspiration
• Correct: Wipe off any armpit perspiration before taking the temperature.

- Taking the temperature after being under the blankets for a long time.