

Polarizador JYM[®]



Science that takes care of you.



clinic



pulse

Product features



Pulse Polarizer

Weight:	50 grams approx.
Made of:	stainless steel
Lifetime:	more than 20 years

Wristwatch to be used the whole time, whether wearred on the wrist or fastened to another part of the body.




Clinic Polarizer

Weight:	220 grams approx.
Made of:	stainless steel
Lifetime:	more than 20 years

Cylinder used to stir up liquids and rub the affected zones.

Scientific principle


This technology operates based on generating a force acting through the ions or molecules that come into contact with the force fields emitted by the polarizer.



*F: Force field generated by the polarizer

When stirring up liquids with the Clinic Polarizer, and wearing the Pulse Polarizer on the wrist, the force fields emitted by the Polarizer modify the binding energy within atoms and molecules, obtaining less viscous solutions with higher hydration and diffusion energy, optimizing bodily functions such as: blood circulation, digestion, anti-inflammatory and diffusive processes, thus achieving a satisfactory treatment for the diseases attended.

Polarization effects



The great majority of diseases attended with the polarizer have its conception in at least one of the processes mentioned in the figure above. That is why when optimizing blood circulation, digestion, anti-inflammatory and diffusive processes, real and positive effects are achieved in all the affections treated with this technology.

Treated diseases



Diseases treated only with the: Clinic Polarizer

- Sciatica
- Constipation
- Gastritis
- Hemorrhoids
- Migraine
- Pulmonary Emphysema



Diseases treated with both the: Clinic Polarizer + Pulse Polarizer

- Uric acid (gout)
- Allergies
- Atherosclerosis
- Arthritis
- Diabetes
- Hypertensión
- Chemotherapy damages
- Fybromyalgia
- Thrombosis
- Dysautonomia
- Gall and kidney stones

Use instructions



Clinic Polarizer

The Clinic Polarizer has two ways of usage: stirring up fluids for their further intake and rubbing it up against the aching areas.

Stirring

- Heat by any means the liquid you are drinking (water, coffee, tea, etc.).
- Stir up the previously heated liquid with the Clinic Polarizer as fast as possible for one or two minutes (once the chronic disease stage has been overcome, the stirring processes could be decreased to 30 seconds).
- The stirring direction has no influence in the efficiency of the treatment
- After stirring, leave the Clinic Polarizer inside the cup or glass and continue stirring sporadically while drinking.



Use instructions

Amount of liquids


- In order to achieve satisfactory results when treating a disease, the patient should intake two liters of liquids per day.
- This amount of fluids should be distributed throughout the day.
- Every cup o glass of liquid should be stirred up and not the two liters at once.
- Some people (not everyone) will lose weight by following the treatment. If this is your case and you are not interested in this benefit, please reduce the the temperature and/or amount of liquids intake (note that by reducing the amount and/or temperature of liquid intake, the results of the treatment could be delayed).
- Ingesting more than two liters of hot liquids could occur: weakness and excessive weight loss.
- In time of pain or high blood pressure take three glasses of polarized liquids immediately.



Use instructions

Liquid temperature

- The ingested liquid temperature should be the highest that you stand to take them without burning your mouth or tongue.
- The average temperature recommended is similar to that of the liquid heated in a microwave for 1:30 minutes (with reference to the average power of a microwave oven). However, because each body is different, this temperature may vary from one person to another.
- The higher the temperature of the liquids the higher the efficiency of the treatment and thus lower the time taken to see results.
- If you are suffering from hypotension (low pressure) it is recommended to take warm liquids instead of hot, as hot can bring it down (in people who do not present hypotension liquid temperature must be high).
- Some people experience excessive weight loss, in this case you can reduce both, the amount of liquid ingested per day or the temperature of the same (note that by reducing the amount and/or temperature of liquid intake, the results of the treatment could be delayed).




Use instructions

Rubbing

In some diseases which present pain or swelling, in addition to using the Clinic Polarizer for stirring liquids to drink, it is important to rub the affected areas (as long as possible) with such equipment as specified in the table below:

- Arthritis Rub the aching area
- Sciatica Rub the aching area
- Constipation Rub the belly area
- Gastritis Rub the aching area
- Migraine Rub the neck and the aching area
- Fibromyalgia Rub the aching area
- Thrombosis Rub the aching area
- Diabetes Rub the thighs of both legs, arms and belly

*Note: Rubbing is not recommended while pregnant.



Use instructions

Clinic Polarizer cleansing

Wash the Clinic Polarizer with soap and water (as if it were a fork or spoon) after each use.



Pulse Polarizer

This polarizer complements the treatment for some diseases (see the treatment required for each disease).

Use

The Pulse Polarizer is usually placed on either arm's wrist.

The base of the Pulse Polarizer should be placed making contact with the inside of the wrist, like a watch that has been turned to the other side of the arm (see the photo on the backcover).

Use instructions

In certain diseases where a specific ailment is shown, the Pulse Polarizer can be positioned as follows:

- Hypertension Use on the right wrist
- Swollen leg Attach* to the swollen area
- Swollen gallbladder Attach* to the gallbladder area until inflammation ceases
- Espolón Attach* to the ankle area until the ache disappears
- Swollen kidney Attach* to the kidney zone until inflammation ceases
- Gall and kidney stones Attach* to the gallbladder area until the stones disappear

The usage of the Pulse Polarizer is complementary treatment and does not reduce fluid intake polarized.

* Micropore or bandage can be used to attach the Pulse Polarizer to the affected area.



Use instructions

Pulse Polarizer treatment

The Pulse Polarizer is recommended for a 24/7 use, though it can be unworn when showering or sleeping (except in the treatment of allergy, in this case the Pulse Polarizer must be worn at all times).


When dealing with diseases that involve swelling or aching, the Pulse Polarizer could be worn on any arm's wrist once these discomforts cease.

Pulse Polarizer Cleansing

If the Pulse Polarizer begins to stain (this depends on the skin of each individual), it can be cleaned with creams and liquids for polishing metals.

Important

The Pulse Polarizer is made of stainless steel and leather strap. If you are allergic to any of these materials you may have the following symptoms:

- Redness at the contact site
 - Hives and/or itching
 - Weakness
- 

Use instructions

If you observe this you can use the Pulse Polarizer with a thin cloth between the base of the polarizer and the skin. If prior to purchasing you are aware of an allergic reaction when using these materials, please notify us for a quotation of special a material polarizer.



Clinic Polarizer + Pulse Polarizer

The following observations apply to both the Clinic and Pulse Polarizer.

Length of Treatment

The treatment of polarization comprises 2 phases. The first phase corresponds to the crisis stage of the disease, and has as main purpose the elimination of aches and pains of the disease attended. The second phase starts from ceasing of the aches and pains of the first phase, and continues indefinitely, with the main focus of preventing other diseases.

Length of Treatment


Phase 1: has an average length of 15 to 60 days. During this phase you will notice a decrease in pain, swelling and other discomforts of the disease attended. This time will depend on variables such as:

- Extent of disease progression
- Liquid intake: the higher the amount of liquids the less time taken
- Liquid temperature: the higher the temperature the less time taken
- Rubbing time with the Clinic Polarizer
- Wearing time of the Pulse Polarizer

Phase 2: starts when the disease is controlled and continues indefinitely, preventing diseases such as:

- Circulation: hemophilia, thrombosis, hemorrhoids, among others.
- Digestion: gastritis, colitis and gastric reflux among others.
- Anti-inflammatory processes, Parkinson's and Alzheimer's, among others.

During this phase, if desired, the amount of polarized liquids ingested can be reduced (from 2 liters to 0.5 liter per day). When attending with diseases that causes swelling or pain, after ceasing these discomforts, you can use the Pulse Polarizer on the either arm's wrist.



Preventive polarization

Polarizing healthy people

Polarization is always recommended for healthy people who seek a better performance of the body functions and prevent diseases.

In the prevention stage, high temperature polarized drinking water is not necessary, if the person so chooses, he can polarize warm or even cold liquids.


Polarizing babies

- Severe colics Stir up the milk with the Clinic Polarizer prior to the babies intake.
- Baby chafing Boil water and stir it up with the Clinic Polarizer. Leave the polarizer in the container while the water cools. Once cooled, apply it on the chafing.



Important

- Some diseases may have symptoms similar to others and, therefore, be confused. For example, gastritis has symptoms similar to the lack of enzymes, therefore, if you buy the product for treating gastritis but you actually suffer lack of enzymes, the results obviously are not going to be the expected. For further information, please visit the following address:
www.cienciaquetecuida.com/ingles/tests.html
- Some people (not everyone) get to feel slight dizziness or increased pain for the first days of the treatment. This effect is temporary and disappears after a few days.
- Any wound, diabetic and nondiabetic, washed with polarized water, diabetics and nondiabetics, heals faster (follow the baby chafing procedure of application).
- Polarization can be administered to any person, regardless of age and health status, even during pregnancy.
- Polarization does not oppose with any medical treatment, on the contrary, all drugs efficiency increases when taken with polarized water. If you are taking medicines we recommend to take them with cold polarized water (to facilitate the pills' intake).
- Polarization does not cause any collateral damages.



Important

- In the case of chemotherapy, polarization does not destroy cancer, it prevents or greatly reduces the adverse side effects of drugs.
- DANGER: In case of having a cardiac pacemaker, never use the Pulse Polarizer like a wristwatch in your arms, it should be used in any of your legs. Do not rub the Polarizador Clinic against your heart.
- Never use any of the products near credit cards or electronic devices but the Pulse Polarizer can be used while working on a computer without risk
- Never expose both polarizers to the fire or the direct flame of a stove or any other source of heat like the microwave.
- Both the strap and lugs of the Pulse Polarizer can be replaced by the client.
- Do not leave your current treatment without supervision of your physician
- Consult your physician.

Testimonails

We invite you to see some testimonials from customers who have found in the scientific principles of polarization an efficient alternative to their various ailments. Visit the following website:

www.cienciaquetecuida.com/ingles/index-3.html

Polarization tests

Foam increase

1. Wash your hands twice using a soap bar. Wash your hands a third time with the same bar of soap. This time analyze the amount and quality of foam generated.
2. Take approximately 200 ml. of water from the same place you previously took water for the prior hand washing. Stir up the water with the Clinic Polarizer for about two minutes.
3. Take the same soap bar in your hands and pour the polarized water over them. Rub hands together to make lather and compare the quantity and quality of the foam generated with the untreated and the polarized water.
4. Foam production and quality increase. This is due to the increased hydration and hydrolysis energy in the polarized water.

Tobacco softening

1. Turn on a strong flavor cigarette (Marlboro for example) and taste.
2. Then rub the cigarette against the Clinic Polarizer for about one or two minutes.
3. Try again the flavor of the cigarette and you will notice that has been softened.



Polarization tests


Alcoholic beverages softening

1. Taste a strong alcoholic drink (tequila for example).
2. Stir it up Clinic Polarizer for about one or two minutes.
3. Taste the drink again and you will notice that its aggressiveness has been diminished.

Removal of gas in soft drinks

1. Have a drink (Cola preferably) in a glass and try it.
2. Stir the drink with the Clinic Polarizer for about one or two minutes.
3. Try the polarized soda and you will notice that the CO₂ concentration was reduced.

Water hardness reduction

1. Please test your water hardness.
 2. Stir the water from the same sample with the Clinic Polarizer for about one or two minutes.
 3. Perform a hardness test (with dropper for a more accurate result) and the hardness value will be reduced within a range of 3 to 10%.
- 



Clinic Polarizer



Pulse Polarizer



química de fluidos

Pino 20 Col. Villas del Descanso
Jiutepec, Mor. 62554 MEXICO

+52 (777) 244 32 20

+52 (777) 124 09 89

www.cienciaquetecuida.com