

M-2 Negative-Ion Massage Chair

User's Manual



Contents

Thank you for choosing our product. Please read this operating instructions book thoroughly before operating the massage chair, and after the careful reading, please keep the instructions where you can pick it easily.

CONTENTS	Cautions on Safety	1
	Functions & Specialties	2~3
	Specifications	3
	Names and Functions of Parts	4
	Preparations Before Using	5
	Names & Function of Remote Control Button...	6
	Proper Ways of Using	7~9
	Adjust Back Angle	8
	Adjust Shelves	9
Ways of Storage	9	

Cautions on Safety

- Don't let the baby or child use or play on it.
- If the wire is broken or it doesn't work, please stop using it. Turn off the electricity and pull out the plug at once. Please consult the dealer or go to the maintenance place first.
- If you feel uncomfortable or feel bad at the strong stimulus when you're using, or feel the machine abnormal, please cut off the electricity and stop using it.



- Don't use it when the surface is broken. Please consult the dealer or go to the maintenance place first.
- If you have following problems, please consult your doctors before using:
Accepting treatment at present or feel bad, vertebra abnormal, vertebra crooked, having acute disease, hyperplasia, intervertebral disc protruding, pregnant, having menstruation, having high fever and so on.



- Don't put the hand into the working part of the Chair.
- Before you adjust the angle, make sure there is nothing or nobody behind the machine, and there is enough room behind it.



- Don't use it when bare-backed or wet over.
- Don't use it in the bathroom where there is a lot of moisture. Besides, don't let the remote control drop into the water.
- Please keep away from direct sunlight or heating appliances like heat stove.
- Don't use diluents, volatile oil or bleaching powder will cause cloth change color or damaged.



Functions & Specialties

Inside the back part, there is an exclusive walking massage mechanism with the coverage of the whole back; the mechanism can massage automatically in the 4 mobile scopes from neck to waist; thru the button control, it can also stop at any suitable spot. This mechanism has a strong function; the operation thru the controller can realize the following modes:

1.The kneading and pinching massage

This exclusive mechanism can imitate the specialized massagist to pull, press, knead and pinch the painful and numb muscles on the back with gentle excitation to the muscles to restore their elasticity and strengthen the blood circulation on the massaged spot, and let the body to relax in a limited time.

2.Swaying-massage

Not only does relaxing swaying-massage adjust the position of spinal column and shake each vertebra and ligament, but also rectifies spinal column, stimulates sympathetic nerve, activates the cells and promotes blood-cycle.

3.The finger pressing massage

Inside the mechanism there is a finger pressing wheel that imitates the hand finger pressing and carries on the pressing-relaxing motions repeatedly on the important points on human back with suitable force. The 3 adjustable finger widths help excite the back effectively to allow for smoother blood flow, fully activated human potentials and fortified human immunity.

4.The vibrating massage

The vibrating massage machine at the cushion can produce durative massage vibration wave which forms an efficient excitation to the leg, promote the blood circulation in the capillary vessels and relieve the tiredness quickly.

5.Thigh-massage

Eight massage-heads are equipped at the thigh part of the cushion to give soft but powerful massage, which helps promptly mitigate ache in the thigh.

6.Calf-massage

Kneading-massage: The equipment can do kneading-massage on the calf from ascend to bottom or from bottom to ascend to release the fatigue in the calf. Flexible massage-heads of the kneading-massage system are scattered in the side of the calf. Well-distributed strength of kneading makes it more comfortable.

Functions & Specialties

7. Negative ion healthy oxygen-bar

This equipment is designed to be your negative ion healthy oxygen-bar. It can bring you natural enjoyments as you have experienced in forests, at seashores and cataracts. The negative ion can raise the density of oxygen in the air, so it is helpful to solve the problems such as fatigue because of lack of oxygen in sealed rooms, brain hypoxia and anoxic skin-aging. It makes you energetic every minute.

8. Angle-adjusting function

Two restriction-free pneumatic angle-adjusting instruments are specially designed. Working-angle of the backrest and the foot-shelf can adjust freely according to what your preference. Any massage posture can be satisfied. You can be sitting or even lying to enjoy the massage.

9. Electronic timer

The inserted electronic timer can time automatically after starting up. It can turn off after working a fixed time. This function can avoid the waste of energy and excessive stimulation because of forgetting to turn off.

10. Automatic program include kinds of massage skills, massage time can be adjusted.

Specifications

Rated voltage: 100V-120V AC 220V-240V AC

Rated frequency: 50/60HZ

Power consumption: 130W

Working time: 15minutes

Weight: 48 kg


Volume: high1000 X width700 X depth900 mm

Safety devices: over heat proof devices, electronic timer

IEC 60335-2-32:1993

※ **Note:** Specifications are subject to change without notice.

Names and Functions of Parts



※ **Note:** Specifications are subject to change without notice.

Preparations Before Using




1. Open the pack of the massage-chair and carry the chassis out with 4 people (Don't lift the handrails roughly). Please keep the Instructions and Guarantee safe.


2. Turn over the backrest cautiously. Fasten the orientation sheath at the lower end of the backrest and the pneumatic spring at the lower part of the chair with provided bolts and block-forelocks.

3. Find the power cords under the chair. Connect the plugs of the same color tight.


4. Leave certain room around the massage-chair to avoid colliding other substances when adjusting the angle and turning the chair.



5. Finishing assemblage, insert the power plug into the socket and turn on the power switch at the lower part on the right side of backrest, then you can use the telecontrol panel to select massage functions and enjoy relaxing, healthy and comfortable massage.




Names & Function of Remote Control Button




Proper Ways of Using


I Start the massage-chair

Press the button  on the telecontrol panel of the massage-chair, and it will be ready to work.

II Set the time of massage





After starting the massage-chair, press the time-selecting button  to do time setting. Three alternatives are available; they are 5-minute, 10-minute and 15-minute. From the display you can see the set working time.

III Automatic massage program





When automatic massage is preferred, press the button , and the massage-chair will be ready to massage automatically. Under this condition, the massage equipment will knead the back and the leg automatically according to the computer program complemented inside.

IV Back-massage



1. Select kneading style

Press the button  and the massage-chair begins to knead the back. Usually the implicit route of the massage-chair is the whole body. You can select different kneading scopes using the button , or you can also move the massage equipment to locate the kneading scope through the button  and .


2. Select swaying style

Press the button  and the massage-chair begins to sway the back. Usually the implicit route of the massage-chair is the whole body. You can select different swaying scopes using the button , or you can also move the massage equipment to locate the swaying scope through the button  and .

3. Select finger-pressing style

Press the button , and the massage-chair begins to massage the back with finger-pressing. Usually the implicit route of the massage-chair is the whole body. Different finger-pressing scopes can be selected through the button .

4. Select the speed of massage

When the massage-chair is at work, the rate of kneading or swinging can be adjusted by pressing the button . Three alternative rates are available. From the display, you can see different rates selection from L1 to L3. And different massage rate can be selected as needed.

Proper Ways of Using

5. Choosing **massage areas**


Press **PARTS**, you can choose one of the three massage areas.

When it is chosen, the indicator light relative to it will be on.




the whole course

It will give massage up and down from your neck to your waist.




shoulder part



It will give massage up and down from your neck part to the upper back.



waist part

It will give massage up and down from your lower back to your waist part.



- Press  or , the massage device of backrest will move up or down. On release of the button, the mechanism stops for spot massage. The pressing massage can't spot massage.

V Thigh-massage

Press the button **THIGH**, the instruments at the cushion will take relaxing massage on the thigh. Press this button again to turn off this function.

VI Calf-massage

The equipment imitates the hand to knead and extrude on the calf. And this stimulates the leg blood circulation, relaxes the leg muscles and releases ache and anaesthesia.

- Press the button **CALF** and the massage-chair begins to massage the calf. Press this button again to stop massaging.

- Press the button **MODE** to adjust the kneading style. Three modes are available. From the display, you can see different mode selected, from L1 to L3. And different modes can be selected as needed.

VII Adjust back angle

- If you want to make the chair back down, lean against the back and pull the back adjusting pole, when found a proper site, loosen adjusting of back will be lock.
- When want the chair back up, leave the back and pull the back adjusting pole.
- ※ It will not effect with the massage function when adjust the chair angle.

Proper Ways of Using




VIII Adjust foot-shelves

- Pull the foot-shelves angle adjusting pulling-pole under the left handrail upward and lift feet slightly, then the foot-shelves will lift automatically.
- When you need lower the foot-shelves, put your legs on the foot-shelves slightly and pull the foot-shelves angle adjusting pulling-pole on the left, then the foot-shelves will lower slowly.
- When you need adjust the length of the foot-shelves, please pull the adjusting pole on the right upward, move the foot-shelf forward or backward to convenient position, then release the adjusting pole. And the foot-shelves will be locked there.

IX Finish operation

After operation, please turn off the power switch of the massage-chair, take the plug off the socket and put the telecontrol panel in the hoarding-box under the right handrail of the massage chair.

Ways of storage

	<ul style="list-style-type: none">● Please keep it where it is dry and there's less dust.◎ Damp and dust will shorten the machine's life and may become the causes of moldy or dirty marks.
	<ul style="list-style-type: none">● Please keep away from direct sunlight or heating appliances like heat stove.◎ The surface may change color or lose high quality.
	<ul style="list-style-type: none">● Make sure the power is turned off, and pull out the plug from the outlet when it is dept. Besides don't use wet hands to pull or insert the plug.◎ To avoid electric shock.
	<ul style="list-style-type: none">● Don't use diluents, volatile oil or bleaching powder.◎ will cause cloth change color or damaged.
	<ul style="list-style-type: none">● Use lukewarm water or thinner soapsuds (neuter detergent) drop into the soft cloth and then dry it to wipe.