## Body Fat Analyzer

Item # 03-303-000

Thank you for purchasing the HealthSmart Body Fat Analyzer. The Body Fat Analyzer emits infrared rays into human body (skin, fat and muscles), and estimates the percentage of fat in organic tissue by analyzing the absorption rate of infrared energy by the human body.

### Operation Instructions

1. Remove cap. Press the ON/SCAN key to turn the device on. When the LCD displays 0.0% it is ready to use.
2. Gently press the sensor against the skin surface.  
**Note:** Do not let any light in between the sensor and the skin.
3. Press the ON/SCAN key. The screen will display “- - -”, meaning the measurement has started. After 0.5 second, the screen will display the results.
4. Press the “M” key to enter the memory mode. Continue to press “M” to recall each memory recorded, “M1” being the first and “M9” being the last.
5. If the unit is not being used for approximately 60 seconds it will automatically shut off.
6. If the battery voltage is below  $2.4V \pm 0.2V$ , the device may not give an accurate measurement. Replace the batteries with new ones.


### Cleaning Instructions

To clean the sensor wipe with rubbing alcohol (70% isopropyl alcohol) or soap and water. Do not submerge in any liquid. Wipe the unit with a dry cloth and allow to air dry.

### Caution

- When measuring, DO NOT let any light get between the sensor and the skin surface.
- The results of the measurements is related to the pressure. Try to apply the same pressure to every measurement.
- The results of the measurements may not represent actual body fat ratio, but rather an estimation of relative body fat.

Toll-Free Customer Care Help Line: 800.622.4714  
Monday – Friday 8 a.m. to 6 p.m. EST

### Specifications

Measuring Range	7.0%~50.0%
Measurement Accuracy	$\pm 1.5\%$ when measuring a 13.4% to 20.0% reference object
Digital Increment	0.1%
Measuring	Recommend back of hand, arm, abdomen or thigh
Measuring Time	< 1 second
Out of Range Alarm	Displays “HHH” when result < 50% Displays “LLL” when result > 7.0%
Sound Indication	Beeps when power-on/measuring/ memory checking
Auto Power-off	$60 \pm 20$ seconds
Low Battery Detection	If the voltage is below $2.4V \pm 0.2V$ , the screen displays symbol
Memory Function	Nine
Current	Work current < 5mA, peak current < 50mA, leakage current when power-off < 3uA
Environment for storage	14°F – 140°F (-10°C – 60°C)
Environment for use	50°F – 104°F (10°C – 40°C)
Battery	Use 2 AAA battery
Battery life	Over 5000 measurements
Dimensions (LxWxH)	4.3” x 2” x 1.25”
Weight	1.7 oz.

This is not a medical device and can not be used for medical diagnostic and treatment activities.

Manufactured for:  
Briggs Medical Service Company  
1931 Norman Drive  
Waukegan, IL 60085

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© 2011 Briggs Medical Service Company