# Smilyss®

FR

ES

NL

IT

EN

DE

MODE D'EMPLOI

## JET DENTAIRE «PRESSION EXPERT»

**SPORT-ELEC® GROUP** - 31 Rue du Val Breton - BP 35 - 27520 BOURGTHEROULDE - FRANCE  
TEL : + 33 (0)2 32 96 50 50 - FAX : + 33 (0)2 32 96 50 59 - [info@sport-elec.com](mailto:info@sport-elec.com)

[www.sport-elec.com](http://www.sport-elec.com)

## JET DENTAIRE «PRESSION EXPERT»

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
## ORAL IRRIGATOR «PRESSION EXPERT»

### IMPORTANT SAFEGUARDS RELATED TO THE USE OF ELECTRICAL PRODUCTS .....

- ▶ Before using your oral irrigator read this manual attentively, it contains important information regarding the use, safety and maintenance of your device. Keep this instruction for future reference if necessary.
- ▶ Only connect the appliance to electrical socket protected by Residual Current Device (RCD) with a maximum rating of 30 mA.
- ▶ Use this appliance only with alternating 100 or 240 VAC voltage (50/60Hz),
- ▶ Before plugging in the appliance, you must make sure that the electrical rating (voltage and frequency, indicated on the rating plate of the device) corresponds to the household electrical supply.
- ▶ The unit must be connected to an electrical supply system that complies with safety standards.

### ATTENTION : NEVER IMMERSE EITHER POWER UNIT OR POWER ADAPTER IN WATER. NEVER USE ANY ABRASIVE PRODUCT OR A SPONGE.


- ▶ Oral irrigators can be used by children and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
- ▶ Cleaning and user maintenance shall not be made by children without supervision.
- ▶ Oral irrigators can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge under supervision.
- ▶ Children shall not play with the appliance.
- ▶ Do not use the product if it has a damaged cord or plug, if it is not working properly. Turn off and unplug the appliance.
- ▶ Before any cleaning, maintenance or service, and after use, turn off and unplug the unit.
- ▶ Do not let cord hang over the edge of a table, sharp corners.
- ▶ Keep the device away from heated surfaces
- ▶ Only use an extension cord suited for this device and in good operating condition.
- ▶ When unplugging the unit from the power outlet, never pull on the power cord.

- ▶ Children should be supervised to ensure that they do not play with the appliance.
- ▶ Do not use the product in a rather humid environment (sauna, hammam), close to electrical machinery equipments linked to the earth and to grounded electrical outlets, i.e. grounded piping system
- ▶ Don't immerse in water or other liquid.
- ▶ Do not use the device outdoors.
- ▶ Do not place the appliance when it functions
- ▶ Only use the accessories provided with the device.
- ▶ Always switch off the unit and unplug it before fitting or removing accessories..
- ▶ Do not operate with fumes or explosive substances.
- ▶ Do not attempt to repair the device yourself. All repairs must be performed by an authorized service center: Improper repairs can cause danger for the user
- ▶ The device is for household use only. Improper use that does not conform to instructions will void the warranty and cannot engage the responsibility of the manufacturer.

### CONNECTING THE UNIT .....

- ▶ The unit base can be connected to any AC outlet 100-240V AC 50/60Hz, including sockets for razors.
- ▶ The power supply may be slightly heated during use, but this has no effect on its operation.

### INSTALLATION OF WALL HOLDER.....

- ▶ Choose a place on a wall away from sources of moisture (shower, sink ...) and close to a power outlet.
- ▶ Make the holes for the two screws to be used to attach the provided plate (mind to use screws and plugs fitted

to your wall).

- ▶ Fix the plate by screwing to the wall.
- ▶ Push the back side of the irrigator on it.
- ▶ Check that your oral jet is stable. SEE sketch-cons.

### GETTING STARTED .....

1. Remove water tank from the base. Be careful to hold the base (1) in place with one hand as you lift the tank with the other hand. Fill the watertank with water. Replace the tank over the base and press down firmly.

2. Plug the power cord into a supply outlet.

**Note:** Before first use, rinse out a full tank of water. If the device is used with a «mouthwash», rinse the device with clean water after each use to remove any residue. If the unit has not been used for some time, it is recommended for reasons of hygiene to follow the instructions for initial employment and rinse out a full tank of water

3. Choose a nozzle and put it on the handle and turn it until it stops.

4. Turn on the device by pushing the On /OFF  power switch on the base of the unit. Lean low over sink and place the tip in mouth. Aim the tip toward teeth. Turn the knob « water/stop » at the top of the handle to "ON"(I) position.

**Note:** Keep your elbow up on height of your mouth, slightly close lips to avoid splashing but allow water to flow freely from mouth into the sink.

5. Turn the pressure regulating button (indicated by drawings "water drops") on the base of the unit to a minimum position (one LED lights) at the beginning. You can gradually increase the pressure and set it at desired setting (from LED 1 to LED 4 the lights shows the pressure rise at 5 bars).

6. Direct the jet stream along your gums and between your teeth. Never point the irrigating head directly to your gums. For best results, start in the molar area (back teeth) working toward the front teeth.

7. Once completed, push the water regulating button on the handle to 'off' position (O). Then, switch off the device by pressing the button «On/Off»  on the base of the unit. Remove the nozzle and replace it in the storage compartment. Remove all the liquid from the device. Unplug the unit.

**Note:** Never allow the unit to remain in the position I (On) if this is not used. The water pressure will continue to develop and can damage the device. Always put the water pressure regulator and the power switch of your oral irrigator in the position 'Off'  (O) if the device is not in use.

**WARNING: Oral irrigators shall be used only with cold tap water or warm tap water (less than 40 DegC)**

### RECYCLING .....

Please dispose of the appliance and the battery in accordance with the relevant environmental protection and recycling regulations.

### THE PACK CONTAINS .....

- ▶ 1 oral irrigator
- ▶ 4 nozzles
- ▶ 1 AC power cable
- ▶ 1 kit for wall fixing

### GUARANTEE ..... 2 YEARS

 When the product is at the end of life, you must ensure compliance with national rules. Dispose of used device in an environmentally-friendly by using the existing public collection.

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