Venus Medical

Breast Enhancement Device

User's Guide **Model 2.0**

-  Gebrauchsanweisung
-  Manual de instrucciones
-  Manuel d'instructions
-  Manuale di istruzioni
-  Gebruiksaanwijzing
-  Instructions Manual

 Gebrauchsanweisung.....	2
 Manual de instrucciones	7
 Manuel d'instructions.....	12
 Manuale di istruzioni.....	17
 Gebruiksaanwijzing.....	22
 User manual.....	27

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Note: If you suffer an illness and want to make sure that the device is suitable for you, we recommend that you consult your doctor or specialist. If you suffer from some kind of allergy or rash, stop using the device immediately.

Maintenance

- This device does not require special care.
- In case you want to wash the vibratory units, do so with a slightly damp cloth or towel and dry thoroughly.
- Do not use acetone or any other kind of solvent in any part of the device.
- Do not store the device with the cables tied in a knot.

Time of use

For optimal results it is recommended to use the device between 20 and 30 minutes, 4 to 5 times a week.

Precautions

If you are not going to use the product during a long period of time, it is recommended to remove the batteries to avoid leakage.

- Do not use the device while showering or bathing.
- Do not expose the device to water or drop it into water.
- Do not place the vibratory units onto wet or damp breasts.
- Only use the device according to the instructions provided in the manual, and follow the directions as specified by the manufacturer.

Note: If the device suffers any kind of damage, do not try to repair it. Venus Medical does not have any part that could be repaired by the user.

Contraindications

The use of this device is suitable to all women except:

- Minors (under 18-year-olds).
- If you are pregnant or breastfeeding.
- Those who have suffered from cancer, inflammation or other diseases in the chest area.
- Those who have pacemakers, heart problems, lumps, sores, pains, burns or other chest injuries.
- Those who have suffered from any kind of disease or injury.

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Venus-Medical offers 8 modes or different massage techniques which stimulate your breasts. Each massage is different thus completing the stimulation cycle. The different modes are:

1. Push
2. Pulse
3. Maintenance
4. Acupuncture
5. Vibration
6. Massage
7. Passive care
8. Firming

Our device also has 10 levels of intensity and another 10 to gradually increase the speed as the treatment advances.

USE:

1. Turn on the device by pressing the "ON/▲" button (the default power program). You get 15 minutes in one power and speed level, that can be seen on the LCD screen.
2. You can select the massage mode by pressing the "MODE" button.
3. Press the "ON/▲" button to increase the power. To decrease it just press the "OFF/▼" button.
4. Switch the power program to the speed program by pressing the "S/P" button, then you can increase or decrease the speed with the "ON/▲" and "OFF/▼" buttons.
5. Switch off the device once you have finished the massage.

*VERY IMPORTANT:

You must always apply lubricant onto the vibratory units before starting to use the device.

We recommend using the device for 20 to 30 minutes a day, 4 to 5 times a week. Even though there are no negative side effects related to the use, it has been proven that the best results are obtained by following the instructions here specified.

1. First, remove the device from the box, open the lid of the central unit and insert two batteries (must be type AAA and 1,5 volts each). After correctly inserting the batteries, place the lid back to the starting position.

2. Next, plug in the cable connecting the two vibrating units onto the central unit. Once the cable has been correctly connected, the device is ready for use.

3. At this point, it is important to apply lubricant gel on the vibratory units before starting to use the device.

4. Then place the vibrating units, one on each breast. To achieve a good support and avoid moving or dislocation of the electrodes, we recommend the use of a bra.

5. Finally, activate the micro-computerized central system. After this you will start noticing a gentle massage. The intensity and speed can be regulated according to your initial sensations and more importantly, once you get used to the device, to your own preferences.

These massages help you get a firmer and fuller chest while the sensation is very stimulating and pleasant, perfect for relaxation. Just a few minutes a day will improve the general condition of your chest in volume, shape and firmness.

Note: When the batteries are running low, you may notice some instable or intermittent functions in the sensors. If this happens, you only need to replace the used batteries with new ones and reconnect.

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Description

Venus Medical is a vibrating device designed and developed to improve the health of the female breasts. It holds the CE-certification in accordance with the requirements as specified by the European Union. It consists of a micro-computerized system with a cable attached to the two vibrating units. These electrodes transmit the vibration directly onto the breasts of women performing a soothing massage. The electrodes are made of anti-allergic rubber.

Technical characteristics

- Easy to use micro-computerized system.
- Powered by low voltage batteries (1,5v.).
- Designed to care for the health of the breasts with ergonomic, comfortable and modern design.
- Voltage: 3V DC.
- Electrical power: <2W.

Instructions

Venus Medical stimulates blood circulation in the pectoral area, reactivating the cell renewal and secretion of hormones through vigorous and pleasant vibration massage. It also improves the drainage of the mammary glands and internal secretion. The vibratory effect maintains the ligaments healthy and strong and increases firmness.

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