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(GB)

Congratulations for buying this high quality product. To use all advantages of this thermometer you should familiarize with the manual prior to first measurement. For later use keep the manual stored in a safe place. If you follow the instruction manual you will be happy with your digital fever thermometer for years.

ADVANTAGES

- **Maximum-Thermometer.** The highest measured value remains visible until device switches off.
- **Flexible tip**
- **Acoustic signal** indicates end of measurement
- **Measured value** after approx. 30 sec., depending on method.
- **Fever alarm** from 37,8°C on
- **Extra large digital display with two digits after decimal point.**
- **Waterproof.** Complete cleansing bath possible (max. 30 min.).
- **Storage.** Last measurement memory for comparison with next measurement.
- **Automatic switch off** for battery protection (after approx. 8-10 minutes).
- **Longlife-Battery** Replacement battery available through good specialized dealers.

IMPORTANT INFORMATION

- Do not drop the device. It is not shock-proof.
- Before measurement check intactness of tip and display.
- This thermometer is determined only for measuring human body temperature.
- This thermometer contains small parts which can be swallowed by children. Don't leave children with this device uncontrolled.
- Do not leave the device in direct sunlight or near heat sources.
- Do not open the device (except the battery cover).
- For cleaning use only water or special substances recommended in chapter „Cleaning and Disinfection“.
- Cleansing bath max. 30 min.

INSTRUCTIONS FOR USE

Place the thermometer and press the on-/off-button. After an acoustic signal, all symbols show on the display for a few seconds as a means of checking that the thermometer is functioning correctly. After an acoustic signal the last correctly measured temperature value is displayed indicated by the „M“ symbol in the bottom right corner. Subsequently the device indicates an internal test value of 37°C ±0.02°C. After this the device is ready for measurement.

During the measurement the current temperature is shown indicated by a flashing „C“. When the acoustic signal sounds and the „C“ stops flashing then the temperature increase is less than 0.1°C in between 16 seconds. This means the nearly final temperature if you use the oral or rectal measurement method. The temperature then changes only minimal. With the axillary method you have to keep the specified time of measurement and the acoustic signal can then sound several times.

The measured value remains visible for approx. 8-10 minutes. After that the device switches off itself automatically. To save battery energy switch off the thermometer after use.

The previously measured temperature value remains stored in any case until the next measurement.

MEMORY

Domotherm® Rapid stores the last measured value automatically. The value appears on the display together with a small „M“ mark in the right corner.

MEASUREMENT METHODS

- Oral measurement (in the mouth)**
Be sure to position the sensor tip under the tongue. This measurement takes approx. 1 minute.
- Rectal measurement (in the anus)**
Gently insert the sensor tip 1-2 cm into the anus. This measurement takes approx. 30 seconds.
- Axillary measurement (under the arm)**
Insert the sensor tip under the arm pit so that the sensor tip is in close contact with the skin. This measurement takes approx. 5-10 minutes. From the medical point of view this method always provides inaccurate values. And so it is unsuitable for accurate measuring.

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