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www.braun.com www.oralb.com

4-716-253/00/XII-02 Korea Printed in Germany

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The Braun Oral-B Plak Control ultra power toothbrush has been developed together with leading dental professionals to effectively remove plaque from all areas of your teeth. It has been designed for daily use for the whole family.

Please read the use instructions carefully before first use.

Important Periodically check the cord for damage. Should this be the case, take the charging unit (E) to a Braun Service Centre. A damaged or non-functioning unit should no longer be used.

pti

- Brushhead On switch Off switch A B C
- Handle D
- Е Charging unit Wall holder

fi

Power supply: see bottom of charging unit Handle voltage: 1.2 V

Connecting and charging The toothbrush is electrically safe and can and in the bathroom without hesitation.

- Plug the charging unit into an electrical outlet. Place the handle on the charging unit. Handle must be switched off. Ρlι
- A full charge takes 16 hours and will provide an operation time of approx. 20 minutes.
- For everyday use, the handle can be stored on the plugged-in charging unit to maintain the toothbrush at full power. Overcharging is impossible

Battery maintenance To maintain the maximum capacity of the rechargeable battery, unplug the charging unit and discharge the handle by regular u at least every 6 months. use

ng t

The toothbrush can be used with any leading brand of toothpaste.

- To avoid splashing, guide the brushhead to your teeth before switching the appliance on.
- Slowly guide the brushhead from tooth to tooth, spending a few seconds on each tooth surface.
- Brush the gums as well as the teeth, first the outsides, then the insides and finally the chewing surfaces.
- Do not press too hard or scrub, simply let the brush do all the work.
- To thoroughly remove plaque, brush for at least two minutes.

Timer (D 9525, D 9511 models only) An incorporated timer memorises the elapsed brushing time. After two minutes the timer signals by 2 thin a short stutter.

During the first days of using the toothbrush, your gums may bleed slightly. In general, bleeding should stop after a few days. Should it persist after two weeks, please consult your dentist or hygienist.

Dentists recommend that you replace any toothbrush every 3 months. Therefore, the brushhead is provided with INDICATOR® bristles (I.a). With proper brushing twice a day for two minutes, the colour will disappear half-way down the bristles in a 3-month period, provided you use toothpaste (I.b). If the bristles splay before the colour recedes (I.c), you are exerting too much pressure on teeth and gums.

<u>Replacement brushheads</u> This product can be used with all Braun Oral-B brushheads.

After use, rinse the brushhead for several seconds under running water with the handle switched on (II). Then switch off the handle and take off the brushhead. Rinse both parts separately under running water (III) and wipe them dec them dry.

om tim e to time, also clean the charging

unit with a damp cloth (IV).

For cleaning purposes, the brushhead compartment (not with all models) is remo vable.

This product contains nickelcadmium rechargeable batteries In the interest of protecting the In the interest of protecting the environment, do not dispose of the product in the household waste at the end of its useful life. Disposal ca an take place at a Braun Service Centre.

However, if you prefer to dispose of the battery yourself, please refer to the section «Battery removal».

Subject to change without notice.

This product conforms to the EMC-Requirements as laid down by the Council Directive 89/336/EEC and to the Low Voltage Regulation (73/23 EEC). F

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mmmmm

Open the handle as shown, remove and dispose of it according to local emove the battery environmental regulations.

Caution!: Opening the handle will destroy the appliance and invalidate the guarantee.

Guarantee We grant 1 year guarantee on the product commencing on the date of purchase. Within the guarantee period we will eliminate, free of charge, any defects in the appliance resulting from faults in materials or workmanship, eithe workmanship, either from faults in materials or working isometry by repairing or replacing the complete appliance as we may choose. This guarantee extends to every count this appliance is supplied by Braun or

/ country aun or its here appointed distributor

This guarantee does not cover: damage due to improper use, normal wear or use as well as defects that have a negligible effect on the value or operation of the appliance. The guarantee becomes void if repairs are undertaken by unauthorised persons and if original Braun parts or not used are not used.

To obtain service within the guarantee period, hand in or send the complete appliance with your sales receipt to an authorised Braun Customer Service Centre. Customer