# BRAUN

## Oral-B

D 9525

D 9525 A

D 9525 T

D 9511

D 9513

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## English

The Braun Oral-B Plak Control Ultra has been developed together with leading dental professionals to effectively remove plaque from all areas of your teeth. It has been designed for daily use for the whole family.

Please read the use instructions carefully before first use.

### Important

Periodically check the cord for damage. Should this be the case, take the charging unit (E) to a Braun Service Centre. A damaged or non-functioning unit should no longer be used.

## Description

- A Brushhead
- B On switch
- C Off switch
- D Handle
- E Charging unit
- F Wall holder

## Specifications

Power supply: see bottom of charging unit

Handle voltage: 1.2 V

## Connecting and charging

The appliance is electrically safe and can be used in the bathroom without hesitation.

- Plug the charging unit into an electrical outlet. Place the handle on the charging unit. Handle must be switched off.
- A full charge takes 16 hours and will provide an operation time of approx. 20 minutes.
- For everyday use, the handle can be stored on the plugged-in charging unit to maintain the toothbrush at full power. Overcharging is impossible.

### Battery maintenance

To maintain the maximum capacity of the rechargeable battery, unplug the charging unit and discharge the handle by regular use at least every 6 months.

## Using the toothbrush

The toothbrush can be used with any leading brand of toothpaste.

- To avoid splashing, guide the brushhead to your teeth before switching the appliance on.
- Slowly guide the brushhead from tooth to tooth, spending a few seconds on each tooth surface.
- Brush the gums as well as the teeth, first the outsides, then the insides and finally the chewing surfaces.
- Do not press too hard or scrub, simply let the brush do all the work.
- To thoroughly remove plaque, brush for at least two minutes.

### Timer

An incorporated timer memorises the elapsed brushing time. After two minutes the timer signals by a short stutter.

During the first days of using the toothbrush, your gums may bleed slightly. In general, bleeding should stop after a few days. Should it persist after two weeks, please consult your dentist or hygienist.

Dentists recommend that you replace any toothbrush every 3 months. Therefore, the brushhead is provided with INDICATOR<sup>®</sup> bristles (I.a). With proper brushing twice a day for two minutes, the colour will disappear half-way down the bristles in a 3-month period, provided you use toothpaste (I.b). If the bristles splay before the colour recedes (I.c), you are exerting too much pressure on teeth and gums.

### Replacement brushheads

This product can be used with all Braun Oral-B brushheads.

## Cleaning

After use, rinse the brushhead for several seconds under running water with the handle switched on (II).

Then switch off the handle and take off the brushhead. Rinse both parts separately under running water (III) and wipe them dry.

From time to time, also clean the charging unit with a damp cloth (IV).

For cleaning purposes, the brushhead compartment (D 9525 only) is removable.

## Environmental notice

This product contains nickel-cadmium rechargeable batteries. In the interest of protecting the environment, do not dispose of the product in the household waste at the end of its useful life. Disposal can take place at a Braun Service Centre.

However, if you prefer to dispose of the battery yourself, please refer to the section «Battery removal».

Subject to change without notice.

This product conforms to the EMC-Requirements as laid down by the Council Directive 89/336/EEC and to the Low Voltage Regulation (73/23 EEC).

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## English

### Guarantee

We grant 2 years guarantee on the product commencing on the date of purchase. Within the guarantee period we will eliminate, free of charge, any defects in the appliance resulting from faults in materials or workmanship, either by repairing or replacing the complete appliance as we may choose.

This guarantee extends to every country where this appliance is supplied by Braun or its appointed distributor.

This guarantee does not cover: damage due to improper use, normal wear or use as well as defects that have a negligible effect on the value or operation of the appliance. The guarantee becomes void if repairs are undertaken by unauthorised persons and if original Braun parts are not used.

To obtain service within the guarantee period, hand in or send the complete appliance with your sales receipt to an authorised Braun Customer Service Centre.

### For UK only:

This guarantee in no way affects your rights under statutory law.

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**Australia**

Gillette Australia  
Pty. Ltd.,  
Private Bag 10,  
Scoresby,  
Melbourne,  
Victoria 3179,  
☎ 1 800 641 820

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**Hong Kong**

Audio Supplies  
Company,  
Room 506,  
St. George's Building,  
2 ICE House Street,  
Hong Kong,  
☎ 25 24-93 77

**Hungary**

Kisgep KFT,  
Pannonia u. 23,  
1136 Budapest,  
☎ 13 494 955

**Iceland**

Verzlunin Pfaff hf.,  
Grensasvegur 13,  
Box 714,  
121 Reykjavik,  
☎ 5 33 22 22

**India**

Braun Division,  
c/o Gillette Div.  
Op. Pvt. Ltd.,  
34, Okhla Industrial  
Estate,  
New Delhi 110 020,  
☎ 11 68 30 218

**Iran**

Tehran Bouran Company  
Irtuc Building,  
No 874 Enghelab Ave.,  
P.O. Box 15815-1391,  
Tehran 11318,  
☎ 02 16 70 23 50

**Ireland (Republic of)**

Gavin's Electronics,  
83/84, Lower  
Camden Street,  
Dublin 2,  
☎ 1800 509 448

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Iskra Prins d.d.  
Rozna dolina c. IX/6  
1000 Ljubljana,  
☎ 386 01 476 98 00

**South Africa**

(Republic of)  
Fixnet After Sales  
Service,  
17B Allandale Park,  
P.O. Box 5716,  
Cnr Le Roux and  
Morkels Close,  
Johannesburg 1685,  
Midrand, Gauteng,  
☎ 11 315 9260

**St. Maarten**

Ashoka,  
P.O. Box 79,  
Philipsburg,  
Netherlands Antilles,  
☎ 52 29 31

**St. Thomas**

Boolchand's Ltd.,  
31 Main Street,  
P.O. Box 5667,  
00803 St. Thomas,  
US Virgin Islands,  
☎ 340 776 0302

**Suomi**

Gillette Group Finland  
Oy/Braun,  
P.O. Box 9,  
Niittykatu 8, PL 9,  
02200 Espoo,  
☎ 09-45 28 71

**Sverige**

Gillette Group  
Sverige AB,  
Dept. Sweden,  
Stockholm Gillette  
Räsundavägen 12,  
Box 702,  
16927 Solna,  
☎ 020-21 33 21

**Syria**

c/o Khalil Daoud  
Amro Ibn Kalthoum Str.  
Unisyria,  
P.O. Box 35002,  
Damascus,  
☎ 11 622 28 81

**Taiwan**

Audio & Electr.  
Supplies Ltd.,  
Brothers Bldg.,  
10th Floor,  
85 Chung Shan N Rd.,  
Sec. 1,  
Taipei (104),  
☎ 0 80 221 630

**Thailand**

Gillette Thailand Ltd.,  
175 South Sathorn Road,  
Tungmahamek, Sathorn,  
11/1 Floor, Sathorn  
City Tower  
Bangkok 10520  
☎ 02-344 91 91

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**Yemen (Republic of)**

Saba Stores for Trading,  
26th September Street,  
P.O. Box 5278,  
Taiz,  
☎ 967 4-25 23 80

**Yugoslavia**

BG Elektronik,  
Bulevar kralja  
Aleksandra 34,  
11000 Beograd,  
☎ 11 3240 030

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Driva S.A.,
Venezuela S.A.,