
Panasonic®

Reach Easy Rolling Massager
Masajero Rodante de Fácil Alcance
Vibromasseur pratique

Operating Instructions Instrucciones de funcionamiento Manuel d'utilisation

Model No. EV2510

Before operating this unit, please read these instructions completely.
Antes de utilizar esta unidad, lea completamente estas instrucciones.
Il est recommandé de lire attentivement ce manuel avant d'utiliser l'appareil.

IMPORTANT SAFETY INSTRUCTIONS

When using an electrical appliance, basic precautions should always be followed, including the following:

READ ALL INSTRUCTIONS BEFORE USING THE MASSAGER.


DANGER -To reduce the risk of electric shock:

1. Always unplug this appliance from the electrical outlet immediately after using and before cleaning.
2. Do not reach for an appliance that has fallen into water. Unplug immediately.
3. Do not use while bathing or in a shower.
4. Do not place or store the appliance where it can fall or be pulled into a tub or sink. Do not place in or drop into water or other liquid.

WARNING -To reduce the risk of burns, fire, electric-shock, or injury to persons:

1. An appliance should never be left unattended when plugged in. Unplug from outlet when not in use, and before adding or removing accessories.
2. Do not operate under blanket or pillow. Excessive heat may develop and cause fire, electric shock, or injury to persons.
3. Close supervision is necessary when this appliance is used by invalids or disabled persons.
4. Use this appliance only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
5. Do not use on or near children.
6. Keep hair away from gaps in massager. Hair may become caught causing damage and injury to the scalp.
7. Never operate this appliance if it has a damaged cord or plug, if it is not working properly, or if it has been dropped or damaged, or dropped into water. Call your nearest authorized service center for repair.
8. Do not carry this appliance by power cord or use cord as a handle.
9. Keep the cord away from hot surfaces.
10. Never operate the appliance with the air openings blocked. Keep the air openings free of lint, hair and the like. Never operate on a soft surface such as a bed or couch where the air openings may be blocked.
11. Do not use by placing under the body, such as under the back, with weight on the massager.
12. Do not use outdoors.
13. Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
14. To disconnect, turn all controls to the off position, then remove plug from outlet.
15. Never drop or insert any object into any opening.

SAVE THESE INSTRUCTIONS


Parts Identification

- | | |
|------------------|----------------|
| ① Wavy surface | ⑤ Grip |
| ② Massage roller | ⑥ Power plug |
| ③ Point massager | ⑦ Power cord |
| ④ Power switch | ⑧ Hanging loop |

How to Use

1. Plug the power plug into an AC outlet.

2. Set the power switch to either of the ON positions.
 - Set the power switch to Strong or Gentle according to your preference.

- It is recommended that you use the Gentle setting until you are used to the massager.
 - Be sure to read the section "For a pleasant massage" (Page 3).
3. After use, set the power switch to Off.
 - After use, be sure to unplug the power plug from the AC outlet.
 - When carrying the massager, be sure to hold it by the grip. Never carry it by the power cord.

Note: During use, electrical interference may occur on TVs, etc.

Use the massager on these areas:

Shoulders and back

Neck


Lower back

It can also be used on these areas:


Arms

Thighs and calves

Soles of feet


For a pleasant massage


Following these methods when using the massager will provide a more effective massage.

Concentrate on stiff areas.


1. Relieve muscle tension.
Stroke the area with the massage roller.
 - Stroke all around the areas which feel stiff or tired to relax them.
2. Ease muscle stiffness and pain on the back by the unit.
Use the point massager or wavy surface.

With the point massager, press hard against the center of the stiff or tired area. (Using both hands will provide even stronger vibrations.)

With the wavy surface, move the massager around the stiff or tired area while pressing against it.

3. Finishing the massage
Use the massage roller to finish the massage.
Stroke the massaged areas again.

Note

- As much as possible, avoid direct contact with the skin surface. If pain occurs even with the switch set to Gentle, cover the area with a towel or similar object and massage through the towel or object.

These are the effects:

Massage surface	Pleasant massage	Corresponding skills of a masseuse	Main effect on body	Recommended area of use
	Roll gently around the entire area to be massaged.		- Effective for muscle fatigue and promoting blood circulation. - Relieves muscle tension.	Back, lower back shoulders, soles of feet, wide areas
	Press against the area to be massaged.		Loosens stiffness and pain due to use of the muscle structure.	Particularly for shoulders and areas with concentrated stiffness
	Press against the area to massage and move massager around.		Effective for blood circulation in specific areas.	Arms and legs

Warning

1. This device should not be used on any injured, swollen, or inflamed areas, or wounds, without consulting your physician. Massage should be pleasant and comfortable. If you feel pain or discomfort, discontinue use immediately and consult a physician.
2. Avoid using on the head, abdomen, joints, or thyroid gland (just below Adam's Apple).

Caution

1. Do not hold the massager by the cord, or pull the cord sharply.
2. Do not attempt to make any repairs yourself. Do not touch the internal workings of the massager.
3. If the massager becomes soiled, clean with a cloth dampened with a mild cleanser. Do not use thinner, benzine or similar solvents. Do not immerse in liquid.
4. Avoid storage in any location where there is excessive moisture, such as the bathroom.
5. Do not bend or turn massager head, or expose it to other excessive force. Damage may result.
6. Roller can be removed for easy cleaning by pulling roller straight out. After cleaning with a damp cloth, push roller back in, making sure it firmly snaps in. Do not use massager if roller is not properly in place.

Specifications

Power source: AC 120V, 60Hz
 Power consumption: 3W
 Head vibration distance: Approx. 2mm
 Vibration rate: Strong: 6300 times/min.
 Gentle: 5000 times/min.
 Dimensions: 470 (H) × 85 (W) × 180 (D) mm
 Weight: Approx. 750g

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Division of Matsushita Electric of Puerto Rico, Inc.
Panasonic Canada Inc.
Panasonic de México, S.A de C.V.
Panasonic Latin America, S.A.