## GARANTIE

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## WARRANTY

This appliance is guaranteed for 12 months from the date of purchase. During the period of the guarantee, the distributor will take responsibility for defects occur, the appliance will be repaired or exchanged as the distributor deems necessary. Damage arising from improper use or normal wear and defects, which do not materially affect the performance of the appliance, is not covered under this guarantee. The guarantee is not valid if unauthorized repairs have been attempted or if replacement parts used are not approved by the supplier. The guarantee is valid only if it is stamped and signed by the retailer at the time of the purchase. For all claims made during the period of the guarantee, send the appliance and the warranty card to your retailer. All other claims are excluded from the guarantee unless specifically prohibited by law.

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

**Fig 1****DESCRIPTION**

- 1 - base
- 2 - pressure control
- 3 - tank
- 4 - "water stop" button
- 5 - storage
- 6 - 4 irrigating heads
- 7 - on/off button

**Fig 2****Fig 3****TECHNICAL FEATURES :**

Power sector : 230 V - 50 Hz  
 Power : 46 W  
 Protection index : IPX7  
 CE

**INSTRUCTION MANUAL****PRECAUTIONS OF USE IN ELECTRICAL APPLIANCE CONTEXT :**

During the use of electrical appliance, it is necessary to respect some security instructions in order to avoid any accident. Please read carefully the instruction manual before using the apparatus, and store it somewhere safe for any future need.

It is recommended to use the appliance on a electrical installation using a differential with nominal current loss which doesn't exceeds 20mA.

- Before starting your dental jet, please read carefully the following instructions. They contain important information concerning the use, the safety and the maintenance of your apparatus.

- Please store this instruction manual somewhere safe for any future need.  
 - Please refer to the explanatory figures while reading this manual.

- Only use this with AC current 230 V, check before connecting your device that the tension of your electrical network corresponds to the one mentioned on your product. The device should be connected only on one electrical supply network whose installation is in conformity with the standards of safety.

- Always disconnect and switch off your device in the following cases : in the event of anomaly during use, before any cleaning or maintenance, and after use.

- Do not let hang the cord on the edge of a table, on sharp surfaces or to come into contact with hot surfaces.

- Keep the apparatus away of heat sources.

- Use a extension for your cord only if it is suitable and in perfect state.

- Never disconnect while pulling on the cord.

- Children do not recognize the dangers emanating from an electrical appliance. Do not leave the apparatus under operation within the range of children.

- Never immerse the apparatus in water or in others liquidate.

- Do not use the apparatus outside.

- Do not pose the apparatus when it is under operation.

- Do not use accessories other than those provided with your apparatus, a danger could result from it. Disconnect the socket when you do not use your apparatus, before mounting or to dismantling the accessories or before any cleaning.  
 - Simply clean the apparatus using a dry or slightly wet rag. Do not use metal sponges or abasiv products for cleaning.  
 - Do not make function the apparatus in presence of explosives or of fumes.

**- CAUTION : o avoid any risk of fire or electrocution, not to dismount the apparatus. Any intervention or repair must be carried out by an approved repairer.**

- Do not use your apparatus, disconnect it or do not connect it immediately if :

- it is damaged
- the electrical cord or the socket is damaged
- you noticed an abnormal fonctionning after incident

- Non-professional fix-ups can be at the origin of dangers to the user.

- If the electric cable is damaged, it must be replaced by the Manufacturer, his After-sales service or a person with similar qualification, in order to avoid any danger.

- Your apparatus is planned for domestic use. A non-adequate use and non-conform to the instructions cancels the guarantee and cannot engage the responsibility of the manufacturer.

- Avoid using the apparatus in a damp room.

**- TO ANY DETERIORATION RESULTING FROM THE NON-APPLIANCE OF THESE PRECAUTIONS OF USE AUTOMATICALLY SUPPRESSES THE GUARANTEE SAVE THESE INSTRUCTIONS**

**HOW TO USE THE ORAL IRRIGATOR ? (Figure 1)**

1 - Remove the water tank (3). When removing water tank use one hand to firmly hold the base of the unit (1) in place, and with the other hand, lift the water tank straight up. Fill with lukewarm water and add mouthwash solution if you want. Reinstall the water tank onto the unit.

**Note :** Prior to initial use, rinse the unit by letting a full water tank of water pass through it. If the unit is used with a mouthwash solution, it should be rinsed with clean water after each use to remove any traces of mouthwash. If the unit has not been used for an extended period of time, to maintain hygiene it is recommended that you follow the directions for initial use, letting a full water tank of water pass through.

2 - Select an irrigator tip for your personal use (6). Lift the irrigator handle out of the base and attach the irrigator tip to the irrigator handle. [Figure 2].

3 - Adjust the pressure control (2) on the base unit to minimum pressure (position 1) for first-time use. Increase the pressure over time using position 2 or 3.

4 - Before turning the unit on, check that the 'water stop' (4) switch on the irrigator handle is in the 'off' position (away from the tip). After turning the unit on by pressing the 'on/off' (7) button on the base unit, lean over the sink and position the tip in your mouth. Move the water stop switch to on (towards the tip) to begin the cleaning process.

**Note :** Keep your mouth slightly open and hold your elbow at mouth level so that the water can run back into the sink .

5 - Direct the flow along the gumline (Figure 3) and also between your teeth. Do not point the irrigator tip directly at your gums. For best results, begin with your back teeth and work forward. Consult your dental professional for further recommendations on use.

6 - When finished, move the 'water stop' switch to the off position and turn off the unit. Remove the irrigator tip and replace in the storage. Empty any liquid remaining in the unit. Unplug the unit.

**Note :** NEVER leave the unit on when not in use. Water pressure will continue to build and may cause damage to the unit. You must turn the water stop switch and the power switch to the "O" position when not in use.

**CLEANING AND MAINTENANCE**

1 - After use, empty any liquid remaining in the unit.  
**Note :** If the unit is used with a mouthwash solution, a full water tank of water should be run through after each use to rinse it and remove any traces of mouthwash.

2 - Rinse the irrigator tip and the irrigator handle separately for several seconds under running water.

**Note :** Be sure to point the irrigator handle downward (upside down) when rinsing. Do not allow water to drain into the irrigator handle.

3 - After cleaning, place the irrigator tip back in the storage and the irrigator handle on the base. Unplug the unit.

4 - Clean the base periodically with a soft cloth or paper towel. The base must be dry before plugging it into an electrical outlet.

**MEDICAL SAFEGUARDS**

Please note the following safeguards for your continued safety. Similar to starting any new oral hygiene product, contact your dentist before using Oral Irrigator if you had oral or periodontal surgery.

Contact your dentist if excessive bleeding occurs or if bleeding continues after 2 weeks.

**CONNECTING THE UNIT TO A POWER SOURCE**

Plug the cord into any standard electrical outlet. If a wall switch controls the outlet, be sure it is turned on.

**Caution**

The supply cord cannot be replaced. If the cord is damaged the appliance should be scrapped.

Check the cord occasionally for any damage. Should this be the case, unplug the charging unit immediately and take it to a service centre. Your service centre will replace the complete charging unit. The charging unit cannot be repaired for safety reasons.

A damaged or non-functioning unit may no longer be used.


**CE - (EMC:2004/108/CE)**

**- (LVD:2006/95/CE)**

**RoHs (2002 /95 /EC) :** This product is manufactured without materials which could have a harmful effect on ecology.

When you bin your product at the end of its lifetime, be sure to respect the national ecology regulations. We thus thank you for depositing it in a place provided for that purpose, in order to secure its elimination in an environmentally respectful way.

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