

150mm

# Healthometer® Dial Scale

Model HAP301

## SCALE USE

1. Place scale on a flat, level surface. A hard surface is best for maximum accuracy. Placing scale on a very soft, deep pile rug may affect the weighing accuracy.
2. Locate zero adjustment control on edge of scale and adjust so the indicator reads zero.
3. Step on and off the scale to check the "zero" setting. If pointer does not return to zero, repeat Step #2.
4. Step on scale and position your feet in a manner that will distribute your weight evenly over the entire scale platform. Stand still while reading display.

**NOTE:** It is normal for your weight to vary during a day and from one day to another. Slight changes in your weight are generally normal.

**If service is required, do not return to your retailer.**

For information on service or warranty, you can locate an authorized Sunbeam® service center by visiting our website at [www.healthometer.com](http://www.healthometer.com); emailing us at: [healthometer@consumer.sunbeam.com](mailto:healthometer@consumer.sunbeam.com); or calling Sunbeam Products, Inc. at **1-800-672-5625**.

To assist us in serving you, please have the model number and date of purchase available.

**NOT LEGAL FOR TRADE**

Healthometer® is a registered trademark of Sunbeam Products, Inc.  
Distributed by Sunbeam Products, Inc. Boca Raton, Florida 33431  
Made in China

34741 05/02

130mm

## Limited Lifetime Warranty

### DO NOT RETURN THIS PRODUCT TO THE PLACE OF PURCHASE

Sunbeam Products, Inc. ("Sunbeam") warrants that for the period that you own this product, it will be free from defects in material and workmanship. Sunbeam, at its option, will repair or replace this product or any component of the product found to be defective during the warranty period. Replacement will be made with a new or remanufactured product or component. If the product is no longer available, replacement may be made with a similar product of equal or greater value. This is your exclusive warranty.

This warranty is valid for the original retail purchaser from the date of initial retail purchase and is not transferable. Sunbeam dealers, service centers, or retail stores selling Sunbeam® products do not have the right to alter, modify or any way change the terms and conditions of this warranty.

The warranty does not cover normal wear of parts or damage resulting from negligent misuse of the product. Further, the warranty does not cover Acts of God, such as fire, flood, hurricanes and tornadoes.

Sunbeam shall not be liable for any incidental or consequential damages caused by the breach of any express or implied warranty or condition. Except to the extent prohibited by applicable law, any implied warranty of merchantability or fitness for a particular purpose is limited in duration to the duration of the above warranty. Sunbeam disclaims all other warranties or conditions, express or implied statutory or otherwise. Some provinces or jurisdictions do not allow the exclusion or limitation of incidental or consequential damages or limitations on how long an implied warranty lasts, so the above limitation or exclusion may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state, or province to province.

#### How to Obtain Warranty Service

To receive warranty service information, you can locate an authorized Sunbeam® service center by:

1. Visiting our website at [www.healthometer.com](http://www.healthometer.com)
2. Emailing us at [healthometer@consumer.sunbeam.com](mailto:healthometer@consumer.sunbeam.com)
3. Calling Sunbeam Products, Inc. at 1-800-672-5625.

If your scale is not covered by warranty, or has been damaged, an estimate of repair costs or replacement costs will be provided to you for approval prior to servicing or replacement.

### DO NOT RETURN THIS PRODUCT TO THE PLACE OF PURCHASE

# Healthometer®