Stella

LOW-PROFILE
CERVICAL PLATE

SCIENT'X

YOUR EXPERTISE OUR SUPPORT

Latest innovation in Scier

Stella 1.95mm low profile p

stabilise the cervical sp

Its unique integrat

reduces operating time and

Integrated locking system

➲ The Stella plate integrates a unique locking system to prevent the screw pull-out

nt'x cervical range,
late is designed to
pine during fusion.

ed locking system
increases security.

The Stella plate implanted via anterior approach is designed to stabilise the cervical spine during cervical fusion. Indications include :


- Discs subsidence
- Spine tumors
- Primary or iatrogenic instability
- Cervical disc hernia
- Traumatic lesions
- Cervical corpectomy

Stella plates are available in different lengths making it possible to perform osteosynthesis from one to five levels. Stella plate is precurved thus allowing immediate adaptation to the patient's anatomical curvature. Its 1.95mm thinness facilitates insertion and minimizes the risks of damaging surrounding tissues.The screw fixation system is integrated into the plate and avoids any risk of post-operative screw back-out.

Ultra low profile

- Protection of surrounding tissues thanks to the thinness of the plate

195mm

Safe positioning

- Holes for pre-fixation of the plate with temporary cervical pins to provide plate stabilisation during surgery

Large slot

- Large slot at the intermediate levels to visualise graft vascularisation and control plate alignment
- Multiple possibilities of intermediate screw insertion

Precurved plate


- Pre-shaped plate with both medial and lateral bending to fit with anatomical curvature


Variable or fixed screw trajectory

- Variable or fixed screw insertion for dynamic or static placement

With variable trajectory

With fixed trajectory


1 Surgical protocol / Patient positioning

a The procedure is performed under general anesthesia with endotracheal intubation. The patient is positioned in supine position on a standard table. The head is placed in a neutral position or rotated in the opposite direction to the approach. Rotation should be approximately 30° to limit stress on the sterno-cleido-mastoid muscle. The shoulders are lowered and fixed with self-adhesive strips in order to release the inferior cervical spine.

b The skin incision targets the level of the lesion to be treated. It is marked on the skin after locating the lesion with an image intensifier. It may be horizontal, along a skin fold, extended to the opposite side of the median line in the case of bilateral approach, or beyond according to the number of levels to be treated. The platysma muscle is sectioned on the same plane as the skin, then two upper and lower detachments are created to avoid muscle tension. The aponeurotic planes are sectioned with the finger up to the anterior surface of the spine.

c The muscles along the neck are separated and rugined to thoroughly release the operating area and allow proper retraction using an autostatic distractor.

d Under scopic control the proper level is located, then the Caspar distractor pins are positioned in the middle of the vertebral bodies over- and under-lying the discal space to be treated and parallel to the endplates.


e Then the anterior vertebral ligament is incised, down to the lateral unco-vertebral articular processes, before distracting the space. Discectomy is then carefully performed, if necessary with the help of a microscope.

2 Preparation of the surgical site

The endplate preparation is performed with a curette. This step must be sufficiently thorough in order to extract all discal debris but not so pronounced as not to collapse the cancellous bone, which would lead to secondary graft collapse. The intersomatic height may be restored in different ways : autologous, heterologous graft or interbody cages.

Plate positioning is performed after resecting the anterior osteophytes. This step is necessary to obtain a close contact between the plate and the anterior vertebral body surface. The choice of the plate and its correct positioning are essential steps of the osteosynthesis.

3 Plate bending and holding

All plates are pre-contoured. If necessary, they can be contoured again using the **plate bender** (except for one-level plate). **The plate holder** ensures plate holding and positioning while offering perfect surgical area visibility. The plate must be positioned on the median line of the cervical spine. Removal of the holder will be carried out after prefixing the plate.

INCREASED LORDOSIS :

DECREASED LORDOSIS :


Plate bender
21CIN01


Plate holder
21PRE13

4 Temporary plate fixation

The plate is temporarily fixed by two **cervical pins** positioned at each end. The **holder** ensures the holding of the cervical pins during their positioning. The end of the holder handle allows impaction with the help of a mallet.

Holder for cervical pins
21PPF02



Temporary cervical pins
21FPC01

5 Drill guide positioning

The drill guide is placed in the screw head housings and enables to perform holes alternately. Depending on the case, two types of guide may be used : the first one with a fixed trajectory and the second one with a variable trajectory.

➲ **The drill guide with fixed trajectory** enables to perform holes according to a defined angulation and reproducible within 0° in the coronal plane and 10° in the sagittal plane. Precise positioning is ensured by the pin located on its spherical end which is inserted in the slit of the screw head housings.



Drill guide with fixed trajectory
21GDVF01

✱ *The handle of the drill guide with fixed trajectory can be oriented by simple rotation thus allowing the operation to be carried out regardless of the available space.*

➔ When a different screw angulation is needed, the **drill guide with variable trajectory** can be used. Its spherical end is perfectly adapted to the screw head housings and enables instrument orientation according to the desired angulation.


Drill guide with variable trajectory
21GDVV01

6 Screw hole preparation


The **conic awl** is placed in the drill guide for optimal guidance. It enables the guide hole to be performed, whose depth is limited to 10mm (A). The metal end of its handle allows impaction using a mallet (B).

Conic awl
21PC03

➔ A second **conic awl with stopper** is used without the drill guide. The length of penetration is limited to 10mm thanks to its spherical end which is stopped by the screw head housings.


Conic awl with stopper
21PC04

7 Drilling

Drilling is performed through the drill guide. The drill fixed stop limits its penetration length to 13mm.


Drill Ø 2.4 mm with AO connection
21FOR02AO

Snap-on handle
21MAE02E

8 Screw insertion


Selected screws are positioned using a **hexagonal screwdriver**. Its conical hexagonal end ensures screw holding during the screw insertion. Screws are inserted crosswise : one screw diagonaly opposite to the previous one. Final tightening is carried out when all screws are positioned.
Screw length may be controlled using the measurer.


Hexagonal screwdriver
21TRV04


INSERTION CROSSWISE

9 Removal of the temporary cervical pins

Once fixed, the plate presents optimal bone/implant contact. The cervical pins are removed using the **holder**. The pin head is introduced into the holder end (A) and locked by unscrewing the outer tube (B). Perfectly secured, cervical pins may be then removed safely.


Holder for cervical pins
21PPF02

10 Plate locking

Screws are secured by the locking wings. The central part of the **locking tool** is introduced into the setting hole and its rotation turns the wings in order to ensure final locking. Alignment markings situated on the plate and on the wings indicate that locking is achieved when they are aligned. Unlocking is carried out in the same way, by simple rotation of the device. The shape of the instrument authorizes angulation with the wings.


Locking tool
21TVP01

* *Wings can be turned either clockwise or counter-clockwise*

Final check

A final check is performed before wound closure. The approach is then closed after rinsing. Hemostasis is checked and aspirative drain placed along the anterior surface of the spine. The platysma must be carefully reconstructed and the skin closed by a resorbable intradermal continuous suture.

Implant removal

When it is necessary to remove material, the anterior cervical approach is used, down to the instrumented area. The locking mechanism is then unlocked by rotating the wings. Positioning temporary cervical pins allows to keep the plate in position during removal of the screws. The cervical screws are removed with the help of the **hexagonal screwdriver**. Once all screws have been removed, insertion of a rugine between the cortical bone and the plate facilitates its extraction. After removal of the cervical pins, the plate can be removed using the **plate holder**.

Locking tool
21TVP01


Hexagonal screwdriver
21TRV04

Plate holder
21PRE13

** Wings can be turned either clockwise or counter-clockwise*

Stella

INSTRUCTIONS

IMPLANTS AND INSTRUMENTS

Instructions for use

OBJECTIVE :
Stella implants are designed for orthopedic and/or neurosurgical treatments of the cervical spine.

GENERAL DESCRIPTION :
Stella implants are comprised of plates and fixation screws. The plates are curved to adapt to the lordosis and the vertebra body anatomy. Stella implants have an integrated security system that secures all of the screws. The locking mechanism is activated by rotation. Stella implants are made of implantable titanium alloy parts which conform to ISO 5832-3 and ASTM-F 136-98 standards. They are incompatible with stainless steel and other metals. Stella implants must not be used with parts from other manufacturers.

INDICATIONS :
Stella products are surgical implants for anterior spinal osteosynthesis. They are designed to correct and stabilize the cervical spine, especially in the case of degenerative disease, trauma and tumors. Indications include :
- weakened discs
- spine tumors
- primary or iatrogenic instabilities
- cervical disc hernias
- traumatic lesions
- cervical corpectomies

CONTRA-INDICATIONS :
Contra-indications for Stella implants include :
• local infection or inflammation
• vertebral osteoporosis
• allergy or intolerance to titanium and its alloys
• incompatible patient age or physical state
• any disease not included in the indications

The contra-indications to these materials are identical to those for analogous materials. Stella implants are not designed, destined nor sold for uses other than those indicated.

POSSIBLE ADVERSE EFFECTS :
• infection
• screw displacement due to poor positioning or poor choice of implant size
• lesions to the vertebral bodies adjoining the osteosynthesis
• intolerance to the material

Note : *An additional surgical intervention may prove necessary to correct an adverse effect.*
Warning : *An entirely satisfactory result is not always obtained after every surgical intervention. This is especially true of spinal surgery where numerous extraneous elements can compromise the results. Osteosynthesis performed with Stella implants must be associated with an arthrodesis performed on one or more levels and secured by osteosynthesis.*

SURGICAL PRECAUTIONS :
The surgical operator must be perfectly familiar with Stella products, their application method, instrumental accessories and operating techniques. The choice of plate size is very important : it must be made in function of patient physiology. Plate prehension is performed using an holder adapted to plate profile. The plate must be fixed on the vertebral body. The holes for the screws are prepared using first a conical awl and then a drill if necessary. Screw hole preparation must be performed using the drill guide. Finally, the cervical screws are fixed using the screwdriver. Screw locking is performed after all screws are in place. Batch numbers and references for the implanted Stella plate and associated screws must systematically be recorded in the patient's surgical file. Implants are single use only. Once removed, an implant must never be implanted again.

PACKAGING :
Implants delivered sterilized are conditioned in individual sterile packages. Implants delivered non-sterile and their instruments are delivered in one or more containers, one (or more) "kits", except in the case of non-assorted kits for stocks. Each kit must be closed and sealed. The packaging must be intact upon reception. If the stocked kit system is used, complete kit composition must carefully be checked. The proper state of all implants and instruments must be checked prior to any use. Storage conditions must enable implant, accessories and packaging integrity to be maintained. Damaged products must never be used and must be returned to Scient'x.

DECONTAMINATION, CLEANING AND STERILIZATION :
> **For implants delivered sterile :** implants are sterilized using Gamma radiation at a dose of 25 to 40 kGy. The expiry date is 5 years. The latest utilization date for sterile elements is indicated on the packaging. Re-sterilizing implants which are delivered sterile is prohibited.
> **For implants and instruments delivered non-sterile :** all non-sterile implants and instruments are delivered non-sterile and must be decontaminated, cleaned and sterilized before and after utilization. Decontamination decreases the microbial population and facilitates later cleaning. Furthermore, sterilization will not be effective if the material is not clean. Recommended method :
> **Decontamination :** Immerse implants and instruments in a bactericidal and fongicidal solution such as didecyldimethylammonium chloride diluted to 0.5% (50 mL to 1 litre of warm water). Soaking time : 20 mn. Rinse with distilled water.
> **Cleaning :** Wash implants and instruments in a LANCER type washing machine with the appropriate cleaning products, rinse and dry. Any product which may alter the material is prohibited (bleach, formalin…). The kit must imperatively be sterilized using an autoclave procedure with the following parameters :
> **Sterilization :** We recommend the sterilization mode in an autoclave for implants and instruments :
• preheat 25' at 110°C (230°F, 1 bar)
• vacuum 5' (0.8 bars below the atmospheric pressure)
• heat 5' at 120°C (248°F, 1 bar)
• vacuum 5' (0.8 bar)
• sterilization 18' at 134°C (273°F, 2 bars)
• dry 20' return to room temperature

DISPOSAL :
Product disposal must be performed in accordance with the precautions in effect for operating room wastes, in conformity with NFS94-030 standards.

COMPLAINTS :
Any health professional (for example a surgeon using these products) who has a complaint or has noted a defect which includes implant and/or instrument quality, identification, reliability, security and/or efficacy must inform Scient'x, or the distributor. In the case of a serious incident or the risk of an incident likely to cause or to have caused death or serious degradation to a patient's or a user's state of health, Scient'x (or the distributor) must immediately be informed by telephone or in writing. We request that any complaint be accompanied with the name(s), reference(s) and batch number(s) of the incriminated components, as well as the name and address of the person posing the complaint, its nature and a detailed description of the complaint and an indication of the response requested.

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