



LIBERTY BED ORDER FORM

999-0822-185
REV B 04/11

Model 701

To order NEW beds, please contact a GF Health Products Sales Representative at 1-800-554-9215.

MAIN BED SECTIONS

Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
LBT Series with Standard 3" Locking Caster Carriage					
LBT7174		80" Grid Deck — Full electric with Head, Foot, and Hi/Lo Motor			
Same		76" Grid Deck — Full electric with Head, Foot, and Hi/Lo Motor			
LBTL Series with Optional 4" Total Lock Steel Caster Carriage					
LBTL7174		80" Grid Deck — Full electric with Head, Foot, and Hi/Lo Motor and Total Locks			
Same		76" Grid Deck — Full electric with Head, Foot, and Hi/Lo Motor and Pedal Lock			

NOTE: Footboard with embedded Staff/Nurse Control, Pan Decks, or Wide Decks are not offered on standard Liberty Beds.

TABLE OF CONTENTS:

For new or replacement bed service parts please refer to the following pages.

Main Bed Section (Decking, u-levers, caster carriages, casters, end caps, etc.)	Page 2
36" Bed End Components (Standard boards, mounting kits, etc.)	Page 3
Wallsavers and Mattress Retainers	Page 3
Electronic Components (Cables, motors, control boxes, hand controllers/pendants, etc)	Page 3
3" Twin Replacement Casters and 4" Total Lock Casters	Page 4
Assist Devices (Half rails, foot mounting bracket, assist bars, trapeze adaptors, trapeze supports, I.V. mount, and traction adaptor)	Page 4
Standard 4" Extension Pan Kit	Page 4
Bed Hardware	Page 5
Control Box and Motor Mounting Hardware	Page 6
Miscellaneous Hardware and Literature	Page 6
Standard Matrix II Bed Diagram	Page 7



LIBERTY BED SERVICE PART ORDER FORM

999-0822-185
REV B 04/11

Model 701

The following list of service parts for the Liberty Bed have been provided for your convenience. Please keep this list in a secure location and use as a reference when ordering new or replacement parts. A general bed diagram has been provided on the last page.

To order new or replacement SERVICE PARTS, please contact a Basic American Medical Products Customer Service Representative at 1-800-365-2338 to confirm you are ordering the right part.

MAIN BED SECTIONS		Please see back diagram (page 7) for general part identification.			
Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
SLEEP SURFACE ASSEMBLY PARTS (999-0822-960)					
Same as Part#	999-0822-900G	Main Frame Weldment (Reference-NOT A SERVICE PART)	1		
Same as Part#	999-0822-920G	Standard Head/Back Grid Deck Weldment (76" & 80")	1		
Same as Part#	999-0822-930G	Standard Foot Grid Deck Weldment (80")	1		
	999-0822-936G	Optional Foot Grid Deck Weldment (76")	(1)		
Same as Part#	999-0822-940G	Standard Knee Grid Deck Weldment (76" & 80")	1		
Same as Part#	999-0822-031G	Knee Pivot Channels	2		
Same as Part#	999-0822-036	Leveler Rods	2		
HEAD & FOOT U-LEVERS PARTS					
Same as Part#	999-0822-913G	Head U-Lever Weldment	1		
Same as Part#	999-0822-903G	Foot U-Lever Weldment	1		
Same as Part#	999-0822-012G	Spanner Plates	2		
Same as Part#	999-0822-009	Pivot Bushings (on either side of u-lever cylinder)	4		
Same as Part#	999-0822-008	Pivot Plated Rods (goes through u-lever cylinder)	2		
CASTER CARRIAGE ASSEMBLIES					
Same as Part#	999-0779-955G	Standard Caster Carriage Assembly with four 3" Locking Casters	1		
LBTL83200	-	Optional Locking Caster Carriage Assembly with Pedal Lock at Foot end, Caster Keepers at Head, four non-locking casters, and all mounting hardware.	(1)		
CASTERS (if ordered separately)					
Same as Part#	100-6330-061	Single Non-Locking Casters for Standard Carriage	2		
Same as Part#	100-6330-062	Single Locking Casters for Standard Carriage	2		
Same as Part#	100-6330-061	Single Non-Locking Casters for Locking Carriage	4		
CAPPING (if ordered separately)					
Same as Part#	999-0822-041	Tube Pivot Caps (Black)	8		
Same as Part#	100-4700-015	1.00" Square Tube End Cap (Main Frame)	6		
Same as Part#	100-4700-017	1.00" x 2.00" Rect. Tube End Cap (2 Main Frame & 4 Carriage)	6		
CABLE TIES					
Same as Part#	100-1100-018	Standard Long Cable Ties (#84-1-0811)	6		
Same as Part#	100-1100-060	Wing Style Cable Ties	6		

BED END COMPONENTS			Please see back diagram (page 7) for general part identification.		
Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
*STAFF CONTROL IS NOT INCLUDED IN STANDARD OR CUSTOM BOARD KITS					
AB53600	-	Solar Oak Hd/Ft Boards	1/1		
AB48600	-	Westminster (Med) Oak Hd/Ft Boards	1/1		
AB80600	-	Montana Walnut Hd/Ft Boards	1/1		
AB29600CCR	-	Enter Crown Windsor Mahogany Hd/Ft Boards	1/1		
AB04600CCR	-	Center Crown Wild Cherry Hd/Ft Boards	1/1		
NA	NA	Center Crown Figured Mahogany Hd/Ft Boards	-		
STANDARD MOUNTING KIT 999-0822-000 INCLUDES:					
Same as Part#	999-0822-034G	LH Mounting Bracket (for mounting boards to frame)	2		
Same as Part#	999-0822-035G	RH Mounting Bracket (for mounting boards to frame)	2		
Hardware Bag 999-0822-999 Includes the following . . .					
Same as Part#	100-5425-081	1/4-20 x 20mm Hex Drive JCB-B Bolts	8		
Same as Part#	100-9900-060	M4 4mm Hex Allen Wrench	1		
Same as Part#	100-7931-007	5/16 x 1.50" Clevis Pin with hole (to adjust bracket position)	4		
Same as Part#	100-2001-006	1-9/16 Hair Pin Clip (to adjust bracket position)	4		
Same as Part#	999-0822-015	Board Mounting Instruction Sheet	1		

OTHER COMPONENTS			Please see back diagram (page 7) for general part identification.		
Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
Same as Part#	999-0822-180	Standard Wireform Head-end Wallsaver (Works with 76", 80" & trapeze adaptor bracket positions- also works on Locking Carriage Assembly: LBTL83200)	1		
Same as Part#	999-0822-001	Single Wire Mattress Retainer	1		

ELECTRICAL COMPONENTS			Please see back diagram (page 7) for general part identification.		
Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
Same as Part#	999-0822-052	Head Actuator/Motor	1		
Same as Part#	999-0808-204	Head Actuator/Motor Cable (25.50")	1		
Same as Part#	999-0822-053	Foot/Knee Actuator/Motor	1		
Same as Part#	999-0822-204	Foot/Knee Actuator/Motor Cable (20")	1		
Same as Part#	999-0778-051	Hi/Lo Actuator/Motor	1		
Same as Part#	999-0778-204	Hi/Lo Actuator/Motor Cable (8")	1		
Same as Part#	999-0779-300	110 V J-Box/Control Box	1		
Same as Part#	999-0775-208	110 J-Box/Control Box 3-Prong Power Cable w/Ground	1		
Same as Part#	999-0779-230	220V J-Box/Control Box	1		
Same as Part#	999-0775-305	Pendant/Hand Controller (with coiled cable)	1		
Same as Part#	999-0791-000	Pendant Holster	1		
Same as Part#	999-0822-100	Hi/Lo Motor Shaft Bushings	2		
Same as Part#	999-0822-083	18GA Ground Wires	3		
Same as Part#	100-5710-019	#10-24 x .375" HWHTC Ground Screws TYP23	7		
NOTE:	There is no Staff Control available for the Liberty Bed as of the printing of this form.				

REPLACEMENT 3" CASTERS & OPTIONAL 4" CASTERS

Purchase #	Part #	Description	Qty Needed On Bed	Qty Ordered	Order Date
Same as Part#	100-6330-062	Single Locking 3" Twin Casters with 7/16-14 Stem	4		
Same as Part#	100-6744-001	7/16-14 Full Nylon Lock Nut	4		
Parts for Converting/Replacing) to 4" Total Lock Steel Casters (no kit - order individual parts) . . .					
Same as Part#	100-6340-015	4" Double Lock Caster with .513" Dia. x 1.625" Stem	4		
Same as Part#	999-0806-150	4" Caster Mounting Bracket Adapter	4		
Same as Part#	100-6750-005	1/2-13 High hex Head Nylon Lock Nut	4		

ASSIST DEVICES

NOTE: Rail "Kits" include both Quick Release Lanyard and Locking Bolt/Nut mounting hardware

Purchase #	Part #	Description	Qty Needed On Bed	Qty Ordered	Order Date
Assist Rails					
ZA82800 (both rails & hardware)	-	Left - Head Half Side Rails (76"/80") Right - Head Half Side Rails (76"/80")	1/1/4 for head end		
ZA82900 (includes both rails & one foot bracket kits)	-	Left - Foot Half Side Rails (76"/80") Right - Foot Half Side Rails (76"/80") Foot Bracket Kits with hardware	1/1/4 for foot end		
-	999-0794-907G	Rt. Head (76"/80") or Lft. Foot (76"/80") Half Rail	1		
-	999-0794-908G	Lft. Head (76"/80") or Rt. Foot (76"/80") Half Rail	1		
-	999-0829-905	Foot Bracket Hardware Bag Only (Includes 2 brackets, 2 cross tubes and mounting hardware)	1		
-	999-0829-023G	Single Cross Tube Only	1		
-	999-0829-025G	Single Foot Bracket Only	1		
Assist Bars					
ZA80300 (includes hardware)		One Head Assist Bar with Mounting Hardware (Includes Quick Release Lanyard and Bolt/Nut Hardware)	1/2 for head end		
-	999-0829-905	Foot Bracket Hardware Bag - Must Be Ordered to attach Assist Bar to foot end (Includes 2 brackets, 2 cross tubes and mounting hardware)	1		
-	999-0803-902G	One Assist Bar only (without hardware)	1		
-	999-0784-100	Quick Release Lanyard Assembly only	1		
Trapeze Support					
ZA83000	-	Trapeze Support Adaptor Kit (Includes Trapeze Support and attaching hardware)	1 for head end		
-	999-0830-901G	Trapeze Support Weldment Only	1		
-	999-0830-905	Trapeze Support Hardware Bag Only (Includes mounting hardware and 4mm wrench)	1 Bag		
ZA78100		Trapeze Bar (Patient Helper Bar) Assembly			

STANDARD 4" EXTENSION PAN KIT

Purchase #	Part #	Description	Qty Needed On Bed	Qty Ordered	Order Date
LBT83300	-	4" Standard Extension Foot Pan Kit (Includes Pan, LH & RH brackets, and mounting hardware)	1 kit		
-	999-0833-901G	LH Footboard Extension Bracket only	1		
-	999-0833-902G	RH Footboard Extension Bracket only	1		
-	999-0833-995	Extension Pan Hardware Bag only	1		

WIDE KITS AND WIDE 4" EXPANSION KITS ARE NOT AVAILABLE ON THE LIBERTY BED

BED HARDWARE

Purchase #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
DECK TO FRAME MOUNTING HARDWARE					
Same as Part#	100-5431-002	5/16-18 x .375" Hex Washer Head Machine Screw Secures Head & Knee Decks to Pan Seat and Knee to Foot Deck	6		
Same as Part#	999-0032-010	Short Self-Wrenching Nut Secures Head & Knee Decks to Pan Seat and Knee to Foot Deck	6		
Same as Part#	999-0808-455	12mm x 55mm Clevis Pin with hole Secures Head & Knee Decks to Pan Seat Brackets	2		
Same as Part#	100-7518-001	7/16 - .50" Circle Cotter Secures Clevis Pin for Head & Knee Decks to Pan Seat Brackets	2		
Same as Part#	100-8031-001	5/16" External Retaining Ring Secures leveler Rods to Frame Leveler Mounting Weldments	2		
U-LEVERS TO MAIN FRAME AND RADIAL ARM MOUNTING HARDWARE					
Same as Part#	100-5725-005	1/4-20 x .75" HWH Self-Tapping Screw - Type F Secures Tube Pivot Caps to Vertical Main Frame Tubes	4		
Same as Part#	100-8050-008	1/2" E-Style Retaining Ring Secures Spanner tubes to U-Lever Pivot Rods	4		
ELECTRONIC COMPONENTS TO MAIN FRAME AND DECKING					
Same as Part#	999-0808-455	12mm x 55mm Clevis Pin with hole Secures Motors to Decks, Pivot Channels & Main Frame Brackets	5		
Same as Part#	100-	1/2" Reg .062 Rue Ring Secures Standard Locking Casters to Carriage Mounting Brackets	5		
Same as Part#	100-7518-001	7/16 - .50" Circle Cotter Secures Clevis Pins to Decks, Pivot Channels & Main Frame Brackets	1		
Same as Part#	100-7938-002	3/8-1.75" Clevis Pin with hole Secures Hi/Lo Motor Shaft to Foot Deck Bracket	1		
Same as Part#	100-8431-019	.688OD X .344ID X .055"Thk ID SAE Flat Washer Seats between Hi/Lo Motor Flange Bushings and Pivot Channel Arms	2		
Same as Part#	999-0822-100	Hi/Lo Motor Shaft Flange Bushings Seats inside both sides of Motor Shaft Hole	2		
Same as Part#	100-5710-019	#10-24 x .375" Hex Washer Head Ground Screw Secures Ground Wire to Main Frame Bracket	1		
STANDARD CASTER CARRIAGE MOUNTING HARDWARE					
Same as Part#	100-6744-001	7/16-14 Full Hex Head Nylon Lock Nut Secures Standard Locking Casters to Carriage Mounting Brackets	4		
Same as Part#	100-5725-005	1/4-20 x .75" HWH Self-Tapping Screw - Type F Secures Tube Pivot Caps to Vertical Caster Carriage Tubes	4		
OPTIONAL PEDAL LOCK & CASTER KEEPER HARDWARE					
Same as Part#	100-5431-025	5/16-18 x 2.50" Machine Screw (part of 999-0832-980G) Secures Pedal Lock to Carriage Cross Brace	1		
Same as Part#	100-6731-006	5/16-18 Thin Nylon Lock Nut (part of 999-0832-980G) Secures Pedal Lock Machine Screw to Carriage Cross Brace	1		
Same as Part#	999-0809-150G	Inner Keeper 2 used on Optional Pedal Lock Carriage Caster (head end)	2		
Same as Part#	999-0809-151G	LH Outer Keeper 1 used on Optional Pedal Lock Carriage Caster (head end)	1		
Same as Part#	999-0809-151G	RH Outer Keeper 1 used on Optional Pedal Lock Carriage Caster (head end)	1		
Same as Part#	100-6744-002	7/16-14 Thin Hex Head Nylon Lock Nut Secures Non-Locking Casters to Head End of Pedal Lock Carriage	2		
Same as Part#	100-5425-059	1/4-20 x 2.50" Hex Head Cap Screw Used to secure Inner Keeper to Outer Keeper	2		
Same as Part#	100-6725-003	1/4-20 Thin Nylon Lock Nut Secures Hex Head Cap Screw to Inner Keeper and Outer Keeper	2		

CONTROL BOX AND MOTOR MOUNTING HARDWARE					
Purchase #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
CONTROL BOX MOUNTING HARDWARE					
-	Cannot be ordered separately	Control Box Clip (comes with control box) Secures Control Box to Head Motor Chassis	1		
-	Cannot be ordered separately	Control Box Cord Retainer Clip (comes with control box) Secures Cables to Control Box ports	1		
Same as Part#	100-1100-050	17.5" Nylon Cable Tie Secures Control Box to Head Motor Chassis	1		
Same as Part#	100-5710-013	#10-24 x .375" Phillips Pan Head Self-Tapping Screw Secures Ground Wire to Main Frame Rail	1		
HEAD MOTOR MOUNTING HARDWARE					
Same as Part#	999-0808-455	12mm Dia. x 55mm Clevis Pin with hole Secures Head Motor Housing & Shaft to Head Deck & Frame Brace	2		
Same as Part#	100-7500-015	.50" Reg. .062 Rue Ring Secures Frame Brace & Motor Housing Mounting Clevis Pins	1		
	100-7518-001	7/16 Circle Cotter Secures Head Deck & Motor Shaft Mounting Clevis Pins	1		
FOOT MOTOR MOUNTING HARDWARE					
Same as Part#	999-0808-455	12mm Dia. x 55mm Clevis Pin with hole Secures Foot Motor Housing & Shaft to Knee Deck & Frame Brace	2		
Same as Part#	100-7500-015	.50" Reg. .062 Rue Ring Secures Foot Mounting Clevis Pins to Frame Brace & Deck	2		
HI/LO MOTOR MOUNTING HARDWARE					
Same as Part#	999-0775-445	12mm Dia. x 45mm Clevis Pin with hole Secures Hi/Lo Motor Housing to Center Frame Brace	1		
Same as Part#	999-7938-002	.375" x 1.75" Clevis Pin with hole Secures Hi/Lo Motor Shaft to U-Lever Weldment	1		
Same as Part#	100-7500-015	.50" Reg. .062 Rue Ring Secures Hi/Lo Housing and Shaft Clevis Pins to weldment & deck	2		

MISCELLANEOUS HARDWARE & INSTRUCTION LITERATURE					
Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
Same as Part#	100-0510-001	Pendant Alligator Clip (for pendant cable)	1		
Same as Part#	999-0822-190	Service/Instruction Manual for Liberty Bed	1		



STANDARD LIBERTY BED SERVICE PART DIAGRAM

Model 701; LBT7174 Shown

