

LANDICE

From treating the most acute neurological or post-operative conditions, to training world-class athletes, the Landice Rehabilitation Treadmill is the world's most versatile solution for rehabilitative, medical, and sports conditioning therapy.

Five-Year Commercial Warranty




Rehabilitation Treadmill

PACKED WITH REHABILITATION

Treadmill always starts at 0.0 mph and resets to 0.0 mph.

Reversing belt option is ideal for practitioners who prescribe downhill and retro-walking therapies.

AccuTrack contact heart rate monitoring option complements cardiopulmonary and respiratory therapies and conditioning without the need for wearing a chest strap.



ATION SPECIFIC FEATURES



15% grade electric elevation with variable elevation over-speed protection

Ultra high-torque drive moves users weighing up to 400 pounds at a true 0.1 to 12.0 mph.

Patient-attachable safety lanyard (treadmill stops when lanyard is pulled)

Remote start-stop switch option is ideal for modalities in which a therapist must manipulate the patient's legs while remaining in control of treadmill operation.

L7 Rehabilitation Treadmill Specifications

- Zero starting speed
- 3 HP continuous-duty motor
- Actual 0.1 to 12.0 mph speed range
- Ultra high-torque drive
- 400-pound user weight capacity
- 15% grade electric elevation
- Variable elevation over-speed protection
- 20" x 58" four-ply treadbelt
- Low 5¾" step-up height
- VFX^{Plus} shock absorption system
- Extended parallel medical handrails
- Patient-attachable safety lanyard
- 2½" diameter steel rollers
- 1-inch thick reversible deck
- Rust-free aluminum frame
- Size: 32" x 77"
- Weight: 300 lbs.
- Choice of Titanium or Matte Textured Black finish

Commercial Warranty: 5 years parts / 1 year labor*



Electrical Requirements: 120VAC, 60HZ, 15A (International voltages available)

Certifications: ETL certified to UL-1647 and CAN CSA 22.2 no. 68

*Designed for use in rehabilitation settings. Heavy use pay-for-membership facilities such as health clubs, gyms, and YMCAs are not eligible for this warranty program.

Rehabilitation Treadmill Options

- Reversing belt direction for downhill and retro walking
- Hand-held remote stop-start switch
- AccuTrack contact heart rate monitoring
- L8 frame with 22" x 63" four-ply treadbelt and 500-pound user weight capacity
- 100 micro-amp hospital isolation transformer
- **ElliptiMill ORTHOPEDIC SHOCK ABSORPTION SYSTEM**
Provides the ultimate reduction of impact on joints and knees.



Orthopedic 4-ply, cushioned treadbelt
5X magnification

