

Resources

Helpful information and references



E4D TECHNOLOGIES

Table of Contents

What Comes With The System?.....3

What to Do When You Get Home4

Before Integration.....4

Information Resources6

Customer Support and Materials Contact List8

What to Expect When Your Ivoclar Representative Visits9

CDD Program.....10

CAD CAM Supplies and Documentation 11

What Comes With The System?

The E4D CAD CAM system comes with:

- Design Center (including the scanner)
- Milling Center
- Job Server
- Customer Documentation
 - User Manual - covers the entire system from scanning to milling
 - Online Help - available within the Design Center software
- E4D Starter Kit
 - Mill Coolant
 - Defoaming Solution



Lava Ultimate :

5 A2 LT C14 Blocks
5 A2 HT C14 Blocks

Misc. items:

1 RelyX Ultimate Adhesive (A1)
1 Scotchbond Universal Adhesive
1 3M ESPE Retraction Capsule
1 CoJet Sand Blast Coating Agent
1 Lava Ultimate Guide



1 Diamond Twist Paste Kit
1 Traxodent Sample
2 Sample Prep Burs
1 Milling Tools Sample Pack
- 2 Ellipsoidal
- 2 Conical
- 2 Tapered
1 Sample Cord



Telio CAD:

4 Telio CAD Blocks
Telio CS Link Transparent
Telio CS Desensitizer 5g
OptraPol Test Pack

IPS Empress CAD:

4 IPS Empress CAD Blocks
2 Empress Shades
1 Empress Stain
1 Empress Glaze
1 Empress Glaze Liquid

IPS e.max CAD:

4 e.max CAD Blocks
2 e.max Shades
1 e.max Stain
1 e.max Glaze Paste
1 e.max Glaze Liquid
1 e.max Crystallization Tray

Misc. items:

2 Multilink Primer
1 Monobond Plus
1 Ceramic Etching Gel
1 Multilink Automix Trans
1 OptraStick
1 Optrafine Promo Pack
1 Object Fix Putty
Cementation Navigation DVD

What to Do When You Get Home

- 1 Practice on the take-home model.
- 2 Scan your cases (models and impressions) and evaluate your preps.
- 3 Scan each other. Scan the teeth as pre-ops to get practice scanning in the mouth.
- 4 Check the online customer resources described below.

Before Integration


Before your system is delivered and installed, you can practice with the system online with GoToMyPC.

The E4D Post Training Learning Lab uses GoToMyPC online service to connect a system at E4D University to the computer in your office or at home. A high-speed Internet connection is strongly recommended.

Note: This is a different set of computers, email, and password than what you may have used for training before coming to E4D University.

The E4D Post Training Learning Lab is intended to provide access to the software for practice before your system is installed. Use the online exercises to practice margin marking and design. If you would like to follow specific instructions for these patients, download the Exercise Workbook on the Resources page of our website.

- 1 Open your Internet browser.
- 2 Go to **<http://www.gotomypc.com>**.



- 3 Click **Log In**.
- 4 Type **posttraining@d4dtech.com** in the Email field.
- 5 Type **posttraining1** as the password. The password is case sensitive.
- 6 Click **Log In**.
- 7 Select **D4D Technologies** as the Account and click **Go**.


A list of the Learning Lab computers is displayed.

- 8 Click one of the **Connect** buttons.

Note: When a button says In Session, that computer is in use. If there are no available computers, try again later.

- 9 If this is the first time for GoToMyPC on your computer, you will be prompted to install the GoToMyPC software. Follow the prompts to install the software.

Note: Your browser's security may block the download and require you to click the warning and select Download.




- 10 Type **posttraining1** as the Access Code.

- 11 You may see warnings about the screen size and/or sound. If the sound has been muted on the lab computer, click **Yes** to turn on the speakers. Sound is only necessary when watching the videos.

Note: If the Color Quality Settings on your PC are less than Highest (32 bit), the colors on the E4D software screens will appear streaked and will not match the example pictures.

Before you close your session, please **EXIT** the Design Center software on the Home page. This will reset the data for the next user.



Information Resources

There are many resources available for gathering information.

If you have not already registered on our website, go to **www.e4d.com/eco** and enter your information.

Log In to **www.e4d.com/myaccount** or click **Customer Login** in the top right corner of the website. Click **Go to Member Resources** after you log in.

The Resources page (**www.e4d.com/resources**) includes:

- Documentation available for download. Printed copies are available for \$25 each and can be ordered by emailing educationgroup@e4d.com.
- Create/Edit your Dentist Finder information
- Download Patient Marketing materials
- Chairside Chats (recorded webinars)
- Link to 3C Learning Library for E4D (links to cadcamcan.com - ability to search content by subject)
- Training Videos
- Register for the CAD CAM Dental Designer (CDD) program

Dentist Finder is a tool on the website that enables the general public and potential patients in your area to locate you.

Customer Support - 1.800.537.6070 customersupport@e4d.com is available to assist you from 7am to 7pm Central time (Monday through Thursday) and 7am to 6pm Fridays.

Newsletters, Chairside Chat, update information, and more is usually communicated via email. When you create your ECO Member account, you are automatically added to our email list. You may unsubscribe at any time.



The Online Help is always available within the Design Center software; simply click the red question mark. The Help is updated along with your software. Videos are available on the Milling Center; simply click the movie icon.



RESOURCES

WATCH, READ AND LEARN

Welcome, **Jeremy Hiser**

[Edit My Account](#) [Log out](#)

NEVO SYSTEM

- [What's New in NEVO ▶](#)
- [E4D User Manual - NEVO ▶](#)
- [Workflow Sheet NEVO ▶](#)
- [Mill Quick Reference - NEVO ▶](#)
- [Design Quick Reference - NEVO ▶](#)

EDUCATIONAL DOCUMENTS

- [Pre - Training Booklet NEVO ▶](#)
- [Pre - Training Booklet NEVO Scan Only ▶](#)
- [E4D Elements Resources Booklet ▶](#)
- [Learning the E4D Dentist System ▶](#)
- [Prep Guidelines ▶](#)
- [Mill Maintenance ▶](#)
- [Scanning Patterns Overview ▶](#)

CLASSIC SYSTEM

- [What's New in v4.6 ▶](#)
- [E4D Exercise Workbook - v4.6 ▶](#)
- [Workflow Sheet v4.6 ▶](#)
- [E4D User Manual - v4.6 ▶](#)
- [Mill Quick Reference - v4.6 ▶](#)
- [Design Quick Reference - v4.6 ▶](#)
- If you would like to order a hard copy of these manuals go [here](#)

E4D FACULTY RESOURCES

- [Compare User Manual ▶](#)
- [Compass Manual ▶](#)
- [E4D Elements Resources Booklet ▶](#)
- [Mill Status Poster ▶](#)
- [Pre - Training Booklet NEVO ▶](#)
- [Prep Guidelines ▶](#)
- [Study Guide ▶](#)



DENTIST FINDER UPDATE



PATIENT MARKETING MATERIALS



CHAIRSIDE CHAT



CADCAMCAN LEARNING LIBRARY



TRAINING VIDEO LIBRARY



ENROLL IN CDD PROGRAM



NEED ASSISTANCE?

Contact SOS at 1-800-537-6070

Hours: Monday-Friday 7am-6pm CST (Holiday hours may differ; call for more information).

Customer Support and Materials Contact List

E4D system support

E4D Customer Support
1.800.537.6070
866.361.1333 corporate phone
972.234.3557 corporate fax

customersupport@e4d.com
7am-7pm Central Time Mon-Thurs
7am-6pm Central Time Friday

Customer Support Director: Scott Freeman



Customer Support Specialists:



Kami Barnes



Shawna Bell



Jodi Butler
Mgr: Clinical



Charles Chandler



Thandiwe
Chirunga



Shari Davis



Joe Escobar



Charlotte Feitlin



Todd Guiles



Ivan Lopez




Jeff McClung
Mgr: Technical



Rodney Valentine



Darren Wiggins



Spencer Wittu

Clinical materials and accessories (cements, adhesives, stains & glaze, etc.)

Ivoclar Vivadent for IPS Empress CAD or IPS e.max CAD restorations

Rebecca Ballaron, MS

Ivoclar Vivadent
175 Pineview Drive
Amherst, NY 14228
716.691.2248 phone
716.691.2254 fax
rebecca.ballaron@ivoclarvivadent.us

3M ESPE for Paradigm C and Paradigm Z100 restorations

Kevin Reilly

3M ESPE
3M Center Bldg. 275-02-W-08
St. Paul, MN 55144-1000
651.736.9808 phone
651.733.2481 fax
kjreilly@mmm.com

Premier Dental Products Company

John Bonner

Premier Dental Products Company
1710 Romano Drive
Plymouth Meeting, PA 19462
610.239.6022
888.773.6872 Ex. 1022
jbonner@premusa.com

What to Expect When Your Ivoclar Representative Visits

The Ivoclar reps go over the contents of the IPS Empress CAD Kit that comes with each machine, and this allows them to cover the following:

- When to use the different types of blocks: HT vs. LT vs. Multi
- The adhesive bonding process with Multilink Automix cement
- Polishing ceramic restorations with Optrapol
- Accessories that can make the CAD CAM process easier, such as OptraGate and OptraDam

From there, other topics are possible based on what other materials/equipment your office purchased:

- Use of the Programat CS oven
- Basic glazing or characterization
- The steps for using e.max CAD
- Laser setup and use

CDD Program

The CAD CAM Dental Designer Program (CDD) provides motivated operators with the opportunity to gain professional recognition and establish credibility in proficiency with the latest dental CAD CAM technology.

Requirements

Completion of the two-day Elements of Success class at E4D University plus home study that should take 3-6 months, depending upon your use of the system.

Self Paced Study Elements

- Learn from online resources. Visit the following sites and write a summary in 200 words or more on how the information was implemented in your office.
 - Go to www.e4d.com/chairside-chat and watch at least 4 videos on this site or CadCamCan.com.
 - One of the videos must be Chairside Chat - Workflow Efficiency. There is a [20 minute version](#) and a [Director's cut](#) which is 40 minutes. You may watch either version.
 - One of the videos must be [Chairside Chat – Materials](#).

- Completion of ALL of the Workbook Exercises that are preloaded* on the Design Center, following the instructions in the Workbook. See the **CDD Workbook Requirements** page for information on how to duplicate the exercises for multiple CDD candidates in your office. Your book cover will look like one of the examples to the right.

[Download Version 4.6 Workbook](#) or [Download Version 5.0 NEVO Workbook](#)



- * If the Workbook Exercises are not included in the patients list, email cdd@e4d.com for assistance.
- Complete 20 E4D restorations and fill out the **Doctor Signoff Form** (included in this document) as you complete them. Email the completed form to cdd@e4d.com or fax it to 972-234.3557 Attn: Education Department.
- Take **Before and After** pictures of 4 E4D restorations. Email the photos and **bite wing x-rays** of the seated restorations to CDD@e4d.com. Please combine these attachments into one email if possible.
- For the **Final Clinical Case** we will choose one special case from your E4D system for you to complete. A trainer will evaluate the case when you are finished and provide feedback.
- Satisfactory completion of the Final Exam.

How to register a new team member

Go to www.e4d.com/resources, click **Enroll in CDD Program**, and fill out the form. The new team member will receive the CDD packet by email. If you do not see the Enroll in CDD Program option under Resources, use the Customer Login at the top of the page.

CAD CAM Supplies and Documentation

The materials listed below are all items used at E4D University. They are grouped by item type. For new documentation, go to www.e4d.com/resources and use the Customer Log In to see customer documentation.

Documentation	
Name	Vendor
User Manual	E4D Technologies
Exercise Workbook	E4D Technologies
Design Center Quick Reference	E4D Technologies
Milling Center Quick Reference	E4D Technologies

Infection Control		
Name	Vendor	Item Number
Alcohol Prep Pads	Schein	1048298
Blue Barrier	Schein	1011838
Cidex Plus	Schein	5550498EZ (Qt) 5554973EZ (Gallon)
Cidex Plus Test Strips	Schein	5553646EZ
Gloves	Schein	
X-Small		5654510
Small		5658087
Medium		5657431
Large		5659481
X-Large		5651575
Wand Aseptic Sleeves	Schein	6314114

Preparation Design & Tissue Management		
Name	Vendor	Item Number
Traxodent Starter Pack	Schein	3780177
	Premier	9007093
Traxodent Value Pack	Schein	3780176
	Premier	9007091
Two Striper Full Crown Kit	Schein	3780210
	Premier	2013581
Two Striper Inlay/Onlay Kit	Schein	3780213
	Premier	2013582

Cementation		
Name	Vendor	Item Number
IPS Ceramic Etch-- 5% HF Acid Etch	Schein	6235981

Impression and Model Materials

Name	Vendor	Item Number
Bite Registration	Schein	
Applicator Tips		
Small/Yellow		9000764
Large/Green		9000766
Gun		9000770
Blu-Bite Fast Set		1023597
Virtual Cadbite	Schein	9459137
	Ivoclar	607908AN
Earth Stone - Quick Set Stone	Schein	9662932
Orban 1/2 Perio Blade for trimming bite registration	Premier	1004751
T-Loc Triple Tray	Schein	
Anterior	Schein	3787158
Posterior	Schein	3785488

Scanning Aids

Name	Vendor	Item Number
E4D Accent (Box of 25)	Schein	6310011
E4D Accent Plus (Box of 25)	Schein	6317890
Mouth Mirrors (Plastic)	Schein	1006640
NEVO Scanning Tips (Pack of 3)	Schein	6314915
Optical Wipes - Kimwipes	Schein	1017070
Standoffs (Classic System - Pack of 3)	Schein	6310019

Design Aids

Name	Vendor	Item Number
Accufilm I Single Sided Red Articulating Paper	Schein	1865309

Milling Center

Name	Vendor	Item Number
Coolant	Schein	6311524
Defoaming Solution	Schein	6318999
Two Striper E4D Mill Diamonds (Burs)		
Conical	Schein	3781031
	Premier	2016002
Ellipsoidal	Schein	3780560
	Premier	2016001
Tapered	Schein	3786546
	Premier	2016000
Assorted	Schein	3780206
	Premier	2016004

Blocks		
Name	Vendor	Item Number
Burn Out Blocks (BOB)	Schein	6310022
IPS e.max CAD		
Size C-14 Shade A1 HT	Schein	9454269
	Ivoclar	626519
Size C-14 Shade A1 LT	Schein	9459227
	Ivoclar	605410
IPS e.max CAD (continued)		
Size I-12 Shade A1 HT	Schein	9456183
	Ivoclar	626503
Size I-12 Shade A1 LT	Schein	9456166
	Ivoclar	605400
Size C-14 Shade A2 HT	Schein	9457237
	Ivoclar	626520
Size C-14 Shade A2 LT	Schein	9453240
	Ivoclar	605411
Size I-12 Shade A2 HT	Schein	9458094
	Ivoclar	626504
Size I-12 Shade A2 LT	Schein	9455801
	Ivoclar	605401
Size C-14 Shade A3 HT	Schein	9456143
	Ivoclar	626521
Size C-14 Shade A3 LT	Schein	9450063
	Ivoclar	605412
Size I-12 Shade A3 HT	Schein	9450799
	Ivoclar	626505
Size I-12 Shade A3 LT	Schein	9454707
	Ivoclar	605402
Size C-14 Shade A3.5 HT	Schein	9455748
	Ivoclar	626522
Size C-14 Shade A3.5 LT	Schein	9458085
	Ivoclar	605413
Size I-12 Shade A3.5 HT	Schein	9455292
	Ivoclar	626506
Size I-12 Shade A3.5 LT	Schein	9455753
	Ivoclar	605403
Size C-14 Shade B1 HT	Schein	9454742
	Ivoclar	626524
Size C-14 Shade B1 LT	Schein	9454241
	Ivoclar	605414
Size I-12 Shade B1 HT	Schein	9453498
	Ivoclar	626508
Size I-12 Shade B1 LT	Schein	9450599
	Ivoclar	605404

Blocks		
Name	Vendor	Item Number
IPS e.max CAD (continued)		
Size I-12 Shade B2 HT	Schein	9455419
	Ivoclar	626509
Size I-12 Shade B2 LT	Schein	9454176
	Ivoclar	605405
Size C-14 Shade BL1 HT	Schein	9452232
	Ivoclar	626539
Size C-14 Shade BL2 HT	Schein	9456948
	Ivoclar	626540
Size C-14 Shade B2 LT	Schein	9457717
	Ivoclar	605415
Size C-14 Shade BL3 HT	Schein	9453845
	Ivoclar	626541
Size C-14 Shade B3 LT	Schein	9450092
	Ivoclar	605416
Size I-12 Shade B3 LT	Schein	9453255
	Ivoclar	605406
Size C-14 Shade BL4 HT	Schein	9457992
	Ivoclar	626542
Size I-12 Shade BL4 HT	Schein	9458958
	Ivoclar	626538
Size I-12 Shade BL4 LT	Schein	9451325
	Ivoclar	605423
Size C-14 Shade BL4 LT	Schein	9454146
	Ivoclar	605428
Size C-14 Shade C1 HT	Schein	9453469
	Ivoclar	626528
Size C-14 Shade C2 HT	Schein	9453118
	Ivoclar	626529
Size C-14 Shade C3 HT	Schein	9451021
	Ivoclar	626530
Size C-14 Shade D2 HT	Schein	9455953
	Ivoclar	626531
Size C-14 Shade D2 LT	Schein	9452306
	Ivoclar	626501
Size I-12 Shade D2 HT	Schein	9458468
	Ivoclar	626516
Size I-12 Shade D2 LT	Schein	9456210
	Ivoclar	626494

Blocks		
Name	Vendor	Item Number
IPS e.max CAD Impulse		
Size C-14 Shade 01	Schein	9450411
	Ivoclar	638165
Size C-14 Shade 02	Schein	9450413
	Ivoclar	638166
Size C-14 Shade V1	Schein	9450414
	Ivoclar	638167
Size C-14 Shade V2	Schein	9450416
	Ivoclar	638168
Size C-14 Shade V3	Schein	9450417
	Ivoclar	638169
IPS Empress CAD		
Size C-14 Shade A1 LT	Schein	9458075
	Ivoclar	602670
Size I-12 Shade A2 HT	Schein	9451494
	Ivoclar	602637
Size C-14 Shade A2 LT	Schein	9455106
	Ivoclar	602671
Size I-12 Shade A3 HT	Schein	9451708
	Ivoclar	602638
Size C-14 Shade A3.5 LT	Schein	9456786
	Ivoclar	602673
Size I-12 Shade A3.5 HT	Schein	9454501
	Ivoclar	602639
Size C-14 Shade B1 LT	Schein	9457876
	Ivoclar	602674
Size I-12 Shade B1 HT	Schein	9456028
	Ivoclar	602640
Size I-12 Shade B2 HT	Schein	9454499
	Ivoclar	602641
Size I-12 Shade B3 HT	Schein	9452612
	Ivoclar	602642
Size I-12 Shade C2 HT	Schein	9451081
	Ivoclar	602643
LAVA Ultimate		
Size 12 Shade A1 HT	Schein	7770487
	3M	3112A1HT
Size 12 Shade A2 HT	Schein	7770489
	3M	3112A2HT
Size 12 Shade A3 HT	Schein	7770491
	3M	3112A3HT
Size 12 Shade B1 HT	Schein	7770494
	3M	3112B1HT

Blocks

Name	Vendor	Item Number
------	--------	-------------

LAVA Ultimate (continued)

Size 14 Shade A1 HT	Schein 3M	7770470 3114A1HT
Size 14 Shade A2 HT	Schein 3M	7770472 3114A2HT
Size 14 Shade A3 HT	Schein 3M	7770475 3114A3HT
Size 14 Shade B1 HT	Schein 3M	7770478 3114B1HT
Size 12 Shade A1 LT	Schein 3M	7770486 3112A1LT
Size 12 Shade A2 LT	Schein 3M	7770488 3112A2LT
Size 12 Shade A3 LT	Schein 3M	7770490 3112A3LT
Size 12 Shade A3.5 LT	Schein 3M	7770492 3112A3.5LT
Size 12 Shade B1 LT	Schein 3M	7770493 3112B1LT
Size 12 Shade BL LT	Schein 3M	7770481 3112BLLT
Size 12 Shade C2 LT	Schein 3M	7770495 3112C2LT
Size 12 Shade D2 LT	Schein 3M	7770497 3112D2LT
Size 14 Shade A1 LT	Schein 3M	7770467 3114A1LT
Size 14 Shade A2 LT	Schein 3M	7770471 3114A2LT
Size 14 Shade A3 LT	Schein 3M	7770474 3114A3LT
Size 14 Shade A3.5 LT	Schein 3M	7770476 3114A3.5LT
Size 14 Shade B1 LT	Schein 3M	7770477 3114B1LT
Size 14 Shade BL LT	Schein 3M	7770466 3114BLLT
Size 14 Shade C2 LT	Schein 3M	7770479 3114C2LT
Size 14 Shade D2 LT	Schein 3M	7770480 3114D2LT

Blocks		
Name	Vendor	Item Number
Multi Blocks: IPS Empress CAD MT Blocks		
Size C-14 Shade A1 Multi	Schein	9450141
	Ivoclar	602701
Size C-14L Shade A1 Multi	Schein	9457715
	Ivoclar	602707
Size I-12 Shade A1 Multi	Schein	9450054
	Ivoclar	602695
Size C-14 Shade A2 Multi	Schein	9458174
	Ivoclar	602702
Size C-14L Shade A2 Multi	Schein	9455223
	Ivoclar	602708
Size I-12 Shade A2 Multi	Schein	9451804
	Ivoclar	602696
Size C-14 Shade A3 Multi	Schein	9456118
	Ivoclar	602703
Multi Blocks: IPS Empress CAD MT Blocks (continued)		
Size I-12 Shade A3 Multi	Schein	9453108
	Ivoclar	602697
Size C-14 Shade A3.5 Multi	Schein	9454152
	Ivoclar	602704
Size I-12 Shade A3.5 Multi	Schein	9451432
	Ivoclar	602698
Size C-14 Shade B1 Multi	Schein	9457948
	Ivoclar	602705
Size C-14L Shade B1 Multi	Schein	9457356
	Ivoclar	602711
Size I-12 Shade B1 Multi	Schein	9452086
	Ivoclar	602699
Paradigm C		
Size 12 Shade A2	Schein	7770739
	3M	68878
Size 14 Shade A2	Schein	7774559
	3M	68884
Paradigm MZ 100 Composite		
Size 14 Shade A1	Schein	7779726
	3M	3714A1
Size 14 Shade A2	Schein	7770031
	3M	3714A2

Blocks

Name

Vendor

Item Number

Telio CAD

Starter Kit	Schein	9450327
	Ivoclar	632883AN
Size 40 Shade A1 LT	Schein	9450329 (3/box)
		9450337 (9/box)
	Ivoclar	631671AN (3/box)
		631676AN (9/box)
Size 40 Shade A2 LT	Schein	9450330 (3/box)
		9450339 (9/box)
	Ivoclar	631672AN (3/box)
		631677AN (9/box)
Size 40 Shade A3 LT	Schein	9450331 (3/box)
		9450340 (9/box)
	Ivoclar	631673AN (3/box)
		631678AN (9/box)
Size 40 Shade A3.5 LT	Schein	9450332 (3/box)
		9450342 (9/box)
	Ivoclar	631674AN (3/box)
		631679AN (9/box)
Size 40 Shade B1 LT	Schein	9450333 (3/box)
		9450343 (9/box)
	Ivoclar	631670AN (3/box)
		631680AN (9/box)
Size 40 Shade BL3 LT	Schein	9450328 (3/box)
		9450334 (9/box)
	Ivoclar	631669AN (3/box)
		631675AN (9/box)
Size 55 Shade A1 LT	Schein	9450345 (Ea)
		9450352 (3/box)
	Ivoclar	631683AN (Ea)
		631689AN (3/box)
Size 55 Shade A2 LT	Schein	9450346 (Ea)
		9450353 (3/box)
	Ivoclar	631684AN (Ea)
		631692AN (3/box)
Size 55 Shade A3 LT	Schein	9450347 (Ea)
		9450354 (3/box)
	Ivoclar	631685AN (Ea)
		631693AN (3/box)
Size 55 Shade A3.5 LT	Schein	9450349 (Ea)
		9450355 (3/box)
	Ivoclar	631686AN (Ea)
		631696AN (3/box)

Blocks		
Name	Vendor	Item Number
Telio CAD (continued)		
Size 55 Shade B1 LT	Schein	9450350 (Ea)
		9450356 (3/box)
	Ivoclar	631682AN (Ea)
		631688AN (3/box)
Size 55 Shade BL3 LT	Schein	9450344 (Ea)
		9450351 (3/box)
	Ivoclar	631681AN (Ea)
		631687AN (3/box)

Restoration Finishing		
Name	Vendor	Item Number
#00 Stain and Glaze Brush	Schein	1010454
#1 Stain and Glaze Brush	Schein	1012773
Calipers	Schein	9004799
Carbide Burs H557-FG	Schein	2281041
Cotton Pliers	Schein	1003313
DiaSheen Diamond Polish Paste Fine	Schein	2289542
Hemostats	Schein	1004662
Honeycomb Porcelain Furnace Tray (for Empress)	Schein	2312774
IPS e.max Crystallization Tray	Schein	9451229
IPS e.max Crystallization Pins	Schein	9459564
IPS e.max Stain and Glaze Starter Kit	Schein	9450095
IPS e.max Stains and Glazes		
A2-A3.5	Schein	9454763
	Ivoclar	605349
B2-B4	Schein	9457600
	Ivoclar	605350
Crystall Shade I1 (Incisal)	Schein	9452748
	Ivoclar	605353
Crystall Shade I2 (Incisal)	Schein	9456852
	Ivoclar	605354
Crystall Shade 1	Schein	9454763
	Ivoclar	605349
Crystall Shade 2	Schein	9457600
	Ivoclar	605350
Glaze	Schein	9456632
	Ivoclar	605363
Mahogany	Schein	9459602
	Ivoclar	605362
Stain and Glaze Liquid	Schein	9455331
	Ivoclar	605366

Restoration Finishing

Name	Vendor	Item Number
------	--------	-------------

IPS e.max Stains and Glazes (continued)

Sunset	Schein	9458690
	Ivoclar	605358

White	Schein	9450732
	Ivoclar	605356

IPS Empress Stain and Glaze Starter Kit	Schein	9456513
---	--------	---------

IPS Empress Stains and Glazes

A1	Schein	9456090
	Ivoclar	562472

A2-A3.5	Schein	9456512
	Ivoclar	562473

B1	Schein	9452318
	Ivoclar	562475

B2-B4	Schein	9457905
	Ivoclar	562476

Glaze	Schein	9457639
	Ivoclar	562495

Incisal 1	Schein	9454051
	Ivoclar	602471

Incisal 2	Schein	9459639
	Ivoclar	602472

Mahagony	Schein	945222
	Ivoclar	562483

Orange	Schein	9453834
	Ivoclar	562485

Stain and Glaze Liquid	Schein	9450973
	Ivoclar	562495

White	Schein	9459779
	Ivoclar	562481

Metal Firing Pins	Schein	187200
-------------------	--------	--------

Object Fix Flow	Schein	9458518
	Ivoclar	605525AN

Object Fix Putty	Schein	9451668
	Ivoclar	605526AN

Programat CS Oven	Schein	9459235
	Ivoclar	597374

Stain and Glaze Palette	Schein	1324049
-------------------------	--------	---------

Two Striper Finishing Kit	Schein	3780201
	Premier	2013553

Virtual Putty (Fast)	Schein	9450478
	Ivoclar	562840

Virtual Putty (Regular)	Schein	9454699
	Ivoclar	562839