

# Catalogue

## Hearing solutions





# Peace of mind for you and your loved ones

For people with hearing loss, everyday things that others take for granted, like for instance hearing the doorbell, can become a real challenge. But it doesn't have to be like that. Nowadays, there are many ways to improve your hearing. Bellman & Symfon offers products that alert you to signals in your home and helps you to join the conversation around the dinner table. Here you can learn more about them.

## Alerting products

|                      |    |
|----------------------|----|
| Overview             | 3  |
| Receivers            | 5  |
| Transmitters         | 7  |
| Accessories          | 9  |
| Value packs          | 11 |
| Stand-alone products | 12 |

## Listening products

|                                           |    |
|-------------------------------------------|----|
| Overview                                  | 13 |
| Listening systems and personal amplifiers | 15 |
| Product comparison                        | 16 |
| Accessories                               | 17 |





# The Visit system

## Hear your home

The Visit system helps people with hearing loss to notice signals in their home, like for instance the doorbell, telephone or smoke alarm. Always knowing when something happens provides peace of mind and an increased sense of security. Here you can learn more about how the various Visit transmitters and receivers work in a typical home.

### Transmitters

- 1 **Door transmitter**  
Detects when the doorbell rings  
Art. no. BE1410



- 2 **Telephone transmitter**  
Detects when the telephone rings  
Art. no. BE1430 + BE9250



- 3 **Smoke alarm transmitter**  
Detects both smoke and fire  
Art. no. BE1480



### Receivers

- 4 **Alarm clock with bed shaker**  
Wakes you with sound, flashes and vibrations  
Art. no. BE1580



- 5 **Flash receiver**  
Alerts you with bright flashing lights  
Art. no. BE1441



- 6 **Pager**  
Alerts you with discreet vibrations  
Art. no. BE1470



1. Doorbell detection



4. Night time alerting



2. Landline and mobile phone detection



3. Fire detection



5. Daytime alerting



6. Outdoors alerting

# Visit receivers

## Smart, wireless and user friendly

### Flash receiver

BE1441 / BE1442



#### Description

The flash receiver alerts you with bright flashes when a transmitter is activated. It is wireless and features a rotating top so you can direct the light. It is also available with battery backup (BE1442) for full functionality even during power cuts.

|                         |                             |
|-------------------------|-----------------------------|
| Bright flashes          | 30 Candela Xenon lights     |
| Clear identification    | Through coloured LED lights |
| Rotating top            | Easy to direct the light    |
| Battery backup (BE1442) | Works during power outages  |
| Size and weight         | Ø 88 x 155 mm, 320 g        |



### Portable receiver

BE1450



#### Description

The portable receiver alerts you with a loud sound signal when a transmitter is activated. It is easy to use and wireless, so you can bring it into the garden, garage or greenhouse.

|                      |                              |
|----------------------|------------------------------|
| Loud sound signal    | Adjustable up to 93 dB       |
| Clear identification | Individual sounds and lights |
| Easy to set up       | Wall mount bracket included  |
| Long lasting power   | Up to 3 years battery life   |
| Size and weight      | 130 x 165 x 36 mm, 590 g     |



### Alarm clock

BE1580



#### Description

The alarm clock receiver awakes you with sound, flashes and vibrations when it's time to get up or when a transmitter is activated. It also features a night light that emits a mild blue glow.

|                   |                           |
|-------------------|---------------------------|
| Loud sound signal | Increases to over 100 dB  |
| Bright flashes    | Four flashing LED lights  |
| Bed shaker        | Wakes you with vibrations |
| Battery backup    | Works during power cuts   |
| Size and weight   | 121 x 92 x 108 mm, 390 g  |



Alerting signal

Sound

Flash

Vibration



Pager receiver

BE1470



Description

The pager receiver alerts you with vibrations and coloured lights when a transmitter is activated. It is small, lightweight and can be carried in your pocket or on your belt. During the night, you can connect a bed shaker to the charger accessory.

|                      |                                   |
|----------------------|-----------------------------------|
| Discreet vibrations  | On your belt or in your pocket    |
| Clear identification | Individual patterns and lights    |
| Portable             | Incl. a belt clip and safety cord |
| Optional charger     | Features two bed shaker jacks     |
| Size and weight      | 57 x 86 x 29 mm, 70 g             |



Wrist receiver

BE8102



Description

The wrist receiver alerts you with vibrations and coloured lights when a transmitter is activated. Just wear it around your wrist as a regular watch. A desktop charger with a bed shaker connection is included in the package.

|                      |                                |
|----------------------|--------------------------------|
| Discreet vibrations  | Directly on your wrist         |
| Clear identification | Individual patterns and lights |
| Charger included     | Features one bed shaker jack   |
| Customizable         | Includes a canvas wrist band   |
| Size and weight      | 38 x 12 x 49 mm, 27 g          |



# Visit transmitters

## Discreet, maintenance free and easy to set up

### Telephone transmitter

BE1430



#### Description


The telephone transmitter helps you to detect when the phone rings. It signals the receiver that alerts you with sound, lights or vibrations. You can also connect a wide range of Visit accessories to its external trigger input.

|                       |                              |
|-----------------------|------------------------------|
| Reliable transmission | Up to 200 m operating range  |
| Multiple detection    | Four separate trigger inputs |
| Adjustable signal     | Use a signal for each phone  |
| Long lasting power    | Up to 3 years battery life   |
| Size and weight       | 60 x 98 x 30 mm, 110 g       |



### Door transmitter

BE1410



#### Description


The door transmitter helps you to detect when the doorbell rings. When someone rings your doorbell or intercom, it signals the receiver that alerts you with sound, lights or vibrations.

|                       |                              |
|-----------------------|------------------------------|
| Reliable transmission | Up to 200 m operating range  |
| Flexible              | Can also detect the intercom |
| Adjustable signal     | Use a signal for each door   |
| Long lasting power    | Up to 3 years battery life   |
| Size and weight       | 80 x 145 x 36 mm, 190 g      |



### Push button transmitter

BE1420



#### Description

The push button transmitter works as a wireless doorbell. When someone presses the button, it signals the receiver that alerts you with sound, lights or vibrations. It can also be used as portable personal alarm.

|                       |                                |
|-----------------------|--------------------------------|
| Reliable transmission | Up to 200 m operating range    |
| Flexible              | Also works as a personal alarm |
| Adjustable signal     | Use a signal for each door     |
| Long lasting power    | Up to 2 years battery life     |
| Size and weight       | 48 x 66 x 23 mm, 50 g          |



Detection
 





## Smoke alarm transmitter

BE1480



### Description

The smoke alarm transmitter will instantly signal your receiver in case of fire. It features battery supervision that alerts you when the batteries run low and is equipped with toast mode to avoid false alarms while cooking.

|                     |                                     |
|---------------------|-------------------------------------|
| Opto-thermal        | Detects both smoke and heat         |
| Battery supervision | Automatically signals the receivers |
| Toast mode          | Use in areas with cooking fumes     |
| Broadcast mode      | Overrides the radio key settings    |
| Size and weight     | Ø 100 x 30 mm, 110 g                |



## Baby cry transmitter

BE1490



### Description

The baby cry transmitter will notify you as soon as your baby needs you, even if you are asleep or a couple of rooms apart. It signals the receiver that alerts you with sound, lights or vibrations so you are always close to your baby.

|                       |                                  |
|-----------------------|----------------------------------|
| Reliable transmission | Up to 200 m operating range      |
| Flexible              | Adjustable sensitivity and delay |
| Customizable signal   | Use a signal for each baby       |
| Long lasting power    | Up to 1 year battery life        |
| Size and weight       | 80 x 145 x 36 mm, 190 g          |



# Accessories

## Expand your Visit solution

### Mobile phone sensor

BE9250/BE9251



#### Description

The mobile phone sensor is activated when the mobile device's screen lights up, for example from an incoming call or message. Place it on top of the device and connect it to the telephone transmitter or directly to the flash receiver.

|                            |                                                                 |
|----------------------------|-----------------------------------------------------------------|
| <b>Optical detection</b>   | Activated by the screen light                                   |
| <b>Light sensitivity</b>   | Visible light >3 lux                                            |
| <b>Connectors (BE9250)</b> | Sensor: 2.5 mm tele jack plug<br>Adapter: 3.5 mm tele jack plug |
| <b>Battery (BE9250)</b>    | 1 x AAA 1.5 V alkaline<br>Operating time > 1 year               |



### Bed shaker

BE1270

#### Description

The bed shaker generates powerful vibrations that wake even heavy sleepers. Just place it under your mattress or pillow and connect it to any of the Visit receivers or charger accessories.

|                            |                         |
|----------------------------|-------------------------|
| <b>Powerful vibrations</b> | Wakes heavy sleepers    |
| <b>Anti-slip surface</b>   | Stays in place          |
| <b>Safe to use</b>         | Operates at low voltage |
| <b>Cord length</b>         | 2 m                     |
| <b>Size and weight</b>     | 88 x 88 x 27 mm, 120 g  |



### Pager charger

BE1260

#### Description

The pager charger is used to recharge the battery of your pager receiver. Place it on the bedside table, and connect a bed shaker to be alerted if any of the transmitters are activated while you are asleep.

|                        |                         |
|------------------------|-------------------------|
| <b>Charging time</b>   | ~ 6 h                   |
| <b>Battery backup</b>  | Works during power cuts |
| <b>Connectors</b>      | Two bed shaker jacks    |
| <b>Cord length</b>     | 1.7 m                   |
| <b>Size and weight</b> | 78 x 88 x 43 mm, 385 g  |



### Repeater

BE1510

#### Description

The repeater is used to extend the operating range of your Visit system. It is useful in large areas or facilities fitted with dense walls. The repeater can also be used as multi-function transmitter.

|                              |                             |
|------------------------------|-----------------------------|
| <b>Reliable transmission</b> | Up to 200 m operating range |
| <b>Signal repetition</b>     | Allows for multiple cuts    |
| <b>Mains powered</b>         | Minimizes maintenance       |
| <b>Battery backup</b>        | Works during power cuts     |
| <b>Size and weight</b>       | 60 x 98 x 30 mm, 110 g      |



### Bridge

BE1520

#### Description

The bridge allows you to use Visit products that operate on the 433 MHz frequency range together with those that operate on 868 MHz. It communicates in both directions between the frequencies.

|                              |                             |
|------------------------------|-----------------------------|
| <b>Reliable transmission</b> | Up to 200 m operating range |
| <b>Two-way radio</b>         | Works in both directions    |
| <b>Mains powered</b>         | Minimizes maintenance       |
| <b>Battery backup</b>        | Works during power cuts     |
| <b>Size and weight</b>       | 60 x 98 x 30 mm, 110 g      |



## External microphone

BE9199/BE9200

### Description

The external microphone is used together with the door transmitter. In an apartment you can for instance use it to detect the intercom, while using the door transmitter's internal microphone to detect the doorbell.

|             |                       |
|-------------|-----------------------|
| Cord length | 2.5 m / 0.75 m        |
| Microphone  | Piezoelectric         |
| Connector   | 3.5 mm mono jack plug |
| Weight      | 15 g                  |
| Colour      | White                 |



## External trigger cable

BE9086

### Description

The external trigger cable is used to link the stand-alone alarm clocks BE1350 and BE1370 to the BE1441 flash receiver. Then the alarm clock will wake you if any of the transmitters are activated while you are asleep.

|                       |                         |
|-----------------------|-------------------------|
| Cord length           | 1.5 m                   |
| Alarm clock connector | 3.5 mm stereo jack plug |
| Receiver connector    | 3.5 mm mono jack plug   |
| Weight                | 20 g                    |
| Colour                | White                   |



## Magnetic switch

BE9023

### Description

The magnetic switch is used with the telephone transmitter to detect when for example a door or window is opened. When the magnets are separated the transmitter signals the receiver that alerts you with sound, lights or vibrations.

|               |                       |
|---------------|-----------------------|
| Cord length   | 0.5 m                 |
| Connector     | 3.5 mm mono jack plug |
| Open switch   | < 1 cm from magnet    |
| Closed switch | > 2 cm from magnet    |
| Weight        | 25 g                  |



## Wall mount bracket

BE9075

### Description

The wall mount bracket is used to mount the flash receiver or the Bellman ringer flash on a wall, so the flashing lights can be seen everywhere in the room. It comes complete with screws, plugs and mounting instructions.

|          |                     |
|----------|---------------------|
| Material | Aluminium           |
| Usage    | For indoor use only |
| Size     | 45 x 90 x 50 mm     |
| Weight   | 110 g               |
| Colour   | White               |



## Outdoor telephone ringer

BE9003

### Description

This outdoor speaker helps you to hear the receiver's signal when you are out in the garden. It is mounted on the house exterior, weather resistant and can be connected to the Bellman ringer and Bellman ringer flash.

|             |                       |
|-------------|-----------------------|
| Cord length | 9.5 m                 |
| Connector   | 3.5 mm mono jack plug |
| Size        | 102 x 115 x 115 mm    |
| Weight      | 500 g                 |
| Colour      | Black                 |



## Contact mat

BE9024

### Description

The contact mat is used together with the telephone transmitter to detect when someone enters or leaves a room. When you step on the mat, the transmitter signals the receiver that alerts you with sound, lights and vibrations.

|             |                       |
|-------------|-----------------------|
| Cord length | 1.5 m                 |
| Connector   | 3.5 mm mono jack plug |
| Size        | 745 x 420 x 3 mm      |
| Weight      | 675 g                 |
| Colour      | Black                 |



# Visit value packs

## An easy way to get started

### Safe value pack

BE8037



**Description**

The Safe value pack will immediately wake you with sound, flashes and vibrations if a fire occurs - even if you sleep several rooms away from the smoke detector. It includes an alarm clock, smoke detector and bed shaker.

|                |                                    |
|----------------|------------------------------------|
| Wireless       | Up to 200 m operating range        |
| Easy to set up | Includes everything you need       |
| Battery backup | Works during power cuts            |
| Reliable       | Approved by healthcare authorities |
| In the package | BE1580, BE1480, BE1272             |



### Care value pack

BE8039



**Description**

The Care value pack will alert you with bright flashes when the doorbell or telephone rings. It includes a wireless doorbell and a telephone cord that is connected to a flash receiver.

|                |                                    |
|----------------|------------------------------------|
| Bright flashes | 30 Candela Xenon lights            |
| Easy to set up | Includes everything you need       |
| Rotating top   | Easy to direct the light           |
| Reliable       | Approved by healthcare authorities |
| In the package | BE1441, BE1420, RJ11 phone cord    |



# Stand-alone products

## Alarm clocks and signal amplifiers

Alerting signal   

### Alarm clock Pro

BE1370



#### Description

This alarm clock will wake you with an extra loud sounding alarm, bright flashes and vibrations from the included bed shaker. It is easy to use and packed with clever features like night light, backup battery and a built-in telephone jack.

|                   |                           |
|-------------------|---------------------------|
| Loud sound signal | Increases to over 100 dB  |
| Bright flashes    | Four flashing LED lights  |
| Bed shaker        | Wakes you with vibrations |
| Battery backup    | Works during power cuts   |
| Size and weight   | 121 x 92 x 108 mm, 390 g  |



### Alarm clock Classic

BE1350



#### Description

This alarm clock will wake you with an extra loud sounding alarm and vibrations from the included bed shaker. It features a clear display and large ergonomic dials that make it easy to set the time and alarm.

|                     |                                |
|---------------------|--------------------------------|
| Loud sound signal   | Increases to over 100 dB       |
| Bed shaker          | Wakes you with vibrations      |
| Easy to use         | Large rotating dials           |
| Illuminated display | Adjustable backlight intensity |
| Size and weight     | 121 x 92 x 108 mm, 390 g       |



### Bellman Ringer

BE1010

#### Description

The Bellman Ringer alerts you with a powerful sound signal and a clear light when the telephone rings. It's easy to adjust the volume, tone and sound pattern to fit your needs. You can also connect a speaker for outdoor use.

|                    |                              |
|--------------------|------------------------------|
| Loud sound signal  | Adjustable up to 95 dB       |
| Flexible settings  | Personalize volume and tone  |
| Long lasting power | Up to 3 years battery life   |
| Easy to set up     | Incl. phone cord and adapter |
| Size and weight    | 95 x 155 x 32 mm, 350 g      |



### Bellman Ringer flash

BE1390

#### Description

The Bellman Ringer flash shows with a bright flashing light when the telephone rings. It is easy to set up and can be used both at home and at work. You can also connect a speaker for outdoor use.

|                 |                               |
|-----------------|-------------------------------|
| Bright flashes  | 30 Candela Xenon lights       |
| Rotating top    | Easy to direct the light      |
| Easy to set up  | Incl. telephone cord and plug |
| Robust design   | Withstands tough conditions   |
| Size and weight | Ø 70 x 140 mm, 300 g          |



# Assistive listening devices

## Hear better in busy surroundings

An assistive listening device helps people with hearing loss to hear better in busy surroundings, such as the workplace or school. It filters out the noise and brings out the words so you can hear clearly even when there is a lot going on around you. This means you don't have to make so much effort to hear what people are saying.

Bellman & Symfon's assistive listening devices are fully digital and work with hearing aids and cochlear implants. We offer both personal amplifiers for quiet family dinners and TV-listening as well as advanced listening systems for busy lectures and conferences.


Here you can read more about our listening products.

### Domino

The Domino listening system consists of a transmitter and a receiver. The receiver (1) picks up speech in your vicinity. Wear it around your neck or place it in front of you. The transmitter (2) picks up speech further away. Place it on the other side of the table or ask e.g. your teacher to wear it.

### Maxi and Mino

Maxi and Mino consists of a single unit with a built-in microphone. Just place it on the table in front of you, adjust the volume to your liking, and the words will come out clearly. You can listen with headphones or use a neck loop if you have hearing aids or cochlear implants.




#### At the restaurant

A listening system will help you to distinguish speech in noisy environments.



#### At work

Put the transmitter on the other side of the table.



#### At home


With a personal amplifier, your partner doesn't have to raise her voice or repeat herself.




#### By the TV

Connect it to your own volume.






transmitter and receiver on either side  
to take part in the conversation.



#### In school

Ask your teacher to wear the transmitter to pick  
up the words in class.



#### In the car

Hand over the transmitter to your travel com-  
panion to filter out the traffic noise.



to the TV or radio and listen at your  
while sparing your family's ears.



#### At the pharmacy

Put it on the counter in front of you to hear clearly  
without missing any important information.



#### At the theatre

Capture the sound with the built-in telecoil in  
public venues fitted with a loop system.

2



# Assistive listening devices

## Discover the joy of better hearing

### Domino listening system

BE8005



#### Description

Domino is a fully digital listening system that brings out speech in noisy environments so you hear the words clearly. It consists of a receiver that picks up speech in your vicinity and a transmitter that picks up speech further away.

|                   |                            |
|-------------------|----------------------------|
| Microphones       | Omni- and directional      |
| Signal processing | 16 bit digital stereo      |
| Radio link        | 2.4 GHz broadband          |
| Batteries         | Rechargeable Lithium-ion   |
| Size and weight   | 99 x 48 x 22 mm, 74 g/unit |



### Mino personal amplifier

BE2030



#### Description


Mino is a digital personal amplifier, designed for those who want a small, lightweight device that fits discreetly in the pocket. Listen with headphones or use a neck loop if you have hearing aids or cochlear implants.

|                   |                                     |
|-------------------|-------------------------------------|
| Microphones       | Omnidirectional and directional     |
| Signal processing | 16 bit digital stereo               |
| Built-in telecoil | For venues fitted with loop systems |
| Batteries         | Rechargeable Lithium-ion            |
| Size and weight   | 99 x 48 x 22 mm, 74 g               |



### Maxi personal amplifier

BE2020



#### Description

Maxi is a robust and user-friendly digital personal amplifier, specially designed for people with impaired vision or low dexterity. Listen with headphones or use a neck loop if you have hearing aids or cochlear implants.

|                   |                                     |
|-------------------|-------------------------------------|
| Microphone        | Omnidirectional                     |
| Signal processing | 16 bit digital stereo               |
| Built-in telecoil | For venues fitted with loop systems |
| Batteries         | 2 x AA 1.5V batteries               |
| Size and weight   | 140 x 45 x 27 mm, 133 g             |



## Comparison

### Domino

### Mino

### Maxi



An advanced listening system that deliver clear and crisp speech even in noisy surroundings



A small and versatile personal amplifier that fits discreetly in your pocket or on your lapel



A robust and user-friendly personal amplifier with large tactile controls and clear indicators

#### Ideal listening situations

Easy

TV-listening  
Private conversation  
Family gathering  
Restaurant and café  
Shop and department store  
School and workplace  
Difficult Lecture and conference

TV-listening  
Private conversation  
Family gathering  
Restaurant and café  
Shop and department store

TV-listening  
Private conversation  
Family gathering

#### Radio link

2.4 GHz broadband, 128-bit encryption  
Range: Up to 40 meters indoors

–

–

–

–

#### Microphones

2 Omnidirectional  
2 Directional

1 Omnidirectional  
1 Directional

1 Omnidirectional  
–

#### Signal processing

16 bit digital stereo  
Frequency range: 20 Hz to 10 KHz  
10-channel dynamic compression  
10-channel dynamic noise reduction  
Adaptive feedback manager  
–

16 bit digital stereo  
Frequency range: 20 Hz to 10 KHz  
10-channel dynamic compression  
10-channel dynamic noise reduction  
Adaptive feedback manager  
Built-in telecoil for loops systems

16 bit digital stereo  
Frequency range: 40 Hz to 10 KHz  
10-channel dynamic compression  
10-channel dynamic noise reduction  
Adaptive feedback manager  
Built-in telecoil for loop systems

#### Power and battery

Rechargeable lithium-ion battery  
Operating time: up to 12 h  
Charging time: up to 2.5 h

Rechargeable lithium-ion battery  
Operating time: up to 18 h  
Charging time: up to 2.5 h

2 x AA alkaline batteries  
Operating time: up to 150 h  
–

#### Buttons and controls

Recessed On/Off button  
9 step volume control (5dB/step)  
9 step tone control (3dB/step)  
9 step balance control (5dB/step)  
Clear indication of Volume, Tone, Mic and RF mode  
Advanced settings menu that contains 9 parameters

Recessed On/Off button  
9 step volume control (5dB/step)  
9 step tone control (3dB/step)  
9 step balance control (5dB/step)  
Clear indication of Volume, Tone, Mic, Telecoil and RF mode  
Advanced settings menu that contains 9 parameters

Recessed On/Off button  
9 step volume control (5dB/step)  
5 step tone control  
9 step balance control (5dB/step)  
Clear indication of Volume, Tone, Mic, Telecoil and PTL mode  
Advanced settings menu that contains 3 parameters

#### Connections

2.5 mm stereo jack  
3.5 mm stereo headphone jack  
Universal charger connector

2.5 mm stereo jack  
3.5 mm stereo headphone jack  
Universal charger connector

2.5 mm stereo jack  
3.5 mm stereo headphone jack  
–

#### Weight and dimensions

99 mm x 48 mm x 22 mm (w x h x d)  
Weight: 74 g / unit

99 mm x 48 mm x 22 mm (w x h x d)  
Weight: 74 g

104 mm x 45 mm x 27 mm (w x h x d)  
Weight: 133 g (incl. batteries)

#### In the package

✓ Domino transmitter and receiver  
✓ Belt clip in stainless steel  
✓ Lanyard with safety clasp  
✓ Durable travel case and earphones  
✓ External microphone and cable kit  
✓ Universal charger  
✓ User manual and quick guide

✓ Mino personal amplifier  
✓ Belt clip in stainless steel  
✓ Lanyard with safety clasp  
✓ Universal charger  
✓ User manual and quick guide

✓ Maxi personal amplifier  
✓ Belt clip  
✓ Lanyard with safety clasp  
✓ 2 x AA alkaline batteries  
✓ User manual and quick guide

# Accessories

## Complement your solution

### Neckloop

BE9159/BE9161



#### Description

The neckloop is comfortable, lightweight and generates high output power. Plug it into your speech amplifier or mobile phone and switch your hearing aids to "T" position. Available with a long (BE9159) and short cord (BE9161).

|                   |                          |
|-------------------|--------------------------|
| High output power | Dual channel output      |
| Soft fabric tube  | Comfortable and hygienic |
| Tangle free cable | Easy to handle           |
| Gold plated plugs | Prevents corrosion       |
| Cord length       | 90 cm or 15 cm           |



### Microset™

BE9129



#### Description

MicroSet™ is an accessory specially developed for Mino which consists of a discreet lapel microphone and an earbud. Attach the microphone on your lapel tuck away Mino in your pocket and leave your hands free for other things.

|             |                         |
|-------------|-------------------------|
| Cord length | 95 cm                   |
| Connector   | 3.5 mm stereo jack plug |
| Microphone  | Omnidirectional         |
| Weight      | 15 g                    |
| Colour      | Black                   |



### External microphone

BE9127

#### Description

The external microphone can be used as a remote microphone during meetings to improve the signal to noise ratio. You can also attach it close to the TV's loudspeaker and use it for TV-listening.

|             |                       |
|-------------|-----------------------|
| Cord length | 5 m                   |
| Connector   | 2.5 mm mono jack plug |
| Microphone  | Omnidirectional       |
| Weight      | 35 g                  |
| Colour      | Black                 |



### Cable kit

BE9126

#### Description

The cable kit allows you to connect the assistive listening device to your TV or stereo and enjoy high quality sound in full stereo. The package includes a 5 m long cable with a stereo jack plug and an RCA (phono) adapter.

|                     |                        |
|---------------------|------------------------|
| Cable cord length   | 5 m                    |
| Cable connectors    | 2.5 to 3.5 mm stereo   |
| Adapter cord length | 25 cm                  |
| Adapter connectors  | Phono to 3.5 mm stereo |
| Colour              | Black                  |



## Headphones

BE9122

### Description

The stereo headphones feature semi-closed drivers and deliver clear sound with minimum feedback. They are equipped with a single side chord and are very robust and user friendly.

|                 |                         |
|-----------------|-------------------------|
| Frequency range | 20 Hz - 20 kHz          |
| Impedance       | 32 $\Omega$ $\pm$ 10%   |
| Sensitivity     | 113 dB/mW @ 1 kHz SPL   |
| Connector       | 3.5 mm stereo jack plug |
| Cord length     | 1.2 m                   |



## Stetoclips

BE9125

### Description

The stetoclips are specially selected for our personal amplifiers and are capable of producing a sound level in mono of up to 140dB. They are comfortable, lightweight and very easy to put on and take off.

|                 |                          |
|-----------------|--------------------------|
| Frequency range | 800 - 3000 Hz            |
| Output level    | 124 dB SPL @ 1kHz, 50 mW |
| Connector       | 3.5 mm stereo jack plug  |
| Weight          | 38 g                     |
| Cord length     | 1.5 m                    |



## Earbuds

BE9124

### Description

The stereo in-ear earbuds are lightweight and deliver clear sound with great detail. They are equipped with closed drivers that result in minimum sound leakage and great isolation from ambient noise.

|                 |                         |
|-----------------|-------------------------|
| Frequency range | 20 Hz - 20 kHz          |
| Impedance       | 16 $\Omega$ $\pm$ 10%   |
| Sensitivity     | 104 dB/mW @ 1 kHz SPL   |
| Connector       | 3.5 mm stereo jack plug |
| Cord length     | 1.2 m                   |



## Travel case

BE9147

### Description

The travel case is specially designed for our range of assistive listening devices. It is made of durable black canvas with a reinforced cover that efficiently protects your listening accessories during travel.

|                  |                          |
|------------------|--------------------------|
| Durable canvas   | Withstands heavy use     |
| Reinforced cover | Keeps your products safe |
| Magnetic lock    | Easy to open and close   |
| Colour           | Black w. white trimmings |
| Size and weight  | 19 x 15 x 5 cm, 240 g    |





## We are passionate about hearing.

Bellman & Symfon of Sweden is dedicated to improve the quality of life for people with hearing loss. Our people and partners are devoted to this mission and we work closely with hearing healthcare professionals and leading experts to make better hearing possible for everyone.

Bellman & Symfon  
Södra Långebergsgatan 30  
421 32 Västra Frölunda  
Sweden  
Phone +46 31 68 28 20  
Fax +46 31 68 28 90  
Email [info@bellman.com](mailto:info@bellman.com)

CAT\_001ART001\_UK

© and TM 2014, Bellman & Symfon AB, all rights reserved. Bellman Visit and Bellman & Symfon are registered trademarks owned by Bellman & Symfon AB.



[bellman.com](http://bellman.com)